AW Family

AW VolumeShare 7

The AW Family is GE Healthcare's ecosystem of processing and advanced visualization solutions across platforms and applications.

Advantage Workstation (AW) with VolumeShare 7 enables multi-modality image review, comparison, and processing with simplicity and power. Featuring, 64-bit technology, this workstation is a cornerstone solution to Radiology departments around the world.

This document is intended to overview the features, specifications, use cases, and other key information of AW VolumeShare 7.

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SECRET DE
AFACERI

Contents

Overview	3
Key Features	3
Industry Standards	3
DICOM Conformance Standards:	3
Filming Protocols	3
Indications for Use	4
Regulatory Compliance	4
Product Details	5
Modalities	5
License Management	5
User Interface	5
Standard Applications & Capabilities	6
Workflow Management	9
Security Capabilities	11
System Components	11
Workstation Configuration	11
About GE HealthCare	13

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SECRET DE
AFACERI

Overview

Key Features

- Simple drag and drop action for networking, media interchange, and filming.
- Multi-tasking capabilities.
- Active directory integration allows enterprise level user authentication. Directories supported include Microsoft® Active Directory® (MSAD) and other LDAP authentication services.
- Patient list to manage images from local workstation or remote systems.
- Postfetch feature, designed to optimize your reading workflow by automatically gathering a patient's prior exams according to your pre-set criteria.
- Support for Key Images.
- End Review automates routine filming and networking tasks with just one click.
- Search Advantage for fast and easy search of a patient's exam history on PACS or any other DICOM®-compliant device.
- Enhanced Quick Filters of the Patient List filters studies by Modality, Date, End Review status or Exam Description.
- DICOM CD/DVD/USB Creation Tool.
- 2D Viewer for image display, manipulation, annotation, review.
- Integrated Filmer with enhanced flexibility to perform all filming and data exporting tasks.
- Access to a wide variety of applications for greater diagnostic flexibility.

Industry Standards

AW VolumeShare 7 complies with a wide variety of industry standards to facilitate adoption of features and performance improvements as the computing and medical imaging industry evolves.

DICOM Conformance Standards:

- DICOM 3.0 Storage Service Class for RT, CT, MR, CR, X-ray (Angio and R&F), Digital X-ray (DX), MG, NM, PET, U/S, Secondary Capture, Secondary Capture Color DICOM Image Objects. (Service Class User (SCU) for image send and Service Class Provider (SCP) for image receive).
- DICOM 3.0 Query/Retrieve Service Class (SCU and SCP).
- DICOM 3.0 Storage Commitment Service Class (SCU).
- DICOM Print (Color and B&W).
- DICOM Media Interchange (CD-R, DVD+R(W)).

Filming Protocols

DICOM Print (Color and B&W) and Adobe® Postscript (Color and B&W) for supported printers.

CONFIDENTIAL

SECRET DE AFACERI

Indications for Use

AW VolumeShare 7 is a review workstation, which allows easy selection, review, processing and filming of multi-modality DICOM images from a variety of diagnostic imaging systems. When interpreted by a trained physician, filmed or displayed images on the AW monitor may be used as a basis for diagnosis, except in the case of mammography images.

Regulatory Compliance

This product complies with the Regulation (EU) 2017/745 of the European Parliament and the Council on the medical device (MDR).

This product or its features may not be available in some countries or regions. Please contact your sales associate.

Rx Only

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SECRET DE
AFACERI

Product Details

The Advantage Workstation VolumeShare 7 software technology makes a multi-modality advanced visualization workflow solution that helps to enhance diagnostic confidence and productivity. This solution features software that is optimized for 64-bit technology and multicore processor hardware to provide leading edge performance. It includes AW's premier 3D image analysis package, Volume Viewer, with a suite of volumetric visualization, and analysis tools for CT, MR, 3D X-ray and PET.

In addition, this package includes the multi-modality 2D Viewer with dedicated functions for review of CT, MR, X-Ray, Angio, DX, U/S and PET images. It also includes the Filmer, a multimedia export tool for creating electronic films, filming, and exporting to an internal web server, CD, or DVD.

AW VolumeShare 7 is also available for purchase without Volume Viewer for users who require only 2D image applications. Volume Viewer dependent features will not be supported with this configuration.

Modalities

Supported DICOM modalities include:

- Computed tomography (CT),
- Magnetic resonance imaging (MR),
- Radiofluoroscopy (RF),
- Xray angiography (XA),
- Computed radiography (CR),
- Digital radiography (DX),
- Mammography (MG),

- Nuclear medicine (NM),
- Positron emission tomography (PET),
- Ultrasound (US)
- Structured Reporting (SR),
- Key Objects (KO)

License Management

- Application Usage Monitor tracks how often a particular license was used, when, and by whom. You can compile and view usage reports for a specific period of time. This helps you understand how your systems contribute to departmental productivity and to plan for future license requirements based on actual usage data.
- For floating licenses, a secondary license server can be assigned for redundancy. If the primary license server becomes unavailable, the secondary can be contacted so that work can continue.

User Interface

- Displays on one or two 1280x1024 monitors with scroll or optical mouse.
- Simple drag-and-drop mechanism for networking, media interchange, filming and data export options.
- Automatic help message display briefly describes function currently under the mouse pointer.
- Shortcut keys and programmable window/level function keys for accelerated control.
- Multi-tasking capabilities that allow Patient List, 2D Viewer, Filmer and Volume Viewer to run simultaneously with Fast Switch capabilities between applications.

CONFIDENTIAL

SECRET DE
AFACERI

Standard Applications & Capabilities

Separate application licenses must be purchased to access advanced applications functionality.

2D Viewer

The 2D Viewer is an application used to display, manipulate, annotate and review 2D images by a trained physician for diagnostic interpretations.

Display Customization

Display customization allows the user to manage layouts to display data. Key benefits include:

- Allows modification of number of exams or series displayed to facilitate either single or multi-exam reviews.
- Standard layouts provide flexibility to tailor the image display from 1x1 to 8x8. The current layout is retained if the images are from the same modality and orientation.
- Annotation levels allow selection of the image information fields to display.
- Toolbar customization enables control of which buttons are displayed based on user preference.

Study Navigation

- The navigator lets users assign a series to a view on the fly.
- Cine mode also provides temporal, spatial, or manual playback loops.
- With two exams or series loaded, Cine mode supports a side-by-side display format with synchronized playback loops for more efficient comparisons.

Image Review

The initial image window and level setting is based on the DICOM header. Once displayed, several methods of adjusting image window and level are provided.

Provides routine image manipulation features:

- Flip/Rotate, Zoom, Pan, Magnifying Glass,
- Inverse Video: inverts grayscale color map,
- Display normal: lets you return the image to its default viewing parameters.

You can access the following features using a single mouse click directly on the image. Having these direct manipulation tools easily at hand gives you fast access with less distraction from your review task by eliminating the need to return to the graphical user interface controls: 2D distance, Angle, Report cursor, Box ROI, elliptical ROI and free-hand ROI.

The program continuously updates statistics on the fly.

- The Annotation feature lets you highlight areas of interest by adding text and line/arrow pointers to any image feature.
- The Copy/Paste/Erase feature lets you copy/paste/remove any text or graphic placed on an image.
- Cross-reference indicates the position of the current slice over the scout or localizer or any other non-parallel series.
- Save lets you store a copy of the image as it appears on the screen for future review.
- The Key Images feature allows you to flag images and create key objects.

Advanced X-Ray Analysis

The user can apply shutters to X-ray images to focus on specific areas within the image. Shutters are black, opaque overlays with elliptical or

CONFIDENTIAL

SECRET DE
AFACERI

rectangular cutouts that you can size as desired. Once a shutter is applied, you can move the underlying image to shift the viewable area. Shutters can also be applied to CT and MR images with the Image Matte feature.

The user can apply several levels of edge enhancement filters to an image.

- **Image Subtraction:** A graphical user interface lets you select a mask and subtraction of all images in an associated sequence. Subtracted images may be saved as a sequence with one mouse click.
- **Landscape:** With this tool, you can introduce a percentage of the mask image into the subtracted image for anatomical reference. A graphical user interface lets you specify the percentage of mask to be reintroduced.
- **Pixel Shift:** Auto pixel shift optimizes a mask to image registration. You can also move a mask manually to optimize subtraction quality in a specific region of interest.
- **Split Pixel Shift:** You can split the screen horizontally or vertically for pixel shift.
- **Maximum/Minimum Opacifications:** You can integrate selected images to provide a resulting Max. Op./Min. Op. image.

Batch Filming

The Print Series feature lets users automatically batch film an entire series with a single keystroke.

Volume Viewer

Volume Viewer is a multi-modality AW Family advanced visualization application. It provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-Ray, PET, PET/MR, and PET/CT datasets. It

features a breadth of high-performance analysis tools such as Multiplanar Reformat, Curvilinear Reformatting, Volume Rendering and Lumen Navigation.

Volume Viewer is the pre-requisite application for what are known as Volume Viewer Foundation (VVF) applications within the AW Family Advanced Visualization Application portfolio. VVF applications share a similar user interface and have some interoperability. For detailed description of Volume Viewer features, please refer to the Volume Viewer Product Data Sheet.

Filmer

The integrated Filmer enhances the efficiency of the review station and gives users greater filming and exporting flexibility. This feature supports two modes: the Mini Filmer mode and Full-Screen mode, which provides the ability for customization and film layout template creation. Three key mechanisms give users flexibility:

- Free format filming
- DICOM structured reporting (SR)
- Data export (HTML/PDF and JPEG, PNG, MPEG, AVI, or QTVR)

With the Filmer, users can easily extract significant images from any AW Family AV application (2D Viewer, Volume Viewer, Ready View, etc.).

- You can film images individually by dragging and dropping to the on-screen Filmer, or by a single F1 keystroke.
- Multiple image formatting allows filming multiple images in a single page frame with the F2 keystroke.
- With Film MID users can send multiple images to a single Filmer frame.

CONFIDENTIAL

SECRET DE
AFACERI

- Batch filming is supported by applications that provide that capability (Print Series in the 2D Viewer, Batch Film Protocols in Volume Viewer).

Communication between applications and the Filmer is accomplished in the Mini Filmer mode, which provides a minimized footprint. The Mini Filmer mode provides the following:

- Compatibility with Batch Filming from Volume Viewer (may be optional in your configuration)
- Store and position images transferred from an application.
- Type of export (film, media, database).
- Rapid switching between Full-Screen Filmer and application.

The resultant electronic films become a quick summary of the patient study and radiological interpretation which can then be reviewed by clinicians and physicians. They can also be saved as independent files (DICOM SR and Secondary Captures) for teaching purposes. Electronic film can contain one or several pages with specific layouts for each page or all pages.

A flexible Edit Mode provides the ability to easily add, manipulate, format, or delete images from the film. Images can contain text and graphics from measurements and your annotations, and may be window/leveled, magnified, flipped, rotated, or cine. Additional annotations can be added to the image in edit mode.

Preview Mode displays the film as it will be printed or exported, considering the layout applied to each page and the compression level specified for non-DICOM exporting (JPEG/PNG and MPEG/AVI/QTVR).

Printing

AW VolumeShare 7 includes Network DICOM Print (B&W and color). For supported printers, AW VolumeShare 7 includes Network PostScript capability (B&W and color). For PostScript and DICOM printers, any printing format created in the Filmer is supported (e.g., non-square matrix formats for rectangular images such as CT run-offs).

- AW VolumeShare 7 postscript printing has been tested on the following devices: Codonics 1660M, 1660MD or Horizon, Lexmark Optra 1650N, 1855N, SC1275N, C710N, C720N, T612 or T614, Seiko 1720D, Kodak DMI3600, Quantum GL2101HD, declared with film/thick paper or with plain paper, Tally T8106, HP LaserJet and Xerox Phaser
- Digital cameras and analog cameras using a 3M-952 protocol (including DASM interface) are not supported on AW VolumeShare 7.

Exporting

Data export is integrated in the Filmer, thus providing all image processing tools needed for multimedia image export. The Filmer exports any electronic film to DICOM SR, PDF/HTML, JPEG, PNG, MPEG, AVI or QTVR format. AW VolumeShare 7 enhanced security of the export functionality with the use of SFTP connectivity. Non-DICOM data can also be saved on a multi-session CD/DVD. The Data Export capability is intended only for publishing and communication, not for diagnostic purposes. Its simplicity is reflected in the different export mechanisms available:

- CD/DVD removable media
- USB flash drive
- Network HTTP and FTP protocols

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SECRET DE AFACERI

Cardiac Review and Export

Processing and review of CT, MR and PET cardiac exams with manual oblique reformatted protocols can be exported as a multi-phase Cine movie that allows referring physicians to review exams in a dynamic mode.

Workflow Management

AW VolumeShare 7 provides features designed to optimize your reading workflow:

- The Patient List provides tools that enable you to sort and filter imaging studies stored on the local workstation or on remote systems.
- Postfetch allows you to retrieve a patient's prior related DICOM exams from a remote DICOM host. Retrieval can be triggered by the arrival of a new patient study on the workstation from the network or supported media.
- Support for IHE Key Image Notes profile, allowing you to mark certain images in the 2D or 3D viewer as key images. Key images are displayed as a separate series in the exam list and accessed or sent to an IHE KIN-supporting PACS. Specific tags and Key Image Notes may be created using the 2D Viewer application.
- Support for external USB disks to serve as DICOM storage media. The USB disk can be accessed for reading and writing via the media button, the same way you would access a CD or DVD. The maximum number of images that can be stored on a USB device is restricted to 300,000 across all exams. USB media must be formatted using the FAT32 file system. See the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual for more details on USB storage.

- Preference sharing among users. Sharable preferences include: Filmer layouts, Volume Viewer custom protocols, and viewer preferences. A standard USB flash drive may also be used to share preferences between workstations.
- One-Touch protocols enable you to define an application or preset protocol to launch automatically based on DICOM elements.
- DICOM Query/Retrieve Storage Class User (SCU) and Storage Class Provider (SCP) provide seamless network integration.
- DICOM Storage Commitment SCU lets you know when exams have been archived successfully on DICOM devices, such as PACS, that support Storage Commitment SCP.

End Review

End Review automates the routine tasks required at the end of reviewing each exam. The "End Review" flag in the patient list allows you to mark exams as "Done" after post-processing has been completed. End Review allows you to automatically perform one or several of these actions:

- Print pages prepared in the Filmer to the default printer (DICOM or Postscript filmers supported), with the option to automatically clear the Filmer after printing.
- Save electronic films in the workstation's attached DICOM database.
- Automatically push the entire exam or only the series created on the AW to one or more remote hosts.

Remote Network Host Management

AW VolumeShare 7 Patient List provides a control panel to select a remote host or destination for

network transfer. Icons signify different device types (acquisition, post processing, PACS, etc.) that are accessible to the workstation via the network.

You can initiate a network transfer simply by dragging an exam, patient, series, or image(s) to the remote host icon. You can activate a remote browser to display detailed exam information from the remote host, including filtering information if supported.

Search Advantage

This advanced search engine lets you search a patient's exam history on PACS or other DICOM device with a few mouse clicks.

Worklists and Filters

You can filter the patient list by modality, date, end review status, or exam description. Most filters are available on the remote host patient list as well.

You can filter the patient list further by choosing one or any of these parameters:

- Modality.
- Patient name or patient ID.
- Exam location (hospital name) and exam description.
- Series description
- Date and time of day: today, or specified date or date ranges with a specified time or time ranges.
- Radiologist's or referring physician's name.
- Accession number.

Quick Sort (ascending or descending) and Quick Access (entry field) are available for any of the fields displayed at exam level (patient name, patient ID, exam location and description,

radiologist's and referring physician's name, date and time).

Queues

AW VolumeShare 7 manages three queues: network transfer; media interchange, filming. A menu on the patient list gives you easy access to queue status, and an animated icon gives you continuous network queue status information.

DICOM CD/DVD/USB Creation Tool

- The DICOM media creation tool offers you the flexibility to save and retrieve from supported CD, DVD, or USB storage devices.
- Add selected exam, series, or images by dragging and dropping to the pre-mastering window.
- CD and DVD usage percent is shown before you launch the writing process.
- You can optionally down sample certain 1024²x512² X-ray angiographic images during CD creation.
- CD/DVD composer lets you manage large amounts of data. If data size is larger than medium capacity, you are prompted to supply additional media storage space.
- You can specify the number of copies created during a CD/DVD save session.
- You can embed a DICOM viewer on DICOM CD/DVD so data can be reviewed on a PC running a Windows operating system.
- You get lossless JPEG reading of any DICOM CD/DVD media. Lossless writing is available for X-ray exams.
- CD/DVD drives operate at 16x write speeds and read speeds up to 48x depending on choice of media.

CONFIDENTIAL

SECRET DE
AFACERI

- For limitations on USB storage, refer to the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual.

Database Management

- The database management system classifies data according to patient folder description of the DICOM standard: Patient/Exam; Study/Series; Sequence/Images.
- Automatically deletes images on first-in first-out (FIFO) basis, which can be toggled on/off.
- Lock Exam feature protects specified exams from deletion.
- Permanent display of available disk space facilitates disk management.
- Study anonymization tool modifies DICOM elements to remove protected health information.

Security Capabilities

- Full support of Enterprise Authentication with Microsoft® Active Directory® and Novell® eDirectory™. Users can access their AW Workstation with existing Single Sign On credentials. This may streamline HIPAA compliance.
- If Enterprise Authentication is not used, an administrator-defined UNIX login and password is required per workstation. A lock screen feature enables you to lock the workstation display when leaving it momentarily, preventing unauthorized access to patient data.
- Main actions performed on the AW (save, film, network, export) are flagged and saved to an audit trail which includes username, ID, exam ID, date, and time.
- Multiple users can have different local UNIX logins, use the same AW workstation, and yet keep their user preferences unique.

- AW VolumeShare 7 is secured by Product Network Filters (PNF), which keeps unauthorized users from accessing open ports over the network. You can configure this feature to provide only specific remote devices access to AW ports.
- While not recommended, AW VolumeShare 7 does support ClamAV® anti-virus protection. This program is not activated by default in order to guarantee optimal performance of the AW workstation. If needed, your GE Healthcare field engineer can activate the program. However, it is the customer's responsibility to provide a safe internet connection from the workstation and to keep virus definitions current.
- AW VolumeShare 7 provides de-identification and encryption capabilities to limit privacy risks to sensitive information. The patient data exported during clinical workflow may be encrypted in flight by using either the SFTP or DICOM TLS protocol depending on the required dataflow.

System Components

Workstation Configuration

- HP Z4 G4 Workstation
- SLES 15 operating system
- Intel® Xeon® W-2245 Eight physical core/ sixteen logical cores 3.9 GHz CPU
- 64GB (4x16GB) DDR4 2933 MHz or higher ECC Registered DIMM
- NVIDIA Quadro NVS P620, 2GB Graphics card
- 1 x 512GB M.2 NVMe Drive for OS and Applications
- 1 x 1TB M.2 NVMe Drive for image cache.
- Cache, which is subject to overhead can store approximately:
 - 6,400,000 256*256*16 bits uncompressed images OR

CONFIDENTIAL

SECRET DE
AFACERI

- 1,600,000 512*512*16 bits uncompressed images OR
- 400,000 1024*1024*16 bits uncompressed images OR
- 100,000 2048*2048*16 bits uncompressed images
- Internal DVD Writer drive for read/write of DICOM CD/DVD media, read/write of Data Export CD/DVD data and service use (DVD Install)
- Integrated dual Ethernet 10/100/1000 Mbit/s ports.
- Additional single Ethernet 10/100/1000 Mbit/s port
- 1 USB QWERTY (or regional) Keyboard and Mouse

Footprint

- Height 38.60 cm (15.2 in.)
- Width 16.89 cm (6.65 in.)
- Depth 44.47 cm (17.5 in.)
- Approximate Weight 17 kg (38 lbs.)

Operating Environment

- Temperature: +5°C to +35°C (40° to 95° F)
- Humidity: 10% to 85% (relative non-condensing)
- Altitude: 0 to 5000m (16,404 ft.)
- Acoustics: LWAd less than 4.3 Bels
- Shock: 40 G peak, half-sine, 2-3 ms

Non-Operating Environment

- Temperature: -40°C to +60°C (-40° to 140° F)
- Humidity: 10% to 90% (relative, non-condensing)
- Altitude: 0 to 12192m (37,000 ft.)

Monitors

- (2) 19" color Flat Panel LCD monitors

- DICOM Part 14 factory calibrated
- Native resolution 1280 x 1024 (5:4 aspect ratio)
- Approximate Weight: 5.3 Kg (11.68 lbs.).
- AC 100 - 240 V : 50 / 60 Hz

Image Networking

- Standard 10/100/1000 Base-T Ethernet for DICOM
- 1000 Base-T dedicated network for optimal Direct Connect performance
- Protocols supported:
 - DICOM 3.0 Storage SCU/SCP and Query/Retrieve SCU/SCP
 - RSVP
 - TCP/IP network layer
- SdCNet is no longer supported
- AW VolumeShare 7 does not support the AdvantageNET network protocol.
- AW VolumeShare 7 does not support DICOM images from GE Healthcare Signa™ version 5.x 1.5T MR systems
- AW VolumeShare 7 software is supported on the previous generation AW VolumeShare 5 HP Z820 and Z440 workstations.

CONFIDENTIAL

SECRET DE
AFACERI

GE HealthCare

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](https://www.gehealthcare.com) for more information.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

Adobe is either a registered trademark or trademark of Adobe Systems Incorporated in the United States and/or other countries.

ClamAV is a registered trademark of Cisco Systems.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Intel and Xeon, are either registered trademarks or trademarks of Intel Corporation in the United States and/or other countries.

Microsoft, Active Directory, and Windows, are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Novell and eDirectory are registered trademarks of Novell, Inc. in the United States and other countries.

Oracle and Java are registered trademarks of Oracle and/or its affiliates.

All other product names and logos are trademarks or registered trademarks of their respective companies.

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
©2023 GE HealthCare.

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SECRET DE
AFACERI

VessellIQ™ Xpress with AutoBone™ Xpress

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Overview

VessellIQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

CONFIDENTIAL

SECRET DE AFACERI

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vessellq-xpress>

Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- **Fast tracking detects entire length of vessel and displays the centerline based on a single start point.**
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- **Thrombus detection with measurements and comparison over time.**
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.

- Save State lets you save and restore the state of the processed images at any stage

System Requirements

- AW Workstation AW 4.7 ext 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Indications for Use

VessellQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. VessellQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels. Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing;

density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

CONFIDENTIAL

SECRET DE
AFACERI

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, imagination at work, AutoBone Xpress, VessellQ Xpress are
trademarks of General Electric Company.
All other third-party trademarks are the property of their respective owners
Doc1471990

Features Detail

Anatomy-based protocols

- The program offers you predefined protocols for fast, easy review of all vascular structures:
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail

- Automatic calcium removal for vessel lumen visualization

Auto Aorta

- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.

- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.


Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.

- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
 - Minimum diameter
 - Maximum diameter
 - Mean diameter
 - Cross-sectional area

Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.


Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.


CONFIDENTIAL

SECRET DE
AFACERI

Advanced Vessel Analysis Tools

- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.


- Available measurements include:
 - Distance
 - Cross-sectional area

- Single or dual reference point comparison
- Relative percent stenosis / aneurysm
- Color Identification
- 3D Angle
- Thrombus Extraction
- Volume
- Mean, min, max diameter

- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.

- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.

Summary Table & Report

- Measurements are automatically captured in summary table for easy review.


- Simply click on any measurement in the summary table and you will be taken to that measurement location.

CONFIDENTIAL

SECRET DE AFACERI

- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, film or PDF.

Save State & Protocol

- Save the current state of tracking and measurements for easy reloading with a single mouse click.
Save state captures:
 - 3D images
 - Vessel tracking
 - Measurements
 - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VesselIQ Xpress provides user-friendly tools an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.

Multiple Study Comparison

- With VesselIQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

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
**SECRET DE
AFACERI**

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TAVI Analysis

Automated segmentation & pre-procedure planning software for transcatheter aortic valve implantation.

Transcatheter Aortic Valve Implantation (TAVI) or Replacement (TAVR) is a procedure for select patients with severe symptomatic aortic stenosis (narrowing of the aortic valve opening). TAVI/TAVR is performed on a beating heart and a replacement valve is delivered via a catheter.

CT provides important information that is important for successful TAVI/TAVR procedures. CT is used to help determine aortic annulus size, to guide selection of appropriate valve, provide dimensions of the entire aorta, help determine the access path for the catheter and give guidance for C-arm angulation for deployment of the device.

GE's TAVI Analysis software provides a streamlined, guided workflow to enable efficient, consistent work-ups of your TAVI studies with connectivity directly to the interventional suite.

Overview

TAVI Analysis is a post-processing software package for the Advantage Workstation (AW) and AW Server platforms. It is a planning tool used for TAVI/TAVR procedures. It automatically segments the aorta and displays the aortic valve in multiple views for quick and easy measurements of the annulus. TAVI Analysis provides guided workflow and semi-automated tools to help evaluate appropriate access pathways and can communicate directly with the interventional suite.

Highlights

- Automated Aortic segmentation displayed to highlight calcific areas.
- Guided workflow for complete pre-procedure planning measurements.
- Wizard panel for valve measurements used for sizing requirements.
- Perpendicular view demonstrates working angles for valve deployment.
- Auto bone segmentation
- Fast centerline tracking to subclavian and femoral arteries.
- Lumen and curved views showing entire vascular tree.
- Quantification wizard panel.
- Intuitive and easily accessible tools
- Connectivity to the interventional suite

CONFIDENTIAL

SECRET DE AFACERI

Visit us:

www.gehealthcare.com/aw/applications/

Features

- One comprehensive workflow for valve measurements and access path measurements.
- Automatic segmentation of transparent volume rendered model of the aorta with calcifications highlighted
- Auto detection of valve plane.
- Intuitive workflow for acquiring the necessary measurements for aortic valve sizing
- One click perpendicular view provides guidance for C-arm angulation for deployment of device.
- Automatic bone segmentation for the chest, abdomen and pelvis.
- Real time vessel tracking provides fast centerline detection based on a single start point.
- Guided vessel tracking tools allow for easy access planning for all transcatheter approaches:
 - Subclavian Artery
 - Femoral Artery
 - Apical Approach
- Peripheral artery measurements automatically sent to summary table to allow easy exporting.
- Save State and send to the Cath lab for sharing access to models and measurements.

Indications for Use

TAVI Analysis for CT

TAVI Analysis is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

TAVI Analysis is a post processing application option for the Advantage Workstation (AW) platform, CT Scanner or PACS stations, which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and directional vessel tortuosity visualization.

TAVI Analysis automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel, analyzing calcified and non-calcified plaque to determine the densities of plaque within a coronary artery, measure areas of abnormalities within a vessel.

TAVI Analysis for X-Ray

TAVI Analysis is a software post processing option for the Advantage Workstation (AW) Platform, which can be used in the analysis of the 3D angiography data. It provides a number of display, measurement and batch filming/archive features and will aid Physicians in studying user selected vessels for stenosis analysis, pre & post stent planning, and directional vessel tortuosity visualization.

System Requirements

Minimum platform release:

- AW Workstations VolumeShare 7 and AW Server 3.1 or later.
- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.


CONFIDENTIAL

SECRET DE
AFACERI

Features Detail


Valve Definition

- The TAVI protocol allows users to select single or multiple cardiac phases.
- Automatic segmentation of aorta and auto detection of valve plane with double oblique tool.
- Ability to define minimum and maximum diameters above the annulus along the length of the aorta.


Aorta Analysis


- Automatic bone segmentation and real time vessel tracking.
- Guided vessel tracking tools allow for easy access planning for all transcatheter approaches:
 - Subclavian Artery
 - Femoral Artery
 - Apical Approach

- Analyze, trace and edit annulus contour. Available annulus measurements include:
 - Minimum diameter
 - Maximum diameter
 - Perimeter
 - Area
 - Equiv. diameter based on perimeter and area

- Aortic measurement tools are also available which include:
 - Diameters above annulus
 - Ostia Heights
 - Sinus diameters
- One click perpendicular view demonstrates working angles for valve deployment in the cath lab.

- Peripheral artery measurements automatically sent to summary table to allow easy exporting.

Aortic Measurements	
C17 phase 46%	Min Diameter: 20.0 ± 0.0 mm
Ascending Aorta - Distal	
C17 phase 46%	Min Diameter: 34.1 ± 0.7 mm
Ascending Aorta - Proximal	
C17 phase 46%	Min Diameter: 25.4 ± 0.9 mm
Subclavian Artery - Origin	
C17 phase 46%	Min Diameter: 27.1 ± 0.8 mm
Subclavian Artery - Proximal	
C17 phase 46%	Min Diameter: 24.1 ± 0.9 mm
Transcatheter Aorta - Cross-section	
C17 phase 46%	Min Diameter: 19.8 ± 0.8 mm

- Summary table converts to PDF format for a simple and concise way to share measurements and images with heart team for planning discussions.

Procedure Preparation

- Provides valve measurements, perpendicular view and gantry angulation in the cath lab to prepare for TAVI/TAVR procedure.
- Direct communication with Heart Vision software for easy transition from planning to guidance.
- 3D and calcium overlay VR models to aide in visualization during procedure.

CONFIDENTIAL

SECRET DE AFACERI

Indications for Use

TAVI Analysis for CT

TAVI Analysis is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

TAVI Analysis is a post processing application option for the Advantage Workstation (AW) platform, CT Scanner or PACS stations, which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and directional vessel tortuosity visualization.

TAVI Analysis automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel, analyzing calcified and non-calcified plaque to determine the densities of plaque within a coronary artery, measure areas of abnormalities within a vessel.

TAVI Analysis for X-Ray

TAVI Analysis is a software post processing option for the Advantage Workstation (AW) Platform, which can be used in the analysis of the 3D angiography data. It provides a number of display, measurement and batch filming/archive features and will aid Physicians in studying user selected vessels for stenosis analysis, pre & post stent planning, and directional vessel tortuosity visualization.

System Requirements

Minimum platform release:

- AW Workstations VolumeShare 7 and AW Server 3.1 or later.
- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

Regulatory Compliance


This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

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Vizitați-ne:

www.gehealthcare.com/awi/applications/

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