# Download Instructions for Capillary Electrophoresis Dx Consumables Product Documents

Pub. No. MAN0026540 Rev. D.0

## Download the instructions for use (IFU)

The IFU for this product is provided in a Thermo Fisher™ Connect Platform resource library.

1. Go to [apps.thermofisher.com/apps/spa/#/publiclib/cedx](https://apps.thermofisher.com/apps/spa/#/publiclib/cedx), then sign in to thermofisher.com.  
**Note:** If you do not have a thermofisher.com account, click **Create Account** in the sign in screen, then follow the instructions.
2. In the resource library, locate your product, then download and extract the ZIP folder containing the IFU.
3. Open the instructions for use with a PDF viewer such as Adobe Reader™ (<https://get.adobe.com/reader/>).

[thermofisher.com/support](https://thermofisher.com/support) | [thermofisher.com/askaquestion](https://thermofisher.com/askaquestion)

Life Technologies Holdings Pte Ltd | Block 33 | Marsiling Industrial Estate Road 3 | #07-06, Singapore 739256

Revision history: Pub. No. MAN0026540 D.0

Revision	Date	Description
D.0	7 March 2023	Translations were added.
C.0	6 January 2023	The manufacturing site was changed.
B.0	29 July 2022	Translations were added.
A.0	21 April 2022	New electronic documentation instructions for accessing instructions for use for capillary electrophoresis consumables.

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©2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

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1. Отворете [apps.thermofisher.com/apps/spa/#/publiclib/cedx](https://apps.thermofisher.com/apps/spa/#/publiclib/cedx), след което влезте в thermofisher.com.

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L512JA309

2013-10-05

TEST	SPECIFICATION	RESULT
<b>Capillary Leak Test</b> Each array is tested to assure absence of leaks.		
Leak Test	No Leaks	PASS
<b>Capillary Plug Test</b> Each array is tested to assure absence of plugged capillaries.		
Plug Test	No Plugged Capillaries	PASS
<b>Fluorescence Background Test</b> Each array is tested to assure background fluorescence is within specification.		
Fluorescence Background Test	Conforms to Specification	PASS
<b>Capillary Mapping Test</b> Each array is tested to verify that capillary inlet and outlet spatial positioning correspond.		
Mapping Test	Conforms to Specification	PASS

*When used with the 3500 Dx series systems, this device meets the requirements of Directive 98/79/EC on In Vitro Diagnostic Medical Devices and complies with QSR 21CFR820.*

**ISO13485:2003  
 REGISTERED**

Singapore

**IVD**

Quality Assurance Manager

Life Technologies Holdings Pte Ltd, Block 33, Marsiling Industrial Estate Road 3, #07-06,  
 Singapore 739256 Tel: (65) 6362 9300

For more information go to [www.lifetechnologies.com](http://www.lifetechnologies.com) or email us at [cofarequest@lifetech.com](mailto:cofarequest@lifetech.com).


# NeoNat SCID SMA

---

## REAL TIME PCR KIT

The NeoNat SCID SMA is a Real-Time PCR assay for screening Severe Combined Immunodeficiency (SCID) Syndrome and Spinal Muscular Atrophy (SMA) in newborns DNA from Dried Blood Spot (DBS) samples.

-  Greater flexibility with separate multiplexing of SCID and SMA from one kit
-  Kit contains all the reagents and components to perform DNA extraction from DBS
-  Ready to use reaction mix and calibrators preloaded plates enable convenience for users
-  Four levels of DBS controls with defined target gene copies per  $\mu\text{L}$
-  Compatible with multiple qPCR instruments
-  Short turnaround time of ~120 mins from DNA extraction to result analysis

## Assay Description

- NeoNat SCID SMA Real time PCR kit is based on 5' nuclease technique
- Screens T-cell receptor excision circle (TREC) / Kappa deleting recombination excision circle (KREC) for SCID and Survival Motor Neuron 1 (SMN1) & Survival Motor Neuron 2 (SMN2) for SMA. The  $\beta$ -globin gene serves as an internal control


## Performance

- The assay demonstrates excellent performance with 100% sensitivity and specificity for all the samples tested
- The analytical sensitivity shows a limit of detection of <5.0 copies/ $\mu$ L for TREC & KREC and <2.5 copies/ $\mu$ L for SMN1 & SMN2

## Compatibility

- Kit compatible with most real-time PCR instruments in the market with at least three measurement channels (FAM, VIC/or HEX, & CY5)
- Developed and validated with Biorad CFX 96 and ThermoFisher Quantstudio 5 & 6.

### DNA EXTRACTION FROM DBS & PCR REACTION

## Ordering Information

96 Reactions	8100481
192 Reactions*	8100482
480 Reactions*	8100483

\*Please check the availability

For more information please order your IFU from  
**Labsystems Diagnostics Oy**

Tiilitie 3, FI-01720 VANTAA, Finland Tel: +358 (0) 20 155 7530  
 sales@labsystemsdx.com | www.labsystemsdx.com

## Real-Time PCR Systems Spectral Calibration Kit I, 96-Well

Product No. **4349180**  
Lot No. **2305202**

TEST	SPECIFICATION			RESULT	
<b>Material Test</b>					
Applied Biosystems performs spectrofluorimetric analysis using a calibrated PerkinElmer LS55 Fluorescence Spectrometer to test fluorescence emission wavelength maximum in each lot of component bulk material used in the Spectral Calibration Kit. The PerkinElmer LS55 is calibrated using a mercury arc lamp which verifies the emission monochromator wavelength accuracy and then the excitation wavelength accuracy is verified against the emission wavelength as a reference.					
<b>Component</b>	<b>Part Number</b>	<b>Lot Number</b>			
SYBR® GREEN Dye	4349763	2305278	522 ± 5 nm	<b>520 nm</b>	<b>Pass</b>
ROX™ Dye	4349411	2305328	601 ± 5 nm	<b>603 nm</b>	<b>Pass</b>
TAMRA™ Dye	4349410	2305326	578 ± 4 nm	<b>578 nm</b>	<b>Pass</b>
NED™ Dye	4349408	2305336	573 ± 4nm	<b>572 nm</b>	<b>Pass</b>
JOE™ Dye	4349409	2305329	548 ± 4 nm	<b>550 nm</b>	<b>Pass</b>
VIC® Dye	4349764	2305313	548 ± 4 nm	<b>548 nm</b>	<b>Pass</b>
FAM™ Dye	4349762	2305307	518 ± 4 nm	<b>519 nm</b>	<b>Pass</b>
ROI™ Dye	4349415	2305397	Visual Inspection		<b>Pass</b>
Background	4330124	2305417	Visual Inspection		<b>Pass</b>

*For Research Use Only. Not for use in diagnostic procedures.*

Singapore

*Tan Chuen Fang* **13 JUN 2023**  
**Quality Assurance**  
(Name, Signature & Date)

Life Technologies Holdings Pte Ltd, Block 33, Marsiling Industrial Estate Road 3, #07-06,  
Singapore 739256 Tel: (65) 6362 9300

To obtain a Certificate of Analysis on-line go to [www.thermofisher.com](http://www.thermofisher.com) or email us at  
[QA-SG\\_COA\\_Request@lifetech.com](mailto:QA-SG_COA_Request@lifetech.com)

# 7500 Real-Time PCR Systems Spectral Calibration Kit I

Catalog Number 4349180

Pub. No. 4350071 Rev. D

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

## Contents and storage

Contents	Amount	Storage
Background Plate sealed with an optical cover	1	-25°C to -15°C
Spectral Calibration Plates sealed with optical covers	7	
Region of Interest (ROI) Calibration Plate sealed with an optical cover	1	

## Related Documentation

For detailed information on instrument setup and the calibration process, refer to the *Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System Installation and Maintenance Guide* (Pub. no. 4347828).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

For support visit [thermofisher.com/support](http://thermofisher.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

The information in this guide is subject to change without notice.

### DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

**For Research Use Only. Not for use in diagnostic procedures.**

## CERTIFICATE OF ANALYSIS

**A46109 PowerTrack™ SYBR™ Green Master Mix**

**Packaging Lot: 2746887**

Expiry Date: 30.04.2025 (DD.MM.YYYY)

Storage: at -20±5°C in the dark

Note: **For Research Use Only. Not for use in diagnostic procedures.**

### Filling lots for components in package:

Lot	Quantity	Description
2728794	5 mL	PowerTrack™ SYBR™ Green Master Mix
2737196	1.25 mL	40X Yellow Sample Buffer

### QUALITY CONTROL

Parameter	Method	Requirement	Result
Functional Test	The product is functionally tested by qPCR analysis. It must demonstrate functional performance using a target concentration of plasmid DNA.	Average Ct for the target concentration is between 20 and 23	Pass
dNTP Concentrations	Determined by analytical method.	Within range of target concentration	Pass
Mg <sup>2+</sup> Concentration	Determined by analytical method.	Within range of target concentration	Pass
K <sup>+</sup> Concentration	Determined by analytical method.	Within range of target concentration	Pass
RNase Activity	Determined by analytical method.	No detectable activity level	Pass
DNase Activity	Determined by analytical method.	No detectable activity level	Pass
E. coli DNA Level	Determined by analytical method.	No detectable level	Pass

#### ISO CERTIFICATION

Manufactured by Thermo Fisher Scientific Baltics UAB, in compliance with ISO 9001 and ISO 13485 certified quality management system.

Quality authorized by QC: **J. Žilinskienė**

# PowerTrack™ SYBR™ Green Master Mix

## USER GUIDE

Master mix with a two-dye tracking system for real-time PCR workflows

**Catalog Numbers** A46012, A46109, A46110, A46111, A46112, A46113

**Publication Number** MAN0018825

**Revision** B.0

Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](https://www.thermofisher.com/symbols-definition).

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Revision history:** Pub. No. MAN0018825

Revision	Date	Description
B.0	29 July 2022	The volumes for preparing the PCR reactions were corrected (Table 2 on page 12 and Table 3 on page 13).
A.0	30 January 2020	New document.

**Important Licensing Information:** This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

**TRADEMARKS:** All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2022 Thermo Fisher Scientific Inc. All rights reserved.

# Contents

■	<b>CHAPTER 1</b>	<b>Product information</b>	5
		Product description	5
		Contents and storage	6
		Required materials	7
		Workflow	9
■	<b>CHAPTER 2</b>	<b>Methods</b>	10
		Guidelines	10
		Requirements for input DNA	10
		Guidelines for PCR reactions	10
		Guidelines for no-template control reactions	10
		Before you begin	10
		Set up the plate document or plate file	10
		Prepare the reagents	11
		Prepare the PCR reactions	11
		Set up and run the real-time PCR instrument	14
		Analyze the results	15
■	<b>APPENDIX A</b>	<b>Troubleshooting</b>	16
		General troubleshooting	16
		Troubleshooting baseline settings	18
		Troubleshooting threshold settings	19
■	<b>APPENDIX B</b>	<b>Background information</b>	21
		Master mix components	21
		Overview of the chemistry	22
		Two-step RT-PCR	23
		A typical amplification plot	24
		Baseline and threshold values	24
		Melt curves	25

- **APPENDIX C** Template quality and quantity ..... 26
  - DNA template quality ..... 26
  - RNA guidelines ..... 26
  - Template quantitation using O.D. 260 ..... 27
  - Template storage ..... 27
  
- **APPENDIX D** Primer design, target sequences, and optimizing primer concentration ..... 28
  - Primer design guidelines ..... 28
    - Avoiding primer-dimers ..... 29
  - Identifying target sequence and amplicon size ..... 29
  - Selecting an amplicon site for cDNA ..... 29
  - Optimize primer concentrations for PCR ..... 30
    - Overview ..... 30
    - Quantitate primers ..... 30
    - Determine the optimal primer concentration for PCR ..... 31
    - Confirm the absence of nonspecific amplification ..... 31
  
- **APPENDIX E** Safety ..... 33
  - Chemical safety ..... 34
  - Biological hazard safety ..... 35
  
- **APPENDIX F** Documentation and support ..... 36
  - Related documentation ..... 36
  - Customer and technical support ..... 36
  - Limited product warranty ..... 36

# Product information

■ Product description .....	5
■ Contents and storage .....	6
■ Required materials .....	7
■ Workflow .....	9

---

**IMPORTANT!** Before using this product, read and understand the information in the “Safety” appendix in this document.

---

## Product description

The Applied Biosystems™ PowerTrack™ SYBR™ Green Master Mix is a two-dye tracking system for real-time PCR.

It contains a visual dye-based indicator that allows the user to confirm that the sample was added to the master mix.

The master mix contains an inert blue dye. The optional Yellow Sample Buffer is added to the sample to track the sample through the PCR. The reaction mix turns green due to the inert blue dye and the Yellow Sample Buffer.

The SYBR™ Green dye binds to double-stranded DNA formed during real-time PCR. For more information, see “Overview of the chemistry” on page 22.

The master mix provides flexibility for varying primer melting temperatures ( $T_m$ ). The  $T_m$  can be between 55°C and 65°C. For recommended primer concentrations for different  $T_m$ , see “Guidelines” on page 10.

The master mix contains the following components:

- SYBR™ Green dye
- DNA polymerase, with an antibody-mediated hot start
- Heat-labile Uracil-DNA Glycosylase (UDG)
- ROX™ dye (passive reference dye)
- dNTP blend containing dUTP/dTTP
- Inert blue dye
- Optimized buffer components

For more information about each component, see “Master mix components” on page 21.

The user provides primers, template, and water.

## Contents and storage

Table 1 PowerTrack™ SYBR™ Green Master Mix

Cat. No. <sup>[1]</sup>	Amount	Storage
<a href="#">A46012</a>	<ul style="list-style-type: none"> <li>1 mL of master mix</li> <li>1.25 mL of Yellow Sample Buffer</li> </ul>	<p>Store master mix and Yellow Sample Buffer at –25°C to –15°C.</p> <p>Yellow Sample Buffer can be stored at 2–8°C for up to 1 year.</p>
<a href="#">A46109</a>	<ul style="list-style-type: none"> <li>5 mL of master mix</li> <li>1.25 mL of Yellow Sample Buffer</li> </ul>	
<a href="#">A46110</a> <sup>[2]</sup>	<ul style="list-style-type: none"> <li>2 x 5 mL of master mix</li> <li>2 x 1.25 mL of Yellow Sample Buffer</li> </ul>	
<a href="#">A46111</a> <sup>[2]</sup>	<ul style="list-style-type: none"> <li>5 x 5 mL of master mix</li> <li>5 x 1.25 mL of Yellow Sample Buffer</li> </ul>	
<a href="#">A46112</a> <sup>[2]</sup>	<ul style="list-style-type: none"> <li>10 x 5 mL of master mix</li> <li>10 x 1.25 mL of Yellow Sample Buffer</li> </ul>	
<a href="#">A46113</a>	<ul style="list-style-type: none"> <li>50 mL of master mix</li> <li>4 x 1.25 mL of Yellow Sample Buffer</li> </ul>	

<sup>[1]</sup> Catalog numbers that appear as links open the web pages for those products.

<sup>[2]</sup> Products with multiples of the 5-mL size are shipped as multiple kits with the single 5-mL size.

## Required materials

Unless otherwise indicated, all materials are available through [thermofisher.com](https://www.thermofisher.com). "MLS" indicates that the material is available from [fisherscientific.com](https://www.fisherscientific.com) or another major laboratory supplier.




Catalog numbers that appear as links open the web pages for those products.

**IMPORTANT!** Do not use plastics made of polyethylene terephthalate co-polyester, glycol modified (PTEG) for storage of the master mix or the reaction mixes. SYBR™ Green dye is not compatible with this type of plastic material. Polypropylene, high density polyethylene (HDPE), and polystyrene are recommended for storage.

Item	Source
<b>One of the following Applied Biosystems™ instruments:</b>	
<ul style="list-style-type: none"> <li>QuantStudio™ 6 Pro Real-Time PCR System or QuantStudio™ 7 Pro Real-Time PCR System</li> <li>QuantStudio™ 3 or 5 Real-Time PCR System</li> <li>QuantStudio™ 6 / QuantStudio™ 7 Flex Real-Time PCR System</li> <li>QuantStudio™ 12K Flex Real-Time PCR System</li> <li>ViiA™ 7 Real-Time PCR System</li> <li>StepOnePlus™ Real-Time PCR System</li> <li>StepOne™ Real-Time PCR System</li> <li>7500 Fast Real-Time PCR System</li> <li>7500 Real-Time PCR Instrument</li> </ul> <p>Or use a compatible real-time PCR instrument from another supplier.</p>	Contact your local sales office.
<b>Equipment</b>	
Centrifuge with adapter for 96- or 384-well plates	MLS
Laboratory mixer (Vortex or equivalent)	MLS
Microcentrifuge	MLS
Pipettors	MLS
<b>Plastics and other consumables</b>	
Plates and seals for your instrument	<a href="https://www.thermofisher.com/plastics">thermofisher.com/plastics</a>
Disposable gloves	MLS
Pipette tips with filters	MLS
Polypropylene tubes	MLS

Item	Source
<b>Reagents and kits</b>	
<b>One of the following reverse transcription kits, if performing gene expression analysis:</b>	
SuperScript™ IV VILO™ Master Mix	<a href="#">11756050</a>
SuperScript™ IV VILO™ Master Mix with ezDNase™ Enzyme	<a href="#">11766050</a>
High-Capacity cDNA Reverse Transcription Kit	4368814
High-Capacity cDNA Reverse Transcription Kit with RNase Inhibitor	4374967
High-Capacity RNA-to-cDNA™ Kit	4387406
<b>Other reagents</b>	
Nuclease-Free Water (not DEPC-Treated)	<a href="#">AM9930</a>
TE, pH 8.0, RNase-free	<a href="#">AM9858</a>

# Workflow

PowerTrack™ SYBR™ Green Master Mix	
<p><b>Start with cDNA or gDNA</b></p> <p>See “Requirements for input DNA” on page 10.</p>	
<p><b>Set up the plate document or plate file</b></p>	
<p><b>(Optional) Add the Yellow Sample Buffer to the DNA</b></p> <p>See “Prepare the PCR reactions” on page 11.</p>	
<p><b>Combine the components for PCR.</b></p> <p>See “Prepare the PCR reactions” on page 11.</p> <p>The reaction will turn green due to the Yellow Sample Buffer that was added to the DNA and the inert blue dye in the master mix.</p>	
<p><b>Set up and run the real-time PCR instrument</b></p>	
<p><b>Analyze the results</b></p>	

■ Guidelines .....	10
■ Before you begin .....	10
■ Prepare the PCR reactions .....	11
■ Set up and run the real-time PCR instrument .....	14
■ Analyze the results .....	15

## Guidelines

### Requirements for input DNA

Use 1–10 ng of cDNA or 10–100 ng of gDNA per reaction.

For more information, see “RNA guidelines” on page 26 and “Template storage” on page 27.

### Guidelines for PCR reactions

- Four replicates of each reaction are recommended.
- Reaction mixes can be prepared depending upon experimental requirements. Scale the components according to the number of reactions and include 10% overage.
- If using smaller reaction volumes, scale all components proportionally. Reaction volumes <10  $\mu\text{L}$  are not recommended.
- The recommended final primer concentration for primers with a  $T_m$  of 55°C is 400 nM.

### Guidelines for no-template control reactions

No-template control (NTC) reactions can be used to identify PCR contamination. NTC reactions contain all of the reaction components except for the sample.

## Before you begin

### Set up the plate document or plate file

Configure the plate document or plate file.

See the appropriate instrument user guide for detailed instructions.

## Prepare the reagents

- Thaw the master mix.
- Once the master mix is thawed, swirl it to mix thoroughly.
- Thaw the DNA samples and primers on ice, vortex to mix, then centrifuge briefly.
- Vortex the Yellow Sample Buffer prior to use.

## Prepare the PCR reactions

---

**Note:** The Yellow Sample Buffer is optional for the real-time PCR.

---

The Yellow Sample Buffer is supplied at a 40X concentration. It is added to the DNA template. The concentration of Yellow Sample Buffer in the final PCR must be 1X. It is recommended that the DNA template is 10–20% of the volume of the final PCR.

1. *(Optional)* Add the Yellow Sample Buffer (40X) to the amount of DNA that is used in the PCR.

Final reaction volume	Amount of Yellow Sample Buffer
20 $\mu$ L	0.5 $\mu$ L
10 $\mu$ L	0.25 $\mu$ L

The Yellow Sample Buffer is diluted to 1X in the final reaction. See the following tables.

2. *(Optional)* Vortex, then centrifuge the DNA and Yellow Sample Buffer.
3. Combine the master mix, the primers, and nuclease-free water according to the following tables.

- Combine the master mix, the primers, and nuclease-free water with the DNA and Yellow Sample Buffer according to the following tables.

**Note:** If the Yellow Sample Buffer is not used, add nuclease-free water to achieve the total PCR volume.

**Table 2** 20- $\mu$ L reaction

Component	Stock concentration	Final concentration	Volume for 1 reaction (20- $\mu$ L reaction)	Volume for 4 reactions with 10% overage (20- $\mu$ L reaction) <sup>[1]</sup>
<b>Yellow Sample Buffer and DNA (step 1)</b>				
DNA <sup>[2]</sup>	5 ng/ $\mu$ L	0.5 ng/ $\mu$ L	2 $\mu$ L <sup>[3]</sup>	8.8 $\mu$ L
Yellow Sample Buffer	40X	1X	0.5 $\mu$ L	2.2 $\mu$ L
<b>Master mix, primers, and nuclease-free water (step 3)</b>				
PowerTrack™ SYBR™ Green Master Mix	2X	1X	10 $\mu$ L	44.0 $\mu$ L
Forward and reverse primers <sup>[4]</sup>	8,000 nM	400 nM	1 $\mu$ L	4.4 $\mu$ L
Nuclease-free water	—	—	6.5 $\mu$ L	28.6 $\mu$ L
<b>Total PCR volume</b>	—	—	<b>20 <math>\mu</math>L</b>	<b>88 <math>\mu</math>L</b>

<sup>[1]</sup> 10% overage is recommended for pipetting variations.

<sup>[2]</sup> Use 1–10ng of cDNA.

<sup>[3]</sup> Does not exceed 8.5  $\mu$ L.

<sup>[4]</sup> The final primer concentration can vary from 300–800 nM. A final concentration of 400 nM is recommended for primers with a T<sub>m</sub> of 55°C.

Table 3 10- $\mu$ L reaction

Component	Stock concentration	Final concentration	Volume for 1 reaction (10- $\mu$ L reaction)	Volume for 4 reactions with 10% overage (10- $\mu$ L reaction) <sup>[1]</sup>
<b>Yellow Sample Buffer and DNA (step 1)</b>				
DNA <sup>[2]</sup>	5 ng/ $\mu$ L	0.5 ng/ $\mu$ L	1 $\mu$ L <sup>[3]</sup>	4.4 $\mu$ L
Yellow Sample Buffer	40X	1X	0.25 $\mu$ L	1.1 $\mu$ L
<b>Master mix, primers, and nuclease-free water (step 3)</b>				
PowerTrack™ SYBR™ Green Master Mix	2X	1X	5 $\mu$ L	22.0 $\mu$ L
Forward and reverse primers <sup>[4]</sup>	8,000 nM	400 nM	0.5 $\mu$ L	2.2 $\mu$ L
Nuclease-free water	—	—	3.25 $\mu$ L	14.3 $\mu$ L
<b>Total PCR volume</b>	—	—	<b>10 <math>\mu</math>L</b>	<b>44 <math>\mu</math>L</b>

<sup>[1]</sup> 10% overage is recommended for pipetting variations.

<sup>[2]</sup> Use 1–10ng of cDNA.

<sup>[3]</sup> Does not exceed 4.25  $\mu$ L.

<sup>[4]</sup> The final primer concentration can vary from 300–800 nM. A final concentration of 400 nM is recommended for primers with a T<sub>m</sub> of 55°C.

---

**IMPORTANT!** The reaction turns green due to the Yellow Sample Buffer added to the DNA and the inert blue dye in the master mix.

---

- Mix the components thoroughly, then centrifuge briefly to collect the contents at the bottom of the tube.
- Transfer the appropriate volume of each reaction to each well of an optical plate.
- Seal the plate with an optical adhesive cover, then centrifuge briefly to collect the contents at the bottom of each well and eliminate any air bubbles.

PCR can be performed on the reaction plate up to 8 hours after completing the set-up, when stored at room temperature protected from light.

## Set up and run the real-time PCR instrument

1. Set up the thermal protocol according to one of the following tables.

**Note:** Standard cycling conditions are recommended for genomic DNA templates or long amplicons.

**Table 4 Fast cycling mode**

Step	Temperature	Duration	Cycles
Enzyme activation	95°C	2 minutes	1
Denature	95°C	5 seconds	40
Anneal/extend	60°C	30 seconds	

**Table 5 Standard cycling mode**

Step	Temperature	Duration	Cycles
Enzyme activation	95°C	2 minutes	1
Denature	95°C	15 seconds	40
Anneal/extend	60°C	60 seconds	

2. Set the instrument to perform a default dissociation step, according to one of the following tables.

**Table 6 Fast cycling mode**

Step	Ramp rate <sup>[1]</sup>	Temperature	Time
1	1.99°C/second	95°C	15 seconds
2	1.77°C/second	60°C	1 minute
3 (Dissociation)	0.075°C/second	95°C	15 seconds

<sup>[1]</sup> Use the default ramp rate for the StepOnePlus™ Instrument.

**Table 7 Standard cycling mode**

Step	Ramp rate <sup>[1]</sup>	Temperature	Time
1	1.6°C/second	95°C	15 seconds
2	1.6°C/second	60°C	1 minute
3 (Dissociation)	0.075°C/second	95°C	15 seconds

<sup>[1]</sup> Use the default ramp rate for the StepOnePlus™ Instrument.

**Note:** A dissociation step must be performed immediately after the real-time PCR run with PowerTrack™ SYBR™ Green Master Mix.

3. Set up the options.
  - Experiment type: Standard curve
  - Reagent: SYBR™ Green reagents
  - Reporter: SYBR™ Green
  - Quencher: None
  - Passive reference dye: ROX™ dye
  - Ramp speed: Standard or fast
  - Melt curve ramp increment (all instruments, except StepOnePlus™ instrument): Continuous (*StepOnePlus™ only*): Step and hold
4. Set the reaction volume appropriate for the reaction plate.
5. Load the reaction plate into the real-time PCR instrument.
6. Start the run.

## Analyze the results

1. View the amplification plots.  
For more information, see “A typical amplification plot” on page 24.
2. Determine the baseline and threshold cycles ( $C_q$ ) for the amplification curves using the instrument software.  
For more information, see the following sections:
  - “Baseline and threshold values” on page 24
  - “Troubleshooting baseline settings” on page 18
  - “Troubleshooting threshold settings” on page 19
3. Check for nonspecific amplification using melt curves.  
It is important to check for nonspecific amplification because SYBR™ Green dye detects any double-stranded DNA.  
For more information, see “Melt curves” on page 25.
4. Perform relative or absolute quantitation.

Option	Description
Relative quantitation	The target is compared to an internal standard, using either the standard curve or comparative $C_q$ method.
Absolute quantitation	The $C_q$ of the unknown samples is compared against a standard curve with known copy numbers.

# Troubleshooting

■ General troubleshooting .....	16
■ Troubleshooting baseline settings .....	18
■ Troubleshooting threshold settings .....	19



## General troubleshooting


Observation	Possible cause	Recommended action
The $C_q$ values are high, there is poor precision, or the PCR reactions failed	There is insufficient DNA template.	Use up to 100 ng of DNA template per reaction.  Typically, 1–10 ng cDNA or 10–100 ng genomic DNA per reaction is sufficient.
	The quality of the DNA template is poor.	Quantify the amount of DNA template and ensure the recommended amount is used (see “Template quantitation using O.D. 260” on page 27).
		Test the DNA template for the presence of PCR inhibitors. Repeat the PCR reaction with a DNA template free of PCR inhibitors, if necessary.
	The sample has degraded.	Prepare fresh cDNA or gDNA, then repeat the experiment.
	Incorrect volumes of components were pipetted for the PCR reactions.	Prepare the PCR reactions as described in “Prepare the PCR reactions” on page 11.
	Too few PCR cycles were used.	Increase the number of PCR cycles to the default setting of 40 (see “Set up and run the real-time PCR instrument” on page 14).
	There was primer-dimer formation and residual polymerase activity.	<ul style="list-style-type: none"> <li>• Optimize the thermal cycling temperatures.</li> <li>• Reduce the primer concentration.</li> <li>• Redesign the primers.</li> </ul>
Low $\Delta R_n$ or $R_n$ values are obtained	The extension time was too short.	Use the recommended standard cycling thermal profile settings (see “Set up and run the real-time PCR instrument” on page 14).

Observation	Possible cause	Recommended action
Low $\Delta R_n$ or $R_n$ values are obtained (continued)	There was primer-dimer formation and residual polymerase activity.	Optimize the thermal cycling temperatures.
		Reduce the primer concentration.
		Redesign the primers.
The $R_n$ vs. cycle plot is not displayed	ROX™ dye was not selected as the passive reference when the plate document was set up.	Select ROX™ dye as the passive reference when setting up the plate document.
		Select ROX™ dye as the passive reference, then reanalyze the data. The run does not need to be repeated.
The $\Delta R_n$ or $R_n$ values are extremely high	ROX™ dye was not selected as the passive reference when the plate document was set up.	Select ROX™ dye as the passive reference when setting up the plate document.
		Select ROX™ dye as the passive reference, then reanalyze the data. The run does not need to be repeated.
	There was evaporation from the reaction plate.	Ensure that the reaction plate is sealed completely, especially around the edges.
The $R_n$ values obtained in early cycles are low	The $C_q$ value is less than 15.	Adjust the upper baseline range to a value less than 15.
There is high variability across the reaction plate	ROX™ dye was not selected as the passive reference when the plate document was set up.	Select ROX™ dye as the passive reference when setting up the plate document.
	There was evaporation from the reaction plate.	Ensure that the reaction plate is sealed completely, especially around the edges.
There is high variability between replicates	The reaction plate was not mixed well.	Mix the reaction mix gently by inversion, then centrifuge briefly before aliquoting to the reaction plate.
The <b>BADROX</b> flag is displayed	The were droplets on the sides of the wells of the reaction plate.	Set up and repeat the run. Centrifuge the plate to collect the contents at the bottom of the well.
	There was evaporation of the reagents.	Set up and repeat the run. Ensure that the plate is sealed properly.
	The incorrect concentration of passive reference dye was used.	Ensure that the correct concentration of the master mix is added to the reaction. See “Prepare the PCR reactions” on page 11.
The amplification plot is sub-optimal for primers with a $T_m$ of 55°C	The run protocol is not optimal.	Change the annealing temperature to 55°C.
		Design the primers with a $T_m$ of 60°C.




## Troubleshooting baseline settings

Observation	Possible cause	Recommended action
<p>The baseline is set too low</p>  <p><b>Figure 1</b></p>	<p>The baseline is set too low (cycles 3–5).</p>	<p>Manually adjust the baseline to a higher range of cycles (see Figure 3 on page 19).</p>
<p>The baseline is set too high</p>  <p><b>Figure 2</b></p>	<p>The baseline is set too high (cycles 5–20).</p>	<p>Manually adjust the baseline to a lower range of cycles (see Figure 3 on page 19).</p>

Observation	Possible cause	Recommended action
<p>The baseline is set correctly</p>  <p>Figure 3</p>	<p>The baseline is set correctly (cycles 3–15).</p>	<p>No action is required.</p>

## Troubleshooting threshold settings

Observation	Possible cause	Recommended action
<p>The threshold is set too low</p>  <p>Figure 4</p>	<p>The threshold is set too low.</p>	<p>Manually adjust the threshold to a higher <math>\Delta R_n</math> (see Figure 6 on page 20).</p>

Observation	Possible cause	Recommended action
<p>The threshold is set too high</p>  <p>Figure 5</p>	<p>The threshold is set too high.</p>	<p>Manually adjust the threshold to a lower <math>\Delta R_n</math> (see Figure 6 on page 20).</p>
<p>The threshold is set correctly</p>  <p>Figure 6</p>	<p>The threshold is set correctly.</p>	<p>No action is required.</p>

# Background information

- Master mix components ..... 21
- Overview of the chemistry ..... 22
- Two-step RT-PCR ..... 23
- A typical amplification plot ..... 24
- Melt curves ..... 25

## Master mix components

Table 8 Function of the components of the master mix

Component	Function
Antibody-mediated hot start polymerase	<ul style="list-style-type: none"> <li>• Provides tight control over <i>Taq</i> activation, preventing undesirable early activity of the polymerase at low temperatures that can lead to nonspecific amplification.</li> <li>• Allows flexibility in reaction set-up, including pre-mixing of PCR reagents and storage at room temperature for up to 8 hours prior to cycling.</li> <li>• Allows activation of polymerase after only 2 minutes at 95°C.</li> </ul>
Heat-labile uracil-DNA glycosylase (UDG)	<ul style="list-style-type: none"> <li>• A 26 kDa recombinant enzyme derived from the thermolabile UDG gene isolated from marine bacteria, and expressed in <i>E. coli</i>.</li> <li>• Prevents reamplification of carryover PCR products by removing any uracil incorporated into single- or double-stranded amplicons.</li> <li>• Acts on single- and double-stranded dU-containing DNA by hydrolyzing uracil-glycosidic bonds at dU-containing DNA site, creating an alkali-sensitive apyrimidic site in the DNA.</li> <li>• Prevents reamplification of carryover PCR products in an assay if all previous PCR for the assay was performed using a dUTP-containing master mix.</li> <li>• Allows stability of PCR products for 72 hours post-amplification.</li> <li>• Has no activity on RNA or dT-containing DNA.</li> </ul>
dUTP/dTTP	<ul style="list-style-type: none"> <li>• Enables UDG activity and maintains optimal PCR results.</li> </ul>
SYBR™ Green dye	<ul style="list-style-type: none"> <li>• Detects PCR products by fluorescing upon binding to double-stranded DNA formed during PCR (see “Overview of the chemistry” on page 22).</li> </ul>
ROX™ passive reference	<ul style="list-style-type: none"> <li>• Provides an internal reference to which the reporter-dye signal can be normalized during data analysis.</li> <li>• Normalization is necessary to correct for fluorescence fluctuations due to changes in volume.</li> </ul>

## Overview of the chemistry

The SYBR™ Green dye is used to detect PCR products by binding to double-stranded DNA formed during PCR.

1. When the master mix is added to a sample, the SYBR™ Green dye immediately binds to all double-stranded DNA (dsDNA) present in the sample.  
The SYBR™ Green dye is only fluorescent when bound to dsDNA.
2. During PCR, DNA polymerase amplifies the target sequence which creates the PCR products.
3. The SYBR™ Green dye then binds to each new copy of double-stranded DNA, generating a fluorescent signal.
4. As the PCR progresses, more PCR product is created.  
The SYBR™ Green dye binds to all double-stranded DNA, so the result is an increase in fluorescence intensity proportional to the amount of PCR product produced.


Figure 7 DNA template

Figure 8 SYBR™ Green dye binds to all double-stranded DNA


Figure 9 Denaturation


Figure 10 Polymerization



Figure 11 Completion

## Two-step RT-PCR

For more information, go to [thermofisher.com/qpcr/education](https://www.thermofisher.com/qpcr/education).


When performing a two-step RT-PCR reaction, total RNA or mRNA must first be reverse transcribed into cDNA.

1. In the reverse transcription (RT) step, cDNA is reverse transcribed from total RNA samples using random primers from the reverse transcription kit.
2. In the PCR step, PCR products are synthesized from cDNA samples using the master mix.

## A typical amplification plot


A typical amplification plot is shown below.

- ① Plateau phase
- ② Linear phase
- ③ Exponential phase (geometric phase)
- ④ Background
- ⑤ Baseline

## Baseline and threshold values

- **Baseline**—the initial cycles of PCR in which there is little change in fluorescence signal.
- **$C_q$** —the intersection of the threshold with the amplification plot.
- **Threshold**—set above the background and within the exponential growth phase of the amplification curve.


Automatic calculation of the baseline and threshold can be conducted using the software on your instrument. Alternatively, baseline and threshold can be set manually.

For examples of amplification plots where the baseline values and the threshold values are set too high or too low, see “Troubleshooting baseline settings” on page 18 and “Troubleshooting threshold settings” on page 19.

## Melt curves

A melt curve is a graph that displays dissociation data from the amplicons of quantitative PCR runs. Change in fluorescence, due to a dye interacting with double-stranded DNA, is plotted against temperature. A single peak indicates specific amplification, whereas multiple peaks or shoulders indicate nonspecific amplification or primer-dimer formation.

Primer-dimers are most prevalent in NTC wells and sample wells containing a low concentration of template.

This melt curve shows typical primer-dimer formation. The specific product is shown with a melting temperature ( $T_m$ ) of 80.5°C, but the primer-dimer has a characteristically lower  $T_m$  of 75°C.

- ① Melt curve of a specific product
- ② Melt curve of a primer-dimer

# Template quality and quantity

■ DNA template quality .....	26
■ RNA guidelines .....	26
■ Template quantitation using O.D. 260 .....	27
■ Template storage .....	27

Go to [thermofisher.com/qpcducation](https://www.thermofisher.com/qpcducation) for more information about DNA template quality, RNA guidelines, template quantitation, and template storage.

## DNA template quality

Both agarose gel electrophoresis and spectrophotometry are used to examine DNA quality.

- **Agarose gel electrophoresis** – Purified DNA should run as a single band on an agarose gel. Agarose gels reveal contaminating DNAs and RNAs, but not proteins.
- **Spectrophotometry** – The  $A_{260}/A_{280}$  ratio should be 1.8 to 2.0. Smaller ratios usually indicate contamination by protein or organic chemicals. Spectrophotometry can reveal protein contamination, but not DNA or RNA contamination.

## RNA guidelines

RNA should be reverse transcribed into cDNA prior to use in a PowerTrack™ SYBR™ Green Master Mix reaction. For recommended reverse transcription kits, see “Required materials” on page 7.

For optimal performance prior to reverse transcription, total RNA or mRNA should be:

- Between 0.002 µg/µL and 0.2 µg/µL
- Less than 0.005% of genomic DNA by weight
- Free of inhibitors of reverse transcription and PCR
- Dissolved in PCR-compatible buffer
- Free of RNase activity

---

**IMPORTANT!** If you suspect that the RNA contains RNase activity, add RNase inhibitor to the reverse transcription reaction at a final concentration of 1.0 U/µL.

---

- Nondenatured

---

**IMPORTANT!** It is not necessary to denature the RNA. Denaturation of the RNA may reduce the yield of cDNA for some gene targets.

---

## Template quantitation using O.D. 260

Template quantitation is critical for successful PCR reactions. The most common way to determine DNA quantity is to measure the absorbance (optical density or O.D.) of a sample at 260 nm in a spectrophotometer.

One O.D. unit is the amount of substance dissolved in 1.0 mL that gives an absorbance reading of 1.00 in a spectrophotometer with a 1-cm path length. The wavelength is assumed to be 260 nm unless state otherwise.  $A_{260}$  values can be converted into  $\mu\text{g}/\mu\text{L}$  using Beer's Law:

Absorbance (260 nm) = sum of extinction coefficient contributions  $\times$  cuvette pathlength  $\times$  concentration

The following formulas are derived from Beer's Law:

- Concentration of single-stranded DNA =  $A_{260} \times 33 \mu\text{g}/\mu\text{L}$
- Concentration of double-stranded DNA =  $A_{260} \times 50 \mu\text{g}/\mu\text{L}$
- Concentration of single-stranded RNA =  $A_{260} \times 40 \mu\text{g}/\mu\text{L}$

---

**Note:** Absorbance measurements of highly concentrated (O.D.  $>1.0$ ) or very dilute (O.D.  $<0.05$ ) DNA or RNA samples can be inaccurate. Dilute or concentrate the DNA/RNA to obtain a reading within the acceptable range.

---

## Template storage

- Store purified RNA templates at  $-20^{\circ}\text{C}$  or  $-70^{\circ}\text{C}$  in Nuclease-Free Water.
- Store purified DNA templates at  $-20^{\circ}\text{C}$  or  $-70^{\circ}\text{C}$  in TE, pH 8.0.

# Primer design, target sequences, and optimizing primer concentration

- Primer design guidelines ..... 28
- Identifying target sequence and amplicon size ..... 29
- Selecting an amplicon site for cDNA ..... 29
- Optimize primer concentrations for PCR ..... 30

## Primer design guidelines

Primers should be designed using Primer Express™ Software or similar software. See the *Primer Express™ Software Version 3.0 Getting Started Guide* (Pub. No. 4362460).

- Keep the GC content in the 30–70% range.
- The optimal primer length is 20 bases.
- Avoid runs of identical nucleotides. If repeats are present, there must be fewer than four consecutive G residues.
- Make sure the last five nucleotides at the 3' end contain no more than two G and/or C bases.

Template	Design guideline
DNA	Design the primers as described above.
Plasmid DNA	
Genomic DNA	
cDNA	Design the primers as described above and see “Selecting an amplicon site for cDNA” on page 29.
RNA	Design the primers as described above.

## Avoiding primer-dimers

Use primers that contain dA nucleotides near the 3' ends so that any primer-dimer generated is efficiently degraded by UDG at least as well as any dU-containing PCR products. The farther a dA nucleotide from the 3' end, the more likely partially degraded primer-dimer molecules can serve as a template for a subsequent PCR amplification.

Production of primer-dimers could lower the amplification yield of the desired target region. If primers cannot be selected with dA nucleotides near the ends, consider using primers with 3' terminal dU-nucleotides. Single-stranded DNA with terminal dU nucleotides are not substrates for UDG, and therefore the primers are not degraded. Biotin-dUMP derivatives are not substrates for UDG.

For more information about designing primers, see “Primer design guidelines” on page 28.

Do not use UDG in subsequent amplifications of dU-containing PCR template, such as in nested PCR protocols. The UNG degrades the dU-containing PCR products, preventing further amplification.

## Identifying target sequence and amplicon size

A target template is a DNA sequence, including cDNA, genomic DNA, or plasmid nucleotide sequence that you want to amplify.

Primers are designed to amplify amplicons (segments of DNA) within the target sequence using Primer Express™ Software. Shorter amplicons work best. Consistent results are obtained for amplicon size ranges from 50–150 bp.

## Selecting an amplicon site for cDNA

Selecting a good amplicon site ensures amplification of the target cDNA without co-amplifying the genomic sequence, pseudogenes, and related genes.

- The amplicon should span one or more introns to avoid amplification of the target gene in genomic DNA.
- The primer pair must be specific to the target gene; the primer pair does not amplify pseudogenes or other related genes.
- Primers should be designed according to the guidelines in the Primer Express™ Software.
- Amplicons should be tested and those with the highest signal-to-noise ratio should be selected (low  $C_q$  with cDNA and no amplification with no template control or genomic DNA).
- The sequence may need to be examined and the amplicon redesigned if no good sequence is found. Alternatively, more sites may need to be screened.

If the gene of interest does not have introns, then an amplicon cannot be designed that amplifies the mRNA sequence without amplifying the genomic sequence. RT–minus controls may need to be run.

## Optimize primer concentrations for PCR

### Overview

By independently varying the forward and reverse primer concentrations, you can identify the primer concentrations that provide optimal assay performance. The primer concentrations you select should provide a low  $C_q$  and a high  $\Delta R_n$  when run against the target template, but should not produce nonspecific product formation with NTCs.

### Quantitate primers

1. Measure the absorbance (at 260 nm of a 1:100 dilution) of each primer oligonucleotide in TE buffer.
2. Calculate the sum of extinction coefficient contributions for each primer:
  - Extinction coefficient contribution =  $\Sigma$  (extinction coefficient  $\times$  number of bases in oligonucleotide sequence)
3. Calculate the oligonucleotide concentration in  $\mu\text{M}$  for each primer:
  - Absorbance at 260 nm = sum of extinction coefficient contribution  $\times$  cuvette pathlength  $\times$  concentration / 100
  - Rearrange to solve for concentration:
    - Concentration =  $100 [\text{absorbance at 260 nm} / (\text{sum of extinction coefficient contribution} \times \text{cuvette pathlength})]$

#### Example calculation of primer concentration

In this example, the concentration of a primer (in TE buffer, diluted 1:100), with the sequence 5' - CGTACTCGTTCGTGCTGC - 3' is calculated using the following values:

Chromophore	Extinction coefficient	Number of specific chromophores in example sequence	Extinction coefficient contribution
A	15,200	1	15,200
C	7050	6	42,300
G	12,010	5	60,050
T	8400	6	50,400
Total	—	—	167,950

- Measured absorbance at 260 nm = 0.13
- Sum of extinction coefficient =  $167,950 \text{ M}^{-1}\text{cm}^{-1}$  contributions for probe
- Cuvette pathlength = 0.3 cm
- Absorbance (260 nm) = sum of extinction coefficient contributions  $\times$  cuvette pathlength  $\times$  oligonucleotide concentration / 100

- $0.31 = 167,950 \text{ M}^{-1}\text{cm}^{-1} \times 0.3 \text{ cm} \times C / 100$
- $C = 258 \text{ }\mu\text{M}$

## Determine the optimal primer concentration for PCR

Calibrate your instrument for SYBR™ Green dye, if necessary. See the instrument user guide for calibration instructions. It is recommended to calibrate your instrument every six months.

1. Prepare a 96-well reaction plate.

Use 10–100 ng of gDNA or 1–10 ng of cDNA template. The final concentration of the master mix is 1X.

**Note:** The plate configuration accounts for four replicates of each of the following nine variations of primer concentration applied to both template and NTC wells:

Reverse primer (nM)	Forward primer (nM)		
	300	500	800
300	300 / 300	500 / 300	800 / 300
500	300 / 500	500 / 500	800 / 500
800	300 / 800	500 / 800	800 / 800

2. Set up the thermal protocol (see “Set up and run the real-time PCR instrument” on page 14).
3. Load the plate into the real-time PCR instrument.
4. Start the run.
5. Compile the results for  $\Delta R_n$  and  $C_q$ , then select the minimum forward and reverse primer concentrations that yield the maximum  $\Delta R_n$  values and low  $C_q$  values.

## Confirm the absence of nonspecific amplification

Melt curves help you select the optimal primer concentrations for your quantification assays with SYBR™ Green dye.

1. Review the linear view of the amplification plot in your NTC wells.

**Note:** In Figure 12 on page 32, the strong amplification of the NTC wells indicates that significant nonspecific amplification is occurring.

2. Generate a melt curve with your real-time PCR system.

**Note:** In the example shown in Figure 13 on page 32, the melting temperature of the product generated in the absence of template is lower than the melting temperature of the specific product generated with template. This variation is typical of primer-dimer formation, and it indicates that lower primer concentration may provide optimal results.


Figure 12 Amplification plot (linear view) demonstrating suspected nonspecific amplification in NTC wells


Figure 13 Melt curve analysis confirming that product in NTC wells has a melting temperature different from the specific product

- ① Target amplification
- ② NTC (nonspecific amplification)

# Safety

■ Chemical safety .....	34
■ Biological hazard safety .....	35

**WARNING! GENERAL SAFETY.** Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, and so on). To obtain SDSs, visit [thermofisher.com/support](https://www.thermofisher.com/support).

## Chemical safety

**WARNING! GENERAL CHEMICAL HANDLING.** To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below. Consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the "Documentation and Support" section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with sufficient ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if needed) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

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**WARNING! HAZARDOUS WASTE (from instruments).** Waste produced by the instrument is potentially hazardous. Follow the guidelines noted in the preceding General Chemical Handling warning.

**WARNING! 4L Reagent and Waste Bottle Safety.** Four-liter reagent and waste bottles can crack and leak. Each 4-liter bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position.

## Biological hazard safety

**WARNING! BIOHAZARD.** Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Conduct all work in properly equipped facilities with the appropriate safety equipment (for example, physical containment devices). Safety equipment can also include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 6th Edition, HHS Publication No. (CDC) 300859, Revised June 2020  
<https://www.cdc.gov/labs/pdf/CDC-BiosafetymicrobiologicalBiomedicalLaboratories-2020-P.pdf>
- Laboratory biosafety manual, fourth edition. Geneva: World Health Organization; 2020 (Laboratory biosafety manual, fourth edition and associated monographs)  
[www.who.int/publications/i/item/9789240011311](http://www.who.int/publications/i/item/9789240011311)

# Documentation and support

## Related documentation

Document	Pub. No.
<i>PowerTrack™ SYBR™ Green Master Mix Quick Reference</i>	MAN0018826

## Customer and technical support

Visit [thermofisher.com/support](http://thermofisher.com/support) for the latest service and support information.

- Worldwide contact telephone numbers
- Product support information
  - Product FAQs
  - Software, patches, and updates
  - Training for many applications and instruments
- Order and web support
- Product documentation
  - User guides, manuals, and protocols
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

---

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

---

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

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**UNI EN ISO 9001-2015 (ISO 9001-2015)**

per i seguenti Processi  
*concerning the following kinds of Processes*

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*Management of the manufacturing and placing on the market of sterile tampons for sampling of biological specimens in natural orifice and in surgical field. Design and manufacturing of diagnostic medical devices for laboratories of analysis. Management of the manufacturing and placing on the market of invasive medical devices with respect to body orifices (class I sterile). Manufacturing of invasive medical devices with respect to body orifices (class I sterile). Marketing of medical and diagnostic devices in vitro. Marketing of laboratory articles.*

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*Management of the manufacturing and placing on the market of sterile tampons for sampling of biological specimens in natural orifice and in surgical field.  
Design and manufacturing of diagnostic medical devices for laboratories of analysis. Management of the manufacturing and placing on the market of  
invasive medical devices with respect to body orifices (class I sterile). Manufacturing of invasive medical devices with respect to body orifices (class I sterile).  
Marketing of medical and diagnostic devices in vitro.*

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Dispositivo Latex free

### 5ML CYLINDRICAL TEST TUBES

5ML CYLINDRICAL TEST TUBES WITH "U" BOTTOM SORVALL CW1 TYPE, NOT GRADUATED, WITHOUT RIM.

MANUFACTURED IN MEDICAL POLYPROPYLENE (PP), TRANSLUCENT AND UNBREAKABLE. AUTOCLAVABLE UP TO +121 °C

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CE Marked product - manufactured in compliance with 98/79/CE Directive and D.lgs 332 dtd 08/09/2000

CARATTERISTICHE PRINCIPALI		TECHNICAL FEATURES
Stato microbiologico	<b>NON STERILE / NOT STERILE</b>	<i>Microbiological status</i>
Materiale impiegato	POLIPROPILENE / POLYPROPYLENE	<i>Raw material</i>
Temperature tollerate provetta	MIN -10°C MAX +121°C	<i>Temperature range - test tube</i>
Dimensioni (mm)	Ø 12 x 75	<i>Dimensions (mm)</i>
Volume (ml)	5,0	<i>Volume (ml)</i>
Spessore (mm)	1,0	<i>Thickness (mm)</i>
Peso (gr.)	2,1	<i>Weight (gr.)</i>
Validità del prodotto	5 ANNI / YEARS	<i>Shelf life</i>

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*Intended purpose is "IN VITRO MEDICAL DEVICE" adapted to contain a human biological sample (for example urine, blood, semen, saliva, sputum, pus, etc) in order to perform diagnostic analysis laboratory. **For professional use only.***

**National classification of medical devices (CND - For Italian law) > W050301020102** (Samples analyses, plastic tubes without additives).

**EDMA code > 51091001** - Other containers for samples of human origin

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*Keep out of flame or heat sources which might damage the product*

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*Do not use after expiry date or if packing is opened*

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*Do not re-use: Disposable device*

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*Do not vary the intended purpose of the product*

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*Store in dry place, Temperature range: min -10°C max +50°C*

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*Disposal: use appropriate personal protective equipment and act according to applicable regulations*

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*Before use with particular substance check the resistance / compatibility chart on our catalogue*

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
*For a maximum centrifuge speed of 5,000 r.p.m to be kept for 20" max*

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meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the “final inspection report” of the device. /

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## Gilson Pipette Tips

### Gilson Pipette Tips

Designed for use in a wide variety of pipetting applications, compatibility with Gilson Pipettors.

PP material.

Size: 10ul, 200ul, 1000ul, 5000ul

Various color optional easy to identify

Non sterile.

<b>Cat No.</b>	<b>Description</b>	<b>Qty/Case(bags)</b>
640101	Gilson pipette tips, 10ul, clear, 1000tips/bag	100
640102	Gilson pipette tips, 200ul, yellow, 1000tips/bag	50
640103	Gilson pipette tips, 1000ul, blue, 500tips/bag	50
640104	Gilson pipette tips, 5000ul, clear, 300tips/bag	30

This is a translation of the certificate ES16/20725.01

# DELTALAB, S.L.

Pol.Ind. La Llana, Plaza de la Verneda, 1, 08191 Rubí, Barcelona

Has been assessed under the management system of the certified organisation defined in the main certificate ES16/20725 as meeting the requirements of

## ISO 9001:2015

For the following activities

Design, manufacture and sale of laboratory material for the collection, transport and conservation of samples for microbiological, molecular biology, hematology, biochemistry, histology, microscopy and colorimetric analysis.

Commercialization of equipment for the storage of prepared samples, cryogenic stored samples, general labware and industrial packages.

Commercialization of equipment and instrumentation for the laboratory, diagnostic kits, healthcare products and cosmetics.

This certificate is valid from 11 October 2022 until 11 October 2025 and remains valid subject to satisfactory surveillance audits.

Issue 2.

The validity of this certificate depends on the validity of the main certificate.

SGS International Certification Services Iberica, S.A.U.

C/Trespaderne, 29. 28042 Madrid. España

t +34 91 313 8115 - [www.sgs.com](http://www.sgs.com)

This document is an authentic electronic certificate for Client' business purposes use only. Printed version of the electronic certificate are permitted and will be considered as a copy. This document is issued by the Company subject to SGS General Conditions of certification services available on Terms and Conditions | SGS. Attention is drawn to the limitation of liability, indemnification and jurisdictional clauses contained therein. This document is copyright protected and any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful.

This is a translation of the certificate ES19/86440.01

## DELTALAB, S.L.

Pol.Ind. La Llana, Plaza de la Verneda, 1, 08191 Rubí, Barcelona

Has been assessed under the management system of the certified organisation defined in the main certificate ES19/86440 as meeting the requirements of

### ISO 14001:2015

For the following activities

Design, manufacture and sale of laboratory material for the collection, transport and conservation of samples for microbiological, molecular biology, hematology, biochemistry, histology, microscopy and colorimetric analysis.

Commercialization of equipment for the storage of prepared samples, cryogenic stored samples, general labware and industrial packages.

Commercialization of equipment and instrumentation for the laboratory, diagnostic kits, healthcare products and cosmetics.

This certificate is valid from 31 August 2022 until 29 August 2025 and remains valid subject to satisfactory surveillance audits.

Issue 2.

The validity of this certificate depends on the validity of the main certificate.

SGS International Certification Services Iberica, S.A.U.  
C/Trespaderne, 29. 28042 Madrid. España  
t +34 91 313 8115 - www.sgs.com

This document is an authentic electronic certificate for Client' business purposes use only. Printed version of the electronic certificate are permitted and will be considered as a copy. This document is issued by the Company subject to SGS General Conditions of certification services available on [Terms and Conditions | SGS](#). Attention is drawn to the limitation of liability, indemnification and jurisdictional clauses contained therein. This document is copyright protected and any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful.

Certificate ES10/81671

SGS

The management system of

# DELTALAB, S.L.

Polígono Industrial La Llana, Plaza de la Verneda 1, 08191 Rubi, Barcelona, Spain

has been assessed and certified as meeting the requirements of

**ISO 13485:2016**

**EN ISO 13485:2016**

For the following activities

Design, manufacture and sale of sterile and non-sterile medical devices for the collection, transport and conservation of biological samples for clinical and IVD analysis.

Distribution of non-active medical devices and in vitro diagnostic medical devices.

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This certificate is valid from 12 October 2022 until 11 October 2025 and remains valid subject to satisfactory surveillance audits.

Issue 10. Certified since 12 October 2010.

*Jonathan M. Hall*

Jonathan Hall  
Global Head - Certification Services

SGS United Kingdom Ltd  
Rossmore Business Park, Ellesmere Port, Cheshire, CH65 3EN, UK  
t +44 (0)151 350-6666 - www.sgs.com

This document is an authentic electronic certificate for Client' business purposes use only. Printed version of the electronic certificate are permitted and will be considered as a copy. This document is issued by the Company subject to SGS General Conditions of certification services available on Terms and Conditions | SGS. Attention is drawn to the limitation of liability, indemnification and jurisdictional clauses contained therein. This document is copyright protected and any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful.

Deltalab, S.L. defines and makes public its commitment with the Standard ISO 9001:2015 Quality Management Systems, ISO 14001:2015 Environmental Management Systems and ISO 13485:2016 Medical devices – Quality Management Systems, with the aim to create value and satisfy all its interested parties:

- Shareholders
- Members of the organisation
- Customers and suppliers
- All members of the surrounding community

The development of this Integrated Management System Policy is carried out with the philosophy of Continuous Improvement and with the support of all the processes described in our Integrated Management System, in order to achieve the following objectives:

1. Become leaders in the design and manufacture of single use products for the laboratory.
2. Bring solutions to cover the current and future customer needs, related to:
  - Design, manufacture and sale of laboratory material for the collection, transport and conservation of samples for microbiology, molecular biology, haematology, biochemistry, histology, microscopy and colorimetric analysis.
  - Design, manufacture and sale of sterile and non-sterile medical devices for the collection, transport and conservation of biological samples for clinical and IVD analysis.
  - Commercialization of diagnosis reagents, equipment and instrumentation for laboratory and equipment for the storage of prepared samples, cryogenic stored samples, general labware and industrial packages.
  - Commercialization of personal care, cosmetics and dietetic products
3. Maintain a constant growth, both in local and international markets, by means of mergers, acquisitions and by launching new products.
4. Achieve the full satisfaction of our customers, by means of a strict compliance to the agreements and expectations agreed with them, as well as the excellence in the service.
5. Reach a high level of innovation of our products and processes, in cooperation with universities, research centers, key opinion leaders and experts, both local and international.
6. Fulfil the legislation and regulatory requirements applicable to the activities carried out by the company, including those applicable to the quality of products and the environmental management.
7. Commit ourselves with the environmental protection, including the prevention of pollution.
8. Achieve and keep a high motivation and involvement of all members of the organisation, suppliers, distributors and customers, by fulfilling the highest Quality and environmental protection standards.
9. Improve the working conditions of all employees and ensure the technical capacity of the personnel by giving them the adequate training with the aim to achieve the required competence.
10. Establish a close relationship with the suppliers and guarantee the maximum quality of materials supplied by means of quality agreements.

The Integrated Management System is periodically reviewed to define the required actions to ensure that:

- ✓ The System is efficient, so that it is a tool for the routine of all the members of the organisation.
- ✓ The customer needs and requirements are duly identified, and their expectations are always met.
- ✓ All members of the organisation are familiar with and know the objectives and policy of the Integrated Management System, and that adequate training plans are defined to achieve them.
- ✓ Encourage the Continuous Improvement Philosophy, both related to Quality and Environmental Management.

This Policy is made available for the public and all interested parties.

JOSEP SAEZ  
Managing Director  
January 2019

# Adhesive sealing film

Code: 900300, Deltalab, Spain

Code: 900300

Sealing film for use with microplates, multiwell plates and microtiter plates.

Advantages:

1. Minimises the risk of contamination or reagent spillage during ELISA or PCR processes.
2. Minimises the risk of contamination from tube to tube and from plate to plate.
3. Prevents sample evaporation.

A 5 mm wide strip (opaque white) at the lateral edges of the film helps pull the film from its protective paper and prevents it sticking onto fingers. The film is thermostable and functional from  $-70\text{ }^{\circ}\text{C}$  to  $95\text{ }^{\circ}\text{C}$  at 75% humidity.

RNAse, DNAse, DNA and PCR inhibitors free. DMSO resistant.

We recommend to use the “roller” 900330 to ensure a perfect seal.

## QPCR adhesive film

Code: 900301, Deltalab, Spain

Code: 900301

Sealing film for use with microplates, multiwell plates and microtiter plates.

Advantages:

1. Minimises the risk of contamination or reagent spillage during ELISA or PCR processes.
2. Minimises the risk of contamination from tube to tube and from plate to plate.
3. Prevents sample evaporation.

A 5 mm wide strip (opaque white) at the lateral edges of the film helps pull the film from its protective paper and prevents it sticking onto fingers. The film is thermostable and functional from  $-70\text{ }^{\circ}\text{C}$  to  $95\text{ }^{\circ}\text{C}$  at 75% humidity.

RNAse, DNAse, DNA and PCR inhibitors free. DMSO resistant.

We recommend to use the “roller” 900330 to ensure a perfect seal.

	<b>AO Vector-Best</b>	01-6-6-a1 (version 6)
	EC Declaration of conformity DoC_IVDD_B_ELISA_0522	Page 1 of 3

## EC DECLARATION OF CONFORMITY

AO Vector-Best hereby ensures under own responsibility and declares that the products indicated on pages 2–3 are in conformity with relevant provisions of EC Council Directive 98/79/EC on *in vitro* diagnostic medical devices and fulfill the essential requirements of that directive.

### Classification of products:

Other devices (all devices except Annex II and self-testing devices)

### Conformity assessment procedure:

Annex III (not including section 6).

### Manufacturer:

AO Vector-Best

Address: 630559, Koltsovo, Novosibirsk Region, Research and Production area, building 36, office 211, Russian Federation, tel. +7 383 332 81 34.

### European authorized representative:

BIORON GmbH

Address: In den Rauhweiden 20, 67354 Roemerberg, Germany, tel.: +49 (0) 6232 298 450.

Date: 2022/05/24  
Place: Novosibirsk

Valid until: 2027/05/26

	<b>AO Vector-Best</b>	01-6-6-a1 (version 6)
	EC Declaration of conformity DoC_IVDD_B_ELISA_0522	Page 2 of 3

No.	Product name	Identification data	REF
1.	Vectohep A-IgG	Enzyme immunoassay kit for the detection and quantification of IgG to hepatitis A virus	D-0362
2.	VectoTBEV-IgG	Enzyme immunoassay kit for the detection and quantification of IgG to tick-borne encephalitis virus	D-1156
3.	VectoMeasles-IgG	Enzyme immunoassay kit for the detection and quantification of IgG to measles virus in blood serum (plasma)	D-1356
4.	VectoMeasles-IgM	Enzyme immunoassay kit for the detection of IgM to measles virus in blood serum (plasma)	D-1358
5.	LymeBest-IgG	Enzyme immunoassay kit for the detection of IgG to <i>Borrelia burgdorferi sensu lato</i> complex	D-1452
6.	LymeBest-IgM	Enzyme immunoassay kit for the detection of IgM to <i>Borrelia burgdorferi sensu lato</i> complex	D-1454
7.	Rotavirus-antigen-EIA-BEST	Enzyme immunoassay kit for the detection of human rotavirus antigen	D-1652
8.	Adenovirus-antigen-EIA-BEST	Enzyme immunoassay kit for the detection of human adenovirus antigen	D-1654
9.	VectoEBV-NA-IgG	Enzyme immunoassay kit for the detection of IgG to Epstein-Barr virus nuclear antigen in blood serum (plasma)	D-2170
10.	VectoEBV-EA-IgG	Enzyme immunoassay kit for the detection of IgG to Epstein-Barr virus early antigens in blood serum (plasma)	D-2172
11.	VectoEBV-VCA-IgM	Enzyme immunoassay kit for the detection of IgM to viral capsid antigen of Epstein-Barr virus in blood serum (plasma)	D-2176
12.	VectoMumps-IgG	Enzyme immunoassay kit for the detection of IgG to mumps virus in blood serum (plasma)	D-2602
13.	Toxocara-IgG-EIA-BEST	Enzyme immunoassay kit for the detection of IgG to <i>Toxocara</i> antigens in blood serum (plasma)	D-2752
14.	Trichinella-IgG-EIA-BEST	Enzyme immunoassay kit for the detection of IgG to <i>Trichinella</i> antigens in blood serum (plasma)	D-3152

	<b>AO Vector-Best</b>	01-6-6-a1 (version 6)
	EC Declaration of conformity DoC_IVDD_B_ELISA_0522	Page 3 of 3

15.	Yersinia-IgG-EIA-BEST	Enzyme immunoassay kit for the detection of IgG to causative agents of yersiniosis	D-3202
16.	Yersinia-IgA-EIA-BEST	Enzyme immunoassay kit for the detection of IgA to causative agents of yersiniosis	D-3204
17.	Yersinia-IgM-EIA-BEST	Enzyme immunoassay kit for the detection of IgM to causative agents of yersiniosis	D-3206
18.	Echinococcus-IgG-EIA-BEST	Enzyme immunoassay kit for the detection of IgG to <i>Echinococcus granulosus</i> antigens in blood serum (plasma)	D-3356
19.	Ascaris-IgG-EIA-BEST	Enzyme immunoassay kit for the detection of IgG to <i>Ascaris lumbricoides</i> antigens in blood serum (plasma)	D-3452
20.	IgA-Transglutaminase-EIA-BEST	Enzyme immunoassay kit for the quantification of IgA to tissue transglutaminase in blood serum (plasma)	D-3758
21.	IgG-Transglutaminase-EIA-BEST	Enzyme immunoassay kit for the quantification of IgG to tissue transglutaminase in blood serum (plasma)	D-3760
22.	Pepsinogen 1-EIA-BEST	Enzyme immunoassay kit for the quantification of pepsinogen 1 in blood serum	D-3762
23.	Pepsinogen 2-EIA-BEST	Enzyme immunoassay kit for the quantification of pepsinogen 2 in blood serum	D-3764
24.	VectoNile-IgM	Enzyme immunoassay kit for the detection of IgM to West Nile virus in blood serum (plasma)	D-5150
25.	VectoNile-IgG	Enzyme immunoassay kit for the detection of IgG to West Nile virus in blood serum (plasma)	D-5152
26.	VectoNile-IgG-avidity	Enzyme immunoassay kit for the determination of the avidity index of IgG to West Nile virus in blood serum (plasma)	D-5154
27.	Interleukin-6-EIA-BEST	Enzyme immunoassay kit for the quantification of interleukin-6 in blood serum and urine	A-8768

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**mdc medical device certification GmbH**

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# Certificate

**mdc medical device certification GmbH**  
certifies that

VECTOR

**AO Vector-Best  
Research and Production Area  
Building 36, Office 211, Koltsovo  
630559 Novosibirsk region  
Russian Federation**

with the locations listed in the attachment  
for the scope

**Design and development, production and distribution of  
medical devices for in vitro diagnostics (PCR, ELISA, Biochemistry)**

has introduced and applies a

## Quality Management System

The mdc audit has proven that this quality management system  
meets all requirements of the following standard

### EN ISO 13485

Medical devices – Quality management systems –  
Requirements for regulatory purposes

EN ISO 13485:2016 + AC:2016 - ISO 13485:2016

Valid from	2020-07-04
Valid until	2023-07-03
Registration no.	D1213100019
Report no.	P20-00568-173687
Stuttgart	2020-06-02

  
Head of Certification Body

**Attachment of the certificate**

**No. D1213100019**

date 2020-06-02

Page 1 of 1

<b>Location</b>	<b>Scope</b>
AO Vector-Best Arbuzova str. 1/1, 630117 Novosibirsk Russian Federation	design and development, production and distribution of medical devices for in vitro diagnostics
AO Vector-Best Research and Production area, building 36, Koltsovo, 630559 Novosibirsk region Russian Federation	design and development, production of medical devices for in vitro diagnostics
AO Vector-Best Pasechnaya str, 3, 630117 Novosibirsk Russian Federation	design and development, production of medical devices for in vitro diagnostics

  
Head of Certification Body