Chromosome *kit P*

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Protocol of use

Protocolo de utilización

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Complete medium for culture of peripheral blood cells

Material supplied

EKAMTP, 10 x 5 ml culture tubes
EK AMTB 100, 1 bottle 100 ml
EK AMTB 500, 1 bottle 500 ml

Product stability

Chromosome Kit and Medium P is shipped at +2/+8°C. Once received, the product has to be stored at +2/+8°C. The product shelf life is 12 months from manufacturing date.

English

Protocol of use

Culture preparation

1. Aseptically add up to 0.5 ml of whole blood to one test tube or dispense 5 ml of bottled medium in a T-Flask and add the sample.
2. Replace the cap of the culture tube.
3. Mix the culture tube (or the T-Flask) contents and incubate at $37 \pm 2^\circ\text{C}$ placing the tube horizontally (on its flat surface). Addition of CO_2 is not necessary as the medium composition is designed for optimal lymphocyte proliferation in a closed system.

Chromosome fixation and slides preparation

1. After 72 hours of culture, add 50 μl of Colcemid (10 $\mu\text{g/ml}$) and incubate at $37 \pm 2^\circ\text{C}$ for 90 minutes (approx.). This last incubation time may be reduced to 15 minutes if the aim is obtaining a higher number of prometaphases.
2. After incubation, transfer the culture in a tube, if necessary, and centrifuge for 3 - 4 minutes at 2.000 rpm (higher speeds will not damage the cells).
3. Discard the supernatant, taking care not to lose part of the cell pellet (it is recommended to leave approximately 0.5 ml of material in the tube).
4. Vigorously resuspend the pellet.
5. Add 5 ml of hypotonic solution (distilled H_2O , 0,075M KCl) and mix adequately.
6. Incubate for at least 7 minutes at room temperature.
7. Centrifuge at 2.000 rpm for 3 - 4 minutes.
8. Discard the supernatant according to step 3.
9. Vigorously resuspend the pellet.
10. Add 5 ml of Ibraimov solution (distilled H_2O , 5% acetic acid) and mix adequately.
11. Immediately centrifuge at 2.000 rpm for 3 - 4 minutes.
12. Discard the supernatant according to step 3.
13. Vigorously resuspend the pellet.
14. Add 5 ml of fixative solution (methanol or ethanol: acetic acid, 3 : 1) and mix adequately.
15. Repeat steps 11-14.
16. Immediately centrifuge at 2.000 rpm for 3 - 4 minutes.
17. Carefully discard almost all the supernatant.
18. Add a few drops of freshly prepared fixative, taking into consideration that the number of drops should be appropriate to the cell concentration in the pellet.
19. Drop a small aliquot of the cell suspension onto a clean and wet microscope slide.
20. Evaporate under constant condition of right temperature and relative humidity into Optichrome (Ref. EKAMH950).

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NOTE

This image shows a full page of blank, lined paper. It features approximately 20 horizontal blue lines spaced evenly across the page, typical of notebook paper. The lines are thin and light blue, set against a plain white background. There is no handwriting or other markings on the page.

This image shows a full page of blank, lined paper. It features approximately 20 horizontal blue lines spaced evenly across the page, typical of notebook or primary writing paper. The lines are thin and light blue, set against a plain white background. There are no margins, text, or other markings on the page.

EuroClone®
serving science through innovation

EuroClone®
serving science through innovation

EuroClone S.p.A.

Via Figino 20/22 - 20016 Pero (MI) Italy

☎ +39 02 38195.1 📠 +39 02 38101465

✉ info@euroclone.it 🌐 www.euroclone.it