# STORZ

**KARL STORZ — ENDOSKOPE**

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our products, this product too is the result of years of experience and careful manufacturing. You and your organization have decided in favor of a modern high quality product from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper installation, connection and operation of the ENDOFLATOR® 50. All of the necessary details and actions are clearly explained.

Please read these instructions carefully. Keep this manual in a convenient and conspicuous location and in its protective casing close to the device.

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**RxOnly** Federal (USA) law restricts this device to sale by or on the order of a physician.

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
**It is recommended to check the suitability of the product for the intended procedure prior to use.**

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**2 Controls, displays,  
connectors, and their uses**

- ① Power switch ("O" = off)
- ② Touch screen
- ③ Insufflation connection to patient
- ④ Connection for patient tube heating
- ⑤ Room temperature sensor
- ⑥ Gas connection (American connection)
- ⑦ SCB connectors
- ⑧ Service interface (concealed)
- ⑨ Ethernet interface (concealed)
- ⑩ Potential equalization connector
- ⑪ Power cord socket
- ⑫ Line fuse holder
- ⑬ Holder for CO<sub>2</sub> bottle  
(optional – cat. no. UI 005)

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	Follow instructions for use
	Power on
○	OFF
	Potential equalization connector
	Applied part type CF
~	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).
	Manufacturer

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**Symbols on label and packaging of the device**

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. This can be downloaded from [www.karlstorz.com](http://www.karlstorz.com).

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**3.2 User interface symbols**

	Settings
	Procedure list
	Working area
	Start/Stop
	PEDIATRICS Mode
	HIGH FLOW Mode
	Gas bottle
	House connection
	Flow
	Pressure
	Alarm audio paused (30 s)
	Overpressure valve deactivated
	House connection, input pressure acceptable
	House connection, input pressure too high or too low
	Gas bottle; input pressure above 30 bar
	Gas bottle; input pressure between 20 and 30 bar
	Gas bottle; input pressure between 10 and 20 bar
	Gas bottle; input pressure below 10 bar
	Gas bottle; input pressure at 1 bar
	Gas bottle; input pressure too high
	Display of gas consumption (4 digit) / Reset the display by tapping on the -0- button
	Patient heating plugged in and in operation
	Defective patient heating/heating control

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**3.3 Gas filter label symbols**

	Manufacturer
	Catalogue number
	Do not reuse
	Do not resterilize
	CE mark
	Do not use if package is damaged
	Consult instructions for use
	Federal (USA) law restricts this device to sale by or on the order of a physician.
	Keep dry
	Keep away from sunlight
	Sterilized using ethylene oxide
	Use-by date
	Batch code

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<b>1</b>	<b>Images of the equipment</b> .....	<b>IV</b>
<b>2</b>	<b>Controls, displays, connectors, and their uses</b> .....	<b>V</b>
<b>3</b>	<b>Symbols employed</b> .....	<b>VI</b>
3.1	Symbols on the device .....	VI
3.2	User interface symbols .....	VII
3.3	Gas filter label symbols .....	VIII
<b>4</b>	<b>General information</b> .....	<b>4</b>
4.1	Description of the device .....	4
4.1.1	Property rights .....	5
<b>5</b>	<b>Safety instructions</b> .....	<b>6</b>
5.1	Explanation of warnings and cautions.....	6
5.2	Intended use.....	11
5.2.1	Indication .....	11
5.2.2	Contraindications.....	13
5.2.3	Warnings .....	13
5.3	User qualification .....	18
5.4	User profile of physician and assistants ..	18
5.5	Intended conditions of use.....	19
5.5.1	Use.....	19
5.5.2	Other intended conditions.....	19
5.5.3	User position .....	19
5.6	Safety precautions at the installation site .....	20
5.7	Safety precautions when operating the device.....	20
5.8	Safety features.....	21
5.8.1	Self-test .....	21
5.8.2	Monitoring during operation.....	21
<b>6</b>	<b>Installation and operating instructions</b> ..	<b>23</b>
6.1	Unpacking the equipment .....	23
6.2	Basic equipment.....	23
6.3	Installing and connecting the device .....	23
6.3.1	Connecting the ground line.....	24
6.3.2	Connecting the power cord .....	24
6.3.3	KARL STORZ SCB .....	25
6.3.4	Affixing the CO <sub>2</sub> bottle to the device .....	25
6.3.5	Connecting the CO <sub>2</sub> bottle.....	26
6.3.6	CO <sub>2</sub> bottles equipped with German-standard or ISO connection .....	26

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## 4 General information

### 4.1 Description of the device

The ENDOFLATOR® 50 is an insufflation device for universal application in laparoscopic and thoracoscopic examinations and operations as well as endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy). In addition, the ENDOFLATOR® 50 is designed to expulse ambient air in open and endoscopically assisted cardiac surgery as well as endoscopic vessel harvesting. The very latest technologies for pressure and flow measurement and control allow for various operating modes, including modes tailored to specific situations.

The 'PEDIATRICS' and 'HIGH FLOW' modes enable various applications to be accommodated. Easy handling and clearly arranged controls together with several safety circuits ensure maximum patient safety.

For particularly sensitive applications, the device features a 'PEDIATRICS' mode which is characterized by low insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow ranges (up to max. 15 l/min).

This mode also includes tighter control over the flow rate at low flows.

To quickly compensate for the considerable loss of gas occurring in complex laparoscopic operations, the ENDOFLATOR® 50 has been designed for a high flow rate of up to 50 l/min.

Use of an insufflation tube with integrated gas heater can reduce the risk of hypothermia.

The ENDOFLATOR® 50 is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface.

The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

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## 5 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 5.1 Explanation of warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words Warning, Caution, and Note convey special meanings. Wherever they are used in this instruction manual, the following text must be carefully reviewed to ensure the safe and effective operation of this device. In addition, to make the signal words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** A Note provides special information regarding the operation of the product, or clarifies important information.

**WARNING:** Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk.

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**Safety instructions**  
Warnings and cautions

-  **WARNING:** The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.
-  **WARNING:** Set up the device out of reach of patients.
-  **WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if
  - they are identified as such in the respective instruction manuals or
  - the intended use and interface specifications of the devices used in combination permit this (cf. Clause 16 of IEC 60601-1).
-  **WARNING:** Grounding reliability can only be achieved when the device is connected to a 'Hospital only' or 'Hospital Grade' outlet. Check the plugs and cables regularly and do not use if damaged.
-  **WARNING:** To avoid the risk of an electric shock, this device may only be connected to a power supply network with a protective conductor.
-  **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.
-  **WARNING:** For safety reasons, do not work with a leaky device. There is a risk that the intraabdominal pressure will increase in an uncontrollable manner.
-  **WARNING:** Always rinse out the system with 1 l CO<sub>2</sub> before use.

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## Safety instructions

### Warnings and cautions

-  **WARNING:** Make sure there are no external forces acting on the cavity. This can increase the intracavitary pressure and/or cause pressure fluctuations.
-  **WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.
-  **WARNING:** Always use the lowest possible pressure and flow for insufflation.
-  **WARNING:** Always check which operating mode is set. For safety reasons, use the 'PEDIATRICALS' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.
-  **WARNING:** Any treatment may only be performed if visual observation of the action of the device is ensured.
-  **WARNING:** After 24 h of continuous operation, a self-test must be performed. Take note of the warnings to ensure that the device is fully functional.
-  **WARNING:** The device must not be used if either the touch screen or display is defective.
-  **WARNING:** If the touch screen does not react when pressed, you must turn the device off.
-  **WARNING:** Do not press several points on the touch screen simultaneously.
-  **WARNING:** Only sterile CO<sub>2</sub> gas may be used.
-  **WARNING:** You must secure CO<sub>2</sub> bottles connected up to the device to prevent them from falling over.
-  **WARNING:** Follow the safety instructions when changing the CO<sub>2</sub> bottle.
-  **WARNING:** The CO<sub>2</sub> bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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**Safety instructions**  
Warnings and cautions

-  **WARNING:** Always unplug the device from the power supply before carrying out any cleaning and maintenance work.
-  **WARNING:** Danger of electric shock! Do not open the device. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unauthorized persons will void warranty claims.
-  **WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.
-  **WARNING:** Only operate the device under the specified environmental conditions.
-  **WARNING:** Always keep a spare device to hand in case the first device should fail.
-  **WARNING:** For safety reasons, do not simultaneously touch the device output sockets and the patient.
-  **WARNING:** The device may only be operated with patient tube heating (heatable insufflation tube sets) supplied by KARL STORZ. The connection socket ④ on the front panel is intended solely for connection of the patient tube heating.
-  **WARNING:** If the patient heating symbol lights up red, this indicates that the control unit of the patient tube heating is defective. In this case, do **not** connect the cable for the patient tube heating to the device.

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## Safety instructions

### Warnings and cautions

-  **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.
-  **WARNING:** Always use a hydrophobic bacterial filter to prevent cross-contamination of patients.
-  **WARNING:** You must insert a sterile CO<sub>2</sub> filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. Replace the CO<sub>2</sub> filter after every use.
-  **WARNING:** Ensure that liquid cannot flow back into the device. Position the device higher than the patient and use a hydrophobic filter between the insufflation connection and insufflation tube.
-  **WARNING:** Check the packaging of the tube set for damage! If the packaging is damaged, the sterility of the tube set may be impaired and the tube set must not be used.
-  **WARNING:** Follow the country-specific regulations/laws for the disposal of accessories.
-  **WARNING:** Avoid allowing fluids to enter the housing. Do not store liquids on or directly above the device. If, despite the precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies for the formation of condensation).
-  **WARNING:** The device may only be operated at the voltage stated on the rating plate.

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-  **WARNING:** When replacing fuses, use only fuses of the correct rating.
-  **WARNING:** Applied parts of devices used in combination with this device must be of type BF or CF.
-  **WARNING:** Patient leakage currents may be additive when energized endoscopes are used with energized endotherapy devices. This is particularly important if a CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.
-  **WARNING:** The compatibility of this device with any accessories and/or energized endotherapy devices should be checked before each use according to the criteria for safe use defined in the individual instructions for use.

Before using the device on the patient for the first time it is imperative that you be acquainted with how the device operates and is controlled.

## 5.2 Intended use

CO<sub>2</sub> insufflators and their accessories are used to create and maintain a cavity during diagnostic or therapeutic endoscopic interventions as well as to expulse ambient air during open, endoscopically assisted or endoscopic interventions.

### 5.2.1 Indication

The insufflator is used to create a cavity in the following diagnostic and therapeutic interventions:

- Laparoscopy
- Thoracoscopy
- Endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy)
- Endoscopic vessel harvesting

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In addition, the ENDOFLATOR® 50 is also designed to expulse ambient air in open and endoscopically assisted cardiac surgery.

 **WARNING:** Only sterile CO<sub>2</sub> gas may be used.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

The ENDOFLATOR® 50 may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

 **WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.

 **WARNING:** Unauthorized conversions or modifications to the device are not allowed for safety reasons.

 **WARNING:** Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see Clause 16 of IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

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### 5.2.2 Contraindications

The use of an insufflator is contraindicated if, in the opinion of a qualified physician, such an application would pose risks to the patient, e.g. due to the patient's general condition.

The device is contraindicated for hysteroscopic applications. Use of this device for intraabdominal distension is contraindicated whenever laparoscopy is contraindicated. See the operator's manual of your laparoscope for absolute and relative contraindications.

### 5.2.3 Warnings

#### 5.2.3.1 Warnings

##### Excessive Absorption of CO<sub>2</sub> from Excessive Pressure/Flow

Excessive absorption of CO<sub>2</sub> results from either excessive flow and/or excessive pressure. The abdomen can be adequately distended by pressure in the range of 15-20 mmHg. It is seldom necessary to use an abdominal pressure greater than 20 mmHg. Little intravasation should occur at these levels. Pressures over 20 mmHg are virtually never needed and will increase the amount and rapidity of intravasation. Adequate respirations help avoid problems related to CO<sub>2</sub>.

##### Metabolic Acidosis and Resultant Cardiac Irregularity

Prolonged intra-abdominal pressures greater than 20 mmHg should be avoided. This can cause any of the following:

- Decreased respiration with compromised diaphragmatic excursion
- Decreased venous return
- Decreased cardiac output
- Acidosis

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**Embolism**

Incorrect positioning of the insufflation instrument or a high intraabdominal pressure can lead to a CO<sub>2</sub> embolism. Avoid high pressures and verify the correct positioning of the insufflation instrument.

**Hypothermia**

The CO<sub>2</sub> insufflation can lead to a decrease in body temperature. Monitor the body temperature throughout the intervention.

**Pneumoperitoneum**

When a vessel is harvested from the leg of a patient with a perforated groin, it is possible for CO<sub>2</sub> to reach the abdomen and cause a pneumoperitoneum. Make sure the abdomen does not fill with CO<sub>2</sub> during surgery.

**Idiosyncratic reactions**

The risk of a metabolic imbalance due to increased CO<sub>2</sub> absorption is increased for patients with sickle cell anemia or pulmonary insufficiency.

**Dehydration**

For longer operations with high gas levels, insufflation can lead to a drying out of the tissues or cause tissue damage. Unnecessary leaks must be avoided for this reason.

The decisive factor for the incidence of these complications is the length of the pneumoperitoneum, the level of the intraabdominal pressure and the CO<sub>2</sub> absorption associated with it. Increased pressure and flow values must thus be avoided. Pressures in excess of 15 mmHg are only required in isolated cases.

**Additional Insufflation Sources**

The use of additional insufflation sources increases the intra-abdominal pressure, while the venting rate of the automatic venting system is limited. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

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**5.2.3.2 Warnings: PEDIATRICS MODE****Hypothermia: PEDIATRICS MODE**

The insufflation gas flow usually drops significantly after the target pressure has been reached and it is then only required to maintain the abdominal pressure.

However, leaks within the abdomen or the instrument can lead to a constant gas flow of above 1 l/min. When operating on children younger than 12, a gas flow of more than 1 l/min poses an increased risk of hypothermia for the patient. Corresponding measures to prevent hypothermia include the use of blankets or prewarmed gas. The patient's body temperature has to be monitored at all times during surgery.

**Pneumolabium/pneumoscotum:  
PEDIATRICS MODE**

Children are at risk of a pneumolabium or pneumoscotum.

**Increased airway pressure:  
PEDIATRICS MODE**

When laparoscopic procedures are performed on children, the increased intraabdominal pressure also increases the risk for higher airway pressures. Always strictly monitor respiration and airway function when performing laparoscopic procedures on children younger than 12 years of age.

**Compression of the vena cava:  
PEDIATRICS MODE**

When insufflating the abdomen of a child with medical CO<sub>2</sub>, an increased risk of compressing the vena cava exists. This risk can be reduced by monitoring the systolic and diastolic blood pressure during the entire surgery.

**Recommended work settings:  
PEDIATRICS MODE**

The flow values listed\* for laparoscopic procedures performed on newborns, infants, and children are only suggested values. The selection of the suitable flow and pressure values is solely the responsibility of the attending physician. However, adhering to the values listed\* ensures an optimal performance of the PEDIATRICS MODE of the insufflator.

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**Gas flow limit: PEDIATRICS MODE**

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg (approx. 55 US pounds).

**Hemodynamic stability:  
PEDIATRICS MODE**

A laparoscopy performed on children younger than 12 years of age can result in the phenomenon of the increased CO<sub>2</sub> content in the blood and with that to problems of the haemodynamic system. It is recommended to increase the breathing rate of the patient and to work with low flow values and pressure values not exceeding 12 mmHg. The patient's circulatory system should be monitored at all times.

**5.2.3.3 Warnings: Endoscopic Vessel  
Harvesting****CO<sub>2</sub> absorption: Endoscopic Vessel  
Harvesting**

Due to the special surgical procedures - start of the heart bypass operation, and the endoscopic removal of the vessel - special care has to be taken as CO<sub>2</sub> is always absorbed through the tissue of the patient during insufflation (intravasation). This means the body absorbs part of the CO<sub>2</sub> gas used for insufflation. CO<sub>2</sub> concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO<sub>2</sub>. High pressure or a high gas flow promotes CO<sub>2</sub> absorption.

**Instruments: Endoscopic Vessel Harvesting**

Before using the insufflator to endoscopic harvest vessels, please check whether the used instrument is intended for CO<sub>2</sub> insufflation.

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**Metabolic and Cardiac Reactions: Endoscopic Vessel Harvesting**

Due to the special surgical conditions - start of the heart bypass surgery and vessel harvesting - it is especially important to remember the existing risk of metabolic acidosis when insufflating with CO<sub>2</sub>. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**Dehydration: Endoscopic Vessel Harvesting**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**Embolism: Endoscopic Vessel Harvesting**

Improper placement of the insufflation instrument could cause insufflation of gas into a vessel, resulting in air or CO<sub>2</sub> embolisms. To reduce the risk of air or CO<sub>2</sub> embolism, perform initial insufflation at a low flow rate and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO<sub>2</sub> embolisms can also be caused by a high pressure. Avoid high-pressure settings and close damaged blood vessels at once.

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**5.3 User qualification**

The ENDOFLATOR® 50 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

**5.4 User profile of physician and assistants**

Recognized medical skills in the field of application (specialist physician, qualified medical staff).

Adequate powers of comprehension to rationally assess the situation in hand.

Familiarity with the use of technical apparatus

Suitable command of one of the languages used by the device and in the instruction manual.

Must have been trained thoroughly at least once in the use of the device

No physical impairments which limit or prevent perception of acoustic or visual alarm signals.

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## 5.5 Intended conditions of use

### 5.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

### 5.5.2 Other intended conditions

- Frequency of use: one or more times a day
- Length of use: from a few minutes to several hours a day.
- Place of installation: positioning on a level, vibration-free surface
- Mobility: can be moved if positioned on a cart
- Combination: can be used on the patient at the same time as other equipment required for the operations
- Control: can be controlled via the KARL STORZ SCB.
- The device is placed outside the sterile area.

### 5.5.3 User position

- To operate the device, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of  $\pm 45^\circ$ , at a distance of approx. 30-70 cm to the front panel.
- For observation of the actual values during the application, a visual distance from the device of 2 m is assumed.

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
**5.6 Safety precautions at the installation site**

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

**! WARNING:** **DO NOT** use the device in the presence of flammable anesthetics.

**! WARNING:** The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

**! WARNING:** The device must not be used in oxygenated environments.

This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

**5.7 Safety precautions when operating the device**

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

This device has been tested and found to comply with the EMC limits of the Medical Device Directive 93/42/EEC (EN 55011 Class B and IEC 60601-1-2).

Please observe the Electromagnetic Compatibility information in section 10.

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### 5.8 Safety features

The ENDOFLATOR® 50 has the following safety features:

#### 5.8.1 Self-test

A self-check procedure is performed each time the device is switched on (for this the device must be supplied with gas). If the test is negative, a corresponding error message is shown and the device can no longer be used as normal.

The following components are tested:

- Leakproofness of the system (valves)
- Checking pressure measurement
- Checking flow measurement
- Functionality of the patient tube heating control

#### 5.8.2 Monitoring during operation

- **Overpressure:** If, during insufflation, the intracorporeal pressure rises above the set setpoint value by a limit dependent on the setpoint value for more than 3 seconds, an acoustic signal sounds (see chapter 8.2.2) and the gas supply is stopped. After a time period of 5 seconds, the overpressure is actively reduced via a blow-off valve or shown by a symbol that it is deactivated.
- **Mode-specific pressure monitoring**
- **Monitoring of gas supply**  
If the system recognizes a lack of CO<sub>2</sub> gas when switched on, a warning signal is given. If a missing gas supply is detected and insufflation is stopped, because the gas bottle is empty, a visible and audible alarm is monitored.
- **Checking gas preheating**  
The internal gas preheating is monitored by a second circuit, independently of the control. If a fault arises, the heating is switched off.

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- **Checking the gas temperature to the patient**  
The gas temperature and the energy output are recorded in the device. If a system error occurs, the heating will be switched off. The color of the symbol for the patient tube heating changes from white to red. Insufflation with (cold) CO<sub>2</sub> is, however, still possible.  
The device thus ensures that the gas temperature does not exceed 41 °C.
- **Passive safety measures**
  - The pressure and flow measurement systems have a redundant stand-by.
  - Both pressure control stages and the flow measurement are equipped with safety pressure release valves.
  - The function of the mechanical high pressure regulator is monitored.

**i NOTE:** If the monitoring system recognizes an error which renders safe functioning seemingly impossible, the device is transferred into a "safe state". From this state, the device can only be reset by turning off the power supply.

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**6 Installation and operating  
instructions**

**6.1 Unpacking the equipment**

Carefully remove the ENDOFLATOR® 50 and the accessories from their packaging. Check the delivery for missing items and possible damage.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

**6.2 Basic equipment**

**KARL STORZ Set no.: UI500S1**

- 1 ENDOFLATOR® 50 UI 500
- 1 Power cord 400 A
- 3 Heatable insufflation tubing sets with gas filter, single use, sterile (1x 031210-03, package of 3)
- 1 Universal wrench 20 4000 30
- 1 SCB connecting cable 20 0901 70, length 100 cm
- 1 Instruction manual 96116045D
- 1 CO<sub>2</sub> supply tube (customized)

**6.3 Installing and connecting the  
device**

 **NOTE:** The ENDOFLATOR® 50 including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

 **WARNING:** Position the device higher than the patient.

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1. Place the device on a flat surface (Fig. 1).

**!** **WARNING:** Ensure that there is sufficient air circulation for the device.

**!** **WARNING:** The device must be set up in such a way that the power cord can be simply pulled out at any time.

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**6.3.1 Connecting the ground line**

**!** **WARNING:** The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.

2. Connect the potential equalization cable to the potential equalization connector (Fig. 2).  
3. Where necessary, have grounding performed by qualified personnel.

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**6.3.2 Connecting the power cord**

**!** **WARNING:** The device may only be operated at the voltage stated on the rating plate.

4. Push the power cord into the power socket as far as it will go (Fig. 3).

**!** **WARNING:** Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

**!** **WARNING:** The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

**!** **WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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**6.3.3 KARL STORZ SCB**

**i NOTE:** To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Therefore, to unplug the SCB cable, pull out the connector – do not pull directly on the cable!

Connect the SCB cable to an (equivalent) SCB socket on the back of the device (Fig. 4).

Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ SCB control NEO System Instruction Manual).

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**6.3.4 Affixing the CO<sub>2</sub> bottle to the device (optional)**

Attach the optionally available holder to the rear of the device (Fig. 5). Use the enclosed screws of length 6 mm.

Insert the CO<sub>2</sub> bottle vertically into the holder and close the locking clip.

**! WARNING:** Only sterile CO<sub>2</sub> gas may be used.

**i NOTE:** When using the ENDOFLATOR® 50 with a CO<sub>2</sub> bottle, KARL STORZ recommends the additional use of a high pressure inline filter (20 4000 32) between the device and the CO<sub>2</sub> bottle. This filter is attached to the device's gas connection and if necessary, prevents particles from the bottle from entering the device (see section 6.3.5). Contamination due to unclean gases can lead to the failure of the device.

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**6.3.5 Connecting the CO<sub>2</sub> bottle**

Connect the gas outlet of the CO<sub>2</sub> bottle to the device's gas connection using the high-pressure tube supplied (Fig. 6-7).

**WARNING:** CO<sub>2</sub> bottles connected up to the device must be secured to prevent them from falling over.

**WARNING:** Follow the safety instructions when changing the CO<sub>2</sub> bottle.

**WARNING:** The CO<sub>2</sub> bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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**6.3.6 CO<sub>2</sub> bottles equipped with German-standard or ISO connection**

Place the high-pressure tube (e.g. 20 4002 22; see the section 'Accessories') onto the outlet of the CO<sub>2</sub> bottle and screw tight using the universal wrench provided (Fig. 8).

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**6.3.7 CO<sub>2</sub> bottles equipped with PIN-Index connection**

Place the connecting piece of the high-pressure tube (e.g., 20 4000 22; see the section 'Accessories') onto the CO<sub>2</sub> bottle and screw this tight.

Use the wrench provided (see the illustration) to open the CO<sub>2</sub> bottle (Fig. 9).

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**6.3.8 Opening the valve of the CO<sub>2</sub> bottle, if necessary.**

Slowly open the valve stopcock of the CO<sub>2</sub> bottle by turning it approx. one-half turn counterclockwise (Fig. 10-11).

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**6.3.9 Connecting to the central gas  
supply**

Connect the low-pressure tube (e.g., UI001; see the section 'Accessories') to the device's gas connection and connect the other end directly to the wall outlet of the central gas supply (3.3...7 bar) (Fig. 12).

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**6.4 Commissioning**

**i** **NOTE:** The gas supply must be connected when the device is switched on, otherwise an error message is emitted.

**6.4.1 Commissioning for the first time**

The ENDOFLATOR® 50 is delivered as follows:

- (Screen) Language: English
- Overpressure valve on (default setting each time the device is switched on)
- Default Settings  
PEDIATRICS Pressure: 3 mmHg Flow: 0.1 l/min  
HIGH FLOW Pressure: 3 mmHg Flow: 1 l/min

**i** **NOTE:** Changing/adjustment of the settings for language, volume and the gas source are described on the following pages.

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1. Press the power switch ①.

➤ Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo (Fig. 13). After this, a white screen appears.

The device is now running a self-test. Once the test has been successfully completed, an acoustic availability signal sounds and a mode selection screen appears (Fig. 14).

2. Select the desired mode.

- 'PEDIATRICS' (see Fig. 18, page 29)
- 'HIGH FLOW' (see Fig. 19, page 29)
- Procedure list (see Fig. 20, page 29)

**i NOTE:** Selection of the desired mode is confirmed by an audio sound.

**! WARNING:** For safety reasons, the device may only ever be used if, immediately after selecting the mode, the acoustic availability signal sounded!  
See section 8.1.5.

### In the event of an error

In the event of an error, the selected screen appears with an error message (shown here for the HIGH FLOW Mode by way of example – Fig. 15/16).

If no gas is connected or if the input pressure is too low, the warning 'Input pressure too low' appears (Fig. 15).

At that point, check whether the bottle is open or must be replaced.

If the gas supply was secured with a delay, the switch-on test is delayed and the user is notified via message 110 that the tests are still active (Fig. 16).

**i NOTE:** Normal use is only possible after successful conclusion of the tests.

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If the self-test fails, a screen appears with an information message and an acoustic information signal sounds.


**NOTE:** *If restarting (turning the device off and on again) fails, verify that the device is properly connected and the preparations have been made according to the instructions.*

Further information on troubleshooting can be found in the chapter 'Technical description - Troubleshooting'.

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**In normal cases**


The switch-on message is followed by the selected screen PEDIATRICS, HIGH FLOW or Procedure list (Fig. 18-20).

**NOTE:** *On the work screens 'PEDIATRICS' (Fig. 18) and 'HIGH FLOW' (Fig. 19), insufflation can be started directly, device settings can be changed, and the procedure information can be called up. On the 'Procedure list' screen (Fig. 20), a saved procedure can be selected (factory default 'empty'; see section 6.5.1. for details on creating and editing procedures).*

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
**Changing settings**


3. Press button ① ('wrench') (see Fig. 18 or 19) to change to the settings mode.

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**Changing the language and other presets**

- i** **NOTE:** The screen language is preset as English.
- ☛ The **Settings** screen appears (Fig. 21).
4. Press the Device Settings button.

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
- ☛ The **Device Settings** screen appears with the list of available user settings (Fig. 22).

- Language
- Source of gas supply
- Volume
- Overpressure valve activation

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
5. Select your preferred language under **Language** (Fig. 23).
- Confirm the selection with **15** or reject the changes with **14**.

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
6. Select the desired gas supply source under **Source of gas supply** (Fig. 24).

- i** **NOTE:** With 'Bottle', you can choose between a single audio signal and a repetitive audio signal (if the bottle pressure is < 10 bar).
- Confirm the selection with **15** or reject the changes with **14**.

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
7. Select the desired volume for warnings and key tones under **Volume**. (Fig. 25)  
Warnings setting range: 1...4  
Key tones setting range: 0...4

Confirm the selection with (15) or reject the changes with (14).

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8. Select the desired setting for the overpressure valve under **Overpressure valve activation** (On or Off) (Fig. 26).

Confirm the selection with (15) or reject the changes with (14).

**(i) NOTE:** When the overpressure valve is activated, the overpressure is reduced after 5 seconds.

When the overpressure valve is deactivated, a corresponding symbol is shown on the screen (see (16) in Fig. 27).

**(i) NOTE:** The overpressure valve is active again when the device is restarted.

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


#### 6.4.2 Normal commissioning

After a successful start (see section 'Commissioning for the first time'), one of the screens shown will be displayed depending on your selection (Fig. 28-30).

The device saves the operating parameters which were last used and reloads these (with a few exceptions) when the device is restarted (see tables below).

HIGH FLOW mode	Setting	Setting after restart
Pressure	1...20 mmHg	1...20 mmHg
	21...30 mmHg	20 mmHg
Flow	1...20 l/min	1...20 l/min
	≥ 21 l/min	20 l/min

PEDIATRICS mode	Setting	Setting after restart
Pressure	1...12 mmHg	1...12 mmHg
	13...15 mmHg	12 mmHg
Flow	0.1...1 l/min	0.1...1 l/min
	≥ 1.1 l/min	1 l/min

**i NOTE:** On the work screens 'PEDIATRICS' (Fig. 28) and 'HIGH FLOW' (Fig. 29), insufflation can be started directly, device settings can be changed, and the procedure information can be called up. On the 'Procedure list' screen (Fig. 30), if necessary, a saved procedure can be selected or a new one can be created. 30 user procedures can be saved for each device (detailed information on 'Procedures' can be found in section 6.5.1).

**i NOTE:** The user procedure list is arranged according to when the procedures were last used. The most recently used procedure is shown at the top of the list.

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1. Attach the sterile CO<sub>2</sub> gas filter to the insufflation connection ③ (Fig. 31).

**! WARNING:** You must insert a sterile CO<sub>2</sub> filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. The CO<sub>2</sub> filter must be replaced after each use.

2. Attach the insufflation tube to the CO<sub>2</sub> gas filter (Fig. 31) and connect the other end of the insufflation tube to the VERESS needle.

**i NOTE:** In the case of the disposable insufflation tube set, the filter is firmly connected to the tube.

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### 6.4.3 Heating up the CO<sub>2</sub> gas

The ENDOFLATOR® 50 enables the CO<sub>2</sub> gas to be preheated to body temperature. The risk of hypothermia during the insufflation of sizeable quantities of CO<sub>2</sub> can thus be reduced.

Use the heatable insufflation tube, Art. no. 031210-10 (pack of 10).

### Connecting the heatable patient tube

Connect the connecting cable of the heatable patient tube to the socket ④ on the ENDOFLATOR® 50 (Fig. 32).

The heating function is activated by plugging in the connecting cable.

**i NOTE:** The patient tube heating symbol ⑦ is displayed in the **Working screen** (Fig. 33).

The CO<sub>2</sub> gas flowing through is heated to a temperature of approx. 37 °C. The unit adapts to various flow rates automatically.

**i NOTE:** If the connecting cable is disconnected, the ENDOFLATOR® 50 can be operated without heating.

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**6.4.4 Test for proper functioning**

**WARNING:** Test the device prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.

**Pressure stability and leakproofness**

1. Set the mode to 'PEDIATRICS'.
2. Use 12 mmHg and 1.0 l/min as setpoint values.
3. Check the connection of the patient tube to the device.
4. Close the stopcock on the instrument.
5. Start insufflation.

☛ Pressure then builds up in the tube which increases above the setpoint value and which, after a few seconds, is released via the blow off valve.

6. Observe the instrument's conduct over a 30 second period.

**NOTE:** The flow display may only show values greater than zero for short periods. If the set pressure is only reached slowly even when the instrument's stopcock is closed or if you note stationary flow rates of 0.2 l/min or more, there is a leak in the system filter – patient tube – instrument.

7. In this case, check the connections and, if necessary, change the tube set.

**WARNING:** Do not work with a leaky device. There is a risk that the intraabdominal pressure will increase in an uncontrollable manner.

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#### 6.4.5 Preparing insufflation

You have two possibilities for entering the desired operating parameters for insufflation mode\*, pressure and flow.

##### a) Manual selection

1. If desired, change the selected *working screen* ('PEDIATRICS'/'HIGH FLOW' – Fig. 34 or 35) via the respective button.
2. On the selected *working screen* ('PEDIATRICS'/'HIGH FLOW' – Fig. 34 or 35) set the setpoint values for pressure and flow using the ± buttons.

##### b) Selecting a saved procedure

1. If necessary, tap on button ② on *working screen* ('PEDIATRICS'/'HIGH FLOW').
  - The *Procedure list* screen appears (Fig. 36).
2. Select the desired procedure by tapping (e.g., ③). (By tapping on button ④, the device jumps back to the *working screen* without any changes being made).
  - The *Working screen* appears (Fig. 34 or 35).

\*The ENDOFLATOR® 50 is equipped with two different insufflation modes:

- HIGH FLOW ('man' symbol)
- PEDIATRICS ('baby' symbol)

**i NOTE:** To compensate for the considerable loss of gas occurring in complex operations, the HIGH FLOW mode has been fitted with a flow rate of up to 50 l/min.

For particularly sensitive applications (with small cavities), the device offers the 'PEDIATRICS' mode. This mode is characterized by a lower insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow range (up to max. 15 l/min).

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**WARNING:** Always check the mode set. For safety reasons, use the 'PEDIATRICS' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.

**6.4.6 Recommended settings for pediatric subcategories**

Group	Weight	Flow range
Children younger than 1 year	Approx. 1...9 kg	0.1...0.5 l/min
Children between 1 and 3 years	Approx. 10...15 kg	0.5...1.0 l/min
Children between 3 and 4 years	Approx. 16...19 kg	1.0...2.0 l/min
Children between 4 and 14 years	≥ 20 kg	> 2.0 l/min
All Children	< 25 kg	≤ 14.0 l/min

**6.4.7 Recommended settings for endoscopic vessel harvesting**

Group	Weight	Flow range
Adult endoscopic vessel harvesting	N/A	3.0...6.0 l/min

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**6.4.8 Performing CO<sub>2</sub> insufflation**

**! WARNING:** Before starting insufflation rinse the entire system with 1 l CO<sub>2</sub>.

**! WARNING:** Always use the lowest possible pressure and flow for insufflation.

1. Reset the volume display using the '0' button (Fig. 37).
2. Insert the VERESS needle and open the locking lever. Tap the START/STOP button to initiate the insufflation process (Fig. 37).

**i NOTE:** In order to use the maximum device flow in the HIGH FLOW mode, HiCap® instruments with a suitable connection must be used. The LUER-lock connection fitted to the end of the insufflation tube must be removed in this case.

**i NOTE:** In the case of disposable tubing sets, cut the tube off behind the LUER-lock.

**Displays**

The insufflation pressure, flow and volume of the gas used are continuously displayed in the **Working screen** (Fig. 37).

If the patient pressure rises above the preselected value, an audible warning signal will be sounded. The pressure is reduced via the overpressure valve after a delay of 5 s.

**i NOTE:** The tone can be switched off for 30 s using the 'Bell' button (Fig. 38).

**6.4.9 Decommissioning**

Close the valve on the CO<sub>2</sub> bottle on ending the intervention (when using a high-pressure gas supply) and remove the insufflation tube (incl. filter) used from the device.

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## 6.5 Menu explanations

### 6.5.1 Procedure list

Presaved values for mode, pressure and flow can be called up directly via the procedure list.

Tap on button ② (see Fig. 39 and Fig. 40) to call up the *Procedure list* screen (text is language-dependent) (if not only selected when starting the device) (see Fig. 41).

The *Procedure list* screen displays the saved procedures.

**NOTE:** The user procedure list is arranged according to when the procedures were last used. The most recently used procedure is shown at the top of the list.

Up to 30 procedures can be saved. Four procedures are displayed on each page. If there are more than four procedures, the buttons ⑥ are displayed for turning the pages.

Saved procedures can be deleted or modified, or saved as a new procedure.

**NOTE:** By tapping on the button ④, the device jumps back to the Working screen without any changes being made.

### Creating a new procedure

1. Tap on button ② (see Figs. 39 resp. 40) to call up the *Procedure list* screen if necessary (also see Section 7.4.2).

The *Procedure list* screen shows any saved procedures (Fig. 41).

2. Tap on button ③ (Create a new procedure – Fig. 41).

**NOTE:** If 30 procedures have already been created, the following information dialog appears when you press the button ③: "A maximum of 30 procedures can be saved in the procedure list. All memory locations are currently occupied. Before a new procedure can be saved, an existing procedure must be deleted from the list."

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The *Set procedure values* screen with the current values for mode, pressure and flow appears (see Fig. 42).

- Tap on the button 15 to confirm the set values. Alternatively, the desired values for mode, pressure and flow can be changed with the corresponding buttons 8, 9, 10, 11, 12 or 13.

**NOTE:** Tap on the button 4 or 14 to reject the changes and return to the *Working screen* or *Procedure list* screen.

The *Assign procedure name* screen appears (see Fig. 43).

- Enter a name for the procedure.

**NOTE:** The on-screen keyboard can type lower-case or capital letters, numbers and special characters depending on the status of the buttons 16 and 17. By pressing the arrow buttons 18 or 19, the cursor moves one character to the left or right. Use the button 20 to delete the character to the left of the cursor.

- Tap on the button 15 to confirm the input. When a new procedure is created, the device returns to the *Procedure list* screen (see Fig. 44).

**Editing an existing procedure**

- Tap on button 8 in the *Procedure list* screen (see Fig. 42).

The *Edit procedure list* screen appears (see Fig. 45).

- If necessary, select the procedure using the buttons 6 (procedure is marked orange).

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
- The *Set procedure values* screen with the current values appears (see Fig. 46).
- Set the desired values for mode, pressure and flow using the corresponding buttons ⑧, ⑨, ⑩, ⑪, ⑫ or ⑬.
- Confirm the set values by tapping the button ⑮ or reject the changed settings using ④ or ⑭.

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- The *Assign procedure name* screen appears (see Fig. 47).
- Enter a new name for the procedure or overwrite the old one.
- After pressing the button ⑮ to confirm, or canceling using the button ⑭, the *Edit procedure list* screen appears again (see Fig. 48).

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- Deleting an existing procedure**
- Call up the *Edit procedure list* screen in the same way as described for editing an existing procedure (see Fig. 45, page 35).
  - Select the procedure using the buttons ⑥ if necessary.
  - Tap the trash can button next to the procedure to be deleted.

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
- The *Delete procedure* screen appears (see Fig. 49).
- Confirm the deletion of the procedure by tapping the button ㉓ or cancel the process using the button ㉒.
- The *Edit procedure list* screen appears again (see Fig. 46).

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### 6.6 Settings

- Device Settings
- Service
- System Log
- Administration
- Device information

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#### 6.6.1 Device settings

The Device Settings (Fig. 51) are described in the section on commissioning for the first time.

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#### 6.6.2 Service

The *Service* screen is protected by a password (Fig. 52).

**i NOTE:** The service programs are only available to service personnel. The settings are described in the Service Manual.

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#### 6.6.3 System Log

The System Log contains the saved information and alarms (Fig. 53).

Each entry occupies one row and consists of a sequential number, date, time and alarm/information ID.

The most recent entry is allocated the number 000.

In case of voltage drops or shut downs, the system log is saved and contains entries concerning start-up and shut-down times.

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The system supports up to 50 000 entries.

Once the maximum capacity of entries is reached, the oldest ones are overwritten.

All alarms, user and service information is saved in the system log, showing the time of occurrence. They can be exported via the admin menu.

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### 6.6.4 Administration

- Export system log
- Unit pressure
- Date and Time
- Change administration password

**NOTE:** The Administration screen is protected by a password. On delivery the password is '2132'.

It is the responsibility of the user to change this if necessary.

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- Export system log

For service purposes, data can be exported via the service interface ⑧ using the USB stick.

1. Open the cover on the rear of the device with a Torx driver (size T 10 – Fig. 55).
2. Insert a USB stick into the USB interface (Fig. 56).
3. Select Export system log on the Administration screen (Fig. 54).

Exporting of the system log starts (Fig. 57).

**NOTE:** The current status is shown via a progress bar. This progress bar shows that the export is underway, however it does not show how far it has advanced.

4. Screw the cover back on once export is complete (Fig. 55).

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

- Unit pressure  
Select the desired unit (mmHg or hPa) under *Unit pressure* (Fig. 58).  
Confirm the selection with **15** or reject the changes with **14**.

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
- Date and Time  
The setting of date and time is self-explanatory (Fig. 59).  
Confirm the selection with **15** or reject the changes with **14**.

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- 6.6.5 Device information**  
The following elements are displayed (Fig. 60):
- Device serial number
  - Operating hours
  - Software versions
  - Software license information

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
**7 Maintenance**

**7.1 Fuse replacement**

1. Turn the device off and disconnect from the power supply (Fig. 61).

2. Remove the screw bits on the line fuse holder ⑫ using a screwdriver or another suitable tool (Fig. 62).

**! WARNING:** Only use fuses of the correct rating.

3. Insert new fuses (Fig. 63).

	100...240 V
Line fuse	2 x T2.5 AH 250 V

4. Now put the line fuse holder ⑫ on again.
5. Restore the power supply.
6. Switch the device on (Fig. 64) and ensure that the start screen appears.

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## 7.2 Reprocessing

- UI 500 – ENDOFLATOR® 50
- Heatable insufflation tubing set with gas filter, single use, sterile

### General warnings

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

**WARNING:** Risk of infection: The reusable insufflation tube is not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect the reusable insufflation tube for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the reusable insufflation tube before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**WARNING:** Risk of infection/product damage: KARL STORZ has validated the procedures given in the instruction manual. Only use these procedures. The reprocessing method must be validated by the user on site.

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**Maintenance**  
**Reprocessing**

-  **WARNING:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.
-  **CAUTION:** Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from [hygiene@karlstorz.com](mailto:hygiene@karlstorz.com).
-  **WARNING:** National laws and regulations must be observed.
-  **NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

**7.2.1 Reprocessing of the ENDOFLATOR® 50**

-  **WARNING:** Always disconnect the device from the mains before cleaning!
-  **WARNING:** Liquid must not be allowed to enter the housing under any circumstances.

**Manual wipe-down disinfection of device and supply tube**

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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**7.2.2 Reprocessing of heatable/non-heatable insufflation tubing set with gas filter, single use, sterile**

**WARNING:** Disposable (heatable) insufflation tubing sets with gas filter, locks and hydrophobic bacteria filters must not be reprocessed.

Sterile disposable products are not approved for reprocessing.

**7.2.3 Reprocessing of reusable insufflation tube**

**WARNING:** The reusable insufflation tube must be cleaned, disinfected and sterilized prior to use and subsequent reuse using validated reprocessing procedures.

**7.2.3.1 Accessories**

Accessories required for carrying out reprocessing:

Brushes: 27652, 27650 A, 27650 B

**7.2.3.2 Preparation for cleaning and disinfection**

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

**7.2.3.3 Manual precleaning Disassembly**

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

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**Immersion in cold water for 5 minutes**

To dissolve dried organic residues, immerse the medical device in cold water (10-20 °C) for 5 minutes.

**Brushing the surfaces**

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

**Rinsing out with a pressure water gun**

The working channel, lumina and hollow spaces must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

**7.2.3.4 Manual cleaning**

The medical device must be completely immersed in a cleaning solution. To ensure bubble-free wetting, the lumens must be filled. After the necessary exposure time, clean the instrument mechanically with the aid of brushes or a sponge. Finally, it must be rinsed with cold water to ensure neutralization.

**7.2.3.5 Manual disinfection**

The medical device must be completely immersed in a disinfectant solution. To ensure bubble-free wetting, the lumens must be filled. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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### 7.2.3.6 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D):

#### Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A<sub>0</sub> value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

**i NOTE:** If necessary, the instrument must be dried off afterwards by hand.

#### Connecting

In order to ensure effective machine cleaning and disinfection, as well as adequate rinsing out of the instrument, it must be connected up to the washer and disinfectant.

### 7.2.3.7 Assembly, inspection and care

**! WARNUNG:** A visual and functional check of the reusable insufflation tube must be carried out!

**! WARNUNG:** Check the reusable insufflation tube for leaks after cleaning and disinfection. By doing this it is possible to detect a leak early, through which gas could escape.

**i NOTE:** Discolorations impede visual inspections for possible contamination.

The cleaned and disinfected reusable insufflation tube must be visually inspected to ensure that it is clean, complete, dry and has not been damaged:

- If residues or contamination are still present, the reusable insufflation tube must be manually

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cleaned and subjected to a full cleaning and disinfection procedure once more.

- Damaged medical devices must be withdrawn from use.
- There should be no holes, cuts, dents or scratches on the tube surface. If the reusable insufflation tube is damaged as described above, it must be replaced. If the reusable insufflation tube is discolored as a result of the reprocessing process, it must be removed immediately from the sterile supply cycle.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

The integrity of the reusable insufflation tube can be tested in the following way:

Attach a syringe to one end of the tube and kink the other end. Place the insufflation tube in completely demineralized water or microbiologically pure, sterile water. Force air into the insufflation tube using the syringe. If any bubbles appear in the water while forcing air through the reusable insufflation tube, the tube must not be reused.

**i NOTE:** The leakage tester 13242 XL may also be used to check for leaks (Fig. 65).

To check for leakage use leak tester 13242 XL, for example. To do so, connect one tube end to the leakage tester and block off the other end, for example, with your finger. If there is a leak in the tube set, a pressure drop will be observed on the pressure gauge of the leakage tester.

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**7.2.3.9 Sterilization**

**WARNING:** Prior to the sterilization of the reusable insufflation tube, particularly in cases of manual reprocessing, it must be ensured that all chemical residue has been rinsed off or removed.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

**KARL STORZ has validated and approved steam sterilization using the fractionated prevacuum procedure for this medical device**

**i NOTE:** Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C - 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes

**i NOTE:** This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

- or 134°C - 137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes

**7.2.3.10 Limits of reprocessing**

The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

**i NOTE:** The storage period for the CO<sub>2</sub> supply tubes\* is 5 years. After this, the supply tube must be replaced.

\* Art. no. Low-pressure tubing:  
UI 001, UI 002 and UI 003  
Art. no. High-pressure tubing:  
20 4000 21, 20 4000 22, 20 4000 27,  
20 4000 28 and 20 4002 22)

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### 7.3 **Maintenance and safety check**

#### 7.3.1 **Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

#### 7.3.2 **Safety check/repeat inspection according to IEC 62353**

**WARNING:** *Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.*

#### **Visual inspection**

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

#### **Electric measurements**

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

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**Test for proper functioning**

A functional inspection must be carried out in accordance with the instruction manual (see section 6.4.4) as well as checking the alarm functions in sections 8.3.1/8.3.2.

**Documentation**

The safety check/repeat inspection and results must be documented.

**WARNING:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

**NOTE:** For detailed information on the range and implementation of the safety check/ repeat inspection, please refer to the latest issue of the relevant Service Manual.

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**7.4 Servicing and repair**

**WARNING:** Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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**7.5 Disposal**

This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).

**WARNING:** At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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### 7.6 Repair program

During the repair period, you will generally be loaned a device, which must then be returned to KARL STORZ as soon as the repaired device is returned to you.

For repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany, you can refer the repairs directly to:

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Take-off Gewerbepark 83  
78579 Neuhausen

Service hotline: +49 7461/708 980  
E-mail: technicalsupport@karlstorz.com

### 7.7 Infection prevention

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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**7.8 Limitation of liability**

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the device has been used in accordance with its operating instructions at all times.

**7.9 Warranty**

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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## 8 Technical description

### 8.1 Alarm specification

The device issues two alarms – ‘303: High pressure’ and ‘373: CO<sub>2</sub> empty’ (Fig. 66 and 67).

The alarm signals are not self-retaining. This means that the alarm signal is only issued whilst the signal conditions are present. The alarms are displayed for at least 5 s and at least one sequence is played.

All alarm thresholds and delays are permanently programmed. In case of an operating pause or a voltage loss, the most recently set volume is applied.

#### 8.1.1 Alarm thresholds Overpressure alarm

Medium priority overpressure alarm (‘303: High pressure’ – see 8.1) indicates that an overly high cavity pressure was measured at the device output. The overpressure alarm is activated 0.5 s after the cavity pressure measurement.

**PEDIATRICS mode:** The overpressure alarm is triggered if, as a product of the setpoint, the measured cavity pressure is 3 to 4 mmHg above the set setpoint value or if the cavity pressure is the same or larger than 16 mmHg.

**HIGH FLOW mode:** The overpressure alarm is triggered if, as a product of the setpoint, the measured cavity pressure is 4 to 5 mmHg above the set setpoint value or if the cavity pressure is the same or larger than 31 mmHg.

The device features an overpressure valve which reduces the overpressure after 5 s if it wasn’t deactivated and if insufflation is active.

**i NOTE:** The automatic reduction in pressure for setpoint values close to or at the maximum of the setting range is not initiated at 16 mmHg (PEDIATRICS) or 31 mmHg (HIGH FLOW) but only when a certain level of overpressure is reached (see table).

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	Setting	Reduction from an overpressure of
PEDIATRICS mode	> 12mmHg	4 mmHg
HIGH FLOW mode	> 26 mmHg	5 mmHg

If insufflation is turned off, the pressure is not reduced automatically.

 **WARNING:** Observe the functioning of the device when an overpressure alarm is issued. If the insufflation or the blow-off valve is deactivated, suitable measures to reduce pressure must be initiated.

The following table shows the overpressure thresholds with triggering to one decimal place as the threshold is calculated depending on the setpoint value.

	Pressure set-point value	Overpressure alarm threshold
PEDIA-TRICS mode	1 mmHg ... 15 mmHg	2.6 mmHg ... 4 mmHg or pressure $\geq$ 16 mmHg
HIGH FLOW mode	1 mmHg ... 30 mmHg	3.6 mmHg ... 5 mmHg or pressure $\geq$ 31 mmHg

**CO<sub>2</sub> empty alarm**

Low priority CO<sub>2</sub> empty alarm (‘373: CO<sub>2</sub> empty‘ – see 8.1) indicates that the gas supply failed during the application and that insufflation has been deactivated.

The following table shows the activation and deactivation thresholds for the CO<sub>2</sub> empty alarm.

Activation threshold	Input pressure < 1 bar
Deactivation threshold	Input pressure > 1,5 bar

 **WARNING:** Reconnect the gas supply in the event of a CO<sub>2</sub> empty alarm and restart insufflation.

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**8.1.2 Visual signals**

The alarms are shown in the title line (e.g., '373' in Fig. 68).

Medium priority alarms are displayed with blue font on a yellow background; low priority alarms are displayed with blue font on a cyan-colored background.

Medium priority alarms overwrite low priority alarms.

**8.1.3 Acoustic signals**


The acoustic alarm signals are specified below.

If a condition for an acoustic signal is met, at least one complete sequence is played. The volume of the audio signals can be adapted in 'Settings'.

Furthermore, the alarm audio signals can be deactivated for 30 s by tapping on the 'Bell' (alarm audio paused). The bell turns yellow during this time (Fig. 68).

At the maximum setting the alarm volume is 74 dBA and at the minimum setting 49 dBA.

**Medium priority alarm signal**

The medium priority alarm signal is a burst of 3 tones (frequency 254 Hz) each lasting 170 ms and with breaks of 230 ms. The burst is repeated every 10 seconds. Four different harmonics are generated.

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
**Low priority alarm signal**

The low priority alarm signal is a burst of 2 tones (tone 1 frequency 320 Hz, tone 2 frequency 254 Hz) each lasting 170 ms and with a break of 240 ms. The burst is repeated every 16 seconds. Four different harmonics are generated.

**8.1.4 Alarm preset**

The volume of the audio signals is set to maximum in the factory.

**8.1.5 Verification of the functioning of the alarm system**

The functioning of the alarm system is verified with the sounding of the availability signal after the mode has been selected.


The individual alarm conditions can be tested as described in section 8.3.

**8.2 Information signals**

The information signals are self-explanatory messages which explain the device behavior, support users with the individual operating functions and thus make the device easier to use, and assist the service technician during troubleshooting.

The information signals, with the exception of those which display the cause of an inoperable device ('safe state'), are not self-retaining.

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This means that the signal is only issued whilst the signal conditions are present. If the signal conditions are only very brief, the minimum duration of the display is likewise only 5 s.

The signals given are ordered according to priority. The alarms have a higher priority than the information signals. There are five priorities for the information signals.


Any low priority signal will be overwritten by an alarm signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If numerous signal conditions with the same priority occur, the most recently detected condition overwrites earlier ones in the title line of the display.

**8.2.1 Visual signals**

The information signals are also shown in the title line.

Information signals are displayed with blue font on a white background (inversely to normal display).

**8.2.2 Acoustic information signal**

The acoustic information signal is a double tone c-e (263 Hz – 330 Hz) lasting 300 ms.

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

Depending on the priority of the message, the sequence is either repeated or played once. The volume of the maximum setting is 69 dBA.

In the case of messages which display that the device is no longer operable or if the message '371' (bottle pressure < 10 bar) is coupled with a repetitive audio, the sequence is repeated every 20 s. In all other cases, the sequence is played just once.

• **Availability signal**

On completion of the successful self-test after switching on the device, an availability signal with the following characteristics sounds:

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

• **Key click**

Tapping on a key on the touch screen is indicated by a short beep. The volume of this beep can be set in the 'Device Settings' menu independently of all other information signals. The key click can also be permanently deactivated.

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### 8.3 Test conditions for alarm and information signals

#### 8.3.1 Overpressure alarm

In order to test the overpressure alarm, proceed as follows.

1. Close the instrument stopcock on the insufflation tube.
2. Start insufflation.
3. Stop insufflation.
4. When insufflation is stopped, reduce the setpoint value by 6 mmHg so that the alarm threshold is exceeded.

☛ A visual and acoustic overpressure alarm is issued.

Since automatic pressure release does not work when insufflation is deactivated, the overpressure condition continues to apply when the system is leakproof.

#### 8.3.2 CO<sub>2</sub> empty alarm

1. Before turning on the device, open the gas supply briefly and then close again.
2. Turn on the device and after a successful switch-on test, start insufflation with a flow limit of > 1.

☛ Once the input pressure has been reduced, a CO<sub>2</sub> empty alarm sounds.

#### 8.3.3 Information signals

1. Proceed as with CO<sub>2</sub> empty alarm and before starting insufflation, set the gas supply system to 'Repetitive audio signal'.
2. Stop insufflation immediately if the message '371' is displayed when the gas supply is closed.

The leakproofness of the system is sufficient to verify the characteristic of the information signal.

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**8.4 Troubleshooting**

**WARNING:** Always unplug the device  
before all maintenance work.

**Description of error:**

Complete failure of the device.

**Possible causes:**

1. Power supply failure.
2. Defective line fuse.

**Measure:**

1. Check that there is electricity to the wall outlet.
2. Change fuses as described in the instruction  
manual. Use only fuses of the correct rating.

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
Device defective.

**Measure:**

Switch unit off and back on. If the warning  
message still appears after restarting, contact  
Service. Prior to this, open the System Log under  
Settings and note the ID, Word1, Word2 and  
information text of the last 5 entries.

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No gas flow.

**Measure:**

Observe the warnings (device temperature  
exceeded, temperature of high pressure regulator  
too low, gas supply not sufficient).

If the warning message still appears after  
restarting, contact Service.

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No pressure build-up.

**Possible causes:**

1. Excessively large leak.
2. Excessively high instrument resistance.
3. Leaky insufflation system.

**Measure:**

1. Increase the flow setpoint value (for information: if the gas outflow from the cavity with the desired intraabdominal pressure is larger than the set flow setpoint value, the insufflation regulation goes to flow limitation and the cavity pressure setpoint value is not achieved).
2. Check whether it is possible to change the instrument (for information: there are situations in which the insufflation pressure provided is not sufficient to build up the desired cavity pressure due to a high level of instrument resistance and simultaneous gas outflow. Depending on the instrument and leak conditions, the current values for pressure and flow are below the setpoint values).
3.
  - a) Check the connections of the accessories and the leakproofness of the trocar. If the test for proper functioning is passed, the likelihood of a loose connection is very high.
  - b) Change the tube/tube set and perform a test for proper functioning.
  - c) If necessary, change the filter and repeat the test for proper functioning.
  - d) Contact service if the leak remains.

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
**Error description:**

After switching on the device without the heatable patient tube plugged in, the patient heating symbol lights up red (Fig. 70).

**Possible causes:**

The device is defective and the patient heating is permanently switched off.

**Action:**

Work can be continued without the patient tube heating.

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**Error description:**

After plugging in the heatable patient tube, the patient heating symbol lights up red and the message 'Defective patient tube heating' is displayed.

**Possible causes:**

The patient heating is defective.

**Action:**

Replace the heatable patient tube and put the patient heating into operation by switching the device off and on again.

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**Error description:**

After plugging in a heatable patient tube, the patient heating symbol lights up white and then changes its color from white to red, and the message 'Defective heated tube, please replace' is displayed.

**Possible causes:**

The patient heating is defective.

**Action:**

Replace the heatable patient tube. If the patient heating symbol disappears after the heatable patient tube is unplugged, a new heatable patient tube can be used immediately.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

- Type of IP protection: IP21

**According to IEC 60601-1-2, FCC-B:**

Observe the information on electromagnetic compatibility in chapter 10 (pages 71-87).

**8.5.2 Directive compliance  
(for UI 500)**

**According to the Medical Device Directive (MDD):**

This medical product belongs to Class IIb.

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

**NOTE:** The code number after the CE mark indicates the responsible notified body.

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### 8.6 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

**i NOTE:** We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Article	Cat. no.
<b>High-pressure tube</b> , American connector/ German connector, length: 55 cm	20 4000 21
<b>Same</b> , length 102 cm	20 4000 27
<b>High-pressure tube</b> , American connector/ PIN-Index connector, length: 55 cm	20 4000 22
<b>Same</b> , length 102 cm	20 4000 28
<b>High-pressure tube</b> , American connector/ ISO connector, length 102 cm	20 4002 22
<b>CO<sub>2</sub> bottle</b> , 1 l, empty, with German connector	26 4000 90
<b>CO<sub>2</sub> bottle</b> , 1 l, empty, with PIN-Index connector	26 4000 91
<b>Low-pressure tube</b> , for the central CO <sub>2</sub> gas supply, length 150 cm	UI 001
<b>Same</b> , length 300 cm	UI 002
<b>Same</b> , length 600 cm	UI 003
<b>HP-Inline gas filter</b>	20 4000 32
<b>Universal wrench</b>	20 4000 30
<b>Line fuses</b> T2.5AH 250V, pack of 10	ET15-1708090
<b>Power cord</b> (grounded)	400 A
<b>Power cord</b> 'Hospital Grade' (USA)	400 B
<b>Instruction manual</b>	96116045 D

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**9.2 Accessories**

Article	Cat. no.
<b>Insufflation tube</b> , sterilizable, length 2.5 m, ID 9 mm	UI 004
<b>Gas filter</b> , disposable, sterile, pack of 25	031122-25
<b>Insufflation tube set with gas filter</b> for single use, sterile, pack of 10	031200-10
<b>Heatable insufflation tubing set with gas filter, single use</b> , length 3 m, for gas flow up to 50 l/min, sterile, pack of 10	031210-10
<b>VERESS pneumoperitoneum needle</b> with spring-loaded blunt stylet, LUER lock, length 10 cm	26120 J
<b>Same</b> , length 12 cm	26120 JL
<b>Same</b> , length 7 cm	26120 JK
<b>Same</b> , length 15 cm for creating pneumoperitonea by means of DOUGLAS puncture with adipose patients	26120 JF
<b>Bottle holder</b> , folding, holder for KARL STORZ CO <sub>2</sub> gas bottles 1000 ml	UI 005
<b>Leakage tester</b>	13242 XL

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**10 Electromagnetic  
Compatibility (EMC)  
information**

The ENDOFLATOR® 50 UI 500 is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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**Accessories and cables used for EMC Compliance**

Cable type	Shielded	Length [m]	Ferrite	Used for
PE (potential equalization)	No	> 3	No	Potential equalization
SCB Cable	Yes	≥ 1	No	connecting SCB network, other SCB units
Power Cord	No	3	No	Power supply
031210-01	No	3	No	Heatable insufflation tubing set (patient tube heating)

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 **WARNING:** The use of an accessory or cable with the ENDOFLATOR® 50 other than those specified in this manual may result in increased emissions or decreased immunity of the ENDOFLATOR® 50. When using an accessory or cable with the ENDOFLATOR® 50 other than those specified in this manual, it becomes the responsibility of the user of the ENDOFLATOR® 50 to determine compliance with IEC 60601-1-2 when using this item.

 **WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ENDOFLATOR® 50 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

 **NOTE:** The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting equipment.

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<b>Table 1 – Compliance level for Immunity tests</b> <b>Guidance and manufacturer’s declaration – electromagnetic immunity</b>			
The ENDOFLATOR® 50 is intended for use in the electromagnetic environment specified below. The user of ENDOFLATOR® 50 should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ENDOFLATOR® 50 require continued operation during power mains interruptions, it is recommended that the ENDOFLATOR® 50 be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the ENDOFLATOR® 50 further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	6 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80 % AM modulation 6 V <sub>rms</sub> in ISM bands	

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**Table 3 – Test levels for Radiated and conducted Immunity Tests**  
**Guidance and manufacturer’s declaration – electromagnetic immunity**

The ENDOFLATOR® 50 is intended for use in the electromagnetic environment specified below. The user of ENDOFLATOR® 50 should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	6 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the ENDOFLATOR® 50, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> . $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ENDOFLATOR® 50 is used exceeds the applicable RF compliance level above, the ENDOFLATOR® 50 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ENDOFLATOR® 50.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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<b>Table 4 – Emission class and group</b> <b>Guidelines and manufacturer’s declaration – electromagnetic emissions</b>		
<p>The ENDOFLATOR® 50 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 50 should ensure that it is used in such an environment.</p>		
RF emissions as per CISPR 11	Group 1	The ENDOFLATOR® 50 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The ENDOFLATOR® 50 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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<b>Table 5</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the ENDOFLATOR® 50</b>			
The ENDOFLATOR® 50 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the ENDOFLATOR® 50 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ENDOFLATOR® 50 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies. Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

**HF surgical instrument compatibility**

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga, Ontario  
L5N 3R3 Canada  
Phone: +1 905 816-4500, Fax: +1 905 816-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8525  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
1 South Los Carneros Road  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.  
Edificio Atlantic, Oficina 3G  
Calle D e/ 1ra y 3ra  
10400 Vedado, Havana, Cuba  
Phone: +537 836 95 06, Fax: +537 836 97 76  
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.  
Lago Constanza No 326,  
Col. Granada, Del. Miguel Hidalgo  
C.P. 11520 Ciudad de México  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6° Piso, B1638 CAA – Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Stamveien1  
1483 Hagan, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
127 39 Skårholm, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Taivaltie 5  
01610 Vantaa, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)968247755  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG  
Representation Office  
Kęstucio st. 59 / Lenkiojki st. 27  
08124 Vilnius, Lithuania  
Phone: +370 5 272 0448  
Mobile: +370 685 67 000  
E-Mail: info-it-lv@karlstorz.com

KARL STORZ Endoskopi Danmark AS  
Skovlyoftens 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B.V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.  
Phone: +31 (0)33 4545890  
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S.A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH  
Landstraße Hauptstr. 148/11/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771 051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Cesta v Gorice 34b  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gliwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited  
11<sup>th</sup> Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG  
Interchange 21 Tower, Level 33  
399 Sukhumvit Road,  
North Klontoy, Wattana,  
10110 Bangkok, Thailand  
Phone: +84 28 3823 8000  
Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG  
Resident Representative Office  
14<sup>th</sup> Floor, MPlaza Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
100022, Beijing, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch  
Room 701A Building 5 & Room 501 Building 7,  
No. 3000 Longdong Avenue, Pilot Free Trade Zone,  
201203, Shanghai, People's Republic of China  
Phone: +86 21 60339888, Fax: +86 21 60339808  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch  
Room 803-805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
6100414, Chengdu, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
110014, Shenyang, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch  
Room 02B & 03 & 04A, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
510620, Guangzhou, People's Republic of China  
Phone: +86 20 87321281, Fax: +86 20 87321286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonweath Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonweath Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.  
9F Hyowon-Building  
97, Jungdae-ro, Songpa-gu  
05719 Seoul, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.  
15 Orion Road Lane Cove NSW 2066  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: info@karlstorz.au

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)