INgezim® AHSV COMPAC PLUS

Prod Ref: 14.AHS.K3

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$$PB = \frac{1,859 - 0,865}{1,859 - 0,110} \times 100 = 56,8$$

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials with Positive Control, ready to use	1	2,5mL	2	2,5mL
Vials with Negative Control, ready to use	1	2,5mL	2	2,5mL
Vials with Peroxidase Conjugate, ready to use	1	30mL	2	30mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Diluent (DE01-01), ready to use	1	125mL	2	125mL
Bottles with Stop Solution, ready to use	1	65mL	1	65mL
Bottles with Substrate (ABTS), ready to use	1	30mL	1	60mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (405 nm filter)

TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (**Blocking Elisa**). A brief description of the technique is shown below:
The antigen is fixed in a solid support (polystyrene plate). After incubation with sera samples, an AHSV specific monoclonal antibody (Mab peroxidase conjugated) is added. If the sample contains antibodies specific of the virus, they will not allow the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies, the Mab will bind to the antigen coating the plate.

After washing the plate to eliminate all non-fixed material, presence or absence of labelled Mab can be detected by adding the substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

If there has been colour development, the conjugate has bound to the antigen, being the sample negative. On the other hand, if there are antibodies specific of AHSV in the sample, they will block the binding of the conjugate and there will not be colour development.

The antigen used in this kit is VP7 recombinant protein from the AHSV (serotype 4) obtained using the baculovirus expression system. This confers many advantages to the assay. From the safety point of view, it is clear the total absence of infectivity. From the accuracy point of view, VP7 is one of the mayor proteins from AHSV, as well as the most antigenic one and the most conserved within the 9 different serotypes. It has not been found any infected or vaccinated animal without antibodies to this protein.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution can precipitate at refrigeration, in this case, incubate at 37°C before use, to eliminate the precipitation. The stop solution may form foams. Avoid dispensing bubbles on the plate as this may affect the reading of the results.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Sera samples must be diluted 1/5 in the provided diluent. This step can be made directly in the wells of the plate by adding 80 µL of diluent and 20 µL of serum sample and shaking carefully the plate to homogenize it.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the Kit in 24 parts of distilled or deionized water (40 mL of concentrate solution and 960 mL of water). Once this solution is ready, it remains stable between +2°C and +8°C.

Preparation of Controls (+) y (-):

Controls are ready to use and do not need any preparation. Dispense 100 µL of each one.

Preparation of conjugate:

Conjugate is ready to use and does not need any preparation.

Preparation of substrate:

Ready to use.

TEST PROCEDURE

1. All reagents must be brought to room temperature (**20°C to 25°C**) before use.
2. Dispense 100 µL of diluted samples into appropriate wells (dilution 1/5). We recommend running the samples in duplicate. Dispense 100 µL of positive control into two wells and 100 µL of negative control into two wells. **Seal the plate and incubate for 1 hour at 37°C.**
3. Wash 5 times as described before
4. Add 100 µL/well of conjugate. **Seal the plate and incubate for 30 minutes at 37°C.**
5. Wash 5 times as described before.
6. Dispense 100 µL/well of substrate using a multichannel pipette. **Seal the plate and incubate for 10 minutes at room temperature.**
7. Dispense 100 µL/well of stop solution. Take care to avoid dispensing of bubbles.
8. Read at 405 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

VALIDATION OF CRITERIA:

- The OD of the Positive control < 0.2
- The OD of the Negative control > 1.0

INTERPRETATION:

Determine the Blocking Percentage (BP) of each sample applying the following formula:

$$PB = \frac{\text{Abs Neg. control} - \text{Abs sample}}{\text{Abs Neg. control} - \text{Abs Pos. control}} \times 100$$

- Samples showing BP value > **50%** are considered as **positive** for antibodies to AHS virus.
- Samples showing BP value < **45%** are considered **negative** for antibodies to AHS virus.
- Samples with BP value between **45% and 50%** are considered **doubtful** and they must be retested. If the result is the same, another extraction must be made and tested 2 weeks later.

EXAMPLE:

- Average Abs Positive Control = 0.110
- Average Abs Negative Control = 1.859
- Abs sample = 0.865

$$BP = \frac{1.859 - 0.865}{1.859 - 0.110} \times 100 = 56.8$$

Sample is positive.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

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Distributed in

por

by