

APPAVISC PFS

Hydroxypropyl Methyl Cellulose Ophthalmic Solution USP

PRODUCT DESCRIPTION

Appavisc PFS is a viscoelastic solution of optimum viscosity. It contains highly purified grade of HPMC 2% w/v. It is clear isotonic sterile, non-inflammatory, used as a lubricant and protector in IOL surgery.

COMPOSITION

Each mL contains:

Hydroxy Propyl Methyl Cellulose USP ...2%w/v
Sterile Isotonic base q.s.

→ Product Features

- Clear Visualization
- Preservative free sterile solution
- Optimum Viscosity
- Better chamber maintenance
- Friendly with corneal endothelium
- Provide excellent Lubrication
- Simple to remove from the anterior chamber

→ Shelf Life

24 months


→ Packaging

Primary Packing:- 3 ml syringes, plunger stopper, plunger rod & pouch

Secondary Packing:- Carton box with 23G Sterile Cannula .

→ Storage:

Store between 15°C-25°C



INDICATIONS: Appavisc PFS is indicated as a surgical aid (medical device) in anterior segment surgical procedures involving the anterior chamber of the eye, including extraction of the lens and insertion of intraocular lenses. It maintains the depth of the anterior chamber during the whole surgical procedure and permits greater operative precision without the risk of damaging the endothelium of the cornea or other intraocular tissues.

DOSAGE AND ADMINISTRATION: APPAVISC PFS (Hydroxypropyl Methyl Cellulose Ophthalmic Solution USP) should be carefully introduced into the anterior chamber using the sterile cannula. It may be applied into the chamber prior to or following delivery of intraocular lens (IOL). It may be used to coat an intraocular lens as well as tips of surgical instruments prior to implantation which protect the corneal endothelium from possible damage during cataract surgical procedure.

SUPPLY: APPAVISC PFS 2 mL (Hydroxy Propyl Methyl Cellulose Ophthalmic Solution USP 2% w/v) is a sterile, non-pyrogenic, viscoelastic preparation supplied in 3mL syringe with luer tip and luer lock packed in double pouch and terminally sterilized.