Medtronic

TrailBlazer™ Angled

Support Catheter

Angled Support Catheter with Hydrophilic Coating
Поддържащ катетър с извит връх и хидрофилно покритие
Zahnutý podpurný katétr s hydrofilním potahem
Vinklet supportkateter med hydrofil belægning
Abgewinkelter Support-Katheter mit hydrophiler Beschichtung
Γωνιωτός καθετήρας υποστήριξης με υδρόφιλη επίστρωση
Catéter de soporte angulado con recubrimiento hidrófilo
Nurgaga hüdrofiilse pinnakattega tugikateeter
Taivutettu hydrofiilisesti pinnoitettu tukikatetri
Cathéter de soutien coudé avec revêtement hydrophile
Angled Support Catheter s hidrofílnim premazom
Szögben álló segédkatéter hidrófil bevonattal
Catetere di supporto angolato con rivestimento idrofilo
친수성 코팅 처리된 각진형 지지 카테터
Lenktas palaikomasis kateteris su hidrofiline danga
Izliektais atbalsta katetrs ar hidrófilu pārklājumu
Аголен потпорен катетер со хидрофилна обвивка
Gehoeekte ondersteuningskatheter met hydrofiele coating
Vinklet støttekateter med hydrofílt belegg
Zagięty cewnik podtrzymujący z powłoką hydrofilną
Cateter de suporte angular com revestimento hidrófilico
Cateter de suporte angular com revestimento hidrófilico
Cateter curbat de susținere, cu înveliș hidrófil
Изогнутый поддерживающий катетер с гидрофильным покрытием
Zahnutý podporný katéter s hydrofilnou vrstvou
Ukrivljeni podporni kateter s hidrófilno prevleko
Pomoćni zakrivljeni kateter sa hidrófilnim zaštitnim slojem
Vinklad stödkateter med hydrofil beläggning
Hidrófilik Kaplamalı Açılı Destek Kateteri
Зігнутий допоміжний катетер з гідрóфільним покриттям

! USA

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

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Instructions for use

Device description

The TrailBlazer™ angled support catheter is an over-the-wire (OTW), single lumen catheter with a tapered, angled tip. The catheter system is compatible with 0.36 mm (0.014 in), 0.46 mm (0.018 in), and 0.89 mm (0.035 in) guidewires. Depending on the model, the support catheter has a working length of 65 cm, 90 cm, 135 cm, or 150 cm. All models are compatible with a 4 Fr introducer sheath. For specific sizing information, see Table 1 for model specifications.


Figure 1. TrailBlazer™ angled support catheter

As shown in Figure 1, the lumen begins (2) at the proximal end of the manifold (1) and continues through the catheter (3) to the angled distal tip (4). The lumen is used to pass the catheter over a guidewire. The guidewire-compatible diameter and working length are marked on the strain relief (5). The distal catheter shaft has three radiopaque markers (6) that aid in positioning the catheter. The first 40 cm at the distal end of the catheter (7) are coated with a hydrophilic coating.

Table 1. Size specifications

MODEL	ASC-014-090	ASC-014-135	ASC-014-150	ASC-018-090	ASC-018-135	ASC-018-150	ASC-035-065	ASC-035-090	ASC-035-135	ASC-035-150
Maximum guidewire diameter mm (in)	0.36 (0.014)	0.36 (0.014)	0.36 (0.014)	0.46 (0.018)	0.46 (0.018)	0.46 (0.018)	0.89 (0.035)	0.89 (0.035)	0.89 (0.035)	0.89 (0.035)
Working length cm	90	135	150	90	135	150	65	90	135	150
Minimum guidewire length cm	150	180	180	150	180	180	150	150	180	180
Marker band spacing mm	15	15	15	15	15	15	50	50	50	50
Shaft outside diameter mm (in)	0.76 (0.030)	0.76 (0.030)	0.76 (0.030)	0.86 (0.034)	0.86 (0.034)	0.86 (0.034)	1.40 (0.055)	1.40 (0.055)	1.40 (0.055)	1.40 (0.055)
Tip outside diameter mm (in)	0.51 (0.020)	0.51 (0.020)	0.51 (0.020)	0.58 (0.023)	0.58 (0.023)	0.58 (0.023)	1.04 (0.041)	1.04 (0.041)	1.04 (0.041)	1.04 (0.041)
Minimum introducer sheath diameter mm (Fr)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)	1.40 (4)

Indications for use

The TrailBlazer angled support catheters are percutaneous, single lumen catheters designed for use in the peripheral vascular system. The TrailBlazer angled support catheters are intended to guide and support a guidewire during access of the vasculature, to allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

Contraindications

There are no known contraindications.

Warnings

- The TrailBlazer angled support catheter should be used by physicians experienced with interventional techniques in the vascular system.

- This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and compromised device performance.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. Carefully withdraw the support catheter.
- Do not torque the catheter if it is or seems stuck. Torquing with excessive force can cause damage to the device or the vessel.
- Always use fluoroscopy or other visualization when manipulating the catheter.

Precautions

- Federal (USA) Law restricts this device to sale by or on the order of a physician.
- Carefully inspect the package and catheter before use. Do not use the device if the packaging or the device is damaged.
- Do not wipe the catheter surface with dry gauze.
- Do not attempt to pass the TrailBlazer angled support catheter through a smaller-sized introducer sheath than indicated on the label. See Table 1 for size specifications.
- Ensure that the guidewire diameter and length are compatible with the TrailBlazer angled support catheter. See Table 1 for size specifications.
- To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Do not advance the catheter into a vessel with a diameter smaller than the outer diameter of the catheter.
- Do not exceed the maximum recommended infusion pressure of 2068 kPa (300 psi).
- After use, this product is a potential biohazard. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and national laws and regulations.

Potential adverse effects

Potential complications or adverse events include, but are not limited to the following conditions:

- Access site complications
- Arterial dissection
- Arterial spasm
- Arterial thrombosis
- Catheter fracture with tip separation and distal embolization
- Death
- Distal embolization (air, blood clots, or plaque)
- Intimal disruption
- Local or systemic infection including sepsis
- Perforation and vessel rupture
- Stroke/TIA
- Surgical intervention

How supplied

Warning: This device is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risks of patient infection and of compromised device performance.

Storage

Store the catheter away from sunlight at room temperature in a dry place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light.

Recommended items

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile heparinized saline
- Three-way stopcock
- Appropriately sized guidewire (refer to Table 1 for size specifications)
- Appropriately sized introducer sheath (refer to Table 1 for size specifications)

Preparation

Caution: Carefully inspect the package and device before use. Do not use the device if the packaging or the device is damaged or compromised.

- Carefully remove the catheter from the inner package.
- The TrailBlazer angled support catheter is coated with a hydrophilic coating. Before inserting the catheter, activate the coating by immersing the catheter in normal saline, for approximately 30 to 60 seconds, or by wiping the catheter shaft with a saturated gauze sponge.
 - Caution:** Do not wipe the catheter surface with dry gauze.
- To displace air, connect a 10 cc syringe filled with heparinized saline to the luer lock on the manifold and flush out the lumen.

Directions for use

- Introduce the TrailBlazer angled support catheter percutaneously using an appropriately sized introducer sheath.
 - Caution:** Do not attempt to pass the TrailBlazer angled support catheter through a smaller-sized introducer sheath than indicated on the label. Refer to Table 1 for sizing specifications.
- Place the prepared catheter over a previously inserted guidewire.
 - Caution:** Ensure that the guidewire diameter and length are compatible with the TrailBlazer angled support catheter. See Table 1 for specifications.
 - Caution:** To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

- Carefully advance the TrailBlazer angled support catheter to the desired vasculature, using fluoroscopic imaging to visualize the marker bands.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forced passage may damage the device or the vessel lumen. Carefully withdraw the catheter.

Warning: Do not torque the catheter if it is or seems stuck. Torquing with excessive force can cause damage to the device or the vessel.

Warning: Always use fluoroscopy or other visualization when manipulating the catheter.

Caution: Do not advance the catheter into a vessel with a diameter smaller than the outer diameter of the catheter.

- Carefully remove the TrailBlazer angled support catheter, while maintaining guidewire position.

Note: To use the catheter for infusion, withdraw the guidewire and attach an infusion device (syringe, power injector, and so forth). Do not exceed 2068 kPa (300 psi) inlet infusion pressure.

Caution: Do not exceed the maximum recommended infusion pressure of 2068 kPa (300 psi).

Table 2. TrailBlazer angled support catheter flow rates (ml/second)

Model	Guidewire diameter mm (in)	Length cm	Saline		Contrast*	
			1034 kPa 150 psi	2068 kPa 300 psi	1034 kPa 150 psi	2068 kPa 300 psi
ASC-014-090	0.36 (0.014)	90	0.9	1.3	0.2	0.5
ASC-014-135	0.36 (0.014)	135	0.8	1.1	0.2	0.3
ASC-014-150	0.36 (0.014)	150	0.8	1.0	0.1	0.3
ASC-018-090	0.46 (0.018)	90	1.6	2.1	0.5	1.1
ASC-018-135	0.46 (0.018)	135	1.3	1.8	0.3	0.8
ASC-018-150	0.46 (0.018)	150	1.2	1.7	0.3	0.7
ASC-035-065	0.89 (0.035)	65	7.8	11.2	6.1	8.5
ASC-035-090	0.89 (0.035)	90	7.1	10.2	5.4	8.1
ASC-035-135	0.89 (0.035)	135	5.9	8.3	4.2	7.3
ASC-035-150	0.89 (0.035)	150	5.5	7.7	3.8	7.0

*A mix of 75% Isovue 370 contrast media and 25% saline equivalent solution

Disposal

Caution: After use, this product is a potential biohazard. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and national laws and regulations.

Warranty Disclaimer

This product is subject to standard Medtronic warranty terms.

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- © 2017 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
ev3 Inc.
Medtronic Australasia Pty Ltd