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SLIDES FOR URINE SEDIMENTATION

DISPOSABLE SLIDE IN PMMA FOR THE QUICK AND ACCURATE READING OF 10 SIMULTANEOUSLY CELLULAR COUNTS. EACH CELL IS PROVIDED WITH 5 COUNTING GRIDS. LATEX FREE DEVICE

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CE Marked product - manufactured in compliance with 98/79/CE Directive and D.lgs 332 dtd 08/09/2000

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Following to related calculations, by using these slides you may determine the number of figured elements present in a urine sample in ml.

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Surface on which the sample distributes on the slide is divided in spaces within a grid. Inside these spaces (squares), volume of the sample is predetermined and, therefore, by following given procedures you may have a direct rate between figured elements calculated inside the spaces and their number present into the sample of tested urine.

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The causes of error inherent to the analytical method for the counting of the figured elements in urinary sediment through the above mentioned chambers can be summarized as follows.

Error dependent on:

MATERIAL

- a) pipettes not properly calibrated
- b) defective counter chambers
- c) inappropriate or incorrectly positioned coverslips

TECHNIQUE

- d) defective sample collection
- e) inaccurate pipetting
- f) defective mixing
- g) drying phenomena

OPERATOR

- h) difference between one operator and another
- i) visual tiredness of the operator

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Plastic slides for urine sedimentation code 15101 are different if compared with the traditional glass chambers for their being very easy to use. This, beyond helping the laboratory personnel in their daily routine, highly reduces the possibilities of mistakes in using the mentioned chambers. In addition, it is a best product compared to other plastic chambers as the grid marking is clearer, allowing an homogeneous distribution within the grid of figured elements of the urine sample. As a fact, in this slide the volume of sample within the grid is pre-defined and constant in all cells.

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The grid is divided in 5 squares with internal side of 1 mm. Each of these squares having a side of 1 mm, is divided in other 9 small cells having a side of 0,333 mm. Each grid is therefore divided in 45 small squares with a side of 0.333 mm.

This way you can obtain a precise division of the volume of sample within the grid, as follows:

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Slide for urine sedimentation code 15101 presents the following advantages:

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READING AND COUNTING THE FIGURED ELEMENTS IN THE GRID

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Independently from the counting directions given, there could be a percentage of errors due to the random distribution of the elements in the cells.

The error percentage could be reduced by counting a bigger number of elements in the squares.

It is therefore necessary that the operator establish the way of counting depending on the number of figured elements present in the sample under examination.

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to obtain values expressed in number of elements per } \mu\text{l of urine}

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to obtain values expressed in number of elements per ml of urine

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$k = 0,01111$

N = numero cellette osservate / *number of small cells observed*

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Count the number of elements present in 5 different squares, taking care not to count twice the same position (example n. 1).

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*The total number of figured elements for ml of urine can be obtained using the **table 1** following reported.*

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Count the number of elements present in 10 different squares, taking care not to count twice the same position (example no. 2)

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The total number of figured elements for **ml** of urine can be obtained using the **table 2** following reported.

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Codice EDMA > 29011001 - Other Hardware + accessories + consumables + software

Intended purpose is “In Vitro Medical Diagnostic Device”, to content biological fluids for lab testing.

For professional use only in laboratory of analysis. Valid nomenclatures must be used. Further tests must be selected and implemented according to recognized methods.

National classification of medical devices (CND - For Italian law) > W0101060203 (Urine Sediment test Kit)

EDMA code > 29011001 - Other Hardware + accessories + consumables + software

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Keep out of flame or heat sources which might damage the products

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Do not use after expiry date or if packing is opened

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Do not re-use: disposable device

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Do not vary the intended purpose of the product

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Disposal: use appropriate personal protective equipment and act according to applicable regulations

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Before using with particular substances check the resistance/compatibility chart on our catalogue.

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