# REF 01FK10W

**Abbott**

**BioLine**

**HBsAg WB**

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### ENGLISH

**About the test**  
 The presence of hepatitis B is caused by the infectious Hepatitis B virus (HBV). HBV DNA is estimated to have infected 400 million people throughout the globe, making HBV one of the most common human pathogens. It is mainly transmitted between people through contact with infected blood, frequently from mother to child at birth. Serum specimens, also a possible means of transmission. There are several types of transmission. A variety of serological assays may be employed to differentiate the type of viral infection as well as to discriminate between chronic and acute hepatitis B virus (HBV) infection. The sensitive and specific methods used commercially in diagnosis are the immunochromatographic test, radioimmunoassay (RIA) and enzyme-linked immunosorbent assay (ELISA).

**[Principle]** The membrane is pre-coated with monoclonal anti-HBsAg pool on the test line region and Mouse monoclonal anti-hepatitis B core (anti-HBc) on the control line region. During testing, the specimen is allowed to react with the colorimetric antigen-antibody complex well formed in advance depending upon pre-coated on the test strip. The mixture (mouse monoclonal anti-HBsAg + HBsAg in specimen) then moves upward on the membrane chromatography by capillary action. For a reactive result, a purple-colored line appears above the test line region. If the specimen is negative, no purple-colored line appears. An absence of this purple-colored line in the test line region signifies a non-reactive result. Regardless of the presence of HBsAg, a colorimetric IgG conjugated gold colloid, pre-coated on the test strip, will migrate upwards to modified mouse monoclonal anti-hepatitis B core (anti-HBc) on the control line region of the control line region if the result window appears. The presence of a purple-colored line in the control line region signifies a reactive result. The volume of specimen has been added and will that proper form has been obtained.

**[Intended use]** HBsAg WBs are an in situ immunochromatographic, rapid assay designed for the qualitative detection of Hepatitis B surface antigen in human serum, plasma (heparin, EDTA and sodium citrate) or urine and whole (heparin, EDTA and sodium citrate). BioLine™ HBsAg WB is intended only for professional use as an aid to diagnosis. Reactive specimens should be reflexed for additional testing either by Enzyme immunoassay (EIA) or by radioimmunoassay (RIA) in a laboratory setting. The specimen is prepared to be used in a laboratory with high HIV prevalence. This test may not be suitable for diagnosis of early infection or blood donation screening. Because false non-reactive results may also arise due to the presence of interfering substances, the performance of BioLine™ HBsAg WB in infants or children has not been validated.

**Materials provided and active ingredients of main components**

- The BioLine™ HBsAg WB contains the following items to perform the assay:
  - 1. HBsAg positive serum/plasma samples without genotyping information used to test in Sanquin.
  - 2. Active ingredients of main components:
    - Microcapsules: Mouse monoclonal anti-HBsAg conjugated gold colloid (0.054:0.101 µg), Chitinase IgG gold colloid (0.042:0.108 µg), Control IgG gold colloid (0.042:0.128 µg), Check line: Mouse monoclonal anti-hepatitis B core (anti-HBc) (0.482:0.096 µg)

**Materials required but not provided**

- Microcapsules, Protective gloves, Tissue, Biohazard container

**Kit storage and stability**

The kit should be stored at a temperature between 1°C and 40°C. Do not freeze the kit or its components.

The test device is sensitive to both humidity and humidity. Perform the test immediately after removing the test device from the package. Do not use the kit if the package is damaged or the kit is indicated on the test kit.

Do not use the test kit beyond its expiration date. The shelf life of the kit is as indicated on the test kit.

Do not use the test kit if the pouch is damaged or the kit is broken.

**Warnings**

- Test devices are for in vitro diagnostic use only. Do not reuse the test device.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be proficient.
- Reactive results may be false. If a reactive result is obtained, it is important that all components are being handled.
- Use gloves and protective gloves while handling specimens and wash hands thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all reaction, reagent kits and potentially contaminated materials in a biohazard container or other appropriate waste.
- Do not mix or interchange different specimens.
- Do not use the desiccant in the foil pouch.
- Avoid splashing or spraying reagents or aerosol formation.

**Specimen collection and handling**

- Using venipuncture, draw whole blood and insert it into the collection tube (containing anticoagulant including heparin, EDTA and sodium citrate).
- If the specimen is not immediately available, it must be refrigerated at 2 - 8 °C. If stored at 2 - 8 °C, the blood specimen must be tested within 3 days of refrigeration.
- Do not use a blood specimen stored for more than 3 days; otherwise a non-specific reaction may occur.

**Plasma or serum**

[Plasma] Using venipuncture, draw whole blood and insert it into the collection tube (containing anticoagulant including heparin, EDTA and sodium citrate) and centrifuge the blood to generate a plasma sample.

[Serum] Using venipuncture, draw whole blood and insert it into the collection tube (NOT containing anticoagulant) then leave for 30 minutes to allow blood to clot completely.

Centrifuge the tubes to generate the following plasma or serum samples:

- WHO laboratory reference panels
- WHO HBsAg laboratory reference panels
- WHO HBsAg laboratory reference panels
- WHO HBsAg laboratory reference panels

**Reference assay**

Reference assay	Number of specimens	Reactive	Indeterminate	Non-reactive
EIA	Positive	201	201	0
Negative	313	1	310	0
<b>Sensitivity (95% CI)</b>		100% (98.1 - 100%)		
<b>Specificity (95% CI)</b>		99.0% (97.2 - 99.8%)		

**Test procedure (Refer to figure)**

- Bring the test device to room temperature between 15 - 40 °C prior to testing.
- Remove the test device from the foil pouch, and place it in a flat, dry surface. Label the test device with a patient identifier.
- Using a microcapillary dispenser, 100µl of serum, plasma or whole blood specimen into the specimen well.
- As the test begins to run, you will see a purple color move across the result window in the center of the test device.
- Interpret test results at 20 minutes.

**Caution:** Do not read test results after 20 minutes; late readings can yield false results.

**Test interpretation (Refer to figure)**

A colored control line will appear at "C" in the result window to show that the test is working properly. The "T" reaction of the result window indicates the test result.

**Non-reactive:** The presence of only the control line (C) in the result window indicates a non-reactive result.

**Reactive:** A non-reactive result does not preclude the possibility of infection with HBV. Please refer to the test limitations section in the IFU. The control line (C) in the result window, regardless of which line appears first, indicates a reactive result.

**Invalid result:** If the control line (C) is not visible within the result window, the result is considered reactive. If the control line (C) is visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test kit may have deteriorated. It is recommended that the specimen be retested using a new test device.

**Test limitations**

- This product may have performance limitations for diagnosis in clinical specimens that are expected to be in this phase of infection or low concentration of HBsAg in origin. Although this assay has been demonstrated to be able to detect HBsAg genotype variations, false results may occur. These performance limitations can be a factor of false non-reactive results. A non-reactive result does not preclude the possibility of infection with HBV.
- Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician. After clinical and laboratory findings have been evaluated.
- Due to the inherent design of qualitative IVD tests, a faint or absent test line (false non-reactive) may occur in specimens containing high concentrations of HBsAg or antigen. In order to obtain a definitive result, all clinical and laboratory findings should be evaluated in reference panel acceptance criteria.

**Internal quality control**

The BioLine™ HBsAg WB device has "Test Line" and "Control Line" on the surface of the device. The entire test line and control line in the result window are not visible before applying any specimens. The control line is used for procedural control. The control line of the RDT only shows that the active ingredients of the main components are in the test kit and are functional, but does not guarantee that the specimen has been properly applied; it is not a reactive component.

**Performance characteristics**

The BioLine™ HBsAg WB test kit is designed to have 100% (95% CI: 98.3 - 100%) of sensitivity

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**Test procedure (Refer to figure)**

- Bring the test device to room temperature between 15 - 40 °C prior to testing.
- Remove the test device from the foil pouch, and place it in a flat, dry surface. Label the test device with a patient identifier.
- Using a microcapillary dispenser, 100µl of serum, plasma or whole blood specimen into the specimen well.
- As the test begins to run, you will see a purple color move across the result window in the center of the test device.
- Interpret test results at 20 minutes; late readings can yield false results.

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A membrana é pré-revestida com pool de anti-HBsAg monoclonal de rato na região da linha de teste e IgG monoclonal de rato anti-HBc na região da linha de controle. Durante o teste, a amostra de reação (anti-HBsAg + HBsAg in specimen) then moves upward on the membrane chromatography by capillary action. For a reactive result, a purple-colored line appears above the test line region. If the specimen is negative, no purple-colored line appears. An absence of this purple-colored line in the test line region signifies a non-reactive result. Regardless of the presence of HBsAg, a colorimetric IgG conjugated gold colloid, pre-coated on the test strip, will migrate upwards to modified mouse monoclonal anti-hepatitis B core (anti-HBc) on the control line region of the control line region if the result window appears. The presence of a purple-colored line in the control line region signifies a reactive result. The volume of specimen has been added and will that proper form has been obtained.

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
- 2** **EN** Interpret test results at 20 minutes. Do not read test results after 20 minutes; late readings can yield false results.

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
- EN** The presence of only the control line (C) within the result window indicates a non-reactive result.

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
### REACTIVE / RÉACTIF / REACTIVO / REATIVO

- EN** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a reactive result. **⚠ Caution:** The presence of any test line, no matter how faint, the result is considered reactive.

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### INVALID / NON VALIDE / NO VÁLIDO / INVÁLIDA

- EN** If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test kit may have deteriorated. It is recommended that the specimen be retested using a new test device.

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