CHEMICAL FUMEHOOD

NIȘA CHIMICĂ

Metallic freestanding fume cupboard with exhaust and filtration system and integrated under bench storage for chemical reagents according to EN 14175.

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FUMEHOOD • FILTERS • EXHAUST SYSTEMS
REAGENT STORAGE • WORK SAFETY • AIRFLOW

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The chemical fumehood – a protection equipment for laboratory use. It's main roles are:

- **To control** user's exposure to gases, fumes, vapors aerosols and other dangerous airborne particles while operating laboratory procedures involving potentially harmful chemicals.
- **To prevent** the development of explosive clouds within the working area.
- **To protect** users in case of violent chemical reactions, explosions or splashes with dangerous substances

Operating principle: In order to ensure a contained hazardous perimeter within it's casing, negative pressure relative to surrounding environment is applied by connecting the fumehood to a complex exhaust system.

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The fumehoods casing: Is made of electrostatic powder-coated steel. For ease of transportation and installation, the fumehood's body is made of easy to handle parts.

Work area: Is made of a contained protected space, with sealed walls and joints on all sides, accessible only from the front side, through a sliding vertical window - the sash. Opening the sash allows access for placing materials and equipment in the work area, so that the user can safely operate them. The sash opening can be adjusted at any height and has an automatic safety locking mechanism.

Control panel: It is ergonomically placed under the worktop and comprises the controllers for all utilities: electric sockets, lights, ventilation, etc.

Reagent storage body: Located beneath the workspace - for frequently used reagents and consumables, can be connected to the fumehoods main exhaust system.

Filters and ventilation: The exhaust system is powered by an electric single-phase /three-phase motor and centrifugal fan. In order to insure a constant airflow at the opening of the sash, the fumehood system has a controller unit that adjusts the fan speed to the opening of the sash. Filters: The filter system is made of active carbon plates.

INTENDED USE: Designed according to EN 14175, the fumehood is a safety laboratory equipment insuring safe work environment for operating laboratory staff.

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1. Worktop: According to the intended laboratory activity – acid-proof, with drain sink and /or profiled overlapped anti-drip edges.

2. Storage units: Cabinets placed under the operating panel – for everyday reagents and consumables – can be connected to the main fumehood exhaust system.

3. Operating panel: Placed ergonomically right under the work area, holds together all switches and controllers for the fumehood supplies: light, ventilation, water /gas taps, etc.

a. Exhaust system: Controller with digital screen, assures a constant 0,5 m/s sash airflow, regardless if the sash is raised or lowered, according to **EN 14 175**. Equipped with acoustic sensor – audible warning system in case of low pressure /error.

b. AC plugs: IP44 with protection cover.

c. Faucets: Designed for laboratory use, with remote control valves, mounted on the control panel, and taps mounted on the deflector panel. Color code according to **EN 13 792:2002**.

3. Deflector panel: Safely guides the airflow.

4. Safety sash /work area: Work area is contained by acid-proof casing with sealed joints, front side of the fumehood with vertical sliding sash.

5. Body /case: The fumehoods casing is made of galvanized electrostatic powder-coated steel.

6. Filters: Are made of active carbon plates set up in special casings, easy to replace.

7. Ventilation and ducts: 230 /380 V with acid-proof or antiex ducts.

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