



## EW 1 Automated endoscope reprocessor (AER)



# Steelco EW 1 has been developed to be compliant with the standard UNI EN ISO 15883-4 and UNI EN ISO 14937.

EW 1 is an automatic device for the reprocessing of one (1) flexible endoscope or up to two (2) videobroncho/cystoscopes or up to four (4) single-channel fiber broncho/cystoscopes of the major brands.

It also allows the reprocessing of rigid endoscopes with dedicated washing carts.

EW 1 has been validated for low temperature high level disinfection and liquid sterilization processes and it can perform the thermal self-disinfection, as recommended by the standard UNI EN ISO 15883-4.

Highest flexibility is granted by the compact design, which assures low cycle costs reducing consumption of water, chemicals and energy.

#### Models

The AER is available as single door or double door pass through versions.

EW 1 is also available in a *Rack* configuration having two units installed one on top of the other in a unique frame structure, in order to exploit small spaces and enable the asynchronous treatment of two endoscopes. The *Rack* configuration is available both with single and double door.

#### Specifications

#### Single unit dimensions

External dimensions Single door: 600 mm x 630 mm x 945 mm Double door: 710 mm x 570 mm x 1206 mm

Washing chamber 415 mm x 480 mm x 375 mm

Loading height Single door: 683 mm Double door: 950 mm

Door passage 415 mm x 260 mm

#### **Rack configuration dimensions**

External dimensions Single door: 750 mm x 690 mm x 1643 mm Double door: 860 mm x 570 mm x 1663 mm Washing chamber 415 mm x 480 mm x 375 mm Loading height Single door – lower chamber – 712 mm Single door – upper chamber – 1322 mm Double door – lower chamber – 720 mm Double door – upper chamber – 1350 mm Door passage 415 mm x 260 mm





Water consumption

6 liters per chamber fill

Sound level

< 60 dB

#### Certifications

EW 1 is classified class IIb, according to rule 15 of Annex IX, 93/42/EEC Medical Device Directive.

EW 1 is compliant with the standard EN ISO 15883 parts 1, 4 and 5.

EW 1 is compliant with the standard EN ISO 14937.

#### Construction

#### Exterior

- Main frame made of stainless steel AISI 304 (DIN 1.4301).
- Panels made of AISI 304 MS2 finish.

#### **Hinged Door**

- Stainless steel frame hinged door.
- The glass door grants total visual monitoring of the process and it is realized with a double HST (High Shock Tested) tempered glass.
- Interlocked doors in case of double door version.

#### Washing Chamber

- Constructed using AISI 316L BA Ra < 30μin (Ra < 0.8μm).</li>
- Designed and constructed with smooth edges and corners to allow the total chamber draining at the end of the cycle and to avoid areas where dirt can accumulate and allow bacterial growth.

#### Components

 Manufactured using stainless steel and other materials which are resistant against the effects of aggressive detergents.

#### Insulation

• High performance melamine insulation guards for thermal and acoustic insulation.

#### Standard features

#### Loading/unloading doors

- Door locking system during cycle execution to guarantee operator safety and to prevent any interference with the washing cycle.
- The double door model is equipped with an automatic interlock system to keep the dirty/ clean areas separate.
- Open doors also work as a stand for the washing

cart during loading and/or unloading operations and as an anti-dripping surface, thus helping keep the floor dry, clean and safe.



#### Endoscope loading and connection to the AER

Steelco EW 1 is equipped with a removable washing cart made of stainless steel, with guiding rails to ease its pull-out and allow the ergonomic load of one big size flexible endoscope. Such washing cart can be similarly pulled out from the opposite side in double door versions at the end of the cycle.

#### Washing system

- Washing circuits with two (2) dedicated pumps:
  - For the washing arms placed on the top and on the bottom of the washing chamber. They are set to grant an efficient dynamic washing flow which allows a thorough endoscope cleaning.
  - For the basket connection and endoscope channels.
- Spray arms, made of AISI 316L stainless steel, which can be easily disassembled for cleaning and maintenance operations.
- The drain system grants the complete emptying of the washing chamber and the hydraulic circuit.

#### Water connections and filtration

- Two (2) water line connections available for softened cold/mixed and demineralized water.
- Water line is equipped with two (2) flowmeters for the control and validation of water quantity introduced in the chamber.
- Water filtration system through a sequence of two filters (one 0.45 µm and one 0.1 µm filter).
- The dual stage water filtration system not only ensures the presence of a microbiologically adequate final rinse water, but it also filters and ensures the high quality of all the water used in the disinfection phase.
- Filters are automatically decontaminated during the thermal self-disinfection.





#### Channels treatment and monitoring

- The channels treatment is guaranteed by a dedicated pump, a flow sensor and a pressure transducer for the water circulation control.
- Endoscope channels flow and pressure are monitored during the whole cycle.
- Automatic stop of the cycle in case of alarm and automatic rinse (single or double) to remove chemical product residuals, when needed.

#### Channels purging/drying system

- Endoscope channels are purged through 0.2 µm filtered compressed air (built-in air compressor).
- Alternatively, the machine can also be directly connected to external medical grade quality compressed air (air compressor not available).
- An additional setup for the connection to external medical grade quality compressed air is available as an option. In this case, built-in air compressor and connection for external compressed air are both available, but the default functioning will be with built-in air compressor.
- A forced heated air-drying system is available as an option for the elimination of residual humidity inside the washing chamber and on the external surfaces of the endoscope.

#### Leakage test

- The leak test verifies the instruments suitability to be treated into the AER preventing possible problems before they can cause damages.
- Leak test check is active during the whole cycle with automatic cycle stop in case of anomalies.

#### Washing chamber heating elements

- 2.5 kW electrical heating element in the sump of the washing chamber
- Two (2) independent temperature probes (PT 1000) in the washing chamber
- Additional PT1000 probe to check the inlet water temperature and the temperature during the heating phase.

#### Chemical dosing system

#### The correct dosing quantity is essential for trustworthy disinfection results.

- Three (3) automatic membrane dosing pumps precisely add liquid chemical agents.
- Volumetric control of the dosed chemicals by double, high-accuracy flowmeters to control and validate the disinfectant and detergent quantity.

#### **Chemical storage**

## The endoscope washer EW 1 allows the storage of chemical tanks inside its frame.

• Capacity of up to four (4) tanks equipped with

protective caps to avoid any vapors emission.

- Collection tank in case of chemical leakages.
- The feeding system is equipped with low level sensors.

#### Drain pump

#### Control System

The control system allows the definition and the memorization of customized programs dedicated to the treatment of each single instrument according to the chemical products used.

- Three (3) electrical boards: one (1) master and two (2) slaves.
- Standard pre-memorized programs certified by microbiological hygienic reports carried out by certified laboratories.
- All different programs can be selected from the control panel. Cycle start is possible only after the operator recognition.

#### System control panel

- Soft touch control system on glass panel.
- A 3,5" graphic color LCD display constantly informs the operator about the machine status: cycle phase, time left, chamber temperature, chemical and water quantities. The system registers all the sensible parameters to print and archive them.
- At the end of the cycle, the system highlights the correct instruments reprocessing.
- Acoustic and visual alarms provide quality control for each wash cycle.
- Printer connection to monitor and validate the washing cycles.
- USB port on frontal panel for cycle data, machine parameters and washing programs download. It also allows easy firmware update.
- Optional Ethernet RJ 45 connection via gateway for supervisor SW connection.

#### **Process traceability**

The identification of each single endoscope ensures that the cycle control information is closely related to the single instrument, thus granting its complete traceability.

#### EW 1 prints a final report for each cycle with:

- AER serial number identification
- ID of the instruments
- ID of the operator

It also reports all the process parameters as:

Cycle time





- Water consumption
- Chemical products dosing
- Phase details
- Events

At the end, the cycle status is pointed out (pass/fail).

#### **Built-in printer**

Built-in printer for validating washing phases with detailed information

#### **Reprocessing cycles**

# EW 1 is compatible and has been validated for low temperature high level disinfection (35°C) and liquid sterilization processes.

By using Steelco chemical products, EW 1 grants instruments maximum safety conditions as well as **efficient results proven by microbiological tests** together with complete documentation.

#### Validated high level disinfectant

The validation has been performed by a certified European microbiological laboratory to attest the conformity to UNI EN ISO 15883-4.

#### Validated liquid sterilant

The validation has been performed by a certified European microbiological laboratory to attest the conformity to UNI EN ISO 14937.

#### Self-disinfection cycles

- The device performs the thermal self-disinfection cycle: 80°C/176°F washing phase as requested by the standard UNI EN ISO 15883-4.
- Completely automatic and programmable selfdisinfection cycles.
- Self-disinfection of the hydraulic circuit and of the water filtration system.

#### Optional

#### Washing chamber forced warm air-drying system

- A forced HEPA H14 filtered air drying system ensures the drying of the external surface of the endoscopes.
- 0.75 kW electric heating element.

#### Extra heating elements

4 kW heating element to be used where there is only cold water available (to reduce the total cycle time without the installation of boilers).

#### Barcode reader

For the recognition of operators and instruments.

#### Stand on wheels (only for single door machine)

#### SteelcoData Live software

Dedicated software for the traceability of all the information of the endoscope washing cycles.

#### SteelcoData ARES software

Dedicated software for the monitoring of the whole endoscope reprocessing system including the use, the manual prewashing, the automatic reprocessing, the storage, the service and all the waiting times between any of these phases.

#### Light inside the chamber

And others, please verify with your local distributor/ agent.

#### Accessories

- Complete range of connectors for endoscopes brands and models.
- Duodenoscope dedicated reprocessing support.
- Steelco Process Challenge Device (PCD) kit for the control of the cleaning efficacy.
- Trolleys for the transport and short time storage of contaminated and/or disinfected endoscopes and accessories.

#### **Required utilities**

For any connection detail, please refer to installation drawing of the selected model/version.

#### Demineralized and softened cold/mixed water

#### **Drain connection**

#### **Electrical requirements**

- The total power of the machine in the standard configuration is 3,05 kW.
- 230V 50Hz
- Other electrical connections are available to match electrical requirements of installation site.

**Rev.13** 



**CE** 0051



Group

Member

Míele



Via Balegante, 27 31039 Riese Pio X (TV) ITALY

Manufacturer:

### STEELCO S.p.A.

Via Balegante, 27 31039 Riese Pio X (TV) ITALY



### CONTENTS

1.	GE	NERAL RULES	6
	1.1	LIMITS OF MANUFACTURER'S LIABILITY	6
	1.2	MANUAL VALIDITY, CONTENTS AND CONSERVATION	6
	1.3	REGULATIONS	7
2.	SAI		8
	2.1	ACHINE COMPONENTS DESCRIPTION	8
	22		8
	2.2	1 APPLICATION FIELDS	8
	2.2.	2 VALIDATED CYCLES	9
:	2.3	IMPORTANT WARNINGS AND SUGGESTIONS	10
	2.4	SAFETY RECOMMENDATIONS	11
	2.4.	1 INLET WATER	11
	2.4.	2 Rinse water	12
:	2.5	RECOMMENDATIONS TO ENSURE HIGH QUALITY PERFORMANCE	13
:	2.6	Residual risks	14
:	2.7	SAFETY SIGNALS USED.	15
:	2.8	TRAINING	16
	2.8.	1 STAFF QUALIFICATION	16
:	2.9	INDICATION OF SOUND LEVEL	17
:	2.10	TRANSPORT AND STORAGE	17
:	2.11	CHEMICAL PRODUCTS CONNECTIONS	18
	2.11	.1 PRESENCE SENSOR OF CHEMICAL PRODUCT	19
	2.11	.2 METER QUANTITY OF CHEMICAL PRODUCT	19
	2.11	.3 CONTAINER BLOCK	19
	Eac	H CONTAINER MUST BE LOCKED BY THE SPECIAL BRACKET TO ENSURE THE CORRECT DATA READING	19
	2.11	.34Replacement of chemical product container	19
:	2.12	AMBIENT VENTILATION REQUIREMENTS	20
2	2.13 T/	ABLE OF SYMBOLS	21
3.	USI	NG THE MACHINE	22
;	3.1	CHECKS	22
;	3.2	OPENING AND CLOSING THE DOOR	22
	3.2.	1 DOOR UNLOCKING	22
;	3.3	SWITCHING ON	23
;	3.4	PREPARATION OF THE ENDOSCOPES	24
	3.4.	1 CHECK OF ENDOSCOPE INSTRUMENT CONNECTORS	24
;	3.5	THERMAL-DISINFECTION CYCLE IN CONFORMITY WITH EN ISO 15883:4	25
;	3.6	LOAD PREPARATION	27
:	3.7	START THE WASHING CYCLE	33
	3.7.	1 ENDOSCOPE REQUIREMENTS	36
;	3.8	INSTRUCTION IN CASE OF 24H STOP OF THE MACHINE	36
4.	CO	NTROL PANEL AND SYMBOLS USED	. 37
4	4.1	SWITCHES (LOADING SIDE)	37
4	4.2	CONTROL PANEL (LOADING SIDE)	38
4	4.3	SWITCHES (UNLOADING SIDE – IF PRESENT)	43
4	4.4	CONTROL PANEL (UNLOADING SIDE – IF PRESENT)	44
5.	MA	CHINE CONDITIONS	45
!	5.1	STAND-BY	45
1	5.2	ON CYCLE	45
	J.Z		
!	5.3	ALARM	45
6.	5.2 5.3 SPF	ALARM	45 <b>46</b>
6.	5.3 SPE	ALARM	45 . <b>. 46</b> 46



	6.1.1	1 Power failure during "WAIT" status	
	6.1.2	2 Power failure during a cycle or during an Alarm	46
6.	.2	RESET PROCEDURE	
7.	wo	RK PROCEDURES	47
7.	.1	INTRODUCTION	
7.	.2	INSTRUCTIONS TO PERSONNEL	47
7.	.3	STERILIZATION - WARNINGS (CYCLE VALIDATED IN ACCORDANCE WITH UNI EN ISO 14937 STANDARDS)	
8.	CLC	ЭСК	49
9.	HIS	TORICAL DATA AND PRINTOUT	49
10.	U	SB PORT	51
1(	0.1	DATA SAVING	
1(	0.2	DATA SAVING DURING THE CYCLE	
1(	0.3	OPERATOR ARCHIVE MANAGEMENT	
1(	0.4	INSTRUMENT ARCHIVE MANAGEMENT	56
11.	Α	LARMS AND EVENTS LIST	58
1	1.1	LOGICAL DESCRIPTION OF ALARM INTERVENTIONS	
1	1.2	LIST OF ALARM MESSAGES	58
1	1.3	LIST OF WARNINGS MESSAGES	64
1	1.4	LIST OF HISTORICAL EVENTS	64
12.	IN	ISTALLATION AND OPERATION	65
13.	RAC	K VERSION	65
14.	М	AINTENANCE	66
14	4.1	GENERAL RECOMMENDATIONS ON MAINTENANCE	
	14.1	.1 MAINTENANCE REQUEST	66
14	4.2	TABLE OF ROUTINE MAINTENANCE	67
15.	Р	ROBLEMS - CAUSES – SOLUTIONS	83
1	5.1	INTRODUCTION	
1	5.2	PROBLEMS - CAUSES - SOLUTIONS	83
16.	D	ECOMMISSIONING	85
16	6.1	INSTRUCTIONS FOR DISASSEMBLY OF THE MACHINE	
10	6.2	MACHINE DISPOSAL	85
AN	NEX	A – TRAINING CERTIFICATE	
AN	NEX	B – WATER SAMPLING DURING THE FINAL RINSING CYCLE	87
AN	NEX	C – TEST PROCEDURE WHEN INSTALLING	92



Thank you for purchasing this appliance.

The installation, maintenance and operating instructions given in the following pages have been prepared to ensure the long life and good performance of the appliance.

Following the instructions carefully.

The appliance was designed and constructed using the latest technological innovations available. Please take good care of it.

Your satisfaction is our best reward.

#### WARNING:

NON OBSERVANCE, EVEN IN PART, OF THE RULES INDICATED IN THIS MANUAL WILL CAUSE THE PRODUCT GUARANTEE TO BECOME INVALID AND RELIEVES THE MANUFACTURER OF ANY RESPONSIBILITY.

THE MANUAL CONTAINS ALL THE OPTIONALS THE MACHINE CAN BE EQUIPPED WITH.



## 1. GENERAL RULES

### 1.1 Limits of manufacturer's liabili

The manufacturer shall not be held liable for failures or problems which arise due to tampering and/or incorrect applications and/or improper use of the machine.

The purchaser must comply with all instructions set forth in the user's manual, and he must in particular:

- Always work within the allowable limits for the use of the machine;
- Always carry out constant and diligent maintenance;
- Allow use of the machine by persons with proper skills and abilities for their role and purpose who have been
  properly trained and instructed;
- Use only manufacturer original spare parts.

Any modifications, adaptation or the like which may be made to machines which are subsequently placed on the market do not oblige the manufacturer to intervene on previously supplied machines, nor to consider the machine and the related user's manual lacking and inadequate.

The installation, maintenance and operating instructions given in the following pages have been prepared to ensure the long life and outstanding performance of the appliance.

For some especially demanding programming or maintenance operations, this manual serves as a memorandum of the main operations to be carried out.

Education on these topics can be obtained by attending training course held by the manufacturer.

The instructions in this manual do not replace but rather are in addition to employer requirements to adhere to current legislation on standards of prevention and safety.

The machine is guaranteed for 15 months as from the time of shipment.

### **1.2 Manual validity, contents and conservation**

This manual reflects the state of the art at the moment of manufacture and delivery of the appliance and is valid for its entire life cycle.

The manufacturer is at clients' disposal for further information or to receive suggestions for making the manual more compliant with the needs for which it was prepared.

The translation of the contents into the client's language has been carefully prepared.

In order to prevent possible accidents to persons or property due to in correct translation of the instructions, the client must:

- Not perform operations or manoeuvres with the machine if there are any doubts or uncertainties about the operation to be performed;
- Ask technical service for clarification of the instruction.
- If lost, ask for a new copy from the manufacturer.

It is important to keep this instruction manual with the machine for future reference.

If the machine is sold or transferred, the manual must be handed over to the new owners or user in order for them to become acquainted with its functioning and the relative warnings.

#### Read the warnings carefully before installing and using the machine.

This is a translation of the Italian text, which prevails in case of doubts.



### 1.3 Regulations

The purpose of the warnings is to safeguard the user in compliance with following Regulations and "Technical Product Standards":

#### EUROPE:

- 2006/42/CE (Machine directive);
- 93/42/EEC and s.m.i (Medical Devices Directive);
- 2014/35/EU (Low Voltage Directive);
- 2014/30/EU (EMC Electromagnetic compatibility directive);
- 2014/68/EU (PED Directive);
- EN 61010-1 (Safety);
- EN 61010-2-040 (Safety);
- 2011/65/EC (RoHS II);
- 2012/19/EC (WEEE);

#### and recognized international standards:

- IEC 61000 (Electromagnetic compatibility);
- IEC 61326-1 (Electromagnetic compatibility);
- ISO 14971 (Medical devices risk analysis);
- ISO 15883-1 (General requirements, terms, definitions and tests);
- ISO 15883-4 (Requirements and tests for medical washer and disinfectors that carry out the chemical disinfection of flexible endoscopes);
- ISO/TS 15883-5 (Soil test A soil test method to prove the effectiveness of cleaning activities);
- IEC 60529 (IP Grade);
- UNI EN ISO 14937 (Sterilization of healthcare products).

Steelco declares that this product, when it is equipped with a water steam version, is in accordance with PED 2014/68/UE directive art. 4 par. 3 and it has been designed and built in accordance with the correct building procedure.



## 2. SAFETY INFORMATION

Compliance with safety standards allow the operator to work productively and calmly, without the danger of harming himself or others.

Before starting work, the worker must be completely familiar with the functions and proper operation of the machine. He must know the precise function of all command and control devices of the machine.

### 2.1 Machine components description



### 2.2 Intended use, improper use

#### INTENDED USE:

The use of this device is intended for washing, high level disinfection and/or sterilization of thermolabile endoscopes by using only chemical agents approved and validated by the device manufacturer. This device is intended only and exclusively for the medical devices treatment.

#### 2.2.1 Application fields

The use of this device is intended only and exclusively for the treatment of medical devices, such as:

- Flexible endoscopes
- Rigid endoscopes



### 2.2.2 Validated cycles

The cycles of high chemical disinfection and sterilization have been validated by the following chemical products:

HIGHLEVEL DISINFECTION CYCLES	C H E MIC ASL
	SteelcoXide-DT (detergent)
HYDROGEREROXIDE	SteelcoXide-A (component A)
	SteelcoXide-B (component B)
	Neodisher SC (detergent)
PERAETIC ACID	Neodisher Septo PAC (paracetic acid)
	Neodisher SC
GLOTAKALDEHTDE	Neodisher Septo GDA
STIR ILIS ATION CYCLES	C HE MIC ALS
	SteelcoXide-DT (detergent)
HYDR OG ERIER OXIDE	SteelcoXide-A (component A)
	SteelcoXide-B (component B)

The improper use of this device is any use other than that for which the machine is intended.

The device has been validated to achieve high-level disinfection and liquid sterilization.

	WARNING
	Any use other than the one intended is forbidden.
<u>.</u>	Improper use of this unit may be hazardous to the operator and may seriously damage the machine itself.
	If the appliance is used in a manner not specified by the manufacturer, protection of the appliance may be compromised.



### 2.3 Important warnings and suggestions

For proper use of the machine, and in order to safeguard employed staff, carefully comply with the following general and specific standards.

#### THE OPERATOR MUST:

- Carefully adhere to the provisions and instructions provided by the employer, managers and supervisors for individual and group safety.
- Use safety devices appropriately and with care, as well as group and individual safety gear provided by the employer.
- Immediately inform the employer, the manager and the supervisor of deficiencies in the aforementioned devices and means, as well as any hazardous conditions which he may become aware of, taking action directly in urgent cases within their scope of responsibilities and abilities to eliminate or reduce the deficiencies or hazards.

#### THE OPERATOR MUST NEVER:

- Remove or modify, without authorization, the safety devices, nor those for signalling and measuring, nor the individual and group safety gear.
- Undertake on his own initiative operations or manoeuvres which are not his responsibility which may compromise safety.
- Insert foreign objects into the electrical parts. Do not insert foreign bodies into the covers of the electrical motors or into the moving parts of the machine.
- Provide power to the machine by tampering with the main switch and the safety devices.



### 2.4 Safety recommendations

- If the new machine seems damaged, contact the retailer before using it.
- Any modification of electrical and hydraulic systems necessary to install the machine must be carried out by qualified, authorised persons only.
- This machine must be operated by trained persons only.
- This machine has been designed for the reprocessing of flexible and rigid endoscopes and the thermal disinfection of washing chambers;
- Any use other than that for which the machine was intended is forbidden.
- The user is forbidden to carry out any work or repairs on the machine.
- Technical Assistance for this machine should be carried out by qualified and authorised operators only.
- The equipment should be installed by authorised persons only.
- The electrical safety of machine is only guaranteed if it is connected to an efficient earth system.
- Take great care when handling detergents and additives: avoid contact, wear gloves, goggles and mask and act in compliance with the safety recommendations indicated by the manufacturer of the chemical products.
- Do not inhale the fumes produced by chemical products.

### WARNING: The chemical products are an irritant for the eyes, in case of contact rinse thoroughly with plenty of water and consult a doctor.

#### If these products come into contact to the skin, rinse with plenty of water.

- The water in the tank is not drinking water.
- Do not lean on the door and do not use it as a step.
- Do not install the equipment in rooms where there is the risk of explosion (ATEX).
- Do not expose the equipment to intense cold.
- Do not wash the machine using high-pressure jets of water.
- The machine reaches a temperature of 80°C during the self-disinfection: take great care to avoid burns.
- Disconnect the machine from the electrical supply before carrying out maintenance work.
- The acoustic pressure of the machine is < 60 dB(A).
- The operator has always to verify before starting of the cycle the presence of the filters water in the sump and their correct positioning.

#### 2.4.1 Inlet water

#### **Physical Properties**

Min. flow pressure	200 kPa (2,0 bar g)
Max. pressure	300 kPa (3,0 bar g)
Max. temperature	35° C
Max. hardness	7° f (70 ppm CaCO3)
Max. conductivity / Ph:	n.a. / 58 pH

#### **Chemical Properties**

Heavy metal ions	Iron	min 0 mg/l (ppm)	max 2 mg/l (ppm)
	Manganese	min 0 mg/l (ppm)	max 2* mg/l (ppm)
	Copper	min 0 mg/l (ppm)	max 2* mg/l (ppm)
	Total heavy metal ions	min 0 mg/l (ppm)	max 10 mg/l (ppm)
Halides	Chloride	min 0 mg/l (ppm)	max 50 mg/l (ppm)
Others ionic contaminants	Phosphates (P <sub>2</sub> O <sub>5)</sub>	min 0 mg/l (ppm)	max 0,2 mg/l (ppm)
	Nitrates (N <sub>i</sub> )	min 0 mg/l (ppm)	max 20* mg/l (ppm)
	Silicates (SiO <sub>2</sub> )	min 0 mg/l (ppm)	max 2 mg/l (ppm)

#### **Microbiological parameters**

Parameter	Parametric Value
Escherichia coli	0/100 ml
Enterococci	0/100 ml
Pseudomonas aeruginosa	0/250 ml
Colony count 22 °C	100 CFU/ml
Colony count 37 °C	20 CFU/ml
Bacterial endotoxins	max 0,25 EU/ml



#### 2.4.2 Rinse water

#### **Physical Properties**

Min. flow pressure	200 kPa (2,0 bar g)
Max. pressure	300 kPa (3,0 bar g)
Max. temperature	35° C max
Max. hardness	1,5° f (15 ppm CaCO3)
Max. conductivity / Ph:	30 µS/cm / 5…8 pH

#### **Chemical Properties**

Heavy metal ions	Iron	min 0 mg/l (ppm)	max 0,2 mg/l (ppm)
	Manganese	min 0 mg/l (ppm)	max 0,2* mg/l (ppm)
	Copper	min 0 mg/l (ppm)	max 0,2* mg/l (ppm)
	Total heavy metal ions	min 0 mg/l (ppm)	max 10 mg/l (ppm)
Halides	Chloride	min 0 mg/l (ppm)	max 10 mg/l (ppm)
Others ionic contaminants	Phosphates (P <sub>2</sub> O <sub>5)</sub>	min 0 mg/l (ppm)	max 0,2 mg/l (ppm)
	Nitrates (N <sub>i</sub> )	min 0 mg/l (ppm)	max 20* mg/l (ppm)
	Silicates (S <sub>i</sub> O <sub>2</sub> )	min 0 mg/l (ppm)	max 0,2 mg/l (ppm)

#### Microbiological parameters

Parameter	Parametric Value
Escherichia coli	0/100 ml
Enterococci	0/100 ml
Pseudomonas aeruginosa	0/100 ml
Mycobacterium Sp.	0/100 ml
Colony count	< 10 CFU/100 ml
Bacterial endotoxins	max 0,25 EU/ml





# 2.5 Recommendations to ensure high quality performance

- The user must oversee the machine during the cycle.
- The injection tube for washing water must always be connected to the appropriated basket.
- When the machine is running do not interrupt the cycle since this jeopardises disinfection.
- Check periodically the correct disinfection with final rinse water checks.
- Use recommended detergents and chemical additives only.
- The use of other products may damage the machine.
- The use of opportune PPEs is compulsory, in order to avoid contact with infected material and to prevent contamination during the handling procedures of medical devices to be reprocessed.
- The chemical products recommended by the manufacturer are those that have been certified and validated by ISO 15883:4, 15883:5-TS and ISO 14937 standard.
- Check that type of chemical product is suitable for the specific washing program used.
- Comply with the instructions provided by the chemical product manufacturer.
- The machine was designed for use with water and chemical additives.
   Do not use organic or other types of solvent as this may result in the risk of explosion or the rapid
- eterioration of certain machine parts. • Residues of solvents or ahkonicclsa,cipcla?r,ticcaunladralmya g"ehysctreced.
- Contact should be avoided.
- Repairs and servicing of this machine must be carried out by authorised persons only.
- Do never use soap powder.
- Do never use foaming detergent.
- Use original accessories only.
- Under no circumstances should the user attempt to carry out repairs.
- The machine is to be used only with the accessories provided by the manufacturer.
- Accessories which are not approved by the manufacturer may compromise the results achieved as well as user safety.
- Do never use chemical products based on chlorides (bleaches, sodium hypochlorite, hydrochloric acid and so on).
- These kinds of chemical detergents irreparably damage the machine and jeopardise the integrity of materials and objects treated.
- Check at every cycle the integrity of the connexions used to connect the endoscopes.
- Wet location.
- Mains supply voltage fluctuations: +/- 10%.
- Overvoltage category: II.
- Pollution degree: 2.
- IP:00.

#### **ATTENTION:**

The taps of the water must be always turned off, as the safety and diagnosis system will be deactivated, in the following situations:

- If the machine is left unused;
- If the machine is disconnected from the electrical connection.

The Manufacturer cannot be held responsible for damage or injury caused by failure to observe the above rules.

The non-observance of these rules produces the total and prompt cancellation of the guarantee.



### 2.6 Residual risks

The appliance includes a series of fixed guards to prevent access to hazardous internal parts or zones. It is however considered that the **EW 1** includes some residual risks. Hereunder for each phase or significant work intervention are useful measures to be taken:

PHASE	BASKET LOADING	
RISK	<b>Contusions and cuts to the upper limbs,</b> due to accidental contact with due to falling or striking against tools, objects and instruments, mainly while loading and handling the basket.	
MEASURE	Assign staff that is instructed and equipped with work equipment (e.g. basket with protections, transport carts) and appropriate clothing and individual protection gear (e.g. shirts and protective gloves).	

PHASE	OBTAINING DETERGENTS/CHEMICAL ADDITIVES
RISK	Contact with body parts with chemical washing products.
MEASURE	Assign staff that is instructed and equipped with appropriate clothing and individual protection gear. Wear clothing, gloves, goggles and mask and act in compliance with the safety recommendations indicated by the manufacturer of the chemical products.
FIRST AID MEASURE	<ul> <li>Immediately take off/remove clothing which has been contaminated or soaked by the product.</li> <li>If the substances come into contact with the skin, wash off affected skin areas immediately and rinse with water.</li> </ul>
RISK	Inhalation of vapours of chemical wash products.
MEASURE	Assign staff that is instructed and equipped with appropriate clothing and individual protection gear. Comply with the safety instructions provided by the manufacturer of the chemical products and if there are none, wear a mask for the protection of the respiratory airways.
RISK	Accidental release of chemical wash product
MEASURE	Do not flush concentrate into drains, surface or ground waters. Collect spillage with adsorbent material (e.g. sand, earth, vermiculite, diatomaceous earth). Flush away minor amounts with plenty of water.
	IN CASE OF CONTACT WITH BODY OR RELEASE OF CHEMICAL PRODUCT LOOK ALWAYS AT THE SAFETY MEASURES INDICATED IN THE CHEMICAL TECHNICAL DATASHEET.

PHASE	MAINTENANCE OF INTERNAL EQUIPMENT	
RISK	Burns of body parts by hot parts of the appliance.	
MEASURE	Allow maintenance to be performed only by trained personnel, equipped with appropria clothing and individual protection gear. Wear suitable clothing and protective gloves.	

PHASE	EMISSION OF HAZARDOUS GAS	
RISK	Inhalation of vapours of hazardous gas.	
MEASURE	With a correct installation, concurring with the manufacturer prescription, using the authorized chemical product and concurring with the rules in force in your country, the machine doesn't generate hazardous gas. However, the machine is supplied with vapours discharge, that have to be connected concurring with the instruction in chapter 3.	



### 2.7 Safety signals used

To inform personnel operating on the machines of obligations of behaviour and residual risks, adequate safety signals (as set forth by 92/58 EEC) are applied to the machine and near the workplace.

#### GENERIC SAFETY SIGNALS:

In particular, labels with signals of obligation, prohibition and danger contained in this manual and pertinent to this machine and most commonly used are:





Electrical risk

Warning! See annex documentation



**Caution hot surface** 

#### INDIVIDUAL SAFETY WEAR:

The evaluation of risks for the health and safety of workers carried out in the workplace and on any equipment used, as well as the evaluation of residual risks as indicated, allow the employer to evaluate the need to adopt the individual protection gear which is most suitable and appropriate to be provided to workers. Considering the type of machine, it is felt that the individual protection gear should be provided to staff.



### 2.8 Training

Instructions for use of the machine will be provided by the **STEELCO INSTALLATION TECHNICIAN** during the start-up phase to **MACHINE OPERATORS** and **MAINTENANCE TECHNICIANS** for their areas of responsibility, which will be thus instructed and trained. Moreover, an appropriate course certificate is issued (see Annex A). It will be the duty of the **EMPLOYER** to check that the degree of staff training is suitable for assigned duties.

#### 2.8.1 Staff qualification

Depending on the difficulty of certain installation operations, and of the operation and maintenance of the system, professional profiles are identified as follows:

#### **IS** INSTALLATION and REPAIR TE CHNICIAN:

Specialized installation and maintenance staff capable of carrying out all machine positioning and installation operations, connection of various systems and machine start-up at the client's place of business, as well as all routine and special maintenance operations.

This operator is responsible for training staff for machine operation and for testing the machine.

#### **As** RESPONSIBLE AUTHORITY FOR THE MACHINE IN THE WORKPLACE:

Specialized staff assigned to the verification of safety devices and procedures for proper use of the machine in complete absence or hazards.

The *responsible authority* is personally responsible for training courses for staff assigned to machine operation and maintenance.

He must ensure that staff assigned to operation has acquired all information required for use and routine maintenance of the machine, registering attendance and documenting comprehension tests.

The *responsible authority must* have a perfect understanding of all command, control and safety devices of the machine.

He must inform all personnel assigned to machine operation and maintenance of the instructions concerning *safety standards*, the *actions to be avoided* and the *first aid interventions* connected with use of the machine and the chemical wash agents it contains.

The *responsible authority* must be aware of all correct procedures for carrying out in absolute absence of danger all operation and maintenance of the machine, as well as all procedures for disposal of any residual pollutants and manufacturing wastes.

He must always be present during extraordinary or routine maintenance and give his *approval to proceed* to staff assigned to operation or to personnel assigned to routine or special maintenance.

The *responsible authority* will be responsible for operation of all command, control and safety devices in the machines of the system.

He shall carry out scheduled verification of those devices in order to ensure their continued operation over time.

#### **AC** MACHINE OPERATOR:

Skilled personnel assigned to machine operation.

The machine operator must be perfectly aware of all of the machine's command and control devices.

Only after approval by the *safety supervisor, the machine operator* must be capable of using the assigned commands to do the following:

- Commissioning and start-up of the machine;
- Loading and unloading of material to be washed in the baskets;
- Operation of the machine in the various possible working modes, such as the start of various programmed wash cycles.
- Programming and setting data from the operator panel, adjustment of single control devices during working phases, starting or resetting of work functions.
- In addition, the *machine operator* must, by making use of all required individual protection gear and following adequate safety measures, be capable of performing some routine maintenance such as cleaning inside the machine, cleaning clogged filters, and disposing of pollutant waste materials produced during working.



### 2.9 Indication of sound level

The value shown refers to the measurement obtained on a machine of the same type as that covered herein and measured with an instrument at a height of 1,5 m at a distance of 1 m from the machine.

AVERAGE SOUND PRESSURE LEVEL: < 60 dB (A)

### 2.10 Transport and storage

Environment conditions:

- Temperature range +5 ... +50 °C;
- Relative Humidity range 20...90% without condensation;
- Ventilation: Air exchange not required (required only if chemical tanks are installed).



### 2.11 Chemical products connections

The chemical product dosing system consists of two different systems: one for aggressive chemical products (neutralizing detergent etc.) and another for disinfectants (peracetic acid, hydrogen peroxide, aldehydes):

• Dosing pump for chemical products.

Presence sensor chemical product.

- The system can be equipped with meter quantity of dispensed product.

Further dosing pumps and accessories can be ordered as optional. Each pump is combined with a corresponding type of chemical, according with the references on the table below.

PERACETIC		
РАА	DISINFECTANT	YELLOW
SC	DETERGENT	BLUE
3rd	CHEMICAL	BLACK

PEROXIDE		
XIDE B	YELLOW	DISINFECTANT
XIDE A	RED	DISINFECTANT
DT	BLUE	DETERGENT



#### **ATTENTION**

In order to guarantee the right treatment of the objects, we suggest the use of specific products. In the case of necessity, ask for advises to the seller or the producer.



#### 2.11.1 Presence sensor of chemical product

Each dosing pump is combined with a sensor that confirm the presence of chemical product inside the container. If the product is scarce, the electronic control system of the machine sends a message on video of lack of product.

#### 2.11.2 Meter quantity of chemical product

Each individual pump is linked to two volumetric sensors in order to measure the quantity of chemical product dispensed. The dual PLC electronic control unit manages the value of the minimum quantity requested and, if necessary, it interrupts the cycle.

#### 2.11.3Container block

## Each container must be locked by the special bracket to ensure the correct data reading.



#### 2.11.34Replacement of chemical product container

To replace the chemical product container perform the following procedure:

- Take the new product container.
- Open the chemical compartment.



• Replace the chemical product container removing the level sensor from the empty tank and put into the new one.





• Close the topper of the chemical product container and place it in the area for the storage of chemical substances.



• Close the chemical compartment.

<u>.</u>	ATTENTION
	The used chemical product can be dangerous if touched or inhaled. Before the use, read carefully the safety information supplied by the manufacturer of the chemical product and the label on the package.
	During the operations of replacement of chemical product container, use the appropriate tools for individual protection (chemical protective gloves, face masks for breathing, etc.).
	The access to the technical compartment, where are located the chemical product containers, is permitted only with keys and to the authorized personal.
	If the chemical tank is in the vacuum or overpressure condition (tank deformation) use the PPE before doing every operation.

If the machine gave an alarm for chemical lack (Error 17-Error 18), after the replacement of the canister, the machine automatically proposes a special cycle for the circuit filling and for the chamber rinsing.

This cycle called **CHEMICAL LOAD** has to be started imperatively after this alarm: it will be not possible to select other cycles until the machine have filled the detergent circuit.

At the end of this cycle the machine will return on stand-by status waiting for a new command.

NB: the chemical loading cycle will not be saved on the historical and will not be printed.

### 2.12 Ambient ventilation requirements

During the normal operation, the machine warms up itself dispersing heat and hot air increasing the humidity value; in the drying phase, these events increase. Therefore, in order to guarantee a comfortable environment with good temperature and humidity for the operator, it is necessary to prepare an air conditioning or air circulation system capable to balance the emissions reported in the installation plan.

The machines with drying system are equipped with an exhaust which can be connected to an external extraction system.

As for machines designed to use Glutaraldehyde, connect the exhaust air connection to an external air evacuation system or put the machine under an extraction hood.



### 2.13 Table of symbols

Symbols installed on the machine:

Â	Electrical risk	
<u></u>	Warning - hot surface	
***	Manufacturer	
$\sim$	Manufacturing date	
$\triangle$	Attention! See the enclosed documentation for important warnings, such as warnings and precautions.	
ī	See instruction for use	
	Protective conductor terminal	
CE	EC Mark	
	WEEE waste disposal	
MD	Medical device indication	
REF	It indicates the catalogue number of manufacturer.	
REP	Authorized Local Representative.	



## 3. USING THE MACHINE

#### 3.1 Checks

Check the quantity of chemical additives present and replace if necessary, as described below:

- Obtain appropriate individual protection gear (gloves for protection from chemical substances, breathing • protection masks, goggles etc.) and the new detergent container.
- Turn off the machine by pressing the OFF button. •
- Follow the instruction on section 2.11.3. •

#### ATTENTION:

The chemical product which is used may be hazardous if touched or inhaled. Prior to use, carefully read the safety information provided by the detergent supplier and the label applied to the package.

### 3.2 Opening and closing the door

- The machine is fitted with an electric door lock to prevent it being opened when the machine is running.
  - To open the door during a wash cycle, interrupt the cycle and remember that:
  - The items inside the machine may be very hot (ONLY DURING THE SELF-DISINFECTION).
     The objects inside the machine could be contaminated
     The entire wash cycle must be repeated.

#### 3.2.1 Door unlocking

In case of power fail or malfunctioning of door lock, it is possible to unlock and open the door by follow the procedure:

- 1. Identify the hole between the door and the cover panel (see the picture).
- 2. Insert the dedicated instrument.
- 3. Keep pushing the dedicated instrument. In this moment the door is unlocked and it is possible to open it. To close the door, keep pushing the dedicate instrument as described on point 3.



After performing the procedure described previously, remember that:	
The items inside the machine could be very hot and contaminated.	
The entire treatment cycle must be repeat.	
Material contaminated with peri acetic acid could be present inside the mac	hine.



### 3.3 Switching on

Turn on the machine following the procedure:

• Press the ON-OFF switch as shown below:



- ٠
- The control panel starts automatically. Check that no alarm messages are present. In negative case remove it. •



### 3.4 Preparation of the endoscopes

In order to be correctly reprocessed inside the EW 1 system, the endoscopes must undergo pre-cleaning and manual washing procedures.

Follow the current applicable National Guideline instructions as well as any internal protocols in force. Some of the main currently available Guidelines are specified below:

C	OUNTRY	GUIDELINES
	EUROPE	ESGE±ESGENA guideline: Cleaning and disinfection in gastrointestinal endoscopy Update 2008.
	ITALY	ANOTE-ANIGEA - Linee guida Pulizia e disinfezione in endoscopia - Update 2011
	FRANCE	Guide de Bonne Pratique de désinfection des dispositifs médicaux - obligatoire depuis le 14 juin 1998. Conseil Supérieur d'Hygiène Publique de France, section prophylaxie des maladies transmissibles. Comité Technique Nationale des Infections Nosocomiales.
	IRELAND	PART 4: RECOMMENDED PRACTICES FOR ENDOSCOPY UNITS Health Service Executive Code of Practice for Decontamination of Reusable Invasive Medical Devices - Review date August 2008.
	GERMANY	Recommendation of the Commission for Hospital Hygiene and Infection Pr evention at the Robert Koch Institute (RKI). Hygiene Requirements for Reprocessing Flexible Endoscopes and Additional Endoscopic Instrumentation - Published in the Federal Health Gazette in April 2002.
	GREAT BRITAIN	National Endoscopy Programme - Decontamination Standards for Flexible Endoscopes - Updated March 2009.
	UNITED STATES	SGNA Society of Gastroenterology Nurses and Associates, Inc. Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes - Revised in 2012.

#### 3.4.1 Check of endoscope instrument connectors

Connectors for the endoscope instruments supplied with the EW 1 machine must be checked daily by machine users and periodically by technicians in charge of maintenance.

The use of not original Steelco connectors or worn connectors can damage the endoscope instruments, activate machine alarms (ex. "channels disconnected", "channel partially clogged" etc.) and most of all does not guarantee a correct disinfection process, endangering both patients' and sanitary personnel's wealth.



WATE R MANIFOLD

# 3.5 Thermal-disinfection cycle in conformity with EN ISO 15883:4

The Endoscope washer-disinfector must perform a daily thermal self-disinfection (regardless the chemical products used).

It is recommended to run the thermal self-disinfection program daily, to minimize contamination risk. This cycle can be performed automatically at a defined time: in case the endoscope washer-disinfector is not used during the weekend, make sure that the automatic self-disinfection is conducted.

If the endoscope washer-disinfector is not used for more than 24 hours, before restarting the operation:

1. Run 2 thermal self-disinfection processes.

2. Control the shelf life of the reprocessing chemicals and replace if expired.

3. Run 1 reprocessing program without endoscopes.

#### MACHINE WITH ONE DOOR

• Connect the CPC connections of the basket water manifold.



Connection with tubes for auto-disinfection (code 671511)

Basket

- Put the basket inside the machine and connect the connection with tubes for auto-disinfection to the relative connector. Then close the leak test inlet by means of the appropriate cap connector as shown in the illustration.
- The yellow tube must be connected to the basket and the opposite side must be free in order to drain the basket water manifold.





#### MACHINE WITH TWO DOORS

•

• Connect the connection for auto-disinfection inside the washing chamber.





Close the leak test inlet by means of the appropriate cap connector as shown in the illustration.



### 3.6 Load preparation

Once the thermal disinfection cycle has finished, remove the sanitation and leak test connectors and proceed with the loading of the endoscope into the basket by carrying out the following procedure:

PUT THE ENDOSCOPE CAREFULLY INTO THE BASKET USING THE APPROPRIATE SUPPORTS

(SEE PHOTO);





INSERT THE CHANNEL SEPARATOR, ACCORDING TO THE ENDOSCOPE TYPE (SEE PHOTO);



COUPLE THE CONNECTION PIPES TO BE USED FOR CONNECTION TO THE ENDOSCOPE WHICH ARE DIFFERENT DEPENDING ON THEIR BRAND AND MODEL (SEE PHOTO);

EXEMPLE OF BIOPSY CHANNEL CONNECTION;



10/05/2022\_REV.1.27\_COD.671064\_A4





EXEMPLE OF "WATER JLECONNECTHONNNE









EXAMPLE OF A WATER CONNECTON POINT;









EXAMPLE OF AN AIR CONNECTION POINT



EXAMPLE OF A LEAK TEST CONECTION POINT









PUT THE BASKET INTO THE MACHINE, MAKING SURE THAT IT FITS INTO THE APPROPRIATE COUPLING;



CONNECT THE CPC QUICK DISCONNECT LEAK TEST COUPLING TO THE MACHINE;





CLOSE THE DOOR AND START THE DISINFECTION CYCLE.

#### ATTENTION

- If the endoscope has less than 7 channels, keep free the CPC connections of the basket water manifold.
- During the wash the yellow tube must always be connected.
- Verify that all channels are properly connected before starting the cycle and at the end, before unloading the endoscope from the basket.



Below are shown some example of basket's type available for the machine:



	ATTENTION
1	<ul> <li>Do never emptying any solid waste into the machine (excrement etc.). This will block the outlet system with pump and destroy the machine.</li> </ul>
	<ul> <li>The treatment cycle has to be activated only if the basket is present into the machine or if it is used a basket equipped with an injection system.</li> </ul>
	<ul> <li>Non observance, even in part, of the rule here indicated, can cause dangerous leakage of water from the door.</li> </ul>


## 3.7 Start the washing cycle

To start the washing cycle follow the procedure:





Select the instrument or instruments by pressing
 PRG

P+ . Scrolling the instruments list, the LED

button lights up only on the selected instruments.



button to start the

 If the chemical product is about to finish a pop-up appears as well as the number of cycles remaining.
 While the pop up is still on the display, prose the

While the pop-up is still on the display, press the





WITH BARCODE (Optional)

Select the washing cycle and press START





 Select the instruments by reading the instrument codes using the barcode reader. After completing this procedure, press the START







Enter the operator code using the barcode ٠ reader and confirm this information by pressing

the START button.



If the chemical product is about to finish a pop-up ٠ appears as well as the number of cycles remaining.

While the pop-up is still on the display, press the



Press the START

washing cycle.

•



### 3.7.1 Endoscope requirements

The flexible endoscopes treated in the EW 1 system must fulfil the following requirements:

- Resistance to a temperature between +10°C and +60°C;
- Resistance of the endoscope channels to a maximum pressure of 1.5 bar;
- Resistance to a pressure of max. 300 mbar during the leak test;
- Resistance to the detergents and disinfectants used;
- Presence of connections where to attach each of the endoscope channels;
- Waterproof with protection against the effects of immersion and, if applicable, equipped with parts that protect the endoscope from humidity.

In order to facilitate the connection of the endoscopes to the EW 1 system, all connections are identifiable in an unambiguous way by using a color coding system (see following figure).



#### **ENDOSCOPE CHANNELS**

- Brown: bioptic channel 1
- Green: bioptic channel 2
- Blue: water channel
- Yellow: air channel
- Grey: suction channel
- Red: auxiliary water channel
- Black: elevator wire channel
- Trasparent: leak test channel

### 3.8 Instruction in case of 24h stop of the machine

#### It is advisable to execute a sanitation cycle every day.

If the machine has not been used for more than 24 hours, it is advisable to carry out a sanitation cycle prior to using the machine.

It is possible to program on the machine an automatic sanitation cycle every day at the same time, that can be programmed by parameters P3.44 and P3.45. When this function is active (P3.43 = 1) the machine automatically start the sanitation cycle (B4) at the time set form the parameters.

WARNING: programming operations are to be performed by the maintenance technician.



ATTENTION: it is important to follow these instruction to allow the correct execution of the self-disinfection procedure.

- I. The basket have to be connected to the sanitation circuit as explained on section 3.5.3.
- II. The machine have to be left with the door shut.



## 4. CONTROL PANEL AND SYMBOLS USED

## 4.1 Switches (loading side)

There are 6 switches with the following functions:

BUTTON	DESCRIPTION	
$\bigcirc$	The switch <b>STOP</b> interrupts the cycle in progress, will be displayed a message indicating that disinfection did not take place, keeps the door locked and if necessary, indicates a high temperature inside the chamber. To return the machine to normal conditions wait the safety drain procedure and then press once more the button to unlock the door.	
	After having selected the programme to be run, pressing the <b>START</b> button to start the cycle.	
D <sup>P1</sup>	Select " <b>STANDAPRADA</b> öycle.	PRG P+
P2 XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	Select " <b>&amp;LF-DISINFECTION 80</b> "cycle.	
<b>&gt;&gt;&gt;</b>	Pressing this button after having selected a programme and before starting it, it is possible to disable the forced air drying step (if selected by the parameter 3.20).	Steelco
PRG P+	Keep pressed for five seconds during Wait or Shutdown to display the Menu. Press once to display during the cycle information about the water and chemical quantity filled and information about the temperature and pressure measured by the transducers installed.	
USB 🜵	On the control panel board there is an USB port that allows the machine programming and data saving.	



## 4.2 Control panel (loading side)

The display shows the following information:



Initially, while the machine is in the stand-by status, it displays the type of program selected, the temperatures, date and time.

By pressing one of the programme switch (P1, P2, P+), the display shows the program selected at the top and at the bottom in red: "press start" or "door open" or any other warnings.

By pressing the switch "P+/PRG" it is possible to scroll all the programs available.

#### LED

The keys of the glass display are touch-sensitive and backlit.



#### BUZZER

The buzzer sounds each time a key is pressed and intermittently in the case of a machine Shutdown, according to the setting of parameters **P1.7**, **P1.8**, **P1.9** through which is possible to set its volume.





PIC. 2	
02.03.188:3700009tank probe23.0 °Ctank probe 223.3 °Csump probe17.7 °Cpress.sensor LT10.244 BarWash pump0.000 BarWash flow0mL/'purge pressure0.038 Bar	During the cycle, by pressing the switch <b>PRG</b> it will be displayed the various temperatures of the device and the values associated to the pressure transducers (Pic. 2).
02.03.18         8:40         00009           Water 1         0.0L           Water 1         0.0L           Water 2         6.1L           Product 1         0mL           Product 2         0mL           Product 4         0mL	It will be displayed also the water quantity filled and the chemical product quantity for each phase, pressing twice
02.03.18 8:40 € 00009 Product 4 0mL	the <b>PRG</b> switch.



PI	C. 3		
	09.04.14 20:32	100273	
	no chemical 1		
	salt loading		PRG
	open door		Pressing <b>PRG</b> a third time it display the screen with the list of alarms and warnings appears during the
	airfilter obstr.		cycle (Pic.3).





PIC. 5	
09.04.14 20:32 00273 STAND. PAA B1 END CYCLE wait 0 no chemical 1	At the end of the cycle, a special window appears as shown on the picture.
09.04.14 20:32 00273 STAND. PAA B1 MAX TEMP. 36.1 °C end program 0 no chemical 1	After few seconds the final windows will be displayed. Now it is possbile to open the door.



PIC. 6	
09.04.14 20:32 00273 STAND. PAA B1 User interrupt wait 15' no chemical 1	In case of manually stop of the cycle, the machine start the automatic safety drain and purge procedure showing a message as on the picture on the side. <b>NB</b> : if, on the phase stopped, where loaded a chemical product, the machine execute also a rinsing before the door unlocking.
09.04.14 20:32 \$\cdots 00273 STAND. PAA B1 NO DISINFECTION wait 15' no chemical 1	After the safety procedure will be displayed a message that warns about the disinfection lack.
09.04.14 20:32 00273 STAND. PAA B1 NO DISINFECTION 15' no chemical 1	Pressing the <b>STOP</b> switch it is possible to unlock the door. The message is still displayed and will be removed after the door opening.



### 4.3 Switches (unloading side – if present)

In case of machine with pass-through door, if the ID operator parameter is disabled (P1.23=0 or 1), at the end of the OK cycle, the unloading door opens automatically.

If this parameter is enabled, P1.23=2 (Ares) or=3 (Steelcodata not Ares with ID unloading machine operator), the ID unloading machine operator procedure at the end of the cycle is:

-at the end of the OK cycle "ended program" popup alternates with "Press start operator" popup;

-the operator must press one of the doors unloading or loading buttons, in order to enter the operator menu (he has to read alternatively "KEY START" or "KEY STOP" bar codes;

-it displays the "ID OPERATOR": the operator must read the operator bar code.

If it is uncorrect, it displays the timed failure message and it comes out the menu.

If it is correct, it displays the operator name and his list number: at that point he must confirm either with the door unloading button or either alternatively with "KEY START" bar code.

The door unlocks and the operator is registered into the supervisor.

-From the unloading machine operator menu, you can exit (return to the OK end cycle screen), by pressing the door locking button or with the "KEY STOP" bar code;

Anyway, you exit automatically by 60 sec. timeout.







### 4.4 Control panel (unloading side – if present)

The display shows the following information:

#### LED

The keys of the glass display are touch-sensitive and backlit.

#### BUZZER

The buzzer sounds each time a key is pressed and intermittently in the case of a machine Shutdown, according to the setting of parameters **P1.10**, **P1.11**, **P1.12** through which is possible to set its volume.





PAG. 45

## 5. MACHINE CONDITIONS

The machine has three possible statuses.

## 5.1 Stand-by

The machine is ready to work. The diagnostic is active. There will be the possible indication of the door open or warnings: lack of chemical product, full memory (historical) or filter obstructed.

5.2	On	cycle

You can reach this point in the procedure by selecting the desired cycle and by pressing the **START** key to start the machine up.

The cycle is running according to its phases.

The diagnostic and the regulators are active.

The user interface gives all indications about the different phases and about the chamber temperature.

### 5.3 Alarm

The diagnostic found an anomaly which cause a block, the cycle is interrupted keeping the door blocked.

The anomaly is shown on the display; after the failure solution is necessary to unlock the machine using the reset procedure (see chapter 6.2).

After that the machine start the safet	v drain procedu	re the door is unlocke	d and it is nossible t	n start a new operation
Alter that the machine start the salet	y uranı procedu		u, anu it is possible t	o start a new operation.

09.04.14	20:32	¢ 00273
STAND. P	AA B1	
drain	a1	
<b>41.1°</b> sp°	AO 🧯 4	1.3° AO
$\bigcirc$		15

20:32

**BLOCK E46** 

pump

00273

09.04.14

st

°¢ ∩

no chemical 1

STAND. PAA B1







## 6. SPECIAL FEATURES

### 6.1 Power failure

In case of a power failure the machine will remember the state it had before the cut off.

### 6.1.1 Power failure statusing "WAIT"

When tension is restored after a power failure during WAIT status the machine returns to the WAIT status waiting a command.

### 6.1.2 Power failure during a cycle or during an alarm

If the power failure occurs during a CYCLE or during an ALARM, after the restoration the machine will show an alarm (power fail): to return on the WAIT status is necessary to execute the reset procedure. The machine automatically starts the safety drain procedure and goes on the waiting state.



### 6.2 Reset procedure

In the event of a **POWER FAILURE** or an **ALARM**, the door remains locked.

To open the door, the reset procedure must be carried out from the keyboard as follows:



4. The following screen will be shown: pressing switch will be displayed the second screen: the machine is unlocking the door.





- 5. When the control panel shows the first picture again keep the door opened for 5 seconds and wait until the alarm is reset.
- **6.** The machine returns to STAND-BY mode.

**NB**: if the alarm occurs during stand-by status, the procedures descripted on point 4 and 5 are not necessary. After the reset, the machine will automatically return on stand-by status.



### ATTENTION



If the machine shutdown persists due to a fault in one of its components (e.g.: faulty probe, unsuitable levels, etc.), the door remain locked and the machine remains inactive. Seek technical assistance.

Should it be necessary to remove the endoscope, switch the machine off and follow the procedure described in chapter 8.1.2.

## 7. WORK PROCEDURES

### 7.1 Introduction

This machine has been designed only and exclusively for the reprocessing of flexible and rigid endoscopes and the thermal-disinfection of washing chambers. Therefore it is subject to continuous contact with aggressive detergents and contaminated instruments.

For this reason it is necessary to provide some useful instructions for the operators who will be using it.

### 7.2 Instructions to personnel

The machine operator, in normal operating conditions, is not subject to risks if he works safely using suitable means of protection.

In order to work safely the operator must:

- Carefully comply with the instructions set forth in this manual.
- Use safety devices appropriately and with care, as well as group and individual safety gear provided in the workplace.
- Personally take action, or inform appropriate persons in the event of deficiencies in the aforementioned devices and means, as well as any hazardous conditions which he may become aware of, taking action directly in urgent cases within their scope of responsibilities and abilities to eliminate or reduce the deficiencies or hazards.

The maintenance technicians, in normal operating conditions, are not subject to risks if they works safely using suitable means of protection.

In order to work safely the maintenance technician must:

- Carefully comply with the instructions set forth in this manual.
- Use safety devices appropriately and with care, as well as group and individual safety gear provided in the workplace.
- Use special care in making repairs or replacing mechanical parts (e.g. drain pump, etc.) on malfunctioning machines which have not completed the disinfection cycle.



# 7.3 Sterilization - Warnings (cycle validated in accordance with UNI EN ISO 14937 standards)

**1.** The cycle is valid only with the combination of the following chemicals:

- SteelcoXide-DT
- SteelcoXide-A and SteelcoXide-B

2. At the end of the cycle the instrument is sterilized but not sterile, as it has no steril barrier system (SBS).

**3.**The material taken out of the sterilizing chamber must be used immediately and handled using sterile gloves and appropriate PPE.

#### **4.LOAD VALIDATION**

Chemical product presence indicators are available. The indicator has to be insered into the appropriate rack in order to check the right chemical dosing.

Once the cycle has finished, check the colour of the indicator (cod. 9992110):





### WARNING

IT IS NECESSARY TO REPEAT THE CYCLE IN THE EVENT THAT THE TEST RESULTS AS BEING PROCESSED YET IN REALITY IT HAS FAILED.



## 8. CLOCK

The main board is equipped with a real time watch. This clock can be set for the following menu:

SETTING  $\rightarrow$  CLOCK  $\rightarrow$  set date and hour.



The time shown is also to qualify the events historical.

## 9. HISTORICAL DATA and PRINTOUT

During the working cycle, the machine memorizes on a RAM all the working data of the wash cycles that have been performed.

If a printer is installed the machine print step by step all the information shown on the example below.

The card is able to record the fields described below for up to a maximum of about 600 cycles in the permanent memory.

When 90% of the memory is full, a signalling pop-up with the message "5% FREE MEMORY" appears on the display.

Printing out the whole historical data, the message will be deleted.

The fields given in the example below are recorded for each cycle:



**************	* * * * * * * * * * * * * *	* * * * * * * * *	
****	• *********	*****	←Customer/distributor (value set by parameter P2.06)
USEL	•		Coser name (value set by parameter P1.01)
Model	: EW I		←Machine model (value set by parameter P2.01)
Machine	: 12001		←Machine serial number (value set by parameter P2.02)
Station	: 0		←Station number (value set by parameter P2.05)
Software	: 7.00		←SW version on main board
Operator	: 13		←Operator nr (or signature space in case of parameter P1.02=no)
INSTRUMENT 1	: RC001		←Endoscope information
	GD900123		
	GASTROSCO	PE	
	STORZ		
B 6 F. ST. PAA	Record	: 00002	<ul> <li>Cycle name and cycle number on historical</li> </ul>
START:	07/06/14	h: 10:14	←Date and hour of cycle start ( <i>standard cycle 17 min</i> )
-> 1: drain T1= 35.9°C	T2= 36.1°C	h: 10:14	←Phase number, name and ending time ←Temperature chamber probe 1 e 2
11 0000 0	10 00.1 0		
-> leak test		h: 10:14	
-> 2: prewash		h: 10:15	
T1= 35.4°C	T2= 35.3°C		
	programm.	executed	
Water 2	6 T.	6 T.	
matter 2	<u>ч</u>	с <u>т</u>	



-> 3: drain h: 10:16 T1= 35.4°C T2= 35.4°C channel drying	
-> 4: treatment h: 10:21 T1= 37.2°C T2= 37.1°C programm. executed Water 2 6 L 6 L Product 1= OK 30 mL 32 mL T>= 35°C t= 300 s t= 300 s (T1= 35.8°C) (T2= 35.8°C)	←Water consumption ←Chemical consumption
-> 5: drain h: 10:22 T1= 37.2°C T2= 37.1°C channel drying	
-> 6: rinsing 1 h: 10:23 T1= 36.4°C T2= 36.4°C programm. executed Water 2 6 L 6 L T>= 0°C t= 20 s t= 20 s (T1= 36.4°C) (T2= 36.4°C)	
-> 7: drain h: 10:24 T1= 36.4°C T2= 36.4°C channel drying	
-> 8: disinfect. h: 10:29 T1= 37.3°C T2= 37.3°C programm. executed Water 1 6 L 6 L Product 2= 0K 60 mL 61 mL T>= 35°C t= 300 s t= 300 s (T1= 35.9°C) (T2= 35.8°C)	
-> 9: drain h: 10:29 T1= 37.3°C T2= 37.3°C channel drying	
-> 10: rinsing 4 h: 10:30 T1= 36.4°C T2= 36.4°C programm. executed Water 1 6 L 6 L T>= 0°C t= 20 s t= 20 s (T1= 36.5°C) (T2= 36.4°C) -> 11: drain h: 10:31 T1= 36.4°C T2= 36.4°C	
channel drying end program: 07/06/14 h: 10:31	←Date and hour of the cycle end ( <i>standard cycle 17 min</i> )
<><><><><><><><><><><><><><><><><><><>	←Alarm/warning description (if present)
Doctor	←Doctor signature
Patient	<ul> <li>         Patient signature     </li> </ul>
××	

Repetition tickets header and cycle's result (this part of printout is controlled by parameter 1.19.)

*******	* * * * * * * * *	* * * * * * * * * * * *	*****
	******	* * * * * * * * *	←Customer/distributor (value set by parameter P2.06)
User	:	********	**** ←User name (value set by parameter P1.01)
Model	:	EW 1	←Machine model (value set by parameter P2.01)
Machine	:	12001	←Machine serial number (value set by parameter P2.02)



Station Software Operator	: 0 : 7.00 : 13	<ul> <li>←Station number (value set by parameter P2.05)</li> <li>←SW version on main board</li> <li>←Operator nr (or signature space in case of parameter P1.02=no)</li> </ul>
INSTRUMENT 1	: RC001 GD900123 GASTROSCOPE STORZ	←Endoscope information
B 6 F. ST. PAA	Record: 000	<ul> <li>Cycle name and cycle number on historical</li> </ul>
START: end program: <><><><><	07/06/14 h: 10: 07/06/14 h: 10: >>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>	<ul> <li></li> <li>↓ Date and hour of cycle start (standard cycle 17 min)</li> <li>↓ Date and hour of the cycle end (standard cycle 17 min)</li> </ul>
EN <><><><>	D CYCLE OK ><><><><><>	←Alarm/warning description (if present)
Doctor		←Doctor signature
Patient		←Patient signature
><		×

## **10. USB PORT**

On the control panel board there is an USB port that allows the machine programming and data saving.

The card is able to record the fields described below for up to a maximum of about 600 cycles in the permanent memory.



## 10.1 Data saving

Insert the USB key on the dedicated port and enter the menu and select the USB menu.

09.04.14	20:32	00273
MENU		
LIOD		
USB		



►HISTORICAL

It is possible to download from the machine the following information and files:

- Cycles
- Parameters
- Historical
- Maintenance historical

The cycle and parameter files can be used to progra m another mach ine or as back -up of the machine.



### 10.2 Data saving during the cycle

By means of the activation of the cycle save parameters required by the user and carried out by the (STEELCO) authorised installer/technician, it will be possible to save data at the end of each washing cycle by leaving the USB memory stick and by carrying out the following procedure:

• Set the parameters P1.04 at 3 and P1.05 to 1.



- Start the washing program; ٠
- At the end of the cycle the machine creates the file with the samples of temperature and pressure probes with • the information of every washing program phases.

To every washing program are associated two files which contain the data structured as below.

Errore. Non si

possono creare The file \*\*\*\*\*G.TXT contains: oggetti dalla N.B.: the \*\*\*\*\*G.TXT will be saved in automatically at the end of the washing program. modifica di :-: 1:-: 4: \*\*\*\*\*\*\* \*\*\*\*\* Information related to \*\*\*\*\*\*\*\*\* End user machine Model EW1 Machine 01234 operator. work station 0 7.01 Software Operator



and



The file \*\*\*\*\*C.TXT contains:

11 II

00036C.TXT

=

N.B.: To save the \*\*\*\*\*C.TXT, ins ert the key in t he dedicated port, enter the menu, select the USB menu and download from the machine the information of historical





### 10.3 Operator archive management

By means of the activation of the cycle save parameters required by the user and carried out by the (STEELCO) authorised installer/technician, it will be possible to save the file archives of the operators present in the machine memory in a file by carrying out the following procedure:

- Insert the USB key into dedicated port.
- Enter menu: USB → OPERATOR → Insert 3<sup>rd</sup> level password → USB EXPORT → Press START button.



The name of file is "**OPERATxxxxx.CSV**", where "**xxxxx**" stands for an optional generic denomination and it contains the data structured as below:

Defines the type of data store.	Defines the length of operator
Identify operator prefix.	OP;6; 1;AB0001;Nome Operatore01 2:AB0002:Nome Operatore02
	3;AB0003;Nome Operatore03 4;AB0004;Nome Operatore04 5:AB0005:Nome Operatore05
Identify the position in the operator list.	6;AB0006;Nome Operatore06 7;AB0007;Nome Operatore07 8:AB0008:Nome Operatore08
	9;AB0009 Nome Operatore09 10;AB0010;Nome Operatore10

Operator code.

FIELD	FIXED CHARACTERISTICS
Type of data store	Nothing
Operator prefix	Length = 2 Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-"minus sign, "." full stop. <b>ATTENTION:</b> set ' <b>OP</b> ' if it is used the barcode reader while set ' <b>00</b> ' if it is used the keyboard (operator identification).
Lenght of operator code	From 1 to 8
Position in the operator list	Progressive order (Maximum number of operators = 60)
Operator code	Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-"minus sign, "." full stop.
Operator denomination	Length ≤ 16 (can be empty) Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-'minus sign, "." full stop.



It is possible to upload the archive operators by insert the USB key into dedicated port and enter the menu: USB  $\rightarrow$  OPERATOR  $\rightarrow$  Insert 3<sup>rd</sup> level password  $\rightarrow$  USB IMPORT  $\rightarrow$  Press START button.

09.04.14 20:32 v 00273 MENU	09.04.14 20:32 200273 MENU	09.04.14 20:32 ¢ 00273 MENU
USB	PASSWORD	OPERATORS
►OPERATORS	► OPERATORS	►USB IMPORT

ATTENTION
<ul> <li>To modify the operator archive file, it is recommended use a text editor (ex. Notepad).</li> <li>If the file contains an operator with the field "OPERATOR CODE" null (no character), the file is considered valid up to the previous operator. All subsequent elements are ignored.</li> <li>If the file contains values that do not comply with the constraints described above, the file is considered incorrect. During the upload it is displayed the first line where the error is present.</li> </ul>



### 10.4 Instrument archive management

By means of the activation of the cycle save parameters required by the user and carried out by the (STEELCO) authorised installer/technician, it will be possible to save the file archives of the instruments present in the machine memory in a file by carrying out the following procedure:

- Insert the USB key into dedicated port;
- Enter menu: USB → INSTRUMENT → Insert 3<sup>rd</sup> level password → USB EXPORT → Press START button.



The name of file is "**INSTRUXXXX.CSV**", where "**XXXXX**" stands for an optional generic denomination and it contains the data structured as below:

Instru	ument					CHAN	NELS			AUX. C	HANNEL	_s	
RC											-		
1	1	G120432	colonEC-380LKp	pentax	100	30000	10	2500	50	130	900	1100	1
2	2	G124804	gastroEG-290Kp	pentax	100	30000	10	2500	50	130	900	1100	1
3	3	G120387	colonEC-380LKp	pentax	100	30000	10	2500	50	130	900	1100	0
4	4	H122654	colonEC-380LKp	pentax	100	30000	10	2500	50	130	900	1100	0
5	5	generico	colon-gastro	pentax	100	30000	10	2500	50	130	900	1100	1
6	6	generico	colon-gastro	pentax	100	30000	10	2500	0	0	0	0	1
1	2	3	4	5	6	7	8	9	10	11	12	13	14

### **TABLE DESCRIPTION**

- 1. Identify the position in the instrument list
- 2. Instrument code
- 3. Serial number
- 4. Instrument denomination
- 5. Instrument manufacturer
- 6. Parameter value of endoscope lower washing water flow limit (ml/')
- 7. Parameter value of endoscope upper washing water flow limit (ml/)
- 8. Parameter value of endoscope lower limit for washing pressure (mb)
- 9. Parameter value of endoscope upper limit for washing pressure (mb)
- **10.** Parameter value of aux. channel lower washing water flow limit (ml/)
- 11. Parameter value of aux. channel upper washing water flow limit (ml/)
- 12. Parameter value of aux. channel lower limit for washing pressure (mb)
- 13. Parameter value of aux. channel upper limit for washing pressure (mb)
- 14. Channels pump

#### **\*ATTENTION:**

IF THE INSTRUMENT HAS NO AUX. CHANNEL, SET THE VALUES 10, 11, 12, 13 TO 0.

ATTENTION\*



\*\* Channels pump activation;

If it is set to 0=not activated channels pump (When there are small instruments);

If it is set to 1= activated channels pump.

FIELD	FIXED CHARACTERISTICS				
Type of data store	Nothing				
Instrument prefix	Length = 2 Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-" minus sign, "." full stop.				
Position in the instrument list	Progressive order (Maximum number of instruments = 60)				
Instrument code	1 ≤ Length ≤ 8 characters Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-" minus sign, "." full stop.				
Serial number	Length ≤ 16 (can be empty) Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-" minus sign, "." full stop.				
Instrument denomination	Length ≤ 16 (can be empty) Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-" minus sign, "." full stop.				
Instrument manufacturer	Length ≤ 16 (can be empty) Allowable characters: 09 number digits, "AZ" uppercase alphabet, "az" lowercase alphabet, " " space, "-" minus sign, "." full stop.				
CHANNELS					
Parameter value of lower washing water flow limit of endoscope.	$0 \le value \le 9999$ Expressed ml (Es. 115 $\rightarrow$ 11.5 litres)				
Parameter value of upper washing water flow limit of endoscope.	$0 \le value \le 9999$ Expressed ml (Es. 360 $\rightarrow$ 36.0 litres)				
Parameter value of lower limit for washing pressure od endoscope.	-1000 $\leq$ value $\leq$ 2500 Expressed millibar (Es. 200 $\rightarrow$ 0.2 bar)				
Parameter value of upper limit for washing pressure of endoscope.	-1000 $\leq$ value $\leq$ 2500 Expressed millibar (Es. 1500 $\rightarrow$ 1.5 bar)				
Channels pump	If it has been setted to 0= no active channels pump (when small instruments are present);				
	AUX. CHANNELS				
Parameter value of lower washing water flow limit of endoscope.	$0 \le value \le 9999$ Expressed ml (Es. 115 $\rightarrow$ 11.5 litres)				
Parameter value of upper washing water flow limit of endoscope.	$0 \le value \le 9999$ Expressed ml (Es. 360 $\rightarrow$ 36.0 litres)				
Parameter value of lower limit for washing pressure of endoscope.	-1000 $\leq$ value $\leq$ 2500 Expressed millibar (Es. 200 $\rightarrow$ 0.2 bar)				
Parameter value of upper limit for washing pressure of endoscope.	-1000 $\leq$ value $\leq$ 2500 Expressed millibar (Es. 1500 $\rightarrow$ 1.5 bar)				
ATTENION	By setting these four values to zero, it means the instrument is without aux. channel.				



It is possible to upload the archive instruments by insert the USB key into dedicated port and enter the menu: USB  $\rightarrow$  INSTRUMENT  $\rightarrow$  Insert 3<sup>rd</sup> level password  $\rightarrow$  USB IMPORT  $\rightarrow$  Press START button.

09.04.14 20:32 ¢ 00273 MENU	09.04.14 20:32 ¢ 00273	09.04.14 20:32 ¢ 00273
USB	PASSWORD	INSTRUMENT
►INSTRUMENT	►INSTRUMENT	►USB IMPORT

	A٦	TENTION
$\wedge$	•	To modify the instrument archive file, it is recommended use a text editor (ex. Notepad).
	•	If the file contains an instrument with the field " <b>INSTRUMENT CODE</b> " null (no character), the file is considered valid up to the previous instrument. All subsequent element are ignored.
	•	If the file contains values that do not comply with the constraints described above, the file is considered incorrect. During the upload it is displayed the first line where the error is present.

## **11. ALARMS and EVENTS LIST**

### 11.1 Logical description of alarm interventions

During machine operation, the operator is aided by **ALARMS** or **ALARM MESSAGES** which use visual signals on the display to advise him of possible anomalies in progress and machine alarms which have intervened. Intervention of an ALARM during operation of the system is signalled by a message on the operator panel. The alarm which appears on the panel remains active until the cause of intervention is removed, and the alarm reset. The intervention of an alarm stops the wash cycle currently in progress.

## 11.2 List of alarm messages

The message includes the number of the alarm that has intervened and its name.

N° AL.	DISPLAY MESSAGE	ALARM DESCRIPTION
E1	power fail	It shows power failure during cycle (Diagnostic is active if P3.04≠END CYCLE).
E2	open load. door	Loading door open and/or unlocked during cycle.
E3	open unload.door	Unloading door open and/or unlocked during cycle.
E4	load.door fail.	Loading door locked but open (discrepancy).
E 5	unload.door fail	Unloading door blocked but open (discrepancy).
E6	doors problem	Incorrect door position ( both open or unblocked).
Е7	unblock.door 1	<ul><li>Loading door problems (loading side):</li><li>a) Overtime lock door ref. P6.14 (problems on the door lock motor).</li><li>b) During the door lock, the door has been opened.</li></ul>
E 8	unblock.door 2	<ul><li>Loading door problems:</li><li>Overtime lock door ref. P6.14 (problems on the door lock motor).</li><li>During block door, the door has been opened.</li></ul>
E 9	unlocking 1fail.	Overtime unlock loading door ref. P6.13.
E10	unlocking 2fail.	Overtime unlock unloading door ref. P6.13.
E11	water 1 lack	Interruption of water 1 filling (lack of new flowmeter impulse over P6.11). (Alarm for water dosing flowmeters (if P8.26=0)).
E12	water 1	Water 1 filling in the tank is not been completed within the maximum time P6.05. (Diagnostic is active if P6.05≠0)
E13	water 1	During water 1 filling in the tank, the water level is deactivated when the quantity measured by the flowmeter reaches or exceeds the parameter P7.09. (Diagnostic is active if $P7.09 \neq 0$ )



N° AL.	DISPLAY MESSAGE	ALARM DESCRIPTION
E14	water 2 lack	Interruption of water 2 filling (lack of new flowmeter impulse over P6.11)
E15	water 2	Water 2 filling in the tank is not been completed within the maximum time P6.05. (Diagnostic is active if P6.05≠0)
E16	water 2	During water 2 filling in the tank, the water level is deactivated when the quantity measured by the flowmeter reaches or exceeds the parameter P7.09. (Diagnostic is active if P7.09≠0)
E17	no chemical 1	Lack of chemical 1. No impulses within the time defined prom parameter (P6.12) with dosing pump active. (Diagnostic is active if P8.27<2)
E18	no chemical 2	Lack of chemical 2. No impulses within the time defined prom parameter (P6.12) with dosing pump active. (Diagnostic is active if P8.28<2)
E20	no chemical 4	Lack of chemical 4. No impulses within the time defined prom parameter (P6.12) with dosing pump active. (Diagnostic is active if P8.30<2)
E21	heating element 1	Discrepancy between the output control and the feedback input for the heating element (delay on the reading = P6.18)
E22	no tank heating	During the heating phase, the water in the tank is fallen below the level of heating ok.
E23	drain problem	The minimum chamber level during the drain has not been reached within the time defined from parameter P6.03 with drain valve open. (Diagnostic is active if P6.03 $\neq$ 0)
E24	fan problem	Diagnostic enabled if P6.19 ≠ 0, in this cases: - the fan pressure switch is open with the fan running at maximum speed - the fan pressure switch is closed with the fan off
E25	drying min°C	During the dryer the air temperature have not reach the temperature set at parameter P7.13, otherwise the temperature is under the set point if the set point is lower than parameter P7.13.
E26	prewash max°C	Tank temperature over maximum setup (P7.14) during prewashing.
E27	tank probe lim°C	Tank temperature over maximum value P7.31 during the working cycles, or P7.32 during sanitation cycle.
E28	dryingprobelim°C	Dryer temperature over maximum value P7.35.
E29	drying max°C	Tank temperature over maximum value P7.33 during the drying phase of working cycles, or P7.34 during sanitation cycle.
E30	tank probe	Tank temperature probe 1 failure.
E31	tank probe 2	Tank temperature probe 2 failure.
E32	drying probe	Dryer temperature probe 1 failure.
E33	drying probe 2	Dryer temperature probe 2 failure.
E34	check temp.	<ul> <li>Appears when P3.02 = YES, only during the treatment phase and if all of these situations are occurred:</li> <li>a) Tank temperature over value P7.12.</li> <li>b) The temperature between the two probes has a difference higher than P7.11.</li> <li>c) The conditions a) and b) are present by over thirty seconds.</li> </ul>
E35	Serial connect.1	No connection between main board and control panel board.
E36	Serial connect.2	No serial connection between expansion board to the keyboard (unloading side).
E37	CAN serialconnec.	No connection between main and slave board (can bus).
E38	TIME	During a treatment phase, the temperature after reaching the set point had a oscillation upon and under the set point. The delay on this reading is fixed at 30". The alarm can be shown also with wrong setting of the chemical loading temperature: chemical loading temperature higher than the set point temperature.
E39	no tank heating	During tank heating phase, the temperature does not increase of 1°C into the time given by parameter P6.01.
E40	chan.obstructed	<ul> <li>When the washing pump is ON, the water flowmeter (in the washing cycle with instruments) measure a flow lower than the minimum value defined than:</li> <li>parameter for single instrument.</li> <li>parameter P7.36 for selection of 2 instruments.</li> <li>parameter P7.60 for selection of 3(+) instruments.</li> <li>(Diagnostic is active if P3.48≠0)</li> </ul>



N° AL.	DISPLAY MESSAGE	ALARM DESCRIPTION
E41	chan.disconnect.	<ul> <li>When the washing pump is ON, the water flowmeter (in the washing cycle with instruments) measure a flow higher than the maximum value defined than:</li> <li>parameter for single instrument.</li> <li>parameter P7.37 for selection of 2 instruments.</li> <li>parameter P7.61 for selection of 3(+) instruments.</li> <li>(Diagnostic is active if P3.48≠0)</li> </ul>
E44	instrum.leakage	During the cycle, activation number of leak valve (endoscope 1) has reached the maximum limit set by parameter P7.58. (it is excluded the first activation to put under pressure the instrument). (Diagnostic is active if P3.47≠0, P3.50=0 and P7.58>0)
E45	instrum.leakage	During the cycle, activation number of leak valve (endoscope 1) has reached the maximum limit set by parameter P7.58. (it is excluded the first activation to put under pressure the instrument). (Diagnostic is active if P3.47≠0, P3.50=0 and P7.58>0)
E46	washing flow	<ul> <li>Failure of the flowmeter or washing pump:</li> <li>When the washing pump is ON, the water flowmeter doesn't read a new impulse within P6.07.</li> <li>When the washing pump is OFF, the water flowmeter read a number of impulses that exceed P7.20.</li> <li>(Diagnostic is active if P3.48≠0)</li> </ul>
E47	flowmeter fail.1	The chemical flowmeter 1 had count an impulse number higher than set point + P7.21. (Diagnostic is active if P8.27<2)
E48	flowmeter fail.2	The chemical flowmeter 2 had count an impulse number higher than set point + P7.21. (Diagnostic is active if P8.27<2)
E50	flowmeter fail.4	The chemical flowmeter 4 had count an impulse number higher than set point + P7.21. (Diagnostic is active if P8.30<2)
E51	water 1 flowmeter	The water 1 flowmeter (control) had count an impulse number higher than set point + P7.20.
E52	water 1 flowmeter	The redundancy water flowmeter 1 had count an impulse number higher than set point + P7.20. (Diagnostic is active if P8.26=0)
E53	water 1 flowmeter	Discrepancy between control and redundancy flowmeters higher than P7.10. (Diagnostic is active if P8.26=0)
E54	Purge pressure	After the time set by P6.16 from the command to open the drain purge valve has been given (following the purge procedure pressurisation stage), the purge pressure has not dropped below the [P7.54 - P7.61] threshold or below the [P7.03 - P7.61] threshold. Active diagnostic if P3.39 is not =0 and P3.50=0.
E55	Conductivity probe	Fault in conductivity probe
E56	Conductivity	The conductivity value exceeds the maximum P7.60 threshold level, with a 5-second debounce
E57	flowmeter fail.1	Discrepancy between the measured value by time and by flowmeter higher than P7.10. (DETERGENT). (Diagnostic is active if P8.27=0)
E58	flowmeter fail.2	Discrepancy between the measured value by control flowmeter and redundancy flowmeter higher than P7.10. (DISINFECTANT). (Diagnostic is active if P8.28=0)
E60	flowmeter fail.4	Discrepancy between the measured value by control flowmeter and redundancy flowmeter higher than P7.10. (DISINFECTANT – bi-component). (Diagnostic is active if P8.30=0)
E61	no water 1	During water 1 filling in the break tank, the water level is not reaches level within time set by P6.06.
E62	break tank	Emergency level is active.
E63	break tank	Lower level is deactive when work level is active.
E64	drain problem	The minimum break tank level during the drain has not been reached within the time defined from parameter P6.04 with drain valve open.
E65	tank probe lim°C	<ul> <li>The chamber temperature has exceeded of 5°C the set point during the treatment phase (controlled by the control probe), in case of:</li> <li>Temperature set point ≠ 0.</li> <li>Water filling completed;</li> </ul>



N° AL.	DISPLAY MESSAGE	ALARM DESCRIPTION
E66	dryingprobelim°C	The air temperature has exceeded of 5°C the set point during the drying phase (controlled by the control drying probe), when the temperature set point $\neq$ 0.
E67	check temp.	<ul> <li>Appears when P3.33=2, <u>only during the drying phase</u> and if all of these situations are occurred:</li> <li>a) Air temperature is over P7.06 value.</li> <li>b) Discrepancy between the two probes is over the P7.05 value.</li> <li>c) The conditions a) and b) are present by over thirty seconds.</li> </ul>
E68	chan.obstructed	When the washing pump is ON, the water flowmeter (in the washing cycle without instrument) measure a flow lower than P7.44. (Diagnostic is active if P3.48≠0)
E69	chan.disconnect.	When the washing pump is ON, the water flowmeter (in the washing cycle without instrument) measure a flow higher than P7.45. (Diagnostic is active if P3.48≠0)
E70	press.sensor LT1	Failure on the leak test pressure transducer (instrument 1). (Diagnostic is active if P3.47 $\neq$ 0)
E71	press.sensor LT2	Failure on the leak test pressure transducer (instrument 2). (Diagnostic is active if P3.47 $\neq$ 0)
E72	max pressure LT1	The leak test pressure is higher than P7.38. (Diagnostic is active if P3.47≠0)
E73	max pressure LT2	The leak test pressure is higher than P7.38. (Diagnostic is active if P3.47≠0)
E74	low pressure LT1	<ul> <li>During a cycle, the leak test pressure is lower than the minimum value, defined by one of the following conditions:</li> <li>The set point P7.39 is not reached within P6.54</li> <li>Pressure lower than P7.39 – P.7.40 within P6.56</li> <li>Pressure lower than P7.39 –P.7.41 during the washing cycle.</li> <li>This alarm occurs also if the endoscope LT is failed, i.e. if the endoscope channel has a leakage. (Diagnostic is active if 3.47≠0)</li> </ul>
E75	low pressure LT2	<ul> <li>During a cycle, the leak test pressure is lower than the minimum value, defined by one of the following conditions:</li> <li>The set point P7.39 is not reached within P6.54</li> <li>Pressure lower than P7.39 – P.7.40 within P6.56</li> <li>Pressure lower than P7.39 –P.7.41 during the washing cycle.</li> <li>This alarm occurs also if the endoscope LT is failed, i.e. if the endoscope channel hasa leakage. (Diagnostic is active if P3.47≠0)</li> </ul>
E76	LT1 drain failure	The leak test pressure is not fallen under the max pressure value (P7.42) after the drain. (Diagnostic is active if P3.47 $\neq$ 0)
E77	LT2 drain failure	The leak test pressure is not fallen under the max pressure value (P7.42) after the drain. (Diagnostic is active if P3.47≠0)
E78	water 2 flowmeter	The water 2 flowmeter count a number of impulses that exceeds P7.20 with water filling solenoid valve OFF.
E79	water 2 lack	During water 2 filling in the break tank, the water level is not reached level within time set by P6.06.
E80	washing pump	Washing pump pressure transducer failure. (Diagnostic is active if P3.48≠0)
E81	chan.obstructed	<ul> <li>When the washing pump is ON, the transducer measures a pressure higher than:</li> <li>parameter for single instrument.</li> <li>parameter P7.47 for selection of 2 instruments.</li> <li>parameter P7.63 for selection of 3(+) instruments.</li> <li>(Diagnostic is active if P3.48≠0).</li> </ul>
E82	chan.disconnect.	<ul> <li>When the washing pump is ON, the transducer measures a pressure lower than:</li> <li>parameter for single instrument.</li> <li>parameter P7.46 for selection of 2 instruments.</li> <li>parameter P7.62 for selection of 3(+) instruments.</li> <li>(Diagnostic is active if P3.48≠0)</li> </ul>
E83	chan.obstructed	When the washing pump is ON, the transducer measures a pressure higher than P7.49. This alarm occurs during a treatment cycle or during a self-disinfection cycle. Check the filter and if necessary, replace. (Diagnostic is active if P3.48≠0).
E84	chan.disconnect.	When the washing pump is ON, the transducer measures a pressure lower than the P7.48. This alarm occurs during a treatment cycle or during a self-disinfection cycle. (Diagnostic is active if P3.48≠0)



DISPLAY MESSAGE	ALARM DESCRIPTION
lim °C sump probe	Functional limit for the over temperature protection. The temperature measured form the sump probe is higher than the value of parameter P7.31 on working cycle and P7.32 on sanitation cycles. (Diagnostic is active if P3.49≠0)
Sump probe	Failure on the sump probe. (Diagnostic is active if P3.49≠0)
lim °C sump probe	The temperature measured from the sump probe is 5°C higher than the set point set for that phase (only if the set point is different form 0°C) (Diagnostic is active if P3.49≠0)
flowmeter fail.1	The disinfectant redundancy flowmeter measures an excess of impulses higher than P7.21 with dosing pump OFF. (Diagnostic is active if P8.27≠0)
flowmeter fail.2	The disinfectant redundancy flowmeter measures an excess of impulses higher than P7.21 with dosing pump OFF. (Diagnostic is active if P8.28≠0)
flowmeter fail.4	The disinfectant redundancy flowmeter (for bi-component disinfectant) measure an excess of impulses higher than P7.21 with dosing pump OFF. (Diagnostic is active if P8.30 $\neq$ 0)
sanificat. Flow	<ul> <li>The pressure switch of the sanitation circuit is OFF, but the sanitation valve is open and the pump is ON</li> <li>The pressure switch of the sanitation circuit is ON, but the sanitation valve is OFF</li> <li>Delay on the pressure switch status reading = P6.52</li> </ul>
purge pressure	Failure on air pressure sensor for purge. (Diagnostic is active if P3.39≠0)
purge pressure	The setpoint set by P7.54 is not reached within maximum time set by P6.59. (Diagnostic is active if P3.39≠0)
Max. air pressure	The air pressure measured from the purge transducer (when the purge activity is in progress) has exceeds the maximum limit set by P7.57. (Diagnostic is active if $P3.39\neq0$ )
Min. air pressure	The air pressure measured from the purge transducer (when the purge activity is in progress) has exceeds the minimum limit set by P7.59. (Diagnostic is active if P3.50=YES)
CAN serialconnec.	(Diagnostic is active if P1.23≠0) No connection between main board and gateway (CAN/Ethernet) for interface the supervisor (Diagnostic is active when the machine is in wait mode).
AUX chan block	<ul> <li>(Diagnostics are enabled if p3.48 = 2)</li> <li>Possible damage has been caused to the flow meter or to the auxiliary instrument channel pump:</li> <li>When the aux channel pump is active, the aux channel flow meter does not signal any new pulses beyond the time specified in p6.07;</li> <li>When the aux channel pump is not in use, the aux channel flow meter signals an excessive amount of pulses exceeding the threshold specified in p7.20 (diagnostics are disabled during fan operation).</li> </ul>
AUX press probe	(Diagnostics are enabled if p3.48 = 2) Failure of the auxiliary instrument channel pump pressure sensor
AUX chan block	<ul> <li>(Diagnostics enabled if p3.48 = 2 and a single instrument is selected)</li> <li>When the aux channel pump is active, the aux channel water circulation flow meter in cycles with instruments in the chamber has measured a flow rate lower than the min. threshold defined by the parameter of the instrument itself in the instrument file (blocked aux instrument channel) in the following conditions:</li> <li>The calculation of the water flow rate is updated every 5 seconds while diagnostics are updated every 10 seconds after the first minute by the activation of the water circulation system;</li> <li>An alarm is triggered if the flow rate remains continuously below the threshold for a period of time exceeding that specified in p6.26;</li> <li>Diagnostics are disabled in case multiple instruments are selected.</li> </ul>
	DISPLAY MESSAGE lim °C sump probe Sump probe lim °C sump probe flowmeter fail.1 flowmeter fail.2 flowmeter fail.4 sanificat. Flow purge pressure purge pressure Max. air pressure Max. air pressure CAN serialconnec. AUX chan block AUX chan block AUX chan block



N° AL.	DISPLAY MESSAGE	ALARM DESCRIPTION
E104	AUX chan disc	<ul> <li>(Diagnostics are enabled if p3.48 = 2 and a single instrument is selected)</li> <li>When the aux channel pump is active, the aux channel water circulation flow meter in cycles with instruments in the chamber has measured a flow rate exceeding the max. threshold defined by the parameter of the instrument itself in the instrument file (disconnected aux instrument channel) in the following conditions:</li> <li>The calculation of the water flow rate is updated every 5 seconds while the diagnostics are updated every 10 seconds after the first minute by the activation of the water circulation system;</li> <li>An alarm is triggered if the flow rate remains continuously above the threshold for a period of time exceeding that specified in p6.26;</li> <li>Diagnostics are disabled in case multiple instruments are selected.</li> </ul>
E105	AUX chan block	<ul> <li>Diagnostics are enabled if p3.48 = 2 and a single instrument is selected)</li> <li>When the aux channel pump is active, the aux channel pressure transducer in cycles with instruments in the chamber has measured pressure exceeding the max. threshold defined by the parameter of the instrument itself in the instrument file (blocked aux instrument channel) in the following conditions:</li> <li>An alarm is triggered if the pressure remains continuously above the threshold for a period of time exceeding that specified in p6.27;</li> <li>Diagnostics are enabled after the first minute by the activation of the water circulation system;</li> <li>Diagnostics are disabled in case multiple instruments are selected.</li> </ul>
E106	AUX chan disc	<ul> <li>(Diagnostics are enabled if p3.48 = 2 and a single instrument is selected)</li> <li>When the aux channel pump is active, the aux channel pump pressure transducer in cycles with instruments in the chamber has measured pressure lower than the min. threshold defined by the parameter of the instrument itself in the instrument file (disconnected aux instrument channel) in the following conditions:</li> <li>An alarm is triggered if the pressure remains continuously below the threshold for a period of time exceeding that specified in p6.27;</li> <li>Diagnostics are enabled after the first minute by the activation of the water circulation system;</li> <li>Diagnostics are disabled in case multiple instruments are selected.</li> </ul>
E107	AUX press probe	(Diagnostics are enabled if p3.48 = 2) When the water pressure of the aux instrument channel exceeds the maximum tolerance level of the instrument (2000 mbar). Diagnostics are enabled only when cycle is in progress (not while in stand-by status).
E108	Sanif. flow	With sanification valve open and washing pump activated, the water flow is not detected in the hydraulic field of the water flowmeters 1 (lack of new flowmeter pulse from beyond P6.11). Alarm active when for at least one of the two water dosing flowmeters this condition is verified (if P8.26=0).



## 11.3 List of warnings messages

DISPLAY MESSAGE	DESCRIPTION
press start	It is possible to start a cycle during a stand-by state.
no chemical 1	It informs that the product (DETERGENT) is: • run out if P8.27=2. • ending if P8.27<2.
no chemical 2	It informs that the product (DISINFECTANT) is: • run out if P8.28=2. • ending if P8.28<2.
no chemical 4	It informs that the product (DISINFECTANT – BI-COMPONENT) is: • run out if P8.30=2. • ending if P8.30<2.
salt loading	Warning to carry out loading procedure increases following the performance of a certain number of softener resin re-generations defined in p7.27 (as for configured regeneration, see $p7.26 > 10$ ).
- open door -	Inform that the door is open.
wait	Generic warn that inform to wait before to do a new action.
close door!	It warns to close the door open to allow proper initialization of the door.
ON PRINT	It appears when you try to start a cycle during the printing of the previous cycle.
UNBLOCKING wait	Automatic unlocking procedure after the unblock switch sequence after an alarm that interrupted a cycle.
ALARM	The safety procedure after an alarm is completed. It requests the manual unlock of the door pressing the STOP switch.
User interrupt wait	It warns that is in progress the safety procedure after a manual stop of working cycle.
NO DISINFECTION	Message shown at the end of the safety procedure after user interruption of the cycle. It requests the manual unlock of the door pressing the STOP switch.
END CYCLE	It is in progress the end cycle procedure after the last phase of cycle (the door remains locked).
airfilter obstr.	The air filter is obstructed.
Supervisor wait	It is in progress the data download of the performed cycle. It is necessary to wait to start new washing cycle.
NO Authorization	Response from Ares relating to denial of required instrument authorisation
Unidentified instrument	Response from Ares relating to required instrument not recognised by the supervisor's archive of reference
X REPETITION	(x = 1, 2, 3) states that repetition number "x" of the stage is in progress (as for stages which can be set with repetitions)

## 11.4 List of historical events

EVENT	DISPLAY MESSAGE	DESCRIPTION
FROM 1 TO 107	SAME LIST OF ALARMS	(SEE ALARMS LIST)
100	ОК	CYCLE ENDS WITH SUCCESS
101	NO DISINFECTION	CYCLE HAS BEEN INTERRUPTED



## **12. INSTALLATION AND OPERATION**

As highlighted in Annex A.3 of the ISO 15883:4 2009 standard, from the moment in which the EW 1 system has been installed on the site where it will subsequently be used, the customer/user shall be held responsible for its proper installation and operation.

This includes:

- Checks and installation tests;
- Operational tests;
- Performance qualification tests;
- Periodic tests;
- The use of process chemicals recommended by the manufacturer;
- The use of the EW 1 system in accordance with the manufacturer's instructions (by limiting the use of reprocessable devices in the EW 1 system in favor of those recommended by the manufacturer).

## **13. RACK VERSION**

The main device can be in ventcal RAC K configuration.

It refers to a configuration dealing with the dsposing of twoidentical E W1, which arevertically placed, one on the machinetop and the otherone on themachinebottom.

Both these machinesare equipped with:

- > 2 screens(one corresponding to themachinetop and the other one corresponding to themachine bottom);
- 2 printers (one corresponding to themachinetop and the other one corresponding table machine bottom);
- a barcodereader system (one corresponding to themachinetop and the other one corresponding to the machinebottom);
- > a systemof filters dedicated to eab machine
- > a singlechemical compartmentn common forthe two machines
- > a single electrical cabinestituated on thetop of the machine. That system has beenestigned so that each machine is able to work on its own;



## **14. MAINTENANCE**

### 14.1 General recommendations on maintenance

This machine has been designed only and exclusively for the reprocessing of flexible and rigid endoscopes and the thermal-disinfection of washing chambers. Therefore it is subject to continuous contact with aggressive detergents and contaminated instruments.

For this reason it is necessary to provide some useful instructions for the operators who will be performing maintenance on it.

The maintenance technicians, in normal operating conditions, are not subject to risks if they work safely using suitable means of protection.

In order to work safely the maintenance technician must:

- Carefully comply with the instructions set forth in this manual.
- Use safety devices appropriately and with care, as well as group and individual safety gear provided in the workplace.
- Use special care in making repairs or replacing mechanical parts (e.g. drain pump, etc.) on malfunctioning machines which have not completed the thermal disinfection cycle.

Maintenance operations for the machine described in this manual can be divided into "Routine Maintenance" and "Special Maintenance".

### **GENERAL GUIDELINES:**

#### **MACHINE STATUS**

The machine must not be powered electrically and the dedicated safety device must be in the OFF position. The person performing the task must ensure that there is no-one around the machine during this operation.

#### SAFETY SYSTEMS TO BE ADOPTED

The operation must be carried out in compliance with standards governing the use of disinfectant substances used (see technical information for the product being used), in compliance with standards concerning contact with parts of the machine which may be contaminated by pathogenic materials and with use of individual protection gear.

### 14.1.1 Maintenance request

The machine displays the "**MAINTENANCE**" warning after a specified time or after a specified numbers of working hours according to parameter **P6.48**. This warning doesn't affect the normal use of machine. The service technician must do the maintenance operations in the shortest possible time.



To clear the "MAINTENANCE" warning, follow the procedure:

- 1. Do the general maintenance of machine;
- 2. Enter the menu:

UTILITY  $\rightarrow$  MAINTENANIQUE del passworrsde r-to MSAINTEENGAINSCTEERR  $\rightarrow$  PRTe button.





### **14.2 Table of routine maintenance**

The following table shows the various routine maintenance tasks, their frequency, who is to perform them and the reference to the specific intervention form.

Each single task is more fully explained in the single reference forms.

Even if the water supply is relatively soft, the high temperature can cause the formation of residues which may create problems with the heating element, compromising the correct wash cycle and the reaching of the disinfection temperature.

For these reasons it is advisable to carry out regular cleaning as described below.



Cool Co		ĺ		Z	BO	SC	0	PES WASHER		SCE
				Prot	gran	amu	d n	naintenance scheme		RENE
Community	Step			mont	hs			A sets cites of		SIFE
Components	make every	3	<mark>6</mark> 9	12	15	18	24	Activity		ł
Chamber filters	make every day							rake off filters and cleaning.	10'	M1
Water solenoid filter	make every		×			×	0	Check, clean and if necessary replace.	10'	M4
Pre filter dryer F5	make every 100 hours						<u> </u>	Replace.	2	M5
Dryer HEPA filter	make every 300 hours	2					<u></u>	Replace.	2'	M5
Temperature probes	make every			×			×	During periodic validation, check the sensor status.	.09	M2
Safety thermostat	make every			×			×	/erify the sensor.	5	
Chemical dosing pump	make every		×	×		×	×	Check the presence of leakage.	2	
Chemical tank level sensor	make every		×	×		×	×	Check and clean the suction filter.	4	
Connection pipe of dosing pump	make every		×	×		×	×	Check of crashing, any leakage or hardening.	10'	
Washing arms	every week						500	Check for free rotation. Dpen the cleaning caps and wash inside: check and in case cleaning the nozzle.	30'	M3
Door gasket	make every		×	×		×	×	/erify the gasket and replace after 1000 cycles.	20'	
Washing pumps	make every			×			×	Check for water leakage from the arm seal.	2	1
Water heating element	make every			×			×	Check for water leakage fro the gasket.	ŀ	
Dryer heating elements (OPTIONAL)	none						0	Dperation is checked by the control system.		
Water solenoid valves	make every			×			×	Check for any leaks, if necessary remove and clean the membrane seat.	õ	
Drain pump	make every			×			×	Check for any leaks, if necessary remove and clean the membrane seat.	3	
Pressure switches	make every			×			× s tt	Dperation is checked by the control system. In case of defect of control system of water levels, go on by empting he tank, olowing inside the black pipe connected to the pressure swtch, in oder to free from obstructions.	10'	
Pipe of unloading water	make every			×			×	Check the situation of pipe and the seal.	õ	
Pipes of loading water	make every			×			×	Check the situation of pipe and the seal.	3	

### TABLE OF ROUTINE MAINTENANCE TASKS

#### N.B.:

Routine maintenance tasks must be performed at the intervals set forth in the table. It is however advisable to carry out single cleaning tasks anytime you feel they may be necessary.


In cas	e the	machine	requires	the	replacement	of	one	or	more	components,	please	refer	to	the
manufa	icturer	's spare p	art list.											

	It is advisable to carry out a general check-up and to clean the appliance regularly, particularly if the supply water is very hard.			
$\square$	Particular attention should be paid to heating element and the probe of thermostats.			

### WARNING:

- Do not clean the machine outside with high pressure water.
- Please contact the retailer that supplies your cleaning products for details of recommended methods and products for sanitizing the machine regularly.
- The machine has a safety thermostat that shuts down the power supply to the heating elements in the event of overheating.
- To re-start the appliance the fault that caused overheating must be corrected.
- Empty the tank of compressor before performing the maintenance.

### **Every 12 months**

- Clean the diaphragms of solenoid valves and replace if necessary.
- Clean the thermostat probe.

Even if the supply water is soft, the high working temperatures may cause limescale to build-up. Apart from damaging the resistors, limescale can also clog the nozzles in which case the correct tank temperature for thermodisinfection may not be reached.





### **CLEANING OF WASHING CHAMBER DRAIN FILTERS**

Μ1

Worker: Ac

Frequency of Intervention: every day

**METHOD OF INTERVENTION:** clean the washing chamber drain filters in the following manner:

- Open the washing chamber door and extract the basket.
- Extract the drain water filtering assembly from the chamber.







Unscrew the threaded pin and remove the cover of the drain water filter basket.





- Clean the drain water filter basket. Remove residues deposited during various wash cycles. •
- Remove and clean any deposits and incrustations from the wash chamber drain. •
- Replace the clean filter on the wash chamber drain. •
- Put the cover for the drain water filter back in place. Lock it in position with the threaded pin. •
- Put the drain water filter group back in the washing chamber.



	CLEANING OF WASHING CHAMBER TEMPERATURE PROBE							
M2	Worker: Is Frequency of Intervention: 6 months							
METH	<b>METHOD OF INTERVENTION:</b> clean the washing chamber temperature probe in the following manner:							
O     C     ar	pen the washi heck the wash nd an appropri	ing chamber ning chambe iate deterger	door and extract the basket. r thermostat probe and clean it of any deposits or lime incrustations using a damp cloth nt.					
Take	care not to d	lamage or m	nove the probe					
	CLEANING OF WASHING ARMS							







- Put the plugs back in place at the ends of the wash arms. Make sure the gasket is properly positioned and in good
- condition. Replace it if necessary.
- Put the rotors back on the machine. Lock them in place with the previously removed fastening pin.

### CHECK OF WASHING ARMS ROTATION

 Worker: Ac
 Frequency of Intervention: every day

 METHOD OF INTERVENTION: check the washing arms rotation as follows:

- Open the washing chamber door and extract the basket.
- Check the free washing arm rotation for both upper and lower arms checking that no slowdown is shown from the washing arms.

### CLEANING THE EXTERNAL BODY OF THE MACHINE

Worker: Ac Frequency of Intervention: every day

### METHOD OF CLEANING OUTER BODY

Use a damp cloth to clean the outer body of the machine.

Use only neutral detergents.

Do not use abrasive detergents or solvents and/or thinners of any kind.

### METHOD OF CLEANING MARKING LABEL

Use a damp cloth to clean the marking label surface. Use only water or isopropyl alcohol. Do not use abrasive detergents or solvents and/or thinners of any kind.

### METHOD OF CLEANING CONTROL PANEL

Clean the control panel using only a soft cloth dampened with a product for the cleaning of plastic materials.



### LIMESCALE REMOVAL TREATMENT

Worker: Ac

Frequency of Intervention: whenever necessary

### **METHOD OF INTERVENTION:**

Use a descaling agent (we recommend vinegar) during an empty washing cycle with cold water (this is usually carried out every week unless the quality of the water requires a daily treatment in order to prevent the build-up of limescale and the blockage of the water jets).

As regards the quantity of the product to use, please comply with the instructions given on the technical data sheet of the product itself. In case vinegar is used, use 0.5 litres.

The descaling product must be poured into a container of the same size, positioned on an empty loading basket. Use a washing programme with water at room temperature, without activating the drying cycle.

### ATTENTION



Even if the feed water only contains a small amount of limescale, high temperatures can generate the formation of limescale residues. This, as well as problems that could be caused to the heating element, may cause the blockage of the nozzles, jeopardising the correct washing process and preventing the ideal disinfection temperature in the tank to be reached.

### 14.4 Procedure for special maintenance work

All special maintenance work is to be performed only by qualified, skilled personnel. A table is shown below which includes possible special maintenance work that may be required. If your machine should require special maintenance, please contact your retailer/distributor.



### 14.5 Table of special maintenance

See scheduled maintenance table.













### **REPLACEMENT OF WATER FILTERS** M7 Worker: Is Frequency of Intervention: every year **METHOD OF INTERVENTION:** replace water filters (0,1 micron and 0,45 micron) as described below: Remove water filters by unscrew pin of clamp (yellow arrow) and remove gasket between filters and machine pipe ٠ fitting. Replace water filters. Pay attention to flow direction (green arrow). Every filter has an arrow on the body to indicate • the correct flow of water. Insert the gasket. ٠ Check the absence of leak. • STEELCO66003





1 10. IT2460



	CHECKING OF SELF-DISINFECTION CONNECTOR				
M8 Worker: I	Frequency of Intervention: every 3 months				
METHOD OF INTERVENT	<b>ON:</b> check the self-disinfection connector into the washing chamber as described below:				
<ul> <li>Check that self-disinference leakage (check during</li> <li>Check that male condition disconnect male connect male connect</li> <li>ATTENTION:</li> <li>Do not force the connection female connector.</li> </ul>	tion female connector (yellow arrow) is fixed in its place and that the gasket doesn't have elf-disinfection cycle). ector (green arrow) can be freely connected/disconnected to/from female connector. To ctor, press grey button placed on female connector (red arrow).				



CHECKING/REPLACEMENT OF BASKET CONNECTION GASKET (only for one door version)

Worker: Is Frequency of Intervention: every 3 months

**METHOD OF INTERVENTION:** check of integrity of gasket for the connection between machine and basket, as described below:

• Check and/or replace the OR gasket on the connector machine and basket (yellow arrow).

### ATTENTION:

M9

In case the OR-gasket is broken or too worn, the machine will show a disconnection alarm of endoscope channels.







	REPLACEMENT OF AIR FILTER FOR LEAK TEST							
M10	Worker: <b>Is</b>	Frequency of Intervention: every year						
METHOD OF INTE	<b>IETHOD OF INTERVENTION:</b> replace leak test filter (0,2 micron) as described below:							
<ul><li>Remove cover</li><li>The filter is pla</li><li>Replace the filter</li></ul>	Remove cover panels from the machine. The filter is placed on the left side of machine, looking from unloading side. Replace the filter (yellow arrow). Pay attention to the suction direction indicated on the filter body (green arrow).							



CHE	ECKING THE S	UCTION LANCE OF CHEMICAL PRODUCT 2 (DESINFECTANT)					
M11	Worker: Is	Frequency of Intervention: every 3 months					
METHOD O	<b>METHOD OF INTERVENTION:</b> check the suction lance of chemical product 2 as described below:						
<ul> <li>Check p chemica 360°) an</li> <li>Check th</li> </ul>	pipe integrity in the al 2. In case the pip ad cut the damaged he suction pipe is n	e point indicated with yellow arrow, without remove the suction lance from container of le is damaged or crushed, remove pipe from the connector (Inox fast coupling, rotating at l component. ot damaged throughout its length.					



### CLEANING OF SAFETY SIGNALS SURFACES

Worker: **Is** Frequency of Intervention: **1 year** 

METHOD OF INTERVENTION:

Clean the safety signals surfaces with water or isopropyl alcohol, using a cloth.



### WARNING

Use a specific product to remove the lime scale. Avoid using products highly corrosive.



### ASSISTANCE

Should your machine not work properly even after ordinary maintenance has been carried out, contact the Technical Support Centre of reference, describing the fault and giving the machine model and serial numbers.



### **15. PROBLEMS - CAUSES – SOLUTIONS**

### **15.1 Introduction**

This chapter includes possible problems which may occur during machine operation, along with their cause and solution. Should the inconveniences continue or take place frequently even after having carried out all the instructions stated in this chapter, please contact the Technical Support Centre of reference.

### 15.2 Problems - Causes - Solutions

### P. MACHINE WILL NOT START:

- C. Circuit breaker de-activated.
- **S.** Place it in the "ON" working position.
- **C.** Machine start switch de-activated.
- S. Press the start button.

### P. UPON GIVING START-UP COMMAND, WASHING CYCLE DOES NOT START:

- C. The door is not correctly closed or locked.
- **S.** Check door closure. Check that the door micro-switch is properly activated.
- C. Micro-switch failure.
- **S.** Check operation and replace as necessary.
- C. No detergent in tank.
- **S.** Turn the machine off and fill the tank.

### P. MACHINE DOES NOT REACH THE SET TEMPERATURE FOR THE SELECTED WASHING CYCLE:

- C. The thermostat probe of the washing chamber is dirty or covered with lime.
- **S.** Clean the thermostat probe of the wash chamber, performing the routine maintenance described in chapter 13.

### P. MACHINE DOES NOT PROPERLY RUN WASHING CYCLE:

- **C.** The nozzles of the washing rotors are clogged by deposits or lime.
- S. Clean the washing arms by carrying out the routine maintenance set forth in chapter 13.
- C. The correct amount of water required for correct washing cycle does not arrive.
- S. Ensure that the water is supplied at the correct pressure and that there are no obstructions.
- C. The correct amount of water required for correct washing cycle does not arrive.
- **R.** Completely close the tap for connection to the plumbing system located upstream from the machine and clean the filter as described in chapter 13.

### P. DETERGENT FILLING PHASE DOES NOT OCCUR CORRECTLY:

- **C.** Chemical dispensing pump not very efficient.
- S. Perform the routine maintenance set forth in chapter 13.
- C. Chemical dispensing pump failed.
- **S.** Contact the Technical Support Centre of reference and ask for the assistance of an **authorized workshop technician** for the repair or replacement of the pump.

### I. MACHINE DOES NOT PERFORM DRYING PHASE:



- **C.** Air filter of drying system is dirty or clogged.
- **R.** Clean the filter by carrying out the routine maintenance set forth in chapter 13.
- **C.** The fan of the drying system does not work.
- **R.** Check the electrical connections of the drying system.
- **R.** Contact the Technical Support Centre of reference and ask for the assistance of an **authorized workshop technician** for the repair or replacement of the motor.



### **16. DECOMMISSIONING**

### 16.1 Instructions for disassembly of the machine

For demolition and subsequent disposal of your machine, proceed as follows:

- Disconnect the machine from the electrical power and water supply, and from the drain. With the machine disconnected, check that the water circuit is not pressurized.
- Contact the organization responsible for reporting and certifying machine demolition, in accordance with the laws in the country where the machine is installed.
- Carry out draining, storage and subsequent disposal of substances such as oils and grease which may be in the lubrication tanks in accordance with the law.
- When disassembling the machine, make sure to divide the materials it is made of according to their chemical makeup (iron, aluminium, bronze, plastic, etc.).
- Ensure that the floor where the machine or any parts of it are placed is made of washable materials, non-absorbent, and provided with adequate drainage to protect against accidental oil leaks or rust. These drains must carry any leakage to watertight collection containers.
- Cover the machine or parts of it with insulating covers to prevent rain or humidity from damaging the structure through oxidation or rust.

Following the legal requirements where the machine is installed and used, dispose of all materials and substances resulting from its disassembly.

### 16.2 Machine disposal



- For the dispose of the equipment get through to the manufacturer or distributor.
- Do not dispose of this equipment as miscellaneous solid municipal waste but arrange to have it collected separately.
- The re-use or correct recycling of the electronic and electrical equipment (EEE) is important in order to protect the environment and the well-being of humans.
- In accordance with European Directive WEEE 2012/19/EC, special collection points are available to which to
  deliver waste electrical and electronic equipment and the equipment can also be handed over to a distributor at
  the moment of purchasing a new equivalent type.
- The public administration and producers of electrical and electronic equipment are involved in facilitating the processes of the re-use and recovery of waste electrical and electronic equipment through the organisation of collection activities and the use of appropriate planning arrangements.
- Unauthorized disposal of waste electrical and electronic equipment is punishable by law with the appropriate penalties.



### **ANNEX A – TRAINING CERTIFICATE**





Certificate N° 569

### **Training Certificate**

We certify that Mr. XXXXXXXX of the company XXXXXXX has attended the training course for EW1 Person in charge for training Riese Pio X, 17/10/2013 D. Aless Steelco S.p.A. Via Balegante, 27 - 31039 - Riese Pio X (TV) - Italy Tel. +39 0423.7561 Fax +39 0423.755528 Website: www.steelcospa.com

"For the Environmentally conscious"



### ANNEX B – WATER SAMPLING DURING THE FINAL RINSING CYCLE

Annex B provides both the instructions and the methods to take the microbiological sample by using the Steelco - Q water BSK (Professional Sampling Kit for Bacterial Check) as well as the relative key for the decoding and the reading of the results obtained.

In any case, please refer to current Legislation, National Guidelines and/or internal protocols.



Q water BSK Professional Sampling Kit for Bacterial Check (code 99911268)

Endoscopes can be re-processed in the cycle dedicated to biological sampling. The only difference compared to a normal disinfection cycle refers to the interruptions during the final rinsing phase.

Water samples are taken during the final rinsing stage following the disinfection stage. The final rinsing cycle consists of two stages.

The sample must be taken during the final rinsing stage.

It is possible to take a water sample from the washing chamber by carrying out the following procedure:

- Set the parameter P3.42 = YES.
- Note: at the end of the water sample-taking cycle, restore the value to NO.

3

09.04.14	20:32	00273
STAND. P	AA B1	
° <b>41.1°</b> ∎ sp°	A0	41.3° AO
press start		

Start the cycle to be tested.

• The machine will stop at the end of every stage and the door lock will be released.



- The START button will be flashing.
- Press the **START** button again to continue the cycle up to the desired point in order to take the water sample from the washing chamber.
- Open the door of the washing chamber and pick up the water sampling from the sump by using the Steelco Q water BSK (Professional Sampling Kit for Bacterial Check). Be sure to perform the sampling in an aseptic way.









Close the door and press the **START** button again to continue the cycle.

It is possible to activate the validation cycle by using the barcode reader, instead of changing the parameters. The modification is valid for a cycle, then it will automatically reset.



	ATTENTION						
	•	Adequate measures must be implemented to avoid contamination of the sample during the water sample-taking stage. It is recommended to use the Steelco - Q water BSK (Professional Sampling Kit for Bacterial Check).					
	•	The water samples for the microbiological test must be taken only during the final rinsing stage. Samples must not be taken during other stages of the process.					
	•	Correct, periodic maintenance of the equipment prevents the risk of contamination. Follow t he manufacturer's instructiologicalgiblersdes wetl as any pre-filters.					
	•	It is important to schedule appropriate thermal disinfection and chemical self-disinfection cycles in accordance with the manufacturer's					



### SAMPLE TAKING

In order to take water samples during the final rinsing stage use the Steelco - Q water BSK (Professional Sampling Kit for Bacterial Check - code 99911268) and follow the instructions below:



### REPROCESSING CYCLE EFFECTIVENESS TEST CARRIED OUT IN THE SUCTION CHANNEL

It is also possible to test the effectiveness of the entire re-conditioning cycle of the endoscope suction channel in the following way:

- Put 20 ml of sterile water in the suction channel (by using a sterile syringe) collecting it from a sterile container situated at the distal part of the tool;
- Clean the same channel using a special cleaning brush by cutting it directly inside the sterile water container. WARNING: the cleaning brush MUST be sterile.
- Remember to close all air and water channel valves before carrying out this procedure.
- Every single container must be accompanied not only by the request form but also by the endoscope type and serial number, the type of channel on which the brushing has been carried out, the name of the operator performing the sampling procedure as well as the report relating to the re-conditioning cycle carried out.

### INTERPRETATION OF THE RESULTS

As for the interpretation of the results, please refer to the following table and criteria.

NUMBER OF COLONIES/100 ml	INTERPRETATION	ACTION		
0	SATISFACTORY			
1-9 (carried out on a regular basis)	ACCEPTABLE	Count under reasonable control		
10-100	UNSATISFACTORY	Investigate by implementing opportune thermal disinfection and chemical self-disinfection cycles		
>100	UNACCEPTABLE	Stop the endoscope washer. Start two thermal disinfection and two chemical self-disinfection cycles. Repeat the test of control and take necessary actions as regards decontamination.		



After interpreting the results, take opportune measures in accordance with the 2 attached tables.

NUMBER OF COLONIES (NOT PSEUDOMONAS)					
TVC LEVE	ïL	ACTION			
SATISFACTORY	< 1cfu/100 ml	No action required.			
ACCEPTABLE	1 – 9 cfu/100 ml	<ul> <li>The personnel must have carried out a self-disinfection cycle in the morning in accordance with the instructions provided by the endoscope washer manufacturer as well as internal protocols.</li> <li>The carrying out of a daily self-disinfection cycle.</li> </ul>			
	10 – 50 cfu/ ml	<ul> <li>The personnel must carry out two thermal disinfection cycles (one after another - 80°C for 10 min) and a special chemical disinfection cycles (35°C for 10 min double chemical).</li> </ul>			
UNSATISFACTORY	51 – 100 cfu/ ml	<ul> <li>The personnel must carry out a special chemical disinfection cycle as well as a thermal disinfection cycle (one after another 80°C for 10 min one after another 55°C for 10 min double chemical).</li> <li>Take another water sample.</li> </ul>			
UNACCEPTABLE	> 100 cfu/100 ml	<ul> <li>STOP USING THE ENDOSCOPE WASHER.</li> <li>The person in charge of the endoscope washer must: <ul> <li>Carry out a special chemical disinfection cycle (55° double chemical for 10 min).</li> <li>Carry out 2 thermal disinfections cycles (one after another to 80°C for 10 min).</li> <li>Take another water sample.</li> </ul> </li> <li>The personnel must: <ul> <li>Do not use the endoscope washer until it has been confirmed that the water sample has a contamination value of &lt; 100 cfu/100 ml.</li> <li>Continue to carry out a self-disinfection cycle on a daily basis.</li> <li>Carry out a special thermal disinfection cycle on a daily basis.</li> </ul> </li> <li>Note: Advice will be obtained from the Lead Doctor for Infection Prevention and Control if there are recurring unacceptable TVC levels.</li> </ul>			



NUMBER OF COLONIES (PSEUDOMONAS)						
TVC LEVE	iL	ACTION				
SATISFACTORY	< 1cfu/100 ml	No action required.				
ACCEPTABLE	1 – 9 cfu/100 ml	<ul> <li>The personnel must have carried out a self-disinfection cycle in the morning in accordance with the instructions provided by the endoscope washer manufacturer as well as internal protocols.</li> <li>The carrying out of a daily self-disinfection cycle.</li> <li>Do not re-process endoscopes until it has been confirmed that the water sample has a contamination value of &lt; 6 cfu/100 ml.</li> </ul>				
	10 – 50 cfu/ ml	<ul> <li>The personnel must carry out two thermal disinfection cycles (one after another - 80°C for 10 min) and a special chemical disinfection cycles (35°C for 10 min double chemical).</li> <li>Do not re-process endoscopes until it has been confirmed that the water sample has a contamination value of &lt; 6 cfu/100 ml. Do not re-process cystoscopes or bronchoscopes until it has been confirmed that water sample has a contamination value of &lt; 10 cfu/100 ml.</li> </ul>				
UNSATISFACTORY	51 – 100 cfu/ ml	<ul> <li>The personnel must carry out two special chemical disinfection cycles (55°C for 10 min double chemical) and two thermal disinfection cycle (one after another 80°C for 10 min).</li> <li>Take another water sample.</li> <li>Do not re-process endoscopes until it has been confirmed that the water sample has a contamination value of &lt; 6 cfu/100 ml. Do not re-process cystoscopes or bronchoscopes until it has been confirmed that water sample has a contamination value of &lt; 10 cfu/100 ml.</li> <li>Note: Advice will be obtained from the Lead Doctor for Infection Prevention and Control if there are recurring unacceptable TVC levels.</li> </ul>				
UNACCEPTABLE	> 100 cfu/100 ml	<ul> <li>STOP USING THE ENDOSCOPE WASHER. The person in charge of the endoscope washer must:</li> <li>Replace all internal filters and clean the filter cases. Carry out 2 special chemical disinfection cycles (55° double chemical for 10 min).</li> <li>Carry out 2 thermal disinfection cycles (one after another to 80°C for 10 min).</li> <li>Take another water sample.</li> <li>The personnel must:</li> <li>Do not use the endoscope washer until it has been confirmed that the water sample has a contamination value of &lt; 100 cfu/100 ml.</li> <li>Carry out a chemical self-disinfection cycle on a daily basis.</li> <li>Carry out a chemical self-disinfection cycle on a daily basis.</li> <li>Do not re-process endoscopes until it has been confirmed that the water sample has a contamination value of &lt; 6 cfu/100 ml. Do not reprocess cystoscopes or bronchoscopes until it has been confirmed that that water sample has a contamination value of &lt; 10 cfu/100 ml.</li> <li>Note: Advice will be obtained from the Lead Doctor for Infection Prevention and Control if there are recurring unacceptable TVC levels.</li> </ul>				



### ANNEX C – TEST PROCEDURE WHEN INSTALLING



### ATTENTION

During the first installation of the equipment upon the site of use, with the aim of eliminating the risk of contamination of the endoscope washer, following maintenance operations that influence critical parts of the endoscope washer (wash pumps, etc.) or after transferring the endoscope washer to another operational centre, the following procedure **MUST** be carried out:

- 3 empty cycles carried out using only detergent;
- 2 complete cycles with detergent + disinfectant;
- 1 thermal disinfection cycle carried out at 80°C.



### ATTENTION

USE ONLY CHEMICAL PRODUCTS THAT HAVE BEEN TESTED AND APPROVED BY THE MANUFACTURER AND IN USE WITH THIS SYSTEM.



# **CLOSED SYSTEM REQUIREMENTS RATIONALE**







National Standard Requirements, Guidelines, Manufaturer Technical Features and Safety Rules: Main reasons why it's important to keep ISO 15883:4 compliance according with international and

- 1. ISO 15883-4 validation
- N Level 2 Script of the wfhss (World Federation for Hospital Sterilization Science) education group Reprocessing Endoscopes
- ω Reprocessing of flexible endoscopes and endoscopic accessories used in gastrointestinal endoscopy: of Gastroenterology Nurses and Associates (ESGENA) – Update 2018 Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society
- 4 Legal implications linked to the Endoscopy Department Risk Management
- 5. Infections & MDRs
- 6. Patient Safety







## ISO/FDIS 15883-4:2018(E) requirements:

Clause 8 Information to be supplied by the manufacturer

4.4.5); n) the detergent(s) and disinfectant(s) type tested with the WD, (see 4.3.3 and



### ISO/FDIS 15883-4:2018(E)

ANNEX A

A.2 Before installation of the WD

It is expected that the WD manufacturer will:

--- state with which endoscopes the WD can be used;

by the WD manufacturer; this is to be done in the light of information provided by the process chemical manufacturer; - provide test data demonstrating the performance of the WD with respect to both cleaning and disinfection using process chemicals specified

**ANNEX 3** 

A.3 Installation and operation of the WD

It is the responsibility of the purchaser/user facility to ensure that the WD is installed and operated correctly.

This includes: The responsibility for ensuring that the WD is correctly installed and functions correctly usually falls to the purchaser/user.

- performance qualification tests;

use of the recommended process chemicals;



ð compliance with ISO 15883:4. **Qualification**) recomended test and it's not be intended PQ (performance enclose as മ full 0

Table C.	2 — Summary	of tests in	addition	to ISO 1	5883-1:2006-	+Amd 1:2014	
Brief description	Requirement	Test	Type	Works	Operational	Performance	Routinea
oftest	subclause	sub- clause	test	test	qualification test	qualification test	test
1. Leak test failure alarm	4.2.3	6.5.3.3	Х	Х	Х	в	X (Q)
2. Leak test non-connection test	4.2.4 4.2.5	<u>6.5.3.4</u>	X	X	В	Х	R(0)X
3. Leak test over-pressurization prevention test	4.2.4 d)	6.5.3.2	X	X	В	в	В
4. Cleaning efficacy	4.3.5	6.11	X	B	X	X	X (Q)
5. in vitro disinfect- ant efficacy	4.4.2.2	6.12.2	Х	В	В	В	B
6. Disinfection effi- cacy — Type test	4.4.2.5	6.12.6.1	X	в	В	в	в
7. Complete process — Operational and performance quali- fication	<u>4.1.3</u>	<u>6.12.6.2</u> 6.12.6.3	0 0	88	B	x v	X (Q)
8. Drying	4.7	6.8	X	0	X	в	Х
9. Disinfection of liquid transport sys- tem — Type test	4.8.5	6.12.5.1	X	в	В	в	в
10. Disinfection of liquid transport sys- tem — Operational qualification and routine test	4.8.5	6.12.5.2	0	В	Х	в	Х
11. Self-disinfection test — Type test	4.8.7	<u>6.12.3.1</u>	X	в	В	в	в
12. Self-disinfection test — Operational qualification and routine test	4.8.7	6.12.3.2	0	в	Х	В	Х
13. Disinfection of water treatment equipment	4.9.2	6.12.4.1	X	в	0	в	в
14. Final rinse water treatment — Micro- bial quality	4.5.2	6.12.4.2	0	в	Х	в	в
15. Channels non-obstruction test	5.2.2.1	6.6	х	Х	в	Х	X(Q)
16. Channels non- connection test	5.2.2.2	<u>6.7</u>	Х	X	В	X	(Q) X
X - Recommended B - Not recommended							
0 – Optional							
Q - Quarterly, W - Wee Frequency of rout!	kly ne testing may be	established b	ased on rish	assessmer	nt using trend an	alysis.	
					c		

		Tabl	e C.2 (con	tinued)			
Brief description	Requirement	Test	Type	Works	Operational	Performance	Routine
oftest	subclause	sub- clause	test	test	qualification test	qualification test	test
17. Temperature	4.4.3	6.9	X	Х	X	X	X (Q)
throughout process	<u>5.4.2</u> 5.4.3						
18. Minimum pro- cess temperature test	5.4.4	<u>6.9.2</u>	X	Х	X	X	X (Q)
19. Water quality	<u>4.5.2</u> <u>4.9.2.3</u> a) <u>4.9.2.3</u> b)	<u>6.3</u>	в	в	Х	Х	X (Q) X (W)
20. Chemical dosing test (single-dose container)	5.7	6.10	Х	Х	Х	В	0
X - Recommended							
D Natura and ad							

Q - Quarterly, W - Weekly - Optional

Frequency of routine testing may be



C
-
•
0
-
Aemb
× '











worldwide harmonization of sterilization departments and of decontamination practices especially by providing: The WFHSS dedicates itself to the promotion of the

- മ sterilization societies, thus stimulating cooperation and the exchange of information and best practices; meeting place for national and regional non-profit
- Ņ information via its website to all our stakeholders interested parties. and



## Endoscopes is stated the following: In the Level 2 Script of the wfhss education group Reprocessing

together, this can give rise to massive chemical reactions (e.g. white deposits on the endoscope). Ask disinfectant) is an important factor for successful reprocessing. If incompatible substances are used the supplier to give a written confirmation of safety and compatibility! Compatibility of the chemical substances used for reprocessing (i.e. compatible detergent and

disinfectors, which sets out the requirements for an EWD as well as its performance criteria. There is a European (and international) standard (EN ISO 15883-4) regulating endoscope washer-

chemistries by the AER manufacturer The above statement shows that the ISO 15883-4 compliance means the use of certified and validated





SOCIETY OF GASTROINTESTINAL ENDOSCOPY REPROCESSING OF FLEXIBLE ENDOSCOPES **POSITION STATEMENT OF THE EUROPEAN** AND ENDOSCOPIC ACCESSORIES USED IN ASSOCIATES (ESGENA) – UPDATE 2018 GASTROENTEROLOGY NURSES AND (ESGE) AND EUROPEAN SOCIETY OF GASTROINTESTINAL ENDOSCOPY:

standard. All Steelco AERs are fully compliant to EN ISO 15883:4

Certified Labs: Certifications have been released by European



### BIOTECH-GERMANDE

### RECOMMENDATION

reeloo

Miele Group

disinfection, in order to: should be the first choice for endoscope cleaning and EWDs compliant with EN ISO 15883 standard series

- Provide a standardized and validated reprocessing cycle in a closed environment;
- Document the process steps automatically (via a printer or electronically);
- Provide reliable and reproducible reprocessing;
- Minimize staff contact with chemicals and contaminated equipment;
- Minimize contamination of the environment;
- Facilitate the work involved for personnel;
- Lower the risk of damage to endoscopes





### LEGAL IMPLICATIONS LINKED TO THE ENDOSCOPY DEPARTMENT RISK MANAGEMENT


## **RATIONALE 4**



Risk Management/Legal issues

health economy Patients claims against physicians and health institutions are increasing and becoming a serious problem in

<u>health plan.</u> Strategies for decreasing these claims and reducing financial losses are an important part of every

should be an integral part of daily routine daily work in gastroenterology units and endoscopy departments Reporting adverse events, complications and fulfilling protocols, guidelines, standards and norms

Being prepared in case of patients potential claims is important for several reasons

- recent protocols, standard, guidelines, norms and are fully traced & documented To prove that all the endoscopes reprocessing procedures have been put in place according with the most
- To prove that all the reprocessing devices have been installed, served and used according to manufacturer instructions (use of approved & validated chemistry is part of it)  $\rightarrow$  IQ, OQ, PQ
- batch number, exipiration date and residual quantity to run cycles demonstrate the safer traceability of the chemistry by the AER by the adoption of RFID labels able to read To prove endoscope compatibility with the chemistry in use, with the AER and connectors and when possible

## **RATIONALE 5**



## **INFECTIONS & MDRS**



### **RATIONAL 5**



## Prevention of multidrug resistant infections from contaminated duodenoscopes **Position Statement of ESGENA**

safety in endoscopy. Noncompliance with guidelines as well as processing protocols how narrow the margin of safety is, despite compliance with redeviations from standardized and validated reprocessing proto-Safe and effective endoscope reprocessing is crucial to patient patient-to-patient transmission. Recent MDRO outbreaks show cols may lead to ineffective reprocessing, with the possibility of



### **RATIONAL 5**





## **RATIONALE 6**



## PATIENT SAFETY



## **RATIONALE 6**



# Patients Safety & Nurses Piece of mind means that:

- Every gastroenterology unit should adopt guality indicators and guality improvement plans for advancing patient satety
- reprocessing, documentation and management of endoscopy units. Endoscopy staff is responsible for individual, comprehensive patient care, technical assistance including

### Why?:

Several studies reported a lack of compliance with established guidelines for endoscope reprocessing.

# An increased risk of infection transmission has been associated with:

- endoscope and AER use
- human error
- inadequate or delayed reprocessing
- failure to sterilize accessory equipment
- incorrect selection of disinfectants
- improper drying



## Why should we take the risk to «swim against the tide»? Question ourselves then ...







Mod. 4606/0



### **CERTIFICATO CE**

Certificato n. 909/MDD Dichiarazione di approvazione del sistema qualità

(Sistema completo di garanzia qualità)

Visto l'esito delle verifiche condotte in conformità all'Allegato II, con l'esclusione del punto 4, della direttiva 93/42/CEE e s.m.i., si dichiara che la ditta:

### **STEELCO SPA**

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

mantiene nello stabilimento di:

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 3 (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 9/A (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 12 (ITA) - Italy

un sistema qualità che assicura la conformità dei seguenti prodotti:

Lavapadelle a termodisinfezione per uso medicale

Termodisinfettori per uso medicale

Lavaendoscopi per uso medicale

Lavastrumenti per decontaminazione ad uso medicale

Lava disinfettatrice-sterilizzatrice chimica per endoscopi

Sterilizzatore al perossido di idrogeno

Soluzione sterilizzante al perossido di idrogeno

Disinfettori termici e/o chimici per uso medicale

Soluzione sterilizzante alla formaldeide

serie e modelli indicati in Allegato

ai requisiti essenziali della direttiva suddetta ad essi applicabili (in tutte le fasi dalla progettazione al controllo finale) ed è sottoposta alla sorveglianza prevista dal punto 5 dell'Allegato II. Per i dispositivi in classe III questo certificato è valido solamente con il relativo certificato di esame CE della progettazione di Allegato II.4.

Emesso il: Data aggiornamento:	2006-03-10	#.Q.f.	
Sostituisce:	2020-07-03		
Dała scadenza:	2024-05-26		
Questa Dichiarazione di a la certificazione CE dei di	pprovazione è sogge sposilivi medici - Marc	ta alle condizioni previste dall'IMQ nel "Regolamento per atura CE - Direttiva 93/42/CEE".	IMQ S.p.A.   I-20138 Milano  Via Quintiliano 43   www.ima.it

Mod. 4606/



### **CERTIFICATO CE**

Certificato n. 909/MDD Dichiarazione di approvazione del sistema qualità

(Sistema completo di garanzia qualità)

Riferimento pratiche IMQ:

DM17-0017730-01; DM17-0016713-01; DM17-0020271-01; DM18-0027447-01; DM18-0030229-01; DM18-0032276-01; DM19-0034875-01; DM18-0033656-01; DM19-0040534-01; DM19-0044781-01; DM20-0048677-01; DM20-0048683-01; DM20-0052515-01; DM20-0052521-01; DM20-0048679-01; DM20-0055248-01.

Questa Dichiarazione di approvazione è rilasciata dall'IMQ S.p.A. quale organismo notificato per la direttiva 93/42/CEE e s.m.i. Il numero identificativo dell'IMQ S.p.A. quale organismo notificato è: 0051.

Emesso il: Data aggiornamento: Sostituisce: Data scadenza:	2006-03-10 2020-09-11 2020-07-03 2024-05-26	Giergi Fulvio IMQ	DocuSign.	
Questa Dichiarazione di c la cerlificazione CE dei di	pprovozione è soggeti spositivi medici - Marca	ta alle condizioni previsto atura CE - Direttiva 93/42	e dal'IMQ nel "Regolamento per /CEE".	IMQ S.p.A.   I-20138 Milano   Via Quintiliano 43   www.ima.it



### **CERTIFICATO CE**

Certificato n. 909/MDD

### Allegato

Lavapadelle a termodisinfezione per uso medicale

Termodisinfettori per uso medicale

Lavaendoscopi per uso medicale

Lavastrumenti per decontaminazione ad uso medicale

Lava disinfettatrice-sterilizzatrice chimica per endoscopi

Sterilizzatore al perossido di idrogeno

Soluzione sterilizzante al perossido di idrogeno

Disinfettori termici e/o chimici per uso medicale

### Soluzione sterilizzante alla formaldeide

Modd. come da documento allegato "ELENCO PRODOTTI" rev. 02 del 02/09/2020, valido solo se provvisto del timbro IMQ; tale allegato costituisce parte integrante e sostanziale del presente certificato.

Emesso il:	2006-03-10	00		
Data aggiornamento:	2020-09-11	A.C	5	
Sostituisce:	2020-07-03	Giorgi Fulvio	Y	
Data scadenza:	2024-05-26	IMQ	DocuSign	
State State				





### EC CERTIFICATE

### Certificate No 909/MDD Full Quality Assurance System Approval Certificate

On the basis of our examination carried out according to Annex II, excluding section 4, of the Directive 93/42/EEC and its revised version, we hereby certify that:

### **STEELCO SPA**

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

manages in the factory of:

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 3 (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 9/A (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 12 (ITA) - Italy

a quality assurance system ensuring the conformity of the following products:

Flusher disinfectors / Bedpan washers

Washer disinfectors

### **Endoscope** washers

**Decontamination clenears** 

### Chemical washer disinfector and sterilizer for endoscopes

Hydrogen peroxide sterilizer

Hydrogen peroxide sterilizing agent

Thermal and/or chemical washer disinfectors for medical purpose

### Formaldehyde sterilizing solution

series and type refs in the Annex

with the relevant essential requirements of the aforementioned directive (from design to final inspection and testing) and it is subject to surveillance as specified in section 5 of Annex II. For class III devices, this certificate is valid only with the relevant EC Design-Examination Certificate of Annex II.4.

Date: Updated: Substitution Date: Expiry Date:	2006-03-10 2020-09-11 2020-07-03 2024-05-26	Giorgi Fulvio IMQ DocuSign	
This Approval Certificate	e is subjected to the provi	sions laid down in the "IMQ regulation for the	IMQ S.p.A.   1-20138 Milano
certification of Medical	Devices - CE Marking - D	rective 93/42/EEC".	Via Quintiliano 43
This is a	a translation of the Italian	text, which prevails in case of doubts	www.imq.it

Mod. 4606/0



### EC CERTIFICATE

### Certificate No 909/MDD Full Quality Assurance System Approval Certificate

Reference to IMQ files Nos:

DM17-0017730-01; DM17-0016713-01; DM17-0020271-01; DM18-0027447-01; DM18-0030229-01; DM18-0032276-01; DM19-0034875-01; DM18-0033656-01; DM19-0040534-01; DM19-0044781-01; DM20-0048677-01; DM20-0048683-01; DM20-0052515-01; DM20-0052521-01; DM20-0048679-01; DM20-0055248-01.

This Approval Certificate is issued by IMQ S.p.A. as Notified Body for the Directive 93/42/EEC and its revised version. Notified Body notified to European Commission under number: 0051.

Date: Updated: Substitution Date: Expiry Date:	2006-03-10 2020-09-11 2020-07-03 2024-05-26	Giorgi Fulvio IMQ DocuSign	
This Approval Certificate	e is subjected to the provisio	ns laid down in the "IMQ regulation for the	IMQ S.p.A.   1-20138 Milano
certification of Medical	Devices - CE Matking - Direc	clive 93/42/EEC".	Via Quintiliano 43
This is	a translation of the Italian tex	t, which prevails in case of doubts	www.imq.it



### EC CERTIFICATE

Certificate No 909/MDD

### Annex

Flusher disinfectors / Bedpan washers

Washer disinfectors

Endoscope washers

**Decontamination clenears** 

### Chemical washer disinfector and sterilizer for endoscopes

Hydrogen peroxide sterilizer

Hydrogen peroxide sterilizing agent

Thermal and/or chemical washer disinfectors for medical purpose

### Formaldehyde sterilizing solution

Type ref. as to attached document "LIST OF PRODUCTS" rev. 02 dated 2020/09/02, valid only if provided with IMQ stamp; this annex is integral and substantial part of this certificate.

Date: Updated: Substitution Date: Expiry Date:	2006-03-10 2020-09-11 2020-07-03 2024-05-26	Giorgi Fulvio IMQ DocuSign	

This is a translation of the Italian text, which prevails in case of doubts

										Revision	e:	02
S	elco	CERT	TIFIC.	AZIONI	ED	I PRODO	)TT(	) – 909/N	IDD	Data:		02/09/2020
				ELE	NCO	) PRODO	ITT			Pagina:	3	1 di 6
Lavapadelle a termodisinfezione per uso medicale della serie BP comprende le seguenti varianti e denominazioni:												
denomina	ALIMENTAZI	ONE PO	OTENZ	Α	P	ORTE		CARR	OZZEF	RIA	AL	LESTIMENTO
Serie: <b>BP</b> Marca:	Monofase 110V / 60Hz 220-240V / 50-6i Trifase 200-230V / 50-6i 380-480V / 50-6i	)Hz I )Hz a )Hz	z Da 750W 1 o 2 MOVIMENTO LAROHEZZA ALTEZZA ALTEZZA ZA ZALTEZZA MANUALE STANDARD / / STANDARD STANDA ZA									STANDARD / E
Nomi comm BP 1 BP 1 BP 1	erciali modelli: 00 M, BP 100 M 00 HA5, BP 100 00 HLA, BP 100	9, BP 100 HE, BP 1 HL2MM	0 A, BP 100 HS 1, BP 1	2 100 A9, I E, BP 100 00 HL2AN	BP 10 HAI M, BI	00 H, BP 100 E, BP 100 H P 100 HXL	) H5, SER,	BP 100 HS, 1 BP 100 HL,	BP 100 BP 100	HA, ) HLS,	Fas I	scicolo tecnico: DT-FT-01_M
Termodis	infettori per 1	iso med	licale	della seri	e D	S comprend	de le	seguenti va	arianti	c denomi	nazi	oni:
	ALIMENTAZIC	NE POT	ENZA		POR	ТЕ	CAI	RROZZERIA	OP	FIONAL		
Serie: DS Marca:	Monofase 110V / 60Hz 220-240V / 50-60 Trifase 200-230V / 50-60 380-480V / 50-60	Hz Da 2 Hz a 37 Hz	2200W 7000W	NUMERO	MC N AU	IANUALE / TOMATICA	s	TANDARD / ALTA	ASCI 3 o 4 I DEI	UGATURA / DOSATORI / PURATA		
DS 55 DS 55 DS 66 DS 66 DS 66 DS 66 DS 66 DS 86 DS 99 DS 99 DS 91	DS 50/2 DRS, DS 50/2 DRSD       DS 50/2 DRSD       DS 50/2 DRSD       DS 50/2 DRSD         DS 50/2 DRS, DS 50/2 DRSD       DS 500 DRSD, DS 500 DRSD, DS 500 DRSD, DS 500 DRS4, DS 500 LED, DS 500 LEDD,       DT-FT-03_M         DS 500 SC, DS 500 SCD, DS 500 SCL, DS 500 SCDL, DS 500 CL, DS 500 CD, DS 600/2 D, DS 600/2 D, DS 600 C, DS 600 CD, DS 600/1 1S, DS 600/2 1S,       DS 610/1, DS 610/1 D, DS 610/2, DS 610/2 D, DS 610/1 1S, DS 610/1 2S, DS 610/2 1S,         DS 610/2 2S, DS 610/1 SL, DS 610/2 D, DS 610/2 SL, DS 610/2 SLD, DS 610/1 SL1S,       DS 610/1 SL2S, DS 610/2 SL1S, DS 610/2 SL2S         DS 650, DS 700, DS 750, DS 750 1S, DS 750 2S, DS 750 3S, DS 800, DS 800 1S, DS 800 2S,       DT-FT-06_M         DS 800 3S,       DS 900, DS 900 1S, DS 900 2S, DS 900 3S, DS 1000, DS 1000 1S, DS 1000 2S, DS 1000 3S,       DT-FT-05_M         DS 900 N, DS 900 N1S, DS 900 N2S, DS 900 N3S, DS 1000 N, DS 1000 N1S, DS 1000 N2S,       DT-FT-05_M									DT-FT-03_M DT-FT-06_M DT-FT-05_M		
Termodis	infettori per 1 ALIMENTAZIO	ISO MED	icale ENZA	della seri	e W	D comprei	ide l	e seguenti v	variant	i e denom	inaz	zioni:
Serie: WD Marca:	Monofase 110V / 60Hz 220-240V / 50-60 Trifase 200-230V / 50-60 380-480V / 50-60	Hz Da 13 a 25 Hz Hz	3500 W 6000W									
Nomi comm WD WD	erciali modelli: 1110 1115										Fas I I	cicolo tecnico: DT-FT-06_M DT-FT-05_M
										6	Ø 2	020-09-11
Redatto da:	PC		V	/erificato	da:	QA		A	pprova	ato da: TD	)	

						Revisione:	02			
517	eelco	CERTIFI	CAZIOI	NE DI PROI	DOTTO – 909/MDD	Data:	02/09/2020			
			EL	ENCO PROI	DOTTI	Pagina:	2 di 6			
Termodis	sinfettori per	orende le seguenti variant	i e denomina	azioni:						
	ALIMENTAZIO	ONE POTENZ	4	STAZIONI NUMERO	$\neg$					
Serie: TW Marca:	ie: TW Trifase rca: 200-230V / 50-60Hz 380-480V / 50-60Hz 380-480V / 50-60Hz Da 3000W a 5									
Nomi comm TW	nerciali modelli: 3000/2, TW 300	0/3, TW 3000/4	4, TW 300	)0/5		F	ascicolo tecnico: DT-FT-09_M			
Townodis	infottori nor u	o modicale d	lla serie	I C comprende	le sequenti varianti e denor	ninazioni				
1 ci mouis	ALIMENTAZIO	NE POTENZA	NUMERO	PORTE	DIMENSIONI LONGITUDINALI					
Serie: LC Marca:	<i>Trifase</i> 200-230V / 50-601 380-480V / 50-601	Hz Hz Hz	1 o 2	MANUALE / AUTOMATICA	/2 = Da 1500 mm a 2999 mm /3 = Da 3000 mm a 3999 mm /4 = Da 4000 mm a 4999 mm					
Nomi com	nerciali modelli:					F	ascicolo tecnico: DT-FT-12 M			
Disinfetto	ri termici e/o c ALIMENTAZIO	ne POTENZA	o medica	PORTE	DIMENSIONI	/arianti e den	ominazioni:			
			МО	VIMENTO	LONGITUDINALI	-				
Marca:	<i>Trifase</i> 200-230V / 50-601 380-480V / 50-601	Hz Hz Hz	M AUT	ANUALE / FOMATICA	/2 = Da 1500 mm a 2999 mm /3 = Da 3000 mm a 3999 mm /4 = Da 4000 mm a 4999 mm					
Nomi comm LC	nerciali modelli: 70/2, LC 70/3, L	C 80/2, LC 80/	3, LC 80/4	4,		F	ascicolo tecnico: DT-FT-08_M			
*Nota: Per	la Serie LC 80 e	LC 70 esistono	due tipoli	gie di disinfezion	ne: termica o ibrida (termica o	chimica).				
I ovostrur	nonti nor dogo	ntominozion	ad uso	medicale con a	usilio di ultrasuoni della se	erie US com	vrende le sequenti			
varianti e o	lenominazioni:	mammazion	au uso	incurcate con a			Jenue le seguena			
	ALIMENTAZI	ONE POTENZ	A	PORTE NUMERO						
Serie: US Marca:	Monofase 110V / 60Hz 220-240V / 50-60 Trifase 200-230V / 50-60 380-480V / 50-60	0Hz Fino a 14000W 0Hz 0Hz		1 o 2						
Nomi comr US 8 US 1 US 3 US 3	nerciali modelli: 30 100, US 200/1, U 300 200, US 1000	S 200/2, US 20	0/3, US 20	00/4, US 100 XL		F	ascicolo tecnico: DT-FT-20_M DT-FT-13_M DT-FT-21_M DT-FT-07_M			
Nomi com US 8 US 1 US 2 US 9	11938e 200-230V / 50-61 380-480V / 50-61 merciali modelli: 30 100, US 200/1, U 300 200, US 1000	S 200/2, US 20	0/3, US 20	00/4, US 100 XL		F	ascicolo tecnico DT-FT-20_M DT-FT-13_M DT-FT-21_M DT-FT-07_M			
						Ċ	2020-09-11			

					Revisione:	02
		CERTIFIC	TTO – 909/MDD	Data:	02/09/2020	
<b>U</b> E	euco		ELENCO PRODOT		0210312020	
			Pagina:	3 di 6		
Lavastrume	nti per decon	itaminazione	ad uso medicale della serie	DC comprende le segu	enti varianti	e denominazioni:
	ALIMENTAL	IONE FOIEN2	NUMERO			
Serie: DC Marca:	Monofase 110V / 60H 220-240V / 50-( Trifase 200-230V / 50-( 380-480V / 50-(	z 50Hz Fino a 10000W 50Hz 50Hz	1 0 2			
Nomi commer DC 100 DC 100	rciali modelli: 0, DC 200/1, 1 00/1, DC 1000	DC 200/2 //2			F	ascicolo tecnico: DT-FT-17_M DT-FT-19_M
Lavaendose	opi per uso n	nedicale della	serie EW 2 comprende le se	guenti varianti e denon	ninazioni:	
	ALIMENTAZI E	ION POTENZA	NUMERO			
Serie: EW 2 Marca:	<i>Trifase</i> 200-230V / 50-60 380-480V / 50-60	0Hz Fino a 9100W	1 0 2			
Nomi commer EW 2/	rciali modelli: 1, EW 2/2, EW	/ 2/1 2S, EW 2/	2 2S, EW 2/1 3S, EW 2/2 3S		F	ascicolo tecnico: DT-FT-10_M
Lava disinf denominazio	ettatrice-ster ni:	ilizzatrice chi	mica per endoscopi della se	erie EW 1 comprende l	e seguenti va	rianti e
	ALIMENTAZI E	ION POTENZA	NUMERO			
Serie: EW 1 Marca:	Monofase 110V / 60Hz 220-240V / 50-60 Trifase 200-230V / 50-60 380-480V / 50-60	0Hz Fino a 6100 W 0Hz 0Hz	Da 1 a 4			
Nomi commer EW 1,	rciali modelli: EW 1 S, EW 1	I DUAL H			F	ascicolo tecnico: DT-FT-18_M
Stavilizzato	o al norossid	o di idrogono	della serie PL comprende le	e sequenți varianți e de	nominazioni.	
Stermzzator	ALIMENTAZ	IONE POTEN	ZA PORTE NUMERO		ionnin 210111.	
Serie: PL Marca:	<i>Trifase</i> 200-230V / 50-0 380-480V / 50-0	60Hz Fino a 5000 W 60Hz	1 o 2			
Nomi commer	rciali modelli: 1, PL 70/2, PL	, 130/1. PL 130	/2. PL 40/1		F	ascicolo tecnico: DT-FT-27_M
					œ	2020-09-11

			Revisione:	02
Nicelco	CERTIFICAZIO	ONE DI PRODOTTO – 909/MDD	Data:	02/09/2020
	E	ELENCO PRODOTTI	Pagina:	4 di 6
Soluzione sterilizzante denominazioni: Serie: SteelcoPro PL Marca:	al perossido di id CONCENTRAZIONE 58% di H2O2	rogeno della serie SteelcoPro PL comp	rende le seș	guenti varianti e
Nomi commerciali modelli: SteelcoPro PL			F	ascicolo tecnico: DT-FT-17_C
Soluzione sterilizzante a Serie: SteelcoPro LTSF Marca:	Ila formaldeide della CONCENTRAZIONE Dal 10% al 20%	a serie SteelcoPro LTSF comprende le segue	enti varianti e	denominazioni:
Nomi commerciali modelli:	StaalcoPro I TSE 6. S	teeleoPro I TSE 8	F	ascicolo tecnico: DT-FT-13 C
			œ	2020-09-11

							Revisione	02		
		CER	TIFICAZI	ONE DI	PRODOTTO	0 – 909/MDD	Data:	02/09/2020		
	Geleo		I	ELENCO	PRODOTTI		Pagina:	5 di 6		
Lavapadelle :	a termodisinfe	ezione po	er uso medica	le della seri	e BP comprende	le seguenti varianti	e denominazio	ni:		
ALIMENTAZIONE POTENZA PORTE CARROZZERIA ALLESTIMENTO										
Serie: <b>PWD</b> Marca:	Monofa 110V / 60 220-240V / 50 Trifase 200-230V / 50	se Hz )-60Hz ? )-60Hz	Da 750W a 6900W	1 0 2	MANUALE / SEMI- AUTOMATICA /	STANDARD / RIDOTTA / MAGGIORATA	STANDARD	STANDARD / E		
AUTOMATICA       MAGGIORATA         Nomi commerciali modelli:       Fascicolo tecnico:         PWD 8541 MD, PWD 8541 AD, PWD 8545 MD, PWD 8545 SAD, PWD 8545 AD, PWD 8546,       DT-FT-01_M         PWD 8546 AD, PWD 8549, PWD 8543       DT-FT-01_M										
<b>Fermodisinfe</b>	ttori per uso 1	nedicale	e della serie co	omprende le	seguenti varianti	e denominazioni:				
	ALIMENTA	ZIONE	POTENZA	P	ORTE	CARROZZER	IA A	LLESTIMENTO		
Serie: PWD	<i>Monofas</i> 110V / 60F	e Iz	Da 2220W	1	MANUALE /	STANDARD		ASCIUGATUIRA		
	220-240V / 50- <i>Trifase</i> 200-230V / 50- 380-480V / 50-	60Hz 60Hz 60Hz	a 37000W	102	ΑυτοΜΑΤΙCA	ALTA		/ DEPURATA		
Nomi commer PWD 8 PWD 6	ciali modelli: 531, PWD 85 5121, PWD 61	31 WS, J 22, PWI	PWD 8532, P D 6121 TH1, J	WD 8532 W PWD 6122 '	VS, PWD 8534, I TH1, PWD 6121	PWD 8534 WS TH2, PWD 6122 7	1 T <b>H2</b>	Fascicolo tecnico: DT-FT-04_M DT-FT-03_M		
						2				
							(i)	2020 00 14		
							4	2020-09-11		

Frequencies         PRODUCT CERTIFICATION - 909/MDD LIST OF PRODUCTS         Date: Da									Revision:	02	
LIST OF PRODUCTS         Page: 1 di 6           Figure place         1 di 6           Figure place         DOOIS         DOOY         Page: 1 di 6           Line: BP         DOOIS         DOOY         SET UP           Tade Mark:         State of the set of the	Ne	260	PROD	TOU	CERT	FICATION	- 909	/MDD	Date:	02/09/2020	
Flusher Disinfectors / Bedpan washers for medical purpose BP line includes the following variants and names:       IF REDING     POWER     DOORS     POOPY       IF IN 6011       Line: BP     Single-phase       Trace phase       Trace phase       Trace phase       1010     101     101     101     101     101     101     101     101       Trace phase       1010     100     <		Page:	1 di 6								
FEEDING         POWER         DOORS         DOOR         BOON         SET UP           Line: BP         Single-phase 1187 (4014)         France phase 200-200 (54-001)         France phase 200-200 (54-001)         France phase 200-200 (54-001)         France phase 200-200 (54-001)         STANDARD 1 or 2         STANDARD 200-200 (54-001)         Technical file: DD STANDARD 200-200 (54-001)         Technical file: DT-FT-01_M         Technical file: DT-FT-01_M         DT-FT-01_M         Technical file: DT-FT-04_M           Line: DS         Single-phase 1180 (401)         From 22000 1 or 2         NUMMER         MOVEMENT         WITH PRUNC 3 or 4 DOSING         Technical file: DT-FT-04_M         DT-FT-04_M         DT-FT-04_M           Trade Mark: 200-200 (54-001)         Single-phase 1180 (401)         To 2         MATUAL AUTOMATIC         WITH PRUNC 3 or 4 DOSING         Technical file: DT-FT-04_M         DT-FT-04_M           DS 500, DS 500 DRS, DS 500 DRSD, DS 500 DRSD, DS 500 DRSD, DS 500 DLSD, S500 CL, DS 500 CL, DS	Flusher Disinfectors / Bedpan washers for medical purpose BP line includes the following variants and names:										
Line: BP Trade Mark:         Single-phase 110V / 00Hz 380-380V / 50-00Hz         NUMBER is 990W         MOVEMENT 1 or 2         VIDTH ANNAL AUTOMATIC         VIDTH WEDUCED INCREASED         STANDARD HEDUCKD INCREASED         STANDARD HEDUCKD INCREASED         STANDARD HEDUCKD HEDUCKD INCREASED         STANDARD HEDUCKD HEDUCKD INCREASED         STANDARD HEDUCKD HEDUCKD INCREASED         STANDARD HEDUCK		FEEDING POWER DOORS BODY									
Line: BP Trade Mark: 200-2007 / Sector 200-2007 / Sec					NUMBER	MOVEMEN	Т	WIDTH	HEIGHT		
Trade Mark:     202-3697 (Sector):     Prom 720W     1 or 2     SEAN- AUTOMATIC     REDUCED     SIANJARD     SIANJARD       Trade manes:     Design / Sector):     0 900W     1 or 2     AUTOMATIC     REDUCED     SIANJARD     Technical file:       DT-FT-01_M     BP 100 HAS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HLS, MP 100 HL	Line: BP	Single-phase 110V / 60Hz	e			MANUAL /	S	TANDARD	CTINDADD	CTANDADD	
Discrete         Three-phase 200-200 / 58-001z         to 9900W         AUTOMATE AUTOMATE         INCREASED         HIGH         E           Trade names: BP 100 N, BF 100 MS, BP 100 A, BP 100 A, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HLS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLS, BP 100 HLA, BP 100 HLAM, BP 100 HL2AM, BP 100 HS, BP 100 HS, BP 100 HLAM, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HLAM, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HA, BP 100 HLA, BP 100 HL2MM, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BP 100 HS, BS 00 DRS, DS 500 DRS, DS 500 DRS, DS 500 LED, DS 500 LS, DS 500 DRS, DS 500 DRS, DS 500 CL, DS 500 LED, DS 500 LED, DS 500 LS, DS 6007 IS, DS 1000 DRS, DS 1000 LS, DS 600 LS, DS 6007 IS, DS 6007 IS, DS 1000 ZS, DS 500 DRS, DS 1000 LS, DS 1000 ZS,	Trade Mark:	220-240V / 50-6	OHz From	750W	1 or 2	SEMI-		REDUCED	STANDARD	STANDARD /	
Trade names:       Trade names:       Technical file:         DS 50, DS 50 D, DS 50 DRS, DS 50 DRSD, DS 50 HDRS, DS 50 HDRSD; DS 500 LEDD, DS 500/2 D, DS 500/2 D, DS 500 DRSD, DS 500 HCA, DS 500/2 D, DS 500/2 D, DS 500/2 D, DS 500 DRSD, DS 500 DRSD, DS 500 HCA, DS 500/2 D, DS 500/2 D, DS 500 DRSD, DS 500 DRSD, DS 500 LEDD, DS 500/2 D, DS 500/2 D, DS 500 DRSD, DS 500 DRSD, DS 500 HCA, DS 500/2 D,	Steelco	Three-phase 200-230V / 50-60 380-480V / 50-60	e to 99 OHz DHz	000W				/ CREASED	HIGH	E	
Washer Disinfectors for medical purpose DS line include the following variants and names:       FEEDING     POWER     NUMBER     MOVEMENT     BODY     OPTIONALS       Inter phase       Single-phase       100 / 6011/       210 - 200 / 1 or 2     MUMBER     MUTH DRIVING 3 or 4 DONING PUMPS       Trade Mark: 210 - 200 / 100 / 50 of DRS, DS 50 DRSD, DS 50 HDRS; DS 50 HDRSD; DS 50/2, DS 50/2 D, DS 500, DS 500 D, DS 500 DRSD, DS 500 DRSD, DS 500 DRSD, DS 500 LED, DS 500 LEDD, DS 500, DS 500 D, DS 500 DRSD, DS 500 DRSD, DS 500 DRSD, DS 500 LL, DS 500 CL, DS 500 DL, DS 500 L, DS 500 DRSD, DS 500 DRSD, DS 500 CL, DS 500 CL, DS 500 CL, SS	Trade names: BP 100 BP 100 BP 100	M, BP 100 M9 HA5, BP 100 H HLA, BP 100 I	, BP 100 A, F IE, BP 100 H <u>1L2MM, BP</u>	BP 100 ISE, B 100 H	A9, BP 10 P 100 HAI IL2AM, BI	00 H, BP 100 H E, BP 100 HSI P 100 HXL	15, BP 10 CR, BP 10	00 HS, BP 00 HL, BP	100 HA, 100 HLS,	Technical file: DT-FT-01_M	
FEEDING         POWER         DOORS         BODY         OPTIONALS           Line: DS Trade Mark:         Single-phase 110V / 601t 20-240V / 50-601t         From 2220W 0: 37000W         Image: MOVEMENT 1 or 2         WITH DRYING AUTOMATIC         WITH DRYING HIGH         3 or 4 DOSING PUMPS WITH WATER SOFTNER           Trade Mark:         20-240V / 50-601t         0: 37000W         1 or 2         MANUAL AUTOMATIC         STANDARD HIGH         WITH WATER SOFTNER           DS 500, DS 500, DS 500 DRS 500 SC DDS 500 DRSD, DS 50 HDRSD, DS 500 LED, DS 500/2, DS 600/2, DS 60/2, DS 750/3, DS 7	Washer Disi	nfectors for m	edical purp	ose D	S line incl	lude the follo	wing var	iants and	names:		
Line: DS Trade Mark:     Single-phase 1101/16102     From 2220W to 37000W     I or 2     MANUAL AUTOMATIC     STANDARD HIGH     WITH DRYING 3 or 4 DOSING PUMPS       Trade Mark:     Three-phase 200-2007/150-6011z     From 2220W     1 or 2     MANUAL AUTOMATIC     STANDARD HIGH     3 or 4 DOSING PUMPS       Trade names:     DS 50, DS 50 DRS 50 DRS, DS 50 DRSD, DS 50 HIDRS; DS 50 HDRSD; DS 507, DS 507, DS 507, DS 507, DS 500 DRSD, DS 600 CD, DS 600 CD, DS 600/1 S, DS 600/1 S, DS 600/2 S, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/1 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 900 NS, DS 900 NS, DS 900 NS, DS 900 NS, DS 1000 NS, DS 1000 NS, DS 1000 NS, DS 1000 NSS, DS 1000 NS		FEEDING	POWEI	2	DO	ORS	BO	DY	OPTIONALS		
Line: DS Trade Mark:         Single-phase 110 Y / 60Hz 200-300 Y / 50-60Hz 300-300 Y / 50-60Hz 300-400 Y / 50-60Hz DS 50, DS 50 D, DS 50 DRS, DS 50 DRSD, DS 50 HDRSD; DS 50 HDRSD; DS 50/2, DS 50/2, DS 50/2, D 500 D, DS 500 D, DS 500 DRS, DS 500 DRSD, DS 500 HDRSD, DS 500 HDRSD; DS 500 LED, D 500 SC, DS 500 D, DS 500 DRS, DS 500 DRSD, DS 500 DRS4, DS 500 LED, DS 500 LED, D 500 SC, DS 500 SC, DS 500 SC, DS 500 DRSD, DS 600 / LS, DS 600 / LS, DS 600/2 IS, DS 600/1, DS 600/1 D, DS 6010/2, DS 600/2 D, DS 600/2 D, DS 600/1 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 500 LED, DS 500 LED, DS 500 LED, DS 610/1, SL, DS 610/2, IS, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL, DS 610/2 SL, DS		and the second sec		NU	UMBER N	MOVEMENT			WITH DRVING		
LILE. DS Trade Mark.         223:200//50:00Hz Threephase         From 2200W 0:37000W         1 or 2         MARVAL AUTOMATIC         STANDARD HIGH         Sort JOSING UNINPS WTTH WATER SOFTNER           Trade Mark.         203:20V/194-00Hz         63:000W         1 or 2         MARVAL AUTOMATIC         STANDARD HIGH         9 or JOSING UNINPS WTTH WATER SOFTNER           Trade mames:         DS 500, DS 500 D, DS 500 D, DS 500 DRS, DS 500 DRSD, DS 500 DRSD, DS 500 DRSD, DS 500 DRSD, DS 500 CDL, DS 500, DS 500 D, DS 500 DRS, DS 500 CLD, DS 500 CLD, DS 500 CLD, DS 500 SC, DS 500 SC, D, DS 500 SC, D, DS 500 SC, DL, DS 500 CLD, DS 500 CLD, DS 500 SC, DS 500 D, DS 500 DRS, DS 500 DRSD, DS 500 DRS4, DS 500 CLD, DS 500 LEDD, DS 500 SC, DS 500 D, DS 500 DRS, DS 500 DRSD, DS 500 CLD, DS 500 CLD, DS 500 CL, DS 600/1 SL, DS 610/2 LS, DS 610/1 SL2, DS 610/2 SL3, DS 610/2 LS, DS 610/1 SL2, DS 610/2 SL3, DS 610/2 SL3, DS 610/1 SL3, DS 610/1 SL2, DS 610/2 SL3, DS 610/2 SL3, DS 610/1 SL3, DS 600, DS 700, DS 750, DS 750 1S, DS 750 2S, DS 750 3S, DS 800, DS 800 1S, DS 800 2S, DS 900 N, DS 900 N1S, DS 900 N2S, DS 900 N3S, DS 1000 N, DS 1000 N1S, DS 1000 N2S, DS 1000 N3S.         DT-FT-06_M DS 900 N3S.           Washers Disinfectors for medical purpose WD line include the following variants and names:         Technical file: DT-FT-06_M DT-FT-06_M DT-FT-06_M           Trade mames:         Single-phase 110V / 001kz         From 13500W         to 25000W           Trade mames:         Verified by: QA         Approved by: TD	Line: DS	Single-phase					1.000 million (1997)				
Sector         200/230 / 50-60Hz 380-869 / 50-60Hz         0 3/000 W         AUTOMATIC         ITCH         WTTT/VATER SOFTAER           Trade names:         DS 50, DS 50 D, DS 50 DRS, DS 50 DRSD, DS 50 HDRS; DS 50 HDRS; DS 500 LED, DS 500 LED, S 500 SC, DS 500 D, DS 600 DRS, DS 600 DRSD, DS 600 DRSD, DS 500 DLED, DS 500 SC, DS 500 D, DS 600 DRS, DS 600 D, DS 600 CD, DS 500 LED, DS 600/1, DS 600/1, DS 600/2, DS 600/2, DS 600/2, DS 600/1 SL, DS 610/1, DS 610/1 D, DS 610/2, DS 600/2, DS 610/2 LD, DS 610/1 SL, DS 610/1, DS 610/1 SL, DS 610/2 SL, DS 610/2 LD, DS 610/1 SL, DS 610/1 SL25, DS 610/2 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL25, DS 610/2 SL, DS 610/2 SL, DS 610/2 SL, DS 610/1 SL25, DS 750, DS 750, DS 750 1S, DS 750 2S, DS 750 3S, DS 800, DS 800 1S, DS 800 2S, DS 900, NS 250, DS 700, DS 750, DS 900 3S, DS 1000, DS 1000 1S, DS 1000 NS, DS 1000 N2S, DS 900, NS 900 N1S, DS 900 N2S, DS 900 N3S, DS 1000, N, DS 1000 N1S, DS 1000 N2S, DS 900 N3S.         DT-FT-06_M           Washers Disinfectors for medical purpose WD line include the following variants and names:         Technical file: DT-FT-06_M           Single-phase 110V / 001k 71acde Mark; 7100 - 001k 7100 -	Trade Mark:	220-240V / 50-60H Three-phase	Iz From 2220	W   1	1 or 2	MANUAL /	STAN	DARD /	3 or 4 DOSING PUMPS		
Trade names:         Technical file:           DS 50, DS 50 D, SS 00 DRS, DS 50 DRSD, DS 50 HDRS; DS 50 HDRS; DS 50 HDRSD; DS 50/2, DS 50/2, DRSD, DS 500 DRS, DS 500 DRS4, DS 500 LED, DS 500 LED, DS 500 LED, DS 500 SCD, DS 500 CL, DS 500 CL, DS 500 CLD, DS 500 CLD, DS 500 CL, DS 500 CLD, DS 600/1 LS, DS 600/1 SL, DS 600/1, D, DS 600/2, DS 600/2, DS 600/2, DS 600/2, DS 600/2, DS 610/2 LS, DS 610/1, D, DS 600/2, DS 610/2 SLD, DS 610/2 S	Reelco	200-230V / 50-60H 380-480V / 50-60H	Iz   10 37000 Iz	~		AUTOMATIC	HI	GH	WITH WATER SOFTNER		
Washers Disinfectors for medical purpose WD line include the following variants and names:         FEEDING       POWER         Line: WD       Single-phase         Trade Mark:       200-230V / 50-60Hz         200-230V / 50-60Hz       From 13500W         to 25000W       to 25000W         Trade names:       Technical file:         WD 1115       DT-FT-06_M         WD 1115       DT-FT-05_M	DS 500 DS 500 DS 600 DS 610 DS 610 DS 610 DS 650 DS 800 DS 900 DS 900 DS 100	DS 50/2 DRS, DS 50/2 DRSD DS 500, DS 500 D, DS 500 DRS, DS 500 DRSD, DS 500 DRS4, DS 500 LED, DS 500 LEDD, DT-FT-03_M DS 500 SC, DS 500 SCD, DS 500 SCL, DS 500 SCDL, DS 500 CL, DS 500 CDL, DS 600/1, DS 600/1 D, DS 600/2, DS 600/2 D, DS 600 C, DS 600 CD, DS 600/1 1S, DS 600/2 1S, DS 610/1, DS 610/1 D, DS 610/2, DS 610/2 D, DS 610/1 1S, DS 610/1 2S, DS 610/2 1S, DS 610/2 2S, DS 610/1 SL, DS 610/1 SLD, DS 610/2 SL, DS 610/2 SLD, DS 610/1 SL1S, DS 610/1 SL2S, DS 610/2 SL1S, DS 610/2 SL2S DS 650, DS 700, DS 750, DS 750 1S, DS 750 2S, DS 750 3S, DS 800, DS 800 1S, DS 800 2S, DS 800 3S DS 900, DS 900 1S, DS 900 2S, DS 900 3S, DS 1000, DS 1000 1S, DS 1000 2S, DS 1000 3S DS 900 N, DS 900 N1S, DS 900 N2S, DS 900 N3S, DS 1000 N, DS 1000 N1S, DS 1000 N2S,								DT-FT-03_M DT-FT-06_M DT-FT-05_M	
Line: WD       Single-phase         110V / 60Hz       From 13500W         200-230V / 50-60Hz       to 25000W         Trade names:       DT-FT-06_M         WD 1115       DT-FT-05_M         WD 1115       Verified by: QA	Washers Dis	sinfectors for r FEEDING	nedical pur	pose V	WD line ii	nclude the fol	lowing v	variants an	nd names:		
Trade names:       Technical file:         WD 1110       DT-FT-06_M         WD 1115       DT-FT-05_M	Line: WD Trade Mark:	Single-phase 110V / 60Hz 220-240V / 50-60H Three-phase 200-230V / 50-60H 380-480V / 50-60H	Iz From 1350 Iz Iz Iz	ow w							
Prepared by: PC Verified by: QA Approved by: TD	Trade names: WD 11 WD 11	10 15								Technical file: DT-FT-06_M DT-FT-05_M	
Prepared by: PC Verified by: QA Approved by: TD		(W) 2020-09-11									
	Prepared by:	: PC		Verif	fied by: Q	A		Арр	roved by: TD		

						Revision:	02
NE	260	PRODUC	ICT CERTIFICATION - 909/MDD			Date:	02/09/2020
LIST OF PRODUCTS				'S	Page:	2 di 6	
Vasher Disi	nfectors for med	ical purpose	TW line	includes the foll	owing variants and na	mes:	
FEEDING POWER STATIONS NUMBER							
Line: TW Trade Mark:	<i>Three-phase</i> 200-230V / 50-60Hz 380-480V / 50-60Hz	From 3000W to 140000W					
Trade names:         Technical file:           TW 3000/2, TW 3000/3, TW 3000/4, TW 3000/5         DT-FT-09_N							
hermal Wa	sher Disinfector	s for medica	l purpose	e LC line include	es the following variar	its and name	s:
	FEEDING	POWER	I NUMBER	DOORS MOVEMENT	LONGITUDINAL DIMENSIONS		
ine: LC Trade Mark:	<i>Three-phase</i> 200-230V /50-60Hz 380-480V / 50-60Hz	From 1200W to 201000W	1200W         1 or 2         MANUAL         /2 = From 1500 mm to 2999 mm           AUTOMATIC         /4 = From 4000 mm to 4999 mm				
rade names: LC 20/	1, LC 20/2.						Technical file: DT-FT-12_M
Thermal an	d/or Chemical V	Washer Disii	ifectors f	or medical pur	pose* LC line includ	les the follo	wing variants and
umes.	FEEDING	G POWER DOORS LONGITUDINAL					
Line: LC Trade Mark:	<i>Three-phase</i> 200-230V /50-60Hz 380-480V / 50-60Hz	From 1200W to 201000W	MANUAL         /2 = From 1500 mm to 29           /         /3 = From 3000 mm to 39           AUTOMATIC         /4 = From 4000 mm to 49			9 mm 9 mm 9 mm	
Trade names: LC 70/2, LC 70/3, LC 80/2, LC 80/3, LC 80/4, LC 80/2 BOT LC 80/3 BOT LC 80/4 BOT							Fechnical file: DT-FT-08_M
Note: For the his distinctio	LC 80 and LC 70 n is guaranteed and	Series there ar recognizable	e two type by the finis	s of disinfection: t shed product code	hermal or hybrid (therma shown on the label.	al or chemical	l).
Decontamin	ation cleaners fo	or medical p	urpose w	ith ultrasonic f	unction US line inclu	des the follo	owing variants and
ames.	FEEDING	POWER	]	DOORS			
Line: US Frade Mark:	Single-phase 110V / 60Hz 220-240V / 50-60Hz Three-phase 200-230V /50-60Hz 380-480V / 50-60Hz	Up to 14000W		1 or 2			
Frade names: US 80 US 100, US 200/1, US 200/2, US 200/3, US 200/4, US 100 XL US 300 US 900, US 1000						;	Fechnical file: DT-FT-20_M DT-FT-13_M DT-FT-21_M DT-FT-07_M
						Ŵ	2020-09-11

				Revision:	02	
	PRODU	CT CERTIFICATIO	N - 909/MDD	Date:	02/09/2020	
		LIST OF PRODUC	TS	Page:	3 di 6	
ation cleaners	for medical p	urpose DC line includes	he following variants a	and names:		
FEEDING	POWER	DOORS	]			
Single-phase 110V / 60Hz 220-240V / 50-601 Three-phase 200-230V /50-601 380-480V / 50-601	Hz Up to 10000W Iz Hz	1 or 2				
Trade names:         Technical file:           DC 100, DC 200/1, DC 200/2         DT-FT-17_N           DC 1000/1, DC 1000/2         DT-FT-19_N						
ashers for me	dical nurnose	EW 2 line includes the fo	lowing variants and r	ames.		
FEEDING	POWER	DOORS NUMBER				
<i>Three-phase</i> 200-230V /50-60H 380-480V / 50-60H	z Up to z 9100W	1 or 2				
, EW 2/2, EW 2	2/1 2S, EW 2/2	2S, EW 2/1 3S, EW 2/2 3S		Т	echnical file: DT-FT-10_M	
sher disinfect	or and steriliz	er for endoscopes EW 1	line includes the follow	wing variants	and names:	
FEEDING	POWER	DOORS NUMBER				
Single-phase 110V / 60Hz 220-240V / 50-60H Three-phase 200-230V /50-60Hz 380-480V / 50-60Hz	z Up to 6100 W z	From 1 to 4				
EW 1 S, EW 1 1	DUAL H			Т	echnical file: DT-FT-18_M	
voxide steriliz	ver PL line inc	ludes the following varian	ts and names.			
FEEDING	POWER	DOORS NUMBER	]			
<i>Three-phase</i> 200-230V / 50-60 380-480V / 50-60	Hz Up to 5000 W	1 or 2				
Trade names:         Technical file:           PL 70/1, PL 70/2, PL 130/1, PL 130/2, PL 40/1         DT-FT-27_N					echnical file: DT-FT-27_M	
A. The second seco						
	tion cleaners FEEDING Single-phase 110V / 60Hz 220-240V / 50-60H Three-phase 200-230V / 50-60H 380-480V / 50-60H 380-480V / 50-60H (Contemportane) (C	PRODUtition cleaners for medical pFEEDINGPOWERSingle-phase 110V / 60Hz 2200-230V / 50-60HzUp to 10000W200-230V / 50-60Hz 380-480V / 50-60HzUp to 10000Wop to 1000/2POWER'ashers for medical purposeFEEDINGPOWERThree-phase 200-230V / 50-60HzUp to 9100W'ashers for medical purposeFEEDINGPOWERThree-phase 200-230V / 50-60HzUp to 9100W, EW 2/2, EW 2/1 2S, EW 2/2sher disinfector and sterilizFEEDINGPOWERSingle-phase 110V / 60HzUp to 6100 W200-230V / 50-60HzUp to 6100 WSingle-phase 110V / 60HzUp to 6100 WW 1 S, EW 1 DUAL Hroxide sterilizer PL line inc FEEDINGPOWERW 1 S, EW 1 DUAL Hroxide sterilizer PL line inc 5000 WFEEDINGPOWERJ80-480V / 50-60HzUp to 5000 W	PRODUCT CERTIFICATION LIST OF PRODUC         Automatical purpose DC line includes to FEEDING       POWER       DOORS         Single-phase 110V / 60Hz       Up to 1 or 2       DOORS         200-230V /50-60Hz       Up to 10000W       1 or 2         ,DC 200/1, DC 200/2 0/1, DC 1000/2       DOORS       NUMBER         ,DC 200/1, DC 200/2 0/1, DC 1000/2       Up to 9100W       1 or 2         ,ashers for medical purpose EW 2 line includes the for NUMBER       Three-phase 9100W       1 or 2         ,BC 200/1, DC 200/2 0/1, DC 1000/2       Up to 9100W       1 or 2         ,ashers for medical purpose EW 2 line includes the for NUMBER       NUMBER         Three-phase 100-230V /50-60Hz       Up to 9100W       1 or 2         ,EW 2/2, EW 2/1 2S, EW 2/2 2S, EW 2/1 3S, EW 2/2 3S       Sher disinfector and sterilizer for endoscopes EW 1         FEEDING       POWER       DOORS NUMBER         Single-phase 110V / 60Hz       Up to 6100 W       From 1 to 4         Swa480V / 50-60Hz       Up to 6100 W       From 1 to 4         SW 1 S, EW 1 DUAL H       Foxide sterilizer PL line includes the following variant FEEDING       POWER       DOORS         W1 S, EW 1 DUAL H       Su0-480V / 50-60Hz       Up to 5000 W       1 or 2         Su-480V / 50-60Hz       Up to 5000 W       1 or 2       1	PRODUCT CERTIFICATION - 909/MDD LIST OF PRODUCTS         stingle-phase 100 / 60Hz 230-340V / 50-60Hz         UP to Single-phase 100 / 60Hz 230-340V / 50-60Hz         200 / 2	PRODUCT CERTIFICATION - 909/MDD LIST OF PRODUCTS         Revision: Date: Page:           ition cleaners for medical purpose DC line includes the following variants and names:         FEEDING         POWER         DOORS           Single-phase 100/v60th 202-30v /50-60th         Up to 1000W         1 or 2         T           abelase         100/v60th         1 or 2         T           abelase         1000W         1 or 2         T           abelase         POWER         DOORS         T           abelase         FEEDING         POWER         DOORS           Three-phase         Up to 9100W         1 or 2         T           precephase         0100 W         From 1 to 4         T           precephase         0100 W         From 1 to 4 </td	

			Revision:	02
Neelco	PRODUCT	Date:	02/09/2020	
		Page:	4 di 6	
Hydrogen peroxide ster	ilizing agent Steelcol	Pro PL line includes the following varia	ants and names:	
Line: SteelcoPro PL Trade Mark:	58% of H2O2			
Trade names: SteelcoPro PL			Т	echnical file: DT-FT-17_C
Formaldahuda atavilizin	ng agant StaalaaDra I	TSE line includes the following varian	its and names:	
Formaldenyde sterilizin	CONCENTRATION	SF mie mendes die fonowing varian	ts and names.	
Line: SteelcoPro LTSF				
Trade Mark:	From 10% to 20%			
Trade names:		4	T	Technical file:
SteelcoPro LTSF 4	; SteelcoPro LISF 6; S	teelcopto LISF 8		D1-11-13_C
6				
				2020-09-11

.

Steelco		PRODUCT CERTIFICATION - 909/MDD					Revision:	02
							Date:	02/09/2020
			LIST OF PRODUCTS					5 di 6
Flusher Disi	nfectors / Bo	edpan v	vashers for	medical p	urpose PWD li	ine includes the fo	llowing varia	nts and names:
	FEEDIN	G	POWER	D	OORS	BODY	IEIOTII	SET UP
Line: PWD Trade Mark:	Monofas 110V / 601 220-240V / 50 Trifase 200-230V / 50 380-480V / 50	e Hz -60Hz -60Hz -60Hz	Da 750W a 6900W	1 o 2	MOVEMENT MANUAL / SEMI- AUTOMATIC / AUTOMATIC	WIDTH STANDARD / REDUCED / INCREASED	STANDARD / HIGH	STANDARD / E Technical file:
PWD 8	8541 MD, PW 8546 AD, PW	D 8541 D 8549.	AD, PWD 8 PWD 8543	545 MD, PV	WD 8545 SAD, I	PWD 8545 AD, PW	D 8546,	DT-FT-01_M
		2 00179						
Washer Disi	infectors for	medica	al purpose l	PWD line	include the follo	owing variants and	names:	
	FEEDIN	١G	POWER	I NUMBER	DOORS MOVEMENT	BODY	1	SET UP
Line: PWD Trade Mark:	Single-ph 110V / 60 220-240V / 50 Three-ph 200-230V / 50 380-480V / 50	nase DHz D-60Hz D-60Hz D-60Hz D-60Hz	From 2220W to 37000W	1 OR 2	MANUAL / AUTOMATIC	STANDA / HIGH	STANDARD / HIGH	
Trade name: PWD 8531, PWD 8531 WS, PWD 8532, PWD 8532 WS, PWD 8534, PWD 8534 WS PWD 6121, PWD 6122, PWD 6121 TH1, PWD 6122 TH1, PWD 6121 TH2, PWD 6122 T				TH2	Technical file: DT-FT-04_M DT-FT-03_M			
							6	2020-09-11



### **CERTIFICATO CE**

Certificato n. 1879/MDD

### Dichiarazione di approvazione del sistema qualità

(Sistema completo di garanzia qualità)

Visto l'esito delle verifiche condotte in conformità all'Allegato II, con l'esclusione del punto 4, della direttiva 93/42/CEE e s.m.i., si dichiara che la ditta:

### **STEELCO SPA**

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

mantiene nello stabilimento di:

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 12 (ITA) - Italy

un sistema qualità che assicura la conformità dei seguenti prodotti:

Soluzione disinfettante per dispositivi medici invasivi

Mod. SteelcoXide-A, SteelcoXide-B e SteelcoXide-DT Marca STEELCO

ai requisiti essenziali della direttiva suddetta ad essi applicabili (in tutte le fasi dalla progettazione al controllo finale) ed è sottoposta alla sorveglianza prevista dal punto 5 dell'Allegato II. Per i dispositivi in classe III questo certificato è valido solamente con il relativo certificato di esame CE della progettazione di Allegato II.4.

Riferimento pratiche IMQ:

DM16A0650339-01; DM17-0007338-01; DM19-0040534-01.

Questa Dichiarazione di approvazione è rilasciata dall'IMQ S.p.A. quale organismo notificato per la direttiva 93/42/CEE e s.m.i. Il numero identificativo dell'IMQ S.p.A. quale organismo notificato è: 0051.

Emesso il:	2016-05-16	¥.Q.	
Data aggiornamento:	2019-09-18	DetuSion	
Sostituisce:	2017-02-02	IMQ	
Data scadenza:	2024-05-26		
Questa Dichiarazione di a la certificazione CE dei di	pprovazione è soggett spositivi medici - Marco	ta alle condizioni previste dall'IMQ nel "Regolamento atura CE - Direttiva 93/42/CEE".	per IMQ S.p.A.   I-20138 Milano   Via Quintiliano 43   www.ima.it



### **EC CERTIFICATE**

Certificate No 1879/MDD

Full Quality Assurance System Approval Certificate

On the basis of our examination carried out according to Annex II, excluding section 4, of the Directive 93/42/EEC and its revised version, we hereby certify that:

### **STEELCO SPA**

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

manages in the factory of:

31039 RIESE PIO X (TV) - VIA BALEGANTE 27 (ITA) - Italy

31039 RIESE PIO X (TV) - VIA DEL LAVORO 12 (ITA) - Italy

a quality assurance system ensuring the conformity of the following products:

### Disinfectant solution for invasive medical devices

Type ref. SteelcoXide-A, SteelcoXide-B e SteelcoXide-DT Trade mark STEELCO

with the relevant essential requirements of the aforementioned directive (from design to final inspection and testing) and it is subject to surveillance as specified in section 5 of Annex II. For class III devices, this certificate is valid only with the relevant EC Design-Examination Certificate of Annex II.4.

Reference to IMQ files Nos:

DM16A0650339-01; DM17-0007338-01; DM19-0040534-01.

This Approval Certificate is issued by IMQ S.p.A. as Notified Body for the Directive 93/42/EEC and its revised version. Notified Body notified to European Commission under number: 0051.



This is a translation of the Italian text, which prevails in case of doubts



Tuesday, May 19, 2020

### COMPATIBILITY DECLARATION LETTER

To whom it may concern,

Based on the declaration letter which issued by Fabio Zardini, President & CEO of STEELCO SPA on 15 May 2020.

We undersigned, Aohua Endoscopy Co., Ltd., with registered office in Aohua Endoscopy Building, No.66, Lane133, Guangzhong Road, Minhang District, Shanghai, China, 201100, hereby declare following compatibility between Aohua endoscope and Steelco systems:

VGT-Q30J Gastroscope can be properly connected with Steelco systems EW1-EW1
 S-EW2-EW2 2S-EW2 3 S-ED200 S- EPW100 S by using Steelco connectors:

Biopsy channel code 660538 Air channel code 660128 Channels separator code 660088 Suction channel code 660528 Air Water channel code 660116 Auxiliary Water Jet channel code 660089 Leak test original AOHUA connector has to be used

- VGT-Q30J Gastroscope can be properly connected with Steelco Drying Storage Cabinets models ED100-ED150-ED200-ED250 by using Steelco connector kit: 99911466
- VRL-1T30 Bronchoscope can be properly connected with Steelco systems EW1-EW1
   S-EW2-EW2 2S-EW2 3 S-ED200 S- EPW100 S by using Steelco connectors:

Biopsy channel code 660538 Suction channel code 660510 *Leak test original AOHUA connector has to be used* 

 VRL-1T30 Bronchoscope can be properly connected with Steelco Drying Storage Cabinets models ED100-ED150-ED200-ED250 by using Steelco connector kit:9991658

Yours sincerely

Jiang Suping 1933 200.5.19

Service Development Manager Aohua Endoscopy Co., Ltd.



### DICHIARAZIONE UE DI CONFORMITÁ EU DECLARATION OF CONFORMITY DÉCLARATION UE DE CONFORMITÉ EU KONFORMITÄTS-ERKLÄRUNG DECLARACIÓN UE DE CONFORMIDAD

### Il sottoscritto, come legale rappresentante della azienda sotto indicata, dichiara che il prodotto:

The undersigned, officer of the under-written company, hereby declares that the product:

Le représentant juridique soussigné de l'usine sous indiquée, il déclare que le produit:

Der Unterzeichner, Handlungsbevollmächtigter des oben genannten unter hingewiesen, erklärt hiermit, daß das Produkt:

El firmante, como representante legal de la empresa indicada, declara que el producto:

Nome/Modello: Name/Type: Nom/Modèle: Name/Model: Nombre/Modelo:



N° di Serie/Lotto: Serial/Lot N.: N° de Série/Lot : Serial N./ Reihe-Zahl:

N° de Serie/Lote:

**19XXXXXXXXXX** 

[GMDN: 44835]

in riferimento alla Direttiva 93/42/CEE e s.m.i., è classificato in classe IIb, in accordo alla regola 15 dell'allegato IX, è stato progettato e costruito in conformità ai requisiti essenziali dell'allegato I, applicando le disposizioni delle norme armonizzate.

La presente dichiarazione è redatta sulla base dei requisiti dell'allegato II della direttiva 93/42/CEE e s.m.i..

La persona giuridica autorizzata a costituire il fascicolo tecnico è la Steelco S.p.A. con sede in via Balegante, 27 - Riese Pio X (TV) - Italia.

Il sistema di garanzia di qualità del prodotto, esclusivamente in accordo all'allegato II della Direttiva 93/42/CEE e s.m.i., è mantenuto sotto controllo dall'organismo notificato IMQ S.p.A con numero identificativo 0051 come da certificato IMQ n.1879/MDD con validità fino al 26/05/2024.

referring to 93/42/EEC Medical Device Directive and s.m.i., is classified on IIb class, according to rule 15 of the annex IX, designed and manufactured in conformity with the annex I, under the harmonized rules. This declaration is written in conformity with the annex II of 93/42/EEC Medical Device Directive and s.m.i..

The juridical person authorized to compile the technical file is Steelco S.p.A., at via Balegante, 27 - Riese Pio X (TV) - Italy.

Exclusively in conformity at annex II of 93/42/EEC Medical Device Directive and s.m.i., the product quality system is guaranteed from the notified authority IMQ S.p.A, under number 0051 as IMQ certificate n. 1879/MDD valid till 2024/05/26.

en référence de la Directive 93/42/CEE et s.m.i., est classifiée en classe IIb, en accord à la règle 15 annexe IX, a été projetée et construit en conformité aux qualités essentielles de l'annexe I, en appliquant les dispositions des normes harmonisées.

Cette déclaration est rédigée sur la base des qualités de l'annexe II de la Directive 93/42/CEE et s.m.i..

La personne juridique autorisée à constituer le dossier technique est la Steelco S.p.A. avec siège en via Balegante, 27 - Riese Pio X (TV) - Italie. Le système de garantie de qualité du produit, exclusivement en accord à l'annexe II de la Directive 93/42/CEE et s.m.i., est maintenu sous contrôle de l'organisme déclaré IMQ S.p.A. avec numéro identificateur CE0051 selon le certificat IMQ n.1879/MDD avec validité jusqu'à 26/05/2024.

Welche gemäß der Richtlinie 93/42/EWG und zusätzliche Änderungen und Ergänzungen, als Medizinprodukt der Klasse IIb klassifiziert ist, konform zur Regel 15 der Anhang IX, wurde geplant und hergestellt, gemäß der wesentlichen Anforderungen laut beiliegender Anlage I, unter der harmonisierten Normen.

Diese Erklärung bestätigt auch die Konformität zu Anhang II der Richtlinie 93/42/EWG und zusätzliche Änderungen und Ergänzungen. Die juristische autorisierte Person, die die technischen Unterlagen zusammenzustellen hat ist Steelco S.p.A. an via Balegante, 27 - Riese Pio X (TV) - Italien.

Ausschließlich konform zu Anhang II der Richtlinie 93/42/EWG über Medizinprodukte, hat der oben genannten Hersteller ein Qualitätsmanagementsystem eingeführt, wie durch die benannte Stelle IMQ S.p.A. mit der Nummer 0051, durch das Zertifikat mit der Nummer IMQ 1879/MDD, gültig bis zum 26/05/2024, bestätigt wird..

en referencia a la Directiva 93/42/CEE y s.m.i., es clasificado en clase IIb, en acuerdo a la regla 15 del anexo IX, ha sido diseñado y construida en conformidad con los requisitos esenciales del anexo I, aplicando las disposiciones de las normas armonizadas.

La presente declaración es redactada sobre la base de los requisitos del anexo II de la norma 93/42/CEE y s.m.i..

La persona jurídica lícita a constituir el expediente técnico es la Steelco S.p.A. con sede en via Balegante, 27 - Riese Pio X (TV) - Italia. El sistema de garantía de calidad del producto, exclusivamente en acuerdo al anexo II a la Directiva 93/42/CEE y s.m.i., es controlado por el organismo notificado IMQ S.p.A. con número identificativo 0051 como de certificado IMQ n. 1879/MDD con validez hasta a 26/05/2024.

Direttive applicate: 93/42/EEC (Medical Devices Directive and s.m.i. - 2007/47/EC) Applied directives: Directives appliquées: Angewandte Richtlinien: Directivas aplicadas:

**Direttore Generale** Managing Director Director Général Geschäftsführer Gerente

Fabio Zardini Steelco S.p.A

STEELCO S.p.A.

MQ204-12 Rev.03

Via Balegante, 27 31039 Riese Pio X (TV) ITALIA – ITALY – ITALIE - ITALIEN

Tel. +39 0423 7561 info@steelcogroup.com www.steelcogroup.com

Fax +39 0423 755528



### DICHIARAZIONE UE DI CONFORMITÁ EU DECLARATION OF CONFORMITY DÉCLARATION UE DE CONFORMITÉ EU KONFORMITÄTS-ERKLÄRUNG DECLARACIÓN UE DE CONFORMIDAD

### Il sottoscritto, come legale rappresentante della azienda sotto indicata, dichiara che il prodotto:

The undersigned, officer of the under-written company, hereby declares that the product:

Le représentant juridique soussigné de l'usine sous indiquée, il déclare que le produit:

Der Unterzeichner, Handlungsbevollmächtigter des oben genannten unter hingewiesen, erklärt hiermit, daß das Produkt: El firmante, como representante legal de la empresa indicada, declara que el producto:

Nome/Modello: Name/Type: Nom/Modèle:

Name/Model:

Nombre/Modelo:

SteelcoXide-B

N° di Serie/Lotto: Serial/Lot N.: N° de Série/Lot : Serial N./ Reihe-Zahl:

N° de Serie/Lote:

**19XXXXXXXXXX** 

[GMDN: 44835]

in riferimento alla Direttiva 93/42/CEE e s.m.i., è classificato in classe IIb, in accordo alla regola 15 dell'allegato IX, è stato progettato e costruito in conformità ai requisiti essenziali dell'allegato I, applicando le disposizioni delle norme armonizzate.

La presente dichiarazione è redatta sulla base dei requisiti dell'allegato II della direttiva 93/42/CEE e s.m.i..

La persona giuridica autorizzata a costituire il fascicolo tecnico è la Steelco S.p.A. con sede in via Balegante, 27 - Riese Pio X (TV) - Italia.

Il sistema di garanzia di qualità del prodotto, esclusivamente in accordo all'allegato II della Direttiva 93/42/CEE e s.m.i., è mantenuto sotto controllo dall'organismo notificato IMQ S.p.A con numero identificativo 0051 come da certificato IMQ n.1879/MDD con validità fino al 26/05/2024.

referring to 93/42/EEC Medical Device Directive and s.m.i., is classified on IIb class, according to rule 15 of the annex IX, designed and manufactured in conformity with the annex I, under the harmonized rules. This declaration is written in conformity with the annex II of 93/42/EEC Medical Device Directive and s.m.i..

The juridical person authorized to compile the technical file is Steelco S.p.A., at via Balegante, 27 - Riese Pio X (TV) - Italy.

Exclusively in conformity at annex II of 93/42/EEC Medical Device Directive and s.m.i., the product quality system is guaranteed from the notified authority IMQ S.p.A, under number 0051 as IMQ certificate n. 1879/MDD valid till 2024/05/26.

en référence de la Directive 93/42/CEE et s.m.i., est classifiée en classe IIb, en accord à la règle 15 annexe IX, a été projetée et construit en conformité aux qualités essentielles de l'annexe I, en appliquant les dispositions des normes harmonisées.

Cette déclaration est rédigée sur la base des qualités de l'annexe II de la Directive 93/42/CEE et s.m.i..

La personne juridique autorisée à constituer le dossier technique est la Steelco S.p.A. avec siège en via Balegante, 27 - Riese Pio X (TV) - Italie. Le système de garantie de qualité du produit, exclusivement en accord à l'annexe II de la Directive 93/42/CEE et s.m.i., est maintenu sous contrôle de l'organisme déclaré IMQ S.p.A. avec numéro identificateur CE0051 selon le certificat IMQ n.1879/MDD avec validité jusqu'à 26/05/2024.

Welche gemäß der Richtlinie 93/42/EWG und zusätzliche Änderungen und Ergänzungen, als Medizinprodukt der Klasse IIb klassifiziert ist, konform zur Regel 15 der Anhang IX, wurde geplant und hergestellt, gemäß der wesentlichen Anforderungen laut beiliegender Anlage I, unter der harmonisierten Normen.

Diese Erklärung bestätigt auch die Konformität zu Anhang II der Richtlinie 93/42/EWG und zusätzliche Änderungen und Ergänzungen. Die juristische autorisierte Person, die die technischen Unterlagen zusammenzustellen hat ist Steelco S.p.A. an via Balegante, 27 - Riese Pio X (TV) - Italien.

Ausschließlich konform zu Anhang II der Richtlinie 93/42/EWG über Medizinprodukte, hat der oben genannten Hersteller ein Qualitätsmanagementsystem eingeführt, wie durch die benannte Stelle IMQ S.p.A. mit der Nummer 0051, durch das Zertifikat mit der Nummer IMQ 1879/MDD, gültig bis zum 26/05/2024, bestätigt wird..

en referencia a la Directiva 93/42/CEE y s.m.i., es clasificado en clase IIb, en acuerdo a la regla 15 del anexo IX, ha sido diseñado y construida en conformidad con los requisitos esenciales del anexo I, aplicando las disposiciones de las normas armonizadas.

La presente declaración es redactada sobre la base de los requisitos del anexo II de la norma 93/42/CEE y s.m.i..

La persona jurídica lícita a constituir el expediente técnico es la Steelco S.p.A. con sede en via Balegante, 27 - Riese Pio X (TV) - Italia. El sistema de garantía de calidad del producto, exclusivamente en acuerdo al anexo II a la Directiva 93/42/CEE y s.m.i., es controlado por el organismo notificado IMQ S.p.A. con número identificativo 0051 como de certificado IMQ n. 1879/MDD con validez hasta a 26/05/2024.

Direttive applicate: 93/42/EEC (Medical Devices Directive and s.m.i. - 2007/47/EC) Applied directives: Directives appliquées: Angewandte Richtlinien: Directivas aplicadas:

**Direttore Generale** Managing Director Director Général Geschäftsführer Gerente

Fabio Zardini Steelco S.p.A

STEELCO S.p.A.

MQ204-12 Rev.03

Via Balegante, 27 31039 Riese Pio X (TV) ITALIA – ITALY – ITALIE - ITALIEN

Tel. +39 0423 7561 info@steelcogroup.com www.steelcogroup.com

Fax +39 0423 755528



### DICHIARAZIONE UE DI CONFORMITÁ EU DECLARATION OF CONFORMITY DÉCLARATION UE DE CONFORMITÉ EU KONFORMITÄTS-ERKLÄRUNG DECLARACIÓN UE DE CONFORMIDAD

### Il sottoscritto, come legale rappresentante della azienda sotto indicata, dichiara che il prodotto:

The undersigned, officer of the under-written company, hereby declares that the product:

Le représentant juridique soussigné de l'usine sous indiquée, il déclare que le produit:

Der Unterzeichner, Handlungsbevollmächtigter des oben genannten unter hingewiesen, erklärt hiermit, daß das Produkt: El firmante, como representante legal de la empresa indicada, declara que el producto:

Nome/Modello: Name/Type: Nom/Modèle: Name/Model:

Nombre/Modelo:



N° di Serie/Lotto: Serial/Lot N.: N° de Série/Lot : Serial N./ Reihe-Zahl: N° de Serie/Lote:

**19XXXXXXXXXX** 

[GMDN: 44835]

in riferimento alla Direttiva 93/42/CEE e s.m.i., è classificato in classe IIb, in accordo alla regola 15 dell'allegato IX, è stato progettato e costruito in conformità ai requisiti essenziali dell'allegato I, applicando le disposizioni delle norme armonizzate.

La presente dichiarazione è redatta sulla base dei requisiti dell'allegato II della direttiva 93/42/CEE e s.m.i..

La persona giuridica autorizzata a costituire il fascicolo tecnico è la Steelco S.p.A. con sede in via Balegante, 27 - Riese Pio X (TV) - Italia.

Il sistema di garanzia di qualità del prodotto, esclusivamente in accordo all'allegato II della Direttiva 93/42/CEE e s.m.i., è mantenuto sotto controllo dall'organismo notificato IMQ S.p.A con numero identificativo 0051 come da certificato IMQ n.1879/MDD con validità fino al 26/05/2024.

referring to 93/42/EEC Medical Device Directive and s.m.i., is classified on IIb class, according to rule 15 of the annex IX, designed and manufactured in conformity with the annex I, under the harmonized rules. This declaration is written in conformity with the annex II of 93/42/EEC Medical Device Directive and s.m.i..

The juridical person authorized to compile the technical file is Steelco S.p.A., at via Balegante, 27 - Riese Pio X (TV) - Italy.

Exclusively in conformity at annex II of 93/42/EEC Medical Device Directive and s.m.i., the product quality system is guaranteed from the notified authority IMQ S.p.A, under number 0051 as IMQ certificate n. 1879/MDD valid till 2024/05/26.

en référence de la Directive 93/42/CEE et s.m.i., est classifiée en classe IIb, en accord à la règle 15 annexe IX, a été projetée et construit en conformité aux qualités essentielles de l'annexe I, en appliquant les dispositions des normes harmonisées.

Cette déclaration est rédigée sur la base des qualités de l'annexe II de la Directive 93/42/CEE et s.m.i..

La personne juridique autorisée à constituer le dossier technique est la Steelco S.p.A. avec siège en via Balegante, 27 - Riese Pio X (TV) - Italie. Le système de garantie de qualité du produit, exclusivement en accord à l'annexe II de la Directive 93/42/CEE et s.m.i., est maintenu sous contrôle de l'organisme déclaré IMQ S.p.A. avec numéro identificateur CE0051 selon le certificat IMQ n.1879/MDD avec validité jusqu'à 26/05/2024.

Welche gemäß der Richtlinie 93/42/EWG und zusätzliche Änderungen und Ergänzungen, als Medizinprodukt der Klasse IIb klassifiziert ist, konform zur Regel 15 der Anhang IX, wurde geplant und hergestellt, gemäß der wesentlichen Anforderungen laut beiliegender Anlage I, unter der harmonisierten Normen.

Diese Erklärung bestätigt auch die Konformität zu Anhang II der Richtlinie 93/42/EWG und zusätzliche Änderungen und Ergänzungen. Die juristische autorisierte Person, die die technischen Unterlagen zusammenzustellen hat ist Steelco S.p.A. an via Balegante, 27 - Riese Pio X (TV) - Italien.

Ausschließlich konform zu Anhang II der Richtlinie 93/42/EWG über Medizinprodukte, hat der oben genannten Hersteller ein Qualitätsmanagementsystem eingeführt, wie durch die benannte Stelle IMQ S.p.A. mit der Nummer 0051, durch das Zertifikat mit der Nummer IMQ 1879/MDD, gültig bis zum 26/05/2024, bestätigt wird..

en referencia a la Directiva 93/42/CEE y s.m.i., es clasificado en clase IIb, en acuerdo a la regla 15 del anexo IX, ha sido diseñado y construida en conformidad con los requisitos esenciales del anexo I, aplicando las disposiciones de las normas armonizadas.

La presente declaración es redactada sobre la base de los requisitos del anexo II de la norma 93/42/CEE y s.m.i..

La persona jurídica lícita a constituir el expediente técnico es la Steelco S.p.A. con sede en via Balegante, 27 - Riese Pio X (TV) - Italia. El sistema de garantía de calidad del producto, exclusivamente en acuerdo al anexo II a la Directiva 93/42/CEE y s.m.i., es controlado por el organismo notificado IMQ S.p.A. con número identificativo 0051 como de certificado IMQ n. 1879/MDD con validez hasta a 26/05/2024.

Direttive applicate: 93/42/EEC (Medical Devices Directive and s.m.i. - 2007/47/EC) Applied directives: Directives appliquées: Angewandte Richtlinien: Directivas aplicadas:

**Direttore Generale** Managing Director Director Général Geschäftsführer Gerente

Fabio Zardini Steelco S.p.A

STEELCO S.p.A.

MQ204-12 Rev.03

Via Balegante, 27 31039 Riese Pio X (TV) ITALIA – ITALY – ITALIE - ITALIEN

Tel. +39 0423 7561 info@steelcogroup.com www.steelcogroup.com

Fax +39 0423 755528

### FUJIFILM

### FUJIFILM ITALIA S.p.A.

S.S. n° 11 Padana Superiore, 2/B 20063 Cernusco sul Naviglio (MI) - Itały Tel. 02.92974.1 - Fax 02.92974.591 http://www.fujifilm.it - info@fujifilm.it

### TO WHOM IT MAY CONCERN

Cernusco Sul Naviglio, June 23rd, 2016 Prot. n° **33761**/2016/GV/dm

The undersigned Giovanni Ettore Valtorta, born in Sovico (MB) – Italy – on January  $31^{st}$ , 1959 and living in Sovico (MB) – Italy – Via San Francesco d'Assisi, 8 in his capacity as Special Attorney of FUJIFILM Italia S.p.A., having its register office in Cernusco sul Naviglio (MI) – Italy – S.S. n° 11 Padana Superiore, 2/b – Fiscal Code 09435590154, VAT n° 11025740157, under his responsibility,

### DECLARES

That FUJIFILM endoscopes are fully compatible with Steelco AERs series EW.

Your successfully.

FUJIFILM Italia S.p.A. A Special Attorney Giovanni Ettore Valtorta









FUJIFILM Europe GmbH

Heesenstraße 31, 40549 Düsseldorf, PO Box 270131, 40524 Düsseldorf, Germany Phone +49 211 5089-100, Fax +49 211 5089-344 http://www.fujifilm.eu

Willich, 23.07.2018

To Whom it may concern

Herewith we confirm, that based on the test results of the company Steelco and the submitted report from the 20.06.2018, the use of the FUJIFILM series ED-580XT in reprocessing machines made by Steelco, is basically permitted and no problems are to be expected.

n

Haruhiko Arai Manager Quality & Service Support



### AZIENDA OSPEDALIERA UNIVERSITARIA INTEGRATA VERONA



(D.Lgs. n. 517/1999 - Art. 3 L.R.Veneto n. 18/2009)

### U.O. Clinica di Gastroenterologia - Direttore Prof. Luca Frulloni U.O. Semplice Organizzata di Endoscopia Digestiva Responsabile: Dr. Armando Gabbrielli

Policlinico G.B. Rossi, P.le L.A. Scuro, 10 - 37134 Verona - Tel. 045 8124743 - Fax 045 8124898 e-mail: endoscopia.digestiva.op@azosp.vr.it

Date 22.6.2016

We are glad to communicate that, in our GI Departement - AOUI of Verona - Steelco S.p.A

AERs Systems (EW) installed in 2012 are succesfully used with Olympus endoscopes series

140, 160 and 180.

Sincerly

Armando Gabbrielli, M.D.

IA INTEGRATA DI VERONIA E POLICLINICO 'BORGO ROMA' di Gastroenterologia zzativa Semplice DIGESTIVA GABBRIELLI



UNITÀ OPERATIVA CON SISTEMA QUALITÀ UNI EN ISO 9001:2008 - Certificato n. 194114

Sede Legale Azienda Ospedaliera Universitaria Integrata: P.le A. Stefani, 1 - 37126 VERONA - Tel 045/812 1111 - Fax 045/916735 C.F. c P. Iva 03901420236 - Portale Aziendale: <u>www.ospedaliverona.it</u>



Acqui Terme, June 23th, 2016

We are glad to comunicate that, in our Ambulatory Department identified as "Endoscopia Digestiva" located in Casa di Cura Villa Igea S.p.A., the Firm SO.E.M. sold and installed in 2014 and 2015 two Systems EW1 of Steelco production. Both systems are successfully used with Olympus endoscopes series 180 and 190.

Best regards

The Legal Representative Casa di Cura Villa Igea" S.p.A. Il Consigliere Delegato Poggio Dott. Giovanni

The Bioingeenering Dept. Maurizio La Gamba

### AZIENDA OSPEDALIERA DI VERONA - OSPEDALE CIVILE MAGGIORE SERVIZIO DI CHIRURGIA ENDOSCOPICA Responsabile: Dott. Luca Rodella

Verona, June 23, 2016

We are glad to communicate that, in our GI Department – Piastra Endoscopica / Endoscopy Dept at Polo Chirurgico "P. Confortini" - Maggiore Hospital of Verona - Steelco S.p.A AERs Systems (EW) installed in 2011 are succesfully used with Olympus endoscopes series 180 and 190.

Sincerely

Luca Rodella Head Endoscopy Unit

InoRolle


Servizio Sanitario Nazionale - Regione Veneto <u>Azienda UNITA' LOCALE SOCIO-SANITARIA N. 13</u> DIPARTIMENTO INTERNISTICO OSPEDALE DI DOLO **U.O.S. di GASTROENTEROLOGIA ed ENDOSCOPIA DIGESTIVA** Resp. **Dr. Renato Marin** 

Dolo, 13.07.2016

#### U.O.S. di Gastroenterologia ed Endoscopia Digestiva

Dr. Renato Marin (Responsabile)

Medici

Dr.ssa Erica Cervellin Dr. ssa Sara Antoniazzi Dr. ssa Adriana Sergio

### Infermieri

I.P. Mara Celin (con ruolo di coordinamento) I.P. Silvia Bizzotto I.P. Federica Carlesso I.P. Federico Cipriotto I.P. Samuele Marcato I.P. Lisa Morelli I.P. Grazia Salmaso I.P. Lisa Scantamburlo

0.S.S.

Sig.ra Elisabetta Zara

## Telefoni

Sala A / Segreteria 041 5133423 Sala B (tel. e fax) 041 5133460

e-mail:

endoscopiadigestiva.dolo@ulss13mirano.ven.it

## Prenotazione prestazioni:

CUP: tel. 041 5103520

#### Prestazioni erogate

- EGDS e Colonscopie (diagnostiche ed operative, routinarie ed in Pronta Disponibilità) -Endoscopia con Videocapsula

- Screening Tumore Colon-Retto

- Visite Gastroenterologiche

We are glad to comunicate that, in our GI department -Gastroenterology and Endoscopy Department U. O. S. D. Dolo Hospital - Steelco S.p.A. AER System (model EW2) installed in 2009 are successfully used with Olympus endoscopes series 140 - 160 - 180 - 190.

Sincerely

Renato Marin Head Endoscopy Unit.

Regione Veneto-AXIENDA U.L.S.S. N. 13 U.O.S. DIPARTIMENTALEA/GASTROENTEROLOGIA ed ENDOSCOPIA DIGESTIVA Resp.: Dott. RENATO MARIN C.F. MRN RNP 52R25 F241H VE 4808



# **NHS** Foundation Trust

Lowestoft Road Gorleston Great Yarmouth Norfolk NR31 6LA

Main Switchboard: 01493 452452

Direct Dial: 01493 453051

E mail: judy.dron@jpaget.nhs.uk

Twitter: @jamespagetNHS

04/07/2016

To whom it may concern-

We use the Steelco AER to reprocess Olympus 240 and 260 series endoscopes. Currently we have had no issues with the process, or problems identified by our endoscope service provider that could be attributed to the Steelco. We have been very satisfied with the machine.

Yours sincerely

J Dron Senior Sister Endoscopy suite James Paget University Hospital

	STEELCO S.P.A.	Revision nr. 3 Dated 18/05/2021
Miele Group Member	Steelcoxide B	Printed on 29/06/2021 Page n. 1/16 Replaced revision:2 (Printed on: 05/06/2019)
SECTION 1. Identification of the 1.1. Product identifier Product name Chemical name and synonym INDEX number EC number	Safety Data Sheet According to Annex II to REACH - Regulation 20 substance/mixture and of the co Steelcoxide B Hydrogen Peroxide Solution 008-003-00-9 231-765-0	<sup>115/830</sup> mpany/undertaking
CAS number	7722-84-1	
Registration Number	01-2119485845-22-XXXX	
1.2. Relevant identified uses of the substar	nce or mixture and uses advised against	
Intended use	Solution "B" of the tri-component syst to be used in combination with the act	em "SteelcoXide". Concentrated sterilizing solution ivator "Solution A". PROFESSIONAL USE ONLY.
1.3. Details of the supplier of the safety dat	a sheet	
Name Full address District and Country	STEELCO S.p.A. Via Balegante, 27 31039 Riese Pio X (TV) ITALY	
	tel. +39 0423 7561	
	fax +39 0423 755528	
e-mail address of the competent person		
responsible for the Safety Data Sheet	info@steelcogroup.com	
Product distribution by:	STEELCO S.p.A.	
1.4. Emergency telephone number		
For urgent inquiries refer to		
	Centros de Orientação de Doentes U	rgentes (CODU): 800250250
SECTION 2. Hazards identificati	on	

## 2.1. Classification of the substance or mixture

The product is classified as hazardous pursuant to the provisions set forth in (EC) Regulation 1272/2008 (CLP) (and subsequent amendments and supplements). The product thus requires a safety datasheet that complies with the provisions of (EU) Regulation 2015/830. Any additional information concerning the risks for health and/or the environment are given in sections 11 and 12 of this sheet.

Hazard classification and indication:

	STEELCO S.P.A	۸.	Revision nr. 3 Dated 18/05/2021
Miele Group	er		Printed on 29/06/2021
	Steelcoxide B		Page n. 2/16 Replaced revision:2 (Printed on: 05/06/2019)
Oxidising solid, category 2 Acute toxicity, category 4 Skin Irritation, category 2 Specific target organ toxicity – single exposu Serious eye damage, category 1 Hazardous to the aquatic environment, chron category 3	H272 H302 H315 re, category 3 H335 H318 hic toxicity, H412	May intensify f Harmful if swal Causes skin irr May cause res Causes seriou Harmful to aqu	ire; oxidiser. lowed. itation piratory irritation s eye damage. atic life with long lasting effects.
	Classification note/note	es according to A	Annex VI to the CLP Regulation: B
2.2. Label elements			
Hazard labelling pursuant to EC Regulation 12	72/2008 (CLP) and subsequent amendme	ents and suppler	nents.
Hazard pictograms:			
	!>		
Signal words: Danger			
Hazard statements:			
H272May intensify fireH302Harmful if swallowH315Causes skin irritaH318Causes serious eH335May cause respinH412Harmful to aquati	; oxidiser. ved. tion ye damage. atory irritation c life with long lasting effects.		
Precautionary statements:			
P305+P351+P338IF IN EYES: Rins rinsing.P280Wear protective of P261P261Avoid breathing v P210P210Keep away from P310Contains:Hydrogen Peroxi INDEX	e cautiously with water for several minute: loves/eye protection/face protection. apours/spray. clothing and other combustible materials. heat, hot surfaces, sparks, open flames ar a POISON CENTER. de Solution 35 ≤wt% <47	s. Remove contained other ignition	act lenses, if present and easy to do. Continue sources. No smoking.
INDEX 008-00	3-00-9		
2.3. Other hazards			
On the basis of available data, the product doe	es not contain any PBT or vPvB in percent	age ≥ than 0,1%	).
SECTION 3. Composition/info	mation on ingredients		
3.1. Substances			
Contains:			
Identification X	= Conc. % Classification 1272/200	8 (CLP)	

	Group	STEELCO S.P.A.	Revision nr. 3 Dated 18/05/2021
	Member	Steelcoxide B	Printed on 29/06/2021 Page n. 3/16 Replaced revision:2 (Printed on: 05/06/2019)
HYDROGEN PEROXIDE SOLUTION			
CAS 7722-84-1 EC 231-765-0 INDEX 008-003-00-9 Reg. no. 01-2119485845-22-xxxx	35 ≤ :	$x < 47$ Ox. Liq. 1 H271, Acute Tox. 4 H302, Acu Eye Dam. 1 H318, Aquatic Chronic 3 H4 according to Annex VI to the CLP RegulaSpecific concentration limit (Annex VI Le STOT SE 3; H335; C $\geq$ 35 % Eye Dam. 1; H318: 8 % $\leq$ C $<$ 50 % Eye Irrit. 2; H319: 5 % $\leq$ C $<$ 50 % Ox. Liq. 1; H271: C $\geq$ 70 % Ox. Liq. 2; H272: 50 % $\leq$ C $<$ 70 % Skin Corr. 1A; H314: C $\geq$ 70 % Skin Irrit. 2; H315: 35 % $\leq$ C $<$ 50 %	ite Tox. 4 H332, Skin Corr. 1A H314, 12, Classification note/notes ation: B gislation 1272/2008)

The full wording of hazard (H) phrases is given in section 16 of the sheet.

## **SECTION 4. First aid measures**

## 4.1. Description of first aid measures

EYES: Remove contact lenses, if present. Wash immediately with plenty of water for at least 30-60 minutes, opening the eyelids fully. Get medical advice/attention.

SKIN: Remove contaminated clothing. Rinse skin with a shower immediately. Get medical advice/attention.

INGESTION: Have the subject drink as much water as possible. Get medical advice/attention. Do not induce vomiting unless explicitly authorised by a doctor.

INHALATION: Get medical advice/attention immediately. Remove victim to fresh air, away from the accident scene. If the subject stops breathing, administer artificial respiration. Take suitable precautions for rescue workers.

PROTECTION MEASURES FOR FIRST AID: for the PPE necessary for first aid interventions refer to section 8.2 of this safety data sheet.

#### 4.2. Most important symptoms and effects, both acute and delayed

EYES: from painful irritation to severe chemical burns (conjunctival hyperemia, conjunctivitis, edema, blepharospasm, iritis, corneal turbidity, epithelial defects, permanent corneal damage).

SKIN: solutions starting from approx. 10%: temporary bleaching (oxygen emphysema in interstitial tissue, gas embolism in blood capillaries), erythema, pain; from a concentration of approx. 70%: chemical burns, comparable to third degree burns; contamination of large areas could lead to systemic oxygen embolism.

INHALATION: mucosal irritation, inflammatory tissue reactions, obstruction, glottic and pulmonary edema, dyspnoea up to respiratory failure; extreme cases could lead to systemic effects.

INGESTION: irritation up to corrosion of the mucous membranes in contact especially in the upper digestive tract, distension of the stomach, displacement of the upper respiratory tract due to foaming, gastritis, duodenitis, colitis, acute visceral congestion, formation of vacuoles in the submucosal gastrointestinal tract, in lymphatic channels, mesenteric lymph nodes or mucosal-associated lymphoid tissue, as well as vacuolation in other organs, systemic effects due to air embolism.

SYSTEMIC EFFECTS: shock, acute coronary insufficiency, status epilepticus, cerebrovascular collapse, respiratory failure.

The most frequent cause of death after ingestion of solutions W> 10% is the obstruction of the respiratory tract due to the formation of foam (-> mechanical asphyxia).

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

In the event of an accident or discomfort, consult a doctor immediately (if possible show the instructions for use or the safety data sheet).

## **SECTION 5. Firefighting measures**

5.1. Extinguishing media



Group Member

Printed on 29/06/2021

Page n. 4/16 Replaced revision:2 (Printed on: 05/06/2019)

#### SUITABLE EXTINGUISHING EQUIPMENT

The extinguishing equipment should be of the conventional kind: carbon dioxide, foam, powder and water spray. UNSUITABLE EXTINGUISHING EQUIPMENT None in particular.

#### 5.2. Special hazards arising from the substance or mixture

HAZARDS CAUSED BY EXPOSURE IN THE EVENT OF FIRE Do not breathe combustion products. By decomposition it releases oxygen. The release of oxygen can promote combustion.

## 5.3. Advice for firefighters

#### GENERAL INFORMATION

Use jets of water to cool the containers to prevent product decomposition and the development of substances potentially hazardous for health. Always wear full fire prevention gear. Collect extinguishing water to prevent it from draining into the sewer system. Dispose of contaminated water used for extinction and the remains of the fire according to applicable regulations.

#### SPECIAL PROTECTIVE EQUIPMENT FOR FIRE-FIGHTERS

Normal fire fighting clothing i.e. fire kit (BS EN 469), gloves (BS EN 659) and boots (HO specification A29 and A30) in combination with self-contained open circuit positive pressure compressed air breathing apparatus (BS EN 137).

## **SECTION 6.** Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1 For those who do not intervene directly

Do not take any action involving any personal risk or without proper training. Evacuate the surrounding areas. Do not touch or walk on the spilled material.

Wear suitable protective equipment (including personal protective equipment referred to in section 8 of this Safety Data Sheet) to prevent contamination of skin, eyes and personal clothing. Wear appropriate respirator when ventilation is inadequate.

Do not inhale the mists / vapors / fumes. Avoid the dispersion of the product into the environment. Follow the appropriate internal procedures for personnel not authorized to intervene directly in the event of accidental release.

## 6.1.2 For those who intervene directly

Stop the leak if there is no danger.

Evacuate unauthorized personnel. Wear suitable protective equipment. (see section 8 of this Safety Data Sheet). Follow the appropriate internal procedures for authorized personnel. Isolate the danger area and deny entry. Ventilate enclosed spaces before entering.

## 6.2. Environmental precautions

The product must not penetrate into the sewer system or come into contact with surface water or ground water.

#### 6.3. Methods and material for containment and cleaning up

Collect the leaked product into a suitable container. Evaluate the compatibility of the container to be used, by checking section 10. Absorb the remainder with inert absorbent material.

Make sure the leakage site is well aired. Contaminated material should be disposed of in compliance with the provisions set forth in point 13.

#### 6.4. Reference to other sections

Any information on personal protection and disposal is given in sections 8 and 13.

## **SECTION 7. Handling and storage**

## 7.1. Precautions for safe handling

Ensure that there is an adequate earthing system for the equipment and personnel. Avoid contact with eyes and skin. Do not breathe powders, vapours

Declar	Group	STEELCO S.P.A.	Revision nr. 3 Dated 18/05/2021
Cicerco	Member	Staalaavida P	Printed on 29/06/2021
		SteelCoxide B	Page n. 5/16 Replaced revision:2 (Printed on: 05/06/2019)

or mists. Do not eat, drink or smoke during use. Wash hands after use. Avoid leakage of the product into the environment.

#### 7.2. Conditions for safe storage, including any incompatibilities

Store only in the original container. Store in a ventilated and dry place, far away from sources of ignition. Keep containers well sealed. Keep the product in clearly labelled containers. Avoid overheating. Avoid violent blows. Keep containers away from any incompatible materials, see section 10 for details.

## 7.3. Specific end use(s)

No use other than that indicated in section 1.2 of this safety data sheet

## **SECTION 8. Exposure controls/personal protection**

#### 8.1. Control parameters

Regulatory References:

DEU	Deutschland	Technischen Regeln für Gefahrstoffe (TRGS 900) - Liste der Arbeitsplatzgrenzwerte und Kurzzeitwerte.
		MAK- und BAT-Werte-Liste 2020, Ständige Senatskommission zur Prüfung gesundheitsschädlicher
		Arbeitsstoffe, Mitteilung 56
ESP	España	Límites de exposición profesional para agentes químicos en España 2019
FRA	France	Valeurs limites d'exposition professionnelle aux agents chimiques en France. ED 984 - INRS
GBR	United Kingdom	EH40/2005 Workplace exposure limits (Fourth Edition 2020)
	TLV-ACGIH	ACGIH 2020

## HYDROGEN PEROXIDE SOLUTION

Threshold Limit Value								
Туре	Country	TWA/8h		STEL/15min		Remarks Observa	; / tions	
		mg/m3	ppm	mg/m3	ppm			
MAK	DEU	0,71	0,5	0,71	0,5			
VLA	ESP	1,4	1					
VLEP	FRA	1,5	1					
WEL	GBR	1,4	1	2,8	2			
TLV-ACGIH		1,4	1					
Predicted no-effect concentra	ation - PNEC							
Normal value in fresh water				12,6	μg	/L		
Normal value in marine water	ſ			12,6	μg	/L		
Normal value for fresh water	sediment			47	hà	/kg		
Normal value for marine wate	er sediment			47	μg	/kg		
Normal value for water, interr	nittent release			19,8	hà	/L		
Normal value of STP microor	ganisms			4,66	mç	g/l		
Normal value for the terrestria	al compartment			2,3	hà	/kg		
Health - Derived no-effe	ct level - DNEL / Effects on consumers	DMEL			Effects on workers			
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic systemic	Acute local	Acute systemic	Chronic local	Chronic systemic
Oral					VND	VND	VND	VND
Inhalation	1,93 mg/m3	VND	210 µg/m³	VND	3 mg/m3	VND	1,4 mg/m3	VND
Skin					VND	VND	VND	VND

Legend:



## Steelcoxide B

Page n. 6/16 Replaced revision:2 (Printed on: 05/06/2019)

(C) = CEILING ; INHAL = Inhalable Fraction ; RESP = Respirable Fraction ; THORA = Thoracic Fraction.

VND = hazard identified but no DNEL/PNEC available ; NEA = no exposure expected ; NPI = no hazard identified.

### 8.2. Exposure controls

As the use of adequate technical equipment must always take priority over personal protective equipment, make sure that the workplace is well aired through effective local aspiration.

When choosing personal protective equipment, ask your chemical substance supplier for advice. Personal protective equipment must be CE marked, showing that it complies with applicable standards. Provide an emergency shower with face and eye wash station.

#### HAND PROTECTION

Protect hands with category III work gloves (see standard EN 374).

The following should be considered when choosing work glove material: compatibility, degradation, failure time and permeability. The work gloves' resistance to chemical agents should be checked before use, as it can be unpredictable. The gloves' wear time depends on the duration and type of use.

#### SKIN PROTECTION

Wear category II professional long-sleeved overalls and safety footwear (see Regulation 2016/425 and standard EN ISO 20344). Wash body with soap and water after removing protective clothing.

#### EYE PROTECTION

Wear airtight protective goggles (see standard EN 166).

In the presence of risks of exposure to splashes or squirts during work, adequate mouth, nose and eye protection should be used to prevent accidental absorption.

#### RESPIRATORY PROTECTION

If the threshold value (e.g. TLV-TWA) is exceeded for the substance or one of the substances present in the product, use a mask with a type B filter whose class (1, 2 or 3) must be chosen according to the limit of use concentration. (see standard EN 14387). In the presence of gases or vapours of various kinds and/or gases or vapours containing particulate (aerosol sprays, fumes, mists, etc.) combined filters are required.

Respiratory protection devices must be used if the technical measures adopted are not suitable for restricting the worker's exposure to the threshold values considered. The protection provided by masks is in any case limited.

If the substance considered is odourless or its olfactory threshold is higher than the corresponding TLV-TWA and in the case of an emergency, wear open-circuit compressed air breathing apparatus (in compliance with standard EN 137) or external air-intake breathing apparatus (in compliance with standard EN 138). For a correct choice of respiratory protection device, see standard EN 529.

ENVIRONMENTAL EXPOSURE CONTROLS

The emissions generated by manufacturing processes, including those generated by ventilation equipment, should be checked to ensure compliance with environmental standards.

Product residues must not be indiscriminately disposed of with waste water or by dumping in waterways.

## **SECTION 9.** Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Appearance	liquid	
Colour	transparent	
Odour	pungent	
Odour threshold	Not determined	
pH	3	
Melting point / freezing point Initial boiling point Boiling range	-33 °C 108 °C Not determined	Concentration:30% Concentration:35%
Flash point	> 60 °C	

	STEELCO	S.P.A.	Revision nr. 3 Dated 18/05/2021
Member	Steelcox	ide B	Printed on 29/06/2021 Page n. 7/16 Replaced revision:2 (Printed on: 05/06/2019)
Evaporation rate Flammability (solid, gas) Lower inflammability limit	Not determined Sulla base dello stato fisico Not applicable		
Upper inflammability limit Lower explosive limit Upper explosive limit	Not applicable Not applicable Not applicable		
Vapour pressure Vapour density	13,33 20 °C 1.2	Concentration:35% Temperature:20 °C	
Relative density Solubility	1,05 soluble in water		
Partition coefficient: n-octanol/water	-1,57 Log Pow	Concentration:60% Temperature:20 °C	
Decomposition temperature	>60 °C		
Viscosity	1,249 mPa s	Concentration:60% Temperature:20 °C	
Explosive properties	Non esplosivo		
Oxidising properties	Ossidante con cat.2		
9.2. Other information			
Molecular weight	24,400		

# **SECTION 10. Stability and reactivity**

## 10.1. Reactivity

Decomposes if exposed to: light,heat. Decomposes on contact with: alkaline metals.Possibility of explosion. The product is a strong oxidizing and reactive agent.

## 10.2. Chemical stability

In contact with impurities, decomposition catalysts, incompatible substances, and combustible substances, the substance can give rise to a selfaccelerated exothermic decomposition reaction with formation of oxygen.

## 10.3. Possibility of hazardous reactions

The product may react violently with water. Danger of decomposition if exposed to heat.

The product is a strong oxidizing and reactive agent. Impurities, decomposition catalysts, incompatible substances, combustible substances, can cause self-accelerated, exothermic decomposition and oxygen formation.

Risk of overpressure and bursting due to decomposition in confined spaces and pipes. The release of oxygen can promote combustion.

## 10.4. Conditions to avoid

Avoid overheating. Prevent moisture or water from penetrating inside the containers. Avoid exposure to: light,heat. Avoid contact with: alkaline substances.

## 10.5. Incompatible materials



# Steelcoxide B

Page n. 8/16

Replaced revision:2 (Printed on: 05/06/2019)

Incompatible with: flammable substances, acetone, ethanol, glycerol, organic sulphides, hydrated bases, oxidising substances, iron,copper,bronze, chromium, zinc,I ead, silver, manganese, acetic acid.

Impurities, decomposition catalysts, metal salts, alkalis, hydrochloric acid, reducing agents (Risk of decomposition.). Flammable substances (Fire hazard). Organic solvents (Explosion hazard)

Group

## 10.6. Hazardous decomposition products

Oxygen.

## **SECTION 11. Toxicological information**

In the absence of experimental data for the product itself, health hazards are evaluated according to the properties of the substances it contains, using the criteria specified in the applicable regulation for classification. It is therefore necessary to take into account the concentration of the individual hazardous substances indicated in section 3, to evaluate the toxicological effects of exposure to the product.

## 11.1. Information on toxicological effects

#### Metabolism, toxicokinetics, mechanism of action and other information

It is an endogenous product formed in the cells of the organism. It penetrates through the skin and mucous membranes and decomposes in the underlying tissues. This causes diffuse infiltration of the released oxygen and the formation of emboli. In the mammalian organism, the enzymes that mostly work for the metabolism of the substance are glutathione peroxidase and catalase (INRS, 2007).

#### Information on likely routes of exposure

The primary route of exposure for hydrogen peroxide solution is through the respiratory tract and skin. (GESTIS Substance Database)

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Acute Effects:

Depending on the concentration: irritating to corrosive effects on the skin and mucous membranes, especially those of the eyes; Inflammatory changes of the respiratory tract, in extreme cases lung damage due to higher vapor / aerosol concentrations.

Chronic Effects:

Irritation of the mucous membranes (mainly of the eyes and throat) and gradual discoloration of the hair; possible skin changes. (GESTIS Substance Database)

Interactive effects

Information not available

ACUTE TOXICITY

ATE (Inhalation) of the mixture:> 20 mg/l ATE (Oral) of the mixture: 1043,88 mg/kg

ATE (Dermal) of the mixture: Not classified (no significant component)

LD50 (Oral) 1193 mg/kg Rat at the concentration of 35% LD50 (Dermal) 2000 mg/kg bw Rabbit

Method: equivalent or similar to OECD 401 Reliability (Klimisch score): 1 Species: Rat (Male / Female)

**Steelcoxide B** 

Printed on 29/06/2021

Page n. 9/16 Replaced revision:2 (Printed on: 05/06/2019)

Routes of exposure: oral Results: LD50 = 1026 mg / kg bw Details of the tested material: 70% Hydrogen Peroxide

Miele

Method: equivalent or similar to OECD 403 Reliability (Klimisch score): 1 Species: Rat (Sprague-Dawley; Male / Female) Routes of exposure: inhalation (vapors) Results: LD50 not determinable at concentrations of 170 mg / m3 Details of the tested material: 50% Hydrogen Peroxide

Group

Membe

Method: equivalent or similar to OECD 402 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White, male / female) Routes of exposure: cutaneous Results: LC50 = not determinable at a concentration of 2000 mg / kg bw Details of the tested material: Hydrogen Peroxide 35%

## SKIN CORROSION / IRRITATION

Corrosive for the skin

Method: equivalent or similar to OECD 404 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White, male / female) Routes of exposure: cutaneous Results: Category 2 Irritant (EC1272 / 2008) Details of the tested material: Hydrogen Peroxide 49.2%

## SERIOUS EYE DAMAGE / IRRITATION

Causes serious eye damage

Method: equivalent or similar to OECD 405 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White, male / female) Routes of exposure: ocular Results: Category 2 Irritant (GHS Criteria) Details of the tested material: Hydrogen Peroxide 10%

## RESPIRATORY OR SKIN SENSITISATION

Does not meet the classification criteria for this hazard class

GERM CELL MUTAGENICITY

Does not meet the classification criteria for this hazard class

Method: equivalent or similar to OECD 474 Reliability (Klimisch score): 1 In vivo test Species: Mouse (Swiss OF1 / ICO: OF1, male / female) Routes of exposure: intraperitoneal Results: negative with and without metabolic activation Details of the tested material: Hydrogen Peroxide 35%

## CARCINOGENICITY

Does not meet the classification criteria for this hazard class

	STEELCO S.P.A.	Revision nr. 3 Dated 18/05/2021
Miele Group Member	<b>.</b>	Printed on 29/06/2021
	Steelcoxide B	Page n. 10/16 Replaced revision:2 (Printed on: 05/06/2019)
Reference: publication (1996) Species: hamster (Syrian; Male / Femal Routes of exposure: oral mucosa Results: negative. NOAEL (carcinogeni Details of the tested material: Hydrogen Bibliographic reference: J Am Coll Toxic	e) city): not detected Peroxide 3% col 15 (1), 45-61	
REPRODUCTIVE TOXICITY		
Does not meet the classification criteria for this	hazard class	
Reference: publication (1958) Species: Rat (Osborne-Mendel; Male / I Routes of exposure: Oral Results: no effect observed Details of the tested material: 0.45% Hy Bibliographic reference: Nature 181, 14	Female) drogen Peroxide 53	
Adverse effects on sexual function and fertility		
Female rats that received 0.45% hydrog The fertility of male mice was not chang	gen peroxide with water to drink for 5 weeks gave birth t ed 3 months after administering a 1% product to drinking	o normal litters after mating with untreated males. g water for 4 weeks (INRS, 2007).
Adverse effects on development of the offspring	1	
Reference: publication (1958) Species: Rat (Wistar) Routes of exposure: Oral Results: no effect observed Details of the tested material: Hydrogen Bibliographic reference: Acta Obst Gyna	Peroxide 10% aec Japan 34, 2149-2154	
STOT - SINGLE EXPOSURE		
Does not meet the classification criteria for this	hazard class	
Based on available data, the substance has n hazard class.	o specific target organ toxicity effects for single expos	ure and is not classified under the relevant CLP
Target organ		
Respiratory trait.		
Route of exposure		
Inalation (vapour).		
STOT - REPEATED EXPOSURE		
Does not meet the classification criteria for this	hazard class	
Method: OECD 40 Reliability (Klimisch score): 1 Species: Mouse (C57BL / 6NCrlBR; n Routes of exposure: oral Results: negative. NOAEL = 100 ppm	nales / females)	
Method: OECD 412		

	STEELCO S.P.A.	Revision nr. 3 Dated 18/05/2021	
Miele Group Member		Printed on 29/06/2021	
	Steelcoxide B	Page n. 11/16 Replaced revision:2 (Printed on: 05/06/2019)	
Reliability (Klimisch score): 1 Species: Rat (Alpk: ApfSD; males / fen Routes of exposure: inhalation (aeroso Results: negative. NOAEL = 2.9 mg / n	nales) I) N <sup>3</sup> air		
ASPIRATION HAZARD			
Does not meet the classification criteria for this h	azard class		
No data are available on the hazard in ca	ase of aspiration.		
SECTION 12. Ecological information	ation		
This product is dangerous for the environment ar	nd the aquatic organisms. In the long term, it have	negative effects on aquatic environment.	
12.1. Toxicity			
Activated sludge:			
EC50 (30 min) = 466 mg/l (OECD 209) EC50 (3 h) >1000 mg/l (OECD 209)			
HYDROGEN PEROXIDE SOLUTION			
LC50 - for Fish	16,4 mg/l/96h Pimephales prome Test Guideline)	elas (USEPA Toxic Substances Control Act	
EC50 - for Crustacea	2,4 mg/l/48h Daphnia pulex (USE Guidelines)	EPA Toxic Substances Control Act Test	
EC50 - for Algae / Aquatic Plants	1,38 mg/l/72h Skeletonema costa	atum, soluzione al 35%	
	5 mg/L/96 n Pimephales promeia Test Guideline)		
Chronic NOEC for Algae / Aquatic Plants	1 mg/l Daphnia pulex (USEPA Toxic Substances Control Act Test Guidelines) 0,63 mg/l Skeletonema costatum, soluzione al 35%		
12.2. Persistence and degradability			
Quickly biodegradable, 99% in 30 min. (OECD	209)		
Solubility in water	100000 mg/l		
12.3. Bioaccumulative potential			
Partition coefficient: n-octanol/water	-1,57		
12.4. Mobility in soil			
Soil / sediment: Log KOC = 0.2 evaporation an Air, volatility, Henrym's constant = 0.75 kPa.m3	d adsorption not significant. 3 / mol. At 20 ° C		
12.5. Results of PBT and vPvB assessment			
On the basis of available data, the product does	not contain any PBT or vPvB in percentage ≥ thar	n 0,1%.	
12.6. Other adverse effects			
Information not available			

# STEELCO S.P.A.

Revision nr. 3 Dated 18/05/2021

Steelcoxide B

Printed on 29/06/2021

Page n. 12/16 Replaced revision:2 (Printed on: 05/06/2019)

## **SECTION 13. Disposal considerations**

### 13.1. Waste treatment methods

Reuse if possible. Product residues are to be considered special hazardous waste. The dangerousness of the waste that partially contains this product must be evaluated according to the laws in force. (Ref. Annex D -

Part IV of Legislative Decree n. 152/2006 and subsequent amendments and adjustments).

Group

Disposal must be entrusted to an authorized waste management company, in compliance with national and possibly local regulations.

The legal responsibility for disposal lies with the producer / holder of the waste.

Different CER (European Refusal Code) codes could be applied to this product according to the specific circumstances that generated the waste, any alterations and contaminations.

The product as it is, out of specification in the original packaging, or poured into a suitable container for disposal as waste, or the product in specification but no longer usable (for example following an accidental spill), it is to be classified with a CER code compatible with the description of use indicated in section 1.2.

The suitable final destination of the waste will be evaluated by the manufacturer according to the chemical-physical characteristics of the waste itself, compatible with the plant

authorized to which it will be conferred for recovery, treatment or final disposal according to the procedures provided for by the regulations in force. Disposal via the wastewater drain is not permitted.

#### CONTAMINATED PACKAGING

Contaminated packaging must be sent, properly labeled, for recovery or disposal in compliance with national regulations on the management of waste and is to be classified with the following CER code:

15 01 10 \*: packaging containing residues of dangerous substances or contaminated by these substances.

## **SECTION 14. Transport information**

#### 14.1. UN number

ADR / RID, IMDG, 2014 IATA:

## 14.2. UN proper shipping name

ADR / RID:	HYDROGEN PEROXIDE, AQUEOUS SOLUTION
IMDG:	HYDROGEN PEROXIDE, AQUEOUS SOLUTION
IATA:	HYDROGEN PEROXIDE, AQUEOUS SOLUTION

## 14.3. Transport hazard class(es)

ADR / RID:	Class: 5.1	Label: 5.1 (8)
IMDG:	Class: 5.1	Label: 5.1 (8)
IATA:	Class: 5.1	Label: 5.1 (8)



## 14.4. Packing group

ADR / RID, IMDG, II IATA:

## 14.5. Environmental hazards

NO

ADR / RID:

Steelco Miele Group Member		STEELCO S.P.A.		Revision nr. 3 Dated 18/05/2021	
				Printed on 29/06/	/2021
		Steelcoxide E	Steelcoxide B		n:2 (Printed on: 05/06/2019)
	NO				
IMDG: IATA:	NO				
14.6. Special precaution	ons for user				
ADR / RID:	ИН	N - Kemler: 58	Limited Quantitie	es: 1 L	Tunnel restriction code: (E)
	Sp	ecial provision: -			
IMDG:	EM	IS: F-H, S-Q	Limited Quantitie	es: 1 L	
IATA:	Ca	rgo:	Maximum quanti	ty: 5 L	Packaging instructions: 554
	Pa	ss.:	Maximum quanti	ty: 1 L	Packaging instructions: 550
	Sp	ecial provision:	-		
14.7. Transport in bull	c according to Annex II	of Marpol and the IBC Code			
Information not relevant	t				
SECTION 15. R	legulatory inform	ation			
15.1. Safety, health a	and environmental regu	lations/legislation specific for the s	ubstance or mixt	ure	
Seveso Category - Dire	ctive 2012/18/EC: P8				
Restrictions relating to t	the product or contained	substances pursuant to Annex XVII to	EC Regulation 19	<u>07/2006</u>	
Product Point	3				
	Liquid : catego a) haza catego b) haza other ti c) haza d) haza	substances or mixtures that correspon ries referred to Annex I of Regulation ( ard classes 2.1 to 2.4, 2.6 and 2.7, 2.8 ries 1 and 2, 2.15 types from A to F; ard classes 3.1 to 3.6, 3.7 adverse effe han narcotic effects, 3.9 and 3.10; ard class 4.1; ard class 5.1.	d to the criteria rel EC) No. 1272/200 types A and B, 2.9 cts on sexual func	ating to one o 18: 9, 2.10, 2.12, 2 tion and fertili	f the following hazard classes or 2.13 categories 1 and 2, 2.14 ity or development, 3.8 effects
Contained substance					
Point	75				
		HYE Reg	ROGEN PEROXI . no.: 01-2119485	DE SOLUTIO 845-22-xxxx	N
	Substa a) subs - categ classifi - Repro exposu - skin s - skin c - seriou b) subs Counci c) subs one of d) subs	Inces included in one or more of the fo stances classified in one of the followin ory 1A, 1B or 2 carcinogenicity, catego ed due to effects following exposure by oductive toxicity of category 1A, 1B or ure by inhalation only; corrosion of category 1, 1A or 1B; corrosion of category 1, 1A, 1B or 1C o us category 1 eye damage or category stances listed in Annex II of Regulation if (*); stances listed in Annex IV of Regulation the columns g, h or i of the table of tha stances listed in Appendix 13 of this an	llowing: g classes in Anne. ory 1A, 1B or 2 ger y inhalation only; 2, but excluding su r skin irritation of c 2 eye irritation; (EC) no. 1223/20 n (EC) no. 1223/20 n (EC) no. 1223/20 t annex; nex. The ancillary	x VI, part 3, oi rm cell mutage ubstances clas category 2; 09 of the Euro 209 for which a requirements	f Regulation (EC) no. 1272/2008: enicity, but excluding substances ssified due to effects following opean Parliament and of the a condition is indicated in at least referred to in points 7 and 8 of

Group Group	STEELCO S.P.A.	Revision nr. 3 Dated 18/05/2021			
Miele Member	Staalaavida D	Printed on 29/06/2021			
	Steelcoxide B	Page n. 14/16 Replaced revision:2 (Printed on: 05/06/2019)			
the column 2 of this item apply to all mixtures intended for tattooing practices, regardless of whether contain one of the substances referred to in points a) to d) of this column and entry.					
Regulation (EC) No. 2019/1148 - on the marketir	ng and use of explosives precursors				
HYDROGEN PEROXIDE (CAS 7722-84-1) - Restricted explosives precursor The acquisition, introduction, possession or use of that restricted explosives precursor by members of the general public is subject to a restriction as set out in Article 5(1) and (3). Restricted explosives precursors shall not be made available to, or introduced, possessed or used by members of the general public.					
obligations as set out in Article 9. All suspicious transactions and significant disapp	earances and thefts must be reported to the releva	nt national contact point.			
Italian Ministry of the Interior Tel.: +390646542182 E-Mail: precursori@dcpc.interno.it					
Substances in Candidate List (Art. 59 REACH)					
On the basis of available data, the product does	not contain any SVHC in percentage ≥ than 0,1%.				
Substances subject to authorisation (Annex XIV	REACH)				
None					
Substances subject to exportation reporting purs	uant to (EC) Reg. 649/2012:				
None					
Substances subject to the Rotterdam Convention:					
None					
Substances subject to the Stockholm Convention:					
None					
Healthcare controls					
Workers exposed to this chemical agent must not undergo health checks, provided that available risk-assessment data prove that the risks related to the workers' health and safety are modest and that the 98/24/EC directive is respected.					
Class IIb medical device according to the 93/42 regulation.					
15.2. Chemical safety assessment					
A chemical safety assessment has not been performed for the preparation/for the substances indicated in section 3.					
SECTION 16. Other information					
Text of hazard (H) indications mentioned in section 2-3 of the sheet:					
Ox. Liq. 1Oxidising liquid, catOx. Sol. 2Oxidising solid, cate	egory 1 egory 2				

teelco

Page n. 15/16 Replaced revision:2 (Printed on: 05/06/2019)

Acute Tox. 4	Acute toxicity, category 4
Skin Corr. 1A	Skin corrosion, category 1A
Skin Corr. 1C	Skin corrosion, category 1C
Eye Dam. 1	Serious eye damage, category 1
Aquatic Chronic 3	Hazardous to the aquatic environment, chronic toxicity, category 3
H271	May cause fire or explosion; strong oxidiser.
H272	May intensify fire; oxidiser.
H302	Harmful if swallowed.
H332	Harmful if inhaled.
H314	Causes severe skin burns and eye damage.
H318	Causes serious eye damage.
H412	Harmful to aquatic life with long lasting effects.

Group

Miele

#### LEGEND:

- ADR: European Agreement concerning the carriage of Dangerous goods by Road
- CAS NUMBER: Chemical Abstract Service Number
- CE50: Effective concentration (required to induce a 50% effect)
- CE NUMBER: Identifier in ESIS (European archive of existing substances)
- CLP: EC Regulation 1272/2008
- DNEL: Derived No Effect Level
- EmS: Emergency Schedule
- GHS: Globally Harmonized System of classification and labeling of chemicals
- IATA DGR: International Air Transport Association Dangerous Goods Regulation
- IC50: Immobilization Concentration 50%
- IMDG: International Maritime Code for dangerous goods
- IMO: International Maritime Organization
- INDEX NUMBER: Identifier in Annex VI of CLP
- LC50: Lethal Concentration 50%
- LD50: Lethal dose 50%
- **OEL: Occupational Exposure Level**
- PBT: Persistent bioaccumulative and toxic as REACH Regulation
- PEC: Predicted environmental Concentration
- PEL: Predicted exposure level
- PNEC: Predicted no effect concentration
- REACH: EC Regulation 1907/2006
- RID: Regulation concerning the international transport of dangerous goods by train
- TLV: Threshold Limit Value
- TLV CEILING: Concentration that should not be exceeded during any time of occupational exposure.
- TWA STEL: Short-term exposure limit
- TWA: Time-weighted average exposure limit
- VOC: Volatile organic Compounds
- vPvB: Very Persistent and very Bioaccumulative as for REACH Regulation
- WGK: Water hazard classes (German).

## GENERAL BIBLIOGRAPHY

- 1. Regulation (EC) 1907/2006 (REACH) of the European Parliament 2. Regulation (EC) 1272/2008 (CLP) of the European Parliament
- 3. Regulation (EU) 790/2009 (I Atp. CLP) of the European Parliament
- 4. Regulation (EU) 2015/830 of the European Parliament
- 5. Regulation (EU) 286/2011 (II Atp. CLP) of the European Parliament
- 6. Regulation (EU) 618/2012 (III Atp. CLP) of the European Parliament

- 7. Regulation (EU) 487/2013 (IV Atp. CLP) of the European Parliament 8. Regulation (EU) 944/2013 (V Atp. CLP) of the European Parliament 9. Regulation (EU) 605/2014 (VI Atp. CLP) of the European Parliament
- 10. Regulation (EU) 2015/1221 (VII Atp. CLP) of the European Parliament
- 11. Regulation (EU) 2016/918 (VIII Atp. CLP) of the European Parliament 12. Regulation (EU) 2016/1179 (IX Atp. CLP)
- 13. Regulation (EU) 2017/776 (X Atp. CLP)



#### Page n. 16/16 Replaced revision:2 (Printed on: 05/06/2019)

- 14. Regulation (EU) 2018/669 (XI Atp. CLP) 15. Regulation (EU) 2018/1480 (XIII Atp. CLP)
- 16. Regulation (EU) 2019/521 (XII Atp. CLP)
- 17. Regulation (EU) 2019/1148
- 18. Regulation (EU) 2020/217 (XIV Atp. CLP)
- The Merck Index. 10th Edition
- Handling Chemical Safety
- INRS Fiche Toxicologique (toxicological sheet)
- Patty Industrial Hygiene and Toxicology
- N.I. Sax Dangerous properties of Industrial Materials-7, 1989 Edition
- IFA GESTIS website
- FCHA website

Database of SDS models for chemicals - Ministry of Health and ISS (Istituto Superiore di Sanità) - Italy

#### TRAINING FOR WORKERS:

The training of workers must include contents, updates and duration according to the risk profiles assigned to the working sectors of membership, according to the procedures provided for by Legislative Decree 81/2008.

#### NOTE FOR THE RECIPIENT OF THE SAFETY DATA SHEET (SDS):

It is the recipient of this SDS who must ensure that the information contained is read and understood by all persons who manipulate, store, use, or in any case come into contact in any way with the substance or mixture to which this sheet refers. In particular the recipient must provide adequate training to personnel assigned to the use of hazardous substances or mixtures. The recipient must make sure of the suitability and completeness of the information in relation to the specific use of the substance or mixture. However, the substance or mixture to which this SDS refers must not be used for uses other than those specified in section 1. Do not assume responsibility for improper use. Since the use of the product does not fall under the direct control of the Supplier, it is the user's obligation observe, under their own responsibility, the laws and regulations in force regarding national and community hygiene and safety. The information contained in this SDS is provided in good faith and is based on the current state of scientific and technical knowledge as of the date revision indicated, available from the Supplier indicated in section 1 of this sheet. The SDS should not be construed as a guarantee of any specific property of the substance or mixture. The information refers only to the substance or mixture specifically designated for section 1 and may not be valid for the substance or mixture used in combination with other materials or in other processes not specifically indicated in the text.

This version of the SDS supersedes all previous versions.

CALCULATION METHODS FOR CLASSIFICATION:

Chemical and physical hazards: Product classification derives from criteria established by the CLP Regulation, Annex I, Part 2. The data for evaluation of chemical-physical properties are reported in section 9.

Health hazards: Product classification is based on calculation methods as per Annex I of CLP, Part 3, unless determined otherwise in Section 11. Environmental hazards: Product classification is based on calculation methods as per Annex I of CLP, Part 4, unless determined otherwise in Section 12.

Classification of the mixture according to Regulation (EC) n. 1272/2008 Classification procedure

May intensify fire; oxidiser.	H272	Based on experimental data
Harmful if swallowed.	H302	Based on literature data
Causes skin irritation	H315	Based on literature data
May cause respiratory irritation	H335	Based on literature data
Causes serious eye damage.	H318	Based on literature data
Harmful to aquatic life with long lasting effects.	H412	Based on literature data

Changes to previous review

The following sections were modified:

01/02/03/04/05/06/07/08/09/10/11/12/13/14/15/16.

Room 201 & 202, 12th Building, Shenzhen Software Park Phase II, 1 Keji Middle 2nd Road, Yuehai Subdistrict, Nanshan District, Shenzhen, 518057, Guangdong, China. Tel:86-755-26722890 Fax:86-755-2652 6612 E-mail: sonoscape@sonoscape.net



Sept. 6th, 2019

Elisa Basso STEELCO SPA Via Balegante, 27 - 31039 Riese Pio X (TV) – ITALY

Dear Ms. Basso,

This is to claim that the following Sonoscape gastroscopes and colonoscopes are able to be reprocessed in Steelco's automated endoscope reprocessor(AER) and also be stored in Steelco's endoscopy drying cabinet, which are both listed here below.

Sonoscape's gastroscope & colonoscope models:

- Gastroscope: EG-550, EG-550L, EG-500, EG-500L
- Colonoscope: EC-550, EC-550T, EC-550L, EC-550L/T, EC-500, EC-500T, EC-500L, EC-500L/T

Steelco's AER & Drying Cabinet model:

- AER: Endoscope Washer EW1
- Drying Cabinet ED 200

Also, Steelco shall provide the scope reprocessing connectors for their AER to reprocess above mentioned Sonoscape endoscopes.

Sincerely 6/9/2019 Ronan Lao

Ronan Lao Global Marketing Director Endoscopy Division Sonoscape Medical Corp. www.sonoscape.com

Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021		
Miele Member		Printed on 29/06/2021		
	Steelcoxide A	Page n. 1/22 Replaced revision:3 (Printed on: 18/05/2021)		
	•			
	Safety Data Shee According to Annex II to REACH - Regulation	<b>et</b> n 2015/830		
SECTION 1. Identification of the	e substance/mixture and of the	company/undertaking		
1.1. Product identifier				
Product name	Steelcoxide A			
1.2. Relevant identified uses of the substan	ce or mixture and uses advised against			
Intended Use:	Soluzione "A" of tricomponent system ' sterilizing solution "SteelcoXide B".	SteelcoXide". Alcohol based activator of concetrated		
Uses advised against:	Uses other than those indicated.			
1.3 Details of the supplier of the safety da	ta sheet			
Name	STEELCO S.p.A.			
Full address District and Country	Via Balegante, 27 31039 Riese Pio X (TV) ITALY			
	tel. +39 0423 7561			
e-mail address of the competent person	fax +39 0423 755528			
responsible for the Safety Data Sheet	info@steelcogroup.com			
1.4. Emergency telephone number				
For urgent	Centros de Orientação de Doentes	s Urgentes (CODU): 800250250		
inquiries refer to				
SECTION 2. Hazards identificat	ion			
2.1. Classification of the substance or mixtu	re			
The product is classified as hazardous pursua supplements). The product thus requires a safe Any additional information concerning the risks	ant to the provisions set forth in (EC) Regul ty datasheet that complies with the provisions for health and/or the environment are given in	ation 1272/2008 (CLP) (and subsequent amendments and of (EU) Regulation 2015/830. sections 11 and 12 of this sheet.		
Hazard classification and indication: Flammable liquid, category 2 Skin corrosion, category 1 Specific target organ toxicity - single exposure	H225 Hig H314 Cau e, category 3 H336 May	hly flammable liquid and vapour. ises severe skin burns and eye damage y cause drowsiness or dizziness.		
2.2. Label elements				
Hazard labelling pursuant to EC Regulation 127	2/2008 (CLP) and subsequent amendments a	and supplements.		
Hazard pictograms:				

	Group	s	TEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
teelco	Miele Member		Staalaavida A	Printed on 29/06/2021
			Steelcoxide A	Page n. 2/22 Replaced revision:3 (Printed on: 18/05/2021)
	!			
Signal words:	Danger			
Hazard statements:				
H225 H314 H336	Highly flammable li Causes severe skir May cause drowsin	quid and vapour. h burns and eye o less or dizziness.	damage	
Precautionary statements:				
P210Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.P312Call a POISON CENTRE if you feel unwell.P261Avoid breathing vapours / spray.P280Wear eye protection / face protection.P305+P351+P338IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.P337+P313If eye irritation persists: Get medical advice.				
Contains:	PROPAN-2-OL			
2.3. Other hazards On the basis of available da	ata, the product does	not contain any	PBT or vPvB in percentage ≥ than 0,1%	5.
SECTION 3. Com	position/inforn	nation on in	gredients	
3.2. Mixtures				
Contains:				
Identification	x = Co	nc. % Cla	ssification 1272/2008 (CLP)	
PROPAN-2-OL           CAS         67-63-0           EC         200-661-7           INDEX         603-117-00-0           Reg. no.         01-211945755	15 ≤ x · 8-25-XXXX	< 35 Flai	n. Liq. 2 H225, Eye Irrit. 2 H319, STOT	SE 3 H336
ACETIC ACID				
CAS 64-19-7	3 ≤ x <	10 Flan note	n. Liq. 3 H226, Skin Corr. 1A H314, Eye e/notes according to Annex VI to the CL	e Dam. 1 H318, Classification P Regulation: B
EC 200-580-7		Spe Skii 2; H	crific concentration limits (Annex VI of R n Corr. 1A; H314: C ≥90 % Skin Corr. 1 I315: 10 % ≤C <25 % Eye Irrit. 2; H319:	Reg. 1272/2008) B; H314: 25 %≤ C <90 % Skin Irrit. : 10 % ≤ C < 25 %

	Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
And a second sec	Member	Steelcoxide A	Printed on 29/06/2021 Page n. 3/22 Replaced revision:3 (Printed on: 18/05/2021)
INDEX 607-002-00-6 Reg. no. 01-2119475328-30-xxxx			
ETHANOL CAS 64-17-5 EC 200-578-6 INDEX 603-002-00-5 Reg. no. 01-2119457610-43-XXXX	3 ≤x< 8	Flam. Liq. 2 H225, Eye Irrit. 2 H319	

The full wording of hazard (H) phrases is given in section 16 of the sheet.

## **SECTION 4. First aid measures**

#### 4.1. Description of first aid measures

EYES: Remove contact lenses, if present. Wash immediately with plenty of water for at least 30-60 minutes, opening the eyelids fully. Get medical advice/attention.

SKIN: Remove contaminated clothing. Rinse skin with a shower immediately. Get medical advice/attention.

INGESTION: Have the subject drink as much water as possible. Get medical advice/attention. Do not induce vomiting unless explicitly authorised by a doctor.

INHALATION: Get medical advice/attention immediately. Remove victim to fresh air, away from the accident scene. If the subject stops breathing, administer artificial respiration. Take suitable precautions for rescue workers.

PROTECTION MEASURES FOR FIRST AID: for the PPE necessary for first aid interventions refer to section 8.2 of this safety data sheet.

#### 4.2. Most important symptoms and effects, both acute and delayed

Specific information on symptoms and effects caused by the product are unknown.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Information not available

## **SECTION 5. Firefighting measures**

### 5.1. Extinguishing media

SUITABLE EXTINGUISHING EQUIPMENT

Extinguishing substances are: carbon dioxide, foam, chemical powder. For product loss or leakage that has not caught fire, water spray can be used to disperse flammable vapours and protect those trying to stem the leak.

UNSUITABLE EXTINGUISHING EQUIPMENT

Do not use jets of water. Water is not effective for putting out fires but can be used to cool containers exposed to flames to prevent explosions.

## 5.2. Special hazards arising from the substance or mixture

HAZARDS CAUSED BY EXPOSURE IN THE EVENT OF FIRE Excess pressure may form in containers exposed to fire at a risk of explosion. Do not breathe combustion products.

PROPAN-2-OL Carbon oxides.

## 5.3. Advice for firefighters

#### GENERAL INFORMATION

Use jets of water to cool the containers to prevent product decomposition and the development of substances potentially hazardous for health. Always

Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021			
Miele Member		Printed on 29/06/2021			
	Steelcoxide A	Page n. 4/22 Replaced revision:3 (Printed on: 18/05/2021)			
wear full fire prevention gear. Collect extinguishing water to prevent it from draining into the sewer system. Dispose of contaminated water used for extinction and the remains of the fire according to applicable regulations. SPECIAL PROTECTIVE EQUIPMENT FOR FIRE-FIGHTERS Normal fire fighting clothing i.e. fire kit (BS EN 469), gloves (BS EN 659) and boots (HO specification A29 and A30) in combination with self-contained open circuit positive pressure compressed air breathing apparatus (BS EN 137).					
SECTION 6. Accidental release	measures				
6.1. Personal precautions, protective equipm	ent and emergency procedures				
6.1.1 For those who do not intervene directly Do not take any action involving any personal material. Wear suitable protective equipment (including p of skin, eyes and personal clothing. Wear appro Do not inhale the mists / vapors / fumes. Ave personnel not authorized to intervene directly in	I risk or without proper training. Evacuate the surroun personal protective equipment referred to in section 8 of priate respirator when ventilation is inadequate. Did the dispersion of the product into the environment the event of accidental release.	ding areas. Do not touch or walk on the spilled this Safety Data Sheet) to prevent contamination . Follow the appropriate internal procedures for			
6.1.2 For those who intervene directly Stop the leak if there is no danger. Evacuate unauthorized personnel. Wear suita procedures for authorized personnel. Isolate the	6.1.2 For those who intervene directly Stop the leak if there is no danger. Evacuate unauthorized personnel. Wear suitable protective equipment. (see section 8 of this Safety Data Sheet). Follow the appropriate internal procedures for authorized personnel. Isolate the danger area and deny entry. Ventilate enclosed spaces before entering.				
6.2. Environmental precautions					
The product must not penetrate into the sewer system or come into contact with surface water or ground water.					
6.3. Methods and material for containment and cleaning up					
Collect the leaked product into a suitable container. Evaluate the compatibility of the container to be used, by checking section 10. Absorb the remainder with inert absorbent material. Make sure the leakage site is well aired. Contaminated material should be disposed of in compliance with the provisions set forth in point 13.					
6.4. Reference to other sections					
Any information on personal protection and disp	Any information on personal protection and disposal is given in sections 8 and 13.				
SECTION 7. Handling and stora	ge				
7.1. Precautions for safe handling					
Keep away from heat, sparks and naked flames; do not smoke or use matches or lighters. Without adequate ventilation, vapours may accumulate at ground level and, if ignited, catch fire even at a distance, with the danger of backfire. Avoid bunching of electrostatic charges. When performing transfer operations involving large containers, connect to an earthing system and wear antistatic footwear. Vigorous stirring and flow through the tubes and equipment may cause the formation and accumulation of electrostatic charges. In order to avoid the risk of fires and explosions, never use compressed air when handling. Open containers with caution as they may be pressurised. Do not eat, drink or smoke during use. Avoid leakage of the product into the environment.					
7.2. Conditions for safe storage, including any incompatibilities					
Store only in the original container. Store the containers sealed, in a well ventilated place, away from direct sunlight. Store in a cool and well ventilated place, keep far away from sources of heat, naked flames and sparks and other sources of ignition. Keep containers away from any incompatible materials, see section 10 for details.					
7.3. Specific end use(s)					

Information not available



Steelcoxide A

Page n. 5/22 Replaced revision:3 (Printed on: 18/05/2021)

# SECTION 8. Exposure controls/personal protection

## 8.1. Control parameters

Regulatory References:

DEU	Deutschland	Technischen Regeln für Gefahrstoffe (TRGS 900) - Liste der Arbeitsplatzgrenzwerte und Kurzzeitwerte.
		MAK- und BAT-Werte-Liste 2020, Ständige Senatskommission zur Prüfung gesundheitsschädlicher
		Arbeitsstoffe, Mitteilung 56
ESP	España	Límites de exposición profesional para agentes químicos en España 2019
FRA	France	Valeurs limites d'exposition professionnelle aux agents chimiques en France. ED 984 - INRS
GBR	United Kingdom	EH40/2005 Workplace exposure limits (Fourth Edition 2020)
EU	OEL EU	Directive (EU) 2019/1831; Directive (EU) 2019/130; Directive (EU) 2019/983; Directive (EU) 2017/2398;
		Directive (EU) 2017/164; Directive 2009/161/EU; Directive 2006/15/EC; Directive 2004/37/EC; Directive
		2000/39/EC; Directive 98/24/EC; Directive 91/322/EEC.
	TLV-ACGIH	ACGIH 2020

## PROPAN-2-OL

Туре	Country	TWA/8h		STEL/15min		Remarks / Observatio	ons	
		mg/m3	ppm	mg/m3	ppm			
AGW	DEU	500	200	1000	400			
MAK	DEU	500	200	1000	400			
VLA	ESP	500	200	1000	400			
VLEP	FRA			980	400			
WEL	GBR	999	400	1250	500			
TLV-ACGIH		492	200	983	400			
Predicted no-effect concentration	- PNEC							
Normal value in fresh water				140,9	mg	/1		
Normal value in marine water				140,9	mg	/I		
Normal value for fresh water sed	iment			552	mg	/kg		
Normal value for marine water se	ediment			552	mg	/kg		
Normal value for water, intermitte	ent release			140,9	mg	/1		
Normal value of STP microorgan	isms			2,251	g/l			
Normal value for the food chain (	secondary poisoni	ng)		160	mg	/kg		
Normal value for the terrestrial co	ompartment			28	mg	/kg		
Health - Derived no-effect I	evel - DNEL / D	MEL						
	Effects on				Effects on			
	consumers				workers			
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic systemic	Acute local	Acute systemic	Chronic local	Chronic systemic
Oral					VND	VND	VND	26 mg/kg bw/d
Inhalation	VND	VND	VND	89 mg/m3	VND	VND	VND	500 mg/m3
Skin	VND	VND	VND	319 mg/kg bw/d	VND	VND	VND	888 mg/kg bw/d
ACETIC ACID								
Threshold Limit Value								
Туре	Country	TWA/8h		STEL/15min		Remarks / Observation	ons	
		mg/m3	ppm	mg/m3	ppm	C DOOI Valle		
AGW	DEU	25	10	50 (C)	20 (C)			
MAK	DEU	25	10	50	20			

•			STEELCO	S.P.A.	Re Da	evision nr. 4 ated 24/06/2021		
teelco	Miele Member		Steelcov	ido A	Pr	inted on 29/06/202	21	
			Sleeicux	iue A	Pa Re	age n. 6/22 eplaced revision:3	(Printed on: 18/05/202	1)
VLA	ESP	25	10	50	20			
VLEP	FRA	25	10	50	20			
WEL	GBR	25	10	50	20			
OEL	EU	25	10	50	20			
TLV-ACGIH		25	10	37	15			
Predicted no-effect concentra	ation - PNEC							
Normal value in fresh water				3,058	m	g/l		
Normal value in marine water	r			305,8	μί	g/L		
Normal value for fresh water	sediment			11,36	m	g/kg		
Normal value for marine wate	er sediment			1,136	m	g/kg		
Normal value for water, interr	mittent release			30,58	m	g/l		
Normal value of STP microor	ganisms			85	m	g/l		
Normal value for the terrestria	al compartment			470	hố	g/kg soil dw		
Health - Derived no-effe	ct level - DNEL / DN Effects on	/EL			Effects on			
Route of exposure	consumers Acute local	Acute systemic	Chronic local	Chronic	workers Acute local	Acute	Chronic local	Chronic
	NDI	Addie Systemie		systemic	Addie 100ai	systemic	officine local	systemic
Oral	NPI		NPI					
Inhalation Skin	25 mg/m3 VND	NPI NPI	25 mg/m3 VND	NPI NPI	25 mg/m3 VND	NPI NPI	25 mg/m3 VND	NPI NPI
ETHANOL Threshold Limit Value								
Type	Country	TWA/8h				Remarks	o /	
	Country			STEL/15min		Remarks	57	
	Country	mg/m3	ppm	mg/m3	ppm	Observa	itions	
AGW	DEU	mg/m3	ppm 200	mg/m3	ppm 800	Observa	ations	
AGW	DEU	mg/m3 380 380	ppm 200 200	mg/m3 1520	ppm 800 800	Observa	tions	
AGW MAK VLA	DEU DEU ESP	mg/m3 380 380	ppm 200 200	mg/m3 1520 1520 1910	ppm 800 800 1000	Observa	tions	
AGW MAK VLA VLEP	DEU DEU ESP FRA	mg/m3 380 380	ppm 200 200	s r EL/15min mg/m3 1520 1520 1910 9500	ppm 800 800 1000 5000	Observa	s /	
AGW MAK VLA VLEP	DEU DEU ESP FRA GBR	mg/m3 380 380 1900	ppm 200 200 1000	STEL/15min           mg/m3           1520           1520           1910           9500	ppm 800 800 1000 5000	Observa	s / . itions	
AGW MAK VLA VLEP WEL TLV-ACGIH	DEU DEU ESP FRA GBR	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500	ppm 800 800 1000 5000	Observa	s /	
AGW MAK VLA VLEP WEL TLV-ACGIH	DEU DEU ESP FRA GBR	mg/m3 380 380 1900 1920	ррт 200 200 1000 1000	STEL/15min           mg/m3           1520           1520           1910           9500           1884	ppm 800 800 1000 5000 1000	Observa	s /	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra	DEU DEU ESP FRA GBR ation - PNEC	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884	ppm 800 800 1000 5000 1000	Observa	itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water	DEU DEU ESP FRA GBR ation - PNEC	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 700	ррт 800 800 1000 5000 1000	Gbserva g/L	s / itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water	DEU DEU ESP FRA GBR ation - PNEC	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 2.0	ррт 800 800 1000 5000 1000 µg т	g/L g/L	s / .	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value for fresh water	DEU DEU ESP FRA GBR ation - PNEC	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6	ррт 800 800 1000 5000 1000 µg т	g/L g/kg	s	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value for fresh water	DEU DEU ESP FRA GBR ation - PNEC	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 3,7	ppm 800 800 1000 5000 1000 µg m m m	g/L g/kg g/kg	itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value for fresh water Normal value for marine water	DEU DEU ESP FRA GBR ation - PNEC	mg/m3 380 380 1900 1920	ppm 200 200 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 2,75	ррт 800 800 1000 5000 1000 ис т т т т т	g/L g/kg g/l	s / itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value for fresh water Normal value for marine water Normal value for marine water Normal value of STP microor	DEU DEU ESP FRA GBR ation - PNEC r sediment er sediment mittent release ganisms	mg/m3 380 380 1900 1920	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 2,75 580	ррт 800 800 1000 5000 1000 ис т т т т т т	g/L g/l g/kg g/l g/l g/l	s / itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in fresh water Normal value for fresh water Normal value for marine water Normal value for marine water Normal value for marine water Normal value for marine water Normal value for the food cha	DEU DEU ESP FRA GBR ation - PNEC r sediment er sediment mittent release ganisms ain (secondary poisonin	mg/m3 380 380 1900 1920 9)	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 2,75 580 380	ррт 800 800 1000 5000 1000 44 т т т т т т т т	g/L g/kg g/kg g/l g/kg g/l g/kg	s / itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value for fresh water Normal value for marine water Normal value for marine water Normal value for marine mater Normal value for the stater Normal value for the food cha Normal value for the food cha	DEU DEU ESP FRA GBR ation - PNEC r sediment er sediment nittent release ganisms ain (secondary poisonin al compartment	mg/m3 380 380 1900 1920 9)	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 2,75 580 380 630	ррт 800 800 1000 5000 1000 ис т т т т т т т	g/L g/l g/kg g/kg g/l g/kg soil dw	s / itions	
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in fresh water Normal value for fresh water Normal value for marine water Normal value for marine water Normal value for marine water Normal value for the stress Normal value for the terrestria Health - Derived no-effe	DEU DEU ESP FRA GBR ation - PNEC r sediment er sediment nittent release ganisms ain (secondary poisonin al compartment ct level - DNEL / DN Effects on consumers	mg/m3 380 380 1900 1920 9) MEL	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 2,75 580 380 630	ppm 800 800 1000 5000 1000 0 0 0 0 0 0 0 0 0 0 0 0	g/L g/L g/l g/kg g/l g/l g/kg soil dw		
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value for fresh water Normal value for marine water Normal value for marine water Normal value for marine water Normal value for the terrestria Normal value for the terrestria Health - Derived no-effe Route of exposure	DEU DEU ESP FRA GBR ation - PNEC r sediment er sediment mittent release ganisms ain (secondary poisonin, al compartment <b>ct level - DNEL / DN</b> Effects on consumers Acute local	mg/m3 380 380 1900 1920 1920 9) MEL Acute systemic	ppm 200 200 1000 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 1884 960 790 3,6 2,9 2,75 580 3,6 2,9 2,75 580 380 630	ppm 800 800 1000 5000 1000 0 0 0 0 0 0 0 0 0 0 0 0	g/L g/L g/l g/kg g/l g/kg g/l g/kg food g/kg soil dw Acute	Chronic local	Chronic
AGW MAK VLA VLEP WEL TLV-ACGIH Predicted no-effect concentra Normal value in fresh water Normal value in marine water Normal value in marine water Normal value for fresh water Normal value for marine water Normal value for marine water Normal value for marine water Normal value for the terrestria Health - Derived no-effe Route of exposure Oral	DEU DEU ESP FRA GBR ation - PNEC r sediment er sediment mittent release ganisms ain (secondary poisonin, al compartment <b>ct level - DNEL / DN</b> Effects on consumers Acute local	mg/m3 380 380 1900 1920 9) MEL Acute systemic	ppm 200 200 1000 1000	STEL/15min mg/m3 1520 1520 1910 9500 1884 960 790 3,6 2,9 2,75 580 380 630 630 Chronic systemic	ррт 800 800 1000 5000 1000 400 00 00 00 00 00 00 00 00 00 00 00	g/L g/L g/l g/kg g/kg g/l g/kg food g/kg food g/kg soil dw Acute systemic NPI	Chronic local	Chronic systemic 87 mg/kg bw/d

	Crown		STEELCO	) S.P.A.		Revision nr. 4 Dated 24/06/2021		
teelco	Miele Member			Printed on 29/06/2021	1			
			Steelcox	cide A		Page n. 7/22 Replaced revision:3 (I	Printed on: 18/05	5/2021)
Skin	NPI	NPI	NPI	206 mg/kg bw/d	NPI	NPI	NPI	343 mg/kg bw/d
Legend:								
(C) = CEILING ; INHAL	(C) = CEILING ; INHAL = Inhalable Fraction ; RESP = Respirable Fraction ; THORA = Thoracic Fraction.							
VND = hazard identified b	out no DNEL/PNEC ava	ailable ; NE	EA = no exposure	expected ; N	IPI = no ha	zard identified.		
8.2. Exposure controls	5							
As the use of adequate t through effective local as When choosing personal Personal protective equip	technical equipment m piration. protective equipment, a ment must be CE mark	ust always ta ask your che ked, showing	ake priority over p mical substance s that it complies w	personal protec supplier for advi vith applicable s	tive equipn ce. .tandards.	nent, make sure	that the wor	kplace is well aired
Provide an emergency sh	ower with face and eye	e wash statio	n.					
HAND PROTECTION Protect hands with category III work gloves (see standard EN 374). The following should be considered when choosing work glove material: compatibility, degradation, failure time and permeability. The work gloves' resistance to chemical agents should be checked before use, as it can be unpredictable. The gloves' wear time depends on the duration and type of use.								
SKIN PROTECTION Wear category I professional long-sleeved overalls and safety footwear (see Regulation 2016/425 and standard EN ISO 20344). Wash body with soap and water after removing protective clothing.								
Consider the appropriater	ness of providing antist	atic clothing i	in the case of wor	king environme	ents in whic	h there is a risk o	f explosion.	
EYE PROTECTION Wear airtight protective g	oggles (see standard E	N 166).						
RESPIRATORY PROTECTION If the threshold value (e.g. TLV-TWA) is exceeded for the substance or one of the substances present in the product, wear a mask with a type AX filter, whose limit of use will be defined by the manufacturer (see standard EN 14387). In the presence of gases or vapours of various kinds and/or gases or vapours containing particulate (aerosol sprays, fumes, mists, etc.) combined filters are required. Respiratory protection devices must be used if the technical measures adopted are not suitable for restricting the worker's exposure to the threshold values considered. The protection provided by masks is in any case limited. If the substance considered is odourless or its olfactory threshold is higher than the corresponding TLV-TWA and in the case of an emergency, wear open-circuit compressed air breathing apparatus (in compliance with standard EN 137) or external air-intake breathing apparatus (in compliance with standard EN 138). For a correct choice of respiratory protection device, see standard EN 529.								
ENVIRONMENTAL EXPOSURE CONTROLS The emissions generated by manufacturing processes, including those generated by ventilation equipment, should be checked to ensure compliance with environmental standards.								
PROPAN-2-OL IBE (Biological Indicators of Exposure - ACGIH 2020): acetone in urine = 40 mg / L (end of shift)								
SECTION 9. Phy	sical and chem	ical prop	erties					
9.1. Information on ba	sic physical and cher	nical proper	ties					
Appearance		liquid						
Colour		colourless						

Steeles Miele Group	STEELCO S	.P.A.	Revision nr. 4 Dated 24/06/2021
Menuber	Steelcoxid	e A	Printed on 29/06/2021 Page n. 8/22 Replaced revision:3 (Printed on: 18/05/2021)
Odour Odour threshold pH Melting point / freezing point Initial boiling point Boiling range Flash point Evaporation rate Flammability (solid, gas) Lower inflammability limit Upper inflammability limit Upper inflammability limit Upper explosive limit Upper explosive limit Vapour pressure Vapour density Relative density Solubility Partition coefficient: n-octanol/water Auto-ignition temperature Decomposition temperature Viscosity Explosive properties Oxidising properties	pungent Not determined 1,6 Not available 105 °C Not available 15 °C Not determined not applicable Not determined Not determined Not determined Not determined Not determined Not determined 0,986 g/L soluble in water 0,05 (Propan-2-olo) Not available Not determined Not determined	Reason for missing data	a:Sulla base dello
9.2. Other information			
VOC (Directive 2010/75/EC) :	25,21 % - 240,23 g/litre		

# **SECTION 10. Stability and reactivity**

## 10.1. Reactivity

There are no particular risks of reaction with other substances in normal conditions of use.

ACETIC ACID It is highly corrosive: it attacks the most common metals with the development of hydrogen.

## 10.2. Chemical stability

The product is stable in normal conditions of use and storage.

ACETIC ACID Hygroscopic. Polymerises in contact with acetic aldehyde.

## 10.3. Possibility of hazardous reactions

The vapours may also form explosive mixtures with the air.

|--|

Steelcoxide A

Revision nr. 4 Dated 24/06/2021

Printed on 29/06/2021

Page n. 9/22 Replaced revision:3 (Printed on: 18/05/2021

#### ACETIC ACID

Vicela

Risk of explosion on contact with: chromium (VI) oxide,potassium permanganate,sodium peroxide,perchloric acid,phosphorus chloride,hydrogen peroxide.May react dangerously with: alcohols,bromine pentafluoride,chlorosulphuric acid,dichromate-sulphuric acid,ethane diamine,ethylene glycol,potassiun hydroxide,strong bases,sodium hydroxide,strong oxidising agents,nitric acid,ammonium nitrate,potassium tert-butoxide,oleum.Forms explosive mixtures with: air. Violent reactions in case of contact with strong bases, oxidants and other reactive compounds.

Formation of hydrogen in case of contact with carbon steel and non-noble metals.

Group

Miele

#### ETHANOL

Risk of explosion on contact with: alkaline metals,alkaline oxides,calcium hypochlorite,sulphur monofluoride,acetic anhydride,acids,concentrated hydrogen peroxide,perchlorates,perchloric acid,perchloronitrile,mercury nitrate,nitric acid,silver,silver nitrate,ammonia,silver oxide,ammonia,strong oxidising agents,nitrogen dioxide.May react dangerously with: bromoacetylene,chlorine acetylene,bromine trifluoride,chromium trioxide,chromyl chloride,fluorine,potassium tert-butoxide,lithium hydride,phosphorus trioxide,black platinum,zirconium (IV) chloride,zirconium (IV) iodide.Forms explosive mixtures with: air.

#### 10.4. Conditions to avoid

Avoid overheating. Avoid bunching of electrostatic charges. Avoid all sources of ignition.

PROPAN-2-OL

Heat, flames and sparks. Extreme temperatures and direct sunlight.

ACETIC ACID Avoid exposure to: sources of heat,naked flames. No ventilation. Open flames, heating and sparks. Humidity.

ETHANOL

Avoid exposure to: sources of heat,naked flames. Avoid high temperatures and proximity to ignition sources

#### 10.5. Incompatible materials

PROPAN-2-OL Oxidizing agents, acid anhydrides, aluminum, halogenated compounds, acids.

#### ACETIC ACID

Incompatible with: carbonates,hydroxides,phosphates,oxidising substances,bases. Bases, acetic anhydride, nitric acid, alcohol, halogens, halogen compounds and oxidizing materials.

#### ETHANOL

Strong mineral acids, oxidizing agents. High temperature aluminum.

## 10.6. Hazardous decomposition products

In the event of thermal decomposition or fire, gases and vapours that are potentially dangerous to health may be released.

## ACETIC ACID

When heated to decomposition, it develops irritating fumes. Carbon monoxide, carbon dioxide (CO2)

## **SECTION 11. Toxicological information**

In the absence of experimental data for the product itself, health hazards are evaluated according to the properties of the substances it contains, using the criteria specified in the applicable regulation for classification.

It is therefore necessary to take into account the concentration of the individual hazardous substances indicated in section 3, to evaluate the toxicological effects of exposure to the product.

Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021					
Miele Member	<b>.</b>	Printed on 29/06/2021					
	Steelcoxide A	Page n. 10/22 Replaced revision:3 (Printed on: 18/05/2021)					
11.1. Information on toxicological effects							
Metabolism, toxicokinetics, mechanism of action	Metabolism, toxicokinetics, mechanism of action and other information						
PROPAN-2-OL It is readily absorbed following inhalation exposure and rapidly spreads to tissues. However, it is also readily excreted in the urine, essentially in the form of the 2-methoxyacetic acid metabolite. (Arch Toxicol, 68, -588-94 - Johanson G, 1994)							
ACETIC ACID The substance is absorbed from the gastrointestinal tract and the lungs. It is rapidly distributed throughout the body. It is almost completely metabolized at the cellular level. After reaction with acetyl coenzyme A, the acetic acid is transformed by the Krebs cycle and is incorporated into the lipids and proteins, part of it is transformed into formic acid. Only a small part is found in the urine in unchanged form							
Information on likely routes of exposure							
ACETIC ACID The main routes of potential exposure substance.	are expected to be inhalation and skin contact in wo	rkers exposed to the manufacture and use of the					
Delayed and immediate effects as well as chroni	c effects from short and long-term exposure						
Information not available							
Interactive effects							
Information not available							
ACUTE TOXICITY							
ATE (Inhalation) of the mixture: Not classified (n ATE (Oral) of the mixture: Not classified (no sigr ATE (Dermal) of the mixture: Not classified (no s	o significant component) ificant component) ignificant component)						
ETHANOL LD50 (Oral) 1187 mg/kg Ratto LC50 (Inhalation) 115,9 mg/l/4h							
ACETIC ACID LD50 (Oral) 3310 mg/kg Ratto LD50 (Dermal) 1060 mg/kg Coniglio LC50 (Inhalation) 11,4 mg/l/4h Ratto							
PROPAN-2-OL LD50 (Oral) 4710 mg/kg Rat LD50 (Dermal) 12800 mg/kg Rat LC50 (Inhalation) 72,6 mg/l/4h Rat							
PROPAN-2-OL							
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 2 Species: Rat (Sherman) Routes of exposure: oral Results: LD50 = 5840 mg / kg	01						
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 2 Species: Rabbit Routes of exposure: cutaneous	02						

	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Miles Member	Steelcoxide A	Printed on 29/06/2021
		Page n. 11/22 Replaced revision:3 (Printed on: 18/05/2021)
Results: LC50 = 16.4 ml / kg		
Method: equivalent or similar to OECD 40 Reliability (Klimisch score): 1 Species: Rat (Fischer 344; Male / Female Routes of exposure: inhalation (vapors) Results: LD50> 10000 ppm / 6h	)3 ?)	
ACETIC ACID		
Bibliographic reference: "The acute oral t 82 (1941))" Reliability (Klimisch score): 2 Species: rat (Male / Female) Routes of exposure: oral Results: LD50: 3310 mg / kg	oxicity of acetic, chloroacetic, dichloroacetic and tr	richloroacetic acids (J Ind Hyg Toxicol, Vol 23, PP 78-
Method: equivalent or similar to OECD 40 Reliability (Klimisch score): 2 Species: rat (Sprague-Dawley Male / Fen Routes of exposure: inhalation (vapors) Results: LC50: 11.4 mg / I 4h	)3 nale)	
Reference: [Lewis, R.J. Sr. (ed) Sax's Da NJ. 2004., p. 16] Species: Rabbit Routes of exposure: cutaneous Results: LD50: 1060 mg / kg	ngerous Properties of Industrial Materials. 11th Ed	lition. Wiley-Interscience, Wiley & Sons, Inc. Hoboken,
ETHANOL		
Method: OECD 401 Reliability (Klimish score): 1 Species: rat (Cox CD; Male / Female) Route of exposure: oral Results: LD50: 10470 mg / kg		
Method: OECD 403 Reliability (Klimish score): 2 Species: rat (Sprague-Dawley; Male / Fer Route of Exposure: inhalation (vapors) LC50 results (male): 116.9 mg / I 4h	nale)	
Reference: Schechter, M. et al, Pharmace Reliability (Klimisch score): 2 Species: Mouse (HS; male / female) Routes of exposure: intraperitoneal Results: LD50 = 9450 mg / kg body weigh	ol Biochem Behav 52 (1): 245-248, 1995 nt	
SKIN CORROSION / IRRITATION		
Does not meet the classification criteria for this ha	azard class	
PROPAN-2-OL		
Reliability (Klimisch score): 2 Species: Rabbit Routes of exposure: cutaneous		
Results: Not irritating Reference: Nixon G et al, Toxicology and	Applied Pharmacology 31, 481-490 (1975)	
ACETIC ACID		

	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Member	Steelcoxide A	Printed on 29/06/2021
		Page n. 12/22 Replaced revision:3 (Printed on: 18/05/2021)
Method: equivalent or similar to OECD 4 Reliability (Klimisch score: 2 Species: rabbit Routes of exposure: cutaneous Results: Based on the strength of eviden skin.	04 ce of the available data as determined by expert judgme	ent, the substance is classified as corrosive to the
ETHANOL		
Method: OECD 404 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White) Routes of exposure: cutaneous Results: non-irritating.		
SERIOUS EYE DAMAGE / IRRITATION		
Causes serious eye irritation		
PROPAN-2-OL		
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White) Routes of exposure: ocular Results: irritating	05	
ACETIC ACID		
Method: equivalent or similar to OECD 44 Reliability (Klimisch score): 2 Species: rabbit (Rsk: NZW) Routes of exposure: ocular Results: Based on the strength of evider eyes.	05 nce of the available data as determined by expert judg	ment, the substance is classified as corrosive to
ETHANOL		
Method: OECD 405 Reliability (Klimisch score): 2 Species: Rabbit Routes of exposure: ocular Results: irritating.		
RESPIRATORY OR SKIN SENSITISATION		
Does not meet the classification criteria for this h	azard class	
PROPAN-2-OL		
Method: OECD 406 Reliability (Klimisch score): 1 Species: Guinea pig (Dunkin-Hurtley; Ma Routes of exposure: cutaneous Results: not sensitizing	le / Female)	
ACETIC ACID		

Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Miele Member		Printed on 29/06/2021
	Steelcoxide A	Page n. 13/22
		Replaced revision:3 (Printed on: 18/05/2021)
Based on the strength of evidence of the ava skin sensitization hazard class.	ailable data as determined by expert judgm	ent, the substance is not classified for the respiratory o
ETHANOL		
Method: equivalent or similar OECD 406 Reliability (Klimisch score): 2 Species: Guinea pig (Pirbright White; Female Routes of exposure: cutaneous Results: not sensitizing.	)	
RM CELL MUTAGENICITY		
es not meet the classification criteria for this hazard	d class	
PROPAN-2-OL		
Based on available data, the substance has n	o mutagenic effects and is not classified un	der the relevant hazard class CLP.
ACETIC ACID		
Method: equivalent or similar to OECD 471 - I Reliability (Klimisch score): 2 Species: Salmonella Typhimurium Results: negative with and without metabolic a	In vitro test activation.	
Method: equivalent or similar to EU B.12 - live Reliability (Klimisch score): 1 Species: rat (CD (Sprague-Dawley) Male / Fe Routes of exposure: inhalation (vapors) Results: negative.	e test male)	
ETHANOL		
Method: equivalent or similar to OECD 471 - I Reliability (Klimisch score): 1 Species: S. typhimurium Results: negative with and without metabolic a Method: equivalent or similar to OECD 474 - I Species: mouse (NMRI; Male / Female) Routes of exposure: intraperitoneal Results: negative.	In vitro test activation In vivo test	
RCINOGENICITY		
es not meet the classification criteria for this hazard	d class	
PROPAN-2-OL Based on available data, the substance has n	o carcinogenic effects and is not classified	under the relevant hazard class CLP.
ACETIC ACID		
Bibliographic reference: "The stimulating effe induced by n-nitrososarcosin ethyl ester in rat Reliability (Klimisch score): 2 Species: rat (outbred white; male) Routes of exposure: oral	ect of acetic acid, alcohol and thermal burn s" (Cancer Letters Vol 47, pp179-185 (1989	n injury on oesophagus and forestomach carcinogenes )))

	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Miele Member		Printed on 29/06/2021
	Steelcoxide A	Page n. 14/22 Replaced revision:3 (Printed on: 18/05/2021)
ETHANOL		
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: rat (Fischer 344 / DuCrj; Male / Routes of exposure: inhalation (vapors) Results: negative.	53 Female)	
REPRODUCTIVE TOXICITY		
Does not meet the classification criteria for this h	azard class	
PROPAN-2-OL		
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: Rat (Sprague-Dawley; Male / F Routes of exposure: Oral Results: negative. NOAEL = 1000 mg / F	:16 iemale) xg bw / day.	
Adverse effects on sexual function and fertility		
PROPAN-2-OL		
Method: equivalent or similar to OECD 4 Reliability (Klimsch score): 1 Species: rat (Sprague-Dawley Male / Fe Routes of exposure: oral Results: negative.	n16 male)	
ETHANOL		
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: mouse (CD-1; Male / Female) Routes of exposure: oral Results: No effect on fertility at doses eq	16 juivalent to 20.7 g / kg / day	
Adverse effects on development of the offspring		
PROPAN-2-OL		
Method: equivalent or similar to OECD 4 Reliability (Klimsch score): 1 Species: rat (Sprague-Dawley) Routes of exposure: oral Results: negative.	14	
ACETIC ACID		
Method: equivalent or similar to EU B.31 Reliability (Klimisch score): 2 Species: rat (Wistar) Routes of exposure: oral Results: negative. NOAEL (developmental toxicity): 1 600 r	ng / kg body weight / day.	
ETHANOL		

Member       Steelcoxide A       Printed on 2906/2021         Page n. 15/22       Page n. 15/22         Replaced revision 3 (Printed on: 18/05/2021)         Method: equivalent or similar to OECD 414         Reliability (Klimisch score): 2         Species: rat (Sprague-Dawley)         Routes of exposure: inhalation         Results: negative. NOAEL (maternal) = 16000 ppm.         NOAEL (fetus)> = 20,000 ppm         STOT - SINGLE EXPOSURE         May cause drowsiness or dizziness         PROPAN-2-OL         Metod: OECD 426         Affidabilità (Klimisch score): 1         Specie: Ratio (Sprague-Dawley; Femmina)         Vie d'esposizione: orale.         Risultati: Può provocare sonnolenza o vertigini.         In base ai dati disponibili, la sostanza presenta effetti di tossicità specifica per organi bersaglio per esposizione singola ed è classificata sotto		STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Steelcoxide A       Page n. 15/22 Replaced revision:3 (Printed on: 18/05/2021)         Method: equivalent or similar to OECD 414 Reliability (Klimisch score): 2 Species: rat (Sprague-Dawley) Routes of exposure: inhalation Results: negative. NOAEL (maternal) = 16000 ppm. NOAEL (fetus)> = 20,000 ppm         STOT - SINGLE EXPOSURE         May cause drowsiness or dizziness         PROPAN-2-OL         Metodo: OECD 426 Affidabilită (Klimisch score): 1 Specie: Ratio (Sprague-Dawley; Femmina) Vie d'esposizione: orale. Risultati: Può provocare sonnolenza o vertigini. In base ai dati disponibili, la sostanza presenta effetti di tossicità specifica per organi bersaglio per esposizione singola ed è classificata sotto	Miele Member	r	Printed on 29/06/2021
Method: equivalent or similar to OECD 414         Reliability (Klimisch score): 2         Species: rat (Sprague-Dawley)         Routes of exposure: inhalation         Results: negative. NOAEL (maternal) = 16000 ppm.         NOAEL (fetus) > = 20,000 ppm         STOT - SINGLE EXPOSURE         May cause drowsiness or dizziness         PROPAN-2-OL         Metodo: OECD 426         Affidabilità (Klimisch score): 1         Specie: Ratio (Sprague-Dawley; Femmina)         Vie d'esposizione: orale.         Risultati: Può provocare sonnolenza o vertigini.         In base ai dati disponibili, la sostanza presenta effetti di tossicità specifica per organi bersaglio per esposizione singola ed è classificata sotto		Steelcoxide A	Page n. 15/22 Replaced revision:3 (Printed on: 18/05/2021)
relativa classe di pericolo CLP. ACETIC ACID In humans, serious effects are reported following accidental single exposures by any route, mainly due to the local corrosive action of the substance with consequent systemic effects (INRS, 2011). Inhalation exposure to vapors or aerosols immediately causes irritation symptoms. ETHANOL Based on available data, the substance has no specific target organ toxicity effects for single exposure and is not classified under the releva CLP hazard class. STOT - REPEATED EXPOSURE Does not meet the classification criteria for this hazard class PROPAN-2-OL Accordi to available data, this substance does not present specific organ toxicity for repeated exposition and it is not classified under the releva CLP hazard class. ACETIC ACID Based on available data, the substance has no specific target organ toxicity effects on repeated exposure and is not classified under the releva CLP hazard class. ACETIC ACID Based on available data, the substance has no specific target organ toxicity effects on repeated exposure and is not classified under the releva CLP hazard class. Bibliographic reference: Antihypertensive effects of acetic acid and vinegar on spontaneously hypertensive rats (Biosci. Biotechnol. Biochem V 65, (12), pp 2690-2694 (2001)) Reliability (Klmisch score): 2 Species: rat (Male) Routes of exposure: oral Bibliographic reference: Acetic acid, a potent stimulator of mouse epidermal macromolecular synthesis and hyperplasia but with weak turco promoting abilit (Nat. Cancer Ints., Vol 55, pp 983-987 (1975)) Reliability (Klmisch score): 2 Species: mouse (CD I Female) Routes of exposure: oral and classified of mouse epidermal macromolecular synthesis and hyperplasia but with weak turco promoting ability (Nat. Cancer Ints., Vol 55, pp 983-987 (1975)) Reliability reguitive. NOAEL: 20 mg / kg body weight / day Bibliographic reference: Acetic acid, a potent stimulator of mouse epidermal macromolecular synthesis and hyperplasia but with weak turco promoting ability	Method: equivalent or similar to OECD Reliability (Klimisch score): 2 Species: rat (Sprague-Dawley) Routes of exposure: inhalation Results: negative. NOAEL (maternal) = NOAEL (fetus)> = 20,000 ppm STOT - SINGLE EXPOSURE May cause drowsiness or dizziness PROPAN-2-OL Metodo: OECD 426 Affidabilità (Klimisch score): 1 Specie: Ratto (Sprague-Dawley; Femm Vie d'esposizione: orale. Risultati: Può provocare sonnolenza o v In base ai dati disponibili, la sostanza p relativa classe di pericolo CLP. ACETIC ACID In humans, serious effects are reports substance with consequent systemic ef ETHANOL Based on available data, the substanc CLP hazard class. STOT - REPEATED EXPOSURE Does not meet the classification criteria for this PROPAN-2-OL Accordi to available data, this substanc CLP hazard class. ACETIC ACID Based on available data, this substanc CLP hazard class. ACETIC ACID Based on available data, the substanc CLP hazard class. ACETIC ACID Based on available data, the substanc CLP hazard class. ACETIC ACID Based on available data, the substance CLP hazard class. Bibliographic reference: Antihypertensii 65, (12), pp 2690-2694 (2001)) Reliability (Klimisch score): 2 Species: rat (Male) Routes of exposure: oral Results: negative. NOAEL: 290 mg / kg Bibliographic reference: Acetic acid, a promoting ability (NAE. Cancer Inst., Vo Reliability (Klimisch score): 2 Species: mouse (CD-1 Female) Routes of exposure: cutaneous Results: negative. NOAEL: 30 mg / anii	Steelcoxide A 414 16000 ppm. ina) vertigini. presenta effetti di tossicità specifica per organi bersaglia ed following accidental single exposures by any route fects (INRS, 2011). Inhalation exposure to vapors or aei e has no specific target organ toxicity effects for single hazard class e does not present specific organ toxicity for repeated e a has no specific target organ toxicity effects on repeate ve effects of acetic acid and vinegar on spontaneously g body weight / day potent stimulator of mouse epidermal macromolecula 15, pp 983-987 (1975)) mal	Page n. 15/22 Replaced revision:3 (Printed on: 18/05/2021)

Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Miele Member		Printed on 29/06/2021
	Steelcoxide A	Page n. 16/22
		Replaced revision:3 (Printed on: 18/05/2021)
ETHANOL		
Method: equivalent or similar OECD 408 Reliability (Klimisch score): 2 Species: Rat (Sprague-Dawley; Male / F Routes of exposure: oral Results: negative. NOAEL: 1730 mg / kg	3 Female) g body weight / day	
ASPIRATION HAZARD		
Does not meet the classification criteria for this h	nazard class	
PROPAN-2-OL		
No data are available on the hazard in case of a	spiration.	
ACETIC ACID		
No data are available on the hazard in case of a	spiration.	
ETHANOL		
No data are available on the hazard in case of a	spiration.	
	- 41	
SECTION 12. Ecological information	ation	
Use this product according to good working contaminate soil or vegetation.	practices. Avoid littering. Inform the competent aut	horities, should the product reach waterways or
12.1. Toxicity		
ETHANOL		
LC50 - for Fish	14,2 g/l/96h Pimephales promelas (L	JS EPA E03-05)
EC50 - for Crustacea	5012 mg/l/48h Ceriodaphnia dubia (A	STM E729-80)
EC50 - for Algae / Aquatic Plants	275 mg/l/72h Chlorella vulgaris (OEC	D 201)
Chronic NOEC for Fish	250 mg/L/5 d	
ACETIC ACID		
LC50 - for Fish	> 300,82 mg/l/96h Oncorhynchus my	kiss (equivalent or similar to OECD 203)
EC50 - for Crustacea	> 300,82 mg/l/48h Daphnia magna (C	DECD 202)
EC50 - for Algae / Aquatic Plants	> 300,82 mg/l/72h Skeletonema costa	atum (equivalent or similar to ISO 10253)
Chronic NOEC for Algae / Aquatic Plants	> 300,82 mg/l/72h Skeletonema costa	atum (equivalent or similar to ISO 10253)
PROPAN-2-OL		
EC50 - for Crustacea	> 10000 mg/l/48h Daphnia magna (ed	quivalent or similar to OECD 203)
12.2. Persistence and degradability		
PROPAN-2-OL Easily biodegradable, 53% in 5 giorni (equivale	ent or similar to EU C.5)	

Steelco Miele Group Member	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
		Printed on 29/06/2021
	Steelcoxide A	Page n. 17/22 Replaced revision:3 (Printed on: 18/05/2021)
ACETIC ACID Easily biodegradable, 96% in 20 days (Price, K.S., Waggy, G.T. And Conway, R.A. 1974, J. Water pollut. Contr. Fed. Vol 46 PP 46-77)		
Solubility in water > 10000 mg/l		
ETHANOL Readily biodegradable, 60% in 10 days (BOD - Standard methods for the examination of water and waste water 1971. 13th ed, American Public Health Assoc, NY)		
12.3. Bioaccumulative potential		
ETHANOL		
Partition coefficient: n-octanol/water	-0,35 Log Kow 24°C (OECD 107)	
ACETIC ACID		
Partition coefficient: n-octanol/water	-0,17 Log Kow (CRC Press Inc. Boca Raton. USA.)	
BCF	3,16 (Q)sar (Meylan,WM, Howard,PH, Boethling,RS et al. 1999)	
PROPAN-2-OL		
Partition coefficient: n-octanol/water	0,05	
12.4. Mobility in soil		
ACETIC ACID Based on an estimated Koc value of approx. 1, very high mobility in soil is expected. Acetic acid volatilizes from the dry soil surface according to the vapor pressure value. Based on the pKa value of 4.76, it mainly exists in anionic form in water. The anionic form adsorbs less than the corresponding neutral form in soils containing organic carbon and clay (HSDB, 2018).		
ACETIC ACID		
Partition coefficient: soil/water	1,153	
12.5. Results of PBT and vPvB assessment		
On the basis of available data, the product does not contain any PBT or vPvB in percentage ≥ than 0,1%.		
12.6. Other adverse effects		
Information not available		
SECTION 13. Disposal considerations		
13.1. Waste treatment methods		
Reuse if possible. Product residues are to be considered special hazardous waste. The dangerousness of the waste that partially contains this product must be evaluated according to the laws in force. (Ref. Annex D - Part IV of Legislative Decree no. 152/2006 and subsequent amendments and adjustments).		

Disposal must be entrusted to an authorized waste management company, in compliance with national and possibly local regulations. The legal responsibility for disposal lies with the producer / holder of the waste. Different CER (European Waste Code) codes could be applied to this mixture according to the specific circumstances that generated the waste, any alterations and contaminations.

The product as it is, out of specification in the original packaging, or poured into a suitable container for disposal as waste, or the product in specification but no longer usable (for example following an accidental spill), is to be classified with a code CER compatible with the description of use indicated in section 1.2.
Declar	Minto Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Ciceico	Member	Steelcoxide A	Printed on 29/06/2021
			Page n. 18/22 Replaced revision:3 (Printed on: 18/05/2021)

The appropriate final destination of the waste will be assessed by the manufacturer according to the chemical-physical characteristics of the waste itself compatible with the authorized plant to which it will be conferred for recovery, treatment or final disposal in the manner prescribed by current regulations. Disposal via the wastewater drain is not permitted.

For hazardous substances registered according to EC Regulation 1907/2006 (REACH) for which a chemical safety report has been prepared, refer to the specific information contained in the exposure scenarios attached to this SDS.

### CONTAMINATED PACKAGING

Contaminated packaging must be sent, properly labeled, for recovery or disposal in compliance with national regulations on waste management and must be classified with the following EWC code:

15 01 10 \*: packaging containing residues of dangerous substances or contaminated by these substances

# **SECTION 14. Transport information**

#### 14.1. UN number

ADR / RID, IMDG, IATA: 1987

# 14.2. UN proper shipping name

ADR / RID:	ALCOHOLS, N.O.S. MIXTURE (Propan-2-OI, Ethanol)
IMDG:	ALCOHOLS, N.O.S. MIXTURE (Propan-2-OI, Ethanol)
IATA:	ALCOHOLS, N.O.S. MIXTURE (Propan-2-OI, Ethanol)

## 14.3. Transport hazard class(es)

ADR / RID:	Class: 3	Label: 3
IMDG:	Class: 3	Label: 3
IATA:	Class: 3	Label: 3



## 14.4. Packing group

ADR / RID, IMDG, II IATA:

# 14.5. Environmental hazards

ADR / RID:	NO
IMDG:	NO
IATA:	NO

For Air transport, environmentally hazardous mark is only mandatory for UN 3077 and UN 3082.

# 14.6. Special precautions for user

ADR / RID:	HIN - Kemler: 33	Limited Quantities: 1 L	Tunnel restriction code: (D/E)
	Special provision: 640D		
IMDG:	EMS: <u>F-E, S-E</u>	Limited Quantities: 1 L	

		STEELCO S.P.	.A.	Revision nr. 4 Dated 24/06/20	21
teelco	Miele Group Member	• • • •		Printed on 29/00	5/2021
		Steelcoxide	A	Page n. 19/22 Replaced revisi	on:3 (Printed on: 18/05/2021)
IATA:	Car Pas Spe	go: ss.: ecial provision:	Maximum quanti Maximum quanti A3	ity: 60 L ity: 5 L	Packaging instructions: 364 Packaging instructions: 353
14.7. Transport in bulk acc	ording to Annex II	of Marpol and the IBC Code			
Information not relevant					
SECTION 15. Regu	ulatory inform	ation			
15.1. Safety, health and e	nvironmental regu	lations/legislation specific for the s	substance or mix	ture	
Seveso Category - Directive	2012/18/EC: P5c				
Restrictions relating to the p	oduct or contained	substances pursuant to Annex XVII to	EC Regulation 19	07/2006	
Product					
Point	3				
	in Annex I of Regu a) hazard classes 2, 2.15 types A to b) hazard classes narcotic effects, 3. c) hazard class 4.1 d) hazard class 5.1	lation (EC) No. 1272/2008: 2.1 to 2.4, 2.6 and 2.7, 2.8 types A ar F; 3.1 to 3.6, 3.7 adverse effects on 5 9 and 3.10; ; 1.	nd B, 2.9, 2.10, 2.1 sexual function an	2, 2.13 categ	ories 1 and 2, 2.14 categories 1 and development, 3.8 effects other than
Point	40				
	Substances classin category 1 or 2, s category 1 pyroph 1272/2008Contain	fied as flammable gases of category substances and mixtures which, in o pric liquids or category 1 pyrophoric s ed substance	1 or 2, flammable contact with water colids, even if not li	liquids of cat , emit flamm isted in Part 3	egory 1, 2 or 3, flammable solids of able gases of categories 1, 2 or 3, of Annex VI to Regulation (EC) No.
Point	75	ACETIC ACID Reg. no.: 01-211947	5328-30-xxxx		
	Substances includ a) substances clas - category 1A, 1B due to effects follo - Reproductive tox inhalation only; - skin sensitization - skin corrosion of - serious category b) substances liste columns g, h or i o d) substances liste this entry apply to referred to in point	ed in one or more of the following: sified in one of the following classes or 2 carcinogenicity, category 1A, 1 wing exposure by inhalation only; icity of category 1A, 1B or 2, but ex of category 1, 1A or 1B; category 1, 1A, 1B or 1C or skin irrita 1 eye damage or category 2 eye irrite d in Annex II of Regulation (EC) no f the table of this annex; ed in Appendix 13 of this annex. The all mixtures intended for tattooing p s a) to d) of this column and entry.	in Annex VI, part 3 B or 2 germ cell i ccluding substance tion of category 2; ation; 1223/2009 of the E . 1223/2009 for wi ancillary requiren ractices, regardles	, of Regulatio mutagenicity, es classified c furopean Parl hich a conditi nents referrec s of whether	n (EC) no. 1272/2008: but excluding substances classified due to effects following exposure by iament and of the Council (*); on is indicated in at least one of the to in points 7 and 8 of column 2 of they contain one of the substances
Regulation (EC) No. 2019/17	148 - on the marketir	ng and use of explosives precursors			
Not applicable					
Substances in Candidate Lis	t (Art. 59 REACH)				

Steelco Miele Group Member	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021			
	per	Printed on 29/06/2021			
	Steelcoxide A	Page n. 20/22			
		Replaced revision:3 (Printed on: 18/05/2021)			
On the basis of available data, the product do	es not contain any SVHC in percentage ≥ than 0,1% <b>.</b>				
Substances subject to authorisation (Annex )	IV REACH)				
None					
Substances subject to exportation reporting p	ursuant to (EC) Reg. 649/2012:				
None					
Substances subject to the Rotterdam Conver	tion:				
None					
Substances subject to the Stockholm Conver	tion:				
None					
Healthcare controls					
Workers exposed to this chemical agent mus workers' health and safety are modest and the	t not undergo health checks, provided that available risk-a at the 98/24/EC directive is respected.	assessment data prove that the risks related to the			
Medical device, class IIb according to legisla	on 93/42.				
D.Lgs. 152/2006 and following modifications					
Emissions according to Chapter V, Annex I:					
TAB. D Class 3 04.97 %					
TAB. D Class 4 21,22 % TAB. D Class 5 03,73 %					
15.2. Chemical safety assessment					
A chemical safety assessment has not been	performed for the preparation/for the substances indicated	in section 3.			
SECTION 16. Other information	SECTION 16. Other information				
Text of hazard (H) indications mentioned in s	ection 2-3 of the sheet:				
Flam, Liq. 2 Flammable liqui	category 2				
Flam Ling 3 Elammable liqui	1, category 2				
Skin Corr. 1A Skin correction	$a_1 = a_1 = a_2 = a_1 = a_2$				
Eve Irrit 2					
Eye initiation, ca	Crocific terret error terrisity structure estimates				
	H225				
	Thy manimum equid and vapour.				
	s and vapour. skin hurns and eve damage				
Causes severe	skin burns and eye damage.				

- H319 Causes serious eye irritation.
- H336 May cause drowsiness or dizziness.

Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021
Miele Member		Printed on 29/06/2021
	Steelcoxide A	Page n. 21/22 Replaced revision:3 (Printed on: 18/05/2021)
LEGEND:		
<ul> <li>ADR: European Agreement concerning the car</li> <li>CAS NUMBER: Chemical Abstract Service Nui</li> <li>CE50: Effective concentration (required to indu</li> <li>CE NUMBER: Identifier in ESIS (European arc</li> <li>CLP: EC Regulation 1272/2008</li> <li>DNEL: Derived No Effect Level</li> <li>EmS: Emergency Schedule</li> <li>GHS: Globally Harmonized System of classific</li> <li>IATA DGR: International Air Transport Associa</li> <li>IC50: Immobilization Concentration 50%</li> <li>IMDG: International Maritime Code for dangerd</li> <li>IMO: International Maritime Code for dangerd</li> <li>IMO: International Maritime Code for dangerd</li> <li>IMO: International Maritime Organization</li> <li>INDEX NUMBER: Identifier in Annex VI of CLF</li> <li>LC50: Lethal Concentration 50%</li> <li>OEL: Occupational Exposure Level</li> <li>PBT: Persistent bioaccumulative and toxic as F</li> <li>PEC: Predicted environmental Concentration</li> <li>PEL: Predicted no effect concentration</li> <li>REL: EC Regulation 1907/2006</li> <li>RID: Regulation concerning the international transport Associal t</li></ul>	riage of Dangerous goods by Road mber uce a 50% effect) hive of existing substances) ation and labeling of chemicals tion Dangerous Goods Regulation ous goods REACH Regulation ansport of dangerous goods by train e exceeded during any time of occupational exposure. <i>re</i> as for REACH Regulation	
GENERAL BIBLIOGRAPHY		
<ol> <li>Regulation (EC) 1907/2006 (REACH) of the E</li> <li>Regulation (EC) 1272/2008 (CLP) of the Euro</li> <li>Regulation (EU) 790/2009 (I Atp. CLP) of the</li> <li>Regulation (EU) 2015/830 of the European Pa</li> <li>Regulation (EU) 286/2011 (II Atp. CLP) of the</li> <li>Regulation (EU) 286/2011 (II Atp. CLP) of the</li> <li>Regulation (EU) 286/2013 (IV Atp. CLP) of the</li> <li>Regulation (EU) 944/2013 (V Atp. CLP) of the</li> <li>Regulation (EU) 2015/1221 (VII Atp. CLP) of the</li> <li>Regulation (EU) 2015/1221 (VII Atp. CLP) of the</li> <li>Regulation (EU) 2016/1179 (IX Atp. CLP) of</li> <li>Regulation (EU) 2016/1179 (IX Atp. CLP)</li> <li>Regulation (EU) 2018/669 (XI Atp. CLP)</li> <li>Regulation (EU) 2018/1480 (XIII Atp. CLP)</li> <li>Regulation (EU) 2019/521 (XII Atp. CLP)</li> <li>Regulation (EU) 2019/1148</li> <li>Regulation (EU) 2020/217 (XIV Atp. CLP)</li> <li>The Merck Index 10th Edition</li> <li>Handling Chemical Safety</li> <li>INRS - Fiche Toxicologique (toxicological shee</li> <li>Patty - Industrial Hygiene and Toxicology</li> <li>N.I. Sax - Dangerous properties of Industrial M</li> <li>IFA GESTIS website</li> <li>ECHA website</li> <li>Database of SDS models for chemicals - Minis</li> </ol>	European Parliament pean Parliament European Parliament a European Parliament e European Parliament e European Parliament e European Parliament f the European Parliament the European Parliament the European Parliament tru e European Parliament tru of Health and ISS (Istituto Superiore di Sanità) – Italy	,
<b>Training for workers:</b> The training of workers must include contents, u membership, according to the procedures provid	pdates and duration according to the risk profiles assigned for by Legislative Decree 81/2008.	ed to the working sectors of

Group Group	STEELCO S.P.A.	Revision nr. 4 Dated 24/06/2021			
Member Member	Ctaslasvida A	Printed on 29/06/2021			
	SteelCoxide A	Page n. 22/22 Replaced revision:3 (Printed on: 18/05/2021)			
Classification of the mixture according to Degulation					
classification of the mixture according to Regulation	(EC) 11. 1272/2006 Class	incation procedure			
Flammable liquid, category 2H225Based on experimental data.		d on experimental data.			
Skin corrosion, category 1 H314 Based on experimental data.		d on experimental data.			
Specific target organ toxicity - single exposure category 3 H336 Based on literature data		d on literature data			

### Note for the recipient of the Safety Data Sheet (SDS):

It is the recipient of this SDS who must ensure that the information contained is read and understood by all persons who handle, store, use, or otherwise come into contact in any way with the substance or mixture to which this sheet refers. In particular, the recipient must provide adequate training to personnel assigned to the use of dangerous substances or mixtures. The recipient must ensure the suitability and completeness of the information in relation to the specific use of the substance or mixture.

However, the substance or mixture to which this SDS refers must not be used for uses other than those specified in section 1. No responsibility is assumed for improper uses. Since the use of the product does not fall under the direct control of the Supplier, it is the user's obligation to observe, under his own responsibility, the laws and regulations in force regarding national and Community hygiene and safety.

The information contained in this SDS is provided in good faith and is based on the current state of scientific and technical knowledge, at the revision date indicated, available from the Supplier indicated in section 1 of this sheet. The SDS should not be interpreted as a guarantee of any specific property of the substance or mixture. The information refers only to the substance or mixture specifically designated in section 1 and may not be valid for the substance or mixture used in combination with other materials or in other processes not specifically indicated in the text.

This version of the SDS supersedes all previous versions.

Changes from the previous revision.

Changes have been made to the following sections: 01/02/03/04/05/06/07/08/09/10/11/12/13/14/15/16.

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021		
Group Miele Group Member	Steelcoxide-DT	Printed on 31/05/2021 Page n. 1/23 Replaced revision:3 (Printed on: 15/09/2020)		
	Safety Data Sheet			
1.1. Product identifier	e substance/mixture and of the compar	ny/undertaking		
Product name	Steelcoxide-DT			
1.2. Relevant identified uses of the substa	nce or mixture and uses advised against			
Intended use Concentrated generation of	decontaminating disinfectant product. Compatible w btical fibers. FOR PROFESSIONAL USE ONLY.	ith thermolabile medical devices and latest		
Uses advised against: Uses other that	n those indicated.			
1.3 Details of the supplier of the safety da	ta sheet			
Name Full address District and Country	STEELCO S.p.A. Via Balegante, 27 31039 Riese Pio X (TV) ITALY			
	tel. +39 0423 7561			
9 II <b>7</b> 0 II I	fax +39 0423 755528			
responsible for the Safety Data Sheet	info@steelcogroup.com			
1.4. Emergency telephone number				
For urgent inquiries refer to	Centros de Orientação de Doentes Urgentes	(CODU): 800250250		
SECTION 2. Hazards identificat	on			
2.1. Classification of the substance or mixtu	re			
The product is classified as hazardous pursuant to the provisions set forth in (EC) Regulation 1272/2008 (CLP) (and subsequent amendments and supplements). The product thus requires a safety datasheet that complies with the provisions of (EU) Regulation 2015/830. Any additional information concerning the risks for health and/or the environment are given in sections 11 and 12 of this sheet.				
Hazard classification and indication:				
Flammable liquid, category 2	H225 Highly flammat	ole liquid and vapour.		
2.2. Label elements				
-lazard labelling pursuant to EC Regulation 1272/2008 (CLP) and subsequent amendments and supplements.				

		STEELCO S.P.	P.A. Revision nr. 4 Dated 21/05/2021
Steelco	Miele Group Member	Steelcoxide-D	Printed on 31/05/2021 Page n. 2/23 Replaced revision:3 (Printed on: 15/09/2020)
Hazard pictograms:			
Signal words:	Danger		
Hazard statements:			
H225 EUH208	Highly flammable lid Contains: SUBTILIS May produce an alle	uid and vapour. NNE argic reaction.	
Precautionary statements:			
P210 P233	Keep away from he Keep container tigh	at, hot surfaces, sparks, open flames a ly closed.	and other ignition sources. No smoking.
2.3. Other hazards			
On the basis of available d	lata, the product does	not contain any PBT or vPvB in percer	entage ≥ than 0,1%.
SECTION 3. Com	position/inform	ation on ingredients	
3.2. Mixtures			
Contains:			
Identification	x = Cor	c. % Classification 1272/2008	8 (CLP)
ETHANEDIOL CAS 107-21-1 EC 203-473-3 INDEX 603-027-00-1 Reg. no. 01-211945681	6 ≤ x < 6-28-xxx	10 Acute Tox. 4 H302, STOT	T RE 2 H373
PROPAN-2-OL CAS 67-63-0 EC 200-661-7 INDEX 603-117-00-0 Reg. no. 01-211945755	6≤ x < 58-25-XXXX	10 Flam. Liq. 2 H225, Eye Irr	rrit. 2 H319, STOT SE 3 H336
ETHANOL CAS 64-17-5 EC 200-578-6 INDEX 603-002-00-5	1 ≤ x •	5 3 Flam. Liq. 2 H225, Eye Irr	rrit. 2 H319

•		STEELCO S.P.A. Foup lember Steelcoxide-DT		Revision nr. 4 Dated 21/05/2021 Printed on 31/05/2021 Page n. 3/23 Replaced revision:3 (Printed on: 15/09/2020)	
Miele Miele	Group Member				
Reg. no. 01-2119457610-43-XXXX					
SUBTILISINE					
CAS 9014-01-1	0,2 ≤ x	c< 0,6	< 0,6 Acute Tox. 4 H302, Eye Dam. 1 H318, Skin Irrit. 2 H315, STOT SE 3 H335,		
EC 232-752-2			Resp. Sens. 1 H334, Aquatic Acute 1 H400 M	1=1, Aquatic Chronic 2 H411	
INDEX 647-012-00-8					
Reg. no. 01-2119480434-38-xxxx	Reg. no. 01-2119480434-38-xxxx				
Γhe full wording of hazard (H) phrases is given in section 16 of the sheet.					
SECTION 4. First aid measures					

# 4.1. Description of first aid measures

EYES: Remove contact lenses, if present. Wash immediately with plenty of water for at least 15 minutes, opening the eyelids fully. If problem persists, seek medical advice.

SKIN: Remove contaminated clothing. Wash immediately with plenty of water. If irritation persists, get medical advice/attention. Wash contaminated clothing before using it again.

INHALATION: Remove to open air. In the event of breathing difficulties, get medical advice/attention immediately.

INGESTION: Get medical advice/attention. Induce vomiting only if indicated by the doctor. Never give anything by mouth to an unconscious person, unless authorised by a doctor.

PROTECTION MEASURES FOR FIRST AID: for the PPE necessary for first aid interventions refer to section 8.2 of the present safety data sheet.

## 4.2. Most important symptoms and effects, both acute and delayed

Specific information on symptoms and effects caused by the product are unknown.

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically. In the event of an accident or discomfort, consult a doctor immediately (if possible show the instructions for use or the safety data sheet).

# **SECTION 5. Firefighting measures**

### 5.1. Extinguishing media

SUITABLE EXTINGUISHING EQUIPMENT

Extinguishing substances are: carbon dioxide, foam, chemical powder. For product loss or leakage that has not caught fire, water spray can be used to disperse flammable vapours and protect those trying to stem the leak.

UNSUITABLE EXTINGUISHING EQUIPMENT

Do not use jets of water. Water is not effective for putting out fires but can be used to cool containers exposed to flames to prevent explosions.

### 5.2. Special hazards arising from the substance or mixture

HAZARDS CAUSED BY EXPOSURE IN THE EVENT OF FIRE Excess pressure may form in containers exposed to fire at a risk of explosion. Do not breathe combustion products.

PROPAN-2-OL Carbon oxides.

5.3. Advice for firefighters

# STEELCO S.P.A.

Revision nr. 4 Dated 21/05/2021

Printed on 31/05/2021

Meelc

Miele Group Member

Steelcoxide-DT

Page n. 4/23 Replaced revision:3 (Printed on: 15/09/2020)

# GENERAL INFORMATION

Use jets of water to cool the containers to prevent product decomposition and the development of substances potentially hazardous for health. Always wear full fire prevention gear. Collect extinguishing water to prevent it from draining into the sewer system. Dispose of contaminated water used for extinction and the remains of the fire according to applicable regulations.

SPECIAL PROTECTIVE EQUIPMENT FOR FIRE-FIGHTERS

Normal fire fighting clothing i.e. fire kit (BS EN 469), gloves (BS EN 659) and boots (HO specification A29 and A30) in combination with self-contained open circuit positive pressure compressed air breathing apparatus (BS EN 137).

# **SECTION 6.** Accidental release measures

# 6.1. Personal precautions, protective equipment and emergency procedures

## Block the leakage if there is no hazard.

Wear suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing. These indications apply for both processing staff and those involved in emergency procedures. Send away individuals who are not suitably equipped. Use explosion-proof equipment. Eliminate all sources of ignition (cigarettes, flames, sparks, etc.) from the leakage site.

## 6.1.1 For those who do not intervene directly

Do not take any action involving any personal risk or without proper training. Evacuate the surrounding areas. Don't touch either walk on the spilled material.

Wear suitable protective equipment (including personal protective equipment referred to in section 8 of this Safety Data Sheet) to prevent contamination of skin, eyes and personal clothing. Wear an appropriate respirator when ventilation is inadequate. Do not inhale the mists / vapors / fumes. Avoid the dispersion of the product into the environment. Follow the appropriate internal procedures provided for

non-personnel authorized to intervene directly in the event of accidental release.

### 6.1.2 For those who intervene directly

Stop the leak if there is no danger.

Evacuate unauthorized personnel. Wear suitable protective equipment (see section 8 of this Safety Data Sheet).

Follow the appropriate internal procedures for authorized personnel. Isolate the danger area and deny entry. Ventilate enclosed spaces before to come in.

# 6.2. Environmental precautions

The product must not penetrate into the sewer system or come into contact with surface water or ground water.

# 6.3. Methods and material for containment and cleaning up

Collect the leaked product into a suitable container. Evaluate the compatibility of the container to be used, by checking section 10. Absorb the remainder with inert absorbent material. Make sure the leakage site is well aired. Contaminated material should be disposed of in compliance with the provisions set forth in point 13.

### 6.4. Reference to other sections

Any information on personal protection and disposal is given in sections 8 and 13.

# **SECTION 7. Handling and storage**

# 7.1. Precautions for safe handling

Keep away from heat, sparks and naked flames; do not smoke or use matches or lighters. Without adequate ventilation, vapours may accumulate at ground level and, if ignited, catch fire even at a distance, with the danger of backfire. Avoid bunching of electrostatic charges. When performing transfer operations involving large containers, connect to an earthing system and wear antistatic footwear. Vigorous stirring and flow through the tubes and equipment may cause the formation and accumulation of electrostatic charges. In order to avoid the risk of fires and explosions, never use compressed air when handling. Open containers with caution as they may be pressurised. Do not eat, drink or smoke during use. Avoid leakage of the product into the environment.

Steelco Miele M		STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021
	Miele Group Member		Printed on 31/05/2021
		Steelcoxide-DT	Page n. 5/23 Replaced revision:3 (Printed on: 15/09/2020)

# 7.2. Conditions for safe storage, including any incompatibilities

Store only in the original container. Store the containers sealed, in a well ventilated place, away from direct sunlight. Store in a cool and well ventilated place, keep far away from sources of heat, naked flames and sparks and other sources of ignition. Keep containers away from any incompatible materials, see section 10 for details.

# 7.3. Specific end use(s)

No use other than that indicated in section 1.2 of this safety data sheet.

# **SECTION 8. Exposure controls/personal protection**

# 8.1. Control parameters

Regulatory References:

DEU	Deutschland	Technischen Regeln für Gefahrstoffe (TRGS 900) - Liste der Arbeitsplatzgrenzwerte und Kurzzeitwerte. MAK- und BAT-Werte-Liste 2020, Ständige Senatskommission zur Prüfung gesundheitsschädlicher Arbeitsstoffe, Mitteilung 56
ESP	España	Límites de exposición profesional para agentes químicos en España 2019
FRA	France	Valeurs limites d'exposition professionnelle aux agents chimiques en France. ED 984 - INRS
ITA	Italia	Decreto Legislativo 9 Aprile 2008, n.81
GBR	United Kingdom	EH40/2005 Workplace exposure limits (Fourth Edition 2020)
EU	OEL EU	Directive (EU) 2019/1831; Directive (EU) 2019/130; Directive (EU) 2019/983; Directive (EU) 2017/2398;
		Directive (EU) 2017/164; Directive 2009/161/EU; Directive 2006/15/EC; Directive 2004/37/EC; Directive
		2000/39/EC; Directive 98/24/EC; Directive 91/322/EEC.
	TLV-ACGIH	ACGIH 2020

# ETHANEDIOL

Country	TWA/8h		STEL/15min		Remarks /	
	ma/m3	maa	ma/m3	maa	Observations	
		PP		PP		
DEU	26	10	52	20	SKIN	
DEU	26	10	52	20	SKIN	
ESP	52	20	104	40	SKIN	
FRA	52	20	104	40	SKIN	
ITA	52	20	104	40	SKIN	
GBR	52	20	104	40	SKIN	
EU	52	20	104	40	SKIN	
		25		50		
			10		INHAL	
ion - PNEC						
			10	mę	g/I	
			1	mę	g/l	
ittent release			10	mę	g/l	
Normal value of STP microorganisms			199,5	mę	g/l	
Normal value for the terrestrial compartment				mę	g/kg	
t level - DNEL /	DMEL					
Effects on				Effects on		
	Country DEU DEU ESP FRA ITA GBR EU ion - PNEC ittent release anisms compartment tt level - DNEL / Effects on consumers	Country       TWA/8h         mg/m3         DEU       26         DEU       26         ESP       52         FRA       52         ITA       52         GBR       52         EU       52         ion - PNEC	Country       TWA/8h         mg/m3       ppm         DEU       26       10         DEU       26       10         ESP       52       20         FRA       52       20         ITA       52       20         GBR       52       20         EU       52       20         ion - PNEC       25         ittent release	Country         TWA/8h         STEL/15min           mg/m3         ppm         mg/m3           DEU         26         10         52           DEU         26         10         52           ESP         52         20         104           FRA         52         20         104           ITA         52         20         104           GBR         52         20         104           EU         52         20         104           EU         52         20         104           EU         52         20         104           IO         25         10         1           itent release         10         1           ittent release         10         199,5           compartment         1,53         1,53	Country         TWA/8h         STEL/15min           mg/m3         ppm         mg/m3         ppm           DEU         26         10         52         20           DEU         26         10         52         20           ESP         52         20         104         40           FRA         52         20         104         40           ITA         52         20         104         40           GBR         52         20         104         40           EU         52         20         104         40           EU         52         20         104         40           IO         25         50         50         50           IO         10         mg         1         mg           itent release         10         mg         1         mg           anisms         199,5         mg         mg         mg           compartment         1,53         mg         mg         mg	Country         TWA/8h         STEL/15min         Remarks / Observations           DEU         26         10         52         20         SKIN           ESP         52         20         104         40         SKIN           FRA         52         20         104         40         SKIN           ITA         52         20         104         40         SKIN           GBR         52         20         104         40         SKIN           EU         52         20         104         40         SKIN           EU         52         20         104         40         SKIN           ion - PNEC         10         mg/l         Intervet         Intervet           ittent release         10         mg/l         Intervet         Intervet           icompartment         1,53         mg/kg         Intervet         Effects on workers

			STEELCO	) S.P.A.		Revision nr. 4 Dated 21/05/2021		
Steelco Miele Group Member			Steelcoxide-DT			Printed on 31/05/2021 Page n. 6/23 Replaced revision:3 (Printed on: 15/09/2020)		
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic	Acute local	Acute	Chronic local	Chronic
Oral		NPI		NPI		Systemic		Systemic
Inhalation	NPI	NPI	7 mg/m3	NPI	NPI	NPI	35 mg/m3	NPI
Skin	NPI	NPI	NPI	53 mg/kg bw/d	NPI	NPI	NPI	106 mg/kg bw/d
PROPAN-2-OL								
Type	Country	TWA/8h		STEL/15min		Remarks	;/	
		ma/m3	maa	ma/m3	maa	Observa	tions	
AGW	DEU	500	200	1000	400			
MAK	DEU	500	200	1000	400			
VLA	ESP	500	200	1000	400			
VLEP	FRA			980	400			
WEL	GBR	999	400	1250	500			
TI V-ACGIH	02.1	492	200	983	400			
Predicted no-effect concentra	tion - PNFC							
Normal value in fresh water				140.9	r	na/l		
Normal value in marine water				140.9	r	ma/l		
Normal value for fresh water	sediment			552	r	na/ka		
Normal value for marine wate	r sediment			552		ng/kg		
Normal value for water intern	nittent release			140.9	r	ng/l		
Normal value of STP microord				2 251		ng/i		
Normal value for the food cha	in (secondary poisoni	20)		160		ma/ka		
Normal value for the terrestria		19)		28	r	ng/kg		
Health - Derived no-effect	ct level - DNEL / D Effects on consumers	MEL		20	Effects on workers	ng/kg		
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic	Acute local	Acute	Chronic local	Chronic
Oral				Systemic	VND	VND	VND	26 mg/kg bw/d
Inhalation Skin	VND VND	VND VND	VND VND	89 mg/m3 319 mg/ka	VND VND	VND VND	VND VND	500 mg/m3 888 mg/ka
				bw/d				bw/d
ETHANOL								
Threshold Limit Value								
Туре	Country	TWA/8h		STEL/15min		Remarks Observa	i / tions	
		mg/m3	ppm	mg/m3	ppm			
AGW	DEU	380	200	1520	800			
MAK	DEU	380	200	1520	800			
VLA	ESP			1910	1000			
VLEP	FRA	1900	1000	9500	5000			
WEL	GBR	1920	1000					
TLV-ACGIH				1884	1000			
Predicted no-effect concentration - PNEC								

			STEELCO	) S.P.A.		Revision nr. 4 Dated 21/05/2021		
Steelco I	Miele Group Member		Steelcox	ide-DT		Printed on 31/05/2 Page n. 7/23 Replaced revision	2021 :3 (Printed on: 15/09/2	2020)
Normal value in fresh water				960	hố	g/L		
Normal value in marine water				790	m	g/l		
Normal value for fresh water sedi	ment			3,6	m	g/kg		
Normal value for marine water se	diment			2,9	m	g/kg		
Normal value for water, intermitte	nt release			2,75	m	g/l		
Normal value of STP microorgani	sms			580	m	g/l		
Normal value for the food chain (s	secondary poisonin	g)		380	m	g/kg food		
Normal value for the terrestrial co	mpartment			630	hố	g/kg soil dw		
Health - Derived no-effect le	Effects on consumers	1EL			Effects on workers			
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic	Acute local	Acute	Chronic local	Chronic
Oral				systemic		NPI		87 mg/kg
Inhalation	950 mg/m3	NPI	NPI	114 mg/m3	1900 mg/m3	NPI	NPI	950 mg/m3
Skin	NPI	NPI	NPI	206 mg/kg bw/d	NPI	NPI	NPI	343 mg/kg bw/d
SUBTILISINE Threshold Limit Value								
Туре	Country	TWA/8h		STEL/15min		Remarks	; / tions	
		mg/m3	ppm	mg/m3	ppm	0030174	10113	
VLA	ESP			6E-05				
WEL	GBR	4E-05						
Predicted no-effect concentration	- PNEC							
Normal value in fresh water				1,7	μί	g/L		
Normal value in marine water				170	nç	g/L		
Normal value for water, intermitte	nt release			900	nç	g/L		
Normal value of STP microorgani	sms			65	m	g/l		
Normal value for the terrestrial co	mpartment			568	μί	g/kg soil		
Health - Derived no-effect le	Effects on	1EL			Effects on			
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic	Acute local	Acute	Chronic local	Chronic systemic
Oral		3,6 mg/kg bw/d		1,8 mg/kg bw/d		oyotonno		Gyotonno
Inhalation	NPI	NPI	15 ng/m3	NPI	NPI	NPI	60 ng/m3	NPI
Substances that could be relea	sed in case of de	ecomposition:			VIND	NF I		NF I
Туре	Country	TWA/	8h		STEL/15min		Re	marks /
		Mg/m3	3 р	pm	mg/m3	ppm	Obs	
VLEP	ITA	0,37	0	,3	0,74	0,6		
OEL	EU	0,37	0	,3	0,74	0,6		
	-	-,-,		1	0.74	03/0	)	
			0,	•	0,74	0,3 (C	/	

# STEELCO S.P.A.

Revision nr. 4 Dated 21/05/2021

Printed on 31/05/2021

teelco

Steelcoxide-DT

Page n. 8/23 Replaced revision:3 (Printed on: 15/09/2020)

Legend:

(C) = CEILING ; INHAL = Inhalable Fraction ; RESP = Respirable Fraction ; THORA = Thoracic Fraction. VND = hazard identified but no DNEL/PNEC available ; NEA = no exposure expected ; NPI = no hazard identified.

# 8.2. Exposure controls

As the use of adequate technical equipment must always take priority over personal protective equipment, make sure that the workplace is well aired through effective local aspiration.

When choosing personal protective equipment, ask your chemical substance supplier for advice.

Group

Member

Miele

Personal protective equipment must be CE marked, showing that it complies with applicable standards.

### HAND PROTECTION

Protect hands with category III work gloves (see standard EN 374).

The following should be considered when choosing work glove material: compatibility, degradation, failure time and permeability. The work gloves' resistance to chemical agents should be checked before use, as it can be unpredictable. The gloves' wear time depends on the duration and type of use.

## SKIN PROTECTION

Wear category I professional long-sleeved overalls and safety footwear (see Regulation 2016/425 and standard EN ISO 20344). Wash body with soap and water after removing protective clothing.

Consider the appropriateness of providing antistatic clothing in the case of working environments in which there is a risk of explosion.

## EYE PROTECTION

Wear airtight protective goggles (see standard EN 166).

## RESPIRATORY PROTECTION

If the threshold value (e.g. TLV-TWA) is exceeded for the substance or one of the substances present in the product, wear a mask with a type AX filter, whose limit of use will be defined by the manufacturer (see standard EN 14387). In the presence of gases or vapours of various kinds and/or gases or vapours containing particulate (aerosol sprays, fumes, mists, etc.) combined filters are required.

Respiratory protection devices must be used if the technical measures adopted are not suitable for restricting the worker's exposure to the threshold values considered. The protection provided by masks is in any case limited.

If the substance considered is odourless or its olfactory threshold is higher than the corresponding TLV-TWA and in the case of an emergency, wear open-circuit compressed air breathing apparatus (in compliance with standard EN 137) or external air-intake breathing apparatus (in compliance with standard EN 138). For a correct choice of respiratory protection device, see standard EN 529.

### ENVIRONMENTAL EXPOSURE CONTROLS

The emissions generated by manufacturing processes, including those generated by ventilation equipment, should be checked to ensure compliance with environmental standards.

PROPAN-2-OL IBE (Biological Indicators of Exposure - ACGIH 2020): acetone in urine = 40 mg / L (end of shift)

# **SECTION 9.** Physical and chemical properties

# 9.1. Information on basic physical and chemical properties

uid
aw yellow
ld
ot determined
5
ot available
0° 0°

	STEELCO	S.P.A.	Revision nr. 4 Dated 21/05/2021
Miele Group Member			Printed on 31/05/2021
	Steelcoxic	le-DT	Page n. 9/23 Replaced revision:3 (Printed on: 15/09/2020)
Boiling range	Not available		
Flash point	35 °C		
Evaporation rate	Not determined		
Flammability (solid, gas)	not applicable	Reason for missing data:	Based on the
Lower inflammability limit	Not available	physical state	
Upper inflammability limit	Not available		
Lower explosive limit	Not available		
Upper explosive limit	Not available		
Vapour pressure	Not available		
Vapour density	Not determined		
Relative density	1,022 g/cm <sup>3</sup>		
Solubility	soluble in water		
Partition coefficient: n-octanol/water	Not determined		
Auto-ignition temperature	Not available		
Decomposition temperature	Not determined		
Viscosity	Not determined		
Explosive properties	Not Explosive		
Oxidising properties	Not Oxidizing		

#### 9.2. Other information

Information not available

# **SECTION 10. Stability and reactivity**

# 10.1. Reactivity

There are no particular risks of reaction with other substances in normal conditions of use.

#### ETHANEDIOL

In the air absorbs moisture.Decomposes at temperatures above 200°C/392°F. It can absorb atmospheric moisture up to twice its own weight. It decomposes at temperatures above 200 ° C.

# 10.2. Chemical stability

The product is stable in normal conditions of use and storage.

#### 10.3. Possibility of hazardous reactions

The vapours may also form explosive mixtures with the air.

#### ETHANEDIOL

Risk of explosion on contact with: perchloric acid.May react dangerously with: chlorosulphuric acid,sodium hydroxide,sulphuric acid,phosphorus pentasulphide,chromium (III) oxide,chromyl chloride,potassium perchlorate,potassium dichromate,sodium peroxide,aluminium.Forms explosive mixtures with: air.

Risk of explosion on contact with: perchloric acid. It can react dangerously with: chlorosulfuric acid, sodium hydroxide, sulfuric acid, phosphorus pentasulfide, chromium (III) oxide, chromyl chloride, potassium perchlorate, potassium dichromate, sodium peroxide, aluminum. Forms explosive mixtures with air.

		STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021
Steelco	Miele Group Member		Printed on 31/05/2021
		Steelcoxide-DT	Page n. 10/23 Replaced revision:3 (Printed on: 15/09/2020)

# ETHANOL

Risk of explosion on contact with: alkaline metals,alkaline oxides,calcium hypochlorite,sulphur monofluoride,acetic anhydride,acids,concentrated hydrogen peroxide,perchlorates,perchloric acid,perchloronitrile,mercury nitrate,nitric acid,silver,silver nitrate,ammonia,silver oxide,ammonia,strong oxidising agents,nitrogen dioxide.May react dangerously with: bromoacetylene,chlorine acetylene,bromine trifluoride,chromium trioxide,chromyl chloride,fluorine,potassium tert-butoxide,lithium hydride,phosphorus trioxide,black platinum,zirconium (IV) chloride,zirconium (IV) iodide.Forms explosive mixtures with: air.

### 10.4. Conditions to avoid

Avoid overheating. Avoid bunching of electrostatic charges. Avoid all sources of ignition.

## ETHANEDIOL

Avoid exposure to: sources of heat,naked flames. Avoid exposure to heat sources and naked flames.

### PROPAN-2-OL

Heat, flames and sparks. Extreme temperatures and direct sunlight.

### ETHANOL

Avoid exposure to: sources of heat,naked flames. Avoid high temperatures and proximity to ignition sources

#### 10.5. Incompatible materials

PROPAN-2-OL Oxidizing agents, acid anhydrides, aluminum, halogenated compounds, acids.

ETHANOL Strong mineral acids, oxidizing agents. High temperature aluminum.

### 10.6. Hazardous decomposition products

In the event of thermal decomposition or fire, gases and vapours that are potentially dangerous to health may be released.

### ETHANEDIOL

May develop: hydroxyacetaldehyde, glyoxal, acetaldehyde, methane, carbon monoxide, hydrogen. Due to thermal decomposition or in the event of fire, gases and vapors potentially harmful to health can be released.

# **SECTION 11. Toxicological information**

### 11.1. Information on toxicological effects

In the absence of experimental toxicological data on the product itself, the possible health hazards of the product have been assessed on the basis of the properties

of the substances contained, according to the criteria established by the reference legislation for classification. Therefore, consider the concentration of the individual dangerous substances possibly mentioned in sect. 3, to evaluate the resulting toxicological effects from exposure to the product.

### Metabolism, toxicokinetics, mechanism of action and other information

### PROPAN-2-OL

It is readily absorbed following inhalation exposure and rapidly spreads to tissues. However, it is also readily excreted in the urine, essentially in the form of the 2-methoxyacetic acid metabolite. (Arch Toxicol, 68, -588-94 - Johanson G, 1994)

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021	
Miele Group Member		Printed on 31/05/2021	
	Steelcoxide-DT	Page n. 11/23 Replaced revision:3 (Printed on: 15/09/2020)	
Information on likely routes of exposure			
Information not available			
Delayed and immediate effects as well as chroni	c effects from short and long-term exposure		
ETHANEDIOL Ingestion initially stimulates the central n uremia. Over-exposure symptoms are: v	ervous system; later replaced by a phase of depression. omiting, drowsiness, difficulty in breathing, convulsions. T	There may be kidney damage, with anuria and he lethal dose for humans is approx. 1.4 ml/kg.	
Interactive effects			
Information not available			
ACUTE TOXICITY			
ATE (Inhalation) of the mixture: Not classified (no ATE (Oral) of the mixture:>2000 mg/kg ATE (Dermal) of the mixture: Not classified (no s	o significant component) ignificant component)		
ETHANEDIOL LD50 (Oral): 7712 mg/kg Ratto LD50 (Dermal): 9530 mg/kg Coniglio LC50 (Inhalation): 2,5 mg/l/6h Ratto			
ETHANOL LD50 (Oral): 1187 mg/kg Ratto LC50 (Inhalation): 115,9 mg/l/4h			
PROPAN-2-OL LD50 (Oral): 4710 mg/kg Rat LD50 (Dermal): 12800 mg/kg Rat LC50 (Inhalation): 72,6 mg/l/4h Rat			
SUBTILISINE LD50 (Oral): 1800 mg/kg Rat			
ETHANEDIOL Harmful if swallowed (Harmonized classi	fication, Annex VI, CLP Reg.)		
Reference: Evaluation of the Developm (Fundamental and Applied Toxicology 24 Reliability (Klimisch score): 2 Species: rat (Sprague-Dawley; Male/Fen Routes of exposure: inhalation (aerosol) Results: LC50> 2.5 mg / L	ental Toxicity of Ethylene Glycol Aerosol in the CD Ra I: 57-75 (1995)) nale)	at and CD-1 Mouse by Whole-Body Exposure.	
Reference: Assessment of the Develop Toxicology 27: 155-166 (1995)) " Reliability (Klimisch score): 2 Species: mouse (CD-1; Male/Female) Routes of exposure: cutaneous Results: LD50> 3500 mg / kg.	pmental Toxicity of Ethylene Glycol Applied Cutaneou	sly to CD-1 Mice. (Fundamental and Applied	
PROPAN-2-OL Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 2 Species: Rat (Sherman)	01		

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021
Steelco Miele Group Member	Steelcoxide-DT	Printed on 31/05/2021
	Oleeicoxide-Di	Page n. 12/23 Replaced revision:3 (Printed on: 15/09/2020)
Routes of exposure: oral Results: LD50 = 5840 mg / kg		
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 2 Species: Rabbit Routes of exposure: cutaneous Results: LC50 = 16.4 ml / kg	02	
Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: Rat (Fischer 344; Male/Female Routes of exposure: inhalation (vapors) Results: LD50> 10000 ppm / 6h	03 )	
ETHANOL Method: OECD 401 Reliability (Klimish score): 1 Species: rat (Cox CD; Male/Female) Route of exposure: oral Results: LD50: 10470 mg / kg		
Method: OECD 403 Reliability (Klimish score): 2 Species: rat (Sprague-Dawley; Male/Fer Route of Exposure: inhalation (vapors) LC50 results (male): 116.9 mg / I 4h	nale)	
Reference: Schechter, M. et al, Pharmac Reliability (Klimisch score): 2 Species: Mouse (HS; Male/Female) Routes of exposure: intraperitoneal Results: LD50 = 9450 mg / kg body weig	col Biochem Behav 52 (1): 245-248, 1995 ht	
SUBTILISINE Method: OECD 401 Reliability (Klimisch score): 1 Species: rat (Wistar; Male/Female) Routes of exposure: oral Results LD50: 1800 mg / kg body weigh The substance is classified as harmful b	t / day y inhalation.	
SKIN CORROSION / IRRITATION		
Does not meet the classification criteria for this h	azard class	
ETHANEDIOL Method: BASF-internal standards Reliability (Klimisch score): 2 Species: rabbit (Vienna White) Routes of exposure: cutaneous Results: non-irritating		
PROPAN-2-OL Reliability (Klimisch score): 2 Species: Rabbit Routes of exposure: cutaneous Results: Not irritating Reference: Nixon G et al, Toxicology and	d Applied Pharmacology 31, 481-490 (1975)	

ETHANOL

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021
Steelco Miele Group Member	Steelcoxide-DT	Printed on 31/05/2021 Page n. 13/23 Replaced revision:3 (Printed on: 15/09/2020)
Method: OECD 404 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White) Routes of exposure: cutaneous Results: non-irritating.		
SUBTILISINE The substance causes skin irritation Har	monized classification, Annex VI, CLP Reg.)	
SERIOUS EYE DAMAGE / IRRITATION		
Does not meet the classification criteria for this h	nazard class	
ETHANEDIOL Method: BASF-internal standards Reliability (Klimisch score): 2 Species: rabbit (Vienna White) Routes of exposure: ocular Results: non-irritating		
PROPAN-2-OL Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: Rabbit (New Zealand White) Routes of exposure: ocular Results: irritating	05	
ETHANOL Method: OECD 405 Reliability (Klimisch score): 2 Species: Rabbit Routes of exposure: ocular Results: irritating.		
SUBTILISINE The substance causes serious eye dama	age (Harmonized Classification, Annex VI, CLP Reg.)	
RESPIRATORY OR SKIN SENSITISATION		
May produce an allergic reaction.Contains:SUB	<b>FILISINE</b>	
ETHANEDIOL Reference: Evaluation of Skin Irritation Pigs. (Dental Materials Journal 15 (2): 2: Reliability (Klimisch score): 2 Species: Guinea Pig (Dunkin-Hartley; Fe Routes of exposure: cutaneous Results: not sensitizing.	and Sensitization of Two Diol Solutions used as Experi 26-232 (1996)) emale)	mental Dentin Primers in Humans and Guinea
PROPAN-2-OL Method: OECD 406 Reliability (Klimisch score): 1 Species: Guinea pig (Dunkin-Hurtley; Ma Routes of exposure: cutaneous Results: not sensitizing	ale/Female)	
ETHANOL Method: equivalent or similar OECD 406 Reliability (Klimisch score): 2		

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021							
Steelco Miele Group Member	Steelcoxide-DT	Printed on 31/05/2021 Page n. 14/23 Replaced revision:3 (Printed on: 15/09/2020)							
Species: Guinea pig (Pirbright White; Female) Routes of exposure: cutaneous Results: not sensitizing.									
SUBTILISINE The substance causes respiratory sensit	SUBTILISINE The substance causes respiratory sensitization (Harmonized Classification, Annex VI, CLP Reg.)								
GERM CELL MUTAGENICITY									
Does not meet the classification criteria for this h	nazard class								
ETHANEDIOL Method: OECD 471 - In vitro test Reliability (Klimisch score): 1 Species: S. typhimurium TA 1535, TA 15 Results: negative	537, TA 98, TA 100 and E. coli WP2								
Method: Publication 1986 - In vivo testin Reliability (Klimisch score): 2 Species: rat (Fischer 344; Male/Female) Routes of exposure: oral Results: negative.	Method: Publication 1986 - In vivo testing Reliability (Klimisch score): 2 Species: rat (Fischer 344; Male/Female) Routes of exposure: oral Results: negative.								
PROPAN-2-OL Based on available data, the substance	has no mutagenic effects and is not classified under the re	elevant hazard class CLP.							
ETHANOL Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: S. typhimurium Results: negative with and without metal	71 - In vitro test bolic activation								
Method: equivalent or similar to OECD 4 Species: mouse (NMRI; Male/Female) Routes of exposure: intraperitoneal Results: negative.	74 - In vivo test								
SUBTILISINE Method: OECD 473 - Test in vitro Reliability (Klimisch score): 1 Species: human (lymphocytes) Results: negative with and without metal	bolic activation								
CARCINOGENICITY									
Does not meet the classification criteria for this h	nazard class								
ETHANEDIOL Available studies have shown no carcinogenic potential. In a carcinogenicity study lasting two years, carried out by the US National Toxicology Program (NTP), in which ethylene glycol was administered in the feed, "no evidence of carcinogenic activity" in male and female B6C3F1 mice was observed (NTP, 1993).									
PROPAN-2-OL Based on available data, the substance	has no carcinogenic effects and is not classified under the	relevant hazard class CLP.							
ETHANOL Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: rat (Fischer 344 / DuCrj; Male/f	-53 Female)								

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021		
Miele Group Miele Member		Printed on 31/05/2021 Page n. 15/23 Replaced revision:3 (Printed on: 15/09/2020)		
	Steelcoxide-DT			
Routes of exposure: inhalation (vapors) Results: negative.				
SUBTILISINE No data available.				
REPRODUCTIVE TOXICITY				
Does not meet the classification criteria for this h	azard class			
ETHANEDIOL Reference: Chronic Toxicity and Oncogenicity Studies of Ethylene Glycol in Rats and Mice. (Fundamental and Applied Toxicology 7: 547-565 (1986)) Reliability (Klimisch score): 2 Species: mouse (CD-1; Male/Female) Routes of exposure: oral Results: negative.				
PROPAN-2-OL Method: equivalent or similar to OECD 416 Reliability (Klimisch score): 1 Species: Rat (Sprague-Dawley; Male/Female) Routes of exposure: Oral Results: negative. NOAEL = 1000 mg / kg bw / day.				
SUBTILISINE No data available.				
Adverse effects on sexual function and fertility				
PROPAN-2-OL Method: equivalent or similar to OECD 416 Reliability (Klimsch score): 1 Species: rat (Sprague-Dawley Male/Female) Routes of exposure: oral Results: negative.				
ETHANOL Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 1 Species: mouse (CD-1; Male/Female) Routes of exposure: oral Results: No effect on fertility at doses eq	16 uivalent to 20.7 g / kg / day			
Adverse effects on development of the offspring				
PROPAN-2-OL Method: equivalent or similar to OECD 4 Reliability (Klimsch score): 1 Species: rat (Sprague-Dawley) Routes of exposure: oral Results: negative.	14			
ETHANOL Method: equivalent or similar to OECD 4 Reliability (Klimisch score): 2 Species: rat (Sprague-Dawley) Routes of exposure: inhalation Results: negative. NOAEL (maternal) = -	14 16000 ppm. NOAEL (fetus)> = 20,000 ppm			

		STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021	
Steelco	Miele Group Member		Printed on 31/05/2021	
multuer		Steelcoxide-DT	Page n. 16/23 Replaced revision:3 (Printed on: 15/09/2020)	
STOT - SINGLE EXPOSURE				
Does not meet the classification	n criteria for this h	azard class		
ETHANEDIOL Based on available dat CLP hazard class.	a, the substance	has no specific target organ toxicity effects for single ex	posure and is not classified under the relevant	
PROPAN-2-OL Method: OECD 426 Reliability (Klimisch score): 1 Species: Rat (Sprague-Dawley; Female) Routes of exposure: oral. Results: May cause drowsiness or dizziness. Based on available data, the substance exhibits specific target organ toxicity effects from single exposure and is classified under the relevant CLP hazard class				
ETHANOL Based on available dat CLP hazard class.	a, the substance	has no specific target organ toxicity effects for single ex	posure and is not classified under the relevant	
SUBTILISINE Based on available data, the substance exhibits specific target organ toxicity effects from single exposure and is classified under the relevant CLP hazard class. (Harmonized classification, Annex VI, CLP Reg.)				
Target organ				
SUBTILISINE Lungs.				
Route of exposure				
SUBTILISINE Inalation.				
STOT - REPEATED EXPOSURE				
Does not meet the classification	n criteria for this h	azard class		
ETHANEDIOL Method: equivalent or similar to OECD 452 Reliability (Klimisch score): 2 Species: Rat (Wistar; Male) Routes of exposure: oral Results: observed toxicity towards the kidneys and bladder. NOAEL = 150 mg / kg bw / day				
PROPAN-2-OL Based on available data, the substance has no specific target organ toxicity effects on repeated exposure and is not classified under the relevant CLP hazard class.				
ETHANOL Method: equivalent or similar OECD 408 Reliability (Klimisch score): 2 Species: Rat (Sprague-Dawley; Male/Female) Routes of exposure: oral Results: negative. NOAEL: 1730 mg / kg body weight / day				

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021
Steel co Miele Group		Printed on 31/05/2021
Member	Stoolcovido-DT	
	Steelcoxide-D1	Page n. 17/23 Replaced revision:3 (Printed on: 15/09/2020)
SUBTILISINE Based on available data, the substance I CLP hazard class.	has no specific target organ toxicity effects on repeate	d exposure and is not classified under the relevant
Target organ		
ETHANEDIOL Kidneys		
Route of exposure		
ETHANEDIOL Oral		
ASPIRATION HAZARD		
Does not meet the classification criteria for this h	azard class	
ETHANEDIOL No data on aspiration hazard are availab	le	
PROPAN-2-OL No data on aspiration hazard are availab	le	
ETHANOL No data are available on the hazard in ca	ase of aspiration.	
SUBTILISINE No data on aspiration hazard are availab	le	
SECTION 12. Ecological information	ation	
Use this product according to good working p contaminate soil or vegetation.	practices. Avoid littering. Inform the competent auth	norities, should the product reach waterways or
12.1. Toxicity		
ETHANEDIOL		
LC50 - for Fish	> 72,86 g/l/96h Pimephales promelas	(EPA 600/4-90/027)
EC50 - for Crustacea	> 100 mg/l/48h Daphnia magna (OEC	D 202)
EC50 - for Algae / Aquatic Plants	> 10940 mg/l/72h Pseudokirchneriella	subcapitata (EPA/600/4-89/001)
Chronic NOEC for Fish	> 15,38 g/l/7 d	
Chronic NOEC for Crustacea	> 8,59 g/l/7 d	
Chronic NOEC for Algae / Aquatic Plants	100 mg/l/72h	
FTHANOI		

LC50 - for Fish EC50 - for Crustacea EC50 - for Algae / Aquatic Plants Chronic NOEC for Fish 14,2 g/l/96h Pimephales promelas (US EPA E03-05) 5012 mg/l/48h Ceriodaphnia dubia (ASTM E729-80) 275 mg/l/72h Chlorella vulgaris (OECD 201) 250 mg/L/5 d

	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021		
		Printed on 31/05/2021		
Member	Stoolcovido-DT			
	Steelcoxide-Di	Page n. 18/23 Replaced revision:3 (Printed on: 15/09/2020)		
PROPAN-2-OL				
LC50 - for Fish	9640 mg/l/96h Pimephales promelas (eo	quivalent or similar to OECD 203)		
EC50 - for Crustacea	> 10000 mg/l/48h Daphnia magna (equi	valent or similar to OECD 202)		
SUBTILISINE				
LC50 - for Fish	8,2 mg/l/96h Oncorhynchus mykiss; (OE	CD 203(		
EC50 - for Crustacea	0,17 mg/l/48h Daphnia Magna; OECD 2	02		
EC50 - for Algae / Aquatic Plants	0,29 mg/l/72h Pseudokirchnerella subca	pitata (OECD 201)		
Chronic NOEC for Fish	0,042 mg/l/34d Pimphales promelas; (O	ECD 210)		
Chronic NOEC for Crustacea	0,019 mg/l/14d Daphnia magna; (OECD	211)		
Chronic NOEC for Algae / Aquatic Plants	0,041 mg/l/72h Pseudokirchnerella subc	apitata; (OECD 201)		
12.2. Persistence and degradability				
ETHANEDIOL Rapidly degradable, 90-100% in 10 days (OEC Solubility in water 1000 - 10000 mg/l	CD 301 A)			
PROPAN-2-OL Rapidly degradable, 53% in 5 days (equivalent	t or similar to EU C.5)			
ETHANOL Readily biodegradable, 60% in 10 days (BOD - Standard methods for the examination of water and waste water 1971. 13th ed, American Public Health Assoc, NY)				
SUBTILISINE Rapidly degradable, 79% in 28 days (OECD 30	01 B).			
12.3. Bioaccumulative potential				
ETHANEDIOL				
Partition coefficient: n-octanol/water	-1,36 Log Kow ( ACS Professional Ret	ference Book, 1995)		
ETHANOL				
Partition coefficient: n-octanol/water	-0,35 Log Kow 24°C (OECD 107)			
PROPAN-2-OL				
Partition coefficient: n-octanol/water	0,05			
12.4. Mobility in soil				
Information not available				
12.5. Results of PBT and vPvB assessment				
On the basis of available data, the product does not contain any PBT or vPvB in percentage $\geq$ than 0,1%.				
12.6. Other adverse effects				
Information not available				

# STEELCO S.P.A.

Revision nr. 4 Dated 21/05/2021

Printed on 31/05/2021

Steelcoxide-DT

Page n. 19/23 Replaced revision:3 (Printed on: 15/09/2020)

# **SECTION 13.** Disposal considerations

## 13.1. Waste treatment methods

Reuse if possible. Product residues are to be considered special hazardous waste. The dangerousness of the waste that partially contains this product must be evaluated according to the laws in force. (Ref. Annex D - Part IV of Legislative Decree no. 152/2006 and subsequent amendments and adjustments).

Disposal must be entrusted to an authorized waste management company, in compliance with national and possibly local regulations.

The legal responsibility for disposal lies with the producer / holder of the waste.

Group

Member

Different CER (European Waste Code) codes could be applied to this mixture according to the specific circumstances that generated the waste, any alterations and contaminations.

The product as it is, out of specification in the original packaging, or poured into a suitable container for disposal as waste, or the product in specification but no longer usable (for example following an accidental spill), is to be classified with a code CER compatible with the description of use indicated in section 1.2.

The appropriate final destination of the waste will be assessed by the manufacturer according to the chemical-physical characteristics of the waste itself compatible with the authorized plant to which it will be conferred for recovery, treatment or final disposal in the manner prescribed by current regulations. Disposal via the wastewater drain is not permitted.

For hazardous substances registered according to EC Regulation 1907/2006 (REACH) for which a chemical safety report has been prepared, refer to the specific information contained in the exposure scenarios attached to this SDS.

# CONTAMINATED PACKAGING

Contaminated packaging must be sent, properly labeled, for recovery or disposal in compliance with national regulations on waste management and must be classified with the following EWC code:

15 01 10 \*: packaging containing residues of dangerous substances or contaminated by these substances

# **SECTION 14. Transport information**

## 14.1. UN number

ADR / RID, IMDG, 1993 IATA:

### 14.2. UN proper shipping name

ADR / RID:	FLAMMABLE LIQUID, N.O.S. (Propan-2-ol; Ethanol) MIXTURE
IMDG:	FLAMMABLE LIQUID, N.O.S. (Propan-2-ol; Ethanol) MIXTURE
IATA:	FLAMMABLE LIQUID, N.O.S. (Propan-2-ol; Ethanol) MIXTURE

### 14.3. Transport hazard class(es)

ADR / RID:	Class: 3	Label: 3
IMDG:	Class: 3	Label: 3
IATA:	Class: 3	Label: 3

### 14.4. Packing group

ADR / RID, IMDG, III IATA:



		STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021		
Steelco Miele Group Member		Steelcoxide-DT	Printed on 31/05/2021 Page n. 20/23 Replaced revision:3 (Printed	Printed on 31/05/2021 Page n. 20/23 Replaced revision:3 (Printed on: 15/09/2020)	
14.5. Environmental haz	ards				
ADR / RID: MDG: MDG: MATA: MDG	NO NO NO				
14.6. Special precaution	s for user				
ADR / RID:	ИН	l - Kemler: 30	Limited Quantities: 5 L	Tunnel restriction code: (D/E)	
IMDG: IATA:	Spo EM Cai Pat	ecial provision: - IS: F-E, S-E rgo: ss.:	Limited Quantities: 5 L Maximum quantity: 220 L Maximum quantity: 60 L	Packaging instructions: 366 Packaging instructions: 355	
	Spo	ecial provision:	A3		
14.7. Transport in bulk a	according to Annex II	of Marpol and the IBC Code			
Information not relevant					
SECTION 15. Re	gulatory inform	ation			
15.1. Safety, health and	l environmental regul	ations/legislation specific for the substance or m	ixture		
Seveso Category - Directi	ive 2012/18/EC: P5c				
Restrictions relating to the	e product or contained	substances pursuant to Annex XVII to EC Regulation	1907/2006		
Product					
Point	3				
	Liquid substances in Annex I of Regu a) hazard classes 2, 2.15 types A to b) hazard classes narcotic effects, 3. c) hazard class 4. d) hazard class 5.	or mixtures that meet the criteria relating to one of t lation (EC) No. 1272/2008: 2.1 to 2.4, 2.6 and 2.7, 2.8 types A and B, 2.9, 2.10, F; 3.1 to 3.6, 3.7 adverse effects on sexual function 9 and 3.10; 1; 1.	he following hazard classes 2.12, 2.13 categories 1 and and fertility or developmen	or categories as set out 2, 2.14 categories 1 and t, 3.8 effects other than	
Point	40				
	Substances classi category 1 or 2, c category 1 pyroph 1272/2008Contair	fied as flammable gases of category 1 or 2, flamma, substances and mixtures which, in contact with wa oric liquids or category 1 pyrophoric solids, even if no red substance	ble liquids of category 1, 2 c ter, emit flammable gases ot listed in Part 3 of Annex V	or 3, flammable solids of of categories 1, 2 or 3, /l to Regulation (EC) No.	
Point	75	SUBTILISINE Reg. No. 01-2119480434-38-xxxx			
	Substances includ	led in one or more of the following:			

r				
	STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021		
Miele Group Member		Printed on 31/05/2021		
	Steelcoxide-DT	Page n. 21/23 Replaced revision:3 (Printed on: 15/09/2020)		
<ul> <li>a) substances classified in one of the following classes in Annex VI, part 3, of Regulation (EC) no. 1272/2008:</li> <li>category 1A, 1B or 2 carcinogenicity, category 1A, 1B or 2 germ cell mutagenicity, but excluding substances classified due to effects following exposure by inhalation only;</li> <li>Reproductive toxicity of category 1A, 1B or 2, but excluding substances classified due to effects following exposure by inhalation only;</li> <li>skin sensitization of category 1, 1A or 1B;</li> <li>skin corrosion of category 1, 1A, 1B or 1C or skin irritation of category 2;</li> <li>serious category 1 eye damage or category 2 eye irritation;</li> <li>b) substances listed in Annex II of Regulation (EC) no. 1223/2009 of the European Parliament and of the Council (*);</li> <li>c) substances listed in Annex IV of Regulation (EC) no. 1223/2009 for which a condition is indicated in at least one of the columns g, h or i of the table of this annex;</li> <li>d) substances listed in Appendix 13 of this annex. The ancillary requirements referred to in points 7 and 8 of column 2 of this entry apply to all mixtures intended for tattooing practices, regardless of whether they contain one of the substances referred to in points a) to d) of this column and entry.</li> </ul>				
Regulation (EC) No. 2019/1148 - on the marketin	ng and use of explosives precursors			
Not applicable				
Substances in Candidate List (Art. 59 REACH)				
On the basis of available data, the product does	not contain any SVHC in percentage $\geq$ than 0,1%.			
Substances subject to authorisation (Annex XIV REACH)				
None				
Substances subject to exportation reporting pursuant to (EC) Reg. 649/2012:				
None				
Substances subject to the Rotterdam Convention:				
None				
Substances subject to the Stockholm Convention:				
None				
Healthcare controls				
Workers exposed to this chemical agent must not undergo health checks, provided that available risk-assessment data prove that the risks related to the workers' health and safety are modest and that the 98/24/EC directive is respected.				
Medical device, class IIb according to legislation 93/42.				
15.2. Chemical safety assessment				
A chemical safety assessment has not been perf	ormed for the preparation/for the substances indicated in	section 3.		
SECTION 16. Other information				
Text of hazard (H) indications mentioned in section 2-3 of the sheet:				
Flam. Liq. 2 Flammable liquid, c	ategory 2			

			· · · · · · · · · · · · · · · · · · ·
		STEELCO S.P.A.	Revision nr. 4 Dated 21/05/2021
Steelco	Miele Group		Printed on 31/05/2021
	Member	Steelcovide-DT	
		Olecicoxide-D1	Page n. 22/23 Replaced revision:3 (Printed on: 15/09/2020)
Flam. Liq. 3	Flammable liquid, c	ategory 3	
Skin Corr. 1A	Skin corrosion, cate	gory 1A	
Eye Irrit. 2	Eye irritation, categ	ory 2	
STOT SE 3	Specific target orga	n toxicity - single exposure, category 3	
H225	Highly flammable lie	quid and vapour.	
H226	Flammable liquid a	nd vapour.	
H314	Causes severe skir	burns and eye damage.	
H319	Causes serious eye	irritation.	
H336	May cause drowsin	ess or dizziness.	
LEGEND:			
ADR: European Agreement concerning the carriage of Dangerous goods by Road - CAS NUMBER: Chemical Abstract Service Number - CESO: Effective concentration (required to induce a 50% effect) - CE NUMBER: Identifier in ESIS (European archive of existing substances) - CEP. EC Regulation 1272/2008 - DNEL: Derived No Effect Level - EmS: Emergency Schedule - GHS: Globally Harmonized System of classification and labeling of chemicals - HTA DGR: International Air Transport Association Dangerous Goods Regulation - ICSO: Immobilization Concentration 50% - IMDC: International Maritime Organization - ICSO: International Maritime Organization - INDEX NUMBER: Identifier in Annex VI of CLP - LCSO: Lethal concentration 50% - UESO: Lethal concentration 50% - UESO: Lethal concentration 50% - UESO: Lethal concentration 50% - UESO: Lethal concentration 50% - DEL: Predicted consource Level - PBT: Persistent bioaccumulative and toxic as REACH Regulation - PEC: Predicted environmental Concentration - PEC: Predicted environmental Concentration - PEL: Predicted environmental Concentration - REACH: EC Regulation 107/2006 - RID: Regulation concerning the international transport of dangerous goods by train - TLV: Threshold Limit Value - TLV CELING: Concentration that should not be exceeded during any time of occupational exposure. - TWA STEL: Short-term exposure limit - VOC: Volatile organic Compounds - wPW: Ware Presistent and very Bioaccumulative as for REACH Regulation - WCK: Water hazard classes (German).			
<ol> <li>Regulation (EC) 1907/2006 (REACH) of the European Parliament</li> <li>Regulation (EU) 1272/2008 (CLP) of the European Parliament</li> <li>Regulation (EU) 790/2009 (I Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2015/830 of the European Parliament</li> <li>Regulation (EU) 286/2011 (II Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 618/2012 (III Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 487/2013 (IV Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 944/2013 (V Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 944/2013 (V Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2015/121 (VII Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2016/1121 (VII Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2016/1121 (VII Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2016/1121 (VII Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2016/1179 (IX Atp. CLP) of the European Parliament</li> <li>Regulation (EU) 2017/776 (X Atp. CLP)</li> <li>Regulation (EU) 2018/1480 (XIII Atp. CLP)</li> <li>Regulation (EU) 2018/1480 (XIII Atp. CLP)</li> </ol>			

Revision nr. 4 Dated 21/05/2021

Printed on 31/05/2021

Veelco Miele

Steelcoxide-DT

Page n. 23/23 Replaced revision:3 (Printed on: 15/09/2020)

- 16. Regulation (EU) 2019/521 (XII Atp. CLP)
- 17. Regulation (EU) 2019/1148
- 18. Regulation (EU) 2020/217 (XIV Atp. CLP) The Merck Index. 10th Edition
- Handling Chemical Safety
- INRS Fiche Toxicologique (toxicological sheet)
- Patty Industrial Hygiene and Toxicology
- N.I. Sax Dangerous properties of Industrial Materials-7, 1989 Edition

Group

Member

- IFA GESTIS website
- ECHA website

- Database of SDS models for chemicals - Ministry of Health and ISS (Istituto Superiore di Sanità) – Italy

# Training for workers:

The training of workers must include contents, updates and duration according to the risk profiles assigned to the working sectors of membership, according to the procedures provided for by Legislative Decree 81/2008.

### Procedure used to derive the classification according to Regulation (EC) 1272/2008 (CLP) in relation to mixtures:

Classification of the mixture according to Regulation (EC) r	n. 1272/2008	Classification procedure
Flammable liquid, category 2	H225	Based on experimental data.

# Note for the recipient of the Safety Data Sheet (SDS):

It is the recipient of this SDS who must ensure that the information contained is read and understood by all persons who handle, store, use, or otherwise come into contact in any way with the substance or mixture to which this sheet refers. In particular, the recipient must provide adequate training to personnel assigned to the use of dangerous substances or mixtures. The recipient must ensure the suitability and completeness of the information in relation to the specific use of the substance or mixture.

However, the substance or mixture to which this SDS refers must not be used for uses other than those specified in section 1. No responsibility is assumed for improper uses. Since the use of the product does not fall under the direct control of the Supplier, it is the user's obligation to observe, under his own responsibility, the laws and regulations in force regarding national and Community hygiene and safety.

The information contained in this SDS is provided in good faith and is based on the current state of scientific and technical knowledge, at the revision date indicated, available from the Supplier indicated in section 1 of this sheet. The SDS should not be interpreted as a guarantee of any specific property of the substance or mixture. The information refers only to the substance or mixture specifically designated in section 1 and may not be valid for the substance or mixture used in combination with other materials or in other processes not specifically indicated in the text.

This version of the SDS supersedes all previous versions.

Changes from the previous revision.

Changes have been made to the following sections: 01/02/03/04/05/06/07/08/09/10/11/12/13/14/15/16.