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LX-2011 Technical Specification

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The constantly high speed rotating anode and the new bearing method are the essential characteristics of the X-ray tube housing assembly LX-2011.

The anode is carried in a sleeve bearing filled with fluid lubricant. It enables noise free rotation and long life expectation.

The tube features the application of high rated cine operation and grid-controlled pulsed fluoroscopic operation.

The anode has high heat storage capacity of 2130 kJ and the tube is cooled by water-cooling heat exchanger.

This high cooling rate of 3.1 kW provides high patient throughput

GENERAL SPECIFICATION

IEC Classification:

IEC60601-1/1988 Amd.1/1991, Amd.2/1995 Class I

Electrical:

Circuit:

High Voltage Generator..... Constant potential high-voltage generator

Grounding..... Center grounded

Tube voltage:

Radiographic 40 ~ 125 kV

Fluoroscopic 40 ~ 125 kV

Nominal focal spot size (IEC60336/1993):

Large focal spot size1.0

Small focal spot size.....0.6

Input power (0.1s) Refer to the rating charts

Large focal spot..... 100 kW

Small focal spot 48 kW

Anode speed Minimum 9000 min⁻¹

Grid cut-off voltage.....-2900 VDC

(Large / Small focal spot)

Dielectric resistance (Between ground terminal and low voltage terminal) 2 MΩ Min.

Input voltage to heat exchanger.....AC 115 V, 2 A

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Mechanical:

Outline dimension.....	Refer to the drawing
Overall length	562 mm
Maximum diameter	213 mm
Anode target:	
Angle	11°
Material.....	Rhenium-Tungsten
Inherent filtration.....	Aluminum equivalent 1.5 mm (Including external aluminum filter of 0.4 mm thickness)
Test technique factor.....	75 kV (IEC60522/1999)
X-ray leakage:	
Complies with IEC60601-1-3/1994	
Test technique factor.....	125 kV, 18 mA
X-ray coverage:	
SID 900 mm	φ 305 mm
Weight:	
X-ray tube assembly.....	Approximately 39 kg
Heat exchanger	Approximately 17 kg
High voltage plug	CLAYMOUNT Corp. MINI-75 type
Low voltage terminal	Refer to the outline drawing
Cooling method	Water cooling heat exchanger type
Operating position	Non-designated

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ABSOLUTE MAXIMUM RATINGS

(Do not exceed the followings in any case)

Nominal X-ray tube voltage :

Radiographic	125 kV
Fluoroscopic	125 kV
Between anode (or cathode) and ground.....	65 kV

Minimum tube voltage 40 kV

Maximum X-ray tube current (Refer to the rating charts) :

Large focal spot	850 mA
Small focal spot	460 mA

Maximum filament current :

Large focal spot	5.8 A
Small focal spot	5.1 A

Maximum filament voltage :

Large focal spot (5.8 A)	13.2 ~ 17.2 V
Small focal spot (5.1 A)	9.4 ~ 12.6 V

Long time maximum Input

60 minutes average	2200 W (3100 HU/s)
30 minutes or less	2500 W (3525 HU/s)
20 minutes or less	2800 W (3950 HU/s)

Thermal characteristic :

Anode heat content	2130 kJ (3000 kHU)
Anode heat dissipation	5500 W (7700 HU/s)
X-ray tube assembly heat content.....	2050 kJ (2890 kHU)

Continuous heat dissipation:

Including water cooling heat exchanger ⁽¹⁾	3100 W (260 kHU/min)
---	----------------------

Note⁽¹⁾ This is specified with the heat exchanger WHE-01B

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Humidity 30 ~ 85 %
(No condensation)

Pressure 70 ~ 106 kPa

Transporting and storage condition (When no cooling liquid is filled) :

Temperature -20 ~ 70 °C

When cooling liquid is filled..... 2 ~ 60 °C

Humidity 20 ~ 90 %
(No condensation)

Pressure..... 50 ~ 106 kPa

Tube positioning.....Cathode side needs to be lower

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
SERIAL RADIOGRAPHIC RATING CHARTS

Operating condition: Constant potential high-voltage generator

Power supply frequency50/60 Hz

Anode rotation speed9000 min.⁻¹


Note: Tube current is depending on filament current for each tube voltage, therefore, emission characteristic needs to be referred.

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


FILAMENT & EMISSION CHARACTERISTIC

Constant potential high-voltage generator

Nominal focal spot value: 1.0


Nominal focal spot value: 0.6

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


THERMAL CHARACTERISTICS

Housing thermal characteristics (room temperature 25°C)

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



TECHNICAL DATA

ANODE ROTATING STARTER UNIT ST-7012

ST-7012 provides the necessary control signals, interlocks and power outputs to drive the anode rotor of the X-ray tube(LM-Bearing) to near 9000 min^{-1} .

GENERAL SPECIFICATION

1. Electrical

1.1 Power supply

Input power	3 phase AC220V $\pm$ 10%
Capacity	5.3 kVA MIN.
Frequency	50 / 60 Hz
Power	1.2kW / 5A MAX.

1.2 Output

Stator

At boosting	AC 220 V
At running	AC 150 V
Pump mounted on X-ray tube	AC 115 V
Getter	AC 10~14 V
Water-cooling unit	AC220V

2. Mechanical

2.1 Outline dimension, Weight

Outline	437(W) $\times$ 359(D) $\times$ 210(H)
(including base plate)	482(W) $\times$ 359(D) $\times$ 210(H)
Weight	16.6 kg

3. Environmental limits

3.1 Operating condition:


Temperature	18 ~ 40 °C
Humidity	30 ~ 85 %
	(No condensation)
Pressure	70 ~ 106kPa

3.2 Transportation and storage condition

Temperature	-20 ~ 65 °C
Humidity	20 ~ 90 %
	(No condensation)
Pressure	70 ~ 106kPa

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
Unit: mm

SYMBOL

- CNN1:ST-7012 control signal connector
- CNN2:ST-7012 power supply connector
- CNN3:WHE-01B limit switch signal connector
- CNN4:LX-3081 limit switch connector
- CNN5:Stator power supply connector
- CNN6:Getter power supply connector

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TECHNICAL DATA

WATER-COOLING HEAT EXCHANGER WHE-01B

WHE-01B is used with water-cooling type X-ray tube assembly and its high cooling performance provides high patient throughput.

GENERAL SPECIFICATION

IEC Classification

(IEC60601-1) Class I
Complied standard
IEC60601-1(General/1988) Amendment2/1995

Electrical:

Applicable X-ray tube:
Water-cooling type X-ray tube housing assembly
Heat dissipation rate Refer to technical documentation for X-ray tube
Power supply 1 ϕ 220V (50/60 Hz)
Power consuming 280W / 340W

Mechanical:

Outline dimension Refer to outline drawing
450 mm(W) $\times$ 370 mm(D) $\times$ 400 mm(H)
Weight Approximately 30 kg
Operating position Horizontal

Installation requirement:

Installation position
..... 3.5 m lower position than X-ray tube
Pipe length 30 m for one way
(40m Maximum)
Power supply cable 220V (3 conductors)
Ground cable length 20 m
Safety signal cable 20 m with 3 conductors

Environmental limits:

Operating condition
Temperature 5 ~ 40 °C
Humidity 30 ~ 90 %
(No condensation)
Pressure 70 ~ 106 kPa

Transportation and storage condition:


Temperature
At shipping (no cooling liquid is filled)
..... -20 ~ 70 °C
When cooling liquid is filled 2 ~ 60 °C
Humidity 20 ~ 90 %
(No condensation)
Pressure 50 ~ 106 kPa

Accessory:

Cooling liquid (pure water) 10 liter

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
OUTLINE DRAWING

Unit: mm

OVERSEAS SUBSIDIARIES AND AFFILIATES

EU REPRESENTATIVE

•TOSHIBA ELECTRONICS EUROPE GMBH

HANSAALLEE 181 40549 DUSSELDORF, GERMANY

PHONE +49 (0) 211 5296 107

FAX +49 (0) 211 5296 402

For Sales & Technical Services, please contact the following representative:

•TOSHIBA ELECTRONICS EUROPE GMBH

HANSAALLEE 181 40549 DUSSELDORF, GERMANY

PHONE +49 (0) 211 5296 107

FAX +49 (0) 211 5296 402

•TOSHIBA AMERICA ELECTRONIC COMPONENTS, INC.

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Toshiba Electron Tubes & Devices Co., Ltd. meets the Environmental Management System Standard, ISO 14001

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Toshiba Electron Tubes & Devices Co., Ltd. meets internationally recognized Standards for Quality Management System ISO 9001, ISO 13485

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TOSHIBA ELECTRON TUBES & DEVICES CO., LTD.