# User Manual

DSH-270WA

LCD Monitor

# Table of contents

|       |                                     |    |
|-------|-------------------------------------|----|
| 1     | Introduction.....                   | 3  |
| 1.1   | About the product.....              | 3  |
| 1.2   | About the user manual.....          | 3  |
| 1.3   | Packing list.....                   | 4  |
| 1.3.1 | Standard.....                       | 4  |
| 1.3.2 | Optional.....                       | 4  |
| 2     | Safety information.....             | 5  |
| 3     | Parts, controls and connectors..... | 8  |
| 3.1   | Front view.....                     | 8  |
| 3.2   | Back view.....                      | 8  |
| 3.3   | Connector pin assignments.....      | 9  |
| 4     | Connection and Quick start.....     | 11 |
| 4.1   | Connection.....                     | 11 |
| 4.2   | Quick start.....                    | 11 |
| 5     | Operation.....                      | 12 |
| 5.1   | Signal indicator.....               | 12 |
| 5.2   | Shortcuts.....                      | 12 |
| 5.3   | Key function.....                   | 12 |
| 5.4   | OSD MENU.....                       | 13 |
| 6     | Technical specification.....        | 14 |
| 7     | Supported signal format.....        | 17 |
| 8     | Cleaning and Disinfection.....      | 19 |
| 9     | Troubleshooting.....                | 20 |
| 10    | Appendix.....                       | 21 |
| 10.1  | Markings and symbols.....           | 21 |
| 10.2  | Conformity declaration.....         | 22 |
| 10.3  | RoHS information.....               | 23 |

## **1 Introduction**

### **1.1 About the product**

DSH-270WA is a FHD surgical medical monitor. This product adopts advanced display technology, specially designed for use in medical environment. True and accurate image restoration, high brightness, high contrast, wide angle of view display performance, ensure its display effect in minimally invasive and operating room environment, provide reliable support for doctors for detailed observation. A variety of video input and output interfaces that fully meet the requirements of various endoscopes, operating rooms and ICU.

### **1.2 About the user manual**

This document describes the functionality and the reasonable usage of DSH-270WA. For ensuring the clarity, please note that it does not contain all detailed information on this product.

This user manual contains important information that you have to read in order to ensure your personal safety and prevent damages to your property. Please read these safety information and product information so that you can operate the machine safely and effectively.

This user manual is only suitable for DSH-270WA.

Please save this manual for future use.

This user manual is updated or changed without any notices.

### 1.3 Packing list

Please check the following articles when you buy the display. Please note that the pictures are only for references. The following pictures only mean the similar items.

#### 1.3.1 Standard

LCD monitor 1pcs  
 Adapter 1pcs

LCD monitor

Adapter

#### 1.3.2 Optional

|                  |        |                     |        |
|------------------|--------|---------------------|--------|
| Certification    | ■ 1pcs | Bolt                | ■ 4pcs |
| DP signal cable  | □ 1pcs | Power cord          | ■ 1pcs |
| BNC signal cable | ■ 1pcs | CD                  | ■ 1pcs |
| DVI signal cable | ■ 1pcs | Equipotential cable | ■ 1pcs |
| Warranty card    | ■ 1pcs | HDMI signal cable   | □ 1pcs |
| Nylon cable ties | ■ 3pcs | Packing list        | ■ 1pcs |

Certification

DP signal cable

BNC signal cable

DVI signal cable

Warranty card


Bolt

Power cord


CD

Equipotential cable

HDMI signal cable

Nylon cable ties

Packing list

## 2 Safety information

To ensure personal safety and proper maintenance, please read this section and the caution statements on the unit.

### WARNING

**If the unit begins to emit smoke, smells like something is burning, or makes strange noises, disconnect all power connections immediately and contact your dealer for advice.**

Attempting to use a malfunctioning unit may result in fire, electric shock, or equipment damage. Information provided for the OPERATOR to safely terminate operation of ME EQUIPMENT.

**Do not open the cabinet or modify the unit.**

Opening the cabinet or modifying the unit may result in fire, electric shock, or burn.

**Refer all servicing to qualified service personnel.**

Do not attempt to service this product yourself as opening or removing covers may result in fire, electric shock, or equipment damage.


Only those who have been trained can use it.

**small objects or liquids away from the unit.**

Small objects accidentally falling through the ventilation slots into the cabinet or spillage into the cabinet may result in fire, electric shock, or equipment damage.

If an object or liquid falls/spills into the cabinet, unplug the unit immediately.

Have the unit checked by a qualified service engineer before using it again.

**This monitor can be a long boot, but the life will be reduced.**

**Suggest that the device is operated by personnel experiencing training.**

**The device is not serviced or maintained while in use with the patient.**

**Place the device on a solid and stable surface.**


If the unit is placed on an inadequate surface may occur fall and result in personal injury or equipment damage. If the unit falls, disconnect the power immediately and ask dealer for help. Please do not continue to use damaged equipment. Using a damaged unit may result in fire or electric shock.

**Set the unit in an appropriate location.**

Not doing so may result in fire, electric shock, or equipment damage.

- Do not place outdoors.
- Do not place in the transportation system (ship, aircraft, trains, automobiles etc.
- Do not place in a dusty or humid environment.
- Do not place in a location where the steam comes directly on the screen.
- Do not place near heat generating devices or a humidifier.
- Do not place in an inflammable gas environment.

**To avoid danger of suffocation, keep the plastic packing bags away from babies and children..**


**Use the enclosed power cord and connect to the standard power outlet of your country.**

Be sure to remain within the rated voltage of the power cord. Not doing so may result in fire or electric shock.

|   |   |
|---|---|
| <p><b>To disconnect the power cord, grasp the plug firmly and pull.</b><br/>Tugging on the cord may damage and result in fire or electric shock.</p>  |  |
| <p><b>The equipment must be connected to a grounded main outlet.</b><br/>To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.</p>  |  |
| <p><b>Use the correct voltage.</b></p> <ul style="list-style-type: none"> <li>➤ The unit is designed for use with a specific voltage only. Connection to another voltage than specified in this User's Manual may cause fire, electric shock, or equipment damage.</li> <li>➤ Do not overload your power circuit, as this may result in fire or electric shock.</li> </ul>        |   |
| <p><b>Handle the power cord with care.</b></p> <ul style="list-style-type: none"> <li>➤ Do not press the power cord under the equipment or other heavy objects.</li> <li>➤ Do not pull or wrap the power cord. If the power cord is damaged, please stop using it. Using a damaged power cord may cause a fire or electric shock.</li> </ul>                                      |  |
| <p><b>Never touch the plug and power cord if it begins to thunder.</b></p> <p style="text-align: center;">  </p> <p>Touching them may result in electric shock.</p>  |   |
| <p><b>Do not touch a damaged LCD panel directly with bare hands.</b></p> <p>The liquid crystal that may leak from the panel is poisonous if it enters the eyes or mouth. If any part of the skin or body comes in direct contact with the panel, please wash thoroughly. If some physical symptoms result, please consult your doctor.</p>  |   |
| <p><b>Follow local regulation or laws for safe disposal.</b></p> <p>The MANUFACTURER will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist SERVICE PERSONNEL to repair those parts of ME EQUIPMENT that are designated by the MANUFACTURER as repairable by SERVICE PERSONNEL.</p> |   |

 **CAUTION**

|   |   |
|---|---|
| <p><b>Handle with care when carrying the unit.</b><br/>Disconnect the power cord and signal cables and remove the optional unit. Moving the unit with the cord or the option attached is dangerous. It may result in injury.</p>  |   |
| <p><b>When handling the unit, grip the bottom of the unit firmly with both hands ensuring the panel faces outward before lifting.</b><br/>Dropping the unit may result in injury or equipment damage.</p>   |  |
| <p><b>Do not block the ventilation slots on the cabinet.</b></p> <ul style="list-style-type: none"> <li>➤ Do not place any objects on the ventilation slots.</li> <li>➤ Do not install the unit in a closed space.</li> <li>➤ Do not use the unit laid down or upside down.</li> </ul> <p>Blocking the ventilation slots prevents proper airflow and may result in fire, electric shock, or equipment damage.</p> |  |
| <p><b>Do not touch the plug with wet hands.</b><br/>Doing so may result in electric shock.</p>  |   |
| <p><b>Periodically clean the area around the plug.</b><br/>Dust, water, or oil on the plug may result in fire.</p>  |   |

**Unplug the unit before cleaning it.**

Cleaning the unit while it is plugged into a power outlet may result in electric shock.

**If you plan to leave the unit unused for an extended period, disconnect the power cord from the wall socket after turning off the power switch for the safety and the power conservation.**

**The ME EQUIPMENT that are not serviced or maintained while in use with the patient .**

**The device cannot used under OXYGEN RICH ENVIRONMENTS**

**The device cannot used with flammable anaesthetics.**

**The device cannot used in conjunction with flammable agents.**

**Do not modify this equipment without authorization of the manufacturer.**

**Do not position ME EQUIPMENT to make it difficult to operate the disconnection device.**

**LCD PANEL**

When the monitor is cold and brought into a room or the room temperature goes up quickly, dew condensation may occur inside and outside the monitor. In that case, do not turn the monitor on and wait until dew condensation disappears, otherwise it may cause some damages to it.

In order to suppress the luminosity change by long-term use and to maintain the stable luminosity, please use the monitor with the lower brightness.

The screen may have defective pixels. These pixels may appear as slightly light or dark area on the screen. This is due to the characteristics of the panel itself, and not the product.

The backlight of the LCD panel has a fixed life span. When the screen becomes dark or begins to flicker, please contact your dealer.

Do not press on the panel or edge of the frame strongly, as this may result in damage to the screen. There will be prints left on the screen if the pressed image is dark or black. If pressure is repeatedly applied to the screen, it may deteriorate or damage your LCD panel. Leave the screen white to decrease the prints.

Do not scratch or press on the panel with any sharp objects, such as a pencil or pen as this may result in damage to the panel. Do not attempt to brush with tissues as this may scratch the LCD panel.

When the screen image is changed after displaying the same image for extended periods of time, an afterimage may appear. Use the screen saver or timer to avoid displaying the same image for extended periods of time.


The monitor can be at power on for a long time, but its lifespan may be reduced. Suggest it should take a rest for 10 minutes per hour.

Turn on the monitor and computer orderly.


The power indicator of the monitor will become green. If images cannot appear, please refer to "Troubleshooting" for help. Turn off the computer and monitor after every use.

### 3 Parts, controls and connectors

#### 3.1 Front view


#### 3.2 Back view


### 3.3 Connector pin assignments

#### ● DVI Connector


| Pin No. | Signal                 | Pin No. | Signal                  |
|---------|------------------------|---------|-------------------------|
| 1       | T.M.D.S. Data2-        | 16      | Hot Plug Detect         |
| 2       | T.M.D.S. Data2+        | 17      | T.M.D.S. Data0-         |
| 3       | DVI Signal Detect      | 18      | T.M.D.S. Data0+         |
| 4       | T.M.D.S. Data 4-       | 19      | T.M.D.S. Data0/5 Shield |
| 5       | T.M.D.S. Data 4+       | 20      | T.M.D.S. Data 5-        |
| 6       | DDC Clock              | 21      | T.M.D.S. Data 5+        |
| 7       | DDC Data               | 22      | T.M.D.S. Clock Shield   |
| 8       | No connect             | 23      | T.M.D.S. Clock+         |
| 9       | T.M.D.S. Data1-        | 24      | T.M.D.S. Clock-         |
| 10      | T.M.D.S. Data1+        | C1      | No connect              |
| 11      | T.M.D.S Data1/3 Shield | C2      | No connect              |
| 12      | T.M.D.S. Data 3-       | C3      | No connect              |
| 13      | T.M.D.S. Data 3+       | C4      | No connect              |
| 14      | +5V Power              | C5      | No connect              |
| 15      | Ground                 |         |                         |

#### ● HDMI Connector


| Pin No. | Signal | Pin No. | Signal   |
|---------|--------|---------|----------|
| 1       | DATA2+ | 11      | GND      |
| 2       | GND    | 12      | CLK-     |
| 3       | DATA2- | 13      | NC       |
| 4       | DATA1+ | 14      | NC       |
| 5       | GND    | 15      | SCL      |
| 6       | DATA1- | 16      | SDA      |
| 7       | DATA0+ | 17      | GND      |
| 8       | GND    | 18      | +5V      |
| 9       | DATA0- | 19      | HOT_PLUG |
| 10      | CLK+   |         |          |

● DP Connector


| Pin No. | Signal        | Pin No. | Signal          |
|---------|---------------|---------|-----------------|
| 1       | ML_Lane 3(N)  | 11      | GND             |
| 2       | GND           | 12      | ML_Lane 0(P)    |
| 3       | ML_Lane 3 (P) | 13      | GND             |
| 4       | ML_Lane 2(N)  | 14      | GND             |
| 5       | GND           | 15      | AUX CH(P)       |
| 6       | ML_Lane 2(P)  | 16      | GND             |
| 7       | ML_Lane 1(N)  | 17      | AUX CH(N)       |
| 8       | GND           | 18      | Hot Plug Detect |
| 9       | ML_Lane 1(P)  | 19      | Return          |
| 10      | ML_Lane 0(N)  | 20      | DP_PWR          |

● VGA Connector

| Pin | Signal      | Pin | Signal         | Pin | Signal      |
|-----|-------------|-----|----------------|-----|-------------|
| 1   | Red video   | 6   | Red ground     | 11  | NC          |
| 2   | Green video | 7   | Green ground   | 12  | Data (SDA)  |
| 3   | Blue video  | 8   | Blue ground    | 13  | H. Sync     |
| 4   | NC          | 9   | NO Pin         | 14  | V. Sync     |
| 5   | GND         | 10  | Ground Shorted | 15  | Clock (SCL) |

● SDI Connector

| Pin No. | Signal |
|---------|--------|
| 1       | SDI    |
| 2       | GND    |

● Power supply Connector


| Pin No. | Function | Pin No. | Function |
|---------|----------|---------|----------|
| 1       | +V (24V) | 4       | +V (24V) |
| 2       | Earth    | 5       | NC       |
| 3       | -V (GND) | 6       | -V (GND) |

## 4 Connection and Quick start

**Note:** Please make sure that the **power switch** of endoscopic equipment and display are **OFF**

### 4.1 Connection

Please reference the connection figure showing below.

### 4.2 Quick start

Turn on the display power, and the power indicator will turn on, start the endoscope host and format the output signal to:

Resolution : **1920 ×1080** Refresh ratio: **60Hz**.

If no image display on the screen, please refer to the "failure analysis" to get help.

**Please use the enclosed power cord to connect the standard power outlet.**

Please ensure that the monitor's power plug can be fully inserted into the socket. Make sure the power cord used within the rated voltage. Or it may result in fire or electric shock.

Unplug the power cord and disconnect from the power grid.

**Equipment must be connected to a outlet that grounded good.**

Or it may result in fire or electric shock.

**Please turn on the monitor and then computer power.**

The monitor's power indicator will light (green). If the image does not appear, refer to the "failure analysis" for help. After each use, turn off the computer and the monitor.

**Note:**

This monitor can run for a long time, but life will be reduced. Recommended 10-minute break every hour.

## 5 Operation


### 5.1 Signal indicator

| Color  | Display mode                |
|--------|-----------------------------|
| Green  | Operating mode              |
| Orange | Standby mode                |
| None   | Power Off or No power input |

### 5.2 Shortcuts

| Key     | Function                |
|---------|-------------------------|
| < Left  | Adjust contrast values  |
| > Right | Adjust backlight values |

### 5.3 Key function

| Keys    | Function  |
|---------|---|
| Source  | Active Input source menu / switch to next input |
| Menu    | Active the Main menu / confirm the selection    |
| EXIT    | Exit menu / return                              |
| Left    | Press to move, select items or adjust value     |
| Right   |   |
| Down    |   |
| Up      |   |
| Control | Key Locked/Unlock: Press and hold  key          |

## 5.4 OSD MENU

| Main Menu     | Option                     | Description   |       |    |
|---------------|----------------------------|---|-------|----|
| Input         | Main Input                 | DVI /VGA /SDI /HMDI /DP /RGB  |       |    |
|               | PIP Mode                   | Off /Side by Side /Small PIP /Large PIP                                 |       |    |
| Image         | Main/Pip Select            | Main /PIP   |       |    |
|               | Brightness                 | Adjust brightness:(0-100)   |       |    |
|               | Contrast                   | Adjust contrast:(0-100)   |       |    |
|               | Saturation                 | Adjust saturation:(0-100)   |       |    |
|               | Sharpness                  | Adjust sharpness:(0-24)   |       |    |
|               | Aspect Ratio               | Full Screen /Native Mode /Fill  |       |    |
|               | Zoom                       | 80,85,90,95,100,105,110,115,120   |       |    |
| Setup         | Color Temp                 | 10000K/9300K/8500K/7500K/6500K/5500K/User                               |       |    |
|               |                            | R   | 0-100 | 50 |
|               |                            | G   | 0-100 | 50 |
|               |                            | B   | 0-100 | 50 |
|               | Gamma                      | DICOM/2.6/2.4/2.2/2.0/1.8/Video   |       |    |
|               | Mirror/Rotation            | Off/Mirror/Rotation   |       |    |
|               | Auto Scan Input            | On/Off  |       |    |
|               | Video Level                | Auto /Limit /Full   |       |    |
|               | RGB/YPbPr SEL              | Set input channel to YPbPr or RGB                                       |       |    |
| Factory Reset | Restore to factory default |   |       |    |
| OSD           | Language                   | English / 中文 /Español /Português /Français /Deutsche /Italiano /русский |       |    |
|               | OSD Position               | Menu position selection:<br>Bottom Center /Bottom Right                 |       |    |
| VGA           | Auto Adjust                | Auto adjust the image position, clock and phase                         |       |    |
|               | H-Position                 | Adjust the horizontal position of the image                             |       |    |
|               | V-Position                 | Adjust the vertical position of the image                               |       |    |
|               | Clock                      | Adjust the frequency of the display's pixel clock                       |       |    |
|               | Phase                      | Adjust the phase of the display's pixel clock                           |       |    |
| Status        |                            | Display    FW Version<br>Main Input<br>Main Format                      |       |    |

## 6 Technical specification

| <b>Monitor characteristics:</b>    |  |
|------------------------------------|--|
| Panel technology                   | 27", TFT, color, LCD screen, anti-glare, IPS,hard coating<br>Optical adhesive bonding, high transparency and low reflection, anti condensation and anti fog.       |
| Active area (H x V)                | 597.89 (H) x 336.31 (V) mm   |
| Pixel pitch                        | 0.3114(H) x 0.3114(V) mm   |
| Resolution                         | 1920 (RGB) x 1080  |
| Refresh rate                       | 60 Hz  |
| Luminance                          | 1000 cd/m <sup>2</sup> (typical)   |
| Contrast ratio                     | 1000:1 (typical)   |
| Viewing angle (CR > 10)            | Horizontal: 178° (typical)   |
|                                    | Vertical: 178° (typical)   |
| Backlight                          | WLED   |
| Lifetime of backlight              | 30000 hours(minimal)   |
| Response time                      | T <sub>GTG</sub> :14ms(typical)  |
| <b>Power supply:</b>               |  |
| DC IN                              | 24 V, 5.0 A (supplied from AC adapter)   |
| Power consumption                  | < 90 W   |
| Standby                            | < 0.5 W  |
| <b>Control and connection:</b>     |  |
| Front                              | Power indicator<br>Touch keypad  |
| Back                               | DVI -IN *1<br>DVI -OUT *1<br>VGA -IN *1<br>SDI -IN *1/OUT *1<br>DP *1<br>HDMI *1<br>RS232 *1<br>DC 5V/0.8A-OUT*1<br>Ground stud*1<br>DC socket*1<br>Power switch*1 |
| <b>Mechanical characteristics:</b> |  |
| Housing components                 | Plastic  |
| cover glass                        | Glass protection: surface treatment AR+AF,high Translucency, low reflection.   |
| Ventilation openings               | Natural heat dissipation   |
| Protection level                   | IP65 (Front) 、 IP20 (back)   |

| <b>Safety regulations:</b>                       |  |                  |
|--|--|------------------|
| Safety standards                                 | CCC: GB4943.1-2011、GB17625.1-2012、GB/T9254.1-2021<br>CE: EN60601-1-2、EN IEC61000-3-2、EN61000-3-3 |                  |
| Protection class                                 | Protection class I   |                  |
| <b>Dimension:</b>                                |  |                  |
| Without packing(D x W x T)                       | 645.26 mm x 391.86 mm x 62.5 mm  |                  |
| With packing (D x W x H)                         | 783 mm x 245 mm x 541 mm   |                  |
| <b>Weight</b>                                    |  |                  |
| Monitor  | Net weight   | 9.6 kg ± 0.5 kg  |
|  | Gross weight   | 14.2 kg ± 0.5 kg |
| Monitor+ Stand                                   | Net weight   | 13.4 kg ± 0.5 kg |
|  | Gross weight   | 18 kg ± 0.5 kg   |
| <b>Climatic conditions:</b>                      |  |                  |
| Operational temperature                          | 5°C ~ 35°C(Recommended 20°C ~ 30°C)  |                  |
| Operational humidity                             | 15% ~ 85% Relative humidity, no condensation   |                  |
| Operational pressure                             | 70KPa ~106KPa  |                  |
| Transport and storage temperature                | -20°C ~ 60°C   |                  |
| Transport and storage humidity                   | 10% ~ 90% Relative humidity, no condensation   |                  |
| Storage and transport atmospheric pressure range | 70KPa ~106KPa  |                  |

**Electromagnetic Emission:**

The LCD Monitor is intended to be used in professional healthcare facility environments such as clinics and hospitals

The requirements special precautions regarding EMC and during installation. You need to carefully read the EMC Information and the “PRECAUTIONS” section of this document, and observe the following instructions when installing and operating the product.

Anyone who connects additional equipment to the signal input parts or signal output parts when configuring a medical system is responsible for ensuring that the system complies with the requirements of IEC/EN60601-1-2.

The LCD Monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of the LCD Monitor should assure that it is used in such an environment.

## Declaration - Electromagnetic Emission

| Description of Test Item  | Compliance   | Limits Performance Criteria | Results |
|---|--|-----------------------------|---------|
| Conducted Disturbance at Mains Terminals CISPR 11   | Class B  | Normal                      | Pass    |
| Radiated Disturbance CISPR 11   | Class B  | Normal                      | Pass    |
| Harmonic Current Emissions IEC 61000-3-2  | Class D  | Normal                      | N/A     |
| Voltage Fluctuation and Flicker IEC 61000-3-3   | Section 5  | Normal                      | Pass    |
| IMMUNITY (EN 60601-1-2:2015, IEC 60601-1-2:2014)  |  |                             |         |
| Description of Test Item  | Compliance level                                   | Performance Criteria        | Results |
| Electrostatic Discharge (ESD) IEC 61000-4-2   | ± 8 kV contact<br>± 15 kV air                      | Normal                      | Pass    |
| Radio-Frequency, Continuous Radiated Disturbance IEC 61000-4-3  | Level: 3 10 V/m<br>80 MHz to 2.5 GHz               | Normal                      | Pass    |
| EFT/B Immunity IEC 61000-4-4  | Level 3: ± 2 kV for power supply lines             | Normal                      | Pass    |
| Surge Immunity IEC 61000-4-5  | ± 1 kV line(s) to lines<br>± 2 kV line(s) to earth | Normal                      | Pass    |
| Conducted RF Immunity IEC 61000-4-6   | Level 2 0.15-80 MHz @ 3 VRMS                       | Normal                      | Pass    |
| Power Frequency Magnetic Field IEC 61000-4-8  | Level 4, 30 A/m                                    | Normal                      | Pass    |
| Voltage Dips, >95% Reduction,<br>Voltage Dips, 30% Reduction,<br>Voltage Interruptions IEC 61000-4-11 |  | Normal                      | Pass    |

## 7 Supported signal format

- DVI

| Video Formats | Horizontal Frequency (kHz) | Vertical Frequency (Hz) | Pixel Clock (MHz) |
|---------------|----------------------------|-------------------------|-------------------|
| 720×400       | 31.50                      | 70.0                    | 28.322            |
| 640×480       | 31.50                      | 60.0                    | 25.175            |
| 800×600       | 37.90                      | 60.0                    | 40.000            |
| 1024×768      | 48.40                      | 60.0                    | 65.000            |
| 1280×1024     | 64.00                      | 60.0                    | 108.00            |
| 1920×1080     | 56.25                      | 50.0                    | 148.50            |
| 1920×1080     | 67.50                      | 60.0                    | 148.50            |

- VGA

| Video Formats | Horizontal Frequency (kHz) | Vertical Frequency (Hz) | Pixel Clock (MHz) |
|---------------|----------------------------|-------------------------|-------------------|
| 720×400       | 31.50                      | 70.0                    | 28.322            |
| 640×480       | 31.5                       | 60.0                    | 25.175            |
| 800×600       | 37.9                       | 60.0                    | 40.00             |
| 1024×768      | 48.40                      | 60.0                    | 65.00             |
| 1280×1024     | 64.00                      | 60.0                    | 108.00            |
| 1920×1080     | 56.25                      | 50.0                    | 148.50            |
| 1920×1080     | 67.50                      | 60.0                    | 148.50            |

- DP/HDMI

| Video Formats | Horizontal Frequency (kHz) | Vertical Frequency (Hz) | Pixel Clock (MHz) |
|---------------|----------------------------|-------------------------|-------------------|
| 720×400       | 31.50                      | 70.0                    | 28.322            |
| 640×480       | 31.50                      | 60.0                    | 25.175            |
| 800×600       | 37.9                       | 60.0                    | 40.000            |
| 1024 × 768    | 48.40                      | 60.0                    | 65.00             |
| 1280 × 1024   | 64.00                      | 60.0                    | 108.00            |
| 1920 x 1080   | 56.25                      | 50.0                    | 148.50            |
| 1920 × 1080   | 67.50                      | 60.0                    | 148.50            |

- RGB/YPbPr

| Video Formats | Horizontal Frequency (kHz) | Vertical Frequency (Hz) | Pixel Clock (MHz) |
|---------------|----------------------------|-------------------------|-------------------|
| 480i/60Hz     | 15.734                     | 59.940                  | 13.5              |
| 576i/50Hz     | 15.625                     | 50.0                    | 13.5              |
| 480p/60Hz     | 31.469                     | 59.940                  | 27.0              |
| 576p/50Hz     | 31.25                      | 50.0                    | 27.0              |
| 720p/50Hz     | 37.50                      | 50.0                    | 74.250            |
| 720p/60Hz     | 45.00                      | 60.0                    | 74.250            |
| 1080i/50Hz    | 28.125                     | 50.0                    | 74.250            |
| 1080i/60Hz    | 33.750                     | 60.0                    | 74.250            |
| 1080p/30Hz    | 33.75                      | 30.0                    | 74.250            |
| 1080p/50Hz    | 56.250                     | 50.0                    | 148.50            |
| 1080p/60Hz    | 67.500                     | 60.0                    | 148.50            |

- SDI

|                                       |
|---------------------------------------|
| 1) 480/ 60I, 59.94I                   |
| 2) 576/ 50I                           |
| 3) 720/ 50p, 60p, 59.94p              |
| 4) 1080/ 25PsF, 25p                   |
| 5) 1080/ 30p, 29.97p, 30PsF, 29.97PsF |
| 6) 1080/ 50I                          |
| 7) 1080/60I, 59.94I                   |
| 8) 1080/50p                           |
| 9) 1080/60p, 59.94p                   |

## 8 Cleaning and Disinfection

Regular cleaning is recommended to keep the appearance of the monitor clean and prolong the service life.

Gently wipe the dirt on the monitor shell or screen with a piece of microfiber cloth dipped with a small amount of water (tap water, distilled water or mineral water) or chemical reagent (70% - 75% alcohol).

Note:

- Please turn off the power first When cleaning the LCD Monitor.
- It is forbidden to use dry fiber cloth to wipe the monitor glass to prevent static electricity generated on the glass surface, which may cause damage to the monitor.
- Don't often use chemical reagent to wipe the LCD monitor; Frequent use of alcohol, disinfectant and other chemical reagents to wipe the LCD monitor may cause the shell or screen to lose luster and fade, thereby reducing the image quality.
- Don't spray detergent directly on the surface of the LCD monitor to prevent the liquid from penetrating into the equipment, damaging the internal circuit and causing equipment failure.
- The detergent should be tested and inspected before use, the other detergent that have not been tested and inspected may damage the monitor shell or screen.

## 9 Troubleshooting



If a problem persists even after applying the suggested remedies, contact a dealer.

| Problems   | Points to check with Possible Solutions   |
|--|---|
| <b>1. No picture</b><br>Indicator status: Off  | Check that the power cord is correctly connected. If the problem persists, turn off the monitor power for a few minutes, then turn it back on and try again.  |
| Indicator status: Green  | Check the <Backlight> shortcut key setting.   |
| Indicator status: Orange   | Try pressing the keys on the computer keyboard or clicking the mouse.   |
| <b>2. Display position is incorrect.</b><br>        | Adjust the image position by using the <H-Position> and <V-Position>. If the problem persists, use the graphics board's utility software to change the display position if available.   |
| <b>3. Vertical bars of distortion appear.</b><br>   | Decrease the horizontal bars by using the <Clock>.  |
| <b>4. Horizontal bars of distortion appear.</b><br> | Decrease the horizontal bars by using the <Phase>.  |
| <b>5. The screen is too bright or too dark.</b>  | Adjust <Backlight> shortcut key (The backlight of the LCD monitor has a fixed life span. When the screen becomes dark or begins to flicker, please contact your dealer.)  |
| <b>6. After image phenomena appear</b>   | <ul style="list-style-type: none"> <li>- Make sure that whether you have used the screen saver or timer when displaying the same image for extended periods of time or not.</li> <li>- Afterimages are particular to LCD monitors. Avoid displaying the same image for extended periods of time.</li> </ul> |
| <b>7. The screen has defective pixels</b>  | This is due to the characteristics of the panel itself, and not the LCD product.  |
| <b>8. Shiatsu traces left on the screen</b>  | Keep the screen background white can solve this problem.  |

## 10 Appendix

### 10.1 Markings and symbols

The markings and symbols on the monitor, power supply and packing have the following meanings (nonrestrictive list).

|   |   |  |   |
|---|---|--|---|
|    | Caution   |    | Earth (ground)  |
|    | Warning: dangerous voltage  |    | Batch number  |
|    | Date of manufacture   |    | Authorized representative in the European Community               |
|    | Manufacturer  |    | Refer to instruction manual/ booklet                              |
|    | Product code  |    | Product ID  |
|    | Medical device  |    | RoHS mark   |
|   | "WEEE (Waste Electrical and Electronic Equipment)". The waste products should be handled legally.               |   | Fragile, handle with care   |
|  | Class II equipment  |  | This way up   |
|  | Keep dry  |  | Temperature limitation  |
|  | Maximum number of identical packages that may be stacked, where "n" stands for the number of permitted packages |  | Humidity limitation   |
|  | Atmospheric pressure limitation   |  | Indicates complies with Part 15 of the FCC rules.                 |
|  | CE mark   |  | Indicates the device is approved according to the CCC regulations |
|  | TUV approval mark   |  | Marking according to ACPEIP (China-RoHS) (name plate)             |
|  | Energy Conservation Certification   |  | Magnetic Resonance (MR)-Unsafe                                    |
|  | UKCA mark   |  | UK Responsible Person   |

## 10.2 Conformity declaration

This equipment is in conformity with the provisions of the CCC standard.

This equipment satisfies the following two conditions:

(1) The device will not bring dangerous interference;

(2) The device will not be influenced by outer interference and cause unpredictable results.

1. To avoid interference with radio or TV receiving, please use the power adapter and cable matched; Use other cable and power adapter may interfere the other electronic equipment;

2. This equipment has passed the standards of CCC; the device will have a wireless electromagnetic wave, incorrect installation and using will cause interference with wireless communication systems.

3. We don't guarantee that the equipment will never produce interference when it is installed. You can judge whether the equipment can disturb radio or the television reception by the way to switch this equipment power source. Please do the following steps if there is interference to the other devices when running this equipment:

- Adjust the angle and position of the radio or the television reception antenna;
- Increase the distance between this equipment and the radio or the television reception antenna;
- This equipment doesn't share the same the power input circuit with the radio or the television reception antenna;
- Please consult the professional person if any problem;

Other ancillary equipment is connected to our monitor must also be passed CCC standards. When configures a medical system, any device is connected to the input or output of the equipment must comply with medical equipment standards IEC 60601-1. If you are unsure, please ask professional technician for help.

### European Authorised Representative

SUNGO Europe B.V.  
Olympisch Stadion 24,  
1076DE Amsterdam, Netherlands

### 10.3 RoHS information

The following product pollution control information is provided according to SJ/T11364-2006 Marking for Control of Pollution caused by Electronic Information Products.

#### Explanation of Pollution Control Label

This symbol indicates the product contains hazardous materials in excess of the limits established by the Chinese standard SJ/T11363-2006 Requirements for Concentration Limits for Certain Hazardous Substances in Electronic Information Products. The number in the symbol is the Environment-friendly Use Period (EFUP), which indicates the period during which the toxic or hazardous substances or elements contained in electronic information products will not leak or mutate under normal operating conditions so that the use of such electronic information products will not result in any severe environmental pollution, any bodily injury or damage to any assets. The unit of the period is "Year". In order to maintain the declared EFUP, the product shall be operated normally according to the instructions and environmental conditions as defined in the product manual, and periodic maintenance schedules specified in Product Maintenance Procedures shall be followed strictly.

Consumables or certain parts may have their own label with an EFUP value less than the product.

Periodic replacement of those consumables or parts to maintain the declared EFUP shall be done in accordance with the Product Maintenance Procedures.

This product must not be disposed of as unsorted municipal waste, and must be collected separately and handled properly after decommissioning.

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\*: Circuit board assembly includes PCB, components and parts that it consists of, such as resistors, capacitors, ICs and connectors etc.

○: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ/T11363-2006.

X: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in SJ/T11363-2006

Data listed in the table represents best information available at the time of publication.

Applications of hazardous substances in this medical device are required to achieve its intended clinical uses, and/or to provide better protection to human beings and /or to environment, due to lack of reasonably (economically or technically) available substitutes.

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**Declaration - Electromagnetic Emission**

| Description of Test Item  | Compliance   | Limits Performance Criteria | Results |
|---|--|-----------------------------|---------|
| Conducted Disturbance at Mains Terminals CISPR 11   | Class B  | Normal                      | Pass    |
| Radiated Disturbance CISPR 11   | Class B  | Normal                      | Pass    |
| Harmonic Current Emissions IEC 61000-3-2  | Class D  | Normal                      | N/A     |
| Voltage Fluctuation and Flicker IEC 61000-3-3   | Section 5  | Normal                      | Pass    |
| <b>IMMUNITY (EN 60601-1-2:2015, IEC 60601-1-2:2014)</b>   |  |                             |         |
| Description of Test Item  | Compliance level                                   | Performance Criteria        | Results |
| Electrostatic Discharge (ESD) IEC 61000-4-2   | ± 8 kV contact<br>± 15 kV air                      | Normal                      | Pass    |
| Radio-Frequency, Continuous Radiated Disturbance IEC 61000-4-3  | Level: 3 10 V/m<br>80 MHz to 2.5 GHz               | Normal                      | Pass    |
| EFT/B Immunity IEC 61000-4-4  | Level 3: ± 2 kV for power supply lines             | Normal                      | Pass    |
| Surge Immunity IEC 61000-4-5  | ± 1 kV line(s) to lines<br>± 2 kV line(s) to earth | Normal                      | Pass    |
| Conducted RF Immunity IEC 61000-4-6   | Level 2 0.15-80 MHz @ 3 VRMS                       | Normal                      | Pass    |
| Power Frequency Magnetic Field IEC 61000-4-8  | Level 4, 30 A/m                                    | Normal                      | Pass    |
| Voltage Dips, >95% Reduction,<br>Voltage Dips, 30% Reduction,<br>Voltage Interruptions IEC 61000-4-11 |  | Normal                      | Pass    |

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