


Nevia™

DESIGNED FOR 1ST PASS SUCCESS
WITH ALL CLOT TYPES



The logo for Vesalio & Eximia features a stylized blue 'V' on the left, composed of a thick stroke and a smaller circle above it. To the right of the 'V' is the word 'VESALIO' in a bold, black, sans-serif font, followed by an ampersand '&' and the word 'EXIMIA' in the same font style.

are happy to welcome you to this

NeVa Induction Webinar

AGENDA:

Why NeVa?

Clinical Experience

User Tips & Tricks

➤ Question & Answer Session



FIRST PASS SUCCESS GETS BETTER OUTCOMES

First Pass Effect - FPE (*full recanalization in the 1st Pass*) is the most powerful predictor of clinical outcome with best safety results

Zaidat et al., Journal of NeuroInterventional Surgery, 2015 → NASA* registry: 354 patients from 24 US centers

	First Pass Effect Group	Total Patient Population	
<i>n</i>	89 patients (25% of total patient population)	354 patients	FPE: 25.1% of patients
MCA occlusions	64.0%	52.5%	FPE more commonly observed in MCA occlusions
ICA occlusions	10.1%	27.7%	
BGCs use in the subgroup	64%	34.7%	FPE more commonly observed when balloon guide catheters were used
Median time to revascularization	34 min	60 min	<i>P</i> =.0003 Median time to revascularization was significantly faster in FPE group
mRS ≤ 2	61.3%	35.3%	<i>P</i> =.013 Patient outcomes were significantly better in FPE group




WHY DOES FIRST PASS SUCCESS GET BETTER OUTCOMES?

Faster recanalization:
TIME is BRAIN

Lower complication
rates

Impact of complete
recanalization


Each hour delay to treatment is associated with a 5.5% absolute decline in the likelihood of achieving good outcome



Proportion of good outcomes by mTICI grade (P<0.0001 for overall comparison)

- Good outcome (mRS 0 - 2)
- Poor outcome (mRS 3 - 6)

Yoo et al, Stroke. 2013;44:2509-2512



WHY DEVELOP ANOTHER STENT-RETRIEVER

1

TREAT ALL OCCLUSIONS

FROM SOFT CLOTS
THAT EASILY DISINTEGRATE
TO HARD, FIBRIN-RICH CLOTS
THAT ARE IMPENETRABLE

2

IMPROVE PROCEDURAL PERFORMANCE

1ST PASS SUCCESS
TIME TO RECANALIZATION
HIGHER TICI 2C/3 RATES

3


PROVIDE EASE OF USE

REAL TIME FEEDBACK
DURING RETRIEVAL
SYNERGISTIC WITH ALL
ACCESS PHILOSOPHIES

TO ACHIEVE BETTER PATIENT OUTCOMES



CAUSES OF ENDOVASCULAR TREATMENT FAILURE



An organized (*hard, fibrin-rich*) clot is more resilient and less sticky than fresh (*soft, red blood cell-rich*) clots, **causing less engagement with a stent retriever and leading to clot missing during retrieval**, especially in the case of a tortuous arterial tree.

Furthermore, because an organized clot may cause more tension in the stent-retriever deployed segment of the parent artery, this **can also likely induce an arterial spasm**.

CAUSES OF ENDOVASCULAR TREATMENT FAILURE

During experimental evaluation of stent retrievers' mechanical properties and effectiveness, differences were observed in their interaction with white versus red clots:


FINDINGS:

- **Large white clots: could not be removed at all**
- **Smaller white clots:** variable results, clot was engaged - then disengaged and rolled between the device and the vessel wall
- **Red clots:** showed a tendency to fragment causing distal embolization
- Device's ability to maintain a constant radial pressure and wall apposition during retrieval played an important role

*“Devices did not expand when interacting with large white clots (6 mm). They remained constrained between the clot and the vessel wall. During retrieval, **devices slid over the clot without capturing it**. None of the devices tested could penetrate and remove large white thrombi.”*

Solitaire*:	0/5	Trevo:	0/5
Embotrap*:	0/5	Eric:	0/5
Preset*:	0/5	Preset LT:	0/5
Catch*:	0/5	Separator 3D:	0/5
Revive*:	0/5	Mindframe:	0/5


Uniquely designed to **CAPTURE**
ALL TYPES OF CLOT INSIDE THE
DEVICE STRUCTURE



the **DROP ZONES™**
allow the capture of
large, organized thrombi
within the NeVa basket

DROP ZONES™

2 or more Drop Zones offset at 90° work by acting as clot pockets: entry points to capture thrombi inside



CLOSED DISTAL TIP

Clot gets inside,
clot stays inside!


SMART MARKERS

2 per drop zone,
for real-time feedback
during retrieval

BALANCED DESIGN

Optimized radial force
balanced with large
openings & closed ends

Clinical Experience with the Nevia™ thrombectomy device



97% RECANALIZATION SUCCESS WITH 1.2 PASSES ACROSS ALL CLOT TYPES



Clot Type	Soft	Hard	Ultra Hard	All Clots
Clot morphology	Whole Blood "RED" Clot	Plasma Rich "WHITE" Clot	Clot modeled from ONYX 500	RED, WHITE and ONYX 500
N =	19	5	11	35
Length of clots - mm	10-40	6-12	4-12	4-40
1 st Pass TICI 3	84%	60%	55%	71%
Final TICI 3	89%	NR	82%	83%
Final TICI 2b/3	100%	100%	91%	97%
Average # of passes for final recanalization	1,05	1,00	1,63	1,23

IN VITRO RESULTS

The same set up was used to test NeVa with large, hard clots:

2. Machi P, et al., 2018

- **NeVa: 6/10 successful complete removals of white thrombi ≥ 6 mm**
- The main difference versus other devices: Drop Zones, allowing for the lateral integration of clots
- The NeVa devices demonstrated good wall apposition and resisted stretch related deformation during retrieval
- NeVa radial pressure values appear slightly higher in comparison to those of other stent retrievers

1. Machi P, et al., 2016

*““Devices did not expand when interacting with large white clots (6 mm). They remained constrained between the clot and the vessel wall. During retrieval, **devices slid over the clot without capturing it**. None of the devices tested could penetrate and remove large white thrombi.”*

Solitaire*:	0/5	Trevo:	0/5
Embotrap*:	0/5	Eric:	0/5
Preset*:	0/5	Preset LT:	0/5
Catch*:	0/5	Separator 3D:	0/5
Revive*:	0/5	Mindframe:	0/5

1. Machi P, et al., “Experimental evaluation of stent retrievers’ mechanical properties and effectiveness”, Journal of NeuroInterventional Surgery. 2016; 0:1-7

* Some versions of these devices were able to minimally displace clots 1-2 times in 5, without succeeding in removing them

2. Machi P, et al., “Experimental evaluation of the NeVa™ thrombectomy device a novel stent retriever conceived to improve efficacy of organized clot removal”, Journal of Neuroradiology. April 2018

CLINICAL RESULTS
MULTI CENTER EXPERIENCE

30
patients

	2B/3	2C/3	
First Pass	19/30 → 63%	14/30 → 47%	Average # of passes for final recan → 1.7
Final Recanalization	28/30 → 93%	19/30 → 63%	

Patient Outcomes

Mean NIHSS @ 24hr : 7

90 day mRS < 2 : 53%

Zero NeVa related adverse events & sICH

- NeVa was effective with both balloon guide and local aspiration strategies
- In the 40 passes where the info was available: 70% clot incorporation into device basket

118
patients

	2B/3	2C/3	
First Pass	61/80 → 56.8%	46/80 → 44.9%	Median # of passes for final recan → 1 (IQR 1–2)
Final Recanalization	77/80 → 95.8%		

Favorable functional outcome (mRS ≤ 2):

- 53% in the “first-pass” subgroup
- 42.4% in the total patient population.

Procedure related complications:

- Symptomatic ICH: 3.3%
- Asymptomatic ICH: 13.6%
- Embolization into new territory: 1.7%
- Dissection that did not require stenting: 1.7 %

Flow Control strategies preferred:

- 92.4% of cases done with distal aspiration (Solombra: Aspiration catheter + NeVa)
- Balloon Guide Catheter used only in 13.6% of cases

CLINICAL RESULTS
SINGLE CENTER EXPERIENCE

80
patients

2B/3

2C/3

First Pass

61/80 → **76.3%**

46/80 → **57.5%**

Average # of passes
for final recan →

Final Recanalization

77/80 → **96.3%**

72/80 → **90.0%**

1.6

Safety Data

No NeVa related intra-procedural complications

- Difficulty of navigation via the MC due to significant tortuosity: 3.7%
- Mild asymptomatic SAH at 24hr follow up: 2.5%
- Asymptomatic petechial hemorrhage (HT1): 15%
- HT-2: 1.2%
- Parenchymal hemorrhage (PH-1): 1.2%

Patient Outcomes

Mean NIHSS @ 24 hours: **7** (Range: **0-33**)

30-day mRS ≤ 2 **51/80 (65%)**

obtained in 64 of 80 patients


90-day mRS ≤ 2 **55/80 (68.7%)**

obtained in 66 of 80 patients

Nevia™

DESIGNED FOR 1ST PASS SUCCESS
WITH ALL CLOT TYPES


DROP ZONE
TECHNIQUE™



NEVA THROMBECTOMY DEVICE

DROP ZONES™

2 or more Drop Zones offset at 90° work by acting as clot pockets:
for lateral integration of clot inside the device structure



BALANCED DESIGN

Optimized radial force
balanced with large openings
& closed ends

SMART MARKERS

2 per Drop Zone,
for real-time feedback during
retrieval


CLOSED DISTAL TIP

Clot gets inside,
clot stays inside!

1. CHOOSING THE CORRECT NEVA SIZE
2. POSITIONING NEVA
3. RETRIEVING NEVA

CHOOSING THE CORRECT NEVA SIZE

NeVa 4.0 x 22 mm, 2 Drop Zones



Vessel diameters **2.0 – 3.5 mm**
Ideal for **MCA, ACA, PCA**
Compatible with: 0.021" MC
Full length: 39 mm



NeVa 4.5 x 29 mm, 3 Drop Zones



Vessel diameters **2.0 – 4.5 mm**
Ideal for **ICA Tip, Proximal MCA, Basilar**
Compatible with: 0.021" MC
Full length: 46 mm



NeVa 4.5 x 44 mm, 5 Drop Zones




Vessel diameters **2.0 – 4.5 mm**
Ideal for **ICA Tip, Proximal MCA**
Ideal for **long occlusions, high clot burden**
Compatible with: 0.021" MC
Full length: 61 mm



NeVa 6.0 x 44 mm, 3 Drop Zones, Flow Restoration Zone



Vessel diameters **3.5 – 6.0 mm**
Ideal for **Proximal Carotid Artery**
Ideal for **large, proximal arteries**
Compatible with: 0.027" MC
Full length: 63 mm



IDEAL POSITIONING

Ideally we want:

Multiple Drop Zones to interact with clot

And we need to:

Balance the benefit and risk of distal placement

To achieve this:

Deploy NeVa with the 1st or 2nd marker at the edge of the occlusion



IDEAL POSITIONING

EXAMPLE WITH NEVA 4.5 X 29 MM

Ideally we want:
And we need to:
To achieve this:

Multiple Drop Zones to interact with clot
Balance the benefit and risk of distal placement
Deploy NeVa with the 1st or 2nd marker at the edge of the occlusion

Example:
4.3 mm vessel,
1 cm clot




The first or second marker
at the edge of the occlusion


Minimum 2 Drop Zones
interacting with clot



A Drop Zone deployed proximal to
clot: **LESS THAN IDEAL**




EXPECT INITIAL ANCHORING AFTER 1CM OF UNSHEATHING




USING THE DROP ZONES TO INCORPORATE CLOTS

Drop Zone markers will get compressed when NeVa is passing next to a hard, calcified clot in the vascular system

Markers compressed together:
YOU MAY BE ADJACENT TO A HARD CLOT:
SLOW DOWN!



Markers spring open:
YOU MAY NOW BE AT THE PROXIMAL EDGE OF
THE HARD CLOT:
THE DROP ZONE IS ON THE CLOT



USING THE DROP ZONES TO INCORPORATE CLOTS


Slow and alert retrieval is recommended



POSITION CORRECTLY

Deploy the device with the 1st or 2nd marker at the occlusion

You do not need wait



START SLOW PULL

Apply slow & gentle vessel straightening traction



WATCH THE SMART MARKERS

Watch Drop Zone markers, observe if one of the pairs is compressing on one another

GENERAL RECOMMENDATIONS

The following recommendations are not specific to NeVa

MICROCATHETER CONSIDERATIONS

Choose

micro-catheters with
sufficient distal support,
especially in tortuous cases


Flush

NeVa
before insertion

Release tension


on the micro-catheter
just before starting the
deployment (unsheathing)

AT SHARP ANGLES



OPTION 1:


Bring NeVa proximally towards the DAC and align NeVa with the tip of your DAC



OPTION 2:


Use NeVa as an anchor, and drive up your DAC

1. Straightening the anatomy eases retrieval
2. Avoids clot fragmentation
3. Aspiration via DAC will be more efficient




WHEN LARGE CLOT BURDEN IS SUSPECTED: PARTIAL RETRIEVAL TECHNIQUE

After deploying NeVa, bring the DAC tip up to the proximal marker



Remove excess tension from the DAC and slowly retrieve NeVa. If significant resistance is encountered, stop retrieval. Clot is likely partially incorporated and trapped between stent and DAC.



- Tighten the RHV of DAC around the MC and retrieve the whole system together (DAC+MC+NeVa) while gently aspirating

RECAP

1. Choose correct NeVa size
2. Flush before use
3. Position NeVa as recommended
4. Release tension on the microcatheter before unsheathing
5. Apply slight forward pressure on the pusher wire during unsheathing until NeVa anchors in the vessel
6. Take your time in deploying NeVa to take advantage of the Drop Zones
7. Retrieval should be slow and gentle: Watch the Drop Zone markers
8. At sharp angles: use NeVa as anchor and drive up your DAC to straighten the anatomy or pull NeVa down to the DAC
9. When combined with aspiration: If you feel a resistance, take the whole system out together to avoid clot shearing/ tooth paste effect

Question & Answer Session



THANK YOU