



Eurolimus™
Product Overview

pre-market dossier for active customers
confidential



EUROLIMUS™

PRODUCT OVERVIEW




Eurolimus™, Sirolimus Eluting Coronary Stent System

- Treatment of **symptomatic ischemic coronary heart disease** (from *de novo lesions & restenosed lesions of the coronary arteries*)
- Treatment of **acute or suspected occlusions** (for patients with unsuccessful balloon dilatation)
- Prevents restenosis, assures **homogeneous drug distribution** and uniform release kinetics
- **Long term efficacy & extremely low rates of thrombosis** (acute, subacute, late & very late)



EUROLIMUS

Sirolimus Eluting Coronary Stent System



Stent Specifications

Type of design	Open-cell, 3 interlinks
Design detail	9 crowns
Material	Cobalt Chromium L605
Expansion range	2.25 mm – 4.00 mm
Strut thickness	0.0026" (65 µm)
Strut width (main segment)	0.0028" (72 µm)
Strut width (interlink)	0.0023" (58 µm)
Foreshortening	< 2.5 %
Mechanical recoil	< 6 %
Metal coverage	< 13.6 %
Matrix thickness	3 – 5 µm
Drug / Polymer	Sirolimus / PLGA
Drug load	1.4 µg / mm ²

Delivery System Specifications


Usable length	138 cm	
Length guidewire lumen	27 cm	
Recommended guide catheter	5F (min. I.D. 0.056")	
Material	Balloon	Proprietary polyamide compound
	Distal shaft	Polyamide, multilayer tube Hydrophilic coating
	Proximal shaft	Stainless steel, PTFE coated
Shaft size	Proximal	1.9 F
	Distal	2.8 – 3.0 F, depending on balloon-Ø
Folding	3-fold balloon	
Marker bands	Embedded Platinum / Iridium	
Tip entry profile	0.017" (0.43 mm)	
Max. guidewire	0.014" (0.36 mm)	



DESIGN & COATING TECHNOLOGY



Low Strut thickness for smooth interventions




Crimped Profile

Overexpansion capabilities


Asymmetric Coating

Expansion Diameter	Crimped Profile	Stent ID at RBP	Post Dilatation Limit
2.25	0.87 mm	2.49 mm	3.50 mm
2.50	0.91 mm	2.82 mm	3.50 mm
2.75	0.94 mm	3.03 mm	3.50 mm
3.00	0.98 mm	3.25 mm	4.00 mm
3.25	1.00 mm	3.53 mm	4.00 mm
3.50	1.01 mm	3.80 mm	4.50 mm
4.00	1.03 mm	4.31 mm	4.50 mm




Cross Section Diagram

Sirolimus mTOR inactivation




Limus Analogs Comparison




Molecular structure of Limus drugs is not altered in the two essential binding sites that interact with target structures mTOR and FKBP-12
Drug's mode of action is essentially the same

Eurolimus® drug release




Coating Structure



Drug release covers entire restenotic cascade (100 days)

SEM of coating 600 x

Polymer Matrix




Lactic acid glycolic acid

- Initial burst phase followed by second phase of progressive release through thicker drug depleted layer.
- PLGA 85/15 degradable polymer formulation regarded as **non toxic and biocompatible**
- Both substances are **fully decomposed to water and carbon dioxide** between 10-13 weeks

Superior Sterilization for optimal performance & shelf life

Electron Beam sterilization utilizes direct current or linear accelerator to emit high energy electrons.



- Sterilization is performed completely dry
- 0% degradation of the polymer
- 100 % repeatability of drug release & working mechanism
- Guarantees synchronized drug & polymer release



3 Year
Shelf Life



EUROLIMUS™

CLINICAL PROGRAM



Two completed Studies on real world patients undergoing percutaneous coronary intervention in Europe

Milan,
Italy



171 Patients
Follow-up 25.6 ± 5.7 months

Retrospective Single Center

Thessaloniki,
Greece



196 Patients
Follow-up 27.8 ± 5.2 months

Retrospective Single Center

End Points



Target Lesion Failure
MACE
Stent Thrombosis



Study Endpoints

MACE



Composite

- Repeat Revascularization (any, including all target & non-target vessels)
- All-cause-mortality
- Any MI (Myocardial Infarction)

Stent
ThrombosisAs defined by the
Academic Research
Consortium

TLF



Composite

- Cardiac Death
- Clinically Driven TLR (Target Lesion Revascularization)
- Target vessel MI (Q-wave or non Q-wave Myocardial Infarction), CABG (Coronary Artery Bypass Graft)

COHORT DEMOGRAPHICS

PATIENT STATUS & MEDICAL HISTORY

Milan, Italy

Thessaloniki, Greece

Demographics	Total	%
N	171	
Age (mean)	70.0 ± 11.9	
Male	138	80.7
BMI	27.2 ± 5.1	
Diabetic	37	25.0
History of MI	35	20.7
History of stroke TIA	6	3.6
Smoker (current)	35	21.7
Hypertension	127	80.9
Hyperlipoproteinemia	71	42.0
Peripheral Vascular Disease	6	3.6

Demographics	Total	%
N	196	
Age (mean)	66.7 ± 10.7	
Male	159	79.8
BMI	27.6 ± 3.3	
Diabetic	54	27.6
History of MI	37	18.8
History of stroke TIA	12	6.1
Smoker (current)	119	60.7
Hypertension	85	43.4
Hyperlipoproteinemia	107	54.6
Peripheral Vascular Disease	9	4.6

COHORT DEMOGRAPHICS

PATIENT STATUS & MEDICAL HISTORY

Milan, Italy

Cardiac Status	Total	%
Chronic stable angina	59	34.5
Silent Ischemia	26	15.2
Unstable angina	26	30.2
MI (Stemi)	38	63.3
MI (Nstemi)	22	36.7

Cardiac History	Total	%
Prior PCI	45	26.3
Prior CABG	29	16.9
Prior MI (> 72h)	25	14.6

Coronary Vessel Disease	Total	%
1 vessel	134	78.4
2	30	17.5
3	7	4.1

Thessaloniki, Greece

Cardiac Status	Total	%
Stable angina	84	42.9
Unstable angina	27	13.8
MI (Stemi)	42	21.4
MI (Nstemi)	43	21.9

Cardiac History	Total	%
Prior PCI	18	9.2
Prior CABG	8	4.1
Prior MI (> 72h)	16	8.2

Coronary Vessel Disease	Total	%
1 vessel	161	82.2
2	32	16.3
3	3	1.5

Location of target lesions

Milan, Italy


- LAD: 88, 38.4%
- RCA: 92, 40.2%
- LCX: 46, 20.1%
- LM: 3, 1.3%

Total of 229 lesions

Thessaloniki, Greece

- LAD: 116, 54.8%
- RCA: 80, 31.6%
- LCX: 57, 22.5%
- LM: 0, 0.0%

Total of 253 lesions



Total distribution of treated target lesion distribution

COHORT DEMOGRAPHICS

PROCEDURE CHARACTERISTICS

Milan, Italy

Left ventricular ejection	Total	%
Fraction	56.6 ± 9.8	34.5
Fraction ≤ 35 %	11	6.4
Amount of Lesions	Total	%
Average per patient	1.34	
1 lesion	124	72.5
2 lesions	36	21.1
3 lesions	11	6.4
Lesion Characteristics	Total	%
RVD	3.38 ± 0.51	
Small vessels RVD (< 2.75 mm)	25	11.1
Diameter Stenosis	83.01 ± 13.31	
Total Occlusions	27	11.79

Thessaloniki, Greece

Cardiac Status	Total	%
Fraction	59.7 ± 7.6	
Amount of Lesions	Total	%
Average per patient	1.29	
1 lesion	152	77.6
2 lesions	35	17.9
3 lesions	6	1.5
4 lesions	3	1.5
Lesion Characteristics	Total	%
RVD	3.09 ± 0.34	
Small vessels RVD (< 2.75 mm)	16	6.2
Diameter Stenosis	86.6 ± 10.7	
Total Occlusions	28	11.1

Follow up results (25.6 ± 5.7 months)

TLF

MACE

Thrombosis

	Total	%
Target Lesion Failure	2	0.9
Cardiac Death	0	0
MI – Attributed to target vessel	0	0
Q-wave	0	0
Non Q-wave	0	0
Clinical Driven TLR	2	0.9
Solved by PCI	2	0.9

	Total	%
MACE	8	3.9
Cardiac Death	0	0
Non Cardiac Death	3	1.7
MI (all)	0	0
TLR	2	0.9
By PCI	2	0.9
TVR	2	0.9
By PCI	2	0.9
Non TVR	1	0.4
By PCI	1	0.4

	Total	%
Thrombosis (all)	0	0
Acute Stent thrombosis ¹	0	0
Subacute stent thrombosis ²	0	0
Late stent thrombosis ³	0	0
Very late stent thrombosis ⁴	0	0

Events

	Total	%
Ischemic Events	0	0
Hemorrhagic/Vascular Events	1	0.6
Hematological Dyscrasia	0	0
Other Complications⁵	6	3.5

¹ ≤ 24 hours ² > 24 hours ≤ 30 days ³ > 30 days ⁴ > 1 year

⁵ Dialysis (FU 12 months), Pulmonary Oedema, acute heart failure (1 month post-PCI), Intramural hematoma of aortic root from coronary dissection (at discharge), acute kidney failure (FU 12 months), Mitral-clip procedure (FU 12 months)

Follow up results (25.6 ± 5.7 months)

TLF

MACE

Thrombosis

	Total	%
Target Lesion Failure	9	4.6
Cardiac Death	2	1.0
MI – Attributed to target vessel	2	1.0
Q-wave	1	0.5
Non Q-wave	1	0.5
Clinical Driven TLR	5	2.6
Solved by PCI	5	2.6

	Total	%
MACE	13	8.4
Cardiac Death	2	3.0
Non Cardiac Death	4	2.0
MI (all)	3	2.6
MI – target vessel	1	0.9
MI – non-target vessel	2	1.7
TLR	5	2.0
By PCI	5	2.0
TVR	1	0.4
By PCI	1	0.4
Non TVR	1	0.7
By PCI	1	0.7

	Total	%
Thrombosis (all)	1	0.5
Acute Stent thrombosis ¹	0	0
Subacute stent thrombosis ²	1	0.5
Late stent thrombosis ³	0	0
Very late stent thrombosis ⁴	0	0

Events

	Total	%
Ischemic Events	0	0
Hemorrhagic/Vascular Events	1	0.5
Hematological Dyscrasia	0	0
Other Complications	0	0

¹ ≤ 24 hours ² > 24 hours ≤ 30 days ³ > 30 days ⁴ > 1 year

THANK YOU FOR YOUR
ATTENTION!



Eurocor

Eurocor Tech GmbH
In den Dauen 6a
53117 Bonn
Germany

Phone: +49 (0)228/20 15 0-0
Fax: +49 (0)228/20 15 0-5
info@eurocor.de
eurocor.de

Eurocor Tech GmbH does not assume responsibility for editorial or clerical errors. All rights to the mentioned studies lie with the respective authors.