mindray

Multi Sera Calibrator

Package sizes

20×3 mL, 10×3 mL

Intended use

Multi Sera Calibrator is used for calibration of quantitative determination of routine chemistry analytes on Mindray BS measurement system.

Summary

Multi Sera Calibrator contains the following analytes:

Albumin (ALB)	Alkaline Phosphatase (ALP) Aspartate	Alanine Aminotransferase (ALT)
α -Amylase (α -AMY)	Aminotransferase (AST)	Direct Bilirubin (Bil-D)
Total Bilirubin (Bil-T)	Calcium (Ca)	Total Cholesterol (TC)
Creatine Kinase (CK)	Creatinine (CREA)	Glucose (Glu)
Gamma-Glutamyltransferase (γ -GT)	α -Hydroxybutyrate Dehydrogenase (α -HBDH)	Lactate Dehydrogenase (LDH)
Magnesium (Mg)	Phosphorus (P)	Total Protein (TP)
Triglycerides (TG)	Urea (UREA)	Uric Acid (UA)
Cholinesterase (CHE)	Lipase(LIP)	

By calibrating the selected analytes, Mindray BS-series chemistry analyzers can generate valid calibration factors, and then calculate the concentrations or activities of each sample automatically after analysis. Mindray BS measurement system is composed of Mindray BS-series chemistry analyzers, Mindray reagent kits, calibrators and control materials.

Components

Multi Sera Calibrator is a lyophilized calibrator based on human serum. The concentrations or activities of the calibrator components are lot-specific.

Warnings and Precautions

1. For in vitro diagnostic use only.
2. Professional or trained person is required. Please take the necessary precautions for handling all laboratory reagents.
3. The calibrator values are lot-specific with the matched models listed in the value sheet.
4. Calibration is recommended when reagent lot is changed or certain maintenance or troubleshooting procedure is performed.

- Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.
- The calibrator was tested with Elisa methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen¹.
- Material safety data sheet is available on request.
- Disposal of all waste material should be in accordance with local guidelines.

Preparation

- Take the vial out of the freezer and allow reaching room temperature.
- Tap the vertically positioned vial gently and ensure that the lyophilized material is at the bottom of the vial.
- Remove the screw cap and rubber stopper carefully, avoiding any loss of lyophilizate.
- Reconstitute by adding 3.0 mL distilled/deionized water to the side of the vial slowly and exactly.
- Carefully replace the rubber stopper and leave to stand for 30 minutes at room temperature.
- During the standing time, mix contents by inverting several times and swirling gently to ensure that all components of lyophilizate are dissolved. Avoid the formation of foam.
- According to the calibration procedure described in the instrument manual, use the specific value listed in the value sheet to set calibration parameter and execute calibration.

Storage and Stability

The calibrator is stable up to expiration date indicated on the label when stored in unopened vial at 2~8°C and protected from light. Once reconstituted, please store calibrator tightly capped when not in use. Once thawed to room temperature, the calibrator should not be re-frozen. The stability information is as below.

Temperatures of storage	Stability		
	Universal components	Total Bilirubin	Direct Bilirubin
15~25°C	8 hours	6 hours	3 hours
2~8°C	2 days	1 day	8 hours
-25~-15°C	4 weeks	2 weeks	2 weeks

Materials required but not provided

1. Distilled/deionized water.
2. Mindray BS Chemistry analyzer and reagent kits.
3. General laboratory equipment.

Traceability of Calibrator

The calibrator values assigned by Mindray standard transfer procedure and routine method are listed in the value sheet. The traceability process is based on ISO 17511² and ISO18153³, each analyte (measurand) in this calibrator is traceable to the available reference method and (or) reference material, or manufacturer's reference standard. Traceability information is given below. Details of traceability are available on request.

Analytes	Reference methods/materials
UREA, Bil-T, Glu, TP, TG, TC,UA, Ca, Mg	International Reference Measurement Method or/and Standard Reference Material from NIST
LDH, γ -GT, CK, α -AMY	IFCC Primary Reference Procedures
ALB	ERM-DA470K
ALT, AST, α -HBDH, ALP, Bil-D, P,CHE,LIP	Manufacturer's Selected Measurement Procedure
CERA,	NIST definitive method ID-LC/MS/MS (NIST SRM 914a, NIST SRM 967)

Quality control

Use Mindray control as sample to monitor the status of calibration. The control result should be within the defined range as shown in the control value sheet. If the control result falls outside the range, the measurement system should be checked. For instance, position of reagent or sample placed on the analyzer; expiration date or storage condition of the calibrator, reagent and control; the parameter setting in software; or performance of the analyzer.

References

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.
2. ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.
3. ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.

In Vitro Diagnostic
medical device

Batch Code

European
Conformity

Authorized representative in
the European Community

Use By

Consult
Instructions for use

Biological
Risks

Temperature
Limit

Manufacturer

Catalogue
number

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Referências

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

2.ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.

3.ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

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Riferimenti

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing

Office; Washington: 2007

- 2.ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.
- 3.ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.

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Indirizzo: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

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Referanslar

- 1.Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.
- 2.ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.
- 3.ISO 18153:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials.

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Üretici:Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray building, Keji 12th Road South, High-Tech Industrial Park,Nanshan, Shenzhen, 518057 P.R.China

E-posta Adresi:service@mindray.com

Web sitesi:www.mindray.com

Tel:+86-755-81888998

Faks:+86-755-26582680

AT Temsilcisi:Shanghai International Holding Corp.GmbH(Europe)

Adres:Eiffestraße 80, Hamburg 20537, Germany

Tel:0049-40-2513175

Faks:0049-40-255726

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© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır