# *Agarose Gel Kits:*

*SPE*

*SPE-Split-Beta*

*alc.Hb*

*LPE*

*UPE*

*IFX*

*MPS*

*2014-07-01*

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---

The BIOTEC-FISCHER **Agarose Gel Kit SPE** is intended for the separation of human serum proteins into the traditional 5 fractions: Albumin, Alpha 1, Alpha 2, Beta and Gammaglobulins by agarose gel electrophoresis. The separated proteins are fixed within the agarose gel and stained with a special staining solution. The transparency of the agarose gel allows an easy and reliable visual as well as densitometrical analysis.

The agarose gel and reagents are specially designed for best results in serum protein electrophoresis. The BIOTEC-FISCHER Agarose Gel Kit SPE is optimized for use with MAESTRO 101, FLILPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

### The BIOTEC-FISCHER Agarose Gel Kit SPE contains:

- 10 agarose gels for serum protein elektrohoresis
- blotting paper for agarose gels / templates
- sample application templates
- buffer solution, concentrated
- staining solution, concentrated

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---

The BIOTEC-FISCHER **Agarose Gel Kit SPE-Split-Beta** is intended for the separation of human serum proteins into 6 fractions: Albumin, Alpha-1, Alpha-2, Beta-1, Beta-2 and Gammaglobulins by agarose gel electrophoresis. The separated proteins are fixed within the agarose gel and stained with a special staining solution. The transparency of the agarose gel allows an easy and reliable visual as well as densitometrical analysis.

The agarose gel and reagents are specially designed for best results in serum protein electrophoresis. The BIOTEC-FISCHER Agarose Gel Kit SPE-Split-Beta is optimized for use with MAESTRO 101, FILIPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

### The Agarose Gel Kit SPE-Split-Beta contains:

- 10 agarose gels for Split-Beta serum protein electrophoresis
- blotting paper for agarose gels / templates
- sample application templates
- buffer solution, concentrated
- staining solution, concentrated

### Article Codes:

- 370602-001 Kit for 10 SPE ( $\beta 1 + \beta 2$ ) on PHERO-gel-A
- 370602-002 Kit for 100 SPE ( $\beta 1 + \beta 2$ ) on PHERO-gel-A
- 370602-101 Kit for 10 SPE ( $\beta 1 + \beta 2$ ) on PHERO-gel-A for MAESTRO 101
- 370602-102 Kit for 100 SPE ( $\beta 1 + \beta 2$ ) on PHERO-gel-A for MAESTRO 101

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---

The BIOTEC-FISCHER **Agarose Gel Kit alc.Hb** is intended for the separation of the normal hemoglobins (HbA and HbA2) and for the detection of the major hemoglobin variants such as HB F, Hb S and Hb C by alkaline agarose gel electrophoresis. The separated proteins are fixed within the agarose gel and stained with a special staining solution. The transparency of the agarose gel allows an easy and reliable visual as well as densitometrical analysis.

The agarose gel and reagents are specially designed for best results in alc. hemoglobin electrophoresis. The BIOTEC-FISCHER Agarose Gel Kit alc.Hb is optimized for use with MAESTRO 101, FLILPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

The BIOTEC-FISCHER Agarose Gel Kit alc.Hb contains:

- 10 agarose gels for hemoglobin electrophoresis
- blotting paper for agarose gels / templates
- sample application templates
- buffer solution, concentrated
- staining solution, concentrated

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---

The BIOTEC-FISCHER **Agarose Gel Kit LPE** is intended for the for the separation and classification of human serum lipoproteins into the four traditional fractions:  $\beta$ , pre- $\beta$ ,  $\alpha$  und Chylomikronen by agarose gel electrophoresis. The separated lipoproteins are fixed within the agarose gel and stained with a special staining solution. The transparency of the agarose gel allows an easy and reliable visual as well as densitometrical analysis.

The agarose gel and reagents are specially designed for best results in lipoprotein electrophoresis. The BIOTEC-FISCHER Agarose Gel Kit LPE is optimized for use with MAESTRO 101, FLILPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

The BIOTEC-FISCHER Agarose Gel Kit LPE contains:

- 10 agarose gels for lipo protein electrophoresis
- blotting paper for agarose gels / templates
- sample application templates
- buffer solution, concentrated
- staining solution, concentrated

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---

The BIOTEC-FISCHER **Agarose Gel Kit UPE** is intended for the fractionation and classification of main urinary proteins by agarose gel electrophoresis to type renal failure. The separated proteins are fixed within the agarose gel and stained with a special high-sensitive staining solution. The transparency of the agarose gel allows an easy and reliable visual as well as densitometrical analysis.

The agarose gel and reagents are specially designed for best results in urine protein electrophoresis. The BIOTEC-FISCHER Agarose Gel Kit UPE is optimized for use with MAESTRO 101, FILIPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

The BIOTEC-FISCHER Agarose Gel Kit UPE contains:

- 10 agarose gels for urin protein electrophoresis
- blotting paper for agarose gels / templates
- sample application templates
- buffer solution, concentrated
- staining solution, concentrated

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The BIOTEC-FISCHER **Agarose Gel Kit IFE** is intended for detection and identification of gammopathies in human sera or other human biological fluids, by agarose gel electrophoresis, followed by an immunofixation. The transparency of the agarose gel and the use of high quality antisera allows an easy and reliable typing of monoclonal bands by the simultaneous immunofixation of the anti human, IgG, IgA, IgM, Kappa and Lambda tracks and staining with special staining solution.

The agarose gel and reagents are specially designed for best results in immunofixation. The BIOTEC-FISCHER Agarose Gel Kit IFE is optimized for use with MAESTRO 101, FILIPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

### The BIOTEC-FISCHER Agarose Gel Kit IFE contains:

- 10 agarose gels for immunofixation
- blotting paper for agarose gels / templates
- blotting paper for samples application / antisera
- application templates for samples / antisera
- buffer solution, concentrated
- staining solution, concentrated
- antisera

### Article Codes:

- 376601-001 Kit for 1 IFE on agarose PHERO-gel-A
- 376601-002 Kit for 10 IFE on agarose PHERO-gel-A
- 376601-101 Kit for 1 IFE on agarose PHERO-gel-A for MAESTRO 101
- 376601-102 Kit for 10 IFE on agarose PHERO-gel-A for MAESTRO 101

(Special manufacture of gels and gel kits to individual requirements possible. Please contact us).

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---

The BIOTEC-FISCHER **Agarose Gel Kit MPS** is intended for the separation of Glykosaminoglycans (GAGs) from urine into the main fractions Chondroitinsulfat, Dermatansulfat, Heparansulfat and Keratansulfat by agarose gel electrophoresis. The separated GAGs are fixed within the agarose gel and stained with a special staining solution. The transparency of the agarose gel allows an easy and reliable visual as well as densitometrical analysis.

The agarose gel and reagents are specially designed for best results in serum protein electrophoresis. The BIOTEC-FISCHER Agarose Gel Kit MPS is optimized for use with MAESTRO 101, FILIPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

The BIOTEC-FISCHER Agarose Gel Kit MPS contains:

- 10 agarose gels for MPS electrophoresis
- blotting paper for agarose gels / templates
- sample application templates

Article Codes:

379601-001 Kit for 100 MPS on agarose PHERO-gel-A

(Special manufacture of gels and gel kits to individual requirements possible. Please contact us).

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---

The BIOTEC-FISCHER **CAF Kit SPE** is intended for the separation of human serum proteins into the traditional 5 fractions: Albumin, Alpha 1, Alpha 2, Beta and Gammaglobulins by electrophoresis on cellulose acetate foils (CAF). The separated proteins are fixed within the CAF and stained with a special staining solution. The special developed transparency solution guarantees excellent transparency of CAF for easy and reliable visual as well as densitometrical analysis.

The reagents are specially designed for best results in serum protein electrophoresis on CAF. The BIOTEC-FISCHER CAF Kits SPE are optimized for use with MAESTRO 101, FLILPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

### The BIOTEC-FISCHER CAF-Kit SPE contains:

- celluloseacetat foils PHERO-cel for serum protein electrophoresis
- buffer solution
- staining solution
- destainer, concentrated
- transparency solution
- blotting paper
- glass plates

### Article Codes:

- 370101 Kit for 800 SPE on CAF PHERO-cel perforated
- 370102 Kit for 200 SPE on CAF PHERO-cel perforated
- 370103 Kit for 800 SPE on CAF PHERO-cel unperforated
- 370104 Kit for 200 SPE on CAF PHERO-cel unperforated

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---

The BIOTEC-FISCHER **CAF Kit alc.Hb** is intended for the separation of the normal hemoglobins (HbA and HbA2) and for the detection of the major hemoglobin variants such as HB F, Hb S and Hb C by alkaline hemoglobin electrophoresis on cellulose acetate foils (CAF). The separated proteins are fixed within the CAF and stained with a special staining solution. The special developed transparency solution guarantees excellent transparency of CAF for easy and reliable visual as well as densitometrical analysis.

The reagents are specially designed for best results in alc. hemoglobin electrophoresis on CAF. The BIOTEC-FISCHER CAF Kit alc.Hb is optimized for use with MAESTRO 101, FLILPO and electrophoresis tanks BFA and 1000 but can be used also in several electrophoresis systems of other brands.

### The BIOTEC-FISCHER CAF-Kit alc.Hb contains:

- cellulose acetate foils PHERO-cel for hemoglobin electrophoresis
- buffer solution
- staining solution
- destainer, concentrated
- transparency solution
- blotting paper
- glass plates

### Article Codes:

- 373101 Kit for 400/800 alk.Hb on PHERO-cel perfor.
- 373102 Kit for 100/200 alk.Hb on PHERO-cel perfor.
- 373103 Kit for 400/800 alk.Hb on PHERO-cel unperf.
- 373104 Kit for 100/200 alk.Hb on PHERO-cel unperf.

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Adresse **BIOTEC-FISCHER GmbH**