Hemolyzer® 3 NG Reagent Pack

► Product Information ◀

Product name	Product code	Number of tests
H3-Compact 100	HE3701	100
H3-Compact 500	HE3705	500

► Product Content ◀

	H3-Compact 100	H3-Compact 500
Diluent	1100 mL	3550 mL
Lyse	350 mL	750 mL
Cleaner	400 mL	1150 mL
Instructions for use	1 pc	1 pc
QR Code	1 pc	1 pc

► Intended Use and Test Principle ◀

H3-Compact 100 and H3-Compact 500 are dedicated set of reagents intended to be used with Hemolyzer® 3 NG semi-automated hematology analyzers. The reagent pack contains the following reagents: Diluent, Lyse and Cleaner.

Diluent reagent is a buffered, stabilized, and micro-filtered electrolyte solution for automated dilution of human blood samples, quantitative and qualitative determination of red blood cells (RBC), white blood cells (WBC), white blood cells subpopulations (lymphocytes [LYM], mid-size cells [MID], granulocytes [GRA]), platelets (PLT) and measurement of hemoglobin (HGB) concentration.

Lyse reagent is a stabilized and micro-filtered lysing agent for stratomolysis of red blood cells (RBC), for quantitative determination of white blood cells (WBC), for WBC 3-part differentiation (lymphocytes [LYM], mid-size cells [MID], granulocytes [GRA]) and hemoglobin (HGB) concentration measurement of human blood.

Cleaner reagent is a buffered, stabilized, and micro-filtered electrolyte solution used for conditioning and cleaning of specific tubing sections and chambers.

The reagent is intended for laboratory professional use as an in vitro diagnostic medical device.

► Reagent Preparation ◀

The product is ready to use, please follow the installation procedures of this instructions for use.

► Warnings and Precautions ◀

- Read instructions for use carefully before use.
- The product is environmentally friendly, azide-free, and does not contain harmful ingredients. This product is classified as non-hazardous in compliance with Regulation (EC) No 1272/2008.
- Do not use the product beyond the expiry date.
- Do not use the product if the primary package is damaged or there is visible sign of chemical deterioration of the liquid inside. In both cases the product should be replaced.
- Store at the recommended temperature range.
- Do not use the product if it has been frozen or kept at excessive heat. Use the product at room temperature (18-30°C). If it has been stored below 18°C or above 30°C, leave it at normal room temperature for 3 hours before use.
- Users are advised to wear approved protective clothing when handling chemical products: lab coat, gloves, and eye protection.
- Observe standard laboratory precautions for use and follow national or local health and safety guidelines.
- Avoid contact with eyes, skin, and clothing. In case of contact with eyes or skin, rinse immediately with plenty of water.
- Please refer to the product's Safety Data Sheet.
- Please report any serious incident involving the product to the manufacturer and the competent local authority.

► Product Composition ◀

	Diluent	Lyse	Cleaner
Sodium chloride	< 1.5 %	-	< 1.5 %
Surfactants	-	< 3.5 %	-
Buffers	< 1.0 %	< 1.0 %	< 1.0 %
Stabilizers	< 0.5 %	< 0.5 %	< 0.5 %
Preservatives	< 0.5 %	< 0.5 %	< 0.5 %
Ion-free water	> 96.5 %	> 94.5 %	> 96.5 %

► Storage and Shelf Life ◀

- Storage conditions: **2-35°C**.
- Shelf life: **4 years**.
- Open bottle stability: **6 months**.
- Expiry date: **refer to the product's labeling**.

► Waste Management ◀

Dispose waste product, unused product, and contaminated packaging in compliance with local regulations.

► Materials Required but not Provided ◄

- Hemolyzer® 3 NG hematology analyzer.
- Reagent pickup tubes and caps.
- Empty waste container.
- Control and calibrator materials.

► Reagent Installation Precautions ◄

- Person installing the reagent must be trained laboratory professional.
- Always keep the reagent in its original container.
- Do not mix the remains of a reagent with another container.
- Use the reagent pickup tubes and caps supplied with the instrument or remove them from the old reagent.
- Clean the pickup tubes by wiping them from top to bottom, using lint-free tissue dampened with distilled water. Avoid any dust or microbial contamination of the tubing and reagent.
- Do NOT place the product above the instrument or on the ground! Place it on the same level as the analyzer.

► Connecting the Reagent ◄

1. Open the caps of the bottles inside the reagent pack. Cut the sealing foil using a clean, sharp object and carefully remove it from the bottles.
2. Immerse the pickup tube into the reagent and tighten the cap. Make sure you match the color on the tube, the cap, and the back of the analyzer with the correct reagent:
 - Diluent – **GREEN**
 - Lyse – **YELLOW**
 - Cleaner – **WHITE**

Make sure the pickup tubes are not bent, broken, twisted, or blocked.

3. Connect an empty waste container with RED cap. Waste container should be placed on the ground.
4. Turn on the analyzer and wait for the startup process to finish. Tap and hold any point on the screen to access the local menu.
5. Tap the "Read QR" icon. Hold the QR code approximately 10 - 15 centimeters from the QR/barcode reader or camera. Ensure the QR sheet is not curved or bent. In case the sheet is damaged, use manual data entry.
6. When the QR code data is entered, the product is ready to use.
7. Run a Fast Blank measurement.

► Reagent Replacement ◄

1. Turn on the analyzer and wait for the startup process to finish. Tap and hold any point on the screen to access the local menu.
2. Tap the "Reag. Replace" icon. Follow the steps on the screen, tap "Next" when a step is complete.
3. The sequence will end with an automated priming cycle and a Fast Blank measurement. When the Fast Blank is in range, the analyzer will display the result and you can use the new reagent.
4. When a new reagent batch is used, verify performance with control measurements and perform calibration if needed.

► Limitations ◄

- The reagent pack is designed to run the number of tests mentioned in the "Product Information" section under the following conditions:
 - A minimum # of tests are measured daily:
 - **H3-Compact 100: 3-5 samples per day.**
 - **H3-Compact 500: 15 samples per day.**
 - Two Fast Blanks are performed daily.
 - The analyzer is powered off every day.
 - The necessary maintenance actions prompted by the software are performed.
- After the specified number of tests, no new samples can be measured.
- The amount of leftover reagent inside the pack may vary depending on usage.
- Quality control, calibration and patient samples equally reduce the number of tests available.
- Running lower number of daily tests or more Fast Blank measurements than the conditions listed above will result in relatively more reagent consumption and the pack can run out of reagents before the specified number of tests.
- No complaints are accepted regarding leftover reagents or about the number of tests, if the above listed conditions are not met.

► Explanation of Symbols ◄

	Catalogue number		Batch number
	Expiry date		In vitro diagnostic medical device
	Permitted storage temperature		CE mark
	See the instructions for use		Manufacturer
	Content		Date of manufacture
	Distributor		

H3-Check

Hemolyzer® 3 NG Control Set

Order information:

Catalog-No.	Article name	Content
HE3504	H3-Check Complete	6 x 2.5 mL (L, N, H)
HE3505	H3-Check Normal	6 x 2.5 mL (N)

Intended Use:

H3-Check Control Sets are intended for quality control measurements on the Analyticon hematology analyzer Hemolyzer® 3 NG (open and closed mode) (HE3100), and serve to validate the measurement results.

Reagents:

The H3-Check Control Set contains human erythrocytes, simulated leukocytes, and mammalian platelets suspended in a plasma-like fluid with preservatives. For lot-specific target values, please refer to the table below.

Warning and Precautions:

For in vitro diagnostic use only.

For trained laboratory professional only.

All body fluid samples should be considered potentially infectious materials. The product was tested negative against anti-HIV 1 and 2, anti-HCV and HBsAg, but since no test method can rule out the danger of infection with absolute certainty, treat all blood and other potentially infectious materials with appropriate precautions. Use gloves, masks and gowns if blood exposure is anticipated.

Exercise the normal precautions required for handling all laboratory reagents. Disposal of the product should be in accordance with local regulations. The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Storage and Stability:

Unopened tube	Until the expiration date	2 – 8 °C
Opened tube	14 days	2 – 8 °C

Store the tubes in upright position. Do not freeze!

Indications of Deterioration:

After mixing, product should be similar in appearance to fresh whole blood. In unmixed vials, the supernatant may appear cloudy and reddish, which is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. Do not use the product if deterioration is suspected.

Procedure:

1. Allow the tubes to warm to room temperature (15 to 30 °C) for at least 15 minutes.
2. Prior to the first use of the control, make sure to open the cap, and gently tighten it back ensuring no leaks. Do not overtighten the cap, otherwise you may experience needle clogs.
3. Invert the tubes for at least 20–30 seconds to mix the control material thoroughly. The red precipitate at the bottom of the tube has to be completely suspended and the red color of the fluid needs to be evenly distributed. Tubes stored for a long time might require longer mixing.
4. Gently invert the tube 8–10 times immediately before sampling and analyze the sample as instructed in the Quality Control section of the Operator's Manual for your instrument.
5. If the tube has been open for sampling, clean residual material from the cap and tube rim with a lint-free tissue after the sampling. Close the cap tightly.
6. Return the tubes to refrigerator (2–8 °C) within 30 minutes of use.

Results and Expected Values:

Verify that the lot number on the tube matches the lot number on the table of target values.

Target values are determined on well-maintained, properly calibrated instruments using the instrument manufacturer's recommended reagents. The "Expected Range" takes inherent imprecision of the method and expected biological variability of the control material into account.

Minor variations may derive from different reagent batches, maintenance status, operating technique, and calibration and might contribute to inter-laboratory variation.

For greater control sensitivity, it is recommended that each laboratory establishes its own target value and expected range, and periodically re-evaluates the target value. The expected range of the laboratory may include values outside of the expected range given by the manufacturer.

Limitations of the Procedure:

The performance of this product is assured only if it is properly stored and handled as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

Control vials can be pierced maximum 40 times. After that the rubber part of the cap has the risk of dropping into the control, causing clogs in the instrument.

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Hypochlorite Solution

► Product Information ◀

Product name	Product code
H3/5-Hypoclean	HE3413

► Product Content ◀

	HE3413
H3/5-Hypoclean	5x100 mL
Instructions for Use	1 pc

► Intended Use and Test Principle ◀

H3/5-Hypoclean detergent is a stabilized and micro-filtered hypochlorite solution for intensive oxidative rinsing and cleaning of capillaries, tubing, and chambers of Hemolyzer® semi-automated hematology analyzers.

The detergent is intended for laboratory professional use as an in vitro diagnostic medical device.

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► Warnings and Precautions ◀

- Read instructions for use carefully before use.
- Do not use the product beyond the expiry date.
- Do not use the product if the primary package is damaged or there is visible sign of chemical deterioration of the liquid inside. In both cases the product should be replaced.
- Always keep the detergent in its original container.
- Store at the recommended temperature range.
- Do not use the product if it has been frozen or kept at excessive heat. Use the product at room temperature (18-30°C). If it has been stored below 18°C or above 30°C, leave it at normal room temperature for 3 hours before use.
- Do not mix the remains of a detergent with another container.
- Do NOT place the product above the instrument or on the ground! Place it on the same level as the analyzer.
- Users are advised to wear approved protective clothing when handling chemical products: lab coat, gloves, and eye protection.
- Observe standard laboratory precautions for use and follow national or local health and safety guidelines.
- Avoid contact with eyes, skin, and clothing. In case of contact with eyes or skin, rinse immediately with plenty of water.
- Please refer to the product's Safety Data Sheet.
- Please report any serious incident involving the product to the manufacturer and the competent local authority.

► Storage and Shelf Life ◀

- Storage conditions: **2-35°C**.
- Shelf life: **2 years**.
- Open bottle stability: **6 months**.
- Expiry date: **refer to the product's labeling**.

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► How to Use: Hemolyzer® 3 NG ◀

Materials required but not provided

- Hemolyzer® 3 NG hematology analyzer.
- Empty cleaning tube.

Cleaning procedure

The detergent is used to perform weekly maintenance procedures on the Hemolyzer® 3 NG hematology analyzer, such as closed sampling cleaning, open sampling cleaning, rinsing and RBC cleaning. When cleaning is required, the calendar icon will appear on the screen.

To perform these regular cleaning procedures, tap on »Management« in the main menu, then go to the »Cleaning« submenu. Please follow the instructions on the analyzer. The software guides you through the cleaning processes.

► How to Use: Hemolyzer® 5 NG ◀

Materials required but not provided

- Hemolyzer® 5 NG hematology analyzer.
- Reagent pickup tubes and caps.

Detergent installation precautions

- Person installing the detergent must be trained laboratory professional.
- Use the reagent pickup tubes and caps supplied with the instrument.
- Clean the pickup tubes by wiping them from top to bottom, using lint-free tissue dampened with distilled water. Avoid any dust or microbial contamination of the tubing and detergent.

Cleaning procedure

The detergent is used to perform RBC Hard Clean cleaning procedure on the Hemolyzer® 5 NG hematology analyzers. Performing the cleaning is recommended if the analyzer reports 5 consecutive samples with Contamination flag.

To perform RBC Hard Clean, tap on »Management« in the main menu, then go to the »Cleaning« submenu. There are two different cleaning options to choose from depending on which reagent pickup tube set is in use:

- Cleaning with old tubing set – without Quick Lock
- Cleaning with new tubing set – with Quick Lock

Please follow the instructions on the analyzer. The software guides you through the cleaning processes.

► Explanation of Symbols ◀

	Catalogue number		Batch number
	Expiry date		In vitro diagnostic medical device
	Permitted storage temperature		CE mark
	See the instructions for use		Manufacturer
	Corrosion – causes severe skin burns and eye damage.		Content
	Date of manufacture		Distributor