FTCE-CALI2-5-v5 (07/2020)

Français - FR

USAGE PRÉVU

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems quantitative tests listed in the value sheet.
This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Sodium azide < 0.1 % (w/w)
- Concentrations for each analyte to test are lot-specific.

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TRACEABILITY

The traceability is indicated in the value sheet enclosed in the kit.

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved or cleared in compliance with European Directive 98/79/EC, Annex II, List A. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

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PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.
- Add in exactly 3 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particulates would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Stored at 2-8 °C and protected from light.
- Do not use after expiration date indicated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.

- Stability of the components:

Between 15 and 25 °C : 8 hours
Between 2 and 8 °C : 2 days
Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light):
Between 15 and 25 °C : 6 hours
Between 2 and 8 °C : 1 day
Between -25 and -15 °C : 2 weeks (When frozen once)
- Stability of direct bilirubin (when stored protected from light):
Between 15 and 25 °C : 3 hours
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 2 weeks (when frozen once)

Note: Store vials tightly capped and protected from light.

PROCEDURE

To use ELICAL 2, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using ELICAL 2 has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet).

Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.
For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

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English - EN

INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of ELITech Clinical Systems quantitative tests listed in the value sheet. These *in vitro* diagnostic devices are for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Concentrations for each analyte to test are lot-specific.

MATERIALS REQUIRED BUT NOT PROVIDED

- ELITech Clinical Systems reagents listed in the value sheet.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated on the instructions for use of the ELITech Clinical Systems reagent used in combination with the recommended calibrator.

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.
- Carefully open the vial, avoiding the loss of lyophilizate.
- Add exactly 5 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

Important: To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).

- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.

- Do not use after the expiry date stated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.
- Stability of the components:
Between 2 and 8 °C : 5 days
Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light) :
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of direct bilirubin (when stored protected from light) :
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 5 days (when frozen once)

Stability of alkaline phosphatase:

- Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 4 weeks (when frozen once)

Note:

- Store the vials tightly capped and protected from light after reconstitution.
- Mix content thoroughly prior to use.
- Thawing frozen specimen at 37 °C, and vigorous mixing (vortex), should be avoided as this may denature ferritin.

PROCEDURE

To use ELITROL I & II, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using ELITROL I & II has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet). Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit. Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values. For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

ASSISTANCE TECHNIQUE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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☛ BIBLIOGRAFIE

1. Johnson, A.M., *Amino Acids and Proteins*, Tietz Fundamentals of Clinical Chemistry, 6th Ed, Burtis, C.A. & Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds), (2008), 286.
2. Sanhai, W.R., et al., *Cardiac and Muscle Disease*, Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc. eds.), (2010), 677 and appendix.
3. Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 190.
4. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
5. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
7. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the quality control of the performances of ELITech Clinical Systems reagents (Ref. ICRP-0400), used for C-reactive protein (CRP) quantitative determination. This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Liquid product prepared from human plasma spiked with human recombinant CRP.
- Sodium azide < 0.1 % (w/w)
- C-reactive protein (CRP) concentrations are lot-specific.

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- ELITech Clinical Systems reagent, used for C-reactive protein (CRP) quantitative determination (Ref. ICRP-0400).
- General Laboratory equipment.

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PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

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These products are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

- Store at 2-8 °C and protect from light. Do not freeze.
- Do not use after the expiry date stated on the vial label.
- These products are stable 6 weeks at 2-8°C after the first opening.
- These products should be immediately and tightly capped to prevent contamination and evaporation.

PROCEDURE

To use CRP IP CONTROL II, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using CRP IP CONTROL II has been validated with the ELITech systems (Analyzers and reagents used for C-reactive protein (CRP) quantitative determination (Ref. ICRP-0400)). Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit. Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European Union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems reagents, used for C-reactive protein (CRP) quantitative determination. This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Liquid product prepared from human plasma spiked with Human recombinant CRP.
- Sodium azide < 0.1 % (w/w)
- C-reactive protein (CRP) concentrations are lot-specific.

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- ELITech Clinical Systems reagent, used for C-reactive protein (CRP) quantitative determination (Ref. ICRP-XXXX).
- General Laboratory equipment.

TRACEABILITY

Calibrators from CRP IP CALIBRATOR SET are traceable to ERM-DA474/IFCC reference material.

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

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PROCEDURE

To use CRP IP CALIBRATOR SET, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using CRP IP CALIBRATOR SET has been validated with the ELITech systems (Analyzers and reagents used for C-reactive protein (CRP) quantitative determination (Ref. ICRP-XXXX)). Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.

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Coagulation Control Plasmas

REF 5186 Routine Control N
REF 5187 Routine Control A
REF 5183 Routine Control SA
REF 5482 Routine Coagulation Control Set

Helena
Biosciences Europe

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-0482P 2016/01 (16)

Coagulation Control Plasmas

en

INTENDED PURPOSE

The Coagulation Control Plasmas kit is intended for use as a quality control material.
Routine Control N, Routine Control A and Routine Control SA are for use as normal, moderately prolonged and markedly prolonged controls for PT and aPTT assays. They are also assigned for Fibrinogen, TCT and ATIII, and are prepared from normal human plasma.

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Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:
Hepatitis B Antigen (HbsAg)
Hepatitis C Antigen (HcAg)
HCV antibody HCV Z antibody
However they should be handled with the same precautions as a human patient sample.

COMPOSITION

REF	Component	Content	Description
5186	Routine Control - N	10 x 1 mL	Prepared from pooled normal plasma.
5187	Routine Control - A	10 x 1 mL	Prepared from atypical human plasma.
5183	Routine Control - SA	10 x 1 mL	Prepared from atypical human plasma.
5482	Routine Coagulation Control Set:		
	Routine Control - N	4 x 1 mL	
	Routine Control - A	3 x 1 mL	
	Routine Control - SA	3 x 1 mL	

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Each kit contains 1 mL of buffered, lyophilised human plasma.

Preparation: Reconstitute each vial of the appropriate control with 1 mL of distilled or deionised water. Swirl gently. Allow to stand for 10 minutes for complete dissolution and mix well before use.

ITEMS REQUIRED BUT NOT PROVIDED

Coagulation Control Plasmas may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with suitable, commercial reagents.

STORAGE SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. The reconstituted controls are stable for 8 hours when kept at 2 – 8°C or 4 weeks at -20°C when flash frozen. Keep covered.

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Not applicable.

PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS

Routine Control N should give values within the laboratory normal range for PT, aPTT and fibrinogen assays. Routine Control A and Routine Control SA have been standardised to give prolonged and markedly prolonged PT and aPTT times respectively. Lot and instrument specific expected values are provided with each pack of controls.

LIMITATIONS

The results obtained with Coagulation Control Plasmas depend on several factors strongly associated with instrumentation. Types of reagents, deficient substrates and laboratory to laboratory variations ^{1,2,3}. Each laboratory should establish an expected range for the particular instrument/reagent system.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

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PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility			
Sample	n	Intra-assay precision aPTT CV (%)	PT CV (%)
Routine Control N	5	2.83	1.01
Routine Control A	5	2.76	1.71
Routine Control SA	5	1.72	1.03

BIBLIOGRAPHY

1. Kirkwood TBL *et al* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, 37:555-568.
2. Goldfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP* 55:561-564.
3. Pallitt HA and Longbery JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP* 59:231-235.

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BIBLIOGRAFIA

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, *Acta Haematol*, 17:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci.* 8:145-192.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Scully RE *et al.* (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, 302(37):37-48.
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, 38(6):196-201

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CE

IVD

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-1787P 2016/01 (6)

APTT Si L Minus

Instructions for use

en

INTENDED PURPOSE

The APTT Si L Minus kit is intended for carrying out clot based haemostasis assays.

For use in the determination of activated partial thromboplastin times (aPTT), and related coagulation procedures using phospholipid extract and a near-colloidal particle activator. The test system can be used as manual, semi-automated and automated methods. From its origins through the work of Langdell and coworkers¹, and later modified by Proctor and Rapaport², the aPTT is used to detect disorders in the intrinsic coagulation system, which involves coagulation factors VIII, IX, XI, XII, prekallikrein, and high molecular weight kininogen. The aPTT is also used in assays which quantitate these factors and is routinely used for presurgical screening and monitoring of heparin therapy³. Commercially available reagents typically use one of three activators: kaolin, silica, or ellagic acid. In the basic screening test, the aPTT indirectly measures the formation of thrombin by its action on fibrinogen forming the fibrin clot. In the test, citrated test plasma is mixed with aPTT reagent for a specified period of time (typically 5 minutes) at 37°C followed by the addition of pre-warmed (37°C) calcium chloride (0.025 M). Timing is begun from the time of addition of calcium chloride. The time required for clot formation is the aPTT. Clot detection can be by mechanical, manual (tilt tube), or photo-optical measurement.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Component	Content	Description	Preparation
APTT Si L Minus	5 x 5 mL (REF 5562 / REF 5562SLQ) 5 x 10 mL (REF 5558 / REF 5558SLQ) 10 x 10 mL (REF 5559 / REF 5559SLQ)	Reagent contains a near colloidal particle activator (magnesium-aluminium-silicate) for optimum sensitivity to factor deficiencies and to heparin. The reagent also contains phospholipids with buffer and stabilisers.	Bring to room temperature prior to use. Mix well by swirling or inversion prior to use.
Calcium Chloride: 0.025M	5 x 5 mL (REF 5562 / REF 5562SLQ) 5 x 10 mL (REF 5558 / REF 5558SLQ) 10 x 10 mL (REF 5559 / REF 5559SLQ)	The reagent is a 0.025 M solution of calcium chloride.	The reagent is ready for use as packaged.

Each kit contains instructions for use.

ITEMS REQUIRED BUT NOT PROVIDED

Any high quality electro-mechanical or photo-optical coagulation instrument designed for performing activated partial thromboplastin times may be used.

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Store at 2° – -8°C. DO NOT FREEZE. Stable for 30 days after opening. Avoid prolonged heating.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 2° – -8°C or 18 – 24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes. This will minimise the neutralisation of the lupus inhibitor⁴. Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, *et al.*⁵.

PROCEDURE

Manual Method

1.

Prewarm well mixed APTT Si L Minus and 0.025 M Calcium Chloride to 37°C.

2.

Prewarm 0.1 mL test plasma in duplicate to 37°C for 2 minutes.

3.

Forcibly add 0.1 mL prewarmed APTT Si L Minus to plasma and start timer. Incubate for exactly 5 minutes at 37°C.

4.

Add 0.1 mL prewarmed 0.025 M Calcium Chloride.

5.

Note time for clot formation. Report result as aPTT time (seconds).

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application notes.

INTERPRETATION OF RESULTS

The results of the aPTT test should be reported to the nearest 1/10 of a second. The normal range (usually X ± 2 standard deviations) for each individual laboratory should be established. Results greater than the upper limits of the normal range should be considered abnormal and follow-up testing should be performed. Any aPTT values less than the lower limits of the normal range should be repeated on a new blood sample. Short aPTT values may be seen in association with *in vivo* thrombosis (e.g. deep vein thrombosis and disseminated intravascular coagulation).

Heparin monitoring

When monitoring heparin therapy, it is important to construct an *in vitro* reference curve which reflects the average heparin response, since individual patients respond differently to heparin. In general, one can consider the therapeutic range for heparin to be 0.2 to 0.5 units/mL^{6,8}. The following precautions should be considered when monitoring heparin therapy:

1.

Time of collection is important, since heparin has an *in vivo* half-life of only 1.5 hours⁷.

2.

Release of platelet factor 4 (heparin neutralising factor) caused by platelet aggregation or damage during collection, should be avoided. Careful blood collection, proper centrifugation and prompt removal of the platelet poor plasma from the cells will help minimise the release of platelet factor 4.

3.

Baseline data on each patient's aPTT should be established before therapy, to determine the respective patient aPTT as it relates to the normal range established by the laboratory.

4.

Different clot detection systems (mechanical, photo-optical, etc.) show variable sensitivities to heparin. The same test system should be used when monitoring heparinised patients.

5.

Heparin response curves should be reestablished when lot numbers of reagent change and at periodic intervals with the same lot number.

6.

The curve should also be constructed using the same heparin employed in therapy, to eliminate variables connected with heparins from different sources (e.g. porcine mucosa or bovine lung).

LIMITATIONS

Expected values for the aPTT test will vary from one laboratory to another, depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of specimen storage are all very important. Plasma sample collection and storage conditions should be standardised and carefully controlled. Unexpected results should be confirmed by additional tests. Platelet fragments present in a specimen may cause the release of phospholipids, and thus the neutralisation of any lupus inhibitor present in the specimen. The use of specimens with small plasma volumes should be avoided due to possible physiological pH changes. Testing could be affected by several drugs⁸. An increase in the aPTT results may be caused by the administration of diphenylhydantoin, heparin, warfarin and radiographic agents^{8,9}. Decreased aPTT values may be seen during the use of oral contraceptives, or male estrogen therapy^{9,10}. Thus, laboratories should establish their own expected values for patients and well defined performance standards for the control.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA
REF 5189	Heparin Control
REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using a photo-optical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility

	Intra-assay precision			Inter-assay precision		
Sample	n	Clot formation (seconds)	CV (%)	n	Clot formation (seconds)	CV (%)
Routine Control N	10	33.0	0.36	100	32.9	2.41
Routine Control SA	10	77.9	0.31	100	78.3	0.77

% Factor

	Factor VIII (seconds)	Factor IX (seconds)	Factor XI (seconds)
<1	85.7	70.9	92.4
10	45.9	44.4	50.7
40	34.2	33.9	34.4
100	28.8	28.8	28.8

BIBLIOGRAPHY

1.

Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, *J. Lab. Clin. Med.* **41**: 637

2.

Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, *Am. J. Clin. Path.* **36**: 212

3.

Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, *AJCP*, **76**: 530-537

4.

Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5

5.

Young DS *et al.* Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACP Press, Washington, D.C., 1990

6.

Kellon JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, *Semin Hematol*, **17**: 259-91

7.

Estes J and Poulin P (1974) Pharmacokinetics of Heparin, *Thromb Diath Haemorr*, **33**: 26-35

8.

Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, *Neurology*, **22**: 1165-71

9.

Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and coagulation factors during open heart surgery, *J Med*, **2**: 65-81

10.

Crowell EB Jr, Clatanoff DV, Kiekhofor W (1971) The effect of oral contraceptives on factor VIII levels, *Lab. Clin. Med.* **77**: 551-7

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BIBLIOGRAFIA

- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, *J. Lab. Clin. Med.* **41**: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, *Am. J. Clin. Path.* **36**: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, *AJCP*. **76**: 530-537
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS *et al.* Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACPS Press, Washington, D.C.,1990
- Kellon JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, *Semin Hematol*, **17**: 259-91
- Estes J and Poulin P (1974) Pharmacokinetics of Heparin, *Thromb Diath Haemorr*, **33**: 26-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, *Neurology*, **22**: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassarman HB, Moore RH, Melzer J (1971) Effect of antithrombolytic agents and estrogens on blood loss and coagulation factors during open heart surgery, *J Med*, **2**: 65-81
- Crowell EB Jr, Clatanoff DV, Kiekhofe W (1971) The effect of oral contraceptives on factor VIII levels, *Lab. Clin. Med.* **77**: 551-7

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helenabiosciences Europe

REF 5265HL
REF 5265L
REF 5267L

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helenabiosciences.com
Web: www.helena-biosciences.com

HL-2-3035P 2015/10 (1)

Thromboplastin L

Instructions for use

INTENDED PURPOSE

The Thromboplastin L kit is intended for carrying out clot based haemostasis assays.

The first standardised one-stage prothrombin test was developed by Dr. Armand Quick in 1935. It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors II, V, VII and X)^{1,2}. It is also used for the induction and monitoring of oral anticoagulant therapy³⁻⁶ and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders. Thromboplastin L is of rabbit brain origin but resembles human preparations in its low International Sensitivity Index (ISI). The ISI of Thromboplastin L is approximately 1.1 and is calibrated against the WHO international reference preparation⁷. Thromboplastin L is particularly suited to the monitoring of oral anticoagulant therapy and, in conjunction with the appropriate factor deficient plasma, the measurement of factor activity in the extrinsic pathway. Tissue thromboplastin, in the presence of calcium ions, is an activator which initiates the extrinsic pathway of coagulation. When a mixture of tissue thromboplastin and calcium ions is added to normal citrated plasma, the clotting mechanism is activated, leading to a fibrin clot. If a deficiency exists within the extrinsic pathway, the time required for clot formation will be prolonged depending on the severity of the deficiency.

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COMPOSITION

Composition	Content	Description	Preparation
Thromboplastin L	2 x 5 mL (REF 5265HL) 8 x 5 mL (REF 5265L) 10 x 10 mL (REF 5267L)	Liquid Rabbit Brain Thromboplastin containing Calcium Chloride, stabilisers and preservatives.	The liquid, calcified thromboplastin is ready-for-use. No further calcium is required to carry out standard PT Assays. The contents of the vial should be mixed well before use. (5 minutes on roller).

Each kit contains Instructions For Use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

The below products can be used in conjunction with Thromboplastin L:

REF 5519 ISI Calibrant Plasma Set
REF 5490 INR Reference Set

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thromboplastin L Opened vials are stable for 2 months at *2 – *8°C, 5 days at *15°C (on-board Sysmex CA-1500) and 6 hours at *37°C (on-board AC-4 including reagent container and cap). A shift-use stability of 7 days (Sysmex CA-1500) can be achieved. DO NOT FREEZE. Large clumps of particles or changes in expected values may indicate product deterioration.

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Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *18 – *24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes⁸.

PROCEDURE

For accurate INR reporting, it is recommended to determine the laboratory specific ISI of the reagent with the testing system in use. The Helena Biosciences Europe ISI Calibrant Plasma Set (REF 5519) is recommended for this purpose^{1,2}. This should be performed for each new reagent batch. The Helena Biosciences Europe INR Reference Set (REF 5490) should be used to check for shifts in the local system ISI which have been noted with changes in laboratory temperature and post instrument servicing, amongst other local variables.

Manual Method

1. Mix sufficient Thromboplastin L to complete the anticipated testing for the day and incubate at *37°C for no more than 4 hours.

2. Pre-warm 0.1 mL of the test plasma at *37°C for 2 minutes.

3. Add 0.2 mL of freshly mixed thromboplastin reagent to the plasma while simultaneously starting a stopwatch.

4. Note the time for clot formation to the nearest 0.1 seconds.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Results should be reported to the nearest 0.1 seconds and duplicates should agree within 5% of each other. %PT values can be interpolated from the calibration graph (%PT of PT Calibration plasmas versus measured clot time), which should be a straight line when plotted on log-log graph paper.

INR values can be calculated using the following formula: INR = (PT Time Patient / Mean Normal PT Time)^{ISI}

For clear guidance on the indications for and management of patients on warfarin, please refer to The British Society for Haematology, for their most current edition of 'Guidelines on oral anticoagulation with warfarin'. At time of printing this is the 2011 fourth edition⁹.

LIMITATIONS

The use of serial dilutions of a reference plasma for the %PT curve is not recommended as this can lead to discrepancies caused by the low fibrinogen in the reference plasma dilutions which are not reflected in patient samples having predominantly normal fibrinogen levels. Helena Biosciences Europe advise use of the 5504R %PT/Direct INR kit for this purpose.

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REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. This is particularly important for local ISI calibration. Using the Sysmex series of instruments, normal values ranging from 11.50 - 14.60 seconds; 0.930 - 1.160 INR; 79.10 - 112.80 %PT are typical.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using a Sysmex CA-1500 coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility						
Sample	SD	CV (%)	SD	CV (%)	SD	CV (%)
Repeatability	0.07	0.59	0.24	1.09	0.45	1.11
Between-run	0.10	0.83	0.16	0.75	0.49	1.20
Between-day	0.04	0.32	0.06	0.27	0.25	0.62
Within-device / Laboratory	0.12	1.07	0.29	1.35	0.72	1.75

Interferences

Helena Thromboplastin L is insensitive to Heparin levels of up to 2 U/mL. Using a 5% interference threshold, there is no significant interference from Haemoglobin at concentrations up to 10 g/L. Using a 5% interference threshold, there is no significant interference from Bilirubin at concentrations up to 0.5 g/L for Thromboplastin L. Lipid interference testing demonstrates that lipid levels do not directly affect the clot time of the reagent up to 3.75g/L. Lipid concentrations in excess of this prevent clot detection.

Method Comparison

Comparison of clot time in seconds and INR values were determined using Thromboplastin L and Thromboplastin LI on 268 samples. The following correlations were obtained:

Thromboplastin L (Seconds) = 0.9911x + 0.1038 ^r² = 0.9941 n = 268
Thromboplastin L (INR) = 0.9853x + 0.0261 ^r² = 0.9500 n = 268

BIBLIOGRAPHY

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci*, **190**: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 5S-11S.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, **700**: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol*, **103**: 358-365.
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol*, **103**: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, **154**(3): 311-324.

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LITERATURVERZEICHNIS

1. Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci*, **190**: 501.

2. Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.

3. Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 5S-11S.

4. Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.

5. World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, **700**: 19.

6. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.

7. Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol*, **103**: 358-365.

8. Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol*, **103**: 366-371.

9. Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, **154**(3): 311-324.

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BIBLIOGRAFIA

- Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, *Am. J. Med. Sci*, **190**: 501.
- Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
- Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, *Chest*, **95**: 5S-11S.
- Poller L (1986) Laboratory Control of Anticoagulant Therapy, *Sem. Thromb. Haemostasis*, **12**: 13-19.
- World Health Organisation (1984) Expert Committee on Biological Standards, *Technical Series*, **700**: 19.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Poller L., Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, *Amer. J. Clin. Pathol*, **103**: 358-365.
- Poller L., Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, *Amer. J. Clin. Pathol*, **103**: 366-371.
- Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, *British Journal of Haematology*, **154**(3): 311-324.

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helena
Biosciences Europe

HL-7-0135DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5183	Routine Control SA	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0137DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.


Product Code	Description	GMDN Classification Code
5186	Routine Control N	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0138DC DOI 2015/07 (7)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.


Product Code	Description	GMDN Classification Code
5187	Routine Control A	30590

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 28 Jul 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0287DC DOI 2015/08 (1)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.

Product Code	Description	GMDN Classification Code
5559	APTT Si L Minus	55981

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 11 Aug 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0511DC DOI 2015/08 (4)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.


Product Code	Description	GMDN Classification Code
5376	Clauss Fibrinogen 100	55997

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 12 Aug 2015

Tel +44 (0)191 482 8440

Fax +44 (0)191 482 8442

info@helena-biosciences.com

www.helena-biosciences.com

Helena Biosciences Europe

Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,

United Kingdom

Declaration of Conformity

helena
Biosciences Europe

HL-7-0664DC DOI 2015/08 (1)

In Application of the Council Directive 98/79/EC on the approximation of the laws of the Member States relating to *In Vitro* Diagnostic Medical Devices & CE marking.

Declaration of conformance to applicable sections of Annex I - Essential Requirements and Annex III (EC Declaration of Conformity) imposed by sections 2 to 5. The below listed products are not classified under Annex II Lists A or B. Access to the appropriate technical files will be made available to the appropriate body in the event this is required.


Product Code	Description	GMDN Classification Code
5267L	Thromboplastin L	55983

I, the undersigned declare that the devices registered against the above GMDN Classification Code conforms to the said Directives.

Full Name: M.J. Stephenson

Title: Managing Director

Signed:

Date: 06 Aug 2015

Tel +44 (0)191 482 8440
Fax +44 (0)191 482 8442
info@helena-biosciences.com
www.helena-biosciences.com

Helena Biosciences Europe
Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD,
United Kingdom

REF HPC-402	English
-------------	---------

A rapid test for the qualitative detection of Procalcitonin in human's whole blood, serum or plasma.
For professional in vitro diagnostic use only.

【INTENDED USE】

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of Procalcitonin in whole blood, serum or plasma.

【SUMMARY】

Procalcitonin (PCT) is a small protein that comprises 116 amino acid residues with a molecular weight of approximately 13 kDa which was first described by Moulic et al. in 1984.^{1,2} PCT is produced normally in C-cells of the thyroid glands. In 1993, the elevated level of PCT in patients with a system infection of bacterial origin was reported and PCT is now considered to be the main marker of disorders accompanied by systemic inflammation and sepsis.^{3,4} The diagnostic value of PCT is important due to the close correlation between PCT concentration and the severity of inflammation. It was shown that "inflammatory" PCT is not produced in C-cells. Cells of neuroendocrine origin are presumably the source of PCT during inflammation.⁵

【PRINCIPLE】

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) is a qualitative, lateral flow immunoassay for the detection of PCT in whole blood, serum or plasma. The membrane is pre-coated with anti-PCT antibody on the test line region of the strip. During testing, the whole blood, serum or plasma specimen reacts with the particle coated with anti-PCT antibody. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-PCT antibody on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

【REAGENTS】

The test device contains mouse anti-PCT antibody particles and mouse anti-PCT antibody coated on the membrane.

【PRECAUTIONS】

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- Do not eat, drink or smoke in the area where the specimen or kits are handled.
- Handle all the specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Humidity and temperature can adversely affect results.

【STORAGE AND STABILITY】

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

【SPECIMEN COLLECTION AND PREPARATION】

- The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.
- To collect **Fingerstick Whole Blood Specimens**:
 - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
 - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
 - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
 - Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube**:
 - Touch the end of the capillary tube to the blood until filled to approximately 50 µL. Avoid air bubbles.
 - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen well (S) of the Test Cassette.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days, for long-term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents
- EDTA K2, Heparin sodium, Citrate sodium and Potassium Oxalate can be used as the anticoagulant for collecting the specimen.

【MATERIALS】

- | Materials Provided | | | |
|--|---|----------|------------------|
| • Test Cassettes | • Sample Droppers | • Buffer | • Package Insert |
| Materials Required But Not Provided | | | |
| • Specimen collection containers | • Centrifuge (for plasma only) | | |
| • Sterile Lancet(forFingerstick Whole Bloodonly) | • Capillary tube(forFingerstick Whole Blood only) | | |
| • Timer | | | |

【DIRECTIONS FOR USE】

Allow the tests, specimen, buffer and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Remove the test cassette from sealed pouch and used it within one hour. Best results will be obtained if the assay is performed immediately after opening foil pouch.
- Place the test cassette on a clean and level surface.

For **Serum or Plasma** specimens:

Hold the dropper vertically and **transfer 1 drop of serum or plasma** (approximately 25 µL) to the specimen well (S) of the Test Cassette, and **add 1 drop of buffer** (approximately 40 µL), then start the timer. See illustration below.

For **Venipuncture Whole Blood** specimens:


Hold the dropper vertically and **transfer 2 drops of whole blood** (approximately 50 µL) to the specimen well (S) of the Test Cassette, and **add 1 drop of buffer** (approximately 40 µL), then start the timer. See illustration below.

For **Fingerstick Whole Blood** specimens:

To use a capillary tube: Fill the capillary tube and **transfer approximately 50 µL of fingerstick whole blood specimen** to the specimen well (S) of the Test Cassette, then **add 1 drop of buffer** (approximately 40 µL) and start the timer. See illustration below.

- Wait for the colored line(s) to appear. The result should be read at **15 minutes**. Do not interpret the result after **20 minutes**.

Note: It is suggested not to use the buffer, beyond 30 days after opening the vial.

【INTERPRETATION OF RESULTS】

(Please refer to the illustration above)

POSITIVE: * Two colored lines appear. One colored line should be in the control line region (C) and another colored line should be in the test line region (T).

***NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of PCT antigen present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive.

NEGATIVE: One colored line appears in the control line region (C). No apparent colored line appears in the test line region (T).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

An internal procedural control is included in the test. A red line appearing in the control line region (C) is an internal quality procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

【LIMITATIONS】

- The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) is for in vitro diagnostic use only. This test should be used for the detection of PCT in whole blood, serum or plasma specimen.
- The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) cannot detect less than 1ng/mL of PCT in specimens.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- In some instances elevated Procalcitonin levels in due to noninfectious reasons can be observed:
 - During the first days after trauma or surgical intervention burns, release of proinflammatory cytokines, lung cancer (oat cell carcinoma), Medullary Thyroid Carcinoma (C-Cell Carcinoma).
 - New born children, < 48hours.
 - Severe cardiogenic shock.
- The hematocrit of the whole blood should be between 25% and 65%.

【EXPECTED VALUES】

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) has been compared with a leading commercial PCT ELISA test. The correlation between these two systems is 98.8%.

【PERFORMANCE CHARACTERISTICS】

Sensitivity

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) has correctly identified a panel of specimens and has been compared to a leading commercial PCT ELISA test using clinical specimens. The results show that the relative sensitivity of the PCT Rapid Test Cassette (Whole Blood /Serum /Plasma) is 98.7%, and the relative specificity is 98.9%.

Method	ELISA		Total Results
	Results	Positive	Negative
	Positive	231	3
PCT Rapid Test Cassette (Whole Blood/Serum/Plasma)	Negative	3	280
		234	283
Total Results		234	283

Relative Sensitivity: 98.7% (95%CI*: 96.3%-99.7%)

*Confidence Intervals

Relative Specificity: 98.9% (95%CI*: 96.9%-99.8%)

Overall Accuracy: 98.8% (95%CI*: 97.8%-99.7%)

Precision

Intra-Assay

Within-run precision has been determined by using 15 replicates of three specimens containing negative, low positive and high positive. The negative and positive values were correctly identified 99% of the time.

Inter-Assay

Between-run precision has been determined by using the same three specimens of negative, low positive and high positive of PCT in 15 independent assays. Three different lots of the PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested over a 10-days period using negative, low positive and high positive specimens. The specimens were correctly identified 99% of the time.

Cross-Reactivity

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested by anti-HAMA IgG, anti-RF IgG, anti-HAV IgM, anti-Syphilis IgG, anti-HIV IgG, anti-*H. Pylori* IgG, anti-MONO IgM, anti-Rubella IgG, anti-Rubella IgM, anti-CMV IgG, anti-CMV IgM, anti-Toxo IgG and anti-Toxo IgM positive specimens. The results showed no cross-reactivity.

Interfering Substances

The PCT Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested for possible interference from visibly hemolyzed and lipemic specimens. No interference was observed. In addition, no interference was observed in specimens containing up to 2,000 mg/dL Hemoglobin, 1000 mg/dL Bilirubin, and 2000 mg/dL human serum Albumin.

【BIBLIOGRAPHY】

- Le Moulic JM, et al. (1984) The complete sequence of human procalcitonin. FEBS Letters 167(1), 93-97.
- Assicot M, et al. (1993) High serum procalcitonin concentrations in patients with sepsis and infection. Lancet 341(8844), 515-518.
- Meisner M and Reinhart K (2001) Is procalcitonin really a marker of sepsis? Int J Intensive Care 8(1), 15-25.
- Sponholz C, et al. (2006) Diagnostic value and prognostic implications of serum procalcitonin after cardiac surgery: a systematic review of the literature. Critical Care 10, R145.
- Meisner M, (2002) Pathobiochemistry and clinical use of procalcitonin. Clin Chim Acta 323, 17-29.

Index of Symbols

	Attention, see instruction for use		Test per kit		Do not use if package is damaged
	For in vitro diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot number		Catalog #

Inzek International Trading
Vissenstraat 32
7324 AL Apeldoorn
The Netherlands
info@biozek.com
www.biozek.com

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