# Anti-TG Ab

For use on IMMULITE® 2000 systems

**SIEMENS**

# IMMULITE® 2000 Anti-TG Ab

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of autoantibodies to thyroglobulin (TG) in serum, EDTA, and heparinized plasma, as an aid in the clinical diagnosis of thyroid diseases.

Catalog Number: **L2KTG2** (200 tests)  
**L2KTG6** (600 tests)  
Test Code: **ATG** Color: **Orange**

### Summary and Explanation

Thyroglobulin is produced only by the thyroid gland and is a major component of the thyroid follicular colloid. The thyroid hormones 3,5,3',5',-tetraiodothyronine (thyroxine, T4) and 3,5,3'-triiodothyronine (T3) are synthesized from thyroglobulin.<sup>21</sup>

Autoantibodies to thyroglobulin (TG autoantibodies) are often present in patients with autoimmune thyroid disease. Approximately 10 percent of healthy individuals have TG autoantibodies at low levels; higher concentrations are found in 30 and 85 percent of patients with Graves' disease and Hashimoto's thyroiditis, respectively.<sup>10</sup> Elevated levels of antibodies to thyroid peroxidase (TPO autoantibodies) occur more frequently than high anti-TG levels in these diseases, however; anti-TG determinations therefore do not seem to add to the diagnostic information provided by anti-TPO results.<sup>5</sup>

TG autoantibody measurements are most useful for evaluating samples submitted for thyroglobulin measurements because TG autoantibodies can interfere with both competitive immunoassays and immunometric assays for thyroglobulin.<sup>19</sup>

### Principle of the Procedure

IMMULITE 2000 Anti-TG Ab is a solid-phase, enzyme-labeled, chemiluminescent sequential immunometric assay.

**Incubation Cycles:** 2 × 30 minutes  
**Time to first result:** 65 minutes

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Anti-TG Ab has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 5 µL sample

**Storage:** 7 days at 2–8°C, or 6 months at –20°C.

**Automatic Predilution Factor:** 100

Dispense required volume of diluent (L2AAZ) to a suitable test tube with barcode label applied.

### Warnings and Precautions

For *in vitro* diagnostic use.

**CAUTION! POTENTIAL BIOHAZARD**

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>23-25</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Danger!** Toxic in contact with skin. Harmful if swallowed. Toxic to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection.

Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

H311, H302, H411

P280, P273, P301 + P312, P302 + P312, P501

**Contains:** 2-methyl-2H-isothiazol-3-one, sodium azide; Anti-TG Ab Adjustors

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

**Materials Supplied**

Components are a matched set. Labels on the inside box are needed for the assay.

**Anti-TG Ab Bead Pack (L2TG12)**

With barcode. 200 beads, coated with highly purified thyroglobulin. Stable at 2–8°C until expiration date.

**L2KTG2:** 1 pack **L2KTG6:** 3 packs

**Anti-TG Ab Reagent Wedge (L2TGA2)**

With barcode. 11.5 mL buffer matrix. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

**L2KTG2:** 1 wedge **L2KTG6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

**Anti-TG Ab Adjustors (LTGL, LTGH)**

Two vials (Low and High), of lyophilized TG autoantibodies in a human serum/buffer matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion until the lyophilized material is fully dissolved. (No further dilution is required.) Aliquot and freeze in tightly capped plastic containers. Avoid repeated freeze/thaw cycles. Stable at 2–8°C for 7 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KTG2:** 1 set **L2KTG6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

### Thyroid Autoantibody Sample Diluent (L2AAZ)

Test code: **TAD**. For the on-board dilution of patient samples. 50 mL concentrated (ready-to-use) anti-TG/anti-TPO antibody-free buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

**L2KTG2:** 1 vial **L2KTG6:** 3 vials

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2KTG2:** 3 labels **L2KTG6:** 5 labels

### Kit Components Supplied Separately

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

**LAACM:** Bi-level control module

Also Required

Distilled or deionized water; test tubes; controls.

### Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**  
2 weeks

**Quality Control Samples:** Use controls or serum pools with at least two levels (low and high) of TG and TPO autoantibodies.

### Expected Values

Based on its relationship to IMMULITE Anti-TG Ab (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Serum samples with normal TSH and free T4 levels from 117 apparently healthy adults with no history of thyroid disease were processed using the IMMULITE Anti-TG Ab assay. The results showed a 95<sup>th</sup> percentile of approximately 40 IU/mL, suggesting a *preliminary* reference range for adults of

Nondetectable to 40 IU/mL.

Thirty-four serum samples from Graves' disease patients and 36 serum samples from Hashimoto's disease patients were tested using the IMMULITE Anti-TG Ab assay. Fourteen out of 34 Graves' disease samples (41%) and 25 out of 36 Hashimoto's disease samples (69%) had elevated anti-TG antibody levels.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitations

For diagnostic purposes, anti-TG results should be used in conjunction with other test results, the overall clinical presentation to the physician, and all other appropriate information.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Reportable Range:** 20–3000 IU/mL  
(WHO 1st IRP 65/93)

**Analytical Sensitivity:** 2.2 IU/mL

**Hook Effect:** None up to 55,552 IU/mL

**Precision:** Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See “Precision” table.)

**Linearity:** Samples were assayed under various dilutions. (See “Linearity” table for representative data.) Due to the heterogeneity of TG autoantibodies, certain samples may not maintain dilutional parallelism in the IMMULITE 2000 Anti-TG Ab procedure.

**Recovery:** Samples spiked 1 to 19 with three anti-TG Ab solutions (200, 800, and 2000 IU/mL) were assayed. (See “Recovery” table for representative data.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 4000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 21 volunteers into plain, heparinized, EDTA and Becton Dickinson SST<sup>®</sup> vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of Anti TG Ab, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Anti-TG Ab procedure.

(Heparin) = 1.05 (Serum) – 11 IU/mL  
r = 0.988

(EDTA) = 1.02 (Serum) – 2.6 IU/mL  
r = 0.988

(SST) = 1.07 (Plain Tubes) – 12 IU/mL  
r = 0.993

Means:  
752 IU/mL (Serum)  
779 IU/mL (Heparin)  
767 IU/mL (EDTA)  
795 IU/mL (SST)

**Method Comparison:** The assay was compared to IMMULITE Anti-TG Ab on 178 samples. (Concentration range: approximately 20 to 2500 IU/mL. See graph.) By linear regression:

(IML 2000) = 1.1 (IML) – 21 IU/mL  
r = 0.959

Means:  
297 IU/mL (IMMULITE 2000)  
290 IU/mL (IMMULITE)

## References

- 1) Becker W, Reimers C, Börner W. Immunological criteria for the differential diagnosis of thyroid disease. *Therapiewoche* 1985;35:1167-76.
- 2) Czarnocka B, Ruf J, Ferrand M, et al. Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid diseases. *FEBS Lett* 1985;190:147-52.
- 3) Czarnocka B, Ruf J, Ferrand M, et al. Interaction of highly purified thyroid peroxidase with anti-microsomal antibodies in autoimmune thyroid diseases. *J Endocrinol Invest* 1986;9:135-8.
- 4) Falk U, et al. Radioimmunoassay for serum thyroglobulin designed for early detection of metastases and recurrences in the follow-up of patients with differentiated thyroid carcinoma. *J Clin Chem Clin Biochem* 1984;22:661-70.
- 5) Feldt-Rasmussen U. Serum thyroglobulin and thyroglobulin autoantibodies in thyroid diseases. *Allergy* 1983;38:369-87.
- 6) Gons MH, et al. Concentration of plasma thyroglobulin and urinary excretion of iodinated material in the diagnosis of thyroid disorders in congenital hypothyroidism. *Acta Endocrinol* 1983;104:27-34.
- 7) Hashimoto H. Understanding the lymphatic changes of the thyroid. *Arch Klin Chir* 1912;97:219-48.
- 8) Horster FA. The importance of MAK, TAK, TRAK and thyroglobulin in the diagnosis of thyroid disease. *Internist* 1988;29:538-40.
- 9) Ket JL, et al. Serum thyroglobulin levels: physiological decrease in infancy and the absence in athyroidism. *J Clin Endocrinol Metab* 1981;58:1301-3.
- 10) Larsen PR, Ingbar SH. The thyroid gland. In: Wilson JD, Foster DW, editors. *Williams textbook of endocrinology*. 8th ed. Philadelphia: WB Saunders, 1992: 357-487.
- 11) Libert F, Ruel J, Ludgate M, et al. Thyroperoxidase, an auto-antigen with a mosaic structure made of nuclear and mitochondrial gene modules. *EMBO J* 1987;6:4193-6.
- 12) McKenzie JM, Zakarija M, Sato A. Humoral immunity in Graves' disease. *Clin Endocrinol Metab* 1978;7:1-45.
- 13) Ruf J, Czarnocka B, Ferrand M, et al. Thyroid peroxidase is the organ-specific 'microsomal' antigen involved in thyroid autoimmunity. *Acta Endocrinol (Copenh)* 1987; Suppl. 281:49-55.
- 14) Ruf J, Czarnocka B, Ferrand M, et al. Relationship between immunological structure and biochemical

properties of human thyroid peroxidase. *Endocrinol* 1988;125(3):1211-8. 15) Schatz H, Löbig H. The diagnostic and prognostic importance of antibodies against thyroid microsomes and thyroglobulin. *Akt Endokr Stoffw* 1989;10:146-53. 16) Scherbaum WA, Berg PA. The importance of autoantibodies in the diagnosis of endocrine disease. *Dtsch Med Wschr* 1981;106:308-13. 17) Scherbaum WA, Berg PA. Thyroid autoimmune disease – new aspects in development and diagnosis. *Dtsch Med Wschr* 1984;109:1574-81. 18) Spencer CA, Wang C-C. Thyroglobulin measurement: techniques, clinical benefits, and pitfalls. *Endocrinol Metab Clin North Am* 1995;24:841-63. 19) Spencer CA, Takeuchi M, Kazarosyan M. Current status and performance goals for serum thyroglobulin assays. *Clin Chem* 1996;42:164-73. 20) Spencer CA. Recoveries cannot be used to authenticate thyroglobulin (Tg) measurements when sera contain Tg autoantibodies. *Clin Chem* 1996;42:661-3. 21) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994. 22) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 23) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8. 24) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 25) Federal Occupational Safety and Health Administration, *Bloodborne Pathogens Standard*, 29 CFR 1910.1030.

## Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

## Tables and Graphs

### Intraassay Precision (IU/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	40	3.2	8.0%
2	210	11.2	5.3%
3	950	64.5	6.8%
4	2390	199.6	8.4%

### Interassay Precision (IU/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	40	4.8	11.9%
2	210	14.0	6.7%
3	950	89.6	9.4%
4	2390	241.0	10.1%

### Linearity (IU/mL)


	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	130	—	—
	4 in 8	69	65	106%
	2 in 8	33	33	100%
	1 in 8	17	16	106%
2	16 in 16	319	—	—
	8 in 16	166	160	104%
	4 in 16	82	80	103%
	2 in 16	41	40	103%
3	16 in 16	20	20	100%
	8 in 16	172	170	101%
	4 in 16	94	85	111%
	2 in 16	48	43	112%
4	16 in 16	20	21	95%
	8 in 16	430	—	—
	4 in 16	220	215	102%
	2 in 16	109	108	101%
	2 in 16	55	54	102%
	1 in 16	26	27	96%

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
5	16 in 16	754	—	—
	8 in 16	382	377	101%
	4 in 16	201	189	106%
	2 in 16	98	94	104%
	1 in 16	48	47	102%

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2015 Siemens Healthcare Diagnostics. All rights reserved.