

# MUTARS®-Münster

---



implantcast



**Proximal Tibia**  
Surgical Technique



# MUTARS®-Münster

---

## Proximal Tibia Surgical Technique

MUTARS® was developed in co-operation with  
Prof. Dr. W. Winkelmann  
and Prof. Dr. G. Gosheger,  
Clinic and Polyclinic for General Orthopedics  
at the University Hospital of Münster, Germany.  
MUTARS® has been in successful clinical use since 1992.

### Table of Contents

|  |    |
|--|----|
| Table of contents .....                        | 1  |
| Information about the silver-coating .....     | 2  |
| System overview .....                          | 3  |
| Assembling options/reconstruction length ..... | 4  |
| Surgical technique .....                       | 5  |
| Implants .....                                 | 21 |
| Instruments .....                              | 25 |

**Nota Bene:** The described surgical technique is the suggested treatment for the uncomplicated procedure.  
In the final analysis the preferred treatment is that which addresses the needs of the individual patient.

**Copyright note:** MUTARS®, implavit® and implatan® are registered trademarks of implantcast GmbH. The  
use and copy of the content of this brochure are only allowed with prior permit given by the implantcast  
GmbH.



### The silver-coating

Infections represent the most severe complications of tumour arthroplastic treatments. Although local and systemic antibiotic treatments are considered, the scientific literature reports of infection rates from 5 to 35 percent [1]. Reasons for these high rates are, for example, the long surgery time, the large incisions and the immunosuppression due to chemo therapy and radio therapy as well as the increasing resistance of the bacteria against antibiotic drugs.

Silver, in particular free silver ions, is well known for its broad-spectrum antimicrobial activity. The silver-coating has been shown to reduce bacterial colonization on the device surface.

Until now only non-articulating surfaces and surfaces without direct bony contact are coated with silver. In the catalogue information of this surgical technique you can find the supplement \*S indicating which MUTARS® components are available in a silver coated version. The eight digit REF-number receives an addition after the last digit (e.g. 5220-0020S).

### Important intra-operative instructions for the use of silver-coated implants

It is not permitted to flush the wound with antiseptics that contain H<sub>2</sub>O<sub>2</sub>, Iodine or heavy metals (such as Betaisodona®) and acetic acid during surgery since this can lead to a subsequent loss of effectiveness of the silver-coating due to their oxidative properties. Alternatively, solutions such as NaCl or Lavasept® and Prontosan® can be used. The additional use of antibiotic-containing bone cement can be an advantage particular in case of a septic revision.

### The TiN-coating for allergy prophylaxis

All metallic implant components release ions to their environment over time. In some patients such ions can elicit allergic reactions. Nickel, cobalt and chromium, which are elements of the base material CoCrMo of the articulating implant components, are considered the most frequently allergy eliciting metals [2]. The TiN-coating is biocompatible and acts like a barrier; the potential release of allergy eliciting ions of the base material is reduced to a minimum [3]. Also in clinical practice there have never been any evidence of allergic reactions with implants that have been TiN-coated showing an intact surface [5]. Therefore the TiN-coating on implant components is especially suitable for patients with sensitivity to nickel, chromium or cobalt [4][5].

Since almost all components of the MUTARS® tumor system consist of titanium alloy, this only concerns those components, which are made of a CoCrMo alloy. The REF-numbers of the TiN-coated implants have the suffix N after the last digit (e.g. 5720-0005N). Items which are available with silver and TiN-coating have the suffix SN after the last digit (e.g. 5720-0005SN).

**\*S:** Implants are available with silver-coating!

**\*N:** Implants are available with TiN-coating!

**\*SN:** Implants are available with silver and TiN-coating!

[1] Gosheger et al. 2004. Silver-coated megaendoprostheses in a rabbit model – an analysis of the infection rate and toxicological side effects. *Biomaterials* 25, 5547-5556.


[2] Eben R et al. (2009) Implantatallergieregister - ein erster Erfahrungsbericht. *Orthopäde* 38: 557-562

[3] Wisbey et al. (1987) Application of PVD TiN coating to Co-Cr-Mo based surgical implants. *Biomaterials*, 11

[4] Prof. Thomas LMU München Final Report Effect of a TiNbN or TiN surface coating on cobaltchromium- molybdenum and stainless steel test specimens regarding the release of nickel, chromium and cobalt: evaluation via eluate analysis and in-vitro cytokine release from peripheral human blood cells, Data on file

[5] Baumann A. (2001) Keramische Beschichtungen in der KTEP Standardlösung für Allergiker. *JATROS Orthopädie & Rheumatologie* 6: 16-17

## System Overview



### MUTARS® GenuX® stem

length: 160 mm

cementless

diameter: ø 12, 14, 16, 18 mm

cemented

diameter: ø11, 13, 15, 17 mm

### MUTARS® GenuX® femoral component

cementless and cemented

size 3, 4, 5

### MUTARS® coupling 15 mm

### MUTARS® PE-insert xsmall

### MUTARS® proximal tibia

### MUTARS® connecting part for modular proximal tibia

### MUTARS® tibial stem

length: 120 mm


cementless

diameter ø11-16 mm

cemented

diameter ø11, 13, 15 mm

information about loan set content and weight limitation on p.23





# MUTARS® Proximal Tibia

assembling options  
(length in mm)

| reconstruction | components        |                 |                 |       |
|----------------|-------------------|-----------------|-----------------|-------|
|                | femoral component | connecting part | extension piece | screw |
| 115            | x                 | 105             | -               | 25    |
| 135            | x                 | 125             | -               | 45    |
| 155            | x                 | 105             | 40              | 65    |
| 175            | x                 | 105             | 60              | 85    |
| 195            | x                 | 105             | 80              | 105   |
| 215            | x                 | 125             | 80              | 125   |
| 235            | x                 | 125             | 40 + 60         | 145   |

**Note:** Please notice that the amount of implants and instruments send with an individual shipment may differ from the information in the catalogue information of this brochure. Please make sure, during the preoperatively planning, that all necessary implants and instruments are available for the surgery.

## Tumor resection

Resect the tumor and measure the dimension of the explant.  
The minimum bone resection is 115 mm.



**figure 1a**

**figure 1b**

## Femoral preparation

Choose the correct femoral size with the femoral sizing template (fig. 1a). Fix the femoral alignment stylus to the femoral resection block and place the assembly on the distal femoral bone.

Slide the stylus as far as possible under the quadriceps muscle and assure that the stylus stays in contact with the anterior cortex. Open the intramedullary cavity using the 9 mm initiator drill (fig. 1b) and remove the instruments.



**figure 2**

Ream with the 10 mm reamer manual up to a stable fit ( If you use a bigger reamer, make sure that you stay 3-4 mm smaller than the preoperative planed stem) (fig. 2).

Adjust the rotation of the femoral resection block referencing on the posterior femoral condyles. Use two of the 3,2 mm fixation pins to fix the femoral resection block to the bone Remove the alignment stylus.



**figure 3**

## Anterior femoral resection

Check the resection with the resection check (fig. 3)

Place the saw capture to perform the anterior resection by the use of the ACS® saw blade (fig. 3).

## Posterior femoral resection

Change the position of the saw capture to resect the posterior condyles (fig. 4).



**figure 4a**

**figure 4b**

## Distal femoral resection

Mount the distal femoral cutting block<sup>1</sup> 6° facing the „L“ for the left knee or „R“ for the right knee to the femoral alignment guide and lock the resection block in such a way that the block corresponds with the mark on the medial side<sup>2</sup> of the alignment guide to determine the level of the distal bone resection (fig. 5).




figure 5

Place the whole assembly on the femoral resection block (fig. 6).

Leave the initiator drill in the bone for additional stability.




figure 6

The distal femoral cutting block should lie flush with the anterior resection plane<sup>3</sup> and it is attached to the bone using two predrilled pins<sup>4</sup> (fig. 7).






figure 7




To make sure that the distal cut is correct use the resection check (fig. 8).

**figure 8**



Remove all instruments except the distal femoral cutting block; add the saw capture and resect the distal femur with the ACS® saw blade (fig. 9).

**figure 9**



**figure 10**

## Femoral stem preparation

Place the finishing guide to the distal femoral bone. The guide should rest completely on both the distal and the anterior bone surface. Another option is to fix the finishing guide above the last reamer and the long stem sleeve. (fig. 10).

Slide the long stem sleeve (next size) into the guide. Drill with the reamer until the 200 mm mark reaches the top of the sleeve (fig. 12).

Please use the sleeves and reamers of growing diameters in the same way enhancing the diameter in 2 mm steps. For additional stability please slide the femoral reamer sleeve over reamer shaft (fig. 11). Please reference to table 1 and 2 to find out the recommended diameter for the bone preparation when a cementless or cemented stem fixation is planned.

**table 1: cementless implantation**

|                                   |
|-----------------------------------|
| femoral stem 12 mm → reamer 11 mm |
| femoral stem 14 mm → reamer 13 mm |
| femoral stem 16 mm → reamer 15 mm |
| femoral stem 18 mm → reamer 17 mm |

**table 2: cemented implantation**


|                                   |
|-----------------------------------|
| femoral stem 11 mm → reamer 12 mm |
| femoral stem 13 mm → reamer 14 mm |
| femoral stem 15 mm → reamer 16 mm |
| femoral stem 17 mm → reamer 18 mm |

Drill with the final reamer until the 200 mm mark reaches the top the sleeve (fig. 12).


Please leave the reamer and the sleeves in place and insert two fixations pins to stabilise the finishing guide in the correct M/L position (fig. 13).



**figure 11**



**figure 12**



**figure 13**

# MUTARS® Proximal Tibia




figure 14




figure 15



figure 16




figure 17a

Use the 18 mm reamer and the 18 mm sleeve (fig. 14) to remove additional bone, to allow a proper seating of the taper connection of the femoral component.

Please ream deep enough that the reaming part of the reamer will stay app. 1 cm out of the sleeve (fig. 15).



Then remove the 18 mm reamer and the 18 mm sleeve. Slide the box reamer guide into the finishing guide (fig. 16).

Remove the intracondylar bone by using the box reamer until it is stopped by the box reamer guide (fig. 17a and 17b).




figure 17b

Use the narrow ACS® saw blade to perform the chamfer cuts to finalise the femoral bone preparation.



**figure 18**

Start with the anterior chamfer cut (fig. 18) and perform the posterior chamfer (fig. 19).



**figure 19**

Cut out the anterior groove with the osteotome (fig. 20).

The femoral bone preparation is now performed.



**figure 20**




figure 21

## Assembling of the femoral implants

Choose the femoral stem of the correct size and version, cemented or cementless.

Place the stem into the assembling block of the instrument tray and connect the femoral component.

Use the femoral impactor and a mallet to enhance the taper connection (fig. 21).



figure 22

## Implantation of the femoral implants

Insert the femoral component with the assembled stem into the femoral bone and impact the components with the impactor inserted in the notch of the femoral component (fig. 22).

If sufficient seating is achieved, the impactor is removed.




figure 23



figure 24

For adjustment of the femoral component after seating, you could remove the femoral component by the use of the extractor mounted on the slide hammer (fig. 23 and 24).

Assemble the MUTARS® coupling 15 mm and the special MUTARS® instrument for locking mechanism. Therefore turn the attachment part of the lock by 100 degrees until it rests in the sleeve of the locking instrument (fig. 25a and 25b).



**figure 25 a and 25b**

Insert the lock into the intracondylar notch of the femoral joint (fig. 26).




**figure 26**

Use the socket wrench to turn the locking instrument and the lock clockwise by 180 degrees (fig. 26). The lock is correctly positioned when the attachment partly falls out of the sleeve of the locking instrument (fig. 27). Remove the locking instrument.



**figure 27**

## MUTARS® Proximal Tibia



**figure 28**

Fill the 4 suture holes of the Proximal Tibia component with non absorbable sutures (Ethibond is recommended) to allow the fixation of the attachment tube. Insert the PE insert x-small in the Proximal Tibia. Move the PE-insert towards the anterior locking rim and push it down at the posterior part until it is locked securely (fig. 28).




**figure 29**

Attach the coupling to the proximal tibia. Therefore, use the setting instrument or the setting instrument for coupling angled. The screw hole should be placed forward-turned to enable locking (fig. 29).



**figure 30**

The positioner is inserted into the screw hole of the short stem of the coupling mechanism (fig. 30).



**figure 31**

Remove the positioner and insert the screw with the 3.5 mm hex screw driver. To complete the connection, please insert the Multilock security screw also with the 3.5 mm hex driver (fig. 31).

## Tibial bone preparation

Use the medullary cavity reamer to prepare the tibial bone (fig. 32a and 32b).



figure 32a and 32b

## Cemented fixation

Ream the tibial medullary cavity preferably up to a depth of 130 mm with a rigid reamer that is 2 mm larger than the size of the tibial stem (fig. 33a and 33b).

## Cementless fixation

Ream the tibial medullary cavity preferably up to a depth of 130 mm with a rigid reamer that is 1,5 mm smaller than the size of the tibial stem (fig. 33a and 33b). Make sure that at least a 9cm contact between reamer and cortical bone is achieved.

## Remark

The use of a tibial rasp for a **cemented stem** is optional. Generally you can proceed with the trial reduction (see page 16).



figure 33a and 33b

## Cementless preparation

Choose the tibial rasp (fig. 34a and 34b) of the preoperatively planned size.

Assemble the tibial rasp of the appropriated size (see table 3 below), the sleeve and the slide hammer. Lock the rasp on the slide hammer by using the engineers' wrench (fig. 34).

| Stem size | Rasp Size |
|-----------|-----------|
| 12 mm     | 12 mm     |
| 13 mm     | 13 mm     |
| 14 mm     | 14 mm     |
| 15 mm     | 15 mm     |
| 16 mm     | 16 mm     |

table 3




figure 34a



figure 34b

## Optional technique for the use of cemented stems

If you want to prepare for a cemented stem with the tibial rasp, please use the rasp which is 2 mm larger than the preoperatively chosen cemented tibial stem.

That will provide a cement mantle of 1 mm thickness (table 4). Use the 16mm rasp to prepare for the 15 mm stem.

| Stem size | Rasp size |
|-----------|-----------|
| 11 mm     | 13 mm     |
| 13 mm     | 15 mm     |
| 15 mm     | 16 mm     |

**table 4**

Although the tibial stem is not curved it is recommended to mark the anterior aspect of the tibial bone to assure that the rotation of the final stem corresponds to the rotation of the rasp (fig. 35a).

Rasp the medullary cavity with the chosen tibial rasp (fig. 35b). Careful use of the slide hammer is recommended.

To prevent fractures of the cortical bone it is helpful to fix a bone forceps around the tibial bone while rasping.

### Remark

It is recommended to clean the rasp of bone chips during the rasping.

Leave the tibial rasp in the bone for the trialing.



**figure 35a and 35b**



figure 36

## Trial reduction

Attach the MUTARS® connecting part for the Proximal Tibia (length: 105 mm or 125 mm) to the tibial rasp (fig. 36). Mark the rotation of both components with methylene blue.

## Remark

For the **cemented procedure** bone rasps are usually not available. Please insert the cemented stem (without cement) for trialling purposes.

Connect the MUTARS® Proximal Tibia and the connecting part. Perform a trial reduction and check the joint stability and the rotational alignment (fig. 37).



figure 37

Adjust the rotation if necessary. If the joint line could not be restored correctly, it might be necessary to change the length of the tibial reconstruction by a change of the connecting part, or adding of an extension piece in conjunction with an enlarged tibial bone resection.

## Implantation of the tibial stem

Impact the MUTARS® tibial stem (fig. 38).

Insert the stem of the same size as the rasp if a **cementless stem** is used.  
To prevent fractures of the cortical bone it is helpful to fix a bone forceps around the femoral bone during impaction.

If a cemented implantation is planned insert the bone cement and use the **cemented stem** which is 2 mm smaller than the previously used reamer or rasp.

Remove all instruments, especially during the cement hardening to prevent bending moments (fig. 39).



figure 38




figure 39



## Final joint locking

Please attach the connecting part (and the possibly used extension pieces) to the tibial stem. Use the bar screw of the correct length (see table on page 2) to lock the component to the tibial stem (fig. 40a). Lock the screw by using the swing wrench and counter the assembly with the engineers' wrench SW 24 (fig. 40b).

figure 40a and 40b



Slide over the attachment tube. The trevira tube should be turned up inward on the end. If necessary cut the tube to the correct length (fig. 41.)

Combine the Proximal Tibia to the connecting part and insert the two locking screws into the anterior holes and lock them with the swing wrench (fig. 42).


figure 41

figure 42

## Fixation on the attachment tube

Please fix the tube to the upper part of the Proximal Tibia by using the previously inserted 4 sutures.


Fold the tube to achieve a very close covering of the components (fig. 43a and 43b). Insert additional sutures around the attachment pads of the implant components.



**figure 43a and 43b**

To reconstruct the extensor mechanism it is mandatory to perform a gastrocnemius muscle transfer. Release the muscle at its distal insertion (fig. 44a). Suture the muscle to the anterior portion of the attachment tube (fig. 44b).

Reinsert the extensor structures to the gastrocnemius muscle and the tube to restore a reasonable function of the joint.



**figure 44a and 44b**



## IMPLANTS

**\*S:** For anti-infective treatment, silver coated implants are available.

**\*N:** For anti-allergic treatment, TiN coated implants are available.

### MUTARS® Genux® stem cementless

*mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3*

|           |           |
|-----------|-----------|
| 5761-1612 | 12/160 mm |
| 5761-1614 | 14/160 mm |
| 5761-1616 | 16/160 mm |
| 5761-1618 | 18/160 mm |



### MUTARS® GenuX® stem cemented \*N

*mat.: implavit®; CoCrMo acc. to ISO 5832-4*

|           |           |
|-----------|-----------|
| 5762-1611 | 11/160 mm |
| 5762-1613 | 13/160 mm |
| 5762-1615 | 15/160 mm |
| 5762-1617 | 17/160 mm |



### MUTARS® GenuX® femoral component cemented \*N

*mat.: implavit®; CoCrMo acc. to ISO 5832-4*

|           |              |
|-----------|--------------|
| 5720-0310 | right size 3 |
| 5720-0315 | left size 3  |
| 5720-0320 | right size 4 |
| 5720-0325 | left size 4  |
| 5720-0330 | right size 5 |
| 5720-0335 | left size 5  |



### MUTARS® GenuX® femoral component cementless \*N

*mat.: implavit®; CoCrMo acc. to ISO 5832-4*

|           |              |
|-----------|--------------|
| 5720-0210 | right size 3 |
| 5720-0215 | left size 3  |
| 5720-0220 | right size 4 |
| 5720-0225 | left size 4  |
| 5720-0230 | right size 5 |
| 5720-0235 | left size 5  |



## IMPLANTS

### MUTARS® patella replacement cemented

mat.: UHMW-PE acc. to ISO 5834-2  
5720-1000 standard

### MUTARS® PE-insert

mat.: UHMW-PE acc. to ISO 5834-2  
5721-0013 xsmall

### MUTARS® modular proximal tibia incl. coupling 15 mm \*S

incl. screw for coupling, counter screw  
and MUTARS® screw for connecting part  
(2x)

mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to  
ISO 5832-3; mat. of coupling 15 mm:  
CoCrMo acc. to ISO 5832-12; PE-safety  
peg: UHMW-PE acc. to ISO 5834-2  
5750-0003

### MUTARS® connecting part for modular proximal tibia \*S

mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to  
ISO 5832-3

5750-0105 105 mm

5750-0125 125 mm

### MUTARS® extension piece \*S

mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to  
ISO 5832-3

5772-2504 40 mm

5772-2506 60 mm

5772-2508 80 mm

### intramedullary plug

mat.: UHMW-PE acc. to ISO 5834-2

0299-4000 small

0299-4010 large



## IMPLANTS

### MUTARS® screw

mat.: *implatan®*;  $TiAl_6V_4$  acc. to ISO 5832-3


|           |            |
|-----------|------------|
| 5792-1002 | M10x 25 mm |
| 5792-1004 | M10x 45 mm |
| 5792-1006 | M10x 65 mm |
| 5792-1008 | M10x 85 mm |
| 5792-1010 | M10x105 mm |
| 5792-1012 | M10x125 mm |
| 5792-1014 | M10x145 mm |
| 5792-1016 | M10x165 mm |
| 5792-1018 | M10x185 mm |
| 5792-1020 | M10x205 mm |
| 5792-1022 | M10x225 mm |



### MUTARS® tibial stem cemented \*N

mat.: *implavit®*; CoCrMo acc. to ISO 5832-4

|           |          |           |
|-----------|----------|-----------|
| 5750-0511 | 11/120mm | max.75 kg |
| 5750-0513 | 13/120mm |           |
| 5750-0515 | 15/120mm |           |




### MUTARS® tibial stem cementless HA

mat.: *implatan®*;  $TiAl_6V_4$  according to ISO 5832-3 with *implaFix®* HA; HA-coating acc. to ISO 13779-2

|            |          |           |
|------------|----------|-----------|
| 5750-1511* | 11/120mm | max. 60kg |
| 5750-1512  | 12/120mm |           |
| 5750-1513  | 13/120mm |           |
| 5750-1514  | 14/120mm |           |
| 5750-1515  | 15/120mm |           |
| 5750-1516  | 16/120mm |           |


\*stems with this size are not included in loan set and have to be ordered separately.



### MUTARS® attachment tube

mat.: *polyethylenterephtalat*

|           |          |
|-----------|----------|
| 5900-0300 | 300/35mm |
| 5900-0310 | 300/55mm |





### IMPLANTS

#### MUTARS® femoral spacer \*S

distal

mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3

| ll/rm     |             | rl/lm     |
|-----------|-------------|-----------|
| 5722-0530 | size 3/10mm | 5722-1030 |
| 5722-0535 | size 3/5mm  | 5722-1035 |
| 5722-0540 | size 4/10mm | 5722-1040 |
| 5722-0545 | size 4/5mm  | 5722-1045 |
| 5722-0550 | size 5/10mm | 5722-1050 |
| 5722-0555 | size 5/5mm  | 5722-1055 |



#### MUTARS® L- femoral spacer \*S


mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to ISO 5832-3

| ll/rm     |             | rl/lm     |
|-----------|-------------|-----------|
| 5722-1530 | size 3/10mm | 5722-2030 |
| 5722-1535 | size 3/5mm  | 5722-2035 |
| 5722-1540 | size 4/10mm | 5722-2040 |
| 5722-1545 | size 4/5mm  | 5722-2045 |
| 5722-1550 | size 5/10mm | 5722-2050 |
| 5722-1555 | size 5/5mm  | 5722-2055 |




# MUTARS® Proximal Tibia

---



## INSTRUMENTS


**MUTARS® basic container**  
7999-5712




**MUTARS® GenuX® femur container 1**  
7999-5723



**MUTARS® GenuX® femur container 2**  
7999-5722



**MUTARS® proximal tibia container**  
bottom tray  
7999-5734



**MUTARS® proximal tibia container**  
top tray  
7999-5734



**MUTARS® rigid drills container**  
7999-5735



## INSTRUMENTS

### **MUTARS® basic container**

7999-5712

### **MUTARS® universal impactor**

7210-0000



### **MUTARS® impact and extract sleeve**

7230-0000



### **MUTARS® socket wrench**

7420-0000



### **alternatively**

### **MUTARS® socket wrench**

7421-0000



### **MUTARS® swing wrench**

7411-0000




### **MUTARS® engineers wrench SW 24**

7490-0000



### **MUTARS® slide hammer**

7220-0001



### **MUTARS® rasp for femoral stem**

|           |       |
|-----------|-------|
| 7760-0112 | 12 mm |
| 7760-0113 | 13 mm |
| 7760-0114 | 14 mm |
| 7760-0115 | 15 mm |
| 7760-0116 | 16 mm |
| 7760-0117 | 17 mm |
| 7760-0118 | 18 mm |



### **handle for intramedullary plug**

7512-4001



### **MUTARS® medullary cavity reamer**

**cross-hole**  
4220-0000



# MUTARS® Proximal Tibia

## **MUTARS® GenuX® femur container 1** 7999-5723



**distal femoral cutting block 6°**  
4220-0018



**femoral alignment guide**  
4220-0028



**femoral alignment stylus**  
4220-0012



**fixation pins 3,2 mm x 97 mm**  
4223-0008 2x



**modular handle**  
4223-0015 2x



**modular handle "fast fix"**  
4223-0017 4x



**MUTARS® femoral chamfer**  
7630-1035



**slap hammer**  
4223-0005



**fitting block for stem**  
4223-4002



**ACS® SC sleeve for long stem primary**  
4220-4110 – 4220-4118 10 mm – 18 mm



**femoral reamer sleeve**  
4220-5112 – 4220-5118 12 mm – 18 mm



**Mutars® extractor for femoral component**  
7610-0002



**Mutars® impactor for femoral component**  
7610-0000



**initiator drill 9 mm**  
4220-0014



**osteotome size 2-6**  
4223-0060



**ic universal rasp**  
7512-1000



**MUTARS® femoral reamer sleeve**  
7630-1028



**saw capture „fast fix“**  
4221-0102



# MUTARS® Proximal Tibia

## **MUTARS® GenuX® femur container 2** 7999-5722



### **MUTARS® femoral cutting block**

|           |        |
|-----------|--------|
| 7630-1000 | size 3 |
| 7630-1002 | size 4 |
| 7630-1001 | size 5 |



### **femoral sizing template**

|           |          |
|-----------|----------|
| 4220-0010 | size 3/4 |
| 4220-0011 | size 5/6 |



### **MUTARS® femoral resection block**

|           |        |
|-----------|--------|
| 7320-0003 | size 3 |
| 7320-0004 | size 4 |
| 7320-0005 | size 5 |



### **MUTARS® femoral trial component**

|           |              |
|-----------|--------------|
| 7720-0210 | size 3 right |
| 7720-0215 | size 4 right |
| 7720-0220 | size 5 right |
| 7720-0225 | size 3 left  |
| 7720-0230 | size 4 left  |
| 7720-0235 | size 5 left  |



**MUTARS® proximal tibia container**  
**top tray**  
7999-5734

**fixation pin Ø3,2 mm / 77 mm**  
4223-0029 4x



**drill 126mm x 3,2 mm**  
4221-0019 2x



**MUTARS® patella drill**  
7351-0000



**MUTARS® rigid reamer**

|           |         |
|-----------|---------|
| 7700-2110 | 10,0 mm |
| 7700-2210 | 10,5 mm |
| 7700-2111 | 11,0 mm |
| 7700-2211 | 11,5 mm |
| 7700-2112 | 12,0 mm |
| 7700-2212 | 12,5 mm |
| 7700-2113 | 13,0 mm |
| 7700-2213 | 13,5 mm |
| 7700-2114 | 14,0 mm |
| 7700-2214 | 14,5 mm |
| 7700-2115 | 15,0 mm |
| 7700-2116 | 16,0 mm |
| 7700-2117 | 17,0 mm |



# MUTARS® Proximal Tibia

## **MUTARS® proximal tibia container** **bottom tray** 7999-5734



## **MUTARS® patella drill guide** 7350-0000



## **MUTARS® patella clamp** 7352-0000



## *or alternative* **MUTARS® patella - clamp** 7352-0001



## **pin inserter** 4223-0006



## **pin extractor** 4223-0007



## *or alternative* **ic- pin extractor** 7512-0800



## **setting instrument for coupling angled** 7751-1201



## **ic-T-handle Zimmer-Jakobs** 4223-0023



## **resection check** 4223-0009



## **MUTARS® instrument for locking mechanism** 7720-1201



## **extractor universal** 7512-2026



**hexagon screw driver short 3,5 mm**  
0280-1007



**MUTARS® positioner for locking mechanism**  
7610-0003



**setting instrument for locking mechanism**  
7751-1200




**MUTARS® impactor for PE-inlay**  
7210-0001



**MUTARS® trial PE-insert**  
7721-0003 extra small

*or alternative*

**MUTARS® trial PE-insert**  
7721-0013 xsmall



**MUTARS® rasp for tibial stem**

|           |       |
|-----------|-------|
| 7750-0312 | 12 mm |
| 7750-0313 | 13 mm |
| 7750-0314 | 14 mm |
| 7750-0315 | 15 mm |
| 7750-0316 | 16 mm |



## MUTARS® Proximal Tibia

### **MUTARS® rigid drills container**

7999-5735


### **MUTARS® rigid reamer**

|             |         |
|-------------|---------|
| 4220-4010.1 | Ø 10 mm |
| 4220-4011.1 | Ø 11 mm |
| 4220-4012.1 | Ø 12 mm |
| 4220-4013.1 | Ø 13 mm |
| 4220-4014.1 | Ø 14 mm |
| 4220 4015.1 | Ø 15 mm |
| 4220 4016.1 | Ø 16 mm |
| 4220 4017.1 | Ø 17 mm |
| 4220 4018.1 | Ø 18 mm |



### **ic adapter outside A/O, inside ic canulated**

7512-3602





implantcast GmbH  
Lüneburger Schanze 26  
D-21614 Buxtehude  
Germany  
phone: +49 4161 744-0  
fax: +49 4161 744-200  
e-mail: [info@implantcast.de](mailto:info@implantcast.de)  
internet: [www.implantcast.de](http://www.implantcast.de)



Your Local Distributor:

