CEMEX[®] **GENTA**

TECRES[®]
ADVANCING HIGH TECHNOLOGY

CEMEX® GENTA

SURGEON INFORMATION

EN

OVERVIEW

Cemex® Genta bone cements are all two-part (powder and liquid) radiopaque polymethylmethacrylate (PMMA) bone cements which include and release a medicinal substance: gentamicin sulphate. The PMMA bone cement includes also barium sulphate as radiopaque agent and, in some variants, colouring pigments which provide a green colour to the cement.

The green coloring increases the visibility of bone cement and to help the surgeon discriminate between bone tissue and bone cement (E102 Tartrazine Lake or FD&C Yellow N°5 and E133 Brilliant Blue or FD&C Blue N°1).

The cements are all supplied sterile, for single use, available in different viscosity variants.

Cemex® Genta High Viscosity are the high viscosity variant and they are ideal for manual application.

Cemex® Genta Low Viscosity are the low viscosity variant and they are ideal for application with injection device, but they can be used also for manual application.

Cemex® Genta Fast is a high viscosity variant that is suitable for manual application in cases where reduced application times are required.

Cemex® Genta ID green is the medium viscosity variant and it is ideal both for manual and syringe application.

For Cemex® Genta Low Viscosity and Cemex® Genta High Viscosity, the package contains a 40 g sachet of powder sterilized by Ethylene Oxide or Gamma Ray irradiation (according to the variant) and a 14,7 g vial of sterile liquid; the liquid is sterilized by filtration;

For Cemex® Genta Fast, the package contains two 20 g sachets of powder sterilized by Gamma Ray irradiation and a 8,35 g vial of sterile liquid; the liquid is sterilized by filtration;

For Cemex® Genta ID Green, the package contains a 40 g sachet of powder sterilized by Gamma Ray irradiation and a 16,7 g vial of sterile liquid; the liquid is sterilized by filtration;

Composition

	Cemex® Genta Low Viscosity REF: 1400/AG Cemex® Genta LV GAMMA REF: 1400/AG INT	Cemex® Genta High Viscosity REF: 1400/IG Cemex® Genta HV GAMMA REF: 1400/IG INT	Cemex® Genta Fast REF: 12A3100	Cemex® Genta ID Green REF: 12A2420
Liquid component:	14.7 g vial	14.7 g vial	8.35g (x2) vial	16.7 g vial
Methylmethacrylate	98.20% w/w	98.20% w/w	98.20% w/w	98,20 % w/w
N-N Dimethyl-p-Toluidine	1.80% w/w	1.80% w/w	1.80% w/w	1.80 % w/w
Hydroquinone	75 ppm	75 ppm	75 ppm	75 ppm
Powder component:	40 g sachet	40 g sachet	20 g sachet (x2)	40 g sachet
Polymethylmethacrylate	82.78 % w/w	82.78 % w/w	82.78 % w/w	82,68 % w/w
Barium Sulphate	10.00 % w/w	10.00 % w/w	10.00 % w/w	10,00 % w/w
Benzoyl Peroxide	3.00 % w/w	3.00 % w/w	3.00 % w/w	3,00 % w/w
Gentamicin sulphate	4.22 % w/w*	4.22 % w/w*	4.22 % w/w*	4,22 % w/w*
Colouring pigments	na	na	na	0,10 % w/w

*Equivalent to 1.0 g (1.0 M.I.U.), 2.5% Gentamicin base in 40 g unit and 0.5 g (0.5 M.I.U.) Gentamicin base in 20g unit

Fully formed polymethylmethacrylate (PMMA) bone cement loaded with gentamicin sulphate will have the following final chemistry formulation:

	Cemex® Genta Low Viscosity REF: 1400/AG Cemex® Genta LV GAMMA REF: 1400/AG INT	Cemex® Genta High Viscosity REF: 1400/IG Cemex® Genta HV GAMMA REF: 1400/IG INT	Cemex® Genta Fast REF: 12A3100	Cemex® Genta ID Green REF: 12A2420
Polymethylmethacrylate (% w/w)	88,85	88,85	89,17	89,17
Barium Sulphate (% w/w)	7,31	7,31	7,05	7,05
Gentamicin sulphate (% w/w)	3,09	3,09	2,98	2,98
Methylmethacrylate (% w/w)	0,27	0,27	0,27	0,27
N-N Dimethyl-p-Toluidine (% w/w)	0,48	0,48	0,53	0,53

* No other components above 0.1% w/w are present in the final product

INTENDED PURPOSE

CEMEX® GENTA bone cements are intended for the fixation of joint prosthesis (hip, knee) implants to the

host bone

INDICATIONS FOR USE

In particular, CEMEX® GENTA bone cements are indicated where there is the risk or presence of infections caused by organisms sensitive to Gentamicin.

CONTRAINDICATIONS

Cemex® Genta are contraindicated:

- where the loss of musculature or neuromuscular compromise in the affected limb would render the surgical procedure unjustifiable.
- in patients who are allergic or sensitive to any of its components including the coloring pigments (E102 Tartrazine Lake or FD&C Yellow N°5 and E133 Brilliant Blue or FD&C Blue N°1) for Cemex® Genta ID Green
- Severe renal impairment
- Pre-existing calcium metabolism disorder

Cemex® Genta must be considered carefully in the presence of myasthenia gravis.

INTENDED USERS

Cemex® Genta are intended for use by suitably qualified orthopaedic surgeons, experienced in arthroplasty procedures, in an operating theatre environment.

PATIENT TARGET GROUP(S)

Skeletally mature patients undergoing a joint arthroplasty.

CLINICAL BENEFITS & PERFORMANCES

Main:

- long-term fixation of joint replacement prostheses to living bone tissue

Ancillary:

- local antibiotic release
- bacterial anti-adhesivity of susceptible species

SUMMARY OF SAFETY AND CLINICAL PERFORMANCES

The Summary of Safety and Clinical Performances (SSCP) is available on the European Database on medical devices (EUDAMED) at the link <https://ec.europa.eu/tools/eudamed>, within six months after the date of publication of the notice of its full functionality.

ADVERSE EFFECTS

The following adverse effects may be associated with use of bone cement. Although some effects are not directly attributable to the device itself, the surgeon should be aware of these potential complications and be ready to treat the patient accordingly.

Serious:

- Anaphylaxis
- Cardiac arrest
- Cerebrovascular incidents
- Myocardial infarct
- Pulmonary embolism
- Sudden death

More frequent

- Haematoma-haemorrhage
- Infection of surface/deep surgical wound
- Short-term cardiac irregularities
- Thrombophlebitis
- Transitory fall in blood pressure
- Trochanteric bursitis

Other referred reactions:

- Adverse soft tissue reaction
- Allergic pyrexia
- BCIS, bone cement implantation syndrome (hypoxia, hypotension, cardiac arrhythmias, increased pulmonary vascular resistance, and cardiac arrest)
- Bladder fistula
- Breakage of the bone cement
- Delayed sciatic nerve entrapment from extrusion of bone cement beyond the region of its intended application
- Dysuria
- Elevated serum gamma-glutamyl-transpeptidase (GGTP) up to 10 days post-operation
- Haematuria
- Heterotopic new bone formation
- Intestinal obstruction because of adhesion and stricture of the ileum from heat released during polymerization.
- Local neuropathy
- Local vascular erosion and occlusion
- Loosening or detachment of prosthesis
- Pain and/or loss of function
- Trochanteric separation
- Gentamicin-associated toxicity (nephrotoxicity, ototoxicity)

Other adverse reactions associated with systemic use of Gentamicin therapy include hypersensitivity, anaphylactic reactions, nausea, vomiting, urticaria, reversible granulocytopenia, anemia, blood dyscrasia, convulsions, central nervous toxicity, abnormal hepatic function, hypomagnesaemia, stomatitis, purpura, allergic contact sensitization and neuromuscular blockade.

The postoperative follow-up should be closely monitored by the surgeon for patients reporting a previous

experience of sensitivity or allergic reaction to antibiotics, in particular of the aminoglycosides' family.

WARNINGS- PRECAUTIONS FOR USE

The present instructions for use leaflet must be read carefully.

The use of Cemex® Genta should be carefully considered in patients with coagulation disorders and in patients with severe cardio-pulmonary insufficiency. Do not use after the expiration date since the effectiveness of the device may be compromised.

Using the device under conditions other than the indicated use is unlikely to provide benefit to the patient and increases the risk of the development of drug-resistant bacteria

For the use of bone cement, it is required a high level of coordination between the anaesthetist and the surgeon. The surgeon must communicate to the anaesthetist when the cement is about to be introduced during the procedure.

In some cases, events defined as "bone cement implantation syndrome" (BCIS) may occur which are characterized by a number of clinical features that include hypoxia, hypotension, cardiac arrhythmias, increased pulmonary vascular resistance (PVR), and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology. These phenomena are commonly associated with, but are not restricted to, cemented hip arthroplasty which usually occur at one of the five stages in the surgical procedure: femoral reaming, acetabular or femoral cement implantation, insertion of the prosthesis or joint reduction (Donaldson et al., 2009, Br J Anaesth).

Patient blood pressure should be monitored during and immediately after application of bone cement. While inserting the prostheses over-pressurization of the bone cement should be avoided to minimize the possibility of pulmonary embolism.

Do not re-sterilize and/or re-use the device as it is designed for single-use on a single patient.

Avoid re-sterilization as it may cause infection risks for the patient. Re-sterilization may change and affect negatively the product and its performances, including the effectiveness of the antibiotic, causing a malfunction of the same with serious risks for the patient's well-being.

Never divide the product into two or more portions, in order to use it for other clinical applications or at different times. This re-use may lead to an error in the correct proportion of powder-to-liquid mix. It could also cause a sterility loss.

For Cemex® Genta Fast: if only one dose (20 g powder sachet and one 8,35 g liquid vial) of the device has been used during the surgical operation, discard the remaining dose as described in the section Disposal. This would be considered a re-use of the product. It could cause a sterility loss and difficulties in the device use due to the missing of the instructions for correct use.

Cemex® Genta should be used with caution:

- in conjunction with other nephrotoxic or ototoxic drugs (e.g. other aminoglycosides, amphotericin B, cephalosporins, ciclosporins, cisplatin, clindamycin, colistin, polymyxin B, foscarnet, vancomycin, viomycin.) since they can increase gentamicin toxicity (the low serum level of gentamicin makes this interaction very unlikely to occur);
- in patients who are predisposed to or who have pre-existing clinical conditions that would put them at risk for gentamicin toxicity (e.g. renal dysfunction, hearing difficulties, dehydration, advanced age, taking drugs which may affect kidneys, undergoing general anaesthesia, etc).

In such cases, patients should be monitored for toxic blood levels of gentamicin (mainly in the first days following implantation), as well as for kidney function.

Simultaneous use of gentamicin and strong diuretics such as ethacrynic acid or furosemide may intensify the ototoxic effect of gentamicin by altering the antibiotic concentrations in serum and tissue.

The use of Cemex® Genta may result in overgrowth of non-susceptible organisms. If overgrowth of non-susceptible organism occurs, appropriate therapy should be initiated.

In cases of significant obesity gentamicin serum concentrations should be closely monitored and a reduction in dose should be considered.

Co-administration with the following agents should be avoided:

- neuromuscular blocking agents such as succinylcholine, botulinum toxin risk of toxicity due to enhanced neuromuscular block
- anticoagulants such as warfarin and phenindione. as Gentamicin has been known to potentiate them
- Cholinergics: antagonism of effect of neostigmine and pyridostigmine.
- Bisphosphonates: increased risk of hypocalcaemia.
- Sulfites in susceptible people, especially those with a history of asthma or allergy: risk of allergic-type reactions including anaphylactic symptoms and bronchospasm.

In vitro elution studies (microbiological method) have shown that the daily release of gentamicin never exceeds the recommended systemic adult dose for gentamicin (5-7 mg/kg/day, adult with normal renal function) according to Goodman and Gilman's recommendations¹.

It is therefore unlikely that the amount of gentamicin absorbed locally from Cemex® Genta will result in serum levels in the toxic range.

Owing to the presence of gentamicin adverse events typical of these substances may occur, even if occurrence is unlikely because of the very low serum levels reached.

1 *The Pharmacological Basis of Therapeutics 13th Edition (Brunton L, Hilal-Dandan R, Knollmann B) - McGraw Hill, New York. Chapter 58: Aminoglycosides (McDougall C)*

Gentamicin is potentially nephrotoxic and/or ototoxic.

Despite the risk of accumulation is low, care should be taken, and the levels of gentamicin should be monitored in patients with impaired renal function, and in patients with previous neuromuscular disorders (e. g. myasthenia gravis, Parkinson's disease), in particular if they also have a history of renal insufficiency. In case of administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified (using calcium salts it may be reversed). However, this is quite unlikely in patients with healthy kidneys, whereas in patients with impaired renal function, the probability of interaction increases in proportion to serum levels.

Plasma trough concentrations of gentamicin higher than 2µg/ml for periods longer than 10 days have been associated to toxicity.

USE DURING PREGNANCY, BREAST-FEEDING AND IN CHILDREN

There are no available tests which demonstrate the safe use of bone cement during pregnancy or breast-feeding, and in children.

During pregnancy and breast-feeding the surgeon should weight the benefit for the mother against the potential risk to the child before using the cement.

It is advisable not to implant bone cement during pregnancy and breast-feeding. In these patients, bone cement should only be used in critical, life endangering situations.

In the event that the bone cement is used during pregnancy or if the patient becomes pregnant while the bone cement is in use, she should be notified of the possible dangers to the fetus. Women of childbearing potential should consider the benefits and dangers associated with the use of the product.

The bone cement is indicated for use in children (skeletal young patients) only when the surgeon deems that saving the joint through other forms of intervention is not possible.

STERILE PACKAGING INFORMATION

- For ref codes 1400/AG and 1400/IG: the package consists of an aluminium wrapping containing a hermetically sealed blister pack whose contents (40 g sachet of powder and 14.7 g vial of liquid) are sterilized by ethylene oxide. The liquid is sterilized by filtration and the powder by ethylene oxide.
- For codes 1400/AG INT and 1400/IG INT: the package contains a 40g sachet of powder sterilized by gamma ray and a blister pack sterilized by ethylene oxide containing a 14.7g vial of sterile liquid. The liquid is sterilized by filtration.
- For code 12A3100: the package has 2 sachets (each containing 20g of sterile powder) and 2 phials (each containing 8.35g of sterile liquid). The liquid is sterilized by filtration and the powder by Gamma-ray sterilization.
- For code 12A2420: the package contains a 40g sachet of powder sterilized by gamma ray and a blister pack sterilized by ethylene oxide containing a 16.7g vial of sterile liquid. The liquid is sterilized by filtration.
- Please discard Cemex® Genta if the external powder sachet and/or the vial blister are damaged or unintentionally opened before use because the sterility of the device may be compromised.
- Do not re-sterilize.

STORAGE & TRANSPORT CONDITIONS

- Store and transport at a temperature below 25°C and in a cool and dry environment. However, if the temperature during the transport cannot be maintained under 25 °C with a relative humidity < 90%, the medical device can tolerate a temperature between 0 °C and 40 °C with a relative humidity less than 90% for maximum 5 days.
- The device must be shipped and transported in a double wall corrugated box with ECT ≥ 44 lbf/in (or ECT ≥ 7.71 kN/m)
- Prior to using Cemex® Genta it is advised to store the product at a temperature of 23°C ± 1°C for 24 hours before surgery.
- Bone cements are temperature-sensitive. Temperatures of more than 23°C for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. Please see the graphs reporting the different preparation phases timings located at the end of this leaflet.

UTILIZATION and IMPLANTATION

Familiarity with proper use of bone cement, and in particular with Cemex® Genta, is essential for successful use of the device.

It is important to maintain strictly aseptic surgical procedures. Any deep infection of a surgical wound is a serious risk and will affect the successful outcome of the technique. Some infections may appear later, even years after surgery, without clinical signs.

The prosthesis that the surgeon chooses for the implantation must be compatible with the use of bone cement.

It is recommended that the following equipment and accessories be available for the procedure:

- Mixing and delivery accessories suitable for bone cements, ensure that they are compatible with PMMA bone cement (e.g. rubber and elastomer materials should be avoided).
- Equipment for joint arthroplasty surgery

- Cement Restrictor (if needed according to patient anatomical condition)

TEMPERATURE EFFECTS ON PREPARATION AND APPLICATION TIMING

Bone cements are sensitive to temperature. The temperature of the storage and of the operative room affects the preparation and application of the cement.

Temperature's effect on cement setting time was evaluated with a laboratory test according to ISO 5833 - *Implants for surgery — Acrylic resin cements*.

The mixing times, waiting times, working times and setting times of Cemex® Genta are shown on the diagrams at the end of these instructions. Please note that they are stated for guidance only, because different factors can affect them. These include the temperature, humidity, ventilation, mixing technique (speed, use of mixer), the thoroughness of mixing, the inadvertent inclusion of external substances inside the cement (such as saline solution, blood etc.), and the pre-heating of the prosthesis component itself. Any increase in temperature of the working environment / cement components / mixing instrumentation / prosthesis reduces the preparation times. Equally, lower temperatures increase such times.

PREPARATION

- Make sure that the inner packages and components are undamaged. Powder should be smooth (no agglomerations) and not yellow or brown in color. The content in the vial should appear as a low viscosity liquid. If the powder has a yellowish or brownish color or if the liquid is syrupy do not use product. This indicates the product has not been stored properly.
- The preparation of Cemex® Genta (single dose) is done by mixing the entire contents of one sachet of powder with one vial of liquid as per the instructions below. If required, more doses may be mixed. Different lot numbers of Cemex® Genta may be used together when properly mixed.
- The surgeon will determine the amount of cement needed based on the clinical application and needs.

Warnings:

- The liquid component of bone cement is a powerful lipid solvent. Contact with monomer to skin and mucous membranes should be avoided since cases of contact dermatitis have been noted in susceptible subjects.
- In order to reduce the possibility of negative reactions it is recommended to wear gloves to ensure adequate protection against the penetration of the monomer (methyl methacrylate) into the skin.
- PVP (three-layer polyethylene, ethylene-vinyl alcohol-copolymer, polyethylene) and Viton/Butyl gloves have proved to provide good protection over an extended period. Putting on two pairs of gloves – one pair of polyethylene surgical gloves over a pair of standard latex surgical gloves – has also proved to offer adequate protection.
- The use of latex or polystyrene-butadiene (rubber) gloves, however, must be avoided. Please request the confirmation of your glove supplier whether the respective gloves are suitable for the use of this cement.
- The operative area must be correctly ventilated as the liquid component is both flammable and volatile. The liquid monomer and vapors must never directly be exposed to flames or heated items. It has been reported that monomer vapors have ignited by use of electrocautery devices in surgical sites near newly implanted bone cement.
- Use caution while mixing the liquid and powder components of the bone cement to prevent prolonged exposure to the concentrated vapors of liquid monomer. This exposure may cause irritation of the respiratory tract, eyes, and in some cases the liver, sleepiness and drowsiness.
- Personnel with contact lenses should not prepare bone cement or be near the mixing process.

Step 1 – Opening

Open the unit container and place, onto a sterile working surface, the powder sachet and the liquid vial.

Step 2 – Mixing

- Break open the vial and pour all the liquid into a proper mixing container.
- Open the powder sachet and pour all the powder over liquid.
- Mix the cement with a spatula from the outside of the container towards the center to minimize the presence of air bubbles.
- Use the spatula delicately to remove any lumps of unmoistened powder into the overall cement dough.
- Mixing time is between 1-1.5 minutes, but the actual time depends on the temperature, the humidity and the mixing technique, and it is determined by the surgeon's experience.

Warnings:

- Do not open the liquid vial over the mixing container in order to avoid the risk of glass fragments entering the dough.
- The cement should not come into contact with the gloved hand until the consistency of dough has been reached.
- Do not change the ratio between liquid and powder components.
- Do not add foreign bodies or other substances (including other antibiotics) to Cemex® Genta.

Step 3 – Handling

- For manual application: after mixing, continue moving the dough until it no longer sticks to the gloves (see the diagrams). At this point, the mass is ready to be applied.

- For syringe application: immediately after mixing, introduce the cement in an appropriate sterile delivery device.

The temperature and humidity of the storage and operating room, of the mixing/delivery accessories used and of the surgeon's hands may determine differences in the timing for cement preparation and application. These are determined by the surgeon's experience.

Step 4 – Application

- The bone cavity must be properly irrigated with Ringer or saline solutions and dried prior to the application of bone cement. This reduces the risk of inclusion of medullary fat, bone fragments or other foreign bodies in the vascular system.
- During this phase, apply the cement.
- For syringe application: extrusion of the cement must begin in the distal area of the bone cavity and then flow into the proximal area. In order to avoid the inclusion of possible air bubbles, the cement flow must be kept as uniform as possible.

Warnings:

- It is important to avoid the presence of liquid between the bone tissue and the cement. The bone surface must be dried with gauze and/or suction catheters before and during the cementation process.
- Avoid over-pressurization of the bone cement because this may lead to extrusion of the bone cement beyond the site of its intended application and damage to the surrounding tissues. Moreover, the over-pressurization of the bone cement should be avoided to minimize the possibility of pulmonary embolism during the prosthesis insertion.
- The surface of the implant that it is intended to be cemented should be covered by a uniform coating of bone cement. It is important to apply an optimal thickness of bone cement.
- Bone cement reaches temperatures higher than physiological ones during the polymerization reaction. Polymerization of the bone cement is an exothermic reaction that occurs while the bone cement is setting (hardening) in situ. The released heat may damage bone or tissue adjacent to the implant.

Step 5 – Prosthesis insertion

- Once the cement has been applied, the prosthesis can be positioned.
- Hold the prosthesis firmly in place until the cement has hardened. This is essential to ensure optimal implantation results.
- Excess cement should be removed prior to setting.
- Setting time of the cement depends on the temperature and humidity of the storage and operating room, type of cement, and the manipulation.
- The increase of the temperature of the cement indicates the proceeding of the polymerization reaction. The cement cools spontaneously after few minutes, which is the end of the reaction and time when the surgeon can release the prosthesis.

Warnings:

- **Application of the prosthesis should be completed as quickly as possible as the temperature of the host bone cavity accelerates cement polymerization.**
- Inadequate fixation or unanticipated postoperative events may affect the cement-bone interface and lead to micromotion of cement against bone surface. A fibrous tissue layer may develop between the cement and bone, and loosening of the prosthesis may occur.

POST-OPERATIVE FOLLOW-UP

The cement or the prosthesis (or both) may loosen or fracture due to incorrect cement insertion technique, trauma, or latent infection. It is therefore advisable to follow-up all patients regularly and in the long-term following surgery.

The follow up may include drug monitoring tests, renal function monitoring tests, serum electrolytes, urinalysis, and audiograms (in the elderly and/or dehydrated patient in whom there is a higher risk of adverse events associated with Gentamicin use). Elderly patients may have reduced renal function that may not be evident in the results of routine screening tests, such as BUN or serum creatinine. A creatinine clearance determination may be more useful.

BONE CEMENT REMOVAL DUE TO PROSTHESIS REVISION

In case of aseptic/septic prosthesis revision a radical debridement shall be performed to remove all cement and potentially infected and devitalized tissues using common instruments for this type of procedure (e.g. osteotomes, drills, hooks, curettes, forceps, endoscopic instruments, high-energy ultrasounds, lithotripsy, laser, bone cement segmental remover). The procedure to remove the implant is a case by case decision by the responsible surgeon depending on the condition of the implant and the patient.

MRI SAFETY INFORMATION

Cemex® Genta are considered as "MR Safe" in accordance with the definition of the standard ASTM F2503, therefore it does not pose any risk in the Magnetic Resonance Environment.

IMPLANT CARD

The implanted device Cemex® Genta is provided with an implant card to be filled by the Healthcare Institution/provider with patient data, and to be supplied to the patient. This implant card reports the

information allowing the identification of the implanted device. The information to be supplied to the patient about warnings, precautions and measures to be taken is reported on the patient leaflet and on the following website: <https://ifu.tecres.com>.

The patient leaflet included in the device package must be supplied to the patient.

Instruction for Implant Card completion

The following data shall be filled by the HealthCare Institution/provider on the implant card front side:

1. Name of the patient or patient ID.
2. Date of implantation (day/month/year).
3. Name and address of the Healthcare Institution/provider.

The image shows a rectangular implant card with a white background and a black border. On the left side, there is a 3D cube logo with the text "TECRES" and "ADVANCED MEDICAL TECHNOLOGY" below it. To the right of the logo, the text "TECRES" is written in a large, bold font, followed by "IMPLANT CARD" in a smaller font. Below the logo and text, there are three rows of input fields. Each row starts with a small icon: a person with a question mark, a calendar with the number 31, and a person with a plus sign. To the right of each icon is a horizontal line for text entry. Above each line is a small circle containing a number: 1, 2, and 3. At the bottom left of the card, there is a small icon of a laptop with a person, followed by the URL "https://ifu.tecres.com". On the right side of the card, there is a large, stylized black and white graphic that looks like a folded piece of paper or a ribbon.

Once filled the Implant Card shall be detached from its support and provided to the patient.

Explanation of overall symbols present on the implant card (front and back sides) are reported under section "SYMBOLS on Implant Card" of the present leaflet.

DISPOSAL

The surplus bone cement should be allowed to set prior to disposal, which should be performed according to local waste regulations for unused or expired drugs since the device contain antibiotic.

In case of expired shelf-life, left-over unused dose of Cemex® Genta Fast (in case of use of only one single dose of the two provided), damaged or unintentionally opened packaging, please dispose the powder component of Cemex® Genta according to local waste regulations for unused or expired drugs since the device contain antibiotic. Instead the liquid component, due to its volatility and flammability, should be evaporated in a properly ventilated hood or absorbed by an inert material and transferred into a suitable container for disposal in a landfill.

SERIOUS INCIDENT REPORTING

Any serious incident potentially attributable to Cemex® Genta shall be reported to the manufacturer, TECRES S.p.A., and to the competent authority in which the user and/or the patient is established.

The manufacturer can be reached at the following contacts:

- email: info@tecres.it
- tel.: +39 045 9217311

The contact information of the competent authority located in the European Member State in which the patient is established is reported on European Commission website, under section "medical devices" at the following link: https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en

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
The image shows a form titled "TECRES IMPLANT CARD". At the top left is the TECRES logo, which consists of a cube-like icon and the text "TECRES" with "ADVANCED 3D TECHNOLOGY" underneath. To the right of the logo is the title "TECRES IMPLANT CARD". Below the logo are three rows of input fields, each preceded by an icon and a circled number. The first row has a person icon with a question mark and a circled "1". The second row has a calendar icon with "31" and a circled "2". The third row has a person with a plus sign icon and a circled "3". Below these fields are several horizontal lines for writing. At the bottom left is a computer icon and the URL "https://ifu.tecres.com". On the right side of the form is a large, stylized graphic element that looks like a 3D block with a white diagonal line.

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The image shows a rectangular card with a white background and a black border. On the left side, there is a logo consisting of a stylized cube with the word 'TECRES' and 'ADVANCED MEDICAL TECHNOLOGY' below it. To the right of the logo, the text 'TECRES® IMPLANT CARD' is printed. Below the logo, there are three rows of input fields. Each row starts with a small icon: a person with a question mark, a calendar icon with '31', and a person with a plus sign. To the right of each icon is a horizontal line for text. Above each line is a small circle containing a number: 1, 2, and 3 respectively. At the bottom left of the card, there is a small icon of a laptop and the URL 'https://ifu.tecres.com'. On the right side of the card, there is a large, stylized black and white graphic element that looks like a folded piece of paper or a ribbon.

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
The image shows a form titled "TECRES IMPLANT CARD" with the TECRES logo and "ADVANCED BONE TECHNOLOGY" tagline. The form contains several fields for patient information, each preceded by an icon and a numbered circle (1, 2, or 3). The icons are: a person with a question mark, a calendar showing '31', and a person with a plus sign. Below these are three horizontal lines for text. At the bottom left is a laptop icon and the URL "https://ifu.tecres.com". On the right side of the form is a large, stylized black and white graphic element.

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