


# Safari<sup>2</sup><sup>TM</sup>

## Pre-Shaped Guidewire

Boston  
Scientific  
Advancing science for life<sup>TM</sup>

Shaped  
and Sized  
for Safety




**The Safari<sup>2</sup> Pre-Shaped Guidewire enables safe treatment of a broad range of patients with valvular disease**

- Enhanced wire predictability with superior shape retention<sup>\*,†</sup>
- Streamlined device delivery through optimized rail support<sup>‡,§,\*\*</sup>
- Widest guidewire choice with three curve sizes

**AVAILABLE IN  
EXTRA SMALL SIZE**


## Curve Tip Configuration



### Guidewire Specifications

<b>Indications For Use</b>	The Safari <sup>2</sup> Guidewire is intended to facilitate the introduction and placement of interventional devices within the chambers of the heart including those used within transcatheter aortic valve procedures.		
<b>Outer Diameter</b>	0.035" (0.889 mm)		
<b>Length</b>	275 cm		
<b>Core Material</b>	Stainless Steel		
<b>Spring Coil Material</b>	Stainless Steel		
<b>Spring Coil Coating</b>	LUBRIGREEN <sup>TM</sup> PTFE		
<b>Configurations</b>	Extra Small Curve	Small Curve	Large Curve
<b>Distal Grind</b>	16 cm	16 cm	18 cm
<b>Curve Height</b>	3.2 cm	4.2 cm	5.0 cm
<b>Curve Width</b>	2.9 cm	4.2 cm	4.9 cm
<b>Package Contains</b>	Guidewire, J-Straightener		
<b>Use</b>	Single use - not for reuse, reprocessing or resterilization		
<b>Sterilization</b>	Ethylene Oxide		
<b>Shelf Life</b>	3 years		

## Smooth Tapered Core Wire Design



### Order Information

Order Number (GTIN)	Ref/Catalog Number	Description	Outer Diameter (inches)	Outer Diameter (mm)	Length	Quantity per Box
08714729887614	H749 <b>39406XS</b> 1	Safari <sup>2</sup> Guidewire Extra Small Curve	0.035	0.889	275 cm	5 Pack
08714729887591	H749 <b>39406S</b> 1	Safari <sup>2</sup> Guidewire Small Curve	0.035	0.889	275 cm	5 Pack
08714729887577	H749 <b>39406L</b> 1	Safari <sup>2</sup> Guidewire Large Curve	0.035	0.889	275 cm	5 Pack
08714729887638	H749 <b>39407XS</b> 0	Safari <sup>2</sup> Guidewire Extra Small Curve	0.035	0.889	275 cm	1 Single
08714729887652	H749 <b>39407S</b> 0	Safari <sup>2</sup> Guidewire Small Curve	0.035	0.889	275 cm	1 Single
08714729887645	H749 <b>39407L</b> 0	Safari <sup>2</sup> Guidewire Large Curve	0.035	0.889	275 cm	1 Single

The C-code used for this product is C11769 Guidewires. C-codes are used for hospital outpatient device reporting for Medicare and some private payers. Boston Scientific is not responsible for the correct use of codes on submitted claims; this information does not constitute reimbursement or legal advice.

# Safari<sup>2</sup><sup>™</sup>

## Pre-Shaped Guidewire

Bench testing performed by Boston Scientific. Bench test results may not necessarily be indicative of clinical performance. Data on file.

\* Curve Retention test data with Safari<sup>2</sup> Extra Small, Confida<sup>™</sup> Brecker, Amplatz Extra Stiff<sup>™</sup>, and Amplatz Super Stiff<sup>™</sup> Guidewires (n = 1).

† Guidewire Compression test data with Safari<sup>2</sup> Extra Small, Confida Brecker, Amplatz Extra Stiff, and Amplatz Super Stiff Guidewires (n = 1).

‡ Nosecone Tracking test data with Safari<sup>2</sup> Extra Small, Confida Brecker, Amplatz Extra Stiff, and Amplatz Super Stiff Guidewires (n = 1).

§ Simulated Use Static Friction, Dynamic Friction, and Distal Curve Displacement test data with Lotus<sup>™</sup> Valve System with Safari Large and Amplatz Super Stiff Guidewires (n = 4).

\*\* Simulated Use Friction, Dynamic Friction, and Distal Curve Displacement with Evolut<sup>™</sup> R with Safari<sup>2</sup> Extra Small and Confida Brecker Guidewires (n = 3).

The Safari<sup>2</sup> Guidewire is manufactured by Lake Region Medical and distributed by Boston Scientific Corporation.

Illustrations for information purposes – not indicative of actual size or clinical outcome.

All trademarks are the property of their respective owners.

### Safari<sup>2</sup> Preshaped TAVR TAVI Guidewire

**Intended Use/Indications for Use:** The Safari<sup>2</sup> Guidewire is intended to facilitate the introduction and placement of interventional devices within the chambers of the heart including those used within transcatheter aortic valve procedures.

**Contraindications:** This wire is not intended for use in the cerebrovasculature or coronary arteries.

**Warnings:**

- The Safari<sup>2</sup> Guidewire should be used only by physicians trained in the introduction and placement of interventional devices including those used within transcatheter aortic valve procedures.
- Carefully read all instructions prior to use. Observe all warnings and precautions. Failure to do so may result in complications.
- Prior to use, inspect for damage. If damaged, DO NOT USE.
- Monitor wire position throughout the procedure for proper placement of curve and distal tip.
- Do not torque this guidewire.
- Use extreme caution and careful judgment in patients for whom anticoagulation is not indicated.
- Inserting the distal end of the guidewire through a Tuohy-Borst Adapter may result in damage to the tip of the guidewire.
- If contrast agents are used, use extreme caution in patients who have had a severe reaction to contrast agents and who cannot be adequately premedicated.
- When advancing or removing the guidewire, always use fluoroscopic guidance with radiographic equipment that provides high-resolution images. Never position the guidewire blindly, as this may result in misplacement, dissection or perforation.
- Exercise care in handling of the guidewire during a procedure to reduce the possibility of accidental breakage, bending, kinking, or coil separation. Resulting guidewire fractures might require additional percutaneous intervention or surgery.
- Never advance the guidewire against the resistance without first determining the reason for resistance under fluoroscopy. Excessive force against resistance may result in damage to the catheter or vessel/organ. Care should be taken when advancing a guidewire after device deployment.
- The curve of the Safari<sup>2</sup> Guidewire should be constrained within a catheter during insertion into or withdrawal from the body or treatment site.
- The Safari<sup>2</sup> Guidewire is manufactured with a double curve; attempts to modify may alter its performance. Alterations to curve may lead to complications including Perforation or Dissection, Mitral Valve Regurgitation, Pericardial Effusion, Cardiac Tamponade, Cardiac Arrest and Guidewire Replacement.

**Adverse Events:**

- Access site complication
- Additional Surgical Procedure
- Air Embolism/Thromboembolism
- Allergic Reaction
- Amputation
- Aorta Complications
- Arteriovenous (AV) Fistula
- Arrhythmia
- Bleeding
- Cardiac and/or Septal Perforation
- Death
- Embolism
- Hematoma
- Hemorrhage
- Hemoglobinuria
- Hypovolemia
- Infection or Sepsis
- MACCE
- Myocardial Ischemia and/or Infarction
- Pericardial Effusion
- Pseudoaneurysm
- Renal Failure or Injury
- Stroke or other Neurologic event
- Tamponade
- Thrombus
- Valve Complications
- Vascular Complication
- Vessel Occlusion
- Vessel Perforation, Dissection, Trauma, or Damage
- Vessel Spasm
- Wire Entrapment/Entanglement
- Foreign Body/Wire Fracture.

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

Rev AA

**Boston  
Scientific**  
Advancing science for life<sup>™</sup>

### Interventional Cardiology

300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

To order product or for more information  
contact customer service at 1.888.272.1001.

© 2016 Boston Scientific Corporation  
or its affiliates. All rights reserved.

SH-322920-AB-US AUG2016