

# ARCHITECT SYSTEM

## AFP Calibrators

### INTENDED USE

The ARCHITECT AFP Calibrators are for the calibration of the ARCHITECT *i* System when used for the quantitative determination of alpha-fetoprotein (AFP) in human serum, plasma and amniotic fluid. Refer to the ARCHITECT AFP reagent package insert and the ARCHITECT System Operations Manual for additional information.

### CONTENTS

6 Bottles (4.0 mL each) of ARCHITECT AFP Calibrators A-F (CAL A - CAL F). Calibrator A contains buffer solution with protein (bovine) stabilizer. Calibrators B-F contain purified AFP (from human cord serum) prepared in buffer solution with protein (bovine) stabilizer. Preservatives: ProClin 300 and ProClin 950.

The calibrators are at the following concentrations:

| Calibrator | Concentration |       |
|------------|---------------|-------|
|            | ng/mL         | IU/mL |
| CAL A      | 0             | 0     |
| CAL B      | 15            | 12.45 |
| CAL C      | 45            | 37.35 |
| CAL D      | 300           | 249   |
| CAL E      | 1500          | 1245  |
| CAL F      | 2000          | 1660  |


### STANDARDIZATION

The ARCHITECT AFP calibrators are manufactured gravimetrically and are referenced to the World Health Organization (WHO) First International Standard 72/225 for Alpha-fetoprotein at each concentration level. The conversion factor is 0.83 International Units per nanogram of AFP.

### PRECAUTIONS

#### IVD

- For In Vitro Diagnostic Use

-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the CONTENTS section of this insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens<sup>1</sup>. Biosafety Level 2<sup>2</sup> or other appropriate biosafety practices<sup>3,4</sup> should be used for materials that contain or are suspected of containing infectious agents.
- Calibrators B-F contain purified AFP from human cord serum tested and found to be nonreactive for HIV 1 & 2, HCV, and HBV.
- WARNING: SENSITIZER** Warning: May cause an allergic reaction.



en

AFP

REF 3P36-01

S3P360

G4-5528/R02

Read Highlighted Changes  
Revised November 2013

- The following warnings and precautions apply to these components
  - Calibrators A-F



**WARNING:** H317 Contains methylisothiazolones. May cause an allergic skin reaction.

#### Prevention

P261 Avoid breathing mist / vapours / spray.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves / protective clothing / eye protection.

#### Response

P302+P352 IF ON SKIN: Wash with plenty of water.  
P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

~~P362+P364~~ Take off contaminated clothing and wash it before reuse.

This material and its container must be disposed of in a safe way.

- Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

### STORAGE

- ARCHITECT AFP Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- ARCHITECT AFP Calibrators must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

2°C  8°C

### QUALITY CONTROL PROCEDURES

A single sample of each control level must be tested to evaluate the assay calibration. For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.

- Ensure that assay control values are within the ranges specified in the control package insert.
- Once an ARCHITECT AFP calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
- A reagent kit with a new lot number is used
  - Controls are out of range

### PREPARATION FOR USE

- ARCHITECT AFP Calibrators must be mixed by gentle inversion before use.
- To perform a calibration, test the calibrators in duplicate. The calibrators should be priority loaded.
- To obtain the recommended volume requirements for the calibrators, hold the bottles vertically and dispense a minimum of 4 drops of each calibrator into each respective sample cup.
- After each use, tightly close the caps and return the calibrators to 2-8°C storage.
- For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.



**BIBLIOGRAPHY**

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline - Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.  
ProClin is property of its respective owner.



Abbott Ireland  
Diagnostics Division  
Finskillin Business Park  
Sligo  
Ireland  
+353-71-9171712



November 2013  
© 2011, 2013 Abbott Laboratories  
Distributed by Abbott Laboratories  
Abbott Park, IL 60064 USA  
and  
ABBOTT 65205 Wiesbaden, Germany

**Key to symbols used**

**GTIN**

Global Trade Item Number

**PRODUCT OF IRELAND**

Product of Ireland

**INFORMATION FOR USA ONLY**

Information needed for  
United States of America only



# ARCHITECT SYSTEM

## AFP Controls

### INTENDED USE

The ARCHITECT AFP Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT *i* System when used for the quantitative determination of alpha-fetoprotein (AFP) in human serum, plasma and amniotic fluid.

Refer to the ARCHITECT AFP reagent package insert and the ARCHITECT System Operations Manual for additional information.

### CONTENTS

3 Bottles (8.0 mL each) of ARCHITECT AFP Controls. Low (**CONTROL L**), Medium (**CONTROL M**), and High (**CONTROL H**) ARCHITECT AFP Controls contain purified AFP (from human cord serum) prepared in buffer solution with protein (bovine) stabilizer. Preservatives: ProClin 300 and ProClin 950.

The controls are at the following concentrations:

| Control          | Target Conc.<br>(ng/mL) | Control Range<br>(ng/mL) | Target Conc.<br>(IU/mL) | Control Range<br>(IU/mL) |
|------------------|-------------------------|--------------------------|-------------------------|--------------------------|
| <b>CONTROL L</b> | 20                      | 13.50 - 26.50            | 16.6                    | 11.21 - 22.00            |
| <b>CONTROL M</b> | 200                     | 136.00 - 266.00          | 166                     | 112.05 - 219.95          |
| <b>CONTROL H</b> | 1000                    | 675.00 - 1325.00         | 830                     | 560.25 - 1099.75         |


Each laboratory should establish its own concentration ranges for new control lots at each control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days. Sources of variation that can be expected should be included in this study in order to be representative of future system performance. These may include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules
- Data points collected at different times of the day

These results should be applied to your laboratory's quality control practices.

### PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens<sup>1</sup>, Biosafety Level 2<sup>2</sup> or other appropriate biosafety practices<sup>3,4</sup> should be used for materials that contain or are suspected of containing infectious agents.
- Low, Medium, and High Controls contain purified AFP from human cord serum tested and found to be nonreactive for HIV 1 & 2, HCV, and HBV.
- **WARNING: SENSITIZER** Warning. May cause an allergic reaction.



en

AFP

REF 3P36-10

C3P360

G4-5558/R03

Read Highlighted Changes  
Revised November 2013

- The following warnings and precautions apply to these components:

- Controls



**WARNING:** Contains methylisothiazolones.  
H317 May cause an allergic skin reaction.

#### Prevention

P261 Avoid breathing mist / vapours / spray.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves / protective clothing / eye protection.

#### Response

P302+P352 IF ON SKIN: Wash with plenty of water.  
P333+P313 If skin irritation or rash occurs: Get medical advice / attention.  
P362+P364 Take off contaminated clothing and wash it before reuse.

This material and its container must be disposed of in a safe way.

- Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

### STORAGE

- ARCHITECT AFP Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- ARCHITECT AFP Controls must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.



### QUALITY CONTROL PROCEDURES

Refer to the ARCHITECT AFP assay reagent package insert and ARCHITECT System Operations Manual for additional information.

The recommended control requirement for the ARCHITECT AFP assay is that a single sample of each control level be tested once every 24 hours each day of use. If your laboratory quality control procedures require more frequent use of controls to verify test results, follow those procedures. Additional controls may be tested in conformance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

The control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and must be retested. Recalibration may be indicated.

### PREPARATION FOR USE

- ARCHITECT AFP Controls must be mixed by gentle inversion before use.
- To obtain the recommended volume requirements for the controls, hold the bottles vertically and dispense a minimum of 4 drops of each control into each respective sample cup.
- After each use, tightly close the caps and return the controls to 2-8°C storage.
- For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.



**BIBLIOGRAPHY**

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline – Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. ProClin is property of its respective owner.



Abbott Ireland  
Diagnostics Division  
Finisklin Business Park  
Sillgo  
Ireland  
+353-71-9171712



November 2013  
© 2011, 2013 Abbott Laboratories



Distributed by Abbott Laboratories  
Abbott Park, IL 60064 USA  
and  
ABBOTT 65205 Wiesbaden, Germany

**Key to symbols used**

**GTIN**

Global Trade Item Number

**PRODUCT OF IRELAND**

Product of Ireland

**INFORMATION FOR USA ONLY**

Information needed for  
United States of America only



# ARCHITECT

## SYSTEM



# en

AFP

**REF** 3P36

**B3P360**

**G6-2601/R06**

Read Highlighted Changes  
Revised October 2015

# AFP

**Customer Service: Contact your local representative or find country-specific contact information on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)**

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

### Key to symbols used

|               |   |                                 |   |
|---------------|---|---------------------------------|---|
| <b>REF</b>    | List Number                               | <b>CONTROL NO.</b>              | Control Number  |
| <b>IVD</b>    | <i>In Vitro</i> Diagnostic Medical Device | <b>REACTION VESSELS</b>         | Reaction Vessels  |
| <b>LOT</b>    | Lot Number                                | <b>REAGENT LOT</b>              | Reagent Lot   |
|               | Expiration Date                           | <b>REPLACEMENT CAPS</b>         | Replacement Caps  |
| <b>SN</b>     | Serial Number                             | <b>SAMPLE CUPS</b>              | Sample Cups   |
| <b>SEPTUM</b> | Septum                                    | <b>WARNING: SENSITIZER</b>      | Warning: May cause an allergic reaction                             |
|               | Store at 2-8°C                            | <b>CONTAINS: AZIDE</b>          | Contains Sodium Azide. Contact with acids liberates very toxic gas. |
|               | Consult instructions for use              | <b>GTIN</b>                     | Global Trade Item Number  |
|               | Manufacturer                              | <b>PRODUCT OF IRELAND</b>       | Product of Ireland  |
|               |   | <b>INFORMATION FOR USA ONLY</b> | Information needed for United States of America only                |

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

**Abbott**



**WARNING:** The concentration of alpha-fetoprotein (AFP) in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the AFP assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining AFP levels serially is changed, additional sequential testing should be carried out. Prior to changing assays, the laboratory MUST:

1. For Cancer Management - Confirm baseline values for patients being serially monitored.
2. For Prenatal Testing - Establish a range of expected values for the new assay based on serum or plasma and amniotic fluid from pregnant women with confirmed gestational age.

Caution: United States Federal Law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by, or on the order of a physician.

## NAME

ARCHITECT AFP

## INTENDED USE

The ARCHITECT AFP assay is a chemiluminescent microparticle immunoassay (CMIA) for the quantitative determination of alpha-fetoprotein (AFP) in:

1. Human serum or plasma to aid in monitoring disease progression during the course of disease and treatment of patients with nonseminomatous testicular cancer.
2. Human serum, plasma and amniotic fluid at 15 to 21 weeks gestation to aid in the detection of fetal open neural tube defects (NTD). Test results when used in conjunction with ultrasonography or amniography are a safe and effective aid in the detection of fetal open NTD.

## SUMMARY AND EXPLANATION OF TEST

The discovery of alpha-fetoprotein (AFP) in fetal serum was first recorded by Bergstrand and Czar in 1956.<sup>1</sup> Alpha-fetoprotein is a single polypeptide chain glycoprotein with a molecular weight of approximately 70,000 daltons. The physicochemical properties and amino acid composition are similar to those of albumin.<sup>2,3</sup> Synthesis of AFP occurs primarily in the liver and yolk sac of the fetus. It is secreted into fetal serum, reaching a peak at about 13 weeks gestation and gradually declining thereafter. Elevated serum AFP levels subsequently reappear during pregnancy and in conjunction with several malignant diseases.

### Cancer Management

Alpha-fetoprotein (AFP) was first described as a human tumor-associated protein in 1964 by Tatarinov.<sup>4</sup> Since then, it has been shown that elevation of serum AFP above values typically found in healthy individuals occurs in several malignant diseases,<sup>5-8</sup> most notably nonseminomatous testicular cancer and primary hepatocellular carcinoma. In the case of nonseminomatous testicular cancer, a direct relationship has been observed between the incidence of elevated AFP levels and the stage of disease.<sup>9,10</sup> Elevated AFP levels have also been observed in patients diagnosed as having seminoma with nonseminomatous elements but have not been observed in patients with pure seminoma.<sup>9,11,12</sup> Human chorionic gonadotropin (hCG) and AFP are also important prognostic indicators of survival rate among patients with advanced nonseminomatous germ cell testicular tumors.<sup>13</sup>

The usefulness of AFP measurements in the management of patients with nonseminomatous testicular cancers has been well documented.<sup>7,11,14</sup> For patients in clinical remission following treatment, AFP levels generally decrease.<sup>11</sup> Post-operative AFP values which fail to return to normal strongly suggest the presence of residual tumor.<sup>8,7,11</sup> Tumor recurrence is often accompanied by a rise in AFP before progressive disease is clinically evident.<sup>7,9</sup>

Greater than 70% of patients with primary hepatocellular carcinoma have been reported to have elevated levels of serum AFP.<sup>5,6,15</sup> Elevated AFP levels have occasionally been found in association with gastrointestinal tract cancers with and without liver metastases<sup>16</sup> and only rarely in other malignancies.<sup>5,6</sup> Serum AFP has been found to be elevated during pregnancy, in diseases such as ataxia telangiectasia, hereditary tyrosinemia, teratocarcinoma and in benign hepatic conditions such as acute viral hepatitis, chronic active hepatitis and cirrhosis.<sup>8,15,17</sup> Elevation of serum AFP in benign hepatic diseases is usually transient.<sup>5</sup>

AFP testing is not recommended as a screening procedure to detect cancer in the general population.

## Prenatal Testing

Many studies have confirmed the utility of AFP in the early detection of fetal open neural tube defects (NTD),<sup>18-20</sup> in the US. NTD, primarily anencephaly and spina bifida, occur at the rate of between 1 and 2 per 1000 live births and are among the most common major congenital malformations.<sup>21,31</sup> The incidence of NTD varies geographically and across racial groups.<sup>22-26</sup>

Anencephaly is incompatible with life and accounts for one-third to one-half of all NTD. Open spina bifida can vary widely in severity.

Reports from the scientific literature suggest additional factors to be considered when assessing the risk of an NTD being present.<sup>27-28</sup> One is the effect of maternal weight. Maternal blood volume, as reflected by maternal weight, has been reported to affect maternal serum AFP (MSAFP) concentration in maternal circulation; the higher the maternal weight, the lower the MSAFP concentration.<sup>29-29</sup> Another factor to consider is maternal diabetes. Insulin dependent diabetic women reportedly have MSAFP levels significantly lower than non-diabetic women and an increased incidence of NTD.<sup>27,28,30</sup> Maternal serum AFP levels in the black population average about 10% higher than MSAFP values in the non-black population. An adjustment factor or use of an appropriate normative data base have been suggested in the literature.<sup>25,26</sup>

Amniotic fluid AFP (AFAP) levels peak at about 13 weeks gestation after which they rapidly decline until about 22 weeks gestation and then gradually decline until term. Transfer of AFP into maternal circulation is accomplished primarily through diffusion across the placenta.<sup>31</sup> If the fetus has an open neural tube defect, AFP is thought to leak directly into the amniotic fluid (AF) causing unexpectedly high levels of AFAP. Subsequently, the AFAP reaches the maternal circulation, thus producing abnormally elevated levels of MSAFP. Certain fetal abnormalities such as congenital renal disease and esophageal atresia also show AFAP elevations.<sup>32,33</sup> Other fetal distress situations such as omphalocele or gastroschisis, defective kidneys, threatened abortion, prematurity and sometimes fetal demise<sup>34-37</sup> may exhibit abnormally high levels of MSAFP. Increased MSAFP values are also seen in multiple pregnancies<sup>38</sup> and in normal singleton pregnancies in which the gestational age has been underestimated. Low MSAFP values have been associated with molar pregnancy, missed abortion, pseudocyesis, overestimated gestational age and Down Syndrome.<sup>29,39</sup>

In a report on over 18,000 pregnancies, the U.K. Collaborative Study has established multiples of the median (MoM) as the preferred way to express AFP results.<sup>18</sup> The median AFP value for each gestational week is first determined; then individual AFP levels are reported as multiples of this value. This method of expression facilitates comparison of AFP test results across gestational weeks and between laboratories.

AFP testing during pregnancy is recommended as an effective way to determine those women potentially at risk of carrying a fetus affected with an open NTD. Used in conjunction with other confirmatory procedures such as ultrasonography or amniography, measurement of AFP serves as an important tool in the care and management of these patients.

## BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT AFP assay is a two-step immunoassay for the quantitative measurement of AFP in human serum, plasma and amniotic fluid using CMIA technology, with flexible assay protocols, referred to as Chemiflex.

In the first step, sample and anti-AFP coated paramagnetic microparticles are combined. AFP present in the sample binds to the anti-AFP coated microparticles. After washing, anti-AFP acridinium-labeled conjugate is added to create a reaction mixture in the second step. Following another wash cycle, pre-trigger and trigger solutions are added to the reaction mixture. The resulting chemiluminescent reaction is measured as relative light units (RLUs). A direct relationship exists between the amount of AFP in the sample and the RLUs detected by the ARCHITECT / System optics. For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

## REAGENTS

### Reagent Kit, 100 Tests/500 Tests

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT / Systems. Please contact your local distributor.

### ARCHITECT AFP Reagent Kit (3P36)

- **MICROPARTICLES** 1 or 4 Bottle(s) (6.6 mL/27.0 mL) Anti-AFP (mouse, monoclonal) coated microparticles in MES buffer with protein (bovine) stabilizer. Minimum concentration: 0.1% solids. Preservative: ProClin 300.
- **CONJUGATE** 1 or 4 Bottle(s) (5.9 mL/26.3 mL) Anti-AFP (mouse, monoclonal) acridinium-labeled conjugate in MES buffer with protein (bovine) stabilizer. Minimum concentration: 400 ng/mL. Preservatives: antimicrobial agents and sodium azide.



## Assay Diluent

### ARCHITECT i Multi-Assay Manual Diluent (7DRS-50)

- **MULTI-ASSAY MANUAL DILUENT** 1 Bottle (100 mL) ARCHITECT i Multi-Assay Manual Diluent containing phosphate buffered saline solution. Preservative: antimicrobial agent.

## Other Reagents

### ARCHITECT i Pre-Trigger Solution

- **PRE-TRIGGER SOLUTION** Pre-trigger solution containing 1.32% (w/v) hydrogen peroxide.

### ARCHITECT i Trigger Solution

- **TRIGGER SOLUTION** Trigger solution containing 0.35 N sodium hydroxide.

### ARCHITECT i Wash Buffer

- **WASH BUFFER** Wash buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

## WARNINGS AND PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use
- Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

## Safety Precautions

- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and be handled in accordance with the OSHA Standard on Bloodborne Pathogens<sup>40</sup>, Biosafety Level 2<sup>41</sup> or other appropriate biosafety practices<sup>42,43</sup> should be used for materials that contain or are suspected of containing infectious agents.

- The following warnings and precautions apply to this component:

- Conjugate

Contains sodium azide.

**EUH032** Contact with acids liberates very toxic gas.

This material and its container must be disposed of in a safe way.

- The following warnings and precautions apply to this component:

- Microparticles



**WARNING:** Contains methylisothiazolones  
H317 May cause an allergic skin reaction.

### Prevention

P261 Avoid breathing mist/vapours/spray.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves/protective clothing/eye protection.

### Response

P302+P352 IF ON SKIN: Wash with plenty of water.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
P362+P364 Take off contaminated clothing and wash it before reuse.

This material and its container must be disposed of in a safe way.

- Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

- For information on the safe disposal of sodium azide and a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

## Handling Precautions

- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between reagent kits.
- Before loading the ARCHITECT AFP Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the PROCEDURE, Assay Procedure section of this package insert.
- Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the Instructions in this package insert.
  - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
  - Once a septum has been placed on the reagent bottle, do not invert the bottle as this will result in reagent leakage and may compromise assay results.

- Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

- For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

## Storage Instructions

- $2^{\circ}\text{C}$   $-8^{\circ}\text{C}$  The ARCHITECT AFP Reagent Kit must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.
- When stored and handled as directed, the reagents are stable until the expiration date.
- The ARCHITECT AFP Reagent Kit may be stored on board the ARCHITECT i System for a maximum of 30 days. After 30 days, the reagent kit must be discarded. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.
- Reagents may be stored on or off the ARCHITECT i System. If reagents are removed from the system, store them at 2-8°C (with septuma and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded. For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

## Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

- The ARCHITECT AFP assay is designed for use on the ARCHITECT i System.
- The ARCHITECT AFP assay file (assay number 003) must be installed on the ARCHITECT i System before performing the assay.
- For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.
- For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.
- For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.
- The default result unit for the ARCHITECT AFP assay is ng/mL. An alternate result unit, IU/mL, may be selected for reporting results by editing assay parameter "Result concentration units" to IU/mL. The conversion factor used by the system is 0.83 as follows:
  - (Concentration in ng/mL) x (0.83) = IU/mL

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

- Serum, plasma, or amniotic fluid specimens may be used with the ARCHITECT AFP assay.
- The specimen collection tubes listed below were verified for use for serum and plasma with the ARCHITECT AFP assay. Other specimen collection tubes have not been tested with this assay.
  - Human serum, plastic (including serum collected in plastic serum separator tubes).
  - Human plasma collected in:
    - sodium heparin, plastic
    - dipotassium EDTA, plastic
    - lithium heparin, plastic
    - sodium EDTA, glass
- Serum or plasma specimens should be collected aseptically in such a way as to avoid hemolysis.
- For maternal serum or plasma analysis, the blood specimen should be collected prior to the initiation of amniocentesis. It has been demonstrated that increased levels of AFP may occur in maternal serum or plasma following amniocentesis.<sup>44</sup>
- When serial specimens are being evaluated, the same type of specimen should be used.



- Amniotic fluid should be collected aseptically with appropriate precautions relative to both fetal and maternal safety by appropriately trained personnel. Visibly bloodstained specimens should be examined for the presence of fetal blood cells by using the Kleihauer-Betke technique and/or fetal hemoglobin by electrophoresis, immunoelectrophoresis, or other available techniques. Amniotic fluid specimens contaminated with fetal blood may exhibit abnormally high AFP values which may lead to misinterpretation of test results.
- Performance has not been established for the use of cadaveric specimens or body fluids other than human serum, plasma, or amniotic fluid.
- The ARCHITECT *i* System does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the ARCHITECT AFP assay.

#### Specimen Conditions

- Do not use specimens with the following conditions:
  - heat-inactivated
  - pooled
  - grossly hemolyzed
  - obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- Use caution when handling patient specimens to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.
- For optimal results, inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

#### Preparation for Analysis

- Follow the tube manufacturer's processing instructions for serum and plasma collection tubes. Gravity separation is not sufficient for specimen preparation.
- Prepare frozen specimens as follows:
  - Frozen specimens must be completely thawed before mixing.
  - Mix thawed specimens thoroughly by inverting 10 times or by low speed vortexing. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous. If samples are not mixed thoroughly, inconsistent results may be obtained.
  - Centrifuge mixed specimens as described below.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and centrifuged at  $\geq 10,000$  RCF (Relative Centrifugal Force) for 10 minutes before testing if
  - they contain fibrin, red blood cells, or other particulate matter or
  - they were frozen and thawed.
- Centrifuged specimens with a lipid layer on the top must be transferred to a sample cup or secondary tube. Care must be taken to transfer only the clarified specimen without the lipemic material.
- Transfer clarified specimen to a sample cup or secondary tube for testing.

#### Storage

##### Serum or Plasma

- Specimens may be stored on or off the clot, red blood cells, or separator gel for
  - up to 3 days at room temperature or
  - up to 7 days at 2-8°C.
- If testing will be delayed more than 7 days, remove serum or plasma from the clot, red blood cells, or separator gel and store at -20°C or colder.
- Avoid more than 5 freeze/thaw cycles.

##### Amniotic Fluid

- Specimens may be stored for
  - up to 2 days at room temperature or
  - up to 5 days at 2-8°C.
- If testing is delayed more than 5 days, store at -20°C or colder.
- Avoid more than 3 freeze/thaw cycles.

#### Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

#### PROCEDURE

##### Materials Provided

- 3P36 ARCHITECT AFP Reagent Kit

##### Materials Required but not Provided

- ARCHITECT *i* System
- ARCHITECT AFP Assay file, may be obtained from:
  - ARCHITECT *i* System e-Assay CD-ROM found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)
  - ARCHITECT *i* System Assay CD-ROM
- 3P36-01 ARCHITECT AFP Calibrators
- 3P36-10 ARCHITECT AFP Controls or other control material
- 7D82-50 ARCHITECT *i* Multi-Assay Manual Diluent
- ARCHITECT *i* PRE-TRIGGER SOLUTION
- ARCHITECT *i* TRIGGER SOLUTION
- ARCHITECT *i* WASH BUFFER
- ARCHITECT *i* REACTION VESSELS
- ARCHITECT *i* SAMPLE CUPS
- ARCHITECT *i* SEPTUM
- ARCHITECT *i* REPLACEMENT CAPS
- Pipettes or pipette tips (optional) to deliver the specified volumes.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

##### Assay Procedure

- Before loading the ARCHITECT AFP Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
  - Invert the microparticle bottle 30 times.
  - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
  - If the microparticles do not resuspend, DO NOT USE. Contact your Abbott representative.
  - Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the Handling Precautions section of this package insert.
- Load the ARCHITECT AFP Reagent Kit on the ARCHITECT *i* System.
  - Verify that all necessary reagents are present.
  - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
  - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
  - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- The minimum sample cup volume is calculated by the system and is printed on the Orderlist report. No more than 10 replicates may be sampled from the same sample cup. To minimize the effects of evaporation, verify adequate sample cup volume is present before running the test.
  - Priority: 75  $\mu$ L for the first AFP test plus 25  $\mu$ L for each additional AFP test from the same sample cup.
  - $\leq 3$  hours on-board: 150  $\mu$ L for the first AFP test plus 25  $\mu$ L for each additional AFP test from the same sample cup.
  - $> 3$  hours on-board: replace with a fresh sample (patient specimens, controls and calibrators).
  - If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.





- Prepare calibrators and controls.
  - Mix the ARCHITECT AFP Calibrators and Controls by gentle inversion before use.
  - To obtain the recommended volume requirements for the ARCHITECT AFP Calibrators and Controls, hold the bottles vertically and dispense 4 drops of each calibrator or control into each respective sample cup.
- Load samples
  - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

#### Specimen Dilution Procedures

- Specimens with an AFP concentration greater than 2000 ng/mL will be flagged as "> 2000.00 ng/mL" and may be diluted using either the Automated Dilution Protocol or the Manual Dilution Procedure.

#### Automated Dilution Protocol for Serum or Plasma Specimens

- If the Automated Dilution Protocol is chosen, use the 1:10 dilution for serum or plasma. The system automatically calculates the concentration of the sample before dilution and reports the result.
- Dilutions other than the automated 1:10 serum or plasma dilution should be done manually.

#### Automated Dilution Protocol for Amniotic Fluid Specimens

**NOTE: Amniotic fluid specimens must be diluted.**

- If the Automated Dilution Protocol is chosen, amniotic fluid **MUST ONLY USE** the 1:40 dilution. The system automatically calculates the concentration of the sample before dilution and reports the result.
- Dilutions other than the automated 1:40 amniotic fluid dilution should be done manually.

#### Manual Dilution Procedure for All Specimen Types

**NOTE:** The ARCHITECT *i* Multi-Assay Manual Diluent (7D82-50) must be used when performing the manual dilution procedure.

- For a 1:20 dilution, add 50  $\mu$ L of the patient specimen to 950  $\mu$ L of the ARCHITECT *i* Multi-Assay Manual Diluent (7D82-50). For a 1:101 dilution, add 10  $\mu$ L of the patient specimen to 1 mL of the ARCHITECT *i* Multi-Assay Manual Diluent (7D82-50).
- The operator must enter the dilution factor in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution.
- For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

#### Calibration

- To perform an ARCHITECT AFP calibration, test calibrators A, B, C, D, E and F in replicates of 2. The calibrators should be priority loaded.
- Calibration Range: 0 to 2000 ng/mL.
- A single sample of each control level must be tested to evaluate the assay calibration.
  - Order controls as described in the Assay Procedure section.
  - Ensure that assay control values are within the ranges specified in the control package insert.
- Once an ARCHITECT AFP calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
  - A reagent kit with a new lot number is used.
  - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

#### QUALITY CONTROL PROCEDURES

The recommended control requirement for the ARCHITECT AFP assay is that a single sample of each control be tested once every 24 hours each day of use. If your laboratory quality control procedures require more frequent use of controls to verify test results, follow those procedures. Additional controls may be tested in conformance with local, state and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

Each laboratory should establish control means and ranges to monitor the acceptable performance of the assay. If a control is out of its specified range, the associated test results are invalid and samples must be retested. Recalibration may be indicated.

After the median AFP values have been established for maternal serum/plasma and amniotic fluid, the control means should remain within acceptable limits set by the laboratory. The acceptability of each calibration should be closely monitored by the controls using guidance from CLSI (LA25-A2),<sup>45</sup> the National Academy of Clinical Biochemistry (NACB),<sup>46</sup> and/or the laboratories internal operating procedures to detect any shifts which may require assay recalibration or re-evaluation of the maternal serum/plasma and amniotic fluid medians.

#### Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B. The ARCHITECT AFP assay belongs to method group 1.

#### RESULTS

The ARCHITECT AFP assay uses a 4 Parameter Logistic Curve fit (4PLC, Y-weighted) data reduction method to generate a calibration curve.

#### Alternate Result Unit

- The default result unit for the ARCHITECT AFP assay is ng/mL. When the alternate result unit, IU/mL, is selected, the conversion factor used by the system is 0.83.
  - Conversion Formula: (Concentration in ng/mL) x (0.83) = IU/mL
- To convert amniotic fluid values to  $\mu$ g/mL, divide the reported AFP concentration (ng/mL) by 1000, as this calculation is not performed automatically.

#### Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

#### Measuring Interval (Reportable Range)

Measuring interval is defined as the range of values in ng/mL which meets the limits of acceptable performance for both imprecision and bias for an undiluted sample. For the studies described in this package insert, the range was 2.00 ng/mL (Limit of Quantitation - LoQ) to 2000.00 ng/mL.

When using the 1:10 automated dilution protocol, the assay can report values up to 20,000.00 ng/mL. When using the 1:40 automated dilution protocol, the assay can report values up to 80,000.00 ng/mL.

#### LIMITATIONS OF THE PROCEDURE

- If the AFP results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other tests, clinical impressions, etc.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA).<sup>47,48</sup> Specimens containing HAMA may produce anomalous values when tested with assay kits such as ARCHITECT AFP that employ mouse monoclonal antibodies.<sup>47</sup>
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays.<sup>49</sup> Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous results may be observed. Additional information may be required for diagnosis.
- Although the ARCHITECT AFP assay is specifically designed to minimize the effects of HAMA and heterophilic antibodies, assay results that are not consistent with other clinical observations may require additional information for diagnosis.
- The ARCHITECT AFP assay is a valuable aid in the management of nonseminomatous testicular cancer patients when used in conjunction with information available from the clinical evaluation and other diagnostic procedures. Increased serum AFP concentrations have also been observed in ataxia telangiectasia, hereditary tyrosinemia, primary hepatocellular carcinoma, teratocarcinoma, gastrointestinal tract cancers with and without liver metastases and in benign hepatic conditions such as acute viral hepatitis, chronic active hepatitis and cirrhosis.
- The ARCHITECT AFP assay should not be used as a cancer screening test.
- Valid measurements of AFP in maternal serum or plasma CANNOT be made after amniocentesis; therefore, maternal serum or plasma specimens MUST be drawn PRIOR to amniocentesis. For further information, refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section in this package insert.



- Amniotic fluid specimens contaminated with fetal blood may exhibit abnormally high AFP values which may lead to misinterpretation of test results. Visibly bloodstained specimens should be examined for the presence of fetal blood cells by using the Kleihauer-Betke technique and/or fetal hemoglobin by electrophoresis, immunoelectrophoresis, or other available techniques. For further information, refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section in this package insert.
- The reliability of MSAFP evaluation in prenatal testing is dependent upon the accurate determination of gestational age. An inaccurate estimation of gestational age may result in an inaccurate estimation of risk of NTD. When gestational age is uncertain, a reliable ultrasound examination is important.
- While elevated levels of MSAFP indicate increased risk of NTD, they are not diagnostic. Increased serum AFP concentrations have been seen in some cancers and some nonmalignant diseases as described above and, thus, may be indicative of maternal conditions. Other conditions including placental malformations, open fetal malformations such as omphalocele or gastroschisis (ventral wall defects), fetal kidney abnormalities, threatened or imminent abortion and fetal demise are associated with elevated levels of MSAFP. Elevated MSAFP levels have also been associated with premature deliveries and low birth weights and have been seen in multiple births. Rarely, singleton, viable and unaffected pregnancies may exhibit elevated MSAFP levels. Confirmatory testing, such as amniocentesis for AFAFP evaluation, high resolution ultrasonography or amniography is an essential part of the AFP testing process.
- Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section of this package insert for specimen limitations.

#### EXPECTED VALUES

Data in the EXPECTED VALUES section were generated using the ARCHITECT i 2000/i 2000SR Systems.

The distribution of ARCHITECT AFP values was determined in 400 specimens from apparently healthy individuals (200 males and 200 females), in 238 patients with nonmalignant diseases and in 224 patients diagnosed with malignant diseases. The data are summarized in the following tables.

Distribution of ARCHITECT AFP Values

| Group/Category              | n   | Distribution of Values (%)<br>by AFP Concentration Range in ng/mL |               |                  |            |             |              |       |
|-----------------------------|-----|---|---------------|------------------|------------|-------------|--------------|-------|
|                             |     | 0 - 8.78  | >8.78 - 15.00 | >15 - 200        | >200 - 500 | >500 - 1000 | >1000 - 2000 | >2000 |
| Apparently Healthy Subjects | 400 | 97.5  | 2.0           | 0.5 <sup>a</sup> | 0.0        | 0.0         | 0.0          | 0.0   |

<sup>a</sup> These 2 samples had AFP concentrations of 25.16 and 27.81 ng/mL.

The observed nonparametric central 95% of the 400 apparently healthy individuals ranged from 0.89 to 8.78 ng/mL. It is recommended that each laboratory establish its own expected reference range for the population of interest.

Distribution of ARCHITECT AFP Values

| Group/Category                       | n   | Distribution of Values (%)<br>by AFP Concentration Range in ng/mL |               |           |            |             |              |       |
|--------------------------------------|-----|---|---------------|-----------|------------|-------------|--------------|-------|
|                                      |     | 0 - 8.78  | >8.78 - 15.00 | >15 - 200 | >200 - 500 | >500 - 1000 | >1000 - 2000 | >2000 |
| <b>Nonmalignant Disease</b>          |     |   |               |           |            |             |              |       |
| Cirrhosis                            | 49  | 98.0  | 2.0           | 0.0       | 0.0        | 0.0         | 0.0          | 0.0   |
| Genitourinary                        | 26  | 100.0   | 0.0           | 0.0       | 0.0        | 0.0         | 0.0          | 0.0   |
| Hepatitis                            | 149 | 90.6  | 6.7           | 2.0       | 0.0        | 0.0         | 0.0          | 0.7   |
| Pancreatitis                         | 14  | 92.9  | 7.1           | 0.0       | 0.0        | 0.0         | 0.0          | 0.0   |
| <b>Malignant Disease<sup>a</sup></b> |     |   |               |           |            |             |              |       |
| Gastrointestinal                     | 64  | 98.4  | 1.6           | 0.0       | 0.0        | 0.0         | 0.0          | 0.0   |
| Hepatocellular                       | 29  | 69.0  | 0.0           | 17.2      | 0.0        | 0.0         | 0.0          | 13.8  |
| Pancreatic                           | 34  | 65.3  | 0.0           | 5.9       | 2.9        | 0.0         | 2.9          | 2.9   |
| Nonseminoma Testicular               | 72  | 87.5  | 1.4           | 9.7       | 1.4        | 0.0         | 0.0          | 0.0   |
| Seminoma Testicular                  | 25  | 100.0   | 0.0           | 0.0       | 0.0        | 0.0         | 0.0          | 0.0   |

<sup>a</sup> The nonseminoma testicular samples were from treated patients. The disease status was unknown for specimens from the other malignant diseases.

#### AFP Values in Maternal Serum and Amniotic Fluid

Due to variations in populations at different locations, it is important for each laboratory to establish its own medians for each gestational week for maternal serum and amniotic fluid. AFP values should be expressed as multiples of the median (MoM) as shown in the calculation below:

$$\text{MoM} = \frac{\text{Specimen AFP Concentration}}{\text{Median AFP Concentration for Gestational Week}}$$

Each laboratory should attempt to gather approximately 100 or more specimens for each gestational week in order to arrive at median values<sup>60</sup>, then utilize a cutoff value (MoM) which most closely suits its needs for specificity and sensitivity.

A total of 685 maternal serum and 687 amniotic fluid specimens from unaffected or low-risk singleton pregnancies were evaluated with the ARCHITECT AFP assay. The AFP values expressed as the regressed medians and multiples of the regressed medians (MoM) for gestational weeks 15 to 21 are shown in the following tables.

Maternal Serum AFP

| Gestational Week | n   | Regressed Medians <sup>a</sup> (ng/mL) | Multiples of Regressed Medians (ng/mL) |        |        |
|------------------|-----|--|--|--------|--------|
|                  |     |  | 2.0                                    | 2.5    | 3.0    |
| 15               | 101 | 32.17                                  | 64.35                                  | 80.44  | 96.52  |
| 16               | 95  | 36.86                                  | 73.73                                  | 92.16  | 110.59 |
| 17               | 102 | 42.24                                  | 84.48                                  | 105.60 | 126.72 |
| 18               | 103 | 48.40                                  | 96.79                                  | 120.99 | 145.19 |
| 19               | 101 | 55.45                                  | 110.90                                 | 138.63 | 166.35 |
| 20               | 106 | 63.53                                  | 127.07                                 | 158.84 | 190.60 |
| 21               | 77  | 72.80                                  | 145.59                                 | 181.99 | 218.39 |

<sup>a</sup> The regressed median values were determined using a weighted log-linear regression analysis.<sup>16</sup>

Amniotic Fluid AFP

| Gestational Week | n   | Regressed Medians <sup>a</sup> (µg/mL) | Multiples of Regressed Medians (µg/mL) |       |       |
|------------------|-----|--|--|-------|-------|
|                  |     |  | 2.0                                    | 2.5   | 3.0   |
| 15               | 104 | 16.41                                  | 32.82                                  | 41.02 | 49.22 |
| 16               | 108 | 13.38                                  | 26.76                                  | 33.45 | 40.14 |
| 17               | 105 | 10.91                                  | 21.82                                  | 27.27 | 32.72 |
| 18               | 109 | 8.89                                   | 17.79                                  | 22.23 | 26.68 |
| 19               | 102 | 7.25                                   | 14.50                                  | 18.13 | 21.75 |
| 20               | 97  | 5.91                                   | 11.83                                  | 14.78 | 17.74 |
| 21               | 62  | 4.82                                   | 9.64                                   | 12.05 | 14.46 |

<sup>a</sup> The regressed median values were determined using a weighted log-linear regression analysis.<sup>16</sup>

Note: AFP values were assigned on the basis of completed gestational weeks. For example, a specimen obtained on gestational day 132 (week 18, day 6) was assigned week 18, because the gestation had only completed 18 gestational weeks, plus 6 days.

#### Clinical Specificity and Sensitivity

The specificity and sensitivity estimates (and associated 95% confidence intervals) of the ARCHITECT AFP assay were determined for maternal serum and amniotic fluid at various multiples of the median (MoM). As defined here, specificity is the probability that the test will be negative in the absence of open NTD and sensitivity is the probability that the test will be positive in the presence of open NTD.

The specificity table represents data gathered on unaffected singleton pregnancies from 15 to 21 weeks gestation using the ARCHITECT AFP assay. The data are summarized in the following table.

| Specimen Type  | n   | Specificity (95% Confidence Interval) by Multiples of the Median (MoM) |                            |                            |
|----------------|-----|--|----------------------------|----------------------------|
|                |     | 2.0  | 2.5                        | 3.0                        |
| Maternal Serum | 682 | 95.45%<br>(93.61%, 96.89%)   | 98.24%<br>(96.95%, 99.09%) | 99.71%<br>(98.94%, 99.96%) |
| Amniotic Fluid | 222 | 98.65%<br>(96.10%, 99.72%)   | 99.10%<br>(96.78%, 99.69%) | 99.55%<br>(97.52%, 99.92%) |



The sensitivity table represents data gathered on confirmed affected, singleton pregnancies using the ARCHITECT AFP assay. The data are summarized in the following table.

| Specimen Type  | n  | Sensitivity (95% Confidence Interval) by Multiples of the Median (MoM) |                            |                            |
|----------------|----|--|----------------------------|----------------------------|
|                |    | 2.0  | 2.5                        | 3.0                        |
| Maternal Serum | 21 | 95.24%<br>(76.18%, 99.88%)   | 80.95%<br>(58.03%, 94.55%) | 71.43%<br>(47.82%, 88.72%) |
| Amniotic Fluid | 19 | 100.00%<br>(82.35%, 100.00%)   | 94.74%<br>(73.91%, 99.97%) | 94.74%<br>(73.97%, 99.87%) |

#### AFP Serial Monitoring Performance

In conjunction with physical examination, histology/pathology, and other clinical evaluation procedures, changes observed in serial AFP assay values should be evaluated when monitoring nonseminomatous testicular cancer.

The reference change value (RCV) was used to determine if a significant change in AFP occurred.<sup>51</sup> For this calculation, the RCV for each assay (ARCHITECT AFP and the comparator) was derived by taking into account the published biological variation for AFP<sup>52</sup> and the total imprecision of the specific assay. The RCV for the ARCHITECT AFP method was calculated to be 39.22% and that of the comparator to be 39.98%. A minimum of 3 serial samples were obtained from each of 72 subjects and were analyzed to determine the change in AFP concentration per sequential pair (n=207). The data are summarized in the following tables.

| % Change in AFP       | Change in Disease Status |               |               |               |                 |
|-----------------------|--------------------------|---------------|---------------|---------------|-----------------|
|                       | R<br>n (%)               | S<br>n (%)    | NED<br>n (%)  | P<br>n (%)    | Total<br>n (%)  |
| > RCV Increase        | 7<br>(3.38)              | 3<br>(1.45)   | 9<br>(4.35)   | 8<br>(3.86)   | 27<br>(13.04)   |
| No Significant Change | 20<br>(9.66)             | 38<br>(18.35) | 70<br>(33.82) | 18<br>(8.70)  | 146<br>(70.53)  |
| > RCV Decrease        | 8<br>(3.86)              | 12<br>(5.80)  | 5<br>(2.42)   | 9<br>(4.35)   | 34<br>(16.43)   |
| Total                 | 35<br>(16.91)            | 53<br>(25.60) | 84<br>(40.58) | 35<br>(16.91) | 207<br>(100.00) |

R = Responding; S = Stable; NED = No Evidence of Disease, P = Progressing.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

All performance studies were conducted using the ARCHITECT i 2000/i 2000<sub>SR</sub> Systems. Assay results obtained in individual laboratories may vary from data presented.

##### Precision

The ARCHITECT AFP assay is designed to have an imprecision of  $\leq 7.5\%$  within-laboratory (Total) %CV for samples between 10 and 2000 ng/mL and an SD of  $\leq 0.75$  for samples less than 10 ng/mL down to the LoQ (i.e., 2.0 ng/mL).

##### System Reproducibility

A 5-day precision study was performed for the ARCHITECT AFP assay based on guidance from the National Committee for Clinical Laboratory Standards (NCCLS) document EP5-A2<sup>53</sup> and the Clinical and Laboratory Standards Institute (CLSI) document EP15-A2.<sup>54</sup> Testing was conducted at 3 clinical sites using 3 lots each of ARCHITECT AFP Reagents, Calibrators and Controls and 1 ARCHITECT i 2000/i 2000<sub>SR</sub> instrument per site. Three controls and 5 human serum panels were assayed in replicates of 4 at 2 separate times of day for 5 days. The results are summarized in the following table.

| Sample         | n   | Grand Mean (ng/mL) | Within-Run |     | Within-Day |     | Within-Laboratory Precision (Total) <sup>a</sup> |     | Precision with Additional Component of Between-Site |     | Precision with Additional Component of Between-Lot |     | Precision with Additional Components of Site and Lot (Overall) <sup>b</sup> |     |
|----------------|-----|--------------------|------------|-----|------------|-----|--|-----|---|-----|--|-----|---|-----|
|                |     |                    | SD         | %CV | SD         | %CV | SD   | %CV | SD  | %CV | SD   | %CV | SD  | %CV |
| Low Control    | 360 | 19.46              | 0.885      | 3.5 | 0.747      | 3.8 | 0.747  | 3.8 | 0.797   | 4.1 | 0.813  | 4.2 | 0.830   | 4.3 |
| Medium Control | 360 | 203.65             | 7.502      | 3.7 | 7.987      | 3.9 | 7.987  | 3.9 | 13.984  | 6.0 | 17.293   | 8.5 | 17.293  | 8.5 |
| High Control   | 360 | 973.10             | 45.461     | 4.7 | 46.393     | 4.8 | 46.450   | 4.8 | 63.066  | 6.5 | 65.166   | 6.7 | 65.891  | 6.8 |
| Panel 1        | 360 | 2.95               | 0.114      | 3.9 | 0.127      | 4.3 | 0.128  | 4.4 | 0.239   | 8.1 | 0.278  | 9.4 | 0.278   | 9.4 |
| Panel 2        | 360 | 9.47               | 0.355      | 3.7 | 0.359      | 3.8 | 0.377  | 4.0 | 0.506   | 5.3 | 0.486  | 5.1 | 0.527   | 5.6 |
| Panel 3        | 360 | 591.26             | 25.303     | 4.3 | 25.620     | 4.3 | 26.503   | 4.5 | 42.717  | 7.2 | 50.677   | 8.6 | 50.677  | 8.6 |
| Panel 4        | 360 | 1511.80            | 70.406     | 4.7 | 75.992     | 5.0 | 81.050   | 5.4 | 90.629  | 6.0 | 90.629   | 6.0 | 90.629  | 6.0 |
| Panel 5        | 360 | 1743.35            | 79.509     | 4.6 | 86.121     | 4.9 | 89.956   | 5.2 | 100.632   | 5.8 | 100.632  | 5.8 | 100.632   | 5.8 |

<sup>a</sup> Within-Laboratory (Total) variability contains within-run, within-day and between-day variance components.

<sup>b</sup> Overall variability contains within-run, within-day, between-day, between-lot, between-site and lot-site interaction variance components.

##### Within-Laboratory Precision

A study was performed based on guidance from the NCCLS document EP5-A2.<sup>53</sup> Testing was conducted using 3 lots of ARCHITECT AFP Reagents and Calibrators, 1 lot of ARCHITECT AFP Controls and 4 instruments. Three controls and 5 human serum panels were assayed in a minimum of 4 replicates at 2 separate times per day for 20 different days. Each reagent lot used a single calibration curve throughout the study. The data are summarized in the following table.

| % Change in AFP        | Change in Disease Status |                |               |
|------------------------|--------------------------|----------------|---------------|
|                        | Progression              | No Progression | Total         |
| > 39.22% Increase      | 8 (A)                    | 19 (B)         | 27 (A+B)      |
| $\leq$ 39.22% Increase | 27 (C)                   | 153 (D)        | 180 (C+D)     |
| Total                  | 35 (A+C)                 | 172 (B+D)      | 207 (A+B+C+D) |

Specificity =  $D / (B+D) \times 100\% = 88.95\%$ ; 95% CI = 84.35% to 93.55%

Sensitivity =  $A / (A+C) \times 100\% = 22.86\%$ ; 95% CI = 9.38% to 40.00%

Negative Predictive Value =

$D / (C+D) \times 100\% = 85.00\%$ ; 95% CI = 78.34% to 90.68%

Positive Predictive Value =

$A / (A+B) \times 100\% = 29.63\%$ ; 95% CI = 12.00% to 52.17%

In addition, samples were analyzed on a per subject basis. Efficacy is demonstrated when the sum of sensitivity and specificity are greater than one. In this study, the RCV efficacy for monitoring testicular cancer was determined to be 1.12 with a 95% CI of 0.98 to 1.25.

The change in AFP concentration results available on both the ARCHITECT AFP assay and the comparator AFP assay were analyzed for agreement using their respective RCV.

| ARCHITECT AFP          | Comparator AFP    |                        | Total         |
|------------------------|-------------------|------------------------|---------------|
|                        | > 39.98% Increase | $\leq$ 39.98% Increase |               |
| > 39.22% Increase      | 18 (A)            | 9 (B)                  | 27 (A+B)      |
| $\leq$ 39.22% Increase | 12 (C)            | 166 (D)                | 178 (C+D)     |
| Total                  | 30 (A+C)          | 175 (B+D)              | 205 (A+B+C+D) |

Overall Agreement =

$(A+D) / (A+B+C+D) \times 100\% = 89.76\%$ ; 95% CI = 84.77% to 93.55%

Positive Agreement =

$A / (A+C) \times 100\% = 60.00\%$ ; 95% CI = 40.60% to 77.34%

Negative Agreement =

$D / (B+D) \times 100\% = 94.86\%$ ; 95% CI = 90.46% to 97.62%



| Instrument               | Reagent Lot | Sample         | n   | Mean (ng/mL) | Within-Run |     | Within-Laboratory Precision (Total) <sup>a</sup> |     |
|--------------------------|-------------|----------------|-----|--------------|------------|-----|--|-----|
|                          |             |                |     |              | SD         | %CV | SD   | %CV |
| i 2000 <sub>SH</sub> (1) | 1           | Low Control    | 120 | 19.81        | 0.317      | 1.6 | 0.327  | 1.6 |
|                          |             | Medium Control | 120 | 199.11       | 3.105      | 1.6 | 3.263  | 1.6 |
|                          |             | High Control   | 120 | 950.55       | 16.411     | 1.7 | 17.200   | 1.8 |
|                          | 2           | Low Control    | 120 | 20.02        | 0.349      | 1.7 | 0.349  | 1.7 |
|                          |             | Medium Control | 120 | 195.29       | 2.725      | 1.4 | 3.043  | 1.6 |
|                          |             | High Control   | 120 | 928.69       | 16.340     | 1.8 | 17.628   | 1.9 |
|                          | 3           | Low Control    | 120 | 20.23        | 0.248      | 1.2 | 0.286  | 1.4 |
|                          |             | Medium Control | 120 | 198.45       | 2.743      | 1.4 | 3.058  | 1.5 |
|                          |             | High Control   | 120 | 955.96       | 17.389     | 1.8 | 17.389   | 1.8 |
| i 2000 <sub>SH</sub> (2) | 1           | Panel 1        | 120 | 3.01         | 0.070      | 2.3 | 0.082  | 2.7 |
|                          |             | Panel 2        | 120 | 9.54         | 0.191      | 2.0 | 0.201  | 2.1 |
|                          |             | Panel 3        | 120 | 577.98       | 13.137     | 2.3 | 13.977   | 2.4 |
|                          |             | Panel 4        | 120 | 1514.74      | 41.437     | 2.7 | 47.765   | 3.2 |
|                          |             | Panel 5        | 120 | 1763.53      | 43.353     | 2.5 | 51.115   | 2.9 |
|                          | 2           | Panel 1        | 120 | 3.10         | 0.060      | 1.9 | 0.065  | 2.1 |
|                          |             | Panel 2        | 120 | 9.67         | 0.188      | 1.9 | 0.202  | 2.1 |
|                          |             | Panel 3        | 120 | 564.10       | 13.445     | 2.4 | 14.358   | 2.5 |
|                          |             | Panel 4        | 120 | 1489.33      | 43.567     | 2.9 | 44.077   | 3.0 |
|                          |             | Panel 5        | 120 | 1729.13      | 50.297     | 2.9 | 54.344   | 3.1 |
|                          | 3           | Panel 1        | 120 | 3.15         | 0.061      | 1.9 | 0.068  | 2.2 |
|                          |             | Panel 2        | 120 | 9.73         | 0.190      | 2.0 | 0.197  | 2.0 |
|                          |             | Panel 3        | 120 | 559.72       | 12.053     | 2.2 | 12.053   | 2.2 |
|                          |             | Panel 4        | 120 | 1490.94      | 43.967     | 2.9 | 45.619   | 3.1 |
|                          |             | Panel 5        | 120 | 1743.06      | 53.149     | 3.0 | 55.158   | 3.2 |
| i 2000 (1)               | 1           | Low Control    | 120 | 19.53        | 0.403      | 2.1 | 0.419  | 2.1 |
|                          |             | Medium Control | 120 | 192.55       | 3.196      | 2.0 | 4.161  | 2.2 |
|                          |             | High Control   | 120 | 925.44       | 20.138     | 2.2 | 22.571   | 2.4 |
|                          | 2           | Low Control    | 120 | 19.60        | 0.460      | 2.3 | 0.476  | 2.4 |
|                          |             | Medium Control | 120 | 192.92       | 4.090      | 2.1 | 4.233  | 2.2 |
|                          |             | High Control   | 120 | 917.39       | 25.481     | 2.8 | 26.696   | 2.9 |
|                          | 3           | Low Control    | 120 | 19.46        | 0.412      | 2.1 | 0.438  | 2.2 |
|                          |             | Medium Control | 120 | 194.92       | 3.602      | 1.8 | 3.602  | 1.8 |
|                          |             | High Control   | 120 | 942.88       | 25.442     | 2.7 | 26.752   | 2.8 |
| i 2000 (2)               | 1           | Panel 1        | 120 | 3.05         | 0.075      | 2.5 | 0.082  | 2.7 |
|                          |             | Panel 2        | 120 | 9.63         | 0.204      | 2.1 | 0.205  | 2.1 |
|                          |             | Panel 3        | 120 | 569.00       | 15.652     | 2.8 | 16.041   | 2.8 |
|                          |             | Panel 4        | 120 | 1530.79      | 56.561     | 3.7 | 60.769   | 4.0 |
|                          |             | Panel 5        | 120 | 1796.54      | 63.911     | 3.6 | 70.192   | 3.9 |
|                          | 2           | Panel 1        | 120 | 3.10         | 0.070      | 2.3 | 0.074  | 2.4 |
|                          |             | Panel 2        | 120 | 9.58         | 0.186      | 1.9 | 0.212  | 2.2 |
|                          |             | Panel 3        | 120 | 553.43       | 13.334     | 2.4 | 13.972   | 2.5 |
|                          |             | Panel 4        | 120 | 1454.43      | 46.282     | 3.2 | 49.527   | 3.4 |
|                          |             | Panel 5        | 120 | 1693.81      | 53.646     | 3.2 | 56.319   | 3.3 |
|                          | 3           | Panel 1        | 120 | 2.87         | 0.071      | 2.5 | 0.077  | 2.7 |
|                          |             | Panel 2        | 120 | 9.31         | 0.237      | 2.6 | 0.262  | 2.8 |
|                          |             | Panel 3        | 120 | 583.73       | 14.561     | 2.5 | 16.269   | 2.8 |
|                          |             | Panel 4        | 120 | 1456.88      | 44.111     | 3.0 | 45.569   | 3.1 |
|                          |             | Panel 5        | 120 | 1666.40      | 46.687     | 2.8 | 49.738   | 3.0 |

<sup>a</sup> Within-Laboratory (Total) variability contains within-run, within-day and between-day variance components.

### WHO Recovery

The ARCHITECT AFP assay is designed to have a recovery range of 100 ± 10% when analyzing samples spiked with known amounts of AFP using the WHO 1<sup>st</sup> International Standard 72/225.

A study was performed with 16 low-level AFP serum specimens and 14 amniotic fluid specimens. The serum specimens were spiked with the WHO 1<sup>st</sup> International Standard 72/225 to create test samples across the measuring interval of the assay. The amniotic fluid specimens were diluted 1:40 using the ARCHITECT i Multi-Assay Manual Diluent and spiked with the WHO 1<sup>st</sup> International Standard 72/225 to create test samples with AFP concentrations within the range of 312.5 to 1250 ng/mL.\* The samples were tested using the ARCHITECT AFP assay on 1 instrument and the resulting percent recovery was calculated. For serum specimens, the mean percent recovery was 103.1% (range: 99.5% to 108.8%). For amniotic fluid specimens, the mean percent recovery was 101.2% (range: 95.1% to 107.3%).

\* A 1:40 dilution of specimens in the range of 12.5 to 50 µg/mL equates to a range of 312.5 to 1250 ng/mL within the measuring interval of the assay.

### Linearity

The ARCHITECT AFP assay is designed to have a deviation from linearity within ± 1 ng/mL for samples less than 10 ng/mL and within ± 10% for samples between 10 ng/mL and 2000 ng/mL.

A study was performed based on guidance from the NCCLS document EP6-A.<sup>56</sup> Three dilution series were prepared as follows: a high AFP sample (> 2000 ng/mL) was combined in specific ratios with a low AFP sample (< 2.0 ng/mL). The 3 dilution series, including the low-level and high-level samples, were tested using the ARCHITECT AFP assay. The ARCHITECT AFP assay demonstrated linearity from 0.01 ng/mL to 2487.76 ng/mL.

### Sensitivity

#### Limit of Detection and Limit of Quantitation

The ARCHITECT AFP assay is designed to have a limit of Detection (LoD) of ≤ 1.0 ng/mL and a Limit of Quantitation (LoQ) of ≤ 2.0 ng/mL. The LoQ is defined as the lowest amount of analyte in a sample that can be accurately quantitated with a total analytical error of ≤ 2.5 ng/mL.

Based on guidance from the NCCLS document EP17-A,<sup>58</sup> a study was performed with 4 zero-level samples (Calibrator A) and 8 low-level AFP samples (2 samples at each of 4 unique target concentration levels of approximately 0.50, 1.00, 1.50 and 2.00 ng/mL). These samples were tested in 5 separate runs over a minimum of 3 days using 3 reagent lots and 2 instruments. The observed LoD was 0.04 ng/mL and the observed LoQ was 0.5 ng/mL.

#### Limit of Blank

In the same study, the Limit of Blank (LoB) was determined to be 0.0 ng/mL.

### Interference

#### Potentially Interfering Endogenous Substances

The ARCHITECT AFP assay is designed to have a difference in AFP concentration within or equal to ± 10% when comparing samples containing elevated levels of endogenous substances to reference samples.

A study was performed based on guidance from the CLSI document EP7-A2.<sup>57</sup> Potentially interfering endogenous substances were evaluated to determine whether AFP concentrations were affected when using the ARCHITECT AFP assay. The endogenous substances listed below were spiked into samples with 2 levels of AFP (approximately 10 and 1000 ng/mL). The samples were assayed (n = 20) and the AFP concentrations of the spiked samples were compared to reference samples. The data are summarized in the following table.

| Potentially Interfering Endogenous Substance | High Test Level | % Interference <sup>a</sup> |            |
|--|-----------------|-----------------------------|------------|
|  |                 | 10 ng/mL                    | 1000 ng/mL |
| Bilirubin (Unconjugated)                     | 20 mg/dL        | -0.5                        | 0.3        |
| Bilirubin (Conjugated)                       | 20 mg/dL        | -0.9                        | -0.8       |
| Hemoglobin                                   | 500 mg/dL       | -0.2                        | -1.3       |
| Total Protein                                | 12 g/dL         | 2.9                         | -0.2       |
| Triglycerides                                | 3000 mg/dL      | -1.0                        | -1.9       |

Mean/Median Test Result -

Mean/Median Reference Result

$$^a \text{ \% Interference} = \frac{\text{Mean/Median Test Result} - \text{Mean/Median Reference Result}}{\text{Mean/Median Reference Result}} \times 100$$

Mean/Median Reference Result



**Potentially Interfering Substances**

The ARCHITECT AFP assay is designed to have a mean % recovery of 100% ± 10% when analyzing Rheumatoid Factor (RF) and Human Anti-Mouse Antibodies (HAMA) samples spiked with known amounts of AFP.

A study was performed based on guidance from the CLSI document EP7-A2.<sup>57</sup> Potentially interfering substances were evaluated to determine whether AFP concentrations were affected when using the ARCHITECT AFP assay. Specimens from individuals with the substances listed below were divided into 3 samples. Two of the samples were spiked to 2 levels of AFP (approximately 10 and 1000 ng/mL). The samples were assayed and the AFP concentrations of the spiked samples were compared to the samples that were not spiked with AFP. The data are summarized in the following table.

| Potentially Interfering Substances | n  | % Recovery <sup>a</sup> |            |
|------------------------------------|----|-------------------------|------------|
|                                    |    | 10 ng/mL                | 1000 ng/mL |
| Human Anti-Mouse Antibodies        | 13 | 104.7                   | 105.9      |
| Rheumatoid Factor                  | 13 | 104.6                   | 102.3      |

$$^a \text{ \% Recovery} = \frac{\text{Mean/Median Spiked Result} - \text{Mean/Median Unspiked Result}}{\text{Mean/Median Amount AFP Added}} \times 100$$

**Analytical Specificity**

The ARCHITECT AFP assay is designed to have a difference in AFP concentration within or equal to ± 10% when comparing samples containing potential interferents to reference samples.

A study was performed based on guidance from the CLSI document EP7-A2.<sup>57</sup> Potential interferents were evaluated to determine whether AFP concentrations were affected when using the ARCHITECT AFP assay. The potential interferents were spiked into samples with 2 levels of AFP (approximately 10 and 1000 ng/mL). The samples were assayed and the AFP concentrations of the spiked samples were compared to the reference samples. The data are summarized in the following table.

| Potential Interferent     | High Test Level | % Interference <sup>a</sup> |            |
|---------------------------|-----------------|-----------------------------|------------|
|                           |                 | 10 ng/mL                    | 1000 ng/mL |
| 5-Fluorouracil            | 3 mmol/L        | -0.4                        | 0.3        |
| Acetaminophen             | 6.5 mg/mL       | -3.1                        | -3.1       |
| Albumin                   | 160 mg/mL       | 2.4                         | -4.1       |
| Alpha-1-Acid Glycoprotein | 2 mg/mL         | 0.2                         | -1.1       |
| Alpha-1-Antitrypsin       | 5 mg/mL         | 8.1                         | 0.3        |
| Alpha-2-Macroglobulin     | 9 mg/mL         | 0.1                         | -0.1       |
| Aspirin                   | 10 mg/mL        | -4.7                        | -4.9       |
| Bleomycin                 | 1000 µU/mL      | -2.5                        | -4.3       |
| Carboplatin               | 0.432 mg/mL     | 0.3                         | 1.3        |
| Ceruloplasmin             | 2.5 mg/mL       | -0.3                        | -0.6       |
| Chorionic Gonadotropin    | 1000 IU/mL      | -1.1                        | -2.2       |
| Cisplatin                 | 1000 µg/mL      | -0.6                        | -1.2       |
| Cyclophosphamide          | 1437 µmol/L     | 0.3                         | -0.7       |

| Potential Interferent | High Test Level | % Interference <sup>a</sup> |            |
|-----------------------|-----------------|-----------------------------|------------|
|                       |                 | 10 ng/mL                    | 1000 ng/mL |
| Etoposide             | 30 µg/mL        | -0.6                        | 0.3        |
| Gamma-Globulins       | 30 mg/mL        | -2.7                        | -2.4       |
| Haptoglobin           | 6 mg/mL         | 0.7                         | -1.0       |
| Ifosfamide            | 249 µg/mL       | -3.1                        | -2.7       |
| Methotrexate          | 2 mmol/L        | -0.6                        | -0.5       |
| Placental Lactogen    | 100 µg/mL       | -3.3                        | -3.6       |
| Prolactin             | 500 ng/mL       | -4.8                        | -5.0       |
| Transferrin           | 25 mg/mL        | -1.6                        | -3.4       |
| Vinblastine           | 500 µg/mL       | -3.4                        | -3.5       |
| Vincristine           | 1000 ng/mL      | -3.1                        | -4.4       |

$$^a \text{ \% Interference} = \frac{\text{Mean/Median Test Result} - \text{Mean/Median Reference Result}}{\text{Mean/Median Reference Result}} \times 100$$

**Autodilution Verification**

The ARCHITECT AFP assay is designed to have a mean difference in concentration within ± 10% when comparing the autodilution method to the manual dilution method for samples with values > 2000 ng/mL.

Twenty-one serum samples were evaluated with the 1:10 autodilution method versus a 1:10 manual dilution method. Fifteen amniotic fluid samples were evaluated with the 1:40 autodilution method versus a 1:40 manual dilution method.

The manually diluted samples and the undiluted samples designated for autodilution were assayed in replicates of 2 using the ARCHITECT AFP assay. For serum samples, the mean percent difference was 2.9% (range: -5.8% to 10.9%) and for amniotic fluid samples, the mean percent difference was 4.6% (range: -1.1% to 11.2%).

**High Dose Hook**

High dose hook is a phenomenon whereby very high level specimens may read within the measuring interval of the assay. For the ARCHITECT AFP assay, no high dose hook effect was observed when samples containing up to 10,000,000 ng/mL of AFP were assayed.

**ARCHITECT i1000<sub>SR</sub> SYSTEM SPECIFIC STUDIES**

The following studies were conducted using the ARCHITECT i1000<sub>SR</sub> System. Assay results obtained in individual laboratories may vary from data presented.

**Precision**

The ARCHITECT AFP assay is designed to have an imprecision of ≤ 7.5% within-laboratory (Total) %CV for samples between 10 and 2000 ng/mL and an SD of ≤ 0.75 for samples less than 10 ng/mL down to the LoQ (i.e., 2.0 ng/mL).

**System Reproducibility**

A 5-day precision study was performed for the ARCHITECT AFP assay based on guidance from the National Committee for Clinical Laboratory Standards (NCCLS) document EP5-A2<sup>58</sup> and the Clinical and Laboratory Standards Institute (CLSI) document EP15-A2.<sup>54</sup> Testing was conducted at 3 clinical sites using 1 lot of ARCHITECT AFP Reagents, Calibrators, and Controls and 1 ARCHITECT i1000<sub>SR</sub> instrument per site. Three controls and 5 human serum panels were assayed in replicates of 4 at 2 separate times of day for 5 days. The results are summarized in the following table.

| Sample         | n   | Grand Mean (ng/mL) | Within-Run |     | Within-Day |     | Within-Laboratory Precision (Total) <sup>a</sup> |     | Precision with Additional Component of Between-Site (Overall) |     |
|----------------|-----|--------------------|------------|-----|------------|-----|--|-----|---|-----|
|                |     |                    | SD         | %CV | SD         | %CV | SD   | %CV | SD  | %CV |
|                |     |                    |            |     |            |     |  |     |   |     |
| Low Control    | 120 | 20.02              | 0.319      | 1.6 | 0.373      | 1.9 | 0.522  | 2.6 | 0.522   | 2.6 |
| Medium Control | 120 | 200.45             | 2.835      | 1.4 | 2.901      | 1.4 | 4.662  | 2.3 | 5.350   | 2.7 |
| High Control   | 120 | 959.62             | 16.945     | 1.8 | 20.440     | 2.1 | 26.446   | 2.8 | 26.446  | 2.8 |
| Panel 1        | 120 | 3.15               | 0.049      | 1.6 | 0.084      | 2.8 | 0.121  | 4.0 | 0.144   | 4.7 |
| Panel 2        | 120 | 9.59               | 0.149      | 1.5 | 0.198      | 2.0 | 0.264  | 2.7 | 0.290   | 3.3 |
| Panel 3        | 120 | 578.85             | 9.449      | 1.6 | 11.025     | 1.9 | 14.212   | 2.5 | 15.580  | 2.9 |
| Panel 4        | 120 | 1494.96            | 36.675     | 2.5 | 38.393     | 2.6 | 48.435   | 3.2 | 49.885  | 3.5 |
| Panel 5        | 120 | 1723.05            | 37.272     | 2.2 | 44.664     | 2.6 | 50.517   | 2.9 | 49.295  | 3.4 |

<sup>a</sup> Within-Laboratory (Total) variability contains within-run, within-day, and between-day variance components.



### Within-Laboratory Precision

A study was performed based on guidance from the NCCLS document EP5-A2.<sup>53</sup> Testing was conducted using 3 lots of ARCHITECT AFP Reagents and Calibrators, 1 lot of ARCHITECT AFP Controls, and 2 instruments. Three controls and 5 human serum panels were assayed in a minimum of 2 replicates at 2 separate times per day for 20 different days. Each reagent lot used a single calibration curve throughout the study. The data are summarized in the following table.

| Instrument              | Reagent Lot | Sample         | n   | Mean (ng/mL) | Within-Run |     | Within-Laboratory Precision (Total) <sup>a</sup> |     |
|-------------------------|-------------|----------------|-----|--------------|------------|-----|--|-----|
|                         |             |                |     |              | SD         | %CV | SD   | %CV |
| i1000 <sub>SR</sub> (1) | 1           | Low Control    | 120 | 19.91        | 0.431      | 2.2 | 0.495  | 2.5 |
|                         |             | Medium Control | 120 | 195.96       | 3.684      | 1.9 | 4.274  | 2.2 |
|                         |             | High Control   | 120 | 953.56       | 24.91      | 2.6 | 27.613   | 2.9 |
|                         | 2           | Low Control    | 120 | 19.61        | 0.373      | 1.9 | 0.429  | 2.2 |
|                         |             | Medium Control | 120 | 194.04       | 3.712      | 1.9 | 4.630  | 2.4 |
|                         |             | High Control   | 120 | 941.16       | 25.317     | 2.7 | 28.814   | 3.1 |
|                         | 3           | Low Control    | 120 | 19.68        | 0.407      | 2.1 | 0.491  | 2.5 |
|                         |             | Medium Control | 120 | 195.83       | 3.767      | 1.9 | 4.335  | 2.2 |
|                         |             | High Control   | 120 | 939.12       | 23.829     | 2.5 | 24.118   | 2.6 |
| i1000 <sub>SR</sub> (2) | 1           | Panel 1        | 120 | 3.05         | 0.057      | 1.9 | 0.063  | 2.1 |
|                         |             | Panel 2        | 119 | 9.67         | 0.143      | 1.5 | 0.212  | 2.2 |
|                         |             | Panel 3        | 120 | 573.62       | 11.594     | 2.0 | 14.622   | 2.5 |
|                         |             | Panel 4        | 120 | 1492.07      | 35.493     | 2.4 | 45.627   | 3.1 |
|                         |             | Panel 5        | 120 | 1731.53      | 44.680     | 2.6 | 51.398   | 3.0 |
|                         | 2           | Panel 1        | 119 | 3.14         | 0.056      | 1.8 | 0.072  | 2.3 |
|                         |             | Panel 2        | 120 | 9.80         | 0.170      | 1.7 | 0.217  | 2.2 |
|                         |             | Panel 3        | 120 | 557.31       | 10.350     | 1.9 | 12.635   | 2.3 |
|                         |             | Panel 4        | 120 | 1456.72      | 34.841     | 2.4 | 41.001   | 2.8 |
|                         |             | Panel 5        | 120 | 1702.52      | 46.419     | 2.7 | 53.955   | 3.2 |
|                         | 3           | Panel 1        | 120 | 3.09         | 0.051      | 1.9 | 0.070  | 2.3 |
|                         |             | Panel 2        | 120 | 9.72         | 0.155      | 1.6 | 0.190  | 2.0 |
|                         |             | Panel 3        | 120 | 570.03       | 11.540     | 2.0 | 13.176   | 2.3 |
|                         |             | Panel 4        | 119 | 1471.65      | 37.123     | 2.5 | 50.591   | 3.4 |
|                         |             | Panel 5        | 120 | 1698.64      | 39.124     | 2.3 | 48.897   | 2.9 |

<sup>a</sup> Within-Laboratory (Total) variability contains within-run, within-day, and between-day variance components.

### Comparison Between the ARCHITECT i1000<sub>SR</sub> System and the ARCHITECT i2000/i2000<sub>SR</sub> System

The comparison between the ARCHITECT i1000<sub>SR</sub> and the ARCHITECT i2000/i2000<sub>SR</sub> was evaluated by testing 205 serum panel members and 205 manually diluted amniotic fluid panel members using 1 lot each of ARCHITECT AFP Reagents, Calibrators, and Controls. Testing of each sample type was performed on 1 ARCHITECT i1000<sub>SR</sub> instrument at each of 3 clinical testing sites and on 1 ARCHITECT i2000/i2000<sub>SR</sub> instrument at 1 clinical testing site.

### Regression

The panel members were evaluated using the Deming regression method. The data are summarized in the following table.

| Sample Type            | Concentration Range (ng/mL) |                           | Correlation Coefficient (r) |                     | Intercept | 95% CI <sup>a</sup> | Slope | 95% CI <sup>a</sup> |
|------------------------|-----------------------------|---------------------------|-----------------------------|---------------------|-----------|---------------------|-------|---------------------|
|                        | i1000 <sub>SR</sub>         | i2000/i2000 <sub>SR</sub> | r                           | 95% CI <sup>a</sup> |           |                     |       |                     |
| Diluted Amniotic Fluid | 5.29 - 1929.41              | 4.65 - 1918.58            | 0.999                       | (0.998, 0.999)      | -8.16     | (-11.88, -4.49)     | 0.98  | (0.97, 0.99)        |
| Serum                  | 2.93 - 1949.81              | 2.84 - 1955.56            | 0.999                       | (0.999, 0.999)      | 4.71      | (3.00, 6.43)        | 0.96  | (0.95, 0.97)        |

<sup>a</sup> 95% CI = Confidence Interval

### BIBLIOGRAPHY

- Bergstrand CG, Czar B. Demonstration of a new protein fraction in serum from the human fetus. *Scand J Clin Lab Invest* 1956;8:174.
- Ruoslahti E, Engvall E, Kessler MJ. Chemical properties of alpha-fetoprotein. In: Herberman RB, McIntire KR, eds. *Immunodiagnosis of Cancer*. New York: Marcel Dekker, Inc., 1979;101-17.
- Ruoslahti E, Seppälä M. Studies of carcino-fetal proteins: physical and chemical properties of human alpha-fetoprotein. *Int J Cancer* 1971;7:218-25.
- Tatarinov YS. Detection of embryo-specific alpha-globulin in the blood serum of patients with primary liver tumors. *Vopr Med Khim* 1964;10:90-1.
- Silver HKB, Gold P, Feder S, et al. Radioimmunoassay for human alpha-fetoprotein. *Proc Natl Acad Sci USA* 1973;70(2):526-30.
- Waldmann TA, McIntire KR. The use of a radioimmunoassay for alpha-fetoprotein in the diagnosis of malignancy. *Cancer* 1974;34:1510-5.
- Kohn J, Orr AH, McElwain TJ, et al. Serum-alpha-fetoprotein in patients with testicular tumours. *Lancet* 1976;2:433-6.
- Abelev GI. Alpha-fetoprotein as a model for studying reexpression of embryonic antigens in neoplasia. In: Herberman RB, McIntire KR, eds. *Immunodiagnosis of Cancer*. New York: Marcel Dekker, Inc., 1979:76-101.
- Scardino PT, Cox HD, Waldmann TA, et al. The value of serum tumor markers in the staging and prognosis of germ cell tumors of the testis. *J Urol* 1977;118:994-9.
- Bosl GJ, Lange PH, Fraley EE, et al. Human chorionic gonadotropin and alpha-fetoprotein in the staging of nonseminomatous testicular cancer. *Cancer* 1981;47:328-32.
- Lange PH, McIntire KR, Waldmann TA, et al. Serum alpha fetoprotein and human chorionic gonadotropin in the diagnosis and management of nonseminomatous germ-cell testicular cancer. *N Engl J Med* 1976;295(22):1237-40.
- Javadpour N, McIntire KR, Waldmann TA. Human chorionic gonadotropin (hCG) and alpha-fetoprotein (AFP) in sera and tumor cells of patients with testicular seminoma. *Cancer* 1978;42:2768-72.
- Report from the Medical Research Council Working Party on Testicular Tumours. Prognostic factors in advanced non-seminomatous germ-cell testicular tumours: results of a multicentre study. *Lancet* 1985;8:1-1.
- Perlin E, Engeler JE, Edson M, et al. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. *Cancer* 1976;37:215-9.
- Wepsic HT. Alpha-fetoprotein: its quantitation and relationship to neoplastic disease. In: Kirkpatrick AM, Nakamura RM, eds. *Alpha-fetoprotein, laboratory procedures and clinical applications*. New York: Masson Publishing USA Inc., 1981:15-29.
- McIntire KR, Waldmann TA, Moertel CG, et al. Serum alpha-fetoprotein in patients with neoplasms of the gastrointestinal tract. *Cancer Res* 1975;35:991-6.
- Chen DS, Sung JL. Relationship of hepatitis B surface antigen to serum alpha-fetoprotein in nonmalignant diseases of the liver. *Cancer* 1979;44:984-92.
- Report of U.K. Collaborative Study on Alpha-fetoprotein in Relation to Neural-tube Defects. Maternal serum-alpha-fetoprotein measurement in antenatal screening for anencephaly and spina bifida in early pregnancy. *Lancet* 1977:1323-32.



19. Second Report of the U.K. Collaborative Study on Alpha-fetoprotein in Relation to Neural-tube Defects. Amniotic-fluid alpha-fetoprotein measurement in antenatal diagnosis of anencephaly and open spina bifida in early pregnancy. *Lancet* 1979;651-82.
20. Haddow JE, Kloza EM, Smith DE, et al. Data from an alpha-fetoprotein pilot screening program in Maine. *Obstet Gynecol* 1983;62(5):556-60.
21. Brock DJH. The prenatal diagnosis of neural tube defects. *Obstet Gynecol Surv* 1976;31(1):32-40.
22. Main DM, Mennuti MT. Neural tube defects: Issues in prenatal diagnosis and counselling. *Obstet Gynecol* 1986;67(1):1-18.
23. Adams MJ, Windham GC, James LM, et al. Clinical interpretation of maternal serum alpha-fetoprotein concentrations. *Am J Obstet Gynecol* 1984;148(3):241-54.
24. American Society of Human Genetics policy statement for maternal serum alpha-fetoprotein screening programs and quality control for laboratories performing maternal serum and amniotic fluid alpha-fetoprotein assays. *Am J Hum Genet* 1987;40:75-82.
25. Cuckle HS, Nanchahal K, Wald NJ. Maternal serum alpha-fetoprotein and ethnic origin. *Br J Obstet Gynaecol* 1987;94:1111-2.
26. Crandall BF, Leberer TB, Schrott PC, et al. Alpha-fetoprotein concentrations in maternal serum: relation to race and body weight. *Clin Chem* 1983;29(3):531-3.
27. Milunsky A, Alpert E, Kitzmiller JL, et al. Prenatal diagnosis of neural tube defects. VIII. The importance of serum alpha-fetoprotein screening in diabetic pregnant women. *Am J Obstet Gynecol* 1982;142:1030-2.
28. Baumgarten A, Robinson J. Prospective study of an inverse relationship between maternal glycosylated hemoglobin and serum alpha-fetoprotein concentrations in pregnant women with diabetes. *Am J Obstet Gynecol* 1988;159(1):77-81.
29. Palomaki GE, Knight GJ, Kloza EM, et al. Maternal weight adjustment and low serum alpha-fetoprotein values. *Lancet* 1985;468.
30. Wald NJ, Cuckle H, Boreham J, et al. Maternal serum alpha-fetoprotein and diabetes mellitus. *Br J Obstet Gynaecol* 1979;88:101-5.
31. Crandall BF. Second trimester maternal serum screening to identify neural tube defects. In: Kirkpatrick AM, Nakamura RM, eds. *Alpha-fetoprotein, laboratory procedures and clinical applications*. New York: Masson Publishing USA Inc., 1981:93-105.
32. Seppälä M, Rapola J, Huttunen NP, et al. Congenital nephrotic syndrome: prenatal diagnosis and genetic counselling by estimation of amniotic fluid and maternal serum alpha-fetoprotein. *Lancet* 1976; 123-5.
33. Seppälä M. Increased alpha fetoprotein in amniotic fluid associated with a congenital esophageal atresia of the fetus. *Obstet Gynecol* 1973;42(4):613-4.
34. Palomaki GE, Hill LE, Knight GJ, et al. Second-trimester maternal serum alpha-fetoprotein levels in pregnancies associated with gastroschisis and omphalocele. *Obstet Gynecol* 1988;71(6):906-9.
35. Wald N, Barker S, Cuckle H, et al. Maternal serum alpha-fetoprotein and spontaneous abortion. *Br J Obstet Gynaecol* 1977;84:357-62.
36. Brock DJH, Barron L, Duncan P, et al. Significance of elevated mid-trimester maternal plasma-alpha-fetoprotein values. *Lancet* 1979:1281-2.
37. Nelson LH, Bensen J, Burton BK. Outcome in patients with unusually high maternal serum alpha-fetoprotein levels. *Am J Obstet Gynecol* 1987;157(3):572-6.
38. Redford DHA, Whitfield CR. Maternal serum alpha-fetoprotein in twin pregnancies uncomplicated by neural tube defect. *Am J Obstet Gynecol* 1985;152(5):550-3.
39. Davenport DM, Macri JN. The clinical significance of low maternal serum alpha-fetoprotein. *Am J Obstet Gynecol* 1983;146(6):657-61.
40. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
41. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2008.
42. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
43. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
44. Horacek I, Pepperell RJ, Hay DL, et al. Detection of fetomaternal haemorrhage by measurement of maternal serum-alpha-fetoprotein. *Lancet* 1976:200.
45. Clinical and Laboratory Standards Institute. *Maternal Serum Screening: Approved Standard – Second Edition*. CLSI document I/LA25-A2. Wayne, PA: CLSI; 2011.
46. Sherwin JE ed. National Academy of Clinical Biochemistry: Laboratory Medicine Practice Guidelines. *Maternal-fetal risk assessment and reference values in pregnancy*. Washington, DC: AACC Press, 2006:24.
47. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-4.
48. Schrott RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-85.
49. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
50. Bradley LA, Palomaki GE, and McDowell GA. ACMG Standards and Guidelines: technical standards and guidelines: prenatal screening for open neural tube defects. *Genetics in Medicine* 2005;7(5):355-69.
51. Fraser CG. *Biological Variation: From Principles to Practice*. Washington, DC: AACC Press, 2001:71-8.
52. Trapé J, Botargues JM, Porta F, et al. Reference change value for alpha-fetoprotein and its application in early detection of hepatocellular carcinoma in patients with hepatic disease. *Clin Chem* 2003;49:1209-11.
53. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
54. Clinical and Laboratory Standards Institute. *User Verification of Performance for Precision and Trueness; Approved Guideline – Second Edition*. CLSI document EP15-A2. Wayne, PA: CLSI; 2005.
55. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. NCCLS Document EP6-A. Wayne, PA: NCCLS; 2003.
56. National Committee for Clinical Laboratory Standards (NCCLS). *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. NCCLS Document EP17-A. Wayne, PA: NCCLS; 2004.
57. Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP7-A2. Wayne, PA: CLSI; 2005.

The following U.S. Patents are relevant to the ARCHITECT System or its components. There are other such patents and patent applications in the United States and worldwide.

|           |           |           |
|-----------|-----------|-----------|
| 5 468 646 | 5 543 524 | 5 545 739 |
| 5 565 570 | 5 669 819 | 5 783 699 |

ARCHITECT and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions.  
ProClin is property of its respective owner.



Abbott Ireland  
Diagnostics Division  
Finisklin Business Park  
Sligo  
Ireland  
+353-71-9171712



Distributed by Abbott Laboratories  
Abbott Park, IL 60064 USA  
and  
ABBOTT  
65205 Wiesbaden, Germany

October 2015  
© 2011, 2015 Abbott Laboratories





REF 7K68-27  
 REF 7K68-22  
 REF 7K68-35  
 REF 7K68-32



Read Highlighted Changes! Revised November 2015.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

**WARNING:** The concentration of CEA in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CEA assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining CEA levels serially is changed, additional sequential testing should be carried out. Before changing assays, the laboratory MUST confirm baseline values for patients being serially monitored.

**CAUTION:** United States Federal Law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by, or on the order of a physician.

**NAME**

ARCHITECT CEA (carcinoembryonic antigen)

**INTENDED USE**

The ARCHITECT CEA assay is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of Carcinoembryonic Antigen (CEA) in human serum and plasma. The ARCHITECT CEA assay is to be used as an aid in the prognosis and management of cancer patients in whom changing concentrations of CEA are observed.

**SUMMARY AND EXPLANATION OF THE TEST**

Carcinoembryonic antigen (CEA), first described in 1965 by Gold and Freedman,<sup>1</sup> is a tumor associated antigen. CEA was characterized as a glycoprotein of approximately 200,000 molecular weight with a  $\beta$ -electrophoretic mobility.<sup>2,3</sup> Subsequent development of a radioimmunoassay (RIA) by Thomson, et al<sup>4</sup> made it possible to detect the very low concentrations of CEA in blood, other body fluids, and also in normal and diseased tissues.<sup>5-7</sup> Two years later, Hansen, et al<sup>8</sup> developed a modified RIA for CEA.

The result of clinical studies to date indicate that CEA, although originally thought to be specific for digestive tract cancers, may also be elevated in other malignancies and in some nonmalignant disorders.<sup>9-15</sup>

CEA testing can have significant value in the monitoring of patients with diagnosed malignancies in whom changing concentrations of CEA are observed. A persistent elevation in circulating CEA following treatment is strongly indicative of occult metastatic and/or residual disease.<sup>16-20</sup>

A persistently rising CEA value may be associated with progressive malignant disease and a poor therapeutic response.<sup>21-23</sup> A declining CEA value is generally indicative of a favorable prognosis and a good response to treatment.<sup>21, 23, 24</sup> Patients who have low pretherapy CEA levels may later show elevations in the CEA level as an indication of progressive disease.<sup>25</sup>

Clinical relevance of the CEA assay has been shown in the follow-up management of patients with colorectal, gastric, breast, lung, prostatic, pancreatic, and ovarian carcinoma.<sup>18, 24, 26-31</sup> Follow-up studies of patients with colorectal, breast, and lung carcinoma suggest that the preoperative CEA level has prognostic significance.<sup>32-35</sup>

CEA testing is not recommended as a screening procedure to detect cancer in the general population; however, use of the CEA test as an adjunctive test in predicting prognosis and as an aid in the management of cancer patients has been widely accepted.

**BIOLOGICAL PRINCIPLES OF THE PROCEDURE**

The ARCHITECT CEA assay is a two-step immunoassay to determine the presence of CEA in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and anti-CEA coated paramagnetic microparticles are combined. The CEA present in the sample binds to the anti-CEA coated microparticles.
2. After washing, anti-CEA acridinium-labeled conjugate is added to create a reaction mixture.
3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of CEA in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

**REAGENTS**

**Kit Contents**

ARCHITECT CEA 7K68

**NOTE:** Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

| REF                   | 7K68-27  | 7K68-22    | 7K68-35     | 7K68-32     |
|-----------------------|--|------------|-------------|-------------|
| $\Sigma$              | 100  | 400        | 500         | 2000        |
| <b>MICROPARTICLES</b> | 1 x 6.6 mL   | 4 x 6.6 mL | 1 x 27.0 mL | 4 x 27.0 mL |
| <b>CONJUGATE</b>      | 1 x 5.9 mL   | 4 x 5.9 mL | 1 x 26.3 mL | 4 x 26.3 mL |
| <b>MICROPARTICLES</b> | anti-CEA (mouse, monoclonal) coated Microparticles in TRIS buffer with protein (bovine) stabilizer. Minimum concentration: 0.1% solids. Preservative: Antimicrobial Agents.                |            |             |             |
| <b>CONJUGATE</b>      | anti-CEA (mouse, monoclonal) acridinium-labeled Conjugate in phosphate buffer with protein (bovine) stabilizer. Minimum concentration: 0.8 $\mu$ g/mL. Preservative: Antimicrobial Agents. |            |             |             |

**Other Reagents**

**MULTI-ASSAY MANUAL DILUENT** 1 x 100 mL ARCHITECT Multi-Assay Manual Diluent, REF 7D82-50, containing phosphate buffered saline solution. Preservative: antimicrobial agent.

**PRE-TRIGGER SOLUTION** ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

**TRIGGER SOLUTION** ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

**WASH BUFFER** ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

**NOTE:** Bottle and volume varies based on order.





## Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use

### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.

Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

### Reagent Handling

- Do not use reagent kits beyond the expiration date.
- **Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
  - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
  - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
  - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

### Reagent Storage

- When stored and handled as directed, reagents are stable until the expiration date.

|                      | Storage Temperature | Maximum Storage Time  | Additional Storage Instructions  |
|----------------------|---------------------|-----------------------|--|
| Unopened/<br>Opened* | 2-8°C               | Until expiration date | May be used immediately after removal from 2-8°C storage.  |
| On board             | System temperature  | 30 days               | Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5. |

\* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

## Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

The ARCHITECT CEA assay file must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

Verified specimen types to be used with this assay:

| Specimen Types | Collection Tubes                               |
|----------------|--|
| Human serum    | Serum<br>Serum separator tubes                 |
| Human plasma   | Heparin (sodium and lithium)<br>Potassium EDTA |

- Other specimen collection tube types have not been tested with this assay.
- Plasma specimens collected in lithium or sodium heparin have been shown to exhibit an average of 7% to 8% higher results compared to corresponding serum results.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

### Specimen Conditions

- Do not use specimens with the following conditions:
  - grossly hemolyzed
  - obvious microbial contamination
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

### Preparation for Analysis

- Follow the tube manufacturer's processing instructions for specimen collection tubes.
- ARCHITECT CEA Calibrators and Controls should be mixed by gentle inversion prior to use.
  - For optimal results, specimens should be free of fibrin, red blood cells, or other particulate matter. **Centrifuge serum and plasma specimens containing fibrin, red blood cells, or particulate matter prior to use to ensure consistency in the results.**
- Specimens must be mixed THOROUGHLY after thawing, by vortexing. Thawed samples containing red blood cells or particulate matter, **or which are hazy or cloudy in appearance** must be centrifuged prior to use to ensure consistency in the results.



- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

#### Specimen Storage

| Specimen Type | Storage Temperature | Maximum Storage Time |
|---------------|---------------------|----------------------|
| Serum/Plasma  | 2-8°C               | ≤ 7 days             |

If testing will be delayed more than 24 hours, serum or plasma should be removed from the clot, serum separator, or red blood cells. If testing will be delayed more than 7 days, specimens should be stored/frozen at -20°C or colder. Avoid multiple freeze/thaw cycles.

#### Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

### PROCEDURE

#### Materials Provided

7K68 ARCHITECT CEA Reagent Kit

#### Materials Required but not Provided

- ARCHITECT CEA Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com).
- 7K68-02 ARCHITECT CEA Calibrators
- 7D82-50 ARCHITECT Multi-Assay Manual Diluent
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

#### Materials Available but not Provided:

- 7K68-12 ARCHITECT CEA Controls

#### Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
  - Invert the microparticle bottle 30 times.
  - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
  - If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.
  - Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the Reagent Handling section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
  - Verify that all necessary reagents are present.
  - Ensure that septums are present on all reagent bottles.
- Order tests.
  - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.

- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

- Priority:
  - Sample volume for first test: 60 µL
  - Sample volume for each additional test from same sample cup: 10 µL
- ≤ 3 hours on board:
  - Sample volume for first test: 150 µL
  - Sample volume for each additional test from same sample cup: 10 µL
- > 3 hours on board: Additional sample volume is required. Refer to the ARCHITECT System Operations Manual, Section 5 for information on sample evaporation and volumes.
- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT CEA Calibrators and Controls.
  - Mix calibrator(s) and controls by gentle inversion before use.
  - Hold bottles vertically and dispense recommended volumes into each respective sample cup.
  - Recommended volumes:
    - for each calibrator: 5 drops
    - for each control: 4 drops
- Load samples.
  - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

#### Specimen Dilution Procedures

Specimens with a CEA value exceeding 1500 ng/mL are flagged with the code ">1500.00" and may be diluted using either the Automated Dilution Protocol or the Manual Dilution Procedure.

#### Automated Dilution Protocol

The system performs a 1:10 dilution of the specimen and automatically calculates the concentration of the specimen before dilution and reports the result.

#### Manual Dilution Procedure

Suggested dilution: 1:100

An additional 1:10 dilution may be made if needed. It is recommended that dilutions not exceed 1:1000.

- Add 20 µL of the patient specimen to 1980 µL of ARCHITECT Multi-Assay Manual Diluent.
- The operator must enter the dilution factor in the Patient or Control order screen. All assays selected for that order will be diluted. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution and report the result. The dilution should be performed so that the diluted result reads > 4 ng/mL.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

A comparison of the Automated Dilution Protocol to the Manual Dilution Procedure yielded recoveries between 86% and 97%.



### Calibration

- Test Calibrators 1 and 2 in duplicate. The calibrators should be priority loaded.  
A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0 - 500 ng/mL.
- The assay protocol allows for the range to be extended to 1500 ng/mL.
- Once an ARCHITECT CEA calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
  - A reagent kit with a new lot number is used or
  - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

### Quality Control Procedures

The recommended control requirement for the ARCHITECT CEA assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

Ensure that assay control values are within the concentration ranges specified in the package insert.

### Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix 3.

The ARCHITECT CEA assay belongs to method group 1.

## RESULTS

### Calculation

The ARCHITECT CEA assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

### Flags

- The default result unit for the ARCHITECT CEA assay is ng/mL.
- Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

## LIMITATIONS OF THE PROCEDURE

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits such as ARCHITECT CEA that employ mouse monoclonal antibodies. Additional information may be required for diagnosis.<sup>40, 41</sup>  
ARCHITECT CEA reagents contain a component that reduces the effect of HAMA reactive specimens. Additional clinical or diagnostic information may be required to determine patient status.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.<sup>42</sup>
- The ARCHITECT CEA assay should not be used as a cancer screening test.

- Patients with confirmed carcinoma frequently have a pretreatment CEA level in the same range as healthy individuals. Elevations in circulating CEA levels may be observed in smokers as well as patients with nonmalignant disease. For these reasons, a serum or plasma CEA level, regardless of value, should not be interpreted as absolute evidence for the presence or absence of malignant disease. The CEA level should be used in conjunction with information available from clinical evaluation and other diagnostic procedures.

## EXPECTED VALUES

The distribution of ARCHITECT CEA values determined in 1,141 specimens is shown in the following table.\*

Distribution of ARCHITECT CEA Values

|                             | Number of Subjects | Percent (%)   |                 |                  |              |
|-----------------------------|--------------------|---------------|-----------------|------------------|--------------|
|                             |                    | 0 - 3 (ng/mL) | > 3 - 5 (ng/mL) | > 5 - 10 (ng/mL) | > 10 (ng/mL) |
| <u>Healthy Subjects</u>     |                    |               |                 |                  |              |
| Smokers                     | 159                | 74.2          | 18.2            | 6.9              | 0.6          |
| Non-smokers                 | 149                | 83.2          | 11.4            | 5.4              | 0.0          |
| Total                       | 308                | 78.6          | 14.9            | 6.2              | 0.3          |
| <u>Nonmalignant Disease</u> |                    |               |                 |                  |              |
| Ulcerative Colitis          | 50                 | 72.0          | 20.0            | 4.0              | 4.0          |
| Rectal Polyps               | 78                 | 83.3          | 10.3            | 5.1              | 1.3          |
| Pulmonary                   | 60                 | 61.7          | 20.0            | 13.3             | 5.0          |
| Cirrhosis                   | 110                | 47.3          | 30.0            | 15.5             | 7.3          |
| Hepatitis                   | 60                 | 70.0          | 16.7            | 11.7             | 1.7          |
| Renal                       | 20                 | 60.0          | 15.0            | 15.0             | 10.0         |
| <u>Malignant Disease</u>    |                    |               |                 |                  |              |
| Colorectal                  | 150                | 24.0          | 10.7            | 10.7             | 54.7         |
| Gastric                     | 37                 | 62.2          | 5.4             | 10.8             | 21.6         |
| Pulmonary                   | 110                | 47.3          | 19.1            | 9.1              | 24.5         |
| Mammary                     | 117                | 62.4          | 11.1            | 10.3             | 16.2         |
| Ovarian                     | 41                 | 78.0          | 7.3             | 2.4              | 12.2         |

\* Representative data; results in individual laboratories may vary. In this study, 93.5% of healthy subjects (n=308) had CEA values of 5.00 ng/mL or less.

It is expected that each laboratory establish its own expected reference range for the population of interest.

The distribution table above for malignant disease is derived primarily from patients representing both active (clinical evidence of disease progression) and inactive (no clinical evidence of disease progression) disease states. When changing CEA assay methods in the course of monitoring a patient, additional sequential testing should be carried out to confirm baseline values.

## SPECIFIC PERFORMANCE CHARACTERISTICS

### Precision

The Architect CEA assay precision is  $\leq 8\%$ . Precision was determined as described in the Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) Protocol EP5-T2,<sup>43</sup> Five samples, consisting of two serum based panels and three CEA controls, were assayed at three laboratories in replicates of two at two separate times per day for twenty days (n=80 for each sample), using a single lot of reagents and a single calibration. Data from this study are summarized in the following table.\*



**Reproducibility of ARCHITECT CEA**

| Sample         | Lab | Within Run       |        |     | Total  |     |
|----------------|-----|------------------|--------|-----|--------|-----|
|                |     | Mean CEA (ng/mL) | SD     | %CV | SD     | %CV |
| Low Control    | 1   | 5.05             | 0.180  | 3.6 | 0.202  | 4.0 |
|                | 2   | 4.79             | 0.098  | 2.1 | 0.178  | 3.7 |
|                | 3   | 4.86             | 0.110  | 2.3 | 0.162  | 3.3 |
| Medium Control | 1   | 20.17            | 0.512  | 2.5 | 0.641  | 3.2 |
|                | 2   | 19.08            | 0.515  | 2.7 | 0.629  | 3.3 |
|                | 3   | 19.99            | 0.605  | 3.0 | 0.695  | 3.4 |
| High Control   | 1   | 99.45            | 3.074  | 3.1 | 3.182  | 3.2 |
|                | 2   | 93.97            | 2.082  | 2.2 | 2.559  | 2.7 |
|                | 3   | 99.51            | 2.898  | 2.9 | 3.072  | 3.1 |
| Panel 1        | 1   | 417.34           | 9.587  | 2.3 | 10.483 | 2.5 |
|                | 2   | 395.70           | 11.313 | 2.9 | 13.995 | 3.5 |
|                | 3   | 419.93           | 11.591 | 2.8 | 13.870 | 3.3 |
| Panel 2        | 1   | 1294.72          | 40.660 | 3.1 | 46.508 | 3.6 |
|                | 2   | 1185.03          | 21.570 | 1.8 | 26.285 | 2.2 |
|                | 3   | 1309.28          | 29.760 | 2.3 | 38.030 | 2.9 |

\* Representative performance data are shown. Results obtained at individual laboratories may vary.

**Recovery**

Known amounts of CEA were added to normal human serum and plasma samples. The concentration of CEA was determined using the ARCHITECT CEA assay and the resulting percent recovery was calculated.\*

| Sample Type                | Endogenous Level (ng/mL) | Recovery          |                      |                  |
|----------------------------|--------------------------|-------------------|----------------------|------------------|
|                            |                          | CEA Added (ng/mL) | CEA Observed (ng/mL) | Percent Recovery |
| Serum                      | 1                        | 0.86              | 94.76                | 95.1             |
|                            | 2                        | 1.07              | 4.49                 | 100.0            |
|                            | 3                        | 0.94              | 94.76                | 98.8             |
|                            | 4                        | 1.11              | 4.49                 | 105.6            |
| Average % Recovery: 99.9%  |                          |                   |                      |                  |
| EDTA                       | 1                        | 0.81              | 94.76                | 95.77            |
|                            | 2                        | 0.70              | 94.76                | 95.7             |
|                            | 3                        | 1.10              | 4.49                 | 104.0            |
|                            | 4                        | 1.72              | 4.49                 | 100.0            |
| Average % Recovery: 99.2%  |                          |                   |                      |                  |
| Heparin                    | 1                        | 0.93              | 94.76                | 94.59            |
|                            | 2                        | 1.26              | 4.49                 | 6.10             |
|                            | 3                        | 0.92              | 94.76                | 95.24            |
|                            | 4                        | 1.17              | 4.49                 | 5.92             |
| Average % Recovery: 103.0% |                          |                   |                      |                  |

\* Representative data; results in individual laboratories may vary.

$$\% \text{ Recovery} = \frac{\text{Observed (ng/mL)} - \text{Endogenous Level (ng/mL)}}{\text{CEA Added (ng/mL)}} \times 100$$

**Analytical Sensitivity**

The sensitivity of the ARCHITECT CEA assay was calculated to be better than 0.5 ng/mL at the 95% level of confidence (n = 18 runs). Sensitivity is defined as the concentration at two standard deviations above the mean RLU for the ARCHITECT CEA MasterCheck Level 0 and represents the lowest measurable concentration of CEA that can be distinguished from zero.

**Specificity**

The specificity of the ARCHITECT CEA assay was determined by testing sera containing the compounds listed below. These compounds showed less than 10% interference in the ARCHITECT CEA assay at the levels indicated.

| Test Compound | Test Concentration |
|---------------|--------------------|
| Bilirubin     | 22 mg/dL           |
| Hemoglobin    | 550 mg/dL          |
| Total Protein | 1.8 to 13.2 g/dL   |
| Triglycerides | 3300 mg/dL         |

**Carryover**

No detectable carryover (less than 12 PPM) was observed when a sample containing 43,630 ng/mL of CEA was assayed.

**High Dose Hook**

High dose hook is a phenomenon whereby very high level specimens may read within the dynamic range of the assay. For the ARCHITECT CEA assay, no high dose hook effect was observed when samples containing up to approximately 60,000 ng/mL of CEA were assayed.

**BIBLIOGRAPHY**

- Gold P, Freedman SO. Demonstration of Tumor-Specific Antigens in Human Colonic Carcinoma by Immunological Tolerance and Absorption Techniques. *J Exp Med* 1964;121:439.
- Krupey J, Gold P, Freedman SO. Physicochemical Studies of the Carcinoembryonic Antigens of the Human Digestive System. *J Exp Med* 1968;183:387.
- Krupey J, Wilson T, Freedman SO, et al. The Preparation of Purified Carcinoembryonic Antigen of the Human Digestive System from Large Quantities of Tumor Tissue. *Immunochem* 1972;9(6):617-622.
- Thomson DMP, Krupey J, Freedman SO, et al. The Radioimmunoassay of Circulating Carcinoembryonic Antigen of the Human Digestive System. *Proc Natl Acad Sci USA* 1969;64:161-167.
- Zarncheck N. Carcinoembryonic Antigen: Quantitative Variations in Circulating Levels in Benign and Malignant Digestive Tract Diseases. *Adv Intern Med* 1974;19:413.
- Go VLM, Ammon HV, Holtermuller KH, et al. Quantification of Carcinoembryonic Antigen-Like Activities in Normal Human Gastrointestinal Secretions. *Cancer* 1975;36:2346-2350.
- Khoo SK, Warner NL, Lie JT, et al. Carcinoembryonic Antigenic Activity of Tissue Extracts: A Quantitative Study of Malignant and Benign Neoplasms. Cirrhotic Liver, Normal Adult and Fetal Organs. *Int J Cancer* 1973;11(3):681-687.
- Hansen HJ, Lance KP, Krupey J. Demonstration of an Ion Sensitive Antigenic Site on Carcinoembryonic Antigen Using Zirconyl Phosphate Gel. *Clin Res* 1971;19:143.
- Oehr P, Schlosser T, Adolphs HD. Applicability of an Enzymatic Test for the Determination of CEA in Serum and CEA-Like Products in Urine of Patients with Bladder Cancer. *Tumor Diagn* 1980;1:40.
- Ng WW, Tong KJ, Tam TN, et al. Clinical Values of CA 19-9, CA125, and CEA in Malignant Obstructive Jaundice. *Chung Hua I Hsueh Tsa Chih (Taipei)* 1995;55(6):438-446.
- Jerzsek B, Cervek J, Rudolf Z, et al. Clinical Evaluation of Potential Usefulness of CEA, CA 15-3, and MCA in Follow-up Breast Cancer Patients. *Cancer Letters* 1996;110(1-2):137-144.
- Wollenberg B, Jan V, Schmit UM, et al. CYFRA 21-1 is not Superior to SCC Antigen and CEA in Head and Neck Squamous Cell Cancer. *Anticancer Res* 1996;16(5b):3117-3124.
- Cerwonka H, Aigner R, Quehenberger F, et al. Preoperative Differential Diagnosis of Benign and Malignant Pancreatic Lesions - The Value of Pancreatic Secretory Trypsin Inhibitor, Procarboxypeptidase B, CA 19-9 and CEA. *Hepato-Gastroenterology* 1997;44(16):1117-1121.
- Pastor A, Menendez R, Cremades JM, et al. Diagnostic Value of SCC, CEA and CYFRA 21.1 in Lung Cancer: A Bayesian Analysis. *Eur Respir J* 1997;10(3):603-609.
- Maestranzi S, Przemioslo R, Mitchell H, et al. The Effect of Benign and Malignant Liver Disease on the Tumour Markers CA 19-9 and CEA. *Ann Clin Biochem* 1998;35:99-103.
- Duk JM, Aalders JG, Fluere GJ, et al. Tumour Markers CA 125, Squamous Cell Carcinoma Antigen, and Carcinoembryonic Antigen in Patients with Adenocarcinoma of the Uterine Cervix. *Obstet Gynecol* 1989;73(4):661-668.
- Munck-Wikland E, Kuylenstierna R, Lindholm T, et al. Carcinoembryonic Antigen, CA 19-9, and CA 50 in Monitoring Human Squamous Cell Carcinoma of the Esophagus. *Anticancer Res* 1990;10(3):703-708.
- Jager W, Kramer S, Palapelas V, et al. Breast Cancer and Clinical Utility of CA 15-3 and CEA. *Scand J Clin Lab Invest Suppl* 1995;221:87-92.



19. Ebert W, Hoppe M, Muley T, et al. Monitoring of Therapy in Inoperable Lung Cancer Patients by Measurement of CYFRA 21-1, TPA-M, TPS, CEA, and NSE. *Anticancer Res* 1997;17(4B): 2875-2878.
20. King J, Caplehorn JRM, Ross WB, et al. High Serum Carcinoembryonic Antigen Concentrations in Patients with Colorectal Liver Metastases is Associated with Poor Cell-Mediated Immunity, Which is Predictive of Survival. *Br J Surg* 1987;84:1382-1385.
21. Noda M, Kusunoki M, Yanagi H, et al. Serum Carcinoembryonic Antigen (CEA) Correlates with the Survival Time During 5-FU Hepatic Arterial Infusion Chemotherapy for Unresectable Colorectal Hepatic Metastase - Clinical Study. *International Journal of Oncology* 1996;9(4):741-746.
22. Koronaga D, Saeki H, Mawatari K, et al. Serum Carcinoembryonic Antigen Concentration Doubling Time Correlates with Tumor Biology and Life Expectancy in Patients with Recurrent Gastrointestinal Carcinoma. *Archives of Surgery* 1997;132(2):188-194.
23. Nakayama T, Watanabe M, Teramoto T, et al. Slope Analysis of CA 19-9 and CEA for Predicting Recurrence in Colorectal Cancer Patients. *Anticancer Res* 1997;17(2b):1379-1382.
24. Lokich JJ, Zamcheck N, Lowenstein M. Sequential Carcinoembryonic Antigen Levels in the Therapy of Metastatic Breast Cancer. *Ann Intern Med* 1978;89:902.
25. Yasue M, Sakamoto J, Teramukai S, et al. Prognostic Values of Preoperative and Postoperative CEA and CA 19-9 Levels in Pancreatic Cancer. *Pancreas* 1994;9(6):735-740.
26. Zamcheck N, Martin EW. Factors Controlling the Circulating CEA Levels in Pancreatic Cancer. Some Clinical Correlations. *Cancer* 1981;47:1620.
27. Alsabli EAK, Kamel A. Carcinoembryonic Antigen (CEA) in Patients with Malignant and Non-Malignant Disease. *Neoplasma* 1979;26:603.
28. Khoo SK, Whitaker S, Jones I, et al. Predictive Value of Serial Carcinoembryonic Antigen Levels in Long-Term Follow-Up of Ovarian Cancer. *Cancer* 1979;43:2471.
29. Chevinsky AH. CEA in Tumors of Other than Colorectal Origin. *Semin Surg Oncol* 1991;7(3):162-166.
30. Diez M, Gomez A, Hernandez F, et al. Serum CEA, CA125, and SCC Antigens and Tumor Recurrence in Resectable Non-Small Cell Lung Cancer. *Int J Biol Markers* 1995;10(1):5-10.
31. Bast RC, Bates S, Bradt AB, et al. Clinical Practice Guidelines for the Use of Tumor Markers in Breast and Colorectal Cancer. *J Clin Oncol* 1996;14(10):2843-2877.
32. Wanebo HJ, Rao B, Pinsky CM, et al. Preoperative Carcinoembryonic Antigen Level as a Prognostic Indicator in Colorectal Cancer. *N Engl J Med* 1978;299(9):448-451.
33. Concannon JP, Dalbow MH, Hodgson SE, et al. Prognostic Value of Preoperative Carcinoembryonic Antigen (CEA) Plasma Levels in Patients with Bronchogenic Carcinoma. *Cancer* 1978;42:1477-1483.
34. Ikeda Y, Mori M, Kajiyama K, et al. Indicative Value of Carcinoembryonic Antigen (CEA) for Liver Recurrence Following Curative Resection of Stage II and III Gastric Cancer. *Hepatogastroenterol* 1996;43(11):1251-1257.
35. Gebauer G, Muller-Ruchholtz W. Tumor Marker Concentrations in Normal and Malignant Tissues of Colorectal Cancer Patients and Their Prognostic Relevance. *Anticancer Res* 1997;17(4a):2731-2734.
36. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
37. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
38. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
39. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
40. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
41. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
42. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
43. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Clinical Chemistry Devices, Tentative Guideline—Second Edition*. NCCLS Document EP5-T2. Villanova, PA: NCCLS; 1992.

## Key to Symbols

|                                   |  |
|-----------------------------------|--|
|                                   | Consult instructions for use                         |
|                                   | Manufacturer   |
|                                   | Sufficient for                                       |
|                                   | Temperature limitation                               |
|                                   | Use by/Expiration date                               |
| <b>CONJUGATE</b>                  | Conjugate  |
| <b>CONTROL NO.</b>                | Control Number                                       |
| <b>DISTRIBUTED IN THE USA BY</b>  | Distributed in the USA by                            |
| <b>GTIN</b>                       | Global Trade Item Number                             |
| <b>INFORMATION FOR USA ONLY</b>   | Information needed for United States of America only |
| <b>IVD</b>                        | In Vitro Diagnostic Medical Device                   |
| <b>LOT</b>                        | Lot Number   |
| <b>MICROPARTICLES</b>             | Microparticles                                       |
| <b>MULTI-ASSAY MANUAL DILUENT</b> | Multi-Assay Manual Diluent                           |
| <b>PRE-TRIGGER SOLUTION</b>       | Pre-Trigger Solution                                 |
| <b>PRODUCT OF IRELAND</b>         | Product of Ireland                                   |
| <b>REACTION VESSELS</b>           | Reaction Vessels                                     |
| <b>REAGENT LOT</b>                | Reagent Lot  |
| <b>REF</b>                        | List Number  |
| <b>REPLACEMENT CAPS</b>           | Replacement Caps                                     |
| <b>SAMPLE CUPS</b>                | Sample Cups  |
| <b>SEPTUM</b>                     | Septum   |
| <b>SN</b>                         | Serial number  |
| <b>TRIGGER SOLUTION</b>           | Trigger Solution                                     |
| <b>WASH BUFFER</b>                | Wash Buffer  |

ARCHITECT, Chemiflex and MasterCheck are trademarks of Abbott Laboratories in various jurisdictions.



Abbott Ireland  
Diagnostics Division  
Finisklin Business Park  
Sligo  
Ireland  
+353-71-9171712



### DISTRIBUTED IN THE USA BY

Abbott Laboratories  
Abbott Park, IL 60064 USA

**Customer Service: Contact your local representative or find country-specific contact information on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)**

Revised November 2015.

©2005, 2015 Abbott Laboratories

**Abbott**





REF 2K91-24  
REF 2K91-32  
REF 2K91-39



 **en**  
CA 19-9XR  
2K91  
613-031 11/16/R03  
B2K9Y0

Read Highlighted Changes: Revised November 2016.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

**WARNING:** The Abbott ARCHITECT CA 19-9XR CMA assay utilizes an antibody/antigen system based on the 1116-NS-19-9 antibody. The unique reagent formulation employed in the ARCHITECT CA 19-9XR assay may return elevated concentrations when compared to other methods for samples expressing high levels of 1116-NS-19-9 reactive determinants.<sup>1, 2</sup> Additionally, there is no internationally recognized standard for CA 19-9, which can contribute to differences between assay methods. The ARCHITECT CA 19-9XR assay is standardized to a reference standard prepared by Fujirebio Diagnostics, Inc. Performance characteristics of the Abbott ARCHITECT CA 19-9XR assay are NOT transferable to other diagnostic kits.

The concentration of 1116-NS-19-9 reactive determinants obtained with different assay methods cannot be used interchangeably due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CA 19-9 assay used. If, in the course of monitoring a patient, the assay method used for determining serial 1116-NS-19-9 reactive determinant levels is changed, additional sequential testing should be carried out. Prior to changing assays, the laboratory MUST confirm baseline values for patients being serially monitored.

**WARNING:** 1116-NS-19-9 reactive determinants are shed naturally in saliva and other body fluids.<sup>3</sup> Contamination of the samples or the ARCHITECT iSystem disposables with saliva or aerosols (e.g., as a result of sneezing) may cause falsely elevated CA 19-9 assay values. It is recommended that all elevated values be reviewed and testing repeated as appropriate. Gloves should always be worn when handling samples, sample cups, reaction vessels, and septums. Face masks are also recommended.

#### NAME

ARCHITECT CA 19-9XR

#### INTENDED USE

The ARCHITECT CA 19-9XR assay is a chemiluminescent microparticle immunoassay (CMIA) for the quantitative determination of 1116-NS-19-9 reactive determinants in human serum or plasma on the ARCHITECT iSystem. The ARCHITECT CA 19-9XR assay is to be used as an aid in the management of pancreatic cancer patients in conjunction with other clinical methods.

#### SUMMARY AND EXPLANATION OF THE TEST

The ARCHITECT CA 19-9XR assay detects a tumor-associated antigen, which occurs in tissue as a monosialoganglioside and in serum as a high molecular weight, carbohydrate-rich glycoprotein known as a mucin.<sup>4-7</sup>

The ARCHITECT CA 19-9XR assay is based upon a monoclonal antibody, 1116-NS-19-9, which reacts with a carbohydrate antigenic determinant expressed on the circulating antigen.<sup>4-6</sup>

The results of published research studies<sup>8-14</sup> indicate that the CA 19-9 assay value is frequently elevated in the serum of subjects with various gastrointestinal conditions, such as pancreatic, colorectal, gastric, and hepatic carcinomas. No data exist to support the use of CA 19-9 in screening for malignancies.<sup>15, 16</sup> The role of CA 19-9 is to be used as an adjunct with other diagnostic information in the management of patients with pancreatic cancer.<sup>15</sup> Increased serum CA 19-9 assay values have also been observed in patients with metastases and in nonmalignant conditions such as hepatitis, cirrhosis, pancreatitis, and other gastrointestinal disease.<sup>8-11, 17-20</sup> Elevated levels have also been seen in cystic fibrosis.<sup>21-24</sup> Research studies demonstrate that CA 19-9 assay values may have utility in monitoring subjects with the above-mentioned diagnosed gastrointestinal malignancies.<sup>25-28</sup> It has been shown that a persistent elevation in CA 19-9 assay value following treatment may be indicative of occult metastatic and/or residual disease. A persistently rising CA 19-9 assay value may be associated with progressive malignant disease and poor therapeutic response. A declining CA 19-9 assay value may be indicative of a favorable prognosis and a good response to treatment.<sup>29-35</sup>

Testing for 1116-NS-19-9 reactive determinants must not be used as a screening procedure for malignancy. 1116-NS-19-9 reactive determinants are present as a normal constituent in serum and plasma of individuals without gastrointestinal carcinomas or having certain aforementioned non-cancer related conditions.

#### BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT CA 19-9XR assay is a two-step immunoassay for the quantitative determination of 1116-NS-19-9 reactive determinants in human serum or plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and 1116-NS-19-9 coated paramagnetic microparticles are combined. The 1116-NS-19-9 reactive determinants present in the sample bind to the 1116-NS-19-9 coated microparticles.
2. After washing, 1116-NS-19-9 acridinium-labeled conjugate is added to create a reaction mixture.
3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of 1116-NS-19-9 reactive determinants in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.



## REAGENTS

### Kit Contents

ARCHITECT CA 19-9XR 2K91

**NOTE:** Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

| REF                   | 2K91-32  | 2K91-24    | 2K91-39     |
|-----------------------|--|------------|-------------|
|                       | 100  | 400        | 500         |
| <b>MICROPARTICLES</b> | 1 x 6.6 mL   | 4 x 6.6 mL | 1 x 27.0 mL |
| <b>CONJUGATE</b>      | 1 x 5.9 mL   | 4 x 5.9 mL | 1 x 26.3 mL |
| <b>MICROPARTICLES</b> | 1116-NS-19-9 (mouse, monoclonal) coated microparticles in citrate buffer with protein (bovine) stabilizer. Minimum concentration: 0.09% solids. Preservatives: sodium azide and ProClin 300.       |            |             |
| <b>CONJUGATE</b>      | 1116-NS-19-9 (mouse, monoclonal) acridinium-labeled conjugate in phosphate buffer with protein (bovine) stabilizer. Minimum concentration: 0.5 µg/mL. Preservatives: sodium azide and ProClin 300. |            |             |

### Other Reagents

**MULTI-ASSAY MANUAL DILUENT** 1 x 100 mL ARCHITECT Multi-Assay Manual Diluent, **REF** 7D82-50, containing phosphate buffered saline solution. Preservative: antimicrobial agent.

**PRE-TRIGGER SOLUTION** ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

**TRIGGER SOLUTION** ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

**WASH BUFFER** ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

### Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use

### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.<sup>36-39</sup>

| The following warnings and precautions apply to: <b>MICROPARTICLES</b> |  |
|--|--|
|  |  |
| <b>WARNING</b>   | Contains methylisothiazolones and sodium azide.                        |
| H317   | May cause an allergic skin reaction.                                   |
| EUH032   | Contact with acids liberates very toxic gas.                           |
| <b>Prevention</b>  |  |
| P261   | Avoid breathing mist / vapors / spray.                                 |
| P272   | Contaminated work clothing should not be allowed out of the workplace. |
| P280   | Wear protective gloves / protective clothing / eye protection.         |
| <b>Response</b>  |  |
| P302+P352  | IF ON SKIN: Wash with plenty of water.                                 |
| P333+P313  | If skin irritation or rash occurs: Get medical advice / attention.     |
| P362+P364  | Take off contaminated clothing and wash it before reuse.               |

| Disposal  |  |
|---|--|
| P501  | Dispose of contents / container in accordance with local regulations.  |
| The following warnings and precautions apply to: <b>CONJUGATE</b> |  |
|   |  |
| <b>DANGER</b>   | Contains polyethylene glycol cetylphenyl ether, methylisothiazolones and sodium azide.   |
| H317  | May cause an allergic skin reaction.   |
| H318  | Causes serious eye damage.   |
| H412  | Harmful to aquatic life with long lasting effects.   |
| EUH032  | Contact with acids liberates very toxic gas.   |
| <b>Prevention</b>   |  |
| P261  | Avoid breathing mist / vapors / spray.   |
| P280  | Wear protective gloves / protective clothing / eye protection.   |
| P272  | Contaminated work clothing should not be allowed out of the workplace.   |
| P273  | Avoid release to the environment.  |
| <b>Response</b>   |  |
| P302+P352   | IF ON SKIN: Wash with plenty of water.   |
| P333+P313   | If skin irritation or rash occurs: Get medical advice / attention.   |
| P305+P351+P338  | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P310  | Immediately call a POISON CENTER or doctor / physician.  |
| P362+P364   | Take off contaminated clothing and wash it before reuse.   |
| <b>Disposal</b>   |  |
| P501  | Dispose of contents / container in accordance with local regulations.  |

Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

### Reagent Handling

- **1116-NS-19-9 reactive determinants are shed naturally in saliva and other body fluids.<sup>3</sup> Contamination of the samples or the ARCHITECT iSystem disposables with saliva or aerosols (e.g., as a result of sneezing) may cause falsely elevated CA 19-9 assay values. It is recommended that all elevated values be reviewed and testing repeated as appropriate. Gloves should always be worn when handling samples, sample cups, reaction vessels, and septums. Face masks are also recommended.**
- Do not use reagent kits beyond the expiration date.
- **Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**



- To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
- Once a septum has been placed on an open reagent bottle, do not invert the bottle as this will result in reagent leakage and may compromise assay results.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

### Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

|                  | Storage Temperature | Maximum Storage Time  | Additional Storage Instructions   |
|------------------|---------------------|-----------------------|---|
| Unopened/Opened* | 2-8°C               | Until expiration date | May be used immediately after removal from 2-8°C storage.<br>Store in upright position.   |
| On board         | System temperature  | 30 days               | Discard after 30 days.<br>For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5. |

\* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright.

If any reagent bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded. For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

### Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

The ARCHITECT CA 19-9XR assay file must be installed on the ARCHITECT iSystem prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

Verified specimen types to be used with this assay:

| Specimen Types | Collection Tubes      |
|----------------|-----------------------|
| Serum          | Serum                 |
|                | Serum separator tubes |
|                | Tripotassium EDTA     |
| Plasma         | Sodium Heparin        |
|                | Lithium Heparin       |

- Other specimen collection tube types have not been tested with this assay.

- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.
- Follow the tube manufacturer's processing instructions for serum and plasma collection tubes.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.

### Specimen Conditions

- Do not use specimens with the following conditions:
  - grossly hemolyzed
  - obvious microbial contamination
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.
- Performance has not been established using body fluids other than human serum or plasma.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

### Preparation for Analysis

- Follow the tube manufacturer's processing instructions for specimen collection tubes.
- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, centrifuge specimens before testing if
  - they contain fibrin, red blood cells, or other particulate matter,
  - they require repeat testing, or
  - they were frozen and thawed.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

### Specimen Storage

| Specimen Type | Storage Temperature | Maximum Storage Time |
|---------------|---------------------|----------------------|
| Serum/Plasma  | 2-8°C               | ≤ 7 days             |
|               | -20°C or colder     | >7 days              |

- If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator or red blood cells.
- Specimens may be stored for up to 7 days at 2-8°C prior to being tested.
- If testing will be delayed more than 7 days, serum or plasma should be stored frozen at -20°C or colder.
- Avoid multiple freeze/thaw cycles.

### Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.





## PROCEDURE

### Materials Provided

2K91 ARCHITECT CA 19-9XR Reagent Kit

### Materials Required but not Provided

- ARCHITECT CA 19-9XR Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com).
- 2K91 ARCHITECT CA 19-9XR Calibrators
- 2K91 ARCHITECT CA 19-9XR Controls
- 7D82-50 ARCHITECT Multi-Assay Manual Diluent
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

### Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
  - **Invert the microparticle bottle 30 times.**
  - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
  - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
  - Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
  - Verify that all necessary reagents are present.
  - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
  - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
  - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

- Priority:
  - Sample volume for first test: 90  $\mu$ L
  - Sample volume for each additional test from same sample cup: 30  $\mu$ L
- $\leq$  3 hours on board:
  - Sample volume for first test: 150  $\mu$ L
  - Sample volume for each additional test from same sample cup: 30  $\mu$ L
- $>$  3 hours on board: Additional sample volume required

- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT CA 19-9XR Calibrators and Controls.
  - Mix calibrator(s) and controls by gentle inversion before use.
  - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
  - Recommended volumes:
    - for each calibrator: 4 drops
    - for each control: 4 drops
- Load samples.
  - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

### Specimen Dilution Procedures

Specimens with an ARCHITECT CA 19-9XR value exceeding 1200 U/mL are flagged with the code "> 1200.00" and may be diluted using either the Automated Dilution Protocol or the Manual Dilution Procedure.

#### Automated Dilution Protocol

The system performs a 1:10 dilution of the specimen and automatically calculates the concentration of the specimen before dilution and reports the result.

#### Manual Dilution Procedure

Suggested dilution: 1:10

An additional 1:10 dilution may be made if needed.

1. Add 50  $\mu$ L of the patient specimen to 450  $\mu$ L of ARCHITECT Multi-Assay Manual Diluent (7D82).
2. The operator must enter the dilution factor in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution and report the result. The dilution should be performed so that the diluted result reads greater than 30 U/mL.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

### Calibration

- Test Calibrators A-F in duplicate. The calibrators should be priority loaded.

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0 - 1200 U/mL.
- Once an ARCHITECT CA 19-9XR calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
  - A reagent kit with a new lot number is used or
  - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

### Quality Control Procedures

The recommended control requirement for the ARCHITECT CA 19-9XR assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

The ARCHITECT CA 19-9XR values must be within the acceptable ranges specified in the control package insert. If a control is out of the specified range, the associated test results are invalid and samples must be retested. Recalibration may be indicated.



### Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT CA 19-9XR assay belongs to method group 1.

## RESULTS

### Calculation

The ARCHITECT CA 19-9XR assay utilizes a Linear Regression data reduction method to generate a calibration curve.

### Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

## LIMITATIONS OF THE PROCEDURE

- The ARCHITECT CA 19-9XR assay value must be used in conjunction with information available from clinical evaluation and other diagnostic procedures.
- If the ARCHITECT CA 19-9XR results are inconsistent with clinical evidence, additional testing is recommended.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.<sup>40</sup>
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Specimens containing HAMA may produce anomalous values when tested with assay kits that employ mouse monoclonal antibodies. ARCHITECT CA 19-9XR reagents contain a component that reduces the effect of HAMA reactive specimens. Additional clinical or diagnostic information may be required to determine patient status.<sup>41, 42</sup>
- Patients with confirmed carcinoma may have pretreatment CA 19-9 assay values in the same range as healthy individuals. Elevations in circulating 1116-NS-19-9 reactive determinants may be observed in patients with metastases and in nonmalignant conditions such as hepatitis, cirrhosis, pancreatitis, and other gastrointestinal disease. Elevated levels have also been seen in cystic fibrosis.<sup>21</sup> For these reasons, a CA 19-9 assay value, regardless of level, should not be interpreted as absolute evidence for the presence or absence of malignant disease. The ARCHITECT CA 19-9XR assay must not be used as a cancer screening test.
- Patients with the Le<sup>a-b-</sup> phenotype may not express the 1116-NS-19-9 reactive determinant.<sup>43</sup>
- Representative performance data are given in the EXPECTED VALUES and SPECIFIC PERFORMANCE CHARACTERISTICS sections. Results obtained in individual laboratories may vary.

## EXPECTED VALUES

### APPARENTLY HEALTHY SUBJECTS

A study was performed with three hundred sixty (360) serum specimens from apparently healthy individuals. The distribution of ARCHITECT CA 19-9XR assay values from these specimens is shown in the table below.\*

|                             | Number of Subjects | Percent (%)   |                 |                  |                   |              |
|-----------------------------|--------------------|---------------|-----------------|------------------|-------------------|--------------|
|                             |                    | 0-37.0 (U/mL) | 37.1-100 (U/mL) | 100.1-500 (U/mL) | 500.1-1200 (U/mL) | >1200 (U/mL) |
| Apparently Healthy Subjects | 360                | 94.4          | 5.6             | 0.0              | 0.0               | 0.0          |

In this study, 94.4% of the specimens from apparently healthy subjects (n=360) had values of 37 U/mL or less.

\* Representative data; results in individual laboratories may vary from these data.

## NONMALIGNANT DISEASE

A study was performed with four hundred forty one (441) samples from patients with nonmalignant disease to determine the distribution of serum ARCHITECT CA 19-9XR assay values. The distribution of values determined in this study is shown in the table below.\*

| Nonmalignant Disease   | Number of Subjects | Percent (%)   |                 |                  |                   |              |
|------------------------|--------------------|---------------|-----------------|------------------|-------------------|--------------|
|                        |                    | 0-37.0 (U/mL) | 37.1-100 (U/mL) | 100.1-500 (U/mL) | 500.1-1200 (U/mL) | >1200 (U/mL) |
| Rectal Polyps          | 33                 | 97.0          | 3.0             | 0.0              | 0.0               | 0.0          |
| Pancreatitis           | 3                  | 100.0         | 0.0             | 0.0              | 0.0               | 0.0          |
| Gallbladder            | 21                 | 95.2          | 0.0             | 0.0              | 0.0               | 4.8          |
| Diabetes               | 38                 | 94.7          | 5.3             | 0.0              | 0.0               | 0.0          |
| Pulmonary              | 40                 | 100.0         | 0.0             | 0.0              | 0.0               | 0.0          |
| Cirrhosis              | 153                | 92.8          | 4.6             | 0.7              | 0.7               | 1.3          |
| Hepatitis              | 68                 | 92.6          | 7.4             | 0.0              | 0.0               | 0.0          |
| Renal                  | 34                 | 91.2          | 8.8             | 0.0              | 0.0               | 0.0          |
| Other Gastrointestinal | 51                 | 96.1          | 3.9             | 0.0              | 0.0               | 0.0          |

The ARCHITECT CA 19-9XR assay is used in conjunction with other clinical methods in the management of cancer patients. It is recommended that each laboratory establish its own reference value for the population of interest.

\* Representative data; results in individual laboratories may vary from these data.

### Monitoring of Disease State in Patients Diagnosed with Pancreatic Cancer

Changes observed in serial CA 19-9 assay values when monitoring pancreatic cancer patients must be evaluated in conjunction with other clinical methods.

The effectiveness of the ARCHITECT CA 19-9XR assay as an aid in monitoring of disease state in pancreatic cancer patients was determined by assessing changes in levels of 1116-NS-19-9 reactive determinants in serial serum samples from 74 patients compared to changes in disease state. A study involving a total of 261 observations was performed with an average number of 3.5 observations per patient. In this study a significant change in levels of 1116-NS-19-9 reactive determinants was defined as at least a 14.0% increase in assay value (i.e., 2.5 times greater than the average of the assay's observed total %CV [5.6%]). A 14.0% change represents the minimum magnitude change between two serial ARCHITECT CA 19-9XR measurements that could not be attributed to assay variation or noise. Positive concordance between serial samples with at least a 14.0% increase in assay value and disease progression was found to be 48% (16/33). Negative concordance between serial samples with less than a 14.0% increase in assay value and no disease progression was found to be 64% (98/154). The overall concordance was found to be 61% (114/187). The following table presents the data in a 2 x 2 classification scheme\*.

| Change in the Level of 1116-NS-19-9 Reactive Determinants | Change in Disease State per Sequential Pair |                | Total |
|---|---|----------------|-------|
|   | Progression                                 | No Progression |       |
| ≥14.0%  | 16  | 56             | 72    |
| <14.0%  | 17  | 98             | 115   |
| Total   | 33  | 154            | 187   |

The following table provides the per patient distribution\*. Positive concordance between serial samples with at least a 14.0% increase in assay value and disease progression was found to be 68% (15/22). Negative concordance between serial samples with less than a 14.0% increase in assay value and no disease progression was found to be 69% (36/52). The overall concordance was found to be 69% (51/74).





| Change in Disease State per Patient                       |             |                |       |
|---|-------------|----------------|-------|
| Change in the Level of 1116-NS-19-9 Reactive Determinants | Progression | No Progression | Total |
| ≥14.0%  | 15          | 16             | 31    |
| <14.0%  | 7           | 36             | 43    |
| Total   | 22          | 52             | 74    |

\* Representative data; results in individual laboratories may vary from these data.

Below are examples of serial monitoring profiles for two patients with the disease state, ARCHITECT CA 19-9XR assay values, and the CA 19-9 RIA values.\* The disease states are:

- Progression from one collection to the next collection (Progression).
- No Change in disease state (Stable).
- Reduction in the signs and symptoms of the disease from one collection to the next (Responding).



\* Representative data; results in individual laboratories may vary from these data.

## SPECIFIC PERFORMANCE CHARACTERISTICS

### Precision

The ARCHITECT CA 19-9XR assay is designed to have an assay precision of ≤10% total CV.

A study was performed as described per the National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2.<sup>44</sup> Six samples were tested consisting of two panels of pooled serum (panels 1 and 2), one panel of serum to which 1116-NS-19-9 reactive determinants were added (panel 3), and the three ARCHITECT CA 19-9XR Controls. Testing was performed using two lots of reagents, in replicates of two at two separate times per day for 20

days on two separate instruments. Each reagent lot used a single calibration curve throughout the study. Data from this study are summarized below.\*

| Sample         | Reagent |            | n  | Mean Conc. (U/mL) | Within Run |     | Total |     |
|----------------|---------|------------|----|-------------------|------------|-----|-------|-----|
|                | Lot     | Instrument |    |                   | SD         | %CV | SD    | %CV |
| Panel 1        | 1       | 1          | 80 | 56.52             | 1.69       | 3.0 | 2.19  | 3.0 |
|                | 2       | 2          | 80 | 51.20             | 1.80       | 3.5 | 2.10  | 4.1 |
| Panel 2        | 1       | 1          | 80 | 311.49            | 7.22       | 2.3 | 10.72 | 3.4 |
|                | 2       | 2          | 80 | 288.82            | 9.14       | 3.2 | 11.23 | 3.9 |
| Panel 3        | 1       | 1          | 80 | 744.81            | 27.82      | 3.7 | 36.95 | 5.0 |
|                | 2       | 2          | 80 | 728.82            | 42.53      | 5.8 | 47.66 | 6.5 |
| Low Control    | 1       | 1          | 80 | 45.03             | 2.59       | 5.8 | 2.98  | 6.6 |
| Control        | 2       | 2          | 80 | 42.33             | 2.94       | 6.9 | 3.60  | 8.5 |
| Medium Control | 1       | 1          | 80 | 157.66            | 5.99       | 3.8 | 8.52  | 5.4 |
| Control        | 2       | 2          | 80 | 146.93            | 6.26       | 4.3 | 8.14  | 5.5 |
| High Control   | 1       | 1          | 80 | 781.68            | 44.76      | 5.7 | 49.87 | 6.4 |
| Control        | 2       | 2          | 80 | 781.42            | 62.10      | 8.0 | 65.28 | 8.4 |

\* Representative data; results in individual laboratories may vary from these data.

### Recovery

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 15% when 1116-NS-19-9 reactive determinants are added to serum samples.

A study was performed for the ARCHITECT CA 19-9XR assay based on guidance from Tietz Textbook of Clinical Chemistry.<sup>45</sup> Known concentrations of 1116-NS-19-9 reactive determinants were added to human serum samples. The concentration of 1116-NS-19-9 reactive determinants was determined using the ARCHITECT CA 19-9XR assay, and the resulting percent recovery was calculated. Representative data from this study are summarized in the table below.\*

| Sample | Endogenous Assay Value (U/mL) | 1116-NS-19-9 Reactive Determinants Added (U/mL) | Observed ARCHITECT CA 19-9XR Assay Value (U/mL) | % Recovery** |
|--------|-------------------------------|---|---|--------------|
| 1      | 48.50                         | 124.21  | 152.42  | 89           |
|        |                               | 629.91  | 645.00  | 95           |
| 2      | 26.36                         | 124.21  | 145.73  | 96           |
|        |                               | 629.91  | 598.93  | 91           |
| 3      | 38.42                         | 124.21  | 175.18  | 108          |
|        |                               | 629.91  | 652.12  | 98           |

Mean recovery across two separate spiked concentrations shown above = 96 %

$$** \% \text{ Recovery} = \frac{\text{Observed (U/mL)}}{\text{Endogenous Level (U/mL)} + \frac{\text{1116-NS-19-9 Reactive Determinants Added (U/mL)}}{2}} \times 100$$

\* Representative data; results in individual laboratories may vary from these data.

### Dilution Linearity

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 15% of the expected result for diluted specimens. A study was performed for the ARCHITECT CA 19-9XR assay modeled after the National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP6-P2.<sup>46</sup> Samples with known elevated 1116-NS-19-9 reactive determinant concentrations were diluted with ARCHITECT Multi-Assay Manual Diluent. The 1116-NS-19-9 reactive determinants concentration was determined for each dilution and the percent recovery was calculated. Representative data from this study are summarized below.\*



| Sample | Final Dilution Factor | Expected Value (U/mL) | Value Obtained (U/mL) | % Recovery** |
|--------|-----------------------|-----------------------|-----------------------|--------------|
| 1      | Undiluted             | 1024.55               | 1024.55               | -            |
|        | 1:2                   | 512.27                | 472.46                | 92           |
|        | 1:4                   | 256.14                | 264.26                | 103          |
|        | 1:5                   | 204.91                | 208.57                | 102          |
|        | 1:10                  | 102.45                | 108.94                | 106          |
|        | 1:20                  | 51.23                 | 54.33                 | 106          |
| 2      | Undiluted             | 1150.50               | 1150.50               | -            |
|        | 1:2                   | 575.25                | 551.62                | 96           |
|        | 1:4                   | 287.63                | 291.06                | 101          |
|        | 1:5                   | 230.10                | 253.65                | 110          |
|        | 1:10                  | 115.05                | 125.97                | 109          |
|        | 1:20                  | 57.53                 | 62.57                 | 109          |
| 3      | Undiluted             | 1028.25               | 1028.25               | -            |
|        | 1:2                   | 514.12                | 492.39                | 96           |
|        | 1:4                   | 257.06                | 290.24                | 113          |
|        | 1:5                   | 205.65                | 204.03                | 99           |
|        | 1:10                  | 102.82                | 120.76                | 117          |
|        | 1:20                  | 51.41                 | 57.25                 | 111          |

Mean recovery across the three diluted samples shown above = 105%

$$** \% \text{ Recovery} = \frac{\text{Values Obtained} \times \text{Dilution Factor}}{\text{Undiluted Expected Value}} \times 100$$

\* Representative data; results in individual laboratories may vary from these data.

#### Analytical Sensitivity

The analytical sensitivity of the ARCHITECT CA 19-9XR assay was calculated to be better than 2.00 U/mL (n = 18 runs, in replicates of 10).

Analytical sensitivity is defined as the concentration at two standard deviations from the ARCHITECT CA 19-9XR Calibrator A (0 U/mL), and represents the lowest measurable concentration of 1116-NS-19-9 reactive determinants that can be distinguished from zero.

#### Interference

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 12% in the presence of the chemotherapeutic agents listed below and elevated levels of bilirubin, hemoglobin, triglycerides, and total protein at the levels indicated.

A study based on guidance from the NCCLS Protocol EP7-A<sup>47</sup> was performed for the ARCHITECT CA 19-9XR assay. Specimens with 1116-NS-19-9 reactive determinant levels between 49.6 and 509.4 U/mL were supplemented with the following potentially interfering substances and chemotherapeutic agents.

#### POTENTIALLY INTERFERING SUBSTANCES

The average recovery observed during the study ranged from 91% to 102%.\*

| Substance     | Concentration |
|---------------|---------------|
| Bilirubin     | 22 mg/dL      |
| Hemoglobin    | 600 mg/dL     |
| Total Protein | 10 g/dL       |
| Triglycerides | 5100 mg/dL    |

\* Representative data; results in individual laboratories may vary from these data.

#### CHEMOTHERAPEUTIC AGENTS

The average recovery observed during the study ranged from 95% to 104%.\*

| Substance        | Concentration |
|------------------|---------------|
| 5-Fluorouracil   | 0.390 mg/mL   |
| Cisplatin        | 0.057 mg/mL   |
| Cyclophosphamide | 0.375 mg/mL   |
| Cytarabine       | 30 µg/mL      |
| Doxorubicin      | 40 µg/mL      |
| Gemcitabine      | 0.382 mg/mL   |
| Leucovorin       | 0.114 mg/mL   |
| Methotrexate     | 0.909 mg/mL   |
| Paclitaxel       | 0.067 mg/mL   |
| Streptozotocin   | 0.28 mg/mL    |
| Tamoxifen        | 2.28 µg/dL    |

\* Representative data; results in individual laboratories may vary from these data.

#### EVALUATION OF POTENTIALLY INTERFERING CLINICAL CONDITIONS

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 12% in the presence of HAMA and rheumatoid factor (RF).

The ARCHITECT CA 19-9XR assay was evaluated using specimens with HAMA and RF to further assess the clinical specificity. Five specimens positive for HAMA and five specimens positive for RF were evaluated for % recovery with 1116-NS-19-9 reactive determinants spiked into each specimen at 35 and 250 U/mL. Mean percent recovery results are summarized in the following table.\*

| Clinical Condition | Number of Specimens | Mean % Recovery |
|--------------------|---------------------|-----------------|
| HAMA               | 10                  | 93              |
| RF                 | 10                  | 93              |

\* Representative data; results in individual laboratories may vary from these data.

#### Carryover

No significant carryover (less than 2.00 U/mL in CA19-9XR Calibrator A\*) was observed for the ARCHITECT CA 19-9XR assay when a sample containing up to 320,000 U/mL of 1116-NS-19-9 reactive determinants was assayed.

\* Representative data; results in individual laboratories may vary from these data.

#### High Dose Hook

No high dose hook effect was observed for the ARCHITECT CA 19-9XR assay when samples containing up to 1,750,000 U/mL\* of 1116-NS-19-9 reactive determinants were assayed. High dose hook is a phenomenon whereby very high level specimens may falsely read within the dynamic range of the assay.

\* Representative data; results in individual laboratories may vary from these data.



## BIBLIOGRAPHY

- Hotakainen K, Tanner P, Alfthan H, et al. Comparison of three immunoassays for CA 19-9. *Clin Chim Acta* 2009;400:123-127.
- Passerinin R, Cassatella M, Boveri S, et al. Routine Testing and Comparison of Two Automated Immunoassays in a Reference Oncology Center. *Am J Clin Pathol* 2012;136:281-287.
- Uhlenbruck G, Höller U, Heising J, et al. Sialylated Le<sup>a</sup> Blood Group Substances Detected by the Monoclonal Antibody CA 19-9 in Human Seminal Plasma and Other Organs. *Urol Res (Germany West)* 1985;13(15):223-226.
- Herlyn M, Stepiewski Z, Herlyn D, et al. Colorectal Carcinoma-Specific Antigen: Detection by Means of Monoclonal Antibodies. *Proc Natl Acad Sci USA* 1979;76:1438-1442.
- Magnani JL, Brockhaus M, Smith DF, et al. A Monosialoganglioside is a Monoclonal Antibody-Defined Antigen of Colon Carcinoma. *Science* 1981;212:55-56.
- Magnani J, Nilsson B, Brockhaus M, et al. The Antigen of a Tumor-Specific Monoclonal Antibody is a Ganglioside Containing Sialylated Lacto-N-Fucopentaose II. *Fed Proc* 1982;41:998.
- Magnani JL, Stepiewski Z, Koprowski H, et al. Identification of the Gastrointestinal and Pancreatic Cancer-Associated Antigen Detected by Monoclonal Antibody 19-9 in the Sera of Patients as a Mucin. *Cancer Res* 1983;43:5489-5492.
- Del Villano BC, Brennan S, Brock P, et al. Radioimmunoassay for a Monoclonal Antibody-Defined Tumor Marker, CA 19-9. *Clin Chem* 1983;29:549-552.
- Steinberg WM, Gelfand R, Anderson KK, et al. Comparison of the Sensitivity and Specificity of the CA 19-9 and Carcinoembryonic Antigen Assays in Detecting Cancer of the Pancreas. *Gastroent* 1986;90:343-349.
- Ritts RE Jr, Del Villano BC, Go VLW, et al. Initial Clinical Evaluation of an Immunoradiometric Assay for CA 19-9 Using the NCI Serum Bank. *Int J Cancer* 1984;33:339-345.
- Jalanko H, Kuusela P, Roberts P, et al. Comparison of a New Tumor Marker, CA 19-9, with Alpha-Fetoprotein and Carcinoembryonic Antigen in Patients with Upper Gastrointestinal Diseases. *J Clin Pathol* 1984;37:218-222.
- Gupta MK, Arciaga R, Bocci L, et al. Measurement of a Monoclonal-Antibody-Defined Antigen (CA 19-9) in the Sera of Patients with Malignant and Nonmalignant Diseases. Comparison with Carcinoembryonic Antigen. *Cancer* 1985;56(2):277-283.
- Andriulli A, Gindro T, Plantino P, et al. Prospective Evaluation of the Diagnostic Efficacy of CA 19-9 Assay as a Marker for Gastrointestinal Cancers. *Digestion* 1986;33(1):26-33.
- Patali A, Héber S, Dörönte Z, et al. Diagnostic Values of CA 19-9 and CEA in Gastrointestinal Diseases. *Orv Hetil* 1992;133(21):1301-1304, 1307.
- Ritts RE and Pitt HA. CA 19-9 in Pancreatic Cancer. *Surgical Oncology Clinics of North America* 1998;7:93-101.
- Duffy MJ, van Dalen A, Haglund C, et al. Clinical Utility of Biochemical Markers in Colorectal Cancer: European Group on Tumor Markers (EGTM) Guidelines. *Eur J Cancer* 2003;39:718-727.
- Cerwenka H, Aigner R, Quehenberger F, et al. Preoperative Differential Diagnosis of Benign and Malignant Pancreatic Lesions - The Value of Pancreatic Secretory Trypsin Inhibitor, Procarboxypeptidase B, CA 19-9 and CEA. *Hepato-Gastroenterology* 1997;44(10):1117-1121.
- Von Ritter C, Eder M, Stieber P, et al. Biliary Mucin Secreted by Cultured Human Gallbladder Epithelial Cells Carries the Epitope of CA 19-9. *Anticancer Res* 1987;7(4B):2931-2934.
- Adachi Y, Iso Y, Moriyama M, et al. Increased Serum CA 19-9 in Patients with Xanthogranulomatous Cholecystitis. *Hepato-Gastroenterology* 1998;45:77-80.
- Maestranzi S, Przemioslo R, Mitchell H, et al. The Effect of Benign and Malignant Liver Disease on the Tumor Markers CA 19-9 and CEA. *Ann Clin Biochem* 1998;35:99-103.
- Duffy MJ, O'Sullivan F, McDonnell T, et al. Increased Concentrations of the Antigen CA 19-9 in Serum of Cystic Fibrosis Patients. *Clin Chem* 1985;31:1245-1246 (Letter).
- Kane RE, Penny J, Walker K, et al. Changes in the CA 19-9 Antigen and Lewis Blood Group with Pulmonary Disease Severity in Cystic Fibrosis. *Pediatr Pulmonol* 1992;12(4):221-226.
- Wu JT and Chang J. Chromatographic Characterization of CA 19-9 Molecules from Cystic Fibrosis and Pancreatic Carcinoma. *J Clin Lab Anal* 1992;6(4):209-215.
- Wu JT, Olsen J, Walker K. Tumor Markers CA 19-9 and CA 195 Are Also Useful as Markers for Cystic Fibrosis. *J Clin Lab Anal* 1992;5(3):151-161.
- Staab HJ, Brümmendorf T, Hornung A, et al. The Clinical Validity of Circulating Tumor-Associated Antigens CEA and CA 19-9 in Primary Diagnosis and Follow-up of Patients with Gastrointestinal Malignancies. *Klin Wochenschr* 1985;63(3):106-115.
- Ychou M, Tuszinski T, Pignon J-P, et al. Gastric Carcinoma: Comparison between CEA and CA 19-9 for diagnosis of recurrence after gastrectomy. *Gastroenterol Clin Biol* 1992;16(11):848-852.
- Grem J. The Prognostic Importance of Tumor Markers in Adenocarcinomas of the Gastrointestinal Tract. *Curr Opin Oncol* 1997;9(4):380-387.
- Gärtner U, Scheulen ME, Conrad C, et al. Value of Tumor-Associated Antigen CA 72-4 Compared with CEA and CA 19-9 in the Followup of Patients Operated for Gastric Carcinoma. *Dischimed Wschr* 1998;123:69-73.
- Willett CG, Daly WJ, Warshaw AL. CA 19-9 Is an Index of Response to Neoadjuvant Chemoradiation Therapy in Pancreatic Cancer. *Am J Surg* 1996;172(4):350-352.
- Fillela X, Molina R, Grau JJ, et al. Prognostic Value of CA 19.9 Levels in Colorectal Cancer. *Ann Surg* 1992;216(1):55-59.
- Kouri M, Pyrhönen S, Kuusela P. Elevated CA 19-9 as the Most Significant Prognostic Factor in Advanced Colorectal Carcinoma. *J Surg Oncol* 1992;49(2):78-85.
- Gebauer G, Müller-Ruchholtz W. Tumor Marker Concentrations in Normal and Malignant Tissues of Colorectal Cancer Patients and Their Prognostic Relevance. *Anticancer Res* 1997;17(4a):2731-2734.
- Reiter W, Stieber P, Reuter C, et al. Preoperative Serum Levels of CEA and CA 19-9 and Their Prognostic Significance in Colorectal Carcinoma. *Anticancer Res* 1997;17(4B):2935-2938.
- Gogas H, Lolts FJ, Evans TRJ, et al. Are Serial Measurements of CA 19-9 Useful in Predicting Response to Chemotherapy in Patients with Inoperable Adenocarcinoma of the Pancreas? *Br J Cancer* 1998;77:325-328.
- Safi F, Schlosser W, Falkenreck S, et al. Prognostic Value of CA 19-9 Serum Course in Pancreatic Cancer. *Hepato-Gastroenterol* 1998;45:253-259.
- US Department of Labor, Occupational Safety and Health Administration, 28 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Tempero MA, Uchida E, Takasaki H, et al. Relationship of Carbohydrate Antigen 19-9 and Lewis Antigens in Pancreatic Cancer. *Cancer Res* 1987;47:5501-5503.
- National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP6-A7. Wayne, PA: NCCLS; 2004.
- Tietz NW. *Tietz Textbook of Clinical Chemistry*. 3rd Edition. Philadelphia, PA: B. Saunders Company, 1999:330.
- National Committee for Clinical Laboratory Standards. *Evaluation of the Linearity of Quantitative Analytical Methods; Proposed Guideline—Second Edition*. NCCLS Document EP6-P2. Wayne, PA: NCCLS; 2001.
- National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Approved Guidelines*. NCCLS Document EP7-A. Wayne, PA: NCCLS; 2002.



## Key to Symbols

|   |   |
|---|---|
|  | Consult instructions for use  |
|  | Manufacturer  |
|  | Sufficient for  |
|  | Temperature limitation  |
|  | Use by/Expiration date  |
| <b>CONJUGATE</b>  | Conjugate   |
| <b>CONTAINS: AZIDE</b>  | Contains Sodium Azide. Contact with acids liberates very toxic gas. |
| <b>CONTROL NO.</b>  | Control Number  |
| <b>ECC HAZARD</b>   | Ecological hazard   |
| <b>IVD</b>  | <i>In Vitro</i> Diagnostic Medical Device                           |
| <b>LOT</b>  | Lot Number  |
| <b>MICROPARTICLES</b>   | Microparticles  |
| <b>MULTI-ASSAY MANUAL DILUENT</b>   | Multi-Assay Manual Diluent  |
| <b>PRE-TRIGGER SOLUTION</b>   | Pre-Trigger Solution  |
| <b>PRODUCED FOR ABBOTT BY</b>   | Produced for Abbott by  |
| <b>PRODUCT OF USA</b>   | Product of USA  |
| <b>REACTION VESSELS</b>   | Reaction Vessels  |
| <b>REAGENT LOT</b>  | Reagent Lot   |
| <b>REF</b>  | List Number   |
| <b>REPLACEMENT CAPS</b>   | Replacement Caps  |
| <b>SAMPLE CUPS</b>  | Sample Cups   |
| <b>SEPTUM</b>   | Septum  |
| <b>SN</b>   | Serial number   |
| <b>TRIGGER SOLUTION</b>   | Trigger Solution  |
| <b>WARNING: SENSITIZER</b>  | Warning: May cause an allergic reaction.                            |
| <b>WASH BUFFER</b>  | Wash Buffer   |

ARCHITECT and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.



Abbott GmbH & Co. KG  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580



**PRODUCED FOR ABBOTT BY**

Fujirebio Diagnostics Inc., Malvern, PA 19355 USA

**Customer Service:** Contact your local representative or find country-specific contact information on [www.abbottiagnostics.com](http://www.abbottiagnostics.com)

Revised November 2016.

©2005, 2016 Abbott Laboratories





# ARCHITECT CA 125 II Calibrators

REF 2K45-01



**en**  
CA 125 II  
2K45  
602-001 8/14/R11  
S2K450

Read Highlighted Changes: Revised November 2014.

## INTENDED USE

The ARCHITECT CA 125 II Calibrators are for the calibration of the ARCHITECT iSystem when used for the quantitative determination of OC 125 defined antigen in human serum and plasma.

Refer to the ARCHITECT CA 125 II reagent package Insert for additional information.

## CONTENTS

6 Bottles (4 mL each) of ARCHITECT CA 125 II Calibrators. Calibrator A contains TRIS buffer with protein (bovine) stabilizers. Calibrators B - F contain OC125 defined antigen (human) prepared in TRIS buffer with protein (bovine) stabilizers. Preservatives: Sodium Azide and ProClin 300.

The calibrators yield the following concentrations:


| Calibrators | CA 125 II Concentration<br>U/mL |
|-------------|---------------------------------|
| CAL A       | 0                               |
| CAL B       | 20                              |
| CAL C       | 75                              |
| CAL D       | 225                             |
| CAL E       | 500                             |
| CAL F       | 1000                            |


## STANDARDIZATION

CA 125 assay values are expressed as U/mL. A unit is a value related to a Fujirebio Diagnostics, Inc. maintained reference preparation. The calibrators for the ARCHITECT CA 125 II assay are manufactured volumetrically and are referenced to this standard prepared by Fujirebio Diagnostics, Inc. There is no internationally recognized CA 125 standard available at this time.

## PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Calibrators B-F contain antigen derived from a human cell line. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.<sup>1-4</sup>


| The following warnings and precautions apply to: CAL A CAL F                      |  |
|---|--|
|  |  |
| <b>WARNING</b>  | Contains metylisothiazolones and sodium azide.                         |
| H317  | May cause an allergic skin reaction.                                   |
| EUH052  | Contact with acids liberates very toxic gas.                           |
| <b>Prevention</b>   |  |
| P261  | Avoid breathing mist / vapors / spray.                                 |
| P272  | Contaminated work clothing should not be allowed out of the workplace. |
| P280  | Wear protective gloves / protective clothing / eye protection.         |
| <b>Response</b>   |  |
| P302+P352   | IF ON SKIN: Wash with plenty of water.                                 |
| P333+P313   | If skin irritation or rash occurs: Get medical advice / attention.     |
| P362+P364   | Take off contaminated clothing and wash it before reuse.               |
| <b>Disposal</b>   |  |
| P501  | Dispose of contents / container in accordance with local regulations.  |

Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

## STORAGE

- Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.



## PREPARATION FOR ANALYSIS

- Calibrators may be used immediately after removal from 2-8°C storage.
- Prior to use, mix by gentle inversion (5-10 times).
- After each use, tightly close the caps and return the calibrators to 2-8°C storage.

## BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. International and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.



**Key to Symbols**

|                                 |   |
|---------------------------------|---|
|                                 | Caution   |
|                                 | Consult instructions for use  |
|                                 | Manufacturer  |
|                                 | Temperature limitation  |
|                                 | Use by/Expiration date  |
| <b>CAL A</b>                    | Calibrator (A,B,C,D,E or F)   |
| <b>CONTAINS AZIDE</b>           | Contains Sodium Azide. Contact with acids liberates very toxic gas. |
| <b>EC REP</b>                   | Authorized Representative in the European Community                 |
| <b>INFORMATION FOR USA ONLY</b> | Information needed for United States of America only                |
| <b>IVD</b>                      | In Vitro Diagnostic Medical Device                                  |
| <b>LOT</b>                      | Lot Number  |
| <b>PRODUCED FOR ABBOTT BY</b>   | Produced for Abbott by  |
| <b>PRODUCT OF USA</b>           | Product of USA  |
| <b>REF</b>                      | List Number   |
| <b>WARNING: SENSITIZER</b>      | Warning: May cause an allergic reaction                             |

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Laboratories  
Diagnostics Division  
Abbott Park, IL 60064  
USA

**EC REP** ABBOTT  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580



**PRODUCED FOR ABBOTT BY**

Fujirebio Diagnostics Inc., Malvern, PA 19355 USA

**Customer Service: Contact your local representative or find country-specific contact information on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)**

Revised November 2014.

©2003, 2014 Abbott Laboratories

