Control-Test M

REF	CONTENT	SYSTEM
11379194263	▽ 50	cobas u 411 urine analyzer Urisys 1100®

English

Intended use

The Control-Test M calibration strip is used to perform the calibration of the **cobas u 411** urine analyzer and the Urisys 1100®.

For professional use including near-patient testing.

Not for self-testing.

Summary

The ready-to-use calibration strip is made of inert gray plastic material of constant reflectance characteristics. The calibration values were determined with reference to a white standard and are stored in the software of the respective analyzer.

For Urisys 1100® users in near-patient testing environments: The Control-Test M can be used in non-critical healthcare environments such as hospital wards and general practitioner's offices.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

Do not immerse the calibration strip in urine. The strips are designed for single use only.

Reagent handling

Calibration strips are ready for use.

Storage and stability

Store the package at 2-30 °C. The test strips are stable up to the expiration date specified on the box, when stored in the original container.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Store the tube in a horizontal position.

Materials provided

- [REF] 11379194263, package with 50 strips

Materials required (but not provided)

- cobas u 411** urine analyzer or Urisys 1100®
- General laboratory equipment

Assay

Refer to the Operators Manual of the respective instrument for full details on calibration and performance checks.

- Take a calibration strip out of the vial. Do not touch the elevated areas of the strip to avoid contamination and falsification of results.
- Insert the strip into the instrument. Follow the instructions in the Operators Manual for correct positioning of the strip.
- Follow the Operators Manual to execute the calibration. Calibration interval: 1 week with Urisys 1100®, 4 weeks with **cobas u 411** urine analyzer.

Calibration

The instrument measures the elevated areas of the strip by means of reflectance photometry. The results are automatically calculated and printed by the instrument as reflectance values. For interpretation of these values, please refer to the relevant section in the Operators Manual. If the values obtained deviate from the tolerances specified, the calibration has to be repeated. To ensure the correctness of the calibration a fresh calibration strip has to be used for this step. This permits the user to check whether the instrument is malfunctioning or whether the strip used is dirty or damaged.

For further information, please refer to the appropriate Operators Manual for the analyzer concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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