



MAGLUMI[®] Light Check

INTENDED USE

The Light Check is used for completing system test to monitor the status of MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer and Biolumi series Integrated System.

SUMMARY

The ABEI in the Light Check reacts with the Starter 1+2 to generate a certain light intensity, this process can be used to check the system is working properly to avoid false results.

MATERIALS PROVIDED

REF	Description	Specification
130299006M	ABEI (N-(4-Aminobutyl)-N-ethylisoluminol) BSA	5x2 mL

The Light Check is ready-for-use.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- For professional use only.
- Do not use the product beyond the expiration date indicated on the label.
- Exercise the normal precautions required for handling all laboratory reagents.
- Avoid ingestion and contact with skin and eyes, if you unfortunately experience it, please drink or rinse with plenty of water, and seek emergency medical treatment.
- All waste should be considered potentially infectious and should be disposed of in accordance with local guidelines.
- Safety data sheets are available on request.
- Note: If any serious incident has occurred in relation to the device, please report to Shenzhen New Industries Biomedical Engineering Co., Ltd. (Snibe) or our authorized representative and the competent authority of the Member State in which you are established.

PROCEDURE

Materials Provided

Light Check.

Materials Required (But Not Provided)

Fully-auto chemiluminescence immunoassay analyzer Maglumi 600, Maglumi 800, Maglumi 1000, Maglumi 1000 Plus, Maglumi 2000, Maglumi 2000 Plus, Maglumi 4000, Maglumi 4000 Plus, MAGLUMI X3, MAGLUMI X6, MAGLUMI X8, or Integrated System Biolumi 8000 and Biolumi CX8.

Handling

- Carefully open the vial containing the liquid material.
- Prior to use, turn upside down slowly and let it stand for at least 5 minutes at room temperature.
- Foam must be avoided and cleaned before use.
- Place the Light Check vial into the matched patient rack; make the barcode label of the vial face the barcode reader of the sample area (refer to operation instructions of MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer and Biolumi series Integrated System for detailed information).
- Seal the vial and store the rest of the Light Check at 2-8°C immediately after been used.

Target values for Light Check













- RLU: 450,000-650,000.
- Target values may differ due to analyzer setting; contact SNIBE authorized technician for further information

STORAGE AND STABILITY

- Unopened at 2-8°C: until the expiry date as stated on the labels.
- Opened at 2-8°C: 4 weeks.
- Protect from direct sunlight.

SYMBOLS EXPLANATIONS


	Consult instructions for use		Manufacturer
--	------------------------------	--	--------------

	Temperature limit (Store at 2-8°C)		Keep away from sunlight
	Catalogue number		Use-by date
	<i>In vitro</i> diagnostic medical device		Batch code
	Authorized representative in the European Community		Kit component
	CE marking		This way up
	Unique device identifier		Importer

MAGLUMI® and Biolumi® are trademarks of Snibe. All other product names and trademarks are the property of their respective owners.



Shenzhen New Industries Biomedical Engineering Co., Ltd.
No.23, Jinxiu East Road, Pingshan District, 518122 Shenzhen, P.R. China
Tel: +86-755-21536601 Fax: +86-755-28292740



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175 Fax: +49-40-255726