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ALBUMIN (MICROALBUMINURIA)

COD 13324 1 x 50 mL
Only for <i>in vitro</i> use in the clinical laboratory

ALBUMIN (MICROALBUMINURIA)
LATEX

INTENDED USE

Reagents for the quantitative measurement of albumin concentration in human urine as an aid in the diagnosis of the disorders related to the renal glomerular permeability in the general population.

This reagent is for use in the BioSystems A25 and A15 analyzers. For *in vitro* professional use only in the clinical laboratory.

CLINICAL BENEFIT

Urinary albumin concentration values provide a good indicator of changes in glomerular permeability, as occur in a number of renal disease¹.

Diabetic nephropathy is characterized by an early hyper filtration stage resulting in small increases in urinary albumin excretion. That is why the measurement of urinary albumin is considered a clinically important indicator of deteriorating renal function in diabetic subjects.

Urinary albumin excretion is also monitored in hypertensive patients to identify the development of significant nephropathy.

Based on clinical guidelines and textbooks, and when used in conjunction with other diagnostic technologies and options, this medical information is useful for the assessment of albumin imbalance. Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Albumin in the urine sample causes agglutination of the latex particle coated with antibody. The agglutination of the particles is proportional to the albumin concentration and can be measured by turbidimetry^{2,5}.

CONTENTS AND COMPOSITION

- A. Reagent: 1 x 40 mL. Borate buffer 0.1 mol/L, sodium azide 0.95 g/L, pH 10.0.
B. Reagent: 1 x 10 mL. Suspension of latex particles coated with anti-human albumin antibodies, sodium azide 0.95 g/L (Note 1).

STORAGE AND STABILITY

Store at 2-8°C.

Components are stable once opened until the expiry date marked in the label if they are kept at the recommended storage temperature, well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 30 days.

Indications of deterioration: Absorbance of the blank over the limit indicated in the parameterization of the analyzer.

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

S. Albumin Standard: 1 x 1 mL (BioSystems cod. 31130). Human albumin. Albumin concentration is given on the label. The concentration value is traceable to the ERM DA-470 reference standard (Institute for Reference Materials and Measurements, IRMM).

Human serum used in the preparation of the standard has been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, the standard should be handled cautiously as potentially infectious.

Reconstitute with 1.00 mL of distilled water. Stable for 30 days at 2-8°C and 3 months at -20°C. Avoid repeated freeze-thaw cycles.

REAGENT PREPARATION

Reagents are provided ready to use. R1: Use Reagent A, R2: Use Reagent B.

SAMPLES

Urine collected by standard procedures. Urine should be centrifuged before analysis. Albumin in urine is stable for 30 days at 4-8°C⁶.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 30 days, after reagent lot change or as required by quality control procedures. It is recommended to use the calibrators mentioned in the paragraph on Additional Materials Required.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Urine (cod. 18054 and cod. 18066) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Urine, adults⁷: 20 mg/L, <30 mg/24h.

This range is given for orientation only; each laboratory should establish its own reference range.

ANALYTICAL AND CLINICAL PERFORMANCE

The metrological characteristics described below have been obtained using an A25 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI). Results are similar with A15. Details of the evaluation data are available upon request.

- Detection limit: 1.67 mg/L. Quantification limit: 5.22 mg/L.
- Linearity limit: 200 mg/L. Measuring range: 5.22 - 200 mg/L. For samples with higher values, dilute manually or refer to the Test Parameterization for Automatic dilution (note that all these samples will be diluted with the same dilution ratio).
- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
24.4 mg/L	2.3 %	3.3 %
106 mg/L	1.2 %	2.1 %
147 mg/L	1.1 %	1.9 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: Bilirubin (up to 30 mg/dL) and hemolysis (hemoglobin up to 500 mg/dL) do not interfere. Other drugs and substances may interfere⁸.
- Zone effect: This method has no zone effect up to 700 mg/L.

NOTES

- Shake the Reagent B vial gently before using. The dark particles that can be observed in some batches come from the production process and do not affect the functionality of the product.

BIBLIOGRAPHY

- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
- Cambiaso CL, Collet-Cassart D, Lievens M. Immunoassay of low concentrations of albumin in urine by latex particle counting. Clin Chem 1988; 34(2):416-418
- Medcalf EA, Newman DJ, Gorman EG, Price CP. Rapid, robust method for measuring low concentrations of albumin in urine. Clin Chem 1990; 36(3):446-449
- Harmoinen A, Ala-Houhala I, Vuorinen P. Rapid and sensitive immunoassay for albumin determination in urine. Clin Chim Acta 1985 15;149(2-3):269-74
- Bernard A, Lauwers R. Latex immunoassay of urinary albumin. J Clin Chem Clin Biochem 1983; 21(1):25-30
- World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev.2; 2002.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 6th ed. Rifai N, Horvath AR, Wittwer CT. WB Saunders Co, 2017.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

TEST PARAMETERS

	A25	A15
GENERAL		
Name	MICROALB.	MICROALB.
Sample type	URI	URI
Analysis mode	fixed-time bir.	fixed-time bir.
Units	mg/L	mg/L
Turbidimetry test	yes	yes
Active mixing	yes	yes
Decimals	0	0
Type of reaction	increasing	increasing
PROCEDURE		
Reading mode	monoch.	monoch.
Main filter	535	535
Reference filter	-	-
Sample	3	3
Vol. R1	335	335
Vol. R2	85	85
Washing	1.2	1.2
Reading 1 (cycle)	4	3
Reading 2 (cycle)	10	8
Reagent 2 (cycle)	2	2
Predilution factor	-	-
Postdilution reduced factor	1.4	1.4
CALIBRATION AND BLANK		
Calibration type	specific	specific
Number of calibrators	1	1
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	1.500	1.500
Kinetic blank limit	-	-
Linearity limit	200	200
Substrate depletion	-	-

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BIOCHEMISTRY CONTROL URINE I, II

LEVEL I	LEVEL II
COD 18054 1 x 5 mL	COD 18066 1 x 5 mL
Only for <i>in vitro</i> use in the clinical laboratory	

INTENDED USE

The Biochemistry Control Urine is intended for quality control purposes only and is supplied with suggested intervals of acceptable values.

COMPOSITION

Biochemistry Control Urine. For 5 mL. Lyophilized human urine containing component concentrations suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests.

PREPARATION AND USE

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 5 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. Analyze the reconstituted material in the same way as the patient samples.

STORAGE AND STABILITY

Store at 2-8°C.

Biochemistry Control Urine is stable until the expiration date given in the label.

The components of the reconstituted material are stable for at least 7 days at 2-8°C. For longer periods, keep frozen. Avoid repeated freeze-thaw cycles.

SPECIFIED VALUES

The assigned concentration values for components and their traceability are shown in the enclosed value sheets. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

The suggested intervals of acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only. Each laboratory should establish its own precision parameters.

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U: Expanded uncertainty (95%; k=2)
NA: Non applicable.

Note: Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV (I and II) and anti-HCV as well as for HBs antigen. However the product should be handled cautiously as potentially infectious.

This lot meets the quality requirements, and it is released for sale.

Approved by: TECHNICAL RESPONSIBLE / DEPUTY TECHNICAL RESPONSIBLE Date: 2025/11/25

This certificate of analysis is a computer printout and has therefore not been signed by hand.

Notes

Concentration value is traceable to the ERM DA-470/IFCC Human Serum (Proteins). Institute for Reference Materials and Measurements (IRMM).

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