PIT-LLSL-4-v13 (01/2020)

Français - FR

USAGE PRÉVU

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English - EN

INTENDED USE

ELITech Clinical Systems LDH-L SL is an *in vitro* diagnostic reagent intended for the quantitative determination of Lactate dehydrogenase (LDH) in human serum and plasma samples.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾


Lactate dehydrogenase (LDH) can be found in nearly all cells of the body with highest activities in myocardium, liver, kidney and skeletal muscle. Consequently, elevations of LDH in the serum have been considered nonspecific (for any disease or disorder). LDH increases in case of acute myocardial infarction, hepatic disorders (viral hepatitis, cirrhosis), muscular dystrophy, cancer, metastasis, anemia (hemolytic, megaloblastic), kidney diseases and in numerous other diseases involving tissue damage.

☞METHOD ⁽⁴⁾

IFCC method
Kinetic.

PRINCIPLE ⁽⁴⁾

Determination of lactate dehydrogenase by measuring the increase of absorbance at 340 nm.

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- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

WARNINGS AND PRECAUTIONS

- This *in vitro* diagnostic device is for professional use only.
- The reagent R1 contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITIES

Store at 2-8 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent is ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES ^(2,3,5)

- Specimen**
- Serum free from hemolysis.
- Lithium heparinized plasma free from hemolysis.
- Do not use other specimens.

Warnings and precautions

- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.
- Samples must be separated from cells and clot promptly (the presence of cells can increase falsely the result).
- Storage**
Samples are stable at room temperature for up to 7 days, at 2-8 °C for 4 days and in frozen state (-20 °C) for up to 6 weeks.

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English - EN

INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of ELITech Clinical Systems quantitative tests listed in the value sheet. These *in vitro* diagnostic devices are for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Concentrations for each analyte to test are lot-specific.

MATERIALS REQUIRED BUT NOT PROVIDED

- ELITech Clinical Systems reagents listed in the value sheet.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated on the instructions for use of the ELITech Clinical Systems reagents used in combination with the recommended calibrator.

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contaminations.

- Consult Safety Data Sheet (SDS) for a proper handling.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.
- Carefully open the vial, avoiding the loss of lyophilizate.
- Add exactly 5 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

Important: To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).

- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.

- Do not use after the expiry date stated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.

- Stability of the components:

Between 2 and 8 °C : 5 days

Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light):
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of direct bilirubin (when stored protected from light):
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of alkaline phosphatase:

Between 2 and 8 °C : 24 hours

Between -25 and -15 °C : 4 weeks (when frozen once)

Note:

- Store the vials tightly capped and protected from light after reconstitution.

- Mix content thoroughly prior to use.

- Thawing frozen specimen at 37°C, and vigorous mixing (vortex), should be avoided as this may denature ferritin.

PROCEDURE

To use ELITROL I & II, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using ELITROL I & II has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet). Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit. Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values.

For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

ASSISTANCE TECHNIQUE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems reagents, used for C-reactive protein (CRP) quantitative determination. This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Liquid product prepared from human plasma spiked with Human recombinant CRP.
- Sodium azide < 0.1 % (w/w)
- C-reactive protein (CRP) concentrations are lot-specific.

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- ELITech Clinical Systems reagent, used for C-reactive protein (CRP) quantitative determination (Ref. ICRP-XXXX).
- General Laboratory equipment.

TRACEABILITY

Calibrators from CRP IP CALIBRATOR SET are traceable to ERM-DA474/IFCC reference material.

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

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These products are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

- Store at 2-8 °C and protect from light. Do not freeze.
- Do not use after the expiry date stated on the vial label.
- These products are stable 6 weeks at 2-8°C after the first opening.
- These products should be immediately and tightly capped to prevent contamination and evaporation.

PROCEDURE

To use CRP IP CALIBRATOR SET, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using CRP IP CALIBRATOR SET has been validated with the ELITech systems (Analyzers and reagents used for C-reactive protein (CRP) quantitative determination (Ref. ICRP-XXXX)).
Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European Union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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Entre em contato com seu distribuidor local ou com a ELITech Clinical Systems SAS. (CCsupport@elitechgroup.com).