INgezim® BTV COMPAC 2.0

Prod Ref: 12.BTV.K3

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KIT COMPOSITION

REAGENT	2 Plates		5 Plates		10 Plates	
	Uni.	Vol.	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5		10	
Vials with Positive Control	1	600µL	2	600µL	4	600µL
Vials with Negative Control	1	600µL	2	600µL	4	600µL
Vials with Peroxidase Conjugate, ready to use	1	30mL	2	30mL	2	60mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL	2	125mL
Bottles with Diluent (DE01-01), ready to use	1	125mL	1	125mL	2	125mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL	1	125mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL	1	125mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter)

TECHNICAL BASIS

The kit has been designed to detect antibodies specific for BTV in sheep, goats and cattle being able to detect a very low titers of antibodies in sera, and plasma samples of infected animals.

The INgezim® BTV kit is based on the blocking immunoenzymatic assay described below:

The solid phase is plates coated with VP7 protein of BTV. After adding the sample to the well, if it contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the sample does not contain specific antibodies, they will not bind to the antigen.

If we add a specific monoclonal antibody against the viral antigen coated to the plate (conjugated with peroxidase), it will compete with the antibodies of the sample. If the samples contain specific antibodies, they will not permit the binding of the labelled MAb to the antigen whereas if it does not contain specific antibodies the MAb will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled MAb by adding the substrate (TMB) that in presence of the peroxidase will develop a colorimetric reaction.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

Control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

They must be assayed diluted 1/2 in diluent. This dilution could be made directly on the assay plate by adding 50 µL of diluent and 50 µL of sample to each well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (i.e., 40 mL of concentrate and 960 mL of water). Once prepared, this solution remains stable between +2°C and +8°C.

Positive and Negative serum controls:

Controls must be tested as sample serum adding 50 µL/well over 50 µL of diluent.

Preparation of the conjugate:

Conjugate is ready to use and must be used adding 100 µL/well.

TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. **Addition of samples:** Add 50 µL of diluent to each well. Add 50 µL of each sample to be assayed. For controls add 50 µL as it is indicated in point V. We recommend to assay the samples by duplicate.
3. Seal the plate and incubate overnight (18-24 hours) at room temperature or alternatively **2.5 hours at 37°C**.
4. Wash 6 times following the described procedure.
5. Add 100 µL of conjugate ready to use to each well. Seal the plate and incubate for **30 min at 37°C**.
6. Wash 6 times following the described procedure.
7. Add 100 µL of substrate, to each well. Keep the plate for **10 min at room temperature**.
8. Add 100 µL of stop solution to each well.
9. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

In case that the samples had been run in duplicate, it has to be considered the mean of both OD values. In the same way, the mean of the values obtained in the two well of positive and the two well of negative control has to be made.

VALIDATION OF THE RESULTS:

- Positive Control / Negative Control < 0.25

CUT OFF CALCULATION:

- Cut Off (-) = 0.65 x Negative Control
- Cut Off (+) = 0.6 x Negative Control

BLOCKING % CALCULATION:

$$\text{Blocking \% of sample} = 100 - \frac{\text{Sample OD} \times 100}{\text{Negative control OD}}$$

RESULTS INTERPRETATION

- Samples will be considered **POSITIVE** (there are antibodies specific for BTV), when the OD at 450 nm was equal or lower than the positive cut off (60 % of negative control).
- Samples will be considered **NEGATIVE** (there are no antibodies specific for BTV in the sample) when the OD value at 450nm was equal or higher than the negative cut off (65% of negative control).
- Samples with percentages of binding between both cut off have to be considered **DOUBTFUL**. In those cases, a new sample of the animal is recommended to be analyzed in two weeks.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

9191.INGE
9175.INGE

Distribuido en
Distributed in

por
by