M-6LN LYSE

【PRODUCT NAME】

M-6LN LYSE

【MODEL】

M-6LN

【PACKAGE SPECIFICATIONS】

4L×1, 1L×4

【PRINCIPLE】
This product is formulated to lyse red blood cells in blood cell measurement, facilitate the dying effect of NRBCs by M-6FN DYE. It participates in the measurement of NRBC-related parameters which is accomplished by using SF Cube cell analysis technology (3D analysis using information from scatter of laser light at two angles and fluorescence signals).

【INTENDED USE】
This product participates in the measurement of NRBC-related parameters together with M-6FN DYE.

【ACTIVE INGREDIENTS】
Citrate Buffer.....≤0.5%
Sodium Chloride.....0.10%
Surfactant.....0.10%

【PRODUCT SPECIFICATIONS】
1. Appearance: Transparent liquid without deposits, suspended grains or flocks.
2. Blank count results: The blank count results of the reagent tested on the Mindray Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1 Blank Count Requirements	
Parameter	Blank Count Requirements
WBC-N	≤0.2×10 ⁹ /L

【STORAGE AND EXPIRATION DATE】
This product shall be stored in ventilated room free from corrosive gas at 2℃~30℃ with the humidity lower than 90%. The shelf life is 2 years.
It shall be used at 15℃~32℃ and within 60 days after being opened.

【APPLICABLE INSTRUMENT】
This product applies to BC-6000/BC-6000Plus/BC-6100/BC-6100Plus/BC-6200/BC-6600Plus/BC-6700Plus/BC-6800Plus Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

【INSTRUCTIONS FOR USE】
1) Restore the product to usage temperature;
2) Open the external packing and insert the pickup tube into the reagent container based on the matching color of the reagent cap and the connector of analyzer cap assembly;
3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis. See the operator's manual of the analyzer for details.

【WARNINGS】
Do not inhale and/or ingest. In case of inhalation or ingestion, seek medical advice immediately.
Avoid skin contact. In case of skin contact, flush the affected area with plenty of water immediately.
Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and seek medical advice immediately.

【NOTICE】
1) This product shall only be used on matched Mindray Auto Hematology Analyzer. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
2) Do not use the reagent if it is frozen.
3) If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
4) If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
5) Dispose of the waste, residual and contaminated packing based on local regulations.
6) The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
7) This product must be operated by clinical professionals, physician or laboratorian trained by Mindray or Mindray authorized personnel.

【INTELLECTUAL PROPERTY STATEMENT】
SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.
Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.
Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

【COMPANY CONTACT】	
Manufacturer:	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
E-mail Address:	service@mindray.com
Tel:	+86 755 81888998
Fax:	+86 755 26582680
EC-Representative:	Shanghai International Holding Corp. GmbH(Europe)
Address:	Eiffestraße 80, 20537 Hamburg, Germany
Tel:	0049-40-2513175
Fax:	0049-40-255726

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