**Total Protein in Urine/CSF (TPUC) Kit
(Pyrogallol Red-Molybdate Method)****Order Information**

Cat. No.	Package size
105-009168-00	R: 3×18 mL + Cal: 1×1 mL
105-009169-00	R: 3×18 mL + Cal: 1×1 mL
105-009170-00	R: 2×40 mL + Cal: 1×1 mL

Intended Purpose

In vitro test for the quantitative determination of Total Protein concentration in human urine or cerebrospinal fluid (CSF) on Mindray BS series chemistry analyzers. Protein measurements in urine are intended to be used for diagnosing and monitoring the therapeutic effect of proteinuria. And it could also be used for aiding to diagnose renal diseases. CSF protein measurements are intended to be used for aiding to diagnose infectious diseases of the central nervous system.

Summary¹⁻³

Measurement of total protein in urine is important in the diagnosis and treatment of diseases associated with renal, cardiac and thyroid function, which are characterized by proteinuria or albuminuria. CSF protein measurements are used in the diagnosis and treatment of conditions such as meningitis, brain tumors and infections of the central nervous system. Urine is formed by ultrafiltration of plasma across the glomerular capillary wall. Proteins with a relative molecular mass >40,000 are almost completely retained, while smaller substances easily enter the glomerular filtrate. Most CSF protein originates by diffusion from plasma across the blood-CSF barrier. Elevated levels of CSF protein occur as a result of increased permeability of the blood-CSF barrier or with increased local synthesis of immunoglobulins.

Assay Principle

Pyrogallol Red-Molybdate method.

Pyrogallol red-molybdate complex + Protein $\xrightarrow{H^+}$ Red complex

Pyrogallol red is combined with molybdate to form a red complex. The assay is based on the shift in absorption that occurs when the pyrogallol red-molybdate complex binds basic amino groups of protein molecules and a blue-purple complex is formed. The absorbance of this complex is directly proportional to the protein concentration in the sample.

Reagents Components

R:	Succinic Acid	100 mmol/L
	Pluronic F-127	0.10%
	Pyrogallol red	0.05 mmol/L
	Sodium molybdate	0.36 mmol/L
	Sodium benzoate	0.10%
Calibrator:	BSA	0.10%
	Sodium azide	0.09%

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

On board in use, the reagents are stable for 30 days when refrigerated on the analyzer. Once opened in use, the calibrator is stable for 30 days when refrigerated at 2-8°C and protected from light.

Contamination must be avoided.

Do not freeze the reagent and calibrator.

Specimen collection and preparation

■ Specimen types

Urine or cerebrospinal fluid (CSF) are suitable for samples.

■ Preparation for Analysis

1. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
2. Centrifuge samples containing precipitate before performing the assay.
3. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment.

■ Sample Stability⁴

Urine: 1 day at 15-25°C

7 days at 2-8°C

1 month at (-25)-(-15) °C

CSF: 1 day at 15-25°C

6 days at 2-8°C

>1 year at (-25)-(-15) °C

For longer storage periods, samples should be frozen at (-20°C)⁵. Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

Reagent Preparation

R are ready to use.

Calibrator is ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

1. General laboratory materials: NaCl solution 9 g/L (saline), distilled/deionized water.
2. Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.
3. Mindray BS series chemistry analyzers and General laboratory equipment.

Assay procedure

Parameters Item	BS-800 chemistry analyzers
Assay type	Endpoint
Wavelength (Primary/Secondary)	605/700 nm
Reaction direction	Increase
Reagent	200 μ L
Mix, incubate at 37°C for 1-2 min, read the absorbance A1, then add:	
Sample or Calibrator	4 μ L
Mix thoroughly, incubate at 37°C for 3-5 min, read the absorbance A2, Then calculate $\Delta A = (A2 - A1)$	

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

Calibration

1. It is recommended to use the Mindray Calibrator (TPUC Calibrator within the kit or other suitable calibrators) and 9 g/L NaCl (saline) for two-point calibration. The calibrator values are assigned by Mindray standard transfer procedure and routine method. The traceability process is based on ISO 17511⁶, the calibrator is traceable to NIST SRM927d.

2. Calibration frequency

Calibration is stable for approximately 30 days on BS-800 chemistry analyzers. The calibration stability may vary in different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.

Recalibration may be necessary when the following occur:

- As changed reagent lot.
 - As required following quality control procedures or out of control.
 - As executes specific maintenance or troubleshooting procedure of chemistry analyzers.
3. The calibrator values are lot-specific with the matched models listed in the value sheet.

Quality control

1. It is recommended to use the Mindray Control (TPUC Control: 105-009193-00 or other suitable controls) to verify the performance of the measurement procedure; other suitable control material can be used in addition.
2. Two levels of control material are recommended to analyze each batch of samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual.
3. Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control doesn't recover within the acceptable tolerances.

Calculation

$\Delta A = (A_2 - A_1)$ sample or calibrator.

$C \text{ sample} = (\Delta A \text{ sample} / \Delta A \text{ calibrator}) \times C \text{ calibrator}$.

The BS series chemistry analyzer detects the change of absorbance (ΔA) and calculates the TPUC concentration of each sample automatically after calibration.

Dilution

If the value of sample exceeds 2000 mg/L, the sample should be diluted with 9 g/L NaCl solution (saline) (e.g. 1+9) and rerun; the result should be multiplied by 10.

Expected values

Sample Type	Units
Urine ²	24-141 mg/24h
CSF ⁷	150-450 mg/L

The expected value is provided from reference, Mindray has verified the expected value of CSF by 117 CSF samples of people from China.

Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

Performance Characteristics

■ Analytical Sensitivity

The Total Protein in Urine/CSF (TPUC) Kit has an analytical sensitivity of 20 mg/L on BS-800. Analytical sensitivity is defined as the lowest concentration of analyte that can be distinguished from a sample that contains no analyte. It is calculated as the value lying 3 standard deviations above that of the mean from 20 replicates of an analyte-free sample.

■ Measuring range

The Mindray BS series systems provide the following linearity range:

Sample Type	Units
Urine/CSF	20-2000 mg/L

A high concentration sample (approximately 2000 mg/L) is mixed with a low concentration sample (<20 mg/L) at different ratios, generating a series of dilutions. The concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.990$. The reportable range is 20-20000 mg/L.

■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3⁸, each sample was assayed 2 times per run, 2 runs per day, a total of 20 days. The precision data of controls on BS-800 are summarized below*.

Specimen Type (N=80)	Mean (mg/L)	Repeatability		Within-Lab	
		SD (mg/L)	CV %	SD (mg/L)	CV %
Control L	151.89	2.93	1.93	3.23	2.13
Control H	400.15	1.82	0.46	2.09	0.52

*Representative data, results in different instruments, laboratories may vary.

■ Analytical Specificity

The samples with different concentration interfering substance were prepared by addition of interferent to human urine sample pools, and recovers are within $\pm 10\%$ of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The data of interference studies on BS-800 are summarized below.

Interfering Substance	Interferent Concentration (mg/dL)	Analyte Concentration (mg/L)	Relative Deviation (%)*
Ascorbic acid	100	404.60	-1.89
Conjugated bilirubin	20	392.50	+4.67
Glucose	200	406.70	-3.37
Intralipid	10	398.03	+0.72
Nitrite	200	397.97	-3.17
Sodium Citrate	50	402.76	-0.12
Sodium Oxalate	15	402.49	-1.91
Calcium	150	406.16	-7.44
Magnesium	200	397.69	+0.66
Urea	1000	401.67	-3.57
Uric Acid	100	400.73	+0.95
Creatinine	250	400.49	+0.06

*Representative data, results in different instruments, laboratories may vary.

■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP09-A3⁹. The Mindray System (Mindray BS-800/Mindray TPUC Reagent) (y) was compared with comparison system (Mindray BS-800/ Medica System TPUC Reagent) (x) using the same urine and CSF specimens. The statistical data obtained by linear regression are shown in the table below*:

Sample Type	Regression Fit	Correlation Coefficient (r)	Sample (N)	Concentration Range (mg/L)
urine	$y=1.0024x+7.8268$	0.9989	110	20.80-1870.30
CSF	$y=1.0068x+6.7181$	0.9998	110	20.55-1999.45

*Representative data, results in different instruments, laboratories may vary.

Result interpretation

The results could be affected by drugs, disease, or endogenous substances¹⁰. When the reaction curve is abnormal, it is recommended to retest and check the result.

Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The reagents avoid direct exposure to

sunlight and freezing. The results can't be assured when stored at inappropriate condition.

4. If unintentionally opened before used, store the reagents tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix reagents with different lots and bottles.
Do not use the reagents beyond the expiration date and the in-use date.
Do not mix new opened reagents with in-use reagents.
Avoid the formation of foam.
6. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
7. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
8. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
9. When the reagent accidentally enters the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
10. Safety data sheet is available for professional user on request.
11. Disposal of all waste material should be in accordance with local guidelines.
12. All human material should be considered potentially infectious.
13. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
14. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
15. This kit contains components classified as follows in accordance with the Regulation (EC) No 1272/2008:

Warning

H315	Causes skin irritation.
------	-------------------------

H319	Causes serious eye irritation.
------	--------------------------------

Prevention:

P280	Wear protective gloves, protective clothing, eye protection and face protection.
------	--

P264	Wash all exposed external body areas thoroughly after handling.
Response:	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.
P302+P352	IF ON SKIN: Wash with plenty of water.
P332+P313	If skin irritation occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

References

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;556-559, 669-679.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
3. Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019;46:283-286.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
5. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
6. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
7. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.
8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671,3-683.

Graphical symbols

In Vitro Diagnostic
medical device

Unique device
identifier

European
Conformity

Consult Instructions
For use

Use-by
date

Authorized representative in
the European Community

Batch Code

Temperature
limit

Manufacturer

Catalogue
number

Keep away from sunlight

Indicates a medical device that needs protection from light sources

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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- 1.Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;556-559, 669-679.
- 2.Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
- 3.Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019,46:283-286.

- 4.Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
- 5.CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
- 6.ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
- 7.Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.
- 8.CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- 9.CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- 10.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671,3-683.

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Referências

- 1.Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders

2012;556-559, 669-679.

- 2.Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
- 3.Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019,46:283-286.
- 4.Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
- 5.CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
- 6.ISO 17511:2003. In vitro diagnostic medical devices-Metrological traceability of values assigned to calibrators and control materials.
- 7.Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.
- 8.CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- 9.CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- 10.Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671,3-683.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Morada: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, R. P. China

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Referencias

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012; 556-559, 669-679.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
3. Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019, 46:283-286.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
5. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
6. ISO 17511:2003. In vitro diagnostic medical devices-Metrological traceability of values assigned to calibrators and control materials.
7. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.
8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671, 3-683.

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Bibliografia

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;556-559, 669-679.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Francoforte: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
3. Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019;46:283-286.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
5. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
6. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
7. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.
8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671,3-683.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, Repubblica popolare cinese

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Referanslar

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;556-559, 669-679.
2. Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
3. Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019;46:283-286.
4. Ehret W, Heil W, Schmitt Y, Töpfer G, Wissner H, Zawta B, et al. Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
5. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
6. ISO 17511:2003. In vitro diagnostic medical devices- Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and Kontrol materials.
7. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.

8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671,3-683.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. Çin

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Références

- 1.Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;556-559, 669-679.
- 2.Thomas L. Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft 1998; 377-400, 1308-1326.
- 3.Bourque PR, McCudden CR, Warman-Chardon J, Brooks J, Hegen H, et al. A Survey of Cerebrospinal Fluid Total Protein Upper Limits in Canada: Time for an Update? Can J Neurol Sci, 2019;46:283-286.
- 4.Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of

- Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/DIL/LAB/99.1 Rev.2:48-50pp.
5. CLSI. Urinalysis; Approved Guideline-Third Edition. CLSI document GP16-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
 6. ISO 17511:2003. In vitro diagnostic medical devices-Measurement of quantities in biological samples-Metrological traceability of values assigned to calibrators and control materials.
 7. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 916-918.
 8. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
 9. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
 10. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-671,3-683.

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© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)