TERUMOBCT

Trima Accel® Disposable Tubing Sets
Instructions for Use
Rx Only
Jednorázové sety Trima Accel®
Návod k použití
Trima Accel® engangsslangesæt
Brugsanvisning
Trima Accel®-Einwegschlauchsets
Gebrauchsanweisung
Αναλώσιμα σετ σωλήνων Trima Accel®
Οδηγίες χρήσης
Equipos de líneas desechables de Trima Accel®
Instrucciones de uso
Kertakäyttöiset Trima Accel® -letkusarjat
Käyttöohjeet
Kits de collecte à usage unique Trima Accel®
Mode d'emploi
Kompleti cijevi Trima Accel za jednokratnu upotrebu
Upute za upotrebu
Trima Accel® egyszer használatos csőszerelékek
Használati utasítás
Circuiti monouso Trima Accel®
Istruzioni per l'uso
Trima Accel® disposable slangensets
Gebruiksaanwijzing
Trima Accel engangsslangesett
Bruksanvisning
Jednorazowe zestawy drenów Trima Accel®
Instrukcje użytkowania
Conjuntos Descartáveis Trima Accel®
Instruções de Uso
Conjuntos de tubos descartáveis Trima Accel®
Instruções de utilização
Одноразовые комплекты магистралей Trima Accel
Инструкции по применению
Trima Accel® engångsslangset
Bruksanvisning
Trima Accel Tek Kullanımlık Setler
Kullanım Yönergeleri

82300

Trima Accel®LRS® Platelet, Plasma Set
Set Trima Accel®LRS® pro trombocyty a plazmu
Trima Accel®LRS® Trombocyt, Plasmasæt
Trima Accel®LRS®-Set für Thrombozyten + Probennehmer-, Plasma- und Erythrozyten-Verfahren
Σετ Trima Accel®LRS® αιμοπεταλίων, Πλάσμα
Equipo LRS® de Trima Accel® para plaquetas, plasma
Trima Accel®LRS® Platelet-, Plasma-sarja
Kit de collecte Plaquettes, Plasma LRS® Trima Accel®
Komplet za trombocite, plazmu i eritrocite Trima Accel®LRS®
Trima Accel®LRS® trombocita, plazmaszerelék
Set per piastrine, plasma LRS® Trima Accel®
Trima Accel®LRS® bloedplaatjes, plasmaset
Trima Accel®-sett for LRS® blodplate, plasma
System Trima Accel®LRS® z zestawami do pobierania płytek krwi i osocza
Conjunto de Plaquetas, Plasma Trima Accel®LRS®
Sistema Trima Accel®LRS® para plaquetas, Plasma
Комплект Trima Accel®LRS® для заготовки тромбоцитов; плазмы
Set med Trima Accel®Leukocytreducingssystem (LRS®) Trombocyter, Plasma
Trima Accel®LRS® Trombosit, Plazma Seti

82400, 82406

Trima Accel® LRS® Platelet, Plasma, RBC Set
Set Trima Accel® LRS® pro trombocyty, plazmu a erytrocyty
Trima Accel® LRS® Trombocyt, Plasma, RBC-sæt
Trima Accel®-Set für Erythrozyten- und Plasma-Verfahren
Σετ Trima Accel® LRS® αιμοπεταλίων, Πλάσμα, RBC
Equipo LRS® de Trima Accel® para plaquetas, plasma, hematías
Trima Accel® LRS® Platelet-, Plasma-, RBC-sarja
Kit de collecte Plaquettes, Plasma, Hématies LRS® Trima Accel®
Komplet za trombocite, plazmu i eritrocite Trima Accel® LRS®
Trima Accel® LRS® trombocita, plazma, VVT-szerelék
Set per piastrine, plasma, GR LRS® Trima Accel®
Trima Accel® LRS® bloedplaatjes, plasma, RBC-set
Trima Accel® LRS®-sett for blodplater, plasma, RBC
System Trima Accel® LRS® z zestawami do pobierania płytek krwi, osocza i KKCz
Conjunto de Plaquetas, Plasma, Hemácias Trima Accel® LRS®
Sistema Trima Accel® LRS® para plaquetas, Plasma, GV
Комплект Trima Accel® LRS® для заготовки тромбоцитов; плазмы; эритроцитов
Set med Trima Accel® Leukocytreducingssystem (LRS®) Trombocyter, Plasma, RBC
Trima Accel® LRS® Trombosit, Plazma, Eritrosit Seti

82408, 82416

Trima Accel®LRS® Platelet + 35 mL Sampler, MultiPlasma, RBC Set
Set Trima Accel®LRS® – trombocyty + odběr vzorků 35 ml, multiplazma, erytrocyty
Trima Accel®LRS® Trombocyt + 35 ml Prøveudtag, MultiPlasma, RBC-sæt
Trima Accel®-Set für MultiPlasma-Verfahren
Σετ Trima Accel®LRS® αιμοπεταλίων + δειγματολήπτης 35 ml, MultiPlasma, RBC
Equipo LRS® de Trima Accel® para plaquetas + muestreador de 35 ml, multiplasma, hematías
Trima Accel®LRS® Platelet + 35 mL Sampler-, MultiPlasma-, RBC-sarja
Kit de collecte Plaquettes + Échantillonneur 35 ml, Multiplasma, Hématies
LRS® Trima Accel®

Trima Accel®LRS® komplet za trombocite + uređaj za uzimanje uzoraka od 35mL, MultiPlasma, RBC

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Intended Use

Trima Accel® disposable tubing sets are intended for use with the Trima Accel® Automated Blood Collection System. Depending on the tubing set used, platelets, plasma and/or red blood cells may be collected.

- Adequate studies have not been performed to evaluate the quality of gamma irradiation or freezing of ACD-A/AS-3 red blood cell (RBC) products collected with gravity drain leukoreduction process (TLR filter) on the Trima Accel system.
- Studies have not been performed to support gamma irradiation or freezing of ACD-A/ AS-3 RBC products collected with an integrated in-line RBC leukoreduction filter(s) (Auto RBC filter) on the Trima Accel system.

Please refer to the Trima Accel Automated Blood Collection system operator's manual for a complete listing of warnings, cautions, adverse effects and operating instructions.

Warnings

1. If the integrity of the tubing set is compromised for any reason, the tubing set is no longer functionally closed, and the product should not be stored beyond 24 hours.
2. Check that all pump headers have been properly loaded before beginning a procedure. Visually inspect each pump and carefully run your finger around the edge of each pump.
3. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
4. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Exposure to excessive ethylene oxide (EO) residuals
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

5. The blood and fluid pathways of the tubing set are sterilized with ethylene oxide and are nonpyrogenic. The packaging cover of the tubing set serves to protect the set and keep contents in place during shipping. Under any of the following conditions, DO NOT use the set:
 - The end caps are not in place
 - The needle cap plug is not in place
 - Severe kinks in the tubing are apparent
 - The tubing set is incorrectly assembled
 - The tubing set is damaged
 - The clamps are closed on the pre-attached filter tubing assembly

Cautions

1. Each operator should be thoroughly familiar with the Trima Accel system operator's manual. All procedures should be performed by qualified medical personnel under the supervision of a physician.
2. Use aseptic technique throughout all procedures to ensure donor safety and product quality.
3. When loading the channel and the LRS chamber into the centrifuge filler, be careful not to stretch the tubing, as leaks could develop.
4. Be careful to press on tubing only with your fingers when loading the channel. To avoid accidental puncture, never use any sharp object.
5. When using the Trima Accel Automated Blood Collection System, Version 5, do not select a procedure configured with an RBC product when using a disposable tubing set that does not contain an RBC bag. If a procedure with an RBC product is selected in error, a "centrifuge pressure high" alert will be displayed when RBC collection begins. Select for rinseback and continue. Any RBC in the system will be returned to the donor.
6. When using MultiPlasma tubing sets, ensure that at least two plasma line clamps are open for procedures with plasma volumes greater than 600 mL.

Procedures for Use

Refer to the Trima Accel Automated Blood Collection System operator's manual for detailed information regarding the following procedures:

- Setting up the disposable tubing set
- Performing a collection procedure
- Removing products and the disposable tubing set

Supplemental Instructions for Connecting Anticoagulant (AC)

When you are prompted to connect AC, connect the luer on the AC line (orange) to the AC bag. Break the frangible connector on the AC bag to allow the AC to flow. Squeeze the drip chamber in order to fill it approximately halfway. Follow the system prompts to load the AC line into the AC sensor, and continue according to the Trima Accel system operator’s manual.

Note: When you break frangible connectors, bend them in both directions to ensure that you break them completely. Failure to do so may result in restricted flow.

Symbols and Certifications

Symbol	Definition
	Indicates the product quantity when the quantity is placed in the square.
	Indicates the manufacturer of the product when accompanied by the name of the manufacturer.
	Indicates the expiration date of the product when accompanied by a specific date.
	Indicates the date of manufacture (or sterilization date, if the product is sterile) when accompanied by a specific date.
	Indicates the product catalog number when accompanied by a number.
	Indicates the product lot number when accompanied by a number.
	Indicates a sterile fluid pathway using ethylene oxide (EO).
	Indicates that the user should consult the instructions for use.
	Indicates that the product should not be reused.
	Indicates the product was manufactured in accordance with Annex II of the European Council Directive 93/42/EEC, as amended.
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates that the product must be kept away from sunlight.

Symbol	Definition
	Indicates that the product must be kept dry.
	Indicates that the product is made from High Density Polyethylene.
	Indicates the European Authorized Representative of the product when accompanied by the name of the European Authorized Representative.
	Indicates that the product should not be re-sterilized.
	Indicates that the product is non-pyrogenic.

Return of Used Product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department.

IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.

Please contact your local representative for information regarding returned goods and product complaints.

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©2019 Terumo BCT, Inc.