Clareon Vivity™ Extended Vision Hydrophobic Intraocular Lenses (IOLs)
Models: CNWET0, CNWET2, CNWET3, CNWET4, CNWET5, CNWET6

DESCRIPTION

The non-diffractive Clareon Vivity™ Extended Vision Hydrophobic Posterior Chamber Intraocular Lenses (IOLs) Model CNWET0 and Toric Models CNWET2, CNWET3, CNWET4, CNWET5, and CNWET6 are UV-absorbing and blue light filtering foldable IOLs which, compared to a monofocal IOL, provide an extended range of vision from distance to near without increasing the incidence of visual disturbances. The Clareon Vivity™ Toric Extended Vision Hydrophobic IOL models also compensate for corneal astigmatism.

The single-piece design consists of a high refractive index hydrophobic acrylic material with proprietary blue light filtering chromophore, which filters light in a manner that approximates the human crystalline lens in the 400-475 nm blue light wavelength range (Boettner and Wolter, 1962). In addition to standard UV-light filtering, the blue-light filtering chromophore reduces transmittance of blue light wavelengths. The biconvex aspheric optic consists of a high refractive index soft acrylic material capable of being folded prior to insertion, allowing placement through an incision smaller than the optic diameter of the lens. After surgical insertion into the eye, the lens gently unfolds to restore its optical performance. The supporting haptics provide proper positioning and fixation of the IOL optic within the eye.

The depth of field extension and monofocal visual disturbance profile are achieved through the patented wavefront-shaping X-WAVE™ technology located on the anterior surface of the IOL. Like a monofocal IOL, the wavefront-shaping X-WAVE™ technology utilizes essentially all of the transmitted light energy and it does so across the range of vision. The location of the wavefront-shaping X-WAVE™ technology is identical for all lens powers. The anterior surface of the Clareon Vivity™ Extended Vision Hydrophobic IOL is designed with negative spherical aberration to compensate for the positive spherical aberration of the cornea. The posterior surface of the Toric lens is biconic creating a toricity to correct the astigmatism on the cornea. The flat meridian of the Clareon Vivity™ Toric IOLs is identified with indentations (dots) on the posterior surface of the optic. The physical characteristics of the Vivity™ are summarized in Figures 1, 2 and 3 and Table 1.

Figure 1: Physical Characteristics of Clareon Vivity™ IOL
(All dimensions in millimeters)


Figure 2: Physical Characteristics of Clareon Vivity™ Toric IOLs
(All dimensions in millimeters)



Table 1: Physical Characteristics of Clareon Vivity™ IOL

Characteristics	Model					
	CNWET0	CNWET2	CNWET3	CNWET4	CNWET5	CNWET6
Optic Type	Biconvex wavefront-shaping X-WAVE™ technology	Biconvex Toric wavefront-shaping X-WAVE™ technology				
Optic / Haptic Material	Ultraviolet-absorbing and blue light filtering hydrophobic Acrylate/ Methacrylate Copolymer UV cutoff at 10% T: 403 nm (+20.0 diopter lens)					
IOL Powers (spherical equivalent diopters)	+6.0 – +30.0 diopters (in 0.5 diopter increments) +31.0 diopters - +34 in 1.0 diopter increments					
IOL Cylinder Powers (Diopters)	N/A	1.00	1.50	2.25	3.00	3.75
IOL Cylinder Powers – Corneal Plane* (Diopters)	N/A	0.65	0.97	1.46	1.95	2.44
Index of Refraction	1.55					
Haptic Configuration	STABLEFORCE™ Modified-L Haptics					
Optic Diameter (mm)	6.0					
Overall Length (mm)	13.0					
Haptic Angle	0°					

* Based on an average pseudophakic human eye using an SRK/T optical A-constant of 119.2

Figure 3: Spectral Transmittance

**NOTES:**

- The cutoff wavelength and the spectral transmittance curves presented here represent the range of transmittance values of IOLs made from hydrophobic acrylate/methacrylate copolymer with bonded UV-absorber and Alcon's proprietary blue light filtering chromophore.
- Measurements were by direct transmittance using Clareon™ IOLs with center thickness equivalent to the marketed range.
- Human crystalline lens data is from Boettner and Wolter (1962).

INDICATIONS

The non-diffractive Clareon Vivity™ Extended Vision Hydrophobic IOL Model CNWET0 employs a patented wavefront-shaping X-WAVE™ technology that provides an extended range of vision, while maintaining a monofocal visual disturbance profile. It is intended for primary implantation for the visual correction of aphakia in adult patients with <1.00 D of preoperative corneal astigmatism, in whom a cataractous lens has been removed by extracapsular cataract extraction. The lens treats the effects of presbyopia by providing a full range of vision from distance to functional near (40 cm). The lens also provides comparable distance visual acuity and decreased spectacle dependence, relative to an aspheric monofocal IOL. The Clareon Vivity™ IOL is intended for capsular bag placement only.

The non-diffractive Clareon Vivity™ Toric Extended Vision Hydrophobic IOL Models CNWET2, CNWET3, CNWET4, CNWET5, and CNWET6 employ a patented wavefront-shaping X-WAVE™ technology that provides an extended range of vision, while maintaining a monofocal visual disturbance profile. They are intended for primary implantation for the visual correction of aphakia and for reduction of residual refractive cylinder in adult patients with pre-existing corneal astigmatism, in whom a cataractous lens has been removed by extracapsular cataract extraction. The lenses treat the effects of presbyopia by providing a full range of vision from distance to functional near (40 cm). These lenses also provide comparable distance visual acuity and decreased spectacle dependence, relative to an aspheric monofocal IOL. The Clareon Vivity™ Toric IOLs are intended for capsular bag placement only.

IOL IMPLANTATION

During standard cataract surgery and implantation of the Clareon Vivity™ IOL, an Alcon qualified delivery system and viscoelastic combination should be used. The use of an unqualified combination may cause damage to the lens and potential complications during the implantation process. Alcon recommends using the qualified MONARCH™ IOL Delivery System or any other Alcon qualified combination. Currently qualified combinations that can be used with these lenses are listed in Table 2; contact Alcon for further information on qualified combinations.

Table 2: Qualified Combinations of Compatible Products

Lens Model	Diopter Range	Cartridge	Handpiece	Ophthalmic Viscosurgical Device (OVD)
CNWET0 CNWET2	+6.0 to +25.0	MONARCH™ III D (8065977763)	CLAREON MONARCH™ IV (SILVER/GREY) (8065977774)	VISCOAT™ OVD PROVISC™ OVD DISCOVISC™ OVD
CNWET3 CNWET4	+6.0 to +30.0	MONARCH™ III C (8065977762)		
CNWET5 CNWET6	+6.0 to +34.0	MONARCH™ II B (8065977758)		

Contact Alcon for any additional Alcon qualified OVDs, handpieces, and cartridges for use with this lens.

See "INSTRUCTIONS FOR USE" for implantation instructions for the Clareon Vivity™ IOLs, which are similar to those of Alcon's Clareon monofocal IOLs.

CALCULATION OF LENS POWER

Accurate biometry is essential for successful visual outcomes. Preoperative calculation of required lens power for the Clareon Vivity™ IOLs should be determined by the surgeon's experience and preference. A reference SRK/T A-Constant value for optical biometry equipment such as ARGOS[®], IOLMaster[®] or LenStar[®] as well as *contact* ultrasound biometry is listed on the outer label. The reference optical A-Constant anticipates the use of both corneal power and axial length values from optical biometry or contact ultrasound equipment with standard settings for a typical patient population and a spectacle far point at 6 meters. In general, lens constants must be "personalized" to compensate for such things as differences in instrumentation, surgical techniques, and IOL power calculation that may exist between clinical practice. IOL power calculation methods are often included with biometry equipment, and they are also described in the references (Hoffer 1993; Holladay 1997; Olsen 2007; Retzlaff, Sanders & Kraff 1990; Haigis 2014).

WARNINGS

1. This lens should not be implanted if the posterior capsule is ruptured, if the zonules are damaged, or if a primary posterior capsulotomy is planned.
2. Rotation of the Clareon Vivity™ Toric IOL away from its intended axis can reduce the astigmatic correction. Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, lens repositioning should occur as early as possible prior to lens encapsulation. Some clinical cases suggest encapsulation is complete within four weeks of implantation.
3. It is recommended that viscoelastic be removed from the eye at the close of surgery with emphasis on the space between the posterior capsule and lens. This may be accomplished by gently depressing the IOL optic posteriorly with the I/A tip and using standard irrigation/aspiration techniques to remove the viscoelastic agent from the eye. This should force any trapped viscoelastic anteriorly where it can be easily aspirated. Residual viscoelastic may allow the lens to rotate causing misalignment of the Clareon Vivity™ Toric IOL with the intended axis of placement.

PRECAUTIONS

1. Prior to surgery, prospective patients should be informed of the possible risks and benefits of the Clareon Vivity™ IOL Models CNWET0, CNWET2, CNWET3, CNWET4, CNWET5, CNWET6.
2. The safety and effectiveness of this IOL have not been studied in patients with certain pre-existing conditions and/or intraoperative conditions (listed in Tables 3 and 4) as these patients were excluded from the clinical studies. Patients with pre-existing conditions may not achieve the visual acuity of patients without such conditions. As with the implantation of any IOL, careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting a lens in a patient with one or more of these conditions.
3. As with any surgical procedure, there is risk involved. Thus, the patient should be followed postoperatively as medically appropriate. Potential complications and/or residual risks accompanying cataract and/or IOL implantation surgery may include the following: blurred or decreased vision, diplopia, glaucoma, unexpected post-operative refraction, ocular irritation, tissue damage (iris prolapse, iris damage, capsular tissue damage, and corneal endothelial cell damage), pigment dispersion, lens epithelial cell on-growth, posterior capsular opacification, infection (endophthalmitis), inflammatory reactions [e.g., toxic anterior segment syndrome (TASS), hypopyon, toxicity response, vitritis, anterior uveitis, and cyclitic membrane], retinal detachment, cystoid macular edema, corneal edema, corneal decompensation, pupillary block, elevated IOP (transient or persistent), hyphema, reduced contrast sensitivity, color misperception, IOL decentration, tilt or dislocation, visual disturbances or spatial distortions, and secondary surgical interventions. Secondary surgical interventions include: lens repositioning, lens replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.
4. Accurate keratometry and biometry in addition to the use of the Toric Calculator for the Clareon Vivity™ Toric IOL models (<http://www.mylcon-toriccalc.com>) are recommended to achieve optimal visual outcomes.
5. Autorefractors may not provide optimal postoperative refraction of patients with Clareon Vivity™ Toric IOLs. Manual refraction with maximum plus technique is strongly recommended.
6. Clinical data supporting the use of these IOLs were collected with the lens implanted in the capsular bag only. There are no clinical data to demonstrate safety or effectiveness for placement in the ciliary sulcus.
7. Pivotal study results indicate mesopic contrast sensitivity was reduced in patients receiving the parent lens, AcrySof IQ Vivity™ IOL, compared to the Monofocal Control IOL in certain high spatial frequency test conditions.
8. Do not re-sterilize these IOLs by any method (see RETURNED GOODS POLICY).
9. Do not store these IOLs at temperatures over 30°C (86°F).
10. Handle lenses carefully to avoid damage to lens surfaces or haptics.
11. Do not attempt to reshape haptics in any way.
12. DO NOT implant the IOL if the sterility has been compromised or if the sterile package has been unintentionally opened before use.
13. When binocular implantation of the Clareon Vivity™ IOL is planned, both eyes of a subject are not intended to be operated on the same day. Simultaneous binocular implantation has not been studied.
14. A high level of surgical skill is required for any intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.

Table 3: Preexisting Conditions with No Safety and Effectiveness Data

• Clinically severe corneal dystrophy (e.g., epithelial, stromal, or endothelial dystrophy), corneal irregularity (including irregularity due to dry eye syndrome), keratitis, keratoconjunctivitis, keratouveitis, keratopathy, or kerectasia • Irregular corneal astigmatism • Subjects with diagnosed degenerative visual disorders (e.g., macular degeneration or other retinal disorders) that are predicted (by subjective assessment of the retina) to cause future acuity losses to a level worse than 0.3 logMAR • Any inflammation or edema of the cornea • Significant irregular corneal aberration • Previous corneal transplant • Aniridia • Iris neovascularization • Uncontrolled glaucoma • Rubella, congenital, traumatic, or complicated cataracts • Extremely shallow anterior chamber, not due to swollen cataract • Choroidal hemorrhage • Previous retinal detachment • Diabetic retinopathy • Optic nerve atrophy • Recurrent anterior or posterior segment inflammation.	• Amblyopia • Pre-existing ocular conditions which may negatively impact stability of the implant (e.g., diagnosis of pseudoexfoliation syndrome) • Color vision deficiencies - o Studies have shown that color vision discrimination is not adversely affected in individuals implanted with a blue light filtering IOL who have normal color vision. The effect of the blue light filtering IOL in subjects with hereditary color vision defects and acquired color vision defects secondary to ocular disease (e.g., glaucoma, diabetic retinopathy, chronic uveitis, and other retinal or optic nerve diseases) has not been studied • Microphthalmos • Previous refractive surgery • Current or previous usage of alpha-1-selective adrenoceptor blocking agent or antagonist of alpha 1A adrenoceptor [e.g., Flomax[®] (tamsulosin HCL), Hytrin[®], or Cardura[®]] • Pregnancy
---	--

Table 4: Intraoperative Conditions with No Safety and Effectiveness Data

• Other planned ocular surgery procedures, including but not limited to, LASIK, astigmatic keratotomy, and limbal relaxing incisions • Excessive iris mobility or floppy iris syndrome • Mechanical or surgical manipulation required to enlarge the pupil	• Vitreous loss (significant) • Anterior chamber bleeding (significant) • Complications in which the IOL stability could be compromised, including but not limited to zonular dehiscence or capsulorhexis rupture • Uncontrollable IOP
--	---

SELECTION OF TORIC IOL

For selection of the appropriate toric IOL, Alcon provides a web-based tool (www.mylcon-toriccalc.com, Abulafia, Barrett, et al. 2015 and Abulafia, Hill, et al. 2015) that uses pre-operative biometry data, incision location, and surgically induced corneal astigmatism. The astigmatism to be corrected should be determined from biometry data rather than refractive data since the presence of lenticular astigmatism in the crystalline lens to be removed may influence results.

AXIS PLACEMENT OF TORIC IOL

For optimal results, the surgeon must ensure the correct placement and orientation of the lens within the capsular bag. The posterior surface of the IOL is marked with indentations (three at each end) at the haptic/optic junction that identify the flat meridian of the toric IOL optic. These indentations form an imaginary line representing the plus cylinder axis (note: IOL cylinder steep meridian is 90° away). The cylinder axis marks should be aligned with the post-incision steep corneal meridian (intended axis of placement) or as determined by the web-based IOL calculator.

Prior to surgery, the lens placement axis should be marked. One way used in clinical practice to mark the eye is as follows: With the patient sitting upright to prevent cyclotorsion, clearly and precisely mark two reference positions with a surgical skin marker or a marking pencil indicated for ophthalmic use. Using these marks as reference points, an axis marker can be used immediately prior to or during surgery to mark the optimal axis of lens placement identified by the web-based IOL calculator (www.mylcon-toriccalc.com).

Precisely align the axis marking indentations on the toric IOL with the marked intended axis of lens placement. Remove all viscoelastic from both the anterior and posterior sides of the lens since residual viscoelastic may allow the lens to rotate causing misalignment of the toric IOL with the intended axis of placement. This may be accomplished by manipulating the IOL optic with the I/A tip and using standard irrigation/aspiration techniques to remove all viscoelastic from the eye. Bimanual techniques may be used, if preferred, to ensure removal of viscoelastic from behind the lens implant. Special care should be taken to ensure proper positioning of the toric IOL at the intended axis following viscoelastic removal.

Misalignment of the axis of the lens versus the intended axis of placement may compromise its astigmatic correction. Such misalignment can result from inaccurate keratometry or marking of the cornea, inaccurate placement of the toric IOL axis during surgery, an unanticipated surgically induced change in the cornea, or physical rotation of the toric IOL after implantation. In order to minimize this effect, the surgeon should be careful to ensure that preoperative keratometry and biometry is accurate and that the IOL is properly oriented prior to the end of surgery.

INSTRUCTIONS FOR USE

1. Examine the label on the unopened package for model, optical power, cylinder power (if applicable), proper configuration, and expiration date.
2. After the outer box is open, verify that the lens primary label information (e.g., model, powers, and serial number) is consistent with information on the outer box labeling.
3. Inspect the primary package carefully for tears, cuts, punctures, or other signs that the pouch has been opened or damaged. This device is sterile until the inner primary package is opened. DO NOT implant the IOL if the sterility has been compromised or if the sterile package has been unintentionally opened before use (see RETURNED GOODS POLICY).
4. Open the undamaged primary package and transfer the case to a sterile environment. Carefully open the case to expose the lens.
5. Use instrumentation that is scrupulously clean to minimize the occurrence of marks. Any forceps used for lens handling must have round edges and smooth surfaces.
6. Handle the IOL by the haptics only when removing the lens from the case. DO NOT grasp the optical area with forceps. Handle lenses carefully to avoid damage to lens surfaces or haptics. DO NOT attempt to reshape haptics in any way.
7. Examine the IOL carefully prior to insertion to ensure that particles have not adhered during handling.

8. Implant the IOL with the most appropriate surgical procedure for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.

NOTE: During lens loading and insertion, DO NOT allow the Clareon Vivity™ IOL to remain in a folded condition within the selected IOL delivery system for more than 3 minutes prior to completing insertion into the capsular bag.

MAGNETIC RESONANCE COMPATIBILITY

The Clareon Vivity™ IOLs are magnetic resonance (MR) Safe. The IOLs consist of acrylate/methacrylate copolymer material, which is a non-conducting, non-metallic, non-magnetic material that poses no known hazards in all magnetic resonance imaging environments.

PATIENT REGISTRATION

The Patient Identification Card included in the package is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye care practitioner that the patient consults in the future.

Each patient should be registered with Alcon Laboratories, Inc., immediately following implantation of one of these lenses. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens box and mailing it to Alcon Laboratories, Inc. Patient registration is essential for the long-term patient follow-up program and will assist Alcon Laboratories, Inc. in responding to reports of adverse events.

SERIOUS INCIDENT REPORTING

Any serious incident that may reasonably be regarded as device related should be reported to Alcon Canada, Inc.:

Website:
<https://notifeye.alcon.com/external/s/>
By Phone:
In Canada - (800) 613-2245

Each IOL is identified by a serial number which provides traceability, and this information should be given to Alcon.

HOW SUPPLIED

The Clareon Vivity™ IOL is supplied dry, in a package terminally sterilized with ethylene oxide, and must be opened only under aseptic conditions (see INSTRUCTIONS FOR USE).

EXPIRATION DATE

Sterility is guaranteed unless the pouch is damaged or opened. The expiration date is clearly indicated on the outside of the lens package. Any lens held after the expiration date should be returned to Alcon Canada, Inc. (see RETURNED GOODS POLICY).

RETURNED GOODS POLICY

In Canada, returned lenses will only be accepted for credit. All returns must be accompanied by an Alcon Canada, Inc. Returned Goods Number and should be shipped via traceable means. A Returned Goods Number is obtained by contacting Alcon Canada, Inc. Customer Service Department. Issuance of this number does not constitute final acceptance of the returned products. For detailed policy guidelines including exchange, please contact your Sales or Customer Service Representative.

CLINICAL STUDIES

Several prospective clinical studies were performed with parent IOLs that support the safety and effectiveness of the Clareon Vivity™ IOLs (Models CNWET0, CNWET2, CNWET3, CNWET4, CNWET5, CNWET6) and are described in this labeling:

1. Data from the US study of AcrySof IQ Vivity™ IOL Model DFT015 are presented to describe the safety and effectiveness of the Clareon Vivity™ IOL optical design.
2. AcrySof™ Toric clinical study data are presented to describe the safety and effectiveness of the toric optical design of Clareon Vivity™ Toric IOL Models CNWET2, CNWET3, CNWET4, CNWET5 and CNWET6.
3. Data from the US study of Clareon™ IOL demonstrate favorable visual acuity and adverse event outcomes of the Clareon™ IOL. IOL rotational stability of the Clareon™ IOL was also assessed and are shown below to support Clareon Vivity Toric IOL along with safety data to support the Clareon material.

Additionally, a study of the AcrySof IQ Vivity™ IOL was conducted outside the United States (OUS). Data collected on the AcrySof IQ Vivity™ IOL in the OUS study was consistent with that observed in the US AcrySof IQ Vivity™ IOL study, and that of AcrySof™ IOLs generally.

Summaries of these clinical studies are provided below to describe the performance characteristics of the Clareon Vivity™ IOLs. Please use caution when comparing these results with results from similar device studies due to potential differences in patient cohorts, test methods, etc.

NOTE: The Vivity™ optical design on the anterior surface of the Clareon™ IOLs are identical to that of the clinically studied AcrySof IQ Vivity™ IOL. The toric element on the posterior side of the Clareon Vivity™ Toric IOLs is identical to that of the clinically studied AcrySof™ Toric IOLs. Similarly, the Clareon™ IOL material used in the manufacture of the Clareon Vivity™ IOLs is identical to that of the clinically studied Clareon Aspheric Hydrophobic IOL.

1. US AcrySof IQ Vivity™ Model DFT015 Clinical Study

Clinical study ILI875-C002 was a prospective, randomized, controlled study of bilateral implantation of the AcrySof IQ Vivity™ IOL Model DFT015 or the control IOL Model SN60WF. 107 and 113 subjects were implanted with the AcrySof IQ Vivity™ IOL Model DFT015 ("Vivity™ IOL") and the approved monofocal control IOL Model SN60WF, respectively, at 11 investigational sites in the United States. IOL implantation instructions did not change based on randomization assignment.

Primary study effectiveness objectives included demonstrating that the Vivity™ IOL provides non-inferior monocular best-corrected distance visual acuity (BCDVA), superior monocular distance corrected intermediate visual acuity (DCIVA), and a 0.5 D increase in the negative defocus range at 0.2 logMAR compared to the monofocal control IOL. Secondary study effectiveness objectives included demonstrating superior monocular distance corrected near visual acuity (DCNVA) and superiority for percentage of subjects never needing spectacles overall compared to the monofocal control IOL.

Study safety objectives included a comparison of Vivity IOL adverse event rates to ISO Safety and Performance Endpoint (SPE) targets, description of monocular mesopic contrast sensitivity (with and without glare) and estimation of rates of severe and most bothersome visual disturbances.

The clinical study results achieved at 6 months postoperatively demonstrate that the AcrySof IQ Vivity™ IOL is safe and effective for the visual correction of aphakia, provides continuous vision

with superior distance-corrected intermediate and near vision, an increased depth of focus, and decreased spectacle wear when compared to the monofocal control IOL, while demonstrating distance vision comparable to the monofocal control lens, and rates of visual disturbances (e.g. halos, starbursts, glare, etc.) and adverse events similar to the monofocal.

All eyes with successful IOL implantation were considered evaluable for the All Implanted Set (AAS).

Clinical Study Results

Subject Population

A total of 220 subjects were implanted in this study, with 107 subjects receiving the AcrySof IQ Vivity™ IOL and 113 receiving the aspheric Monofocal Control IOL. At the baseline visit, the study consisted of 55.9% females and 44.1% males. Stratifying by race, 97.7% designated White, 0.9% designated Black or African American, 0.5% designated American Indian or Alaska Native, and 0.9% designated "Other". 1.8% of the study population designated their ethnicity as Hispanic or Latino. The mean ($\pm$ SD) age for the study population was 68.8 $\pm$ 7.22 years. Data are reported for the 6 Months visit. One of the subjects in the AcrySof IQ Vivity™ IOL group was not bilaterally implanted and does not contribute to binocular assessments. Two subjects in the monofocal control group discontinued study prior to 6 months.

IOL Power Calculation and Postoperative Manifest Refraction

In this study, surgeons were instructed to select the lens power that targeted emmetropia (closest to 0.0 D) for both the AcrySof IQ Vivity™ and the Monofocal Control IOLs prior to randomization. Table 5 summarizes absolute manifest refraction spherical equivalent (MRSE) data for AcrySof IQ Vivity™ IOL and Monofocal Control IOL subjects at 6 Months. For first eyes, 91.6% implanted with AcrySof IQ Vivity™ IOL and 86.5% implanted with the Monofocal Control IOL achieved MRSE within 0.5 D of emmetropia and mean $\pm$ SD MRSE in first eyes at 6 Months was 0.049 $\pm$ 0.345 D and 0.081 $\pm$ 0.411 D respectively.

Table 5: Mean Absolute MRSE at 6 Months by Treatment Group

MRSE Category	AcrySof IQ Vivity™ IOL		Monofocal Control IOL	
	First Eyes N = 107 n (%)	Second Eyes N = 106 n (%)	First Eyes N = 113 n (%)	Second Eyes N = 113 n (%)
$\leq$ 0.25 D	69 (64.5)	72 (67.9)	63 (56.8)	69 (62.2)
$\leq$ 0.5 D	98 (91.6)	92 (86.8)	96 (86.5)	93 (83.8)
$\leq$ 1.00 D	107 (100.0)	106 (100.0)	108 (97.3)	110 (99.1)
$>$ 1.00 D	0 (0.0)	0 (0.0)	3 (2.7)	1 (0.9)
Total	107	106	111	111

Percentages are calculated as (n/Total)*100

EFFECTIVENESS

Visual Acuity

Visual Acuity was assessed using a computerized test system (CTS, M&S Technologies, Niles, IL). The following subsections present monocular and binocular visual acuity data collected at distance, intermediate, and near. Primary and secondary study analyses included the following assessments of monocular visual acuity:

- Mean best corrected distance visual acuity (BCDVA) and distance corrected intermediate and near visual acuities (DCIVA and DCNVA) by treatment group
- Defocus curves by treatment group

These analyses are discussed in the following subsections.

Monocular Visual Acuity

Table 6 presents monocular, photopic, best corrected distance (4 m) visual acuity (BCDVA), distance corrected intermediate (66 cm) visual acuity (DCIVA), and distance corrected near (40 cm) visual acuity (DCNVA) data for first eyes treated with the AcrySof IQ Vivity™ and monofocal control IOLs at 6 Months.

Table 6: Comparison of Mean Monocular Photopic (First Eyes) BCDVA, DCIVA, and DCNVA Using Least Square Estimates All Implanted, 6 Months Postoperative

Monocular Mean Visual Acuity	AcrySof IQ Vivity™ IOL			Monofocal Control IOL			95% One-Sided UCL	P-Value
	Total	Mean logMAR	Standard Error	Total	Mean logMAR	Standard Error		
BCDVA	107	0.016	0.0091	111	-0.036	0.0089	0.073 ^a	N/A
DCIVA	107	0.148	0.0120	111	0.312	0.0118	N/A	$<$ 0.001 ^b
DCNVA	107	0.359	0.0147	111	0.515	0.0144	N/A	$<$ 0.001 ^b

^aNon-inferiority of AcrySof IQ Vivity™ IOL to the Monofocal Control IOL is demonstrated as the one-sided 95% Upper Confidence Limit (UCL) for treatment difference $<$ 0.10 logMAR

^bSuperiority of AcrySof IQ Vivity™ IOL over the Monofocal Control IOL is demonstrated as the upper bound of two-sided 95% CI for treatment difference is $<$ 0.0 logMAR and p-value is less than alpha = 0.05.

The AcrySof IQ Vivity™ IOL met the clinical performance co-primary endpoint for mean photopic BCDVA non-inferior to the monofocal control IOL (non-inferiority margin of 0.10 logMAR). The upper limit of the one-sided 95% confidence interval (CI) of the mean difference in BCDVA between IOL groups was less than 0.1 logMAR (less than a line), demonstrating that the AcrySof IQ Vivity™ IOL is non-inferior to the control lens in providing BCDVA. Additionally, AcrySof IQ Vivity™ IOL first eyes achieved BCDVA of 0.3 logMAR or better at a rate of 99.1% and 100.0% in the AAS and Best-Case Analysis Set (BAS) set, respectively. This exceeded the ISO BCDVA Safety and Performance Endpoint (SPE) rates of 92.5% and 96.7%, respectively.


The AcrySof IQ Vivity™ IOL provided mean photopic monocular DCIVA superior to the monofocal control IOL with an increase in 1.6 lines read and met the clinical performance co-primary endpoint with $>$ 50% (72.9%) of first eyes achieving DCIVA of 0.2 logMAR or better at 6 Months. When analyzed by pupil size ($<$ 3 mm, $\geq$ 3 - $\leq$ 4 mm, and $>$ 4 mm), mean monocular DCIVA in all pupil size subgroups of the Vivity IOL arm remains greater than 0.2 logMAR, and $>$ 1 line above those of the monofocal control IOL arm.

AcrySof IQ Vivity™ IOL provided mean photopic monocular DCNVA superior to the monofocal control IOL, with 40.2% (95% CI: 30.8%, 50.1%) of AcrySof IQ Vivity™ IOL first eyes achieving DCNVA of 0.3 logMAR or better at 6 Months, at a rate > 25% above that of the monofocal control IOL, and an increase of 1.6 lines read.

Monocular Defocus Curve

Monocular defocus curves were obtained at 6 months for the AcrySof IQ Vivity™ IOL and the monofocal control IOL and are depicted in Figure 4. Data were obtained from 107 and 111 subjects respectively in each arm using a computerized visual acuity test system (CTS, M&S Technologies, Niles, IL). The AcrySof IQ Vivity™ IOL met the clinical performance co-primary endpoint, demonstrating a defocus range of > 0.5 D greater than the monofocal control IOL at 0.2 logMAR at 6 Months.

Figure 4: Mean Monocular (First Eyes) Defocus Curve with 95% Confidence Limits by Lens Model at 6 Months Postoperative All Implanted

Binocular Visual Acuity

Additional descriptive and categorical analyses of binocular uncorrected and best corrected distance visual acuities (UCDVA and BCDVA), uncorrected and distance corrected intermediate visual acuities (UCIVA and DCIVA), and uncorrected and distance corrected near visual acuities (UCNVA and DCNVA), in addition to binocular defocus curves are presented in the following subsections.

Binocular Mean Visual Acuity

Binocular mean BCDVA, UCDVA, DCIVA, UCIVA, DCNVA, and UCNVA are presented by treatment group in Table 7. Consistent with the mean monocular visual acuity data presented in Table 6, binocular visual acuity data for the AcrySof IQ Vivity™ IOL group showed comparable performance at distance, with improvements in intermediate and near visual acuity compared to the monofocal control group. One of the subjects in the AcrySof IQ Vivity™ IOL group was not bilaterally implanted and does not contribute to binocular assessments.

Table 7: Comparison of Binocular Mean Photopic BCDVA, UCDVA, DCIVA, UCIVA, DCNVA, and UCNVA All Implanted, 6 Months Postoperative

Binocular Mean Visual Acuity	AcrySof IQ Vivity™ IOL N = 106			Monofocal Control IOL N = 113		
	n	Mean logMAR	Standard Deviation	n	Mean logMAR	Standard Deviation
BCDVA	106	-0.028	0.084	111	-0.071	0.086
UCDVA	106	0.035	0.102	111	-0.022	0.107
DCIVA	106	0.054	0.093	111	0.196	0.113
UCIVA	106	0.058	0.083	111	0.139	0.122
DCNVA	106	0.253	0.118	111	0.391	0.135
UCNVA	106	0.208	0.104	111	0.339	0.149

Binocular Categorical Visual Acuity

Table 8 presents the percentage of subjects achieving binocular photopic BCDVA and UCDVA by Snellen category at 6 Months. The majority of subjects in both the AcrySof IQ Vivity™ IOL and the monofocal control IOL groups achieved 20/20 or better.

Table 8: Binocular Photopic Distance Best Corrected and Uncorrected Visual Acuity (Snellen) All Implanted, 6 Months Postoperative

Binocular Snellen Visual Acuity	AcrySof IQ Vivity™ IOL N = 106		Monofocal Control IOL N = 113	
	Uncorrected n (%)	Best Corrected n (%)	Uncorrected n (%)	Best Corrected n (%)
20/20 or better	65 (61.3)	93 (87.7)	85 (76.6)	105 (94.6)
20/25 or better	94 (88.7)	103 (97.2)	105 (94.6)	109 (98.2)
20/32 or better	102 (96.2)	104 (98.1)	109 (98.2)	111 (100.0)
20/40 or better	105 (99.1)	106 (100.0)	111 (100.0)	111 (100.0)
Worse than 20/40	1 (0.9)	0 (0.0)	0 (0.0)	0 (0.0)
Total Responses	106	106	111	111

Table 9 presents the percentage of subjects achieving binocular photopic DCIVA and UCIVA by Snellen category at 6 Months. The majority of AcrySof IQ Vivity™ IOL subjects achieved 20/20 or better binocular uncorrected and distance corrected intermediate VA, as compared to 27.0% and 10.8% of control subjects, respectively.

Table 9: Binocular Photopic Intermediate (66 cm) Distance Corrected and Uncorrected Visual Acuity (Snellen)

Binocular Snellen Visual Acuity	AcrySof IQ Vivity™ IOL N = 106		Monofocal Control IOL N = 113	
	Uncorrected n (%)	Distance Corrected n (%)	Uncorrected n (%)	Distance Corrected n (%)
20/20 or better	59 (55.7)	59 (55.7)	30 (27.0)	12 (10.8)
20/25 or better	91 (85.8)	93 (87.7)	70 (63.1)	38 (34.2)
20/32 or better	105 (99.1)	103 (97.2)	92 (82.9)	82 (73.9)
20/40 or better	106 (100.0)	106 (100.0)	106 (95.5)	103 (92.8)
Worse than 20/40	0 (0.0)	0 (0.0)	5 (4.5)	8 (7.2)
Total Responses	106	106	111	111

Table 10 presents the percentage of subjects achieving binocular photopic UCNVA and DCNVA by Snellen category at 6 Months. 67.0% and 57.5% of the AcrySof IQ Vivity™ IOL subjects achieved 20/32 or better binocular uncorrected and distance corrected near VA, compared to 28.8% and 16.2% of control subjects, respectively.


Table 10: Binocular Photopic Near (40 cm) Distance Corrected and Uncorrected Visual Acuity (Snellen) All Implanted, 6 Months Postoperative

Binocular Snellen Visual Acuity	AcrySof IQ Vivity™ IOL N = 106		Monofocal Control IOL N = 113	
	Uncorrected n (%)	Distance Corrected n (%)	Uncorrected n (%)	Distance Corrected n (%)
20/20 or better	3 (2.8)	1 (0.9)	1 (0.9)	0 (0.0)
20/25 or better	37 (34.9)	23 (21.7)	11 (9.9)	2 (1.8)
20/32 or better	71 (67.0)	61 (57.5)	32 (28.8)	18 (16.2)
20/40 or better	96 (90.6)	82 (77.4)	62 (55.9)	48 (43.2)
Worse than 20/40	10 (9.4)	24 (22.6)	49 (44.1)	63 (56.8)
Total Responses	106	106	111	111

Defocus Curves

Binocular defocus curves were obtained at 6 months for the AcrySof IQ Vivity™ IOL and the monofocal control IOL and are depicted in Figure 5. Data were obtained from 106 and 111 subjects in each arm using a computerized visual acuity test system (CTS, M&S Technologies, Niles, IL). Consistent with the monocular defocus curve presented in Figure 4, AcrySof IQ Vivity™ IOL subjects demonstrated a > 0.5 D greater defocus range compared to the monofocal control group.

Figure 5: Mean Monocular (First Eyes) Defocus Curve with 95% Confidence Limits by Lens Model at 6 Months Postoperative All Implanted

Need for Eyeglasses/Contact Lenses

A questionnaire developed by Alcon was used in the study to assess need for eyeglasses/contact lenses following implantation with the IOLs. The AcrySof IQ Vivity™ IOL showed superior performance as 21.6% of subjects (N = 102) versus 3.6% of Monofocal Control IOL subjects (N = 111) reported never needing to wear eyeglasses or contact lenses at 6 Months (difference: 18.0%, 95% CI: 9.65-27.37%).

The questionnaire evaluated subject-reported use of eyeglasses or contact lenses at distance ('far away'), intermediate ('arm's length'), and near ('up close'), in bright and dim lighting conditions. Table 11 summarizes another analysis assessing the proportion of AcrySof IQ Vivity™ IOL and Monofocal Control IOL subjects 'never' or 'rarely' needing eyeglasses or contact lenses. Data for this analysis are presented overall and under each lighting and distance condition.

Table 11: Proportion of Subjects Rarely or Never Needing Eyeglasses or Contact Lenses by Distance and Lighting Condition

Condition		Proportion of Subjects Rarely or Never Needing Eyeglasses or Contact Lenses					
		AcrySof IQ Vivity™ IOL N = 106			Monofocal Control IOL N = 113		
		Total Responses	n	%	Total Responses	n	%
Overall		102	46	45.1	111	19	17.1
Bright Light	Distance ('Far away')	102	96	94.1	111	102	91.9
	Intermediate ('At arm's length')	102	89	87.2	111	64	57.6
	Near ('Up close')	102	47	46.1	111	18	16.2
Dim Light	Distance ('Far away')	102	95	93.2	111	99	89.2
	Intermediate ('At arm's length')	102	86	84.3	111	59	53.1
	Near ('Up close')	102	40	39.2	111	12	10.8

Subjects who reported not using eyeglasses at least some of the time (95.0% and 85.7% in the AcrySof IQ Vivity™ IOL and Monofocal Control IOL groups, respectively) were asked to report their quality of vision without eyeglasses, up close, at arm's length, and far away, in dim and bright light. The proportion of subjects reporting 'good' or 'very good' vision without eyeglasses at 6 Months is presented by treatment group in Table 12.

Table 12: Proportion of Subjects Reporting 'Good' or 'Very Good' Vision Without Eyeglasses by Distance and Lighting Condition

Condition		AcrySof IQ Vivity™ IOL N = 106			Monofocal Control IOL N = 113		
		Total Responses	n	%	Total Responses	n	%
Bright Light	Distance ('Far away')	96	90	93.7	89	82	92.2
	Intermediate ('At arm's length')	96	88	91.6	89	56	62.9
	Near ('Up close')	96	55	57.3	89	22	24.8
Dim Light	Distance ('Far away')	96	84	87.5	90	70	77.8
	Intermediate ('At arm's length')	96	80	83.4	89	45	50.6
	Near ('Up close')	96	36	37.5	89	7	7.8

The results of this questionnaire demonstrate a reduced overall spectacle need in subjects implanted with the AcrySof IQ Vivity™ IOL compared to the monofocal control IOL.

SAFETY

The incidences of cumulative adverse events (AEs) including 1-sided 95% lower confidence limits for the AcrySof IQ Vivity™ IOL as compared to the ISO 11979-7:2014 historical grid rates are provided in Table 13. If the same event occurred multiple times in an eye, only the first occurrence is counted in the table below. The safety and performance endpoint (SPE) rate is considered not exceeded if the 1-sided 95% lower CL for an AE is less than the SPE%. No cumulative adverse events exceeded the established rates according to ISO 11979-7:2014.

No persistent adverse events (adverse events in the ISO grid that are observed at the 6 month postoperative visit) were observed for 107 subjects implanted with AcrySof IQ Vivity™ IOL.

Table 13: Cumulative and Persistent Serious Adverse Events and SPE Rates for AcrySof IQ Vivity™ IOL Group (All Implanted)

Cumulative Serious Adverse Events	First Implanted Eyes N = 107		Second Implanted Eyes N = 106		SPE %
	n (%)	1-sided 95% Lower CL	n (%)	1-sided 95% Lower CL	
Cystoid macular oedema	1 (0.9)	0.05	0 (0.0)	0.00	3.0
Hypopyon	0 (0.0)	0.00	0 (0.0)	0.00	0.3
Endophthalmitis	0 (0.0)	0.00	0 (0.0)	0.00	0.1
Lens dislocated from posterior chamber	0 (0.0)	0.00	0 (0.0)	0.00	0.1
Pupillary block	0 (0.0)	0.00	0 (0.0)	0.00	0.1
Retinal detachment	0 (0.0)	0.00	0 (0.0)	0.00	0.3
Secondary surgical intervention	0 (0.0)	0.00	2 (1.9)	0.34*	0.8
Other					
Visual acuity reduced transiently	1 (0.9)	0.05	0 (0.0)	0.00	N/A
Cataract operation complication	0 (0.0)	0.00	2 (1.9)	0.34	N/A
Persistent Serious Adverse Events					
Corneal stroma oedema	0 (0.0)	0.00	0 (0.0)	0.00	0.3
Cystoid macular oedema	0 (0.0)	0.00	0 (0.0)	0.00	0.5
Iritis	0 (0.0)	0.00	0 (0.0)	0.00	0.3
Raised IOP requiring treatment	0 (0.0)	0.00	0 (0.0)	0.00	0.4

*Both secondary surgical interventions were unrelated to the IOL and were performed to remove retained cataract material from the eye.

Contrast Sensitivity

Mean monocular mesopic contrast sensitivity (with and without glare) was assessed at 6 months for the AcrySof IQ Vivity™ IOL and the monofocal control IOL and results are depicted in Figures 6 and 7. Data were obtained from 107 and 111 subjects, respectively, in each arm using a backlit sine wave grating chart system (CSV1000-HGT, VectorVision, Greenville, OH). The difference in the mean monocular mesopic contrast sensitivity between AcrySof IQ Vivity™ IOL and the Monofocal Control IOL first eyes was < 0.21 log units at each of the spatial frequencies tested for both with and without glare conditions, which is below the threshold considered clinically relevant (< 0.3 log units) as defined in ISO 11979-7:2014 (Section D.8.2.2).

Figure 6: Mean Monocular Mesopic Contrast Sensitivity without Glare with Standard Deviation at 6 Months, First Eye (Safety Analysis Set)



Figure 7: Mean Monocular Mesopic Contrast Sensitivity with Glare with Standard Deviation at 6 Months, First Eye (Safety Analysis Set)

Visual Disturbances

A questionnaire developed by Alcon was used in this clinical study to assess visual disturbances. Subjects were asked to report whether they had experienced a particular visual disturbance (glare, halos, starbursts, hazy vision, blurred vision, double vision in one or both eyes, color distortion, or peripheral dark area). A subject may report multiple symptoms. As depicted in Table 14, the percentage of subjects not experiencing these visual disturbances were similar between the AcrySof IQ Vivity™ IOL group and the Monofocal Control IOL group at 6 Months.

Table 14: Comparison of Percentage of Subjects Not Experiencing Visual Disturbances for AcrySof IQ Vivity™ IOL and Monofocal Control IOL 6 months Postoperative (following second eye implantation)

Visual Disturbance	AcrySof IQ Vivity™ IOL N = 106				Monofocal Control IOL N = 113				Difference (DFT015 - SN60WF)	
	n	Total	%	95% CI	n	Total	%	95% CI	%	95% CI
Starbursts	70	106	66.0	(56.2, 75.0)	68	110	61.8	(52.1, 70.9)	4.2	(-8.8, 17.1)
Halos	78	106	73.6	(64.1, 81.7)	91	110	82.7	(74.3, 89.3)	-9.1	(-20.4, 2.0)
Glare	81	105	77.1	(67.9, 84.8)	81	111	73.0	(63.7, 81.0)	4.2	(-7.6, 15.8)
Hazy Vision	93	105	88.6	(80.9, 94.0)	96	111	86.5	(78.7, 92.2)	2.1	(-7.2, 11.3)
Blurred Vision	96	106	90.6	(83.3, 95.4)	89	111	80.2	(71.5, 87.1)	10.4	(0.6, 20.1)
Double Vision	101	106	95.3	(89.3, 98.5)	108	111	97.3	(92.3, 99.4)	-2.0	(-8.3, 3.6)
Dark Area (Negative Dysphotopsia)	95	106	89.6	(82.2, 94.7)	95	111	85.6	(77.6, 91.5)	4.0	(-5.2, 13.4)

Percentage calculated as (n / Total) * 100

Subjects who reported experiencing a particular visual disturbance (glare, halos, starbursts, hazy vision, blurred vision, double vision in one or both eyes, color distortion, or peripheral dark area) were asked to rate the severity, frequency, and bothersomeness of the disturbance. As demonstrated in Table 15, reports of the bothersomeness and severity of visual disturbances/distortions were similar between the AcrySof IQ Vivity™ IOL and the monofocal control IOL groups at 6 months.

Table 15: Comparison of Visual Disturbance Bothersomeness for AcrySof IQ Vivity™ IOL and Monofocal Control IOL
6 months Postoperative (following second eye implantation)

AcrySof IQ Vivity™ IOL N = 106						
Visual Disturbance	Total	Not at all bothered ^a %	A Little bit %	Some what %	Quite a bit %	Very much %
Starbursts	106	73.6	14.2	9.4	0.9	1.9
Halos	106	83.0	12.3	2.8	0.9	0.9
Glare	105	78.1	11.4	6.7	3.8	0.0
Hazy Vision	105	89.5	4.8	2.9	2.9	0.0
Blurred Vision	106	91.5	0.0	4.7	2.8	0.9
Double Vision	106	98.1	0.9	0.9	0.0	0.0
Dark Area	106	92.5	5.7	0.0	0.9	0.9
Monofocal Control IOL N = 113						
Visual Disturbance	Total	Not at all bothered ^a %	A Little bit %	Some what %	Quite a bit %	Very much %
Starbursts	110	71.8	14.5	10.9	1.8	0.9
Halos	110	89.1	5.5	4.5	0.9	0.0
Glare	111	75.7	9.9	11.7	2.7	0.0
Hazy Vision	111	88.3	3.6	6.3	1.8	0.0
Blurred Vision	111	80.2	9.0	7.2	1.8	1.8
Double Vision	111	99.1	0.0	0.9	0.0	0.0
Dark Area	111	90.1	5.4	2.7	0.9	0.9

^a Includes patients who did not experience the disturbance and those reporting being not at all bothered

2. AcrySof™ TORIC INTRAOCULAR LENS CLINICAL STUDY OVERVIEW

A multicenter, subject-masked, randomized, prospective clinical study was performed in the United States to evaluate the safety and effectiveness of AcrySof™ Toric IOL (Models SA60T3-SA60T5), collectively referred to as SA60TT) compared to Model SA60AT as control lens when implanted into the capsular bag following phacoemulsification. In this study, the results achieved by the subjects successfully followed for 6 months postoperatively (defined as 120 to 180 days postoperative after the second eye implant) provide reasonable assurance that the AcrySof™ Toric IOL is a safe and effective device for the visual correction of aphakia and pre-existing corneal astigmatism following cataract surgery.

Only data from the first operative eye from those subjects who received either a Model SA60TT or Model SA60AT intraocular lens are included. Three different lens models of varying cylinder correction were evaluated in this clinical study. The three different models evaluated and their applicable cylinder powers are listed in Table 16.

Table 16: AcrySof™ Toric IOL: IOL Models and Cylinder Power in Clinical Study

IOL Model™	Cylinder Power		Recommended Corneal Astigmatism Correction Ranges
	at IOL plane	at corneal plane	
SA60T3	1.50	1.03	0.75 - 1.50 D
SA60T4	2.25	1.55	1.50 - 2.00 D
SA60T5	3.00	2.06	2.00 D & up

[†]These IOL models are collectively referred to as SA60TT in the text that follows.

The recommended corneal astigmatism correction ranges are based on 1) the preoperative corneal astigmatism and 2) the predicted effect of 0.5 diopter surgically induced astigmatism for a standardized temporal incision. The combination of these two parameters is used in Alcon provided software to select the appropriate Toric IOL model and recommended axis of placement. As such, the recommended range of corneal astigmatism to be corrected while not identical, is directly related to, the preoperative keratometric cylinder.

Subject Population

Specific requirements for study participation included (1) ≥0.75 D of preoperative With-the-Rule or preoperative oblique regular corneal astigmatism and (2) ≥1.00 D of preoperative Against-the-Rule regular corneal astigmatism. The subject population implanted with a Model SA60TT in the first operative eye consists of 53.3% females and 46.7% males. The subject population implanted with the Model SA60AT (control) intraocular lens consists of 57.2% females and 42.8% males. Stratifying by race for the Model SA60TT population, 97.6% are Caucasian, 2.0% are Black, and 0.4% are other. The control (SA60AT) population is 95.6% Caucasian, 1.6% Black, 1.2% Asian and 1.6% other. The mean age for the population receiving the Model SA60TT and SA60AT were 70.0 and 72.4, respectively.

Uncorrected Distance Visual Acuity (UCDVA)


A summary of uncorrected distance visual acuity achieved for Models SA60TT and SA60AT at six months postoperatively (Form 5 visit) is presented in Table 17.

Table 17: Uncorrected Distance Visual Acuity, Status at Form 5
Lens Models SA60TT and SA60AT, All Implanted

	Sample size	Acuity										20/40 or better	
		20/20 or better		20/25		20/32		20/40		Worse than 20/40			
		N	n	%	n	%	n	%	n	%	n		
Total (SA60TT)	211	81	38.4	59	28.0	36	17.1	22	10.4	13	6.2	198	93.8
Total (SA60AT)	210	40	19.0	46	21.9	37	17.6	39	18.6	48	22.9	162	77.1

At six months postoperatively (Form 5 visit), shown in Figure 8, 93.8% of Model SA60TT subjects achieved 20/40 or better UCDVA (first operative eye of the All Implanted data set) compared to 77.1% of the subjects implanted with the control Model SA60AT. The difference in UCDVA between Models SA60TT and SA60AT was statistically significant (all p-values ≤0.0001) in favor of Model SA60TT.

Figure 8: Cumulative UCDVA, Status at Form 5, Model SA60TT vs. Control

Similar results were noted when data was analyzed by cylinder range:

- The difference in cumulative UCDVA between Models SA60T3 and SA60AT was statistically significant (all p-values ≤0.0115) for each visual acuity category (20/20 or better, 20/25 or better, 20/32 or better and 20/40 or better) in favor of Model SA60T3.
- The difference in cumulative UCDVA between Models SA60T4 and SA60AT was statistically significant (all p-values ≤0.0082) for each visual acuity category (20/25 or better, 20/32 or better and 20/40 or better) in favor of Model SA60T4, with the exception of the 20/20 or better category.
- The difference in cumulative UCDVA between Models SA60T5 and SA60AT was statistically significant (all p-values ≤0.0171) for each visual acuity category (20/20 or better, 20/25 or better, 20/32 or better and 20/40 or better) in favor of Model SA60T5.

Best Spectacle Distance Corrected Visual Acuity (BSCDVA)

A summary of best spectacle corrected distance visual acuity (BSCDVA) achieved at six months postoperatively (Form 5 visit) among subjects who did not have any visually significant preoperative pathology or macular degeneration at any time (Best Case) is presented in Table 18A. Visual acuity achieved by the overall subject population is shown in Table 18C. Control data are found for the same data sets in Tables 18B and 18D, respectively.

Of the first operative eyes implanted with a Model SA60TT and examined at six months postoperatively (Form 5 visit), 100.0% achieved a BSCDVA of 20/40 or better in the Best Case dataset. These rates exceed the FDA grid rates of 96.7%.

Table 18A: BSCDVA by Age Category, Status at Form 5 - Lens Model SA60TT, Best Case

Age Category	Sample size	Acuity										20/40 or better	
		20/20 or better		20/25		20/32		20/40		Worse than 20/40			
		N	n	%	n	%	n	%	n	%	n		
<60	29	27	93.1	1	3.4	1	3.4	0	0	0	0	29	100.0
60-69	51	42	82.4	7	13.7	2	3.9	0	0	0	0	51	100.0
70-79	73	57	78.1	13	17.8	3	4.1	0	0	0	0	73	100.0
≥80	20	14	70.0	4	20.0	1	5.0	1	5.0	0	0	20	100.0
Total	173	140	80.9	25	14.5	7	4.0	1	0.6	0	0	173	100.0

Table 18B: BSCDVA by Age Category, Status at Form 5 - Lens Model SA60AT, Best Case

Age Category	Sample size	Acuity										20/40 or better	
		20/20 or better		20/25		20/32		20/40		Worse than 20/40			
		N	n	%	n	%	n	%	n	%	n	%	n
<60	15	13	86.7	1	6.7	1	6.7	0	0	0	0	15	100.0
60-69	49	38	77.6	11	22.4	0	0	0	0	0	0	49	100.0
70-79	75	48	64.0	21	28.0	6	8.0	0	0	0	0	75	100.0
≥80	32	19	59.4	8	25.0	2	6.3	3	9.4	0	0	32	100.0
Total	171	118	69.0	41	24.0	9	5.3	3	1.8	0	0	171	100.0

Of the first operative eyes implanted with a Model SA60TT and examined at six months postoperatively (Form 5 visit), 100.0% achieved a BSCDVA of 20/40 or better in the All Implanted dataset. These rates exceed the FDA grid rates of 92.5%.

Table 18C: BSCDVA by Age Category, Status at Form 5 - Lens Model SA60TT, All Implanted

Age Category	Sample size	Acuity										20/40 or better	
		20/20 or better		20/25		20/32		20/40		Worse than 20/40			
		N	n	%	n	%	n	%	n	%	n		%
<60	33	30	90.9	2	6.1	1	3.0	0	0	0	0	33	100.0
60-69	56	47	83.9	7	12.5	2	3.6	0	0	0	0	56	100.0
70-79	90	72	80.0	15	16.7	3	3.3	0	0	0	0	90	100.0
≥80	32	22	68.8	5	15.6	4	12.5	1	3.1	0	0	32	100.0
Total	211	171	81.0	29	13.7	10	4.7	1	0.5	0	0	211	100.0

Table 18D: BSCDVA by Age Category, Status at Form 5 - Lens Model SA60AT, All Implanted

Age Category	Sample size	Acuity										20/40 or better	
		20/20 or better		20/25		20/32		20/40		Worse than 20/40			
		N	n	%	n	%	n	%	n	%	n	%	n
<60	15	13	86.7	1	6.7	1	6.7	0	0	0	0	15	100.0
60-69	54	41	75.9	12	22.2	1	1.9	0	0	0	0	54	100.0
70-79	91	59	64.8	22	24.2	10	11.0	0	0	0	0	91	100.0
≥80	49	28	57.1	13	26.5	2	4.1	3	6.1	3	6.1	46	93.9
Total	209	141	67.5	48	23.0	14	6.7	3	1.4	3	1.4	206	98.6

Figures 9A–9C show a summary of cumulative best corrected visual acuities for each Toric model compared to the control subjects in the same cylinder range for the All Implanted dataset.

Figure 9A: Cumulative BSCDVA, Model SA60T3 vs. Control, Form 5, All Implanted


Figure 9B: Cumulative BSCDVA, Model SA60T4 vs. Control, Form 5, All Implanted



Figure 9C: Cumulative BSCDVA, Model SA60T5 vs. Control, Form 5, All Implanted

Absolute Residual Refractive Cylinder

Figures 10A–10C demonstrate that residual refractive cylinder values were statistically significantly lower among those subjects implanted with an AcrySof™ Toric IOL Model SA60T3, SA60T4, or SA60T5 when compared to the corresponding subjects implanted with the control Model SA60AT. Subjects implanted with an AcrySof™ Toric Model SA60T3 showed a 62.4% mean reduction in refractive cylinder from the preoperative visit (keratometric cylinder) as compared to the 10.8% mean reduction for subjects implanted with the concurrent control Model SA60AT. Subjects implanted with an AcrySof™ Toric IOL Model SA60T4 or SA60T5 showed similar results with a mean reduction in refractive cylinder of 54.8 % and 67.8%, respectively, as compared to subjects implanted with the concurrent control model who had a mean reduction in refractive cylinder of 22.1% and 27.7%, respectively. Each of the AcrySof™ Toric Lens Models SA60T3, SA60T4 and SA60T5 had at least a 3-fold increase in the likelihood of achieving residual refractive cylinder of 0.5 D or less as compared to the corresponding control model.

Figure 10A: AcrySof™ Toric IOL: Absolute Residual Refractive Cylinder, Model SA60T3 vs. Control, 6 Month Visit (Form 5), All Implanted


Figure 10B: AcrySof™ Toric IOL: Absolute Residual Refractive Cylinder, Model SA60T4 vs. Control, 6 Month Visit (Form 5), All Implanted



Figure 10C: AcrySof™ Toric IOL: Absolute Residual Refractive Cylinder, Model SA60T5 vs. Control, 6 Month Visit (Form 5), All Implanted

Stability of Cylinder

Subjects implanted with lens Model SA60TT exhibited stability of cylinder at Form 4 (3 months) with greater than 90% of all subjects changing less than or equal to 1.00 diopter at consecutive visits between Form 3 (one month) and Form 6 (twelve months) as described in Tables 19A–19D.

Table 19A: AcrySof™ Toric IOL: Stability of Cylinder
(Eyes that had 2 consecutive exams, but not necessarily every follow-up exam)

Recommended Corneal Astigmatism Correction Ranges	Toric IOL Model	Magnitude of Vector Change in Cylinder	1 and 3 Months n/N,%	3 and 6 Months n/N,%	6 and 12 Months n/N,%
< 1.5 D	SA60T3	≤ 1.00 D	106/107,99.07%	101/105,96.19%	55/55,100.00%
		Mean Change	0.28	0.29	0.20
		SD	0.32	0.33	0.25
≥ 1.5 - < 2.0 D	SA60T4	≤ 1.00 D	54/56,96.43%	53/54,98.15%	25/27,92.59%
		Mean Change	0.40	0.27	0.46
		SD	0.35	0.22	0.45
≥ 2.0 D	SA60T5	≤ 1.00 D	40/45,88.89%	35/40,87.50%	27/30,90.00%
		Mean Change	0.43	0.42	0.41
		SD	0.44	0.45	0.38
Combined	SA60TT	≤ 1.00 D	200/208,96.15% (93.54,98.77)	189/199,94.97% (91.94,98.01)	107/112,95.54% (91.71,99.36)
		Mean Change	0.35	0.31	0.32
		SD	0.36	0.34	0.36
		95% CI	0.30,0.39	0.26,0.36	0.25,0.39

n/N,%,(%)CI are for percent with change between ± 1.00D

Table 19B: AcrySof™ Toric IOL: Stability of Cylinder
(Eyes that had every follow-up exam up to Form 6, 12 months)

Recommended Corneal Astigmatism Correction Ranges	Toric IOL Model	Magnitude of Vector Change in Cylinder	1 and 3 Months n/N,%	3 and 6 Months n/N,%	6 and 12 Months n/N,%
< 1.5 D	SA60T3	≤ 1.00 D	34/34,100.00%	34/34,100.00%	34/34,100.00%
		Mean Change	0.25	0.24	0.21
		SD	0.23	0.22	0.24
≥ 1.5 - < 2.0 D	SA60T4	≤ 1.00 D	17/17,100.00%	16/17,94.12%	16/17,94.12%
		Mean Change	0.27	0.25	0.35
		SD	0.25	0.26	0.33
≥ 2.0 D	SA60T5	≤ 1.00 D	17/19,89.47%	15/19,78.95%	16/19,84.21%
		Mean Change	0.44	0.56	0.52
		SD	0.47	0.50	0.43
Combined	SA60TT	≤ 1.00 D	68/70,97.14% (93.23,100.00)	65/70,92.86% (86.82,98.90)	66/70,94.29% (88.84,99.73)
		Mean Change	0.31	0.33	0.33
		SD	0.32	0.35	0.34
		95% CI	0.23,0.38	0.24,0.41	0.25,0.41

n/N,%,(%)CI are for percent with change between ± 1.00D

Table 19C: AcrySof™ Toric IOL: Stability of Absolute Cylinder
(Eyes that had 2 consecutive exams, but not necessarily every follow-up exam)

Recommended Corneal Astigmatism Correction Ranges	Toric IOL Model	Magnitude of Change in Absolute Cylinder	1 and 3 Months n/N,%	3 and 6 Months n/N,%	6 and 12 Months n/N,%
< 1.5 D	SA60T3	≤ 1.00 D	107/107,100.00%	104/105,99.05%	55/55,100.00%
		Mean Change	0.04	0.02	0.05
		SD	0.32	0.38	0.29
≥ 1.5 - < 2.0 D	SA60T4	≤ 1.00 D	54/56,96.43%	54/54,100.00%	27/27,100.00%
		Mean Change	0.18	0.05	-0.12
		SD	0.42	0.27	0.41
≥ 2.0 D	SA60T5	≤ 1.00 D	44/45,97.78%	37/40,92.50%	29/30,96.67%
		Mean Change	0.09	0.06	0.00
		SD	0.38	0.49	0.45
Combined	SA60TT	≤ 1.00 D	205/208,98.56% (96.93,100.00)	195/199,97.99% (96.04,99.94)	111/112,99.11% (97.36,100.00)
		Mean Change	0.09	0.03	-0.01
		SD	0.37	0.38	0.37
		95% CI	0.04,0.14	-0.02,0.09	-0.08,0.06

n/N,%,(%)CI are for percent with change between ± 1.00D

Table 19D: AcrySof™ Toric IOL: Stability of Absolute Cylinder
(Eyes that had every follow-up exam up to Form 6, 12 months)

Recommended Corneal Astigmatism Correction Ranges	Toric IOL Model	Magnitude of Change in Absolute Cylinder	1 and 3 Months n/N,%	3 and 6 Months n/N,%	6 and 12 Months n/N,%
< 1.5 D	SA60T3	≤ 1.00 D	34/34,100.00%	34/34,100.00%	34/34,100.00%
		Mean Change	0.01	-0.01	0.07
		SD	0.28	0.31	0.28
≥ 1.5 - < 2.0 D	SA60T4	≤ 1.00 D	17/17,100.00%	17/17,100.00%	17/17,100.00%
		Mean Change	0.06	0.19	-0.04
		SD	0.30	0.21	0.42
≥ 2.0 D	SA60T5	≤ 1.00 D	18/19,94.74%	17/19,89.47%	18/19,94.74%
		Mean Change	0.17	0.05	0.01
		SD	0.45	0.54	0.55
Combined	SA60TT	≤ 1.00 D	69/70,98.57% (95.78,100.00)	68/70,97.14% (93.23,100.00)	69/70,98.57% (95.78,100.00)
		Mean Change	0.07	0.05	0.03
		SD	0.34	0.38	0.40
		95% CI	-0.01,0.15	-0.04,0.14	-0.07,0.12

n/N,%,(%)CI are for percent with change between ± 1.00D

Adverse Events

The incidence of cumulative adverse events for the Model SA60TT compared favorably to the FDA historical grid rates. Only the rates for retinal detachment/repair and surgical reintervention exceeded the FDA historical grid (Table 20). However, neither of these rates were statistically significant (p=0.5196 and p=0.1336, respectively). No occurrences of persistent adverse events were observed in any subjects implanted with the AcrySof™ Toric IOL.

Table 20: AcrySof™ Toric IOL: Adverse Events Incidence Rates, First Eye – Safety

Cumulative Adverse Events	Model SA60TT N=244		FDA Grid Rate
	N	%	%
Retinal Detachment/Repair	1	0.4	0.3
Surgical Reintervention	4**	1.6	0.8
IOL Reposition Due to Rotation	1	0.4	NA
IOL Replacement Due to Rotation	1	0.4	NA
Laser Treatment	2	0.8	NA
Paracentesis	1	0.4	NA

The incidence rates in this table are based upon the number of eyes with an event divided by the number of eyes implanted.

Cumulative adverse events are those events that have occurred at any time during the clinical study.

FDA Grid Rate = FDA Grid of Adverse Events with Posterior Chamber Intraocular Lens Historical Controls, FDA Intraocular Lens Guidance Document, Annex B (October 14, 1999)

**There were 5 occurrences of surgical reintervention in 4 eyes for Model SA60TT first eye


The incidence of cumulative adverse events for the Model SA60TT also compared favorably to the concurrent control.

Distance-Vision Spectacle Independence

Spectacle independence was assessed in the study by direct subject responses obtained from a self-reported, binocular subject questionnaire. Since the AcrySof™ Toric IOL study was completed in 2005, the questionnaire utilized was not developed using the US FDA guidance document "Patient-Reported Outcomes Measures: Use in Medical Product Development to Support Labeling Claims" dated December 2009.

Statistically significantly more Model SA60TT subjects reported postoperative distance-vision spectacle independence compared to Model SA60AT subjects (59.9% versus 37.7%, respectively) when unilaterally implanted. Distance-vision spectacle independence was defined as the percentage of subjects who selected the "none of the time" response for distance-vision frequency-of-spectacle-wear. Spectacle dependence was defined as subjects indicating any reliance on glasses for distance-vision and represents the summation of the "some of the time", "half of the time", "most of the time" and "all of the time" frequency-of-spectacle-wear responses. Consequently, fewer Model SA60TT subjects were spectacle dependent at 40.1% compared to 62.3% of the Model SA60AT subjects. **Figure 11** illustrates the distance-vision frequency-of-spectacle-wear distributions between Model SA60TT and Model SA60AT groups. Implantation of an AcrySof™ Toric Intraocular lens in astigmatic subjects provides significantly improved distance-vision spectacle independence relative to a conventional monofocal IOL.

Figure 11: Distance-Vision Spectacle Independence: Frequency of Spectacle Wear, Form 5, All Implanted

3. STUDY RESULTS FOR Clareon ASPHERIC HYDROPHOBIC ACRYLIC IOL

A prospective, multicenter, open-label clinical study in adult subjects requiring cataract surgery with IOL implantation was performed on the Clareon™ IOL in the United States. The objective of this study was to demonstrate favorable visual acuity and adverse event outcomes for the

Clareon™ IOL compared to historical safety and performance endpoint (SPE) rates as reported in EN ISO 11979-7:2014. The results achieved by the subjects successfully followed for 12 Months postoperatively (defined as 330 to 420 days postoperative) provide reasonable assurance that the Clareon™ IOL is a safe and effective device for the visual correction of aphakia following cataract surgery. The subject population in the clinical study consisted of 350 subjects in need of cataract surgery that were unilaterally implanted, the percentage of females (60.9%) exceeded the percentage of males (39.1%) enrolled. The mean age for the total population was 69.7 years.

Visual Acuity

The best corrected distance visual acuity (BCDVA) achieved at one year postoperatively by the Best Case Analysis Set is presented in **Table 21** and visual acuity achieved by All-Implanted Analysis Set is presented in **Table 22**. Similar data for BCDVA expressed as Snellen values are presented in **Tables 23** and **24**. The Best Case analysis included all eyes successfully implanted with the test article, had at least one postoperative visit, did not have any ocular pathology at Visit 0 (preoperative), had no macular degeneration at any visit, and did not have any previous surgery for the correction of refractive errors. All-Implanted Analysis set includes all eyes with successful test article implantation.

The one-sided exact 95% upper confidence limit for the percentage of subjects with monocular BCDVA of 0.3 logMAR or better at 12 months postoperative (Visit 5) is greater than or equal to the SPE rates as reported in EN ISO 11979-7:2014.

Table 21: Best Corrected Visual Acuity in the Best Case Analysis Set at Month 12, Clareon™ IOL

Visual Acuity	n	%
0.0 logMAR or better	271	83.1
0.1 logMAR or better	307	94.2
0.2 logMAR or better	323	99.1
0.3 logMAR or better	325	99.7
Worse than 0.3 logMAR	1	0.3
ISO SPE Rate for % of 0.3 logMAR or better	96.7%	
N	326	
ISO SPE Rate = BS EN ISO 11979-7:2014, Ophthalmic Implants - Intraocular Lenses - Part 7: Clinical Investigations, Table B.4 - Posterior Chamber IOL Adverse Event Safety and Performance Endpoint Rates		

Table 22: Best Corrected Visual Acuity in the All-Implanted Analysis Set at Month 12, Clareon™ IOL

Visual Acuity	n	%
0.0 logMAR or better	277	81.0
0.1 logMAR or better	318	93.0
0.2 logMAR or better	339	99.1
0.3 logMAR or better	341	99.7
Worse than 0.3 logMAR	1	0.3
ISO SPE Rate for % of 0.3 logMAR or better	92.5%	
N	342	
ISO SPE Rate = BS EN ISO 11979-7:2014, Ophthalmic Implants - Intraocular Lenses - Part 7: Clinical Investigations, Table B.3 - Posterior Chamber IOL Adverse Event Safety and Performance Endpoint Rates		

Table 23: Best Corrected Visual Acuity in the Best Case Analysis Set at Month 12 (Snellen), Clareon™ IOL

Visual Acuity	n	%
20/20 or better	287	88.0
20/25 or better	316	96.9
20/32 or better	324	99.4
20/40 or better	325	99.7
Worse than 20/40	1	0.3
N	326	
Note: 20/20 = 0.04 logMAR 20/40 = 0.34 logMAR	20/25 = 0.14 logMAR	20/32 = 0.24 logMAR

Table 24: Best Corrected Visual Acuity in the All-Implanted Analysis Set at Month 12 (Snellen), Clareon™ IOL

Visual Acuity	n	%
20/20 or better	297	86.8
20/25 or better	331	96.8
20/32 or better	340	99.4
20/40 or better	341	99.7
Worse than 20/40	1	0.3
N	342	
Note: 20/20 = 0.04 logMAR 20/40 = 0.34 logMAR	20/25 = 0.14 logMAR	20/32 = 0.24 logMAR

IOL Rotation

The rotational stability of the Clareon™ IOL was evaluated in a subset of the total study population, which included 141 subjects at 6 sites. Results were compared to performance targets in EN ISO 11979-7:2014, defined as the difference between IOL axis orientation on day of surgery and at the Month 6 visit. A total of 98.4% of subjects had IOL rotation <10 degrees and 99.2% had rotation <20 degrees or <30 degrees at Month 6 compared to the day of surgery, which met the pre-specified performance targets. Results shown in **Table 25**.

Table 25: Categorical Statistics for Absolute IOL Rotation (10 Degrees Increment) at Month 6 from Day of Surgery, Rotation Analysis Set, Clareon™ IOL

(N=141)		
Category	n	%
Total	124	-
Less than 10 degrees	122	98.4
Less than 20 degrees	123	99.2
Less than 30 degrees	123	99.2
Greater than 30 degrees	1	0.8
N = Number of eyes in the analysis set n = Number of eyes in specified category Total = Number of eyes with data Percentages are calculated as (n/Total) * 100		

In an additional analysis between Month 1 and Month 6, 98.4% of subjects implanted with the Clareon™ IOL had an absolute rotation of less than or equal to 5 degrees. The remaining subjects (1.6%) had rotation between 5 to 10 degrees.

Absolute rotation of Clareon™ IOL were also described at all post-operative visits (1 Day, 1 Week, 1 Month, and 6 Months) compared to the day of surgery. Descriptive statistics results at each postoperative visit are shown in **Table 26**.

Table 26: Descriptive Statistics for Absolute IOL Rotation (Degrees), Rotation Analysis Set, Clareon™ IOL

Visit	Statistic	(N = 141)*
1 Day	n	127
	Mean (SD)	1.8483 (3.72254)
	Median	0.942
	(Min, Max)	(0.002, 38.232)
1 Week	95% CI	(1.1946, 2.5020)
	n	126
	Mean (SD)	1.9942 (3.86948)
	Median	1.213
1 Month	(Min, Max)	(0.020, 40.547)
	95% CI	(1.3119, 2.6764)
	n	127
	Mean (SD)	2.2335 (3.95863)
6 Months	Median	1.414
	(Min, Max)	(0.008, 41.110)
	95% CI	(1.5384, 2.9287)
	n	124
	Mean (SD)	2.2696 (3.87231)
	Median	1.449
	(Min, Max)	(0.014, 40.033)
	95% CI	(1.5812, 2.9579)
N = Number of eyes in the analysis set n = Number of eyes at visit with data SD = Standard Deviation, CI = Confidence Interval * Of 141, 10 had images at baseline which were not readable and at each visit, individual subjects may have had images which were missing or not readable		

Adverse Events

The cumulative and persistent adverse events rates for the Safety Analysis Set are presented in **Table 27**. Cumulative adverse events are those that occur at any point in the study while persistent adverse events are those that occur or are ongoing at the 12 Month Visit. The Safety Analysis Set included all eyes with attempted implantation with the test article (successful or aborted after contact with the eye). In this study, the number of eyes in the Safety Analysis Set was equivalent to the All-Implanted Analysis Set.

For both cumulative and persistent adverse events, the exact one-sided 95% lower confidence limit of cumulative and persistent adverse events were within limits compared to the SPE rates reported in EN ISO 11979-7:2014.

Table 27: Cumulative and Persistent Adverse Events in the Safety Analysis Set, Clareon™ IOL

	(N = 350)			
	n (%)	2-sided 95% CI	1-sided 95% Lower CL	SPE %
Cumulative Serious Adverse Events				
Cystoid macular oedema	3 (0.9)	(0.18, 2.48)	0.23	3.0
Hypopyon	0 (0.0)	(0.00, 1.05)	0.00	0.3
Endophthalmitis	0 (0.0)	(0.00, 1.05)	0.00	0.1
Lens dislocated from posterior chamber	0 (0.0)	(0.00, 1.05)	0.00	0.1
Pupillary block	0 (0.0)	(0.00, 1.05)	0.00	0.1
Retinal detachment	0 (0.0)	(0.00, 1.05)	0.00	0.3
Secondary surgical intervention*	6 (1.7)	(0.63, 3.69)	0.75	0.8
Other				
Herpes virus infection	1 (0.3)	(0.01, 1.58)	0.01	NA
Macular fibrosis	1 (0.3)	(0.01, 1.58)	0.01	NA
Macular hole	1 (0.3)	(0.01, 1.58)	0.01	NA
Punctate keratitis	1 (0.3)	(0.01, 1.58)	0.01	NA
Retinal tear	2 (0.6)	(0.07, 2.05)	0.10	NA
Persistent Serious Adverse Events				
Corneal stroma oedema	0 (0.0)	(0.00, 1.05)	0.00	0.3
Cystoid macular oedema	1 (0.3)	(0.01, 1.58)	0.01	0.5
Iritis	0 (0.0)	(0.00, 1.05)	0.00	0.3
Raised IOP requiring treatment	0 (0.0)	(0.00, 1.05)	0.00	0.4

CI = Confidence Interval
CL = Confidence Limit
SPE = Safety and Performance Endpoint
Persistent = Present or ongoing at the final scheduled visit
IOP = Intraocular Pressure
SPE rates are from BS EN ISO 11979-7:2014, Ophthalmic Implants - Intraocular Lenses - Part 7: Clinical Investigations, Table B.2 - Posterior Chamber IOL Adverse Event Safety and Performance Endpoint Rates
If an eye has multiple occurrences of an AE, the eye is presented only once in the respective eye count column (n) for the corresponding AE.
"Other" includes the MedDRA Preferred Term for ocular SAEs that do not belong to any predefined SPE categories.
Percentages are calculated as (n/N) * 100.
The SPE rate is considered not exceeded if the 1-sided 95% lower CL for an AE is less than the SPE%.
*None of these secondary surgical interventions were related to the IOL.

IOL Observations

At all five scheduled postoperative visits and at unscheduled visits, slit-lamp examinations were performed. The list of pre-specified slit-lamp findings included IOL observations, which were described as IOL glistenings, scratches/cracks on the IOL, and surface haze on the IOL. Across 1,852 postoperative (including unscheduled) visits, no IOL observations (IOL glistenings, scratches/cracks on the IOL, and surface haze on the IOL) were noted for the Clareon™ IOL.

REFERENCES

- Abulafia A, Barrett GD, Kleinmann G, Ofir S, Levy A, Marcovich AL, *et al.* Prediction of refractive outcome with toric intraocular lens implantation. *J Cataract Refractive Surgery*. 2015; 41: 936-944.
- Abulafia A, Hill WE, Franchina M, Barrett GD. Comparison of methods to predict residual astigmatism after intraocular lens implantation. *J Refractive Surgery*. 2015b; 31(10):699-706.
- Boettner, E.A. and Wolter, J.R. Transmission of the ocular media. *Invest. Ophthalmol*. 1962;1: 776-83.
- Hoffer KJ. The Hoffer Q formula: a comparison of theoretic and regression formulas. *J Cataract Refract Surg*. 1993;19(6):700-12.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry, and intraocular lens power calculations. *J Cataract Refract Surg*. 1997;23(9):1356-70.
- Olsen T. Calculation of intraocular lens power: a review. *Acta Ophthalmol Scand*. 2007;85(5): 472-85.
- Retzlaff JA, Sanders DR, Kraff M. *Lens Implant Power Calculation*. 3rd ed. Thorofare (NJ): Slack, Inc.; 1990.
- Haigis W. ULIB - User group for Laser Interference Biometry. Updated 2014 Jul 17; cited 2015 Mar 13. Available at: <http://ocusoft.de/ulib/>.

STAND ALONE SYMBOLS FROM ISO 7000/ ISO 7001* USED ON LABELING
(ISO 7000/IEC 60417* Title: Graphical Symbols for Use on Equipment)
(*ISO 7001 Title: Graphical symbols – Public information symbols)

Symbol	Reference Number from ISO 7000 / IEC 60417* / ISO 7001*	Symbol Title / Explanatory Text
	1051	Do not re-use
	2608	Do not re-sterilize
	2607	Use-by date
	2501	Sterilized using ethylene oxide
	2498	Serial number
	2493	Catalogue number
	0434A	Caution
	3082	Manufacturer
	0533	Upper Limit of Temperature
	1641	Consult instructions for use or consult electronic instructions for use
	2606	Do not use if package is damaged and consult instructions for use
	2497	Date of Manufacture
	3079	Open Here
	3010	RFID tag, general
	6050*	Model number
	6049*	Country of manufacture
	5662*	Date
	5664*	Patient identification
	PI PF 044*	Health care centre or doctor
	PI PF 002*	Hospital
	3707	Single sterile barrier system
	3705	Patient information website

*This symbol is from ISO 7001.

*This symbol is from IEC 60417.

STAND ALONE SYMBOL FROM ASTM F2503-13 USED ON LABELING
(Title: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment)

Symbol	Reference Number from ASTM F2503-13	Symbol Title / Explanatory Text
	N/A	MR safe

ABBREVIATIONS or SYMBOLS USED ON LABELING

Symbol	Symbol Title / Explanatory Text
	Medical device
	Unique device identifier
IOL 	Intraocular lens
OVD	Ophthalmic viscosurgical device
	UV and Blue Light Filter with wavefront-shaping X-WAVE™ technology (i.e., Extended Vision)
	Eye
	Posterior chamber IOL
UV	Ultraviolet
D	Diopter
$\varnothing_B$	Body diameter (Optic diameter)
$\varnothing_T$	Overall diameter (Overall length)
L	Left
R	Right
CYL	Cylinder Add Power
PWR	Spherical Equivalent Power
	D-size nozzle for MONARCH™ Delivery System cartridge*
	C-size nozzle for MONARCH™ Delivery System cartridge*
	B-size nozzle for MONARCH™ Delivery System cartridge*
	CLAREON™ MONARCH™ IV
	Not made with natural rubber latex
	Does not contain PHT (phthalates)
	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician
	Authorized Representative in the European Community

*The recommendation shown on the labeling is for the smallest qualified cartridge nozzle size per diopter. For all qualified combinations of products, please review Table 2.

†Trademarks are the property of their respective owner

Manufacturer:
Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, Texas 76134-2099 USA

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© 2023, 2024 Alcon Inc.