Medtronic

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1 INTENDED USE

For use in determining the activated clotting time of whole blood.
For in vitro diagnostic use.

2 SUMMARY

The Activated Clotting Time (ACT) test was originally described by Hattersley in 1966.⁷ The Medtronic High Range Activated Clotting Time (HR-ACT) Cartridge utilizes kaolin suspended in a buffer as the activating agent. The test and its variations have been widely used to monitor the anticoagulant effect of heparin.^{2, 3}

3 PRINCIPLES

The activated clotting time test is a functional test that measures the clotting time of fresh whole blood. Activation of Factor XII is accomplished by addition of a contact activator such as silica, diatomaceous earth or kaolin. The test responds to increased heparin concentrations in a linear fashion. The end point of the test is the detection of clot formation. The Medtronic HMS Plus Hemostasis Management System detects clot formation by measuring the rate of fall of the plunger contained in each cartridge channel. The plunger assembly falls rapidly through an unclotted sample. A fibrin web formed during clotting impedes the rate of plunger assembly descent. Clot formation is detected by a photo optical sensor located in the actuator assembly of the instrument. The Medtronic HMS Plus instrument performs simultaneous duplicate tests using a two-channel cartridge. The clotting times for each channel are displayed in seconds following completion of the test.

4 REAGENTS

Each cartridge reagent chamber contains 100 µL of kaolin suspended in HEPES buffer, calcium chloride, and sodium azide as a bacteriostatic agent. The final concentration of kaolin in the cartridge is 12%.

5 STORAGE

Room Temperature

HR-ACT Cartridges are designed for room temperature storage (2°C to 25°C) until the expiration date on the cartridge box. Cartridges should be stored in the original box to reference the expiration date. Cartridges should be disposed of if the expiration date is unknown. Do not use if the cartridge appears discolored, cracked or disfigured; the reagent appears to have evaporated or to have been contaminated; or the expiration date has passed.

Note: Do not use if there is evidence of freezing.

6 INSTRUMENT

Medtronic HR-ACT Cartridges are designed specifically for use in the Medtronic HMS Plus instrument. For operating instructions, consult the HMS Plus instrument operator's manual.

7 SPECIMEN COLLECTION AND PREPARATION

Blood samples that are improperly drawn and handled can yield erroneous results; samples must be collected carefully and handling techniques must be consistent.

If blood is obtained by venipuncture, the venipuncture should be clean and the first several milliliters of blood discarded to avoid contamination of the sample with tissue thromboplastin. Blood should not be taken from heparinized indwelling catheters or other anticoagulated lines without thoroughly flushing the line. Baseline samples should be tested within 60 seconds of collection. Heparinized samples should be tested within 2 minutes of collection.

8 PROCEDURE

This product is intended for use only by trained individuals.

Materials supplied:

- HR-ACT Cartridges for use with the HMS Plus instrument
- 3 mL sterile single-use Monoject™ syringe
- 19 ga, 1 7/16 in (36 mm) sterile single-use blunt needle

Materials required, not supplied:

- Test sample, fresh whole blood
- Medtronic HMS Plus instrument

Steps:

1. Prewarm the cartridge in the HMS Plus instrument actuator heat block for 3 to 5 minutes before collecting the test sample.
2. Tap the cartridge to resuspend the kaolin reagent.
3. Using the supplied syringe, collect at least 1.5 mL of test sample.
4. Attach and prime a clean blunt tip needle.
5. Insert the syringe into the dispenser of the HMS Plus instrument. Lock the syringe in place.
6. Press **[START/STOP]** to initiate the test. The dispenser will automatically deliver 400 µL of test sample into each cartridge channel.
7. If dispenser is not functioning, 400 µL of test sample may be delivered manually into each cartridge channel. Consult the HMS Plus instrument operator's manual for complete instructions.

Note: It is important to accurately deliver the test sample into each cartridge channel.

9 RESULTS

The response of the activated clotting time to heparin varies considerably from individual to individual.^{2,3,6} Various drugs affect the activated clotting time.¹⁴ A number of other factors can affect the response of an individual's activated clotting time to heparin (for example, antithrombin III levels, heparin potency, coagulation factor deficiencies, sample activation, consumptive coagulopathies, excessive sample dilution, sample temperature). These factors should be taken into account when interpreting the results of activated clotting time tests.^{8,10,12}

Duplicate channel results should fall within 10% of each other for baseline (unheparinized) samples and within 12% of each other for extended (anticoagulated) samples.

The **ACT Results** screen displays the clotting time for each cartridge channel and the average result.

10 LIMITATIONS

Medtronic HR-ACT Cartridges are for in vitro diagnostic use only and are to be used only with the Medtronic HMS Plus instrument.

Cartridges that have been used for testing should be considered hazardous. It is the responsibility of the user to dispose of the cartridge in accordance with local regulations and institution procedures.

To obtain accurate and consistent activated clotting times, the following must be observed:

- Both the instrument heat block and cartridge temperature should be at 37°C ± 0.5°C.
- Baseline samples should be tested within 60 seconds of collection. Heparinized samples should be tested within 2 minutes of collection.
- Patient diagnosis and medications should be noted. Medications can alter clotting times. Tests should be repeated if unexplained values are obtained. If the results are abnormal and no cause is apparent, a citrated sample may be collected and sent to the clinical laboratory for further analysis.
- To optimize precision, all technique variables should be consistent from test to test.

4 English Instructions for Use

- A prolonged activated clotting time is not specific for heparin and may be due to heparin, other anticoagulant drugs (warfarin, direct thrombin inhibitors), hemodilution, administration of large volumes of citrated blood products, some antiplatelet drugs, hypothermia, severe thrombocytopenia, markedly abnormal platelet function, or other coagulopathies.¹¹
- Excessive blood dilution during bypass can lead to postsurgically hemodiluted patients. This dilution can be great enough to extend the activated clotting time.¹¹ In severe instances, the dilution may lead to an unmeasurable clotting time.

For HMS Plus instrument precautions and limitations, consult the HMS Plus instrument operator's manual.

11 EXPECTED VALUES

Response to an increased concentration of heparin may vary between individuals. Each laboratory/facility should establish its own reference range based on its local population.

12 SPECIFIC PERFORMANCE CHARACTERISTICS

HR-ACT Cartridges in the clotting time range of 0 to 600 seconds typically do not exceed a variation of 12% between duplicate channels for an individual. In subjects with extreme sensitivity to heparin, an activated clotting time in excess of 600 seconds is generally not considered adequately reliable to predict heparin effect.

13 QUALITY CONTROL

HR-ACT Cartridge liquid controls are used to verify instrument and cartridge performance. Results should fall within the specified ranges listed by the manufacturer, or ranges determined by the user.

Note: When using only liquid controls, 2 levels must be performed every 8 hours of patient testing.

The HEPTrac Electronic Quality Control is used to supplement liquid controls and does not replace them. Refer to the current regulatory agency guidelines on the acceptability and use of electronic controls.

14 DISCLAIMERS OF WARRANTY

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

DISCLAIMER OF WARRANTY

ALTHOUGH THE HR-ACT CARTRIDGES, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

DISCLAIMER OF WARRANTY

ALTHOUGH THE HR-ACT CARTRIDGES, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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

**REFERENCES / RÉFÉRENCES / LITERATUR /
BIBLIOGRAFÍA / BIBLIOGRAFIA**

1. Bowie JE, Kemna GD. Automated management of heparin anticoagulation in cardiovascular surgery. *Proc Amer Acad Perfus.* 1985;6:85.
2. Bull BS, Korpman RA, Huse WM, Briggs BD. Heparin therapy during extra-corporeal circulation. I. Problems inherent in existing heparin protocols. *J Thorac Cardiovasc Surg.* 1975;69:674-684.
3. Bull BS, Huse MH, Brauer FS, Korpman RA. Heparin therapy during extra-corporeal circulation. II. The use of a dose response curve to individualize heparin and protamine dosage. *J Thorac Cardiovasc Surg.* 1975;69:685-689.
4. Clinical and Laboratory Standards Institute: *Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines—Fifth Edition*; H21-A5; 2008.
5. Cohen EJ, Camerlengo LJ, Dearing JP. Activated clotting time and cardiopulmonary bypass. I. The effect of hemodilution and hypothermia upon activated clotting time. *J Extracorp Tech.* 1980;12:139.
6. Effeney DJ, Goldstone J, Chin D, Krupski WC, Ellis RJ. Intraoperative anticoagulation in cardiovascular surgery. *Surgery.* 1981;90:1068-1074.
7. Hattersley PG. Activated coagulation time of whole blood. *JAMA.* 1966;196:436.
8. Henry JB. *Clinical Diagnosis and Management by Laboratory Methods.* Vol 1. 16th ed. Philadelphia, PA: WB Saunders Company; 1979.
9. Kemna GD, Bowie JE. Heparin protamine management for cardiopulmonary bypass patients with abnormal heparin dose response. *Proc Amer Acad Perfus.* 1985;6:89.
10. Moorehead MT, Westengard JC, Bull BS. Platelet involvement in the activated coagulation time of heparinized blood. *Anesth Analg.* 1984;63:394-398.
11. Niinikoski J, Laato M, Laaksonen V, Jalonen J, Ingerg MV. Use of activated clotting time to monitor anticoagulation during cardiac surgery. *Scand J Thor Cardiovasc Surg.* 1984;18:57-61.
12. Reno WJ, Rotman M, Grumbine FC, Dennis LH, Mohler ER. Evaluation of the BART test (a modification of the wholeblood activated recalcification time test) as a means of monitoring heparin therapy. *Am J Clin Path.* 1974;61:78-84.
13. Stenbjerg S, Berg E, Albrechtsen OK. Evaluation of the activated whole blood clotting time (ACT) in vitro. *Scand J Haematol.* 1979;23:239.
14. Young DS, Pastaner LC, Gibberman V. Effects of drugs on clinical laboratory tests. *Clin Chem.* 1975;21:1D.

Europe

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH - 1131 Tolochenaz
Switzerland
Internet: www.medtronic.co.uk
Tel. 41-21-802-7000
Fax 41-21-802-7900

EC REP

Authorized Representative in the European Community

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Asia-Pacific

Japan

Medtronic Japan
Comodio Shiodome 5F
2-14-1 Higashi-Shimbashi, Minato-ku
Tokyo 105-0021
Japan
Tel. 81-3-6430-2011
Fax 81-3-6430-7140

Australia

Medtronic Australasia Pty. Ltd.
97 Waterloo Road
North Ryde NSW 2113
Australia
Tel. 61-2-9857-9000
Fax 61-2-9878-5100

Asia

Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2891-4068
Fax 852-2591-0313

Medtronic

Americas

Latin America

Medtronic Latin America
3750 NW 87th Avenue
Suite 700
Miami, FL 33178
USA
Tel. 305-500-9328
Fax 786-709-4244

Canada

Medtronic of Canada Ltd.
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
Tel. 905-460-3800
Fax 905-826-6620
Toll-free: 1-800-268-5346

United States

Manufacturer:
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
Internet: www.medtronic.com
Tel. 763-514-4000
Fax 763-391-9100
Toll-free: 1-800-328-2518
(24-hour consultation service)

Medtronic Perfusion Systems

7611 Northland Drive
Minneapolis, MN 55428
USA
Tel. 763-391-9000
Fax 763-391-9100
Customer Service and Product Orders
Toll-free: 1-800-854-3570
Visit Medtronic at
www.perfusion.medtronic.com

© 2003, 2011 Medtronic

M947894A001 Rev. 1A

* M 9 4 7 8 9 4 A 0 0 1 *

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Medtronic® and HEPtrac® are registered trademarks of Medtronic, Inc.
Monoject™ is a trademark of COVIDIEN AG.