<b>REF</b>	260177, (1.60 m) 155500, (1.88 m) 260179, (3.1 m)
<b>QTY</b>	Quantity Menge Quantité Cantidad Quantità Quantidade 数量 数量 Количество Miktar
	Follow the Instructions for use. Gebrauchsanweisung beachten. Respecter les instructions d'utilisation. Siga las instrucciones de uso. Siga as Instruções de uso. Seguire le istruzioni per l'uso. 請按照說明使用。 使用説明書の説明に従ってください。 Следуйте инструкциям по эксплуатации. Kullanım talimatlarını izleyin.
	Do not reuse. Nicht wiederverwenden. Ne pas réutiliser. De un solo uso. Non riutilizzare. Não reutilizar. 請勿重複使用。 再使用しないでください。 Не использовать повторно. Tekrar kullanmayın.
	Applicable to neonatal patient group. Für die Patientengruppe „Neonaten“. Applicable au groupe de patients néonataux. Applicable al grupo de pacientes neonatos. Applicable per il gruppo di pazienti neonatale. Aplicável ao grupo de pacientes neonatais. 适用于新生儿病人组。 新生児患者グループに対して適用可能。 Применяется к группе пациентов «Младенцы». Neonatal hasta grubu için uygundur.
	Do not use if the packaging is damaged. Bei beschädigter Verpackung nicht verwenden. Ne pas utiliser si l'emballage est endommagé. No lo utilice si el envase está dañado. Non utilizzare se la confezione è danneggiata. Não utilize se a embalagem estiver danificada. 如果包装损坏，请勿使用。 梱包が破損している場合は使用しないでください。 Не использовать при поврежденной упаковке. Ambalaj zarar görmüşse kullanmayın.
<b>MR</b>	MR Safe MR-sicher MR Safe Compatible con la resonancia magnética Idoneo per RM RM Segura MR 安全 MR対応 Безопасно для МРТ MR Güvenli
<b>MD</b>	Medical device   Medizinprodukt   Dispositif médical Aparato médico   Dispositivo médico   Dispositivo médico   医疗设备   医療機器   Медицинское оборудование   Tibbi cihaz
CAUTION: (USA only) Federal law restricts this device to sale by or on the order of a physician.	

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**EN** Flow sensor  
Instructions for use

**WARNING**

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.
- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.

**CAUTION**

When using active humidification, prevent water accumulation in the flow sensor by ensuring that the flow sensor is positioned at a  $\geq 45^\circ$  angle relative to the floor. Excess water can affect the flow sensor measurements and lead to inaccurate volume delivery, potentially resulting in hypoventilation.

**NOTICE**

- A neonatal flow sensor is required with breathing circuits used on all ventilators for neonatal ventilation. With the HAMILTON-C1/T1/MR1 ventilators, when using the nCPAP or nCPAP-PC modes, remove the flow sensor and use the pressure line with the breathing circuit.
- To prevent inaccurate flow sensor readings, make sure the flow sensor is correctly installed: The flow sensor tubes must not be kinked, and must be secured with the included clamp.
- Replace the flow sensor in accordance with hospital infection control procedures, or depending on patient's secretions and nebulization of medication, if required.
- Visually inspect the flow sensor body, tubing, and internal flap. Discard the flow sensor if there is any sign of damage.
- The flow sensor must be calibrated prior to patient use.

**Intended use**  
The Hamilton Medical neonatal flow sensor is a device to measure patient air flow to be used with the neonatal patient group of Hamilton Medical ventilators.

**Connecting the flow sensor to a breathing circuit**  
Insert the flow sensor between the Y-piece of the patient circuit and the patient connection. The end of the flow sensor which should be towards the patient is marked with „Patient“. It is recommended that a short section of flexible tubing be used between the flow sensor and the actual patient connection.

**Flow sensor calibration**  
The flow sensor must be calibrated prior to patient use. To calibrate, follow the instructions in the appropriate ventilator operator's manual. Connect the supplied calibration connector for the calibration. After successful completion of the procedure, the flow sensor is ready for patient use.

**Delivery**  
The flow sensor is delivered clean, ready for clinical use.

**Disposal**  
A used flow sensor must be handled as contaminated. Follow all local, state and federal regulations with respect to environmental protection when disposing of used flow sensors.

**Environmental requirements and storage conditions**

Operating temperature	-15°C to 40°C (5°F to 104°F)
Storage temperature	-20°C to 60°C (-4°F to 140°F)
Storage and operating humidity	Relative humidity 5% to 95%, noncondensing

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