Vitamin D CONTROL

REF**2130152 - Tri-level Kit: 6 x 2 mL****IVD**

EN EXPLANATION OF SYMBOLS
CS VÝZNAM SYMBOLŮ
DA SYMBOLFORKLARING
DE ERKLÄRUNG DER SYMBOLE
EL ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES SIGNIFICADO DE LOS SÍMBOLOS
ET SÜMBOLITE SELGITUS
FR EXPLICATION DES SYMBOLES
HU JELMAGYARÁZAT
IT SPIEGAZIONE DEI SIMBOLI
JA 記号の説明
LV SIMBOLU SKAIDROJUMS
LT ŽYMĖJIMŲ PAAIŠKINIMAS
NL VERKLARING DER SYMBOLEN
PL OBJAŚNIENIE SYMBOLI
PT EXPLICAÇÃO DOS SÍMBOLOS
RO EXPLICAREA SIMBOLURILOR
SK VÝZNAM SYMBOLOV
SL RAZLAGA SIMBOLOV
SV SYMBOLFÖRKLARING
ZH 符号注解

LOT

Batch code/Číslo šarže/
Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Sarzszzám/
Codice del lotto/バッチ コード/
Partijas kods/Partijos kodas/Lot
nummer/Kod partii/Código do
lote/Lot nr./Číslo šarže/Številka
serije/Satsnummer/批号

Date of manufacture/Datum
výroby/Produktionsdato/
Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha
de fabricación/Tootmiskuuäev/
Date de fabrication/Gyártás
időpontja/Data di fabbricazione/
製造日/Izgatavošanas datums/
Pagaminimo data/Productiedatum/
Data produkcji/Data de fabrico/
Data de fabricație/Datum
výroby/Datum izdelave/
Tilverkningsdatum/生产日期

Use By/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha de
caducidad/Kasutada enne/Utiliser
jusque/Felhasználható/Utilizzare
entro/使用期限/Izlietot līdz/
Sunaudoti iki/Houdbaar tot/Użyć
przed/Prazo de validade/Data
expirării/Použití do/Uporabno
do/Använd före/此日期前使用

REF

Catalogue number/Katalogové číslo/
Katalognummer/Bestellnummer/
Αριθμός καταλόγου/Número de
catálogo/Katalogi number/Référence
du catalogue/Katalógusszám/Numero
di catalogo/カタログ番号/Kataloga
numurs/Katalogo numeris/Catalogus
nummer/Numer katalogowy/Referência
de catálogo/Număr de catalog/
Katalogové číslo/Kataloška številka/
Katalognummer/目録番号

Manufacturer/Výrobce/Producent/
Hersteller/Κατασκευαστής/
Fabricante/Тоотъја/Fabricant/Gyártó/
Fabbicante/製造業者/Izgatavotājs/
Gamintojas/Fabrikant/Producent/
Fabricante/Productor/Výrobca/
Izdevalovca/Tillverkare/生産商

IVD

In Vitro Diagnostic Medical Device/
In Vitro diagnostický zdravotnícký
prostředek/Medicinský udstyr
til in vitro-diagnostik/In-vitro-
Diagnostikum/In Vitro Διαγνωστικό
Ιατροτεχνολογικό προϊόν/Product
sanitario para diagnóstico in vitro/
In vitro diagnostiline meditsiiniseade/
Dispositif médical de diagnostic in
vitro/In vitro diagnostikum/Dispositivo
medico-diagnostico in vitro/体外診断
医療機器/Medicinas ierīce, ko lieto
in vitro diagnostikā/In vitro diagnostikos
medicinos prietaisai/Medisch
hulpmiddel voor in-vitro diagnostiek/
Wyrób do diagnostyki In Vitro/
Dispositivo médico para diagnóstico in
vitro/Dispositiv medical de diagnostic
in-vitro/Zdravotnícka pomůcka in
vitro/In vitro diagnostični pripomoček/
Medicintekniska produkter för in vitro
diagnostik/体外诊断医疗器械

Temperature limitation/
Teplotní rozmezí od do/
Temperaturbegrensning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limite de temperatura/
Temperatuuri piirangud/
Limites de température/
Hőmérséklettartomány/
Limiti di temperatura/温度制限/
Temperatūras ierobežojumi/
Temperatūras apribojimas/
Temperatūrlīmiets/
Przestrzegać zakresu temperatury/
Limites de temperatura/
Limite de temperatură/
Teplotné rozmedie od do/
Omejitve temperature/
Temperaturbegränsning/
温度限制

EC REP

Authorized representative/Pověřený
zástupce/Autoriseret repræsentant/
Bevollmächtigter Vertreter/
Εξουσιοδοτημένος αντιπρόσωπος/
Representante autorizado/
Volitadud esindaja/Märkakeskedö/
Rappresentante autorizzato/
委任された代表/Pilnvarots pārstāvis/
Jgalotasis atstovas/Geautoriseerd
vertegenwoordiger/Autorizovaný
zástupce/Representante autorizado/
Reprezentant autorizat/Autorizovaný
zástupca/Poolbāšēni pārstāvniks/
Auktoriserad representant/授权代表

Consult Instructions for Use/Viz
návod k použití/Se brugsanvisning/
Gebrauchsanweisung beachten/
Συμβουλευτείτε τις οδηγίες χρήσης/
Consulte las instrucciones de uso/
Lugege kasutusjuhendit/Consulte les
instructions d'utilisation/ Nézze meg
a Használati utasítást/Consultare
le istruzioni per l'uso/使用説明書
を参照/lepazīstieties ar lietošanas
instrukciju/Žr. naudojimo instrukciją/
Raadpleeg de gebruiksaanwijzing/
Sprawdź w instrukcji obsługi/Consulte
as instruções de utilização/Consultati
instrucțiunile de utilizare/Vid návod
na použití/Oglejte si navodila za
uporabo/Se brugsanvisningen/
请参考“使用说明”

Biological risks/Biologicky nebezpečné/
Biologisk fare/Biogefährdung/Βιολογικοί
κίνδυνοι/Riesgo biológico/Biologigiline
oht/Risques biologiques/Biológiai
kockázat/Rischio biologico/
生物学的リスク/Biologiskie riski/
Biologinis pavojus/Biologisch risico/
Riziko bioloģicne/Risco biológico/
Risc biologic/Bioloska tveganja/Bioloska
tveganja/Biologisk risk/生物风险

CONT

Contents of kit/Obsah sady/Kittets
inhold/Inhalt des Kits/Περιεχόμενα
του κιτ/Contenido del kit/Komplekti
sisu/Contenu du kit/A készlet
tartalma/Contenuto del kit/キットの
内容/Komplekta saturs/Komplekto
turinys/Inhoud van de set/Zawartość
zestawu/ Conteúdo do kit/Conținutul
kitului/Obsah súpravy/Vsebina
kompleta/Satsinnehåll/试剂盒内容

Reconstitute with/Rozfiedte pomoci/
Rekonstituere med/Rekonstituieren mit/
Ανασύσταση με/Reconstituir con/Lahusti/
Reconstituer avec/Feloldásához/Ricostituire
con/次により調製する:/Atskaidīt ar/Ištirpinti
naudojanti/Reconstituție met/Odtworzyć za
pomocą/Reconstituir com/Reconstituire cu/
Rozried'te pomocoou/Rekonstituiraite z/s/
Rekonstituera med/复溶剂

EN

CS

DA

DE

EL

ES

ET

FR

HU

IT

JA

LT

LV

NL

PL

PT

RO

SK

SL

SV

ZH

INTENDED USE

For In Vitro Diagnostic Use Only.

Fujirebio Diagnostics Vitamin D Control is intended for use as an assayed control serum to monitor the precision of laboratory testing procedures for the analysis of Vitamin D.

SUMMARY AND PRINCIPLE

This quality control product can be used as an objective judgment of the laboratory's procedures and personnel techniques. It is a valuable tool to assess good laboratory practices. Three levels of control are available to compare observations with expected ranges therefore assuring consistent performance of the testing system within the clinical range.

REAGENTS

This product is prepared from human serum, protein (bovine), purified biochemical materials, and chemicals. Contains Proclin 300® / Gentamicin as preservatives. Level 1 of this control contains both 25(OH) Vitamin D₃ only; Level 2 and Level 3 contain both 25(OH) Vitamin D₂ and 25(OH) Vitamin D₃. Each level of control contains clinically relevant concentrations of Total Vitamin D (1). The control is provided in lyophilized form.

STORAGE AND STABILITY

Unopened: This product is stable until the expiration date when stored un-opened at 2-8°C.

Reconstituted: Once the control is reconstituted, the product is stable for 14 days when stored tightly capped at 2-8°C. The product is stable for 60 days when stored at <-10°C. Controls may be frozen and thawed repeatedly for up to 9 cycles. Do not store at ambient temperature.

RECONSTITUTION

Using a calibrated pipette, reconstitute each vial with 2.0 mL of distilled or deionized water. Replace the stopper and lightly tap the vial 5-6 times to loosen the pellet and allow the controls to stand at room temperature, 15-30°C for at least 30 minutes. Invert three times and swirl before use to ensure homogeneity.

PROCEDURE

Follow all instructions provided with the assay or instrument of interest while running the control. The control is to be run as you would a patient sample. Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. Quality control materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements. In the event of damage to packaging, contact your local distributor or Fujirebio Diagnostics at the numbers listed on the last page of this package insert.

WARNING

Human source material. Treat as potentially infectious. Each human donor unit used to manufacture this control was tested by FDA accepted methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV), antibody to HIV-1/HIV-2 and antibody to Treponema Pallidum (Syphilis). This product may also contain other human pathogens for which there are no approved tests. All human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

Never mouth pipette. To avoid cross-contaminations, ensure that the stoppers and caps are not switched during use.

**CONTAINS METHYLISOTHIAZOLONES****Hazard statement:**

H317: May cause an allergic skin reaction.

Precautionary statement:**Prevention:**

P261: Avoid breathing mist or vapor.

P272: Contaminated work clothing must not be allowed out of the workplace.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P302+P352: If on skin: Wash with plenty of water.

P333+P313: If skin irritation or rash occurs: Get medical advice/attention.

P363: Wash contaminated clothing before reuse.

Disposal:

P501: Dispose of contents/container in accordance with local/regional/national/international regulations.

LIMITATIONS

1. Do not use this product beyond the expiration date.
2. Discard the vial if there are indications of microbial contamination or cloudiness in the reconstituted control.

3. This product is not to be used as a standard.

ASSIGNMENT OF VALUES

The assigned ranges shown have been established for the assays listed in the Assigned Values Sheet only. The ranges listed resulted from replicate analyses and are specific for each lot of product. The tests listed were performed by the reagent manufacturer and/or independent laboratories using manufacturer supported reagents (available at the time of assay) and a sampling of this lot of control. Your laboratory means should fall within the printed acceptable range; however, laboratory means may vary from the listed values during the life of this control.

Variations between laboratories and variations over time may be caused by differences in laboratory technique, instrumentation and reagents, or by modifications made to manufacturer's test methods. It is recommended that each laboratory establish its own means and acceptable ranges and use the printed ranges as guidance only. To convert ng/mL to nmol/L use a conversion factor of 2.496 (1).

SPECIFIC PERFORMANCE CHARACTERISTICS

This is a lyophilized product manufactured under strict quality standards. This control must be stored and handled as directed to obtain expected values.

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Fujirebio Diagnostics, Inc.
Fujirebio Diagnostics, AB.