**DIAGNOSTIC KIT FOR HIV
(1+2) ANTIBODY
(COLLOIDAL GOLD)
V2**

Instructions for use

REF R-401-1-C-CE (1 Test)

KHB® Shanghai Kehua Bio-engineering Co., Ltd.

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DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2

REF R-401-1-C-CE (1 Test)

PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE

Intended Use

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB is an *in vitro*, visually read test for the qualitative determination of antibodies against HIV-1 and HIV-2 in human serum, plasma, venous and capillary whole blood. This product is intended for detection of antibodies against HIV from individuals to aid in the diagnosis of HIV infection.

It is intended for professional use in laboratory and at point-of-care.

Summary and Explanation of the Test

AIDS (Acquired Immunodeficiency Syndrome) is an immunosuppressive disease characterized by depletion of the helper T-cell, which leaves the infected individual susceptible to opportunistic infections and some malignancies. Viruses that have been detected to cause AIDS are HIV-1, first isolated in 1983⁽¹⁾ and HIV-2, a distinctive but related type isolated in 1985. Presence of the virus elicits the production of the antibodies specific to either HIV-1 or HIV-2, and this test is designed to detect the antibodies to indicate the presence of HIV infection.

Principle of the Procedure

The Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 from KHB adopts the solid phase colloidal gold immunochromatographic technology for the qualitative detection of antibodies against HIV-1/HIV-2. The gold HIV-1 gp160⁽²⁾ conjugate and the gold HIV-2 gp36⁽³⁾ conjugate are coated to the conjugate pad in advance. The test line (recombinant gp36 and gp41 antigens) and the control line (monoclonal antibody against gp160) are pre-coated on the surface of Nitrocellulose (NC) membrane. When the specimen is added to the sample pad, it migrates through the conjugate pad, reconstitutes and mixes with the colloidal gold-antigen complex and moves on till the end of the strip. A reddish-purple test line will be visible in the strip if there are enough antibodies against HIV-1/HIV-2 in the specimen. If the specific antibodies are absent, or present at a very low level, no test line appears^(4, 5).

Materials Provided

• HIV (1+2) antibody test cassette with desiccant	1 cassette	Cassette			
• Sample diluent	0.5 ml × 1 bottle	Sample Diluent			
• Instructions for use	1 piece				
Accessories		Manufacturer	EC-Representative	CE mark	Classification
• Safety lancet	1 piece	Ningbo Medsun Medical Co., Ltd. No.55, Jinxi Road, Zhenhai, 315221, Ningbo, P.R.China	Shanghai International Holding Corp. GmbH (Europe) Eiffestraße 80 20537 Hamburg, Germany	CE ₀₁₂₃	MDD93/42/EEC, AnnexV, Class IIa
• Alcohol pad (70% Isopropyl)	1 piece	Lights Medical Manufacture Co., Ltd. No.19, Quanda Road, Wuqing	Welkang Ltd., t/a Welkang Tech Consulting Suite B, 29 Harley Street,	CE	MDD93/42/EEC, Annex VII, Class I

KHB [®] DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2			For <i>in vitro</i> diagnostic use EN		
alcohol)		Development Area, Tianjin, 301700, P.R.China	London W1G 9QR, England, United Kingdom		
		Zhejiang Sorfa Life Science Research Co., Ltd.	Sungo Certification Company Limited.		
• Disposable transfer pipette	1 piece	No.148 Longshan Road, Zhongguan Town, Deqing county, Huzhou City, 313220, Zhejiang, P.R.China	RM101, Maple House,118 High Street, Purley, London, England, United Kingdom	CE	98/79/EC, Annex III, Other Devices

Materials Required But Not Provided

- Specimen transfer device (graduated)
- Tips (Pipette)
- Protective gloves
- Specimen and test waste container
- Timer
- Collection tube (with or without anticoagulant, see Specimen Collection and Storage Requirements)

Warnings and Precautions

The test is for *in vitro* diagnostic use only. The test kit is for single use only. The precautions are included, but not limited to the following:

[Safety Precautions]

- Wear gloves during the entire testing process.
- Do not use the accessories if seal is broken or package is damaged.
- All the accessories are for single use. Do not reuse the accessories.
- Do not pipette by mouth.
- Do not eat or smoke while handling specimens.
- Clean and disinfect all the areas that may be contaminated by spills of specimens or reagents with appropriate disinfectant.
- Decontaminate and dispose of all specimens, reagents, accessories and other potentially contaminated materials as if they were infectious wastes, in a biohazard container.
- The gold HIV-1 gp160 conjugate is coated to the conjugate pad in advance. The gp160 antigen is prepared using a mammalian cell expression system. There is no analytical method that can assure that infectious agents in the reagents are completely absent. Therefore, handle reagents and specimens with extreme caution as if capable of transmitting infectious diseases.
- The sample diluent has not to be labelled due to the calculation procedure of Regulation 1272/2008/EC.

[Handling Precautions]

- The accessories are only used for specimen collection, do not store the specimen in the pipette.
- The assay is validated for use with 40µL of specimen, obtained when the specimen reaches the second graduation of the pipette. Use of other volumes may not produce the correct result.
- Do not mix and interchange different specimens.
- The sample diluent is specific for each lot of reagent. Do not use diluents from different lots.
- Do not perform the test under environment which leads to rapid evaporation (e.g. high temperature (> 45 °C) and low humidity, fast air flow)
- Use of hemolytic specimens, rheumatoid factors-contained specimens and lipidemic, icteric specimens may cause invalid results or false results.
- Pooled specimens or specimens other than specified (e.g. saliva, urine) are strictly forbidden to be tested.
- Ensure that the specified volume of specimen is added by either precision pipette or specimen transfer pipette provided within the test kit.
- Intake of air bubbles in the pipette during the collection of fingerstick whole blood may cause false negative result.

[Warning]

- For *in vitro* diagnosis of the infection of HIV only, this assay will not provide a diagnosis for any other disease.
- The cassette, individually enclosed in a foil pouch, is only for "single use".

- Do not use the test cassette if the pouch is damaged or the seal is broken. Contact the supplier when it happens as this is an important indicator for quality and performance of the IVD.
- The sample diluent in this kit contains Sodium Azide(<0.1%) as the preservative, be cautious while using with the kit. In case of contact with sample diluent to eyes and/or skin, flush affected area with copious amount of water. Do not dispose the sample diluent directly in drains.
- Care should be taken to avoid the contact between the sample area of the cassette and the diluent bottle to avoid the contamination of the diluent.

Storage Requirements

- The test cassettes and the sample diluent must be stored at 4-30°C until expiration date.
- Do not use expired devices.
- Test cassette should be used immediately upon opening the pouch.
- Sample diluent should be used within 3 months of opening and be stored at 4-30°C.

Specimen types

Human serum, plasma, and whole blood specimens should be collected through venipuncture (or whole blood through fingerstick).

Plasma specimens and venous whole blood specimens can be prepared with EDTA, heparin or sodium citrate as anticoagulant.

Specimen Collection and Storage Requirements

[Whole Blood]

- Collection by venipuncture

Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture.

Venous whole blood specimens can be stored up to 8 hours at 18-24 °C or up to 7 days at 2-8 °C if not tested at the same day of specimen collection, and should not be frozen.

Whole blood specimens should be brought to room temperature(18-28 °C) prior to use. Haemolytic or turbid whole blood specimens should not be used.

- Collection by fingerstick⁽⁶⁾


Figure a

Figure b

Figure c

Figure d

1. Choose the ring finger of the patient.
2. Open the alcohol pad and clean the fingertip (Figure a). Allow the fingertip dry.
3. Remove the cap of the safety lancet. (Figure b)
4. Place the lancet on the fingertip, press the lancet against the finger and puncture the skin. (Figure c)
5. Discard the lancet into a bio-hazard waste container immediately.
6. Squeeze the fingertip to form a bulb and use the disposable transfer pipette to collect the drop of the blood. (Figure d)
7. Allow the blood to reach the second graduation and the specimen volume will be 40 µL. (Figure e)
Avoid formation of bubbles during the collection of specimen. Disposable transfer pipette is used for fingerstick only.
8. Transfer the specimen to the sample well immediately after collection. Avoid formation of bubbles.
9. Discard the disposable transfer pipette into a bio-hazard waste container immediately.


Figure e

[Plasma or Serum]

- Plasma

Draw the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture and lightly mix, then centrifuge blood to get plasma supernatant.

- Serum

Draw the whole blood into the collection tube by venipuncture, leave to settle for 2 hours for blood coagulation and then centrifuge blood to get serum supernatant.

Serum and plasma specimens may be stored at 2-8 °C for up to 7 days from time of draw, store at -18 °C or below for long time storage. Several freeze-thaw cycles should be avoided.

Plasma or Serum specimens should be brought to room temperature prior to use.


Plasma or serum specimens with visible bio-contamination should not be used.

Test Procedure


Equilibrate all specimens and the devices to room temperature (18-28 °C) before testing.

Sample diluent from different batches cannot be mixed.

1. Take out a test cassette from the foil pouch, and place it on a horizontal surface.
2. Use the specimen of either serum/plasma or whole blood: Add 40µL of specimen to the sample area firstly (Figure 1), then slowly add 1 drop of sample diluent vertically to the same area (Figure 2). Please avoid direct contact of the diluent bottle with the sample area.
3. Incubate the cassette at room temperature (18-28 °C) and read the result after 15 minutes but not more than 25 minutes. Reading the test before 15 minutes or after 25 minutes may cause false result.


Interpretation of Results


Reactive Result/Positive Result (Figure 3)

Clear visible reddish-purple bands appear both at the control line (C-line) and the test line (T-line) of the cassette. The test result is valid even the band of the control line (C-line) appears much lighter or darker than the band of test line (T-line).


Non-reactive Result/Negative Result (Figure 4)

A reddish-purple band appears only at the control line (C-line) of the cassette (Figure 4).

Invalid Result 1 (Figure 5)

Reddish-purple band appears at neither the control line nor the test line of the cassette (Figure 5). Repeat the test of the specimen. If the result remains the same, then the result is invalid and contact KHB supplier.

**Invalid Result 2 (Figure 6)**

A reddish-purple band appears only at the test line (T-line) of the cassette (Figure 6). Repeat the test and contact KHB supplier if the C-line remains invisible.

Limitations of the Procedure

- The kit is designed to detect antibodies against HIV-1 and HIV-2 in human serum, plasma, and whole blood. Specimens other than the specified ones may not supply accurate results.
 - The intensity of the T-line does not necessarily correlate to the titer of antibody in the specimen.
 - The presence of the control line (C-line) only means the flow of the conjugate.
 - The sensitivity on fingerstick whole blood specimen may be lower than the sensitivity on serum or plasma specimen.
 - When specimens contain high concentration of antibody against HIV-1 are tested on the device, the control line could be absent. In this case, please perform further analysis according to the instructions of Figure 6 in Interpretation of Results.
 - As this product is intended to detect antibodies against HIV (1+2) from individuals, clinical diagnosis on AIDS should not be made only based on the test results.
 - A non-reactive result should not exclude the possibility of infection caused by HIV-1 or HIV-2. A false negative result can occur in the following circumstances:
 - Negative result may not have detected recently acquired HIV infection.
 - Low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test.
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration, in exceptional cases this may lead to observation of false negative results.
 - Specimens are not properly stored.
 - False negative result due to very high concentration of antibody(ies) to HIV (which is known as 'hook effect').
- For these reasons, care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.
- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
 - The product is not validated on specimens from infants, children.
 - The product is not specifically designed for detection of antibody to HIV-1 group O.
 - No interference of anti-retroviral treatment is observed on Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2.

Performance Characteristics

- Specificity

A total of 1900 specimens including serum/plasma and whole blood were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 1)

Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Specificity
		Non-reactive	Reactive	
Blood donor	1000	1000	0	100%
Whole blood donation	500	500	0	100%
Hospitalized patient specimens	200	200	0	100%
Pregnant Women	200	200	0	100%
Total	1900	1900	0	100%

A total of 108 potentially cross-reacting specimens were tested, the results are presented in Table 2.

Table 2

Potentially cross-reacting specimens	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2	
		Non-reactive	Reactive
Rheumatoid factor	10	10	0
Cytomegalovirus	5	5	0
Epstein Barr Virus	13	11	2
Malaria	5	5	0
Syphilis	5	5	0
Herpes	5	5	0
HBsAg	15	15	0
Anti-HCV	15	15	0
Anti-HBc	15	15	0
Anti-HTLV1/2	15	15	0
Anti-HEV	5	5	0
Total	108	106	2

- Sensitivity

A total of 500 HIV-1 and HIV-2 antibody positive specimens were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 3)

Table 3

Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Sensitivity
		Reactive	Non-reactive	
HIV-1Ab positive specimens	360	360	0	100%
non-B subtypes HIV-1 Ab positive specimens	40	40	0	100%
HIV-2 Ab positive specimens	100	100	0	100%
Total	500	500	0	100%

A total of 40 HIV-1 Non-B Subtypes positive specimens (in Table 3) were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2, the results are presented in Table 4.

Table 4

HIV Subtype	No. of specimens tested	Reactive on KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2
A	1	1
A1	1	1
A2	1	1
C	3	3
CRF01_AE	3	3
CRF02_AG	3	3
CRF06_cpx	3	3
CRF36_cpx	3	3
D	3	3
F1	2	2
F2	1	1
G	3	3
H	3	3
J	3	3
K	3	3

KHB [®] DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2			For <i>in vitro</i> diagnostic use EN
Group O	3		3
Recombinant Form (recombination of subtype B, D, CRF19 and CRF21_A2D)	1		1
Total	40		40

A total of 100 paired HIV-1 antibody positive whole blood and plasma specimens were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 5)

Table 5				
Population	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Sensitivity
		Reactive	Non-reactive	
HIV-1Ab positive whole blood	100	100	0	100%
HIV-1 Ab positive plasma*	100	100	0	100%

*The 100 HIV-1 Ab positive plasma are included in the HIV-1 Ab positive plasma in Table 3.

Seroconversion panels

The sensitivity was evaluated on 30 commercially available seroconversion panels. The Diagnostic kit for HIV (1+2) Antibody (colloidal gold) V2 detected 115 specimen out of 246 seroconversion panel members and the two CE marked rapid HIV Ab tests detected 112 specimen in total, which is 3 samples more than the CE marked rapid HIV Ab tests among all 30 seroconversion panels.

Among the 18 seroconversion panels from SeraCare Life Sciences tested with the Diagnostic kit for HIV (1+2) Antibody (colloidal gold) V2 by the Institute of Tropical Medicine, there are 27 HIV seroconversion samples and 53 HIV early seroconversion samples. All the 27 seroconversion samples were detected by Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. Among the 53 HIV early seroconversion samples, Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 detected 6 samples and the CE marked rapid HIV Ab tests detected 5 samples.

A total of 30 commercially available seroconversion panels were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2. (Table 6)

Table 6					
	Panel	Panel manufacturer	KHB [®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2		Difference
			CE marked rapid HIV Ab test A	CE marked rapid HIV Ab test B	
Tested in Institute of Tropical Medicine (ITM)	PRB925	SeraCare Life Sciences	2/6	2/6	0
	PRB930		2/4	2/4	0
	PRB934		2/3	2/3	0
	PRB945		2/6	2/6	0
	PRB947		3/4	3/4	0
	PRB950		1/4	0/4	1
	PRB951		1/6	1/6	0
	PRB952		2/6	1/6	1
	PRB953		1/4	0/4	1
	PRB955		2/5	2/5	0
	PRB958		2/6	2/6	0
	PRB959		4/7	4/7	0
	PRB967		2/6	3/6	-1
	PRB968		3/10	3/10	0

KHB® DIAGNOSTIC KIT FOR HIV (1+2) ANTIBODY (COLLOIDAL GOLD) V2				For <i>in vitro</i> diagnostic use EN		
Tested in KHB	PRB969		3/10	3/10	0	
	PRB970		1/4	1/4	0	
	PRB973		1/4	1/4	0	
	PRB977		2/4	1/4	1	
	SCP-HIV-002	Biomex	10/20	9/20	1	
	SCP-HIV-003		2/7	2/7	0	
	SCP-HIV-004		3/7	3/7	0	
	SCP-HIV-005		19/25	19/25	0	
	SCP-HIV-006		15/17	15/17	0	
	SCP-HIV-007		3/9	4/9	-1	
	HIV9075	ZeptoMetrix	2/5	2/5	0	
	HIV9079		15/25	14/25	1	
	HIV9084		1/4	1/4	0	
	HIV9096		2/6	2/6	0	
	HIV12007		4/9	4/9	0	
	HIV12008		3/13	4/13	-1	
Total	30	—	115/246	107/225	5/21	3

Anti-retroviral treatment samples

30 anti HIV positive samples from individuals which are under anti-retroviral treatment (ART) were tested. The 30 samples are collected from individuals of different stages of HIV infection and different combinations of anti-retroviral drugs.

The combination of anti-retroviral drugs listed in table 7 are tested and no interference of anti-retroviral treatment is observed.

Table 7

Combinations of anti-retroviral drugs	No. of specimens tested	KHB [®] Diagnostic Kit for HIV (1+2)Antibody (Colloidal Gold) V2	
		Reactive	Non-reactive
3TC + TDF + EFV	25	25	0
3TC + AZT + EFV	2	2	0
3TC + TDF+ LPV/r	2	2	0
AZT + 3TC + NVP	1	1	0

3TC - Lamivudine, TDF - Tenofovir disoproxil, EFV - Efavirenz, LPV/r - Lopinavir/ Ritonavir, AZT - Zidovudine, NVP – Nevirapine

Fingerstick whole blood


30 paired anti HIV positive fingerstick whole blood and EDTA plasma were tested by KHB[®] Diagnostic Kit for HIV (1+2) Antibody (Colloidal Gold) V2 and all the results are reactive. However, 7 of the 30 fingerstick whole blood specimens give a weaker intensity of T-line than paired EDTA plasma. In the test of anti HIV positive fingerstick whole blood, bubbles were observed in 2 samples during the sample collection and sample addition to the sample well of the cassette, which led to non-reactive results. Fingerstick whole blood samples were collected again from the 2 patients and tested on both Diagnostic kit for HIV (1+2) Antibody (colloidal gold) V2 and reference assay. The test results are both reactive.

References

1. Gallo, R. C., S. Z. Salahuddin, M. Popovic, G. M. Shearer, M. Kaplan, B. F. Haynes, T. J. Palker, R. Redfield, J. Oleske, B. Safai, and et al. 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. Science **224**:500-503

2. **Viscidi, R., E. Ellerbeck, L. Garrison, K. Midthun, M. L. Clements, B. Clayman, B. Fernie, and G. Smith.** 1990. Characterization of serum antibody responses to recombinant HIV-1 gp160 vaccine by enzyme immunoassay. NIAID AIDS Vaccine Clinical Trials Network. AIDS research and human retroviruses **6**:1251-1256.
3. **Chen, Y. H., Y. Xiao, W. Wu, Q. Wang, G. Luo, and M. P. Dierich.**2000. HIV-2 transmembrane protein gp36 binds to the putative cellular receptor proteins P45 and P62. Immunobiology **201**:317-322..
4. **Brandt, J., L. O. Andersson, and J. Porath.** 1975. Covalent attachment of proteins to polysaccharide carriers by means of benzoquinone. Biochimica et biophysica acta **386**:196-202.
5. **Osikowicz, G., M. Beggs, P. Brookhart, D. Caplan, S. Ching, P. Eck, J. Gordon, R. Richerson, S. Sampedro, D. Stimpson, and et al.** 1990. One-step chromatographic immunoassay for qualitative determination of choriogonadotropin in urine. Clinical chemistry **36**:1586.
6. **Spielberg, F., C. M. Kabeya, R. W. Ryder, N. K. Kifuani, J. Harris, T. R. Bender, W. L. Heyward, and T. C. Quinn.** 1989. Field testing and comparative evaluation of rapid, visually read screening assays for antibody to human immunodeficiency virus. Lancet **1**:580-584.

Key to symbols used

	Caution		Temperature limitation (4-30 °C)
	Batch code		Date of manufacture
	In vitro Diagnostic use Medical device		Consult instructions for use
	Catalogue number or order number		Do not use if package is damaged
	Use by		Sufficient for "n" use
	Do not reuse		Biological risks
	Manufacturer		Authorised representative in the European Community
	Identification of notified body		Sterilized using Ethylene Oxide
	Sterilized using irradiation		

Shanghai Kehua Bio-engineering Co., Ltd.


1189 North Qinzhou Road, 200233, Shanghai, P.R.China

Tel: +86-21-64950625

+86-8008203370

Fax: +86-21-64950975

Website: <http://www.skhb.com>

DIAneering® Diagnostics Engineering & Research GmbH

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
Références

1. Gallo, R. C., S. Z. Salahuddin, M. Popovic, G. M. Shearer, M. Kaplan, B. F. Haynes, T. J. Palker, R. Redfield, J. Oleske, B. Safai, and et al. 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* **224**:500-503
2. Viscidi, R., E. Ellerbeck, L. Garrison, K. Midthun, M. L. Clements, B. Clayman, B. Fernie, and G. Smith. 1990. Characterization of serum antibody responses to recombinant HIV-1 gp160 vaccine by enzyme immunoassay. NIAID AIDS Vaccine Clinical Trials Network. *AIDS research and human retroviruses* **6**:1251-1256.
3. Chen, Y. H., Y. Xiao, W. Wu, Q. Wang, G. Luo, and M. P. Dierich. 2000. HIV-2 transmembrane protein gp36 binds to the putative cellular receptor proteins P45 and P62. *Immunobiology* **201**:317-322..
4. Brandt, J., L. O. Andersson, and J. Porath. 1975. Covalent attachment of proteins to polysaccharide carriers by means of benzoquinone. *Biochimica et biophysica acta* **386**:196-202.
5. Osikowicz, G., M. Beggs, P. Brookhart, D. Caplan, S. Ching, P. Eck, J. Gordon, R. Richerson, S. Sampedro, D. Stimpson, and et al. 1990. One-step chromatographic immunoassay for qualitative determination of choriogonadotropin in urine. *Clinical chemistry* **36**:1586.
6. Spielberg, F., C. M. Kabeya, R. W. Ryder, N. K. Kifuani, J. Harris, T. R. Bender, W. L. Heyward, and T. C. Quinn. 1989. Field testing and comparative evaluation of rapid, visually read screening assays for antibody to human immunodeficiency virus. *Lancet* **1**:580-584.

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DIAnceering® Diagnostics Engineering & Research GmbH