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EN : INSTRUCTIONS FOR USE OF
REUSABLE FORCEPS, FOR VESSEL SEALING

THERMOCLAMP / THERMOCISION

1. Scope of application: THERMOCLAMP™ forceps for open surgery are designed for HF bipolar electro-surgical current for vessel sealing: the electrical energy conversion into thermal energy provides melting of tissue collagen and sealing of vascular walls under forceps pressure. THERMOCISION forceps are designed for vessel sealing, haemostasis and cutting of tissues and small vessels. Their use is strictly restricted to qualified surgeons.

2. Meaning of symbols used on labelling

	See instructions for use (IFU) / identification of operating instructions.
	Manufacturer Identification (name and address)
	Manufacturing year (near from the symbol)
	Item Reference Code
	Batch number
	Serial number
	Medical device supplied in non-sterile packaging
	Warning: imperative safety instructions.
	Device complying with the essential requirements of Directive 93/42 / EEC (0459 = Identification number of the notified body LNE / G-MED which issued the EC certificate of conformity)

3. Contraindication: Some physiological factors may affect the vessel sealing result: - The advanced patient age in association with conditions, such as cancer or coronary artery disease, can reduce the healing process effectiveness, which can weaken the fusion area. - Atherosclerosis, aneurysm, radiotherapy may affect hemostasis. It is necessary to limit the vessel sealing to non-affected areas.

- Do not try to seal tissues, or vessels with a diameter over 7 mm. Do not try to seal lymphatic vessels.


4. Checking: Because of their use, the electro-surgical accessories undergo wear by mechanical or thermal or chemical factors. Accidental deterioration can also arise. A running test and a meticulous visual inspection before each procedure must be done. THE FORCEPS PRESENTING A VISIBLE SIGN OF DETERIORATION MUST NOT BE USED.

Deterioration criteria to monitor:

- **Jaws distortion** - The jaws in closed position have to be perfectly parallel, in front and side view. If the jaws converge, they could make contact (short-circuit). If they diverge, the clamping will be insufficient. These distortions compromise the vessel sealing.

- **Insulation deterioration** - Check, under magnification if necessary, the condition of the coating after cleaning. If the coating comes off, there are risks of contamination and burns.

NEVER USE THE FORCEPS WHEN THE COATING COMES OFF OR DEGRADED.

- **Scratch** - The polished surface of the jaws avoids tissue adhesion. Scratches on the jaws internal face will cause cleaning difficulties. Dried tissues on the jaws compromises the vessel sealing on total grasped area.

- **Loosening** - Connector pins and the axis of the forceps should not be able to loosen by hand without tools. It would compromise the vessel sealing.

Functional testing: Before starting the procedure, check the vessel sealing with gauze soaked with saline: connect the forceps to the generator with its cable. Adjust the output power. Place the gauze into the forceps and activate the vessel sealing without releasing. The current crossing causes evaporation. If not, check the adjustment or replace the cable.

5. Electrical compatibilities: bipolar current to apply on forceps must be only the sealing current emitted by generators SURGILEC SEAL, OPTIMA or MCB. The rated voltage of THERMOCLAMP™ forceps is 350 V peak. Using a higher voltage can damage

their insulation and can cause tissue charring. Never activate the bipolar cut, or the bipolar coagulation on these forceps.

6. Patient safety: There is no special recommendation for the patient safety when using bipolar mode. However, if during the surgery, monopolar mode is used, safety instructions must be followed (see the user's manual of generator).

7. Connection to the generator: Connection to the generator: The THERMOCLAMP™ forceps have to be connected to the electrosurgical unit on the dedicated output for vessel sealing with a reusable and sterilizable cable (ref. V11F343). The 2 female plugs with a clamping device have to be turned as indicated on the opposite picture, to enable the device to work correctly. In the case of a bad connection on the plugs, or a break in the cable, an electrical arc can occur, causing an electrical disturbance. In this case, replace the cable. Disconnect the cable by pulling the plug and not the cord.

8. Generator settings - SURGILEC SEAL: Refer and conform to the electrosurgical unit's user manual. - **OPTIMA:** On the Generator screen, in the [Vessel sealing] menu, select the icon of the forceps model used. 3 power levels are available for each forceps (1 to 3 bars). Select the power level according to the result obtained, and the duration of fusion. - **MCB** unit: setting is automatic. It allows you to choose from 3 preset levels.

9. Precautions for use - Inappropriate storage and/or physical stress during transport could lead to a failure or dysfunction/malfunction likely to induce harm. Therefore, it is required for the user to check the proper operation of device before use and to strictly follow the mandatory verification operations as described in IFU.

- Non observance of IFU in general & in particular warnings, such as the strict application of verification and maintenance operations, or for device disposal can lead to misuse, dysfunction/malfunction likely to induce harm (Adverse Events or Serious Adverse Events) or non-effective device (intended use not reached), or hazardous situations for environment.



- Devices shall be used strictly according to IFU: by persons corresponding to the user profile, for indications and intended use as described. All other use are prohibited and contraindicated by the manufacturer because it could lead to Hazardous Situations likely to induce Harm (Adverse Events or Serious Adverse Events) or non-effective device (intended use not reached), or hazardous situations for environment.

- As any device, dysfunction and/or Failure can occur during the device lifecycle. Therefore, it is necessary to ensure proper operation of device before use, as described in IFU. Non-Observance of operations mentioned in IFU could lead to Hazardous Situations likely to induce Harm (Adverse Events or Serious Adverse Events) or non-effective device (intended use not reached), or hazardous situations for environment.

- Failure to observe the instructions for use of the ESU during operation may cause damage or serious damage. Refer to the operating instructions of the ESU, on which the active accessories are connected.


- Benefits outweigh Serious Adverse Events (SAE) occurrence. However, Serious Adverse Events (SAE) is reported in scientific and medical literature and reviewed periodically to ensure adherence to state of art. These Serious Adverse Events need to be addressed immediately by medical staff. SAE reported and known are as follow:

A) **BLEEDING RISK:** Do not try to seal a vessel with the forceps tip. Forceps has to fully grasp the vessel to activate an effective seal process.

B) Grab tissues to be sealed between the forceps jaws and grip until the clamp locks. Activate the vessel sealing current with the single pedal until the audible beep signal indicates the complete sealing of the tissue. Open the forceps and start a second vessel sealing near to the first one. **TISSUE SHOULD BE CUT ONLY BETWEEN THE 2 SEALED AREAS WITHOUT CUTTING OUTSIDE OF THE SEALED AREA.**

C) At a junction between two vessels, maintain a safe distance of some mm to preserve the uncut vessel from thermal spread and maintain the possibility to suture if needed.

D) **BURN RISKS:** DURING THE VESSEL SEALING PROCESS, AVOID CONTACT BETWEEN JAWS AND COLLATERAL TISSUES. IF IT IS POSSIBLE, LIFT THE FORCEPS OR INTERPOSE A RETRACTOR (OR A VAGINAL RETRACTOR FOR VAGINAL HYSTERECTOMY) DO NOT LET THE FORCEPS COME INTO CONTACT WITH THE PATIENT AFTER USE: THE JAWS MAY STAY REALLY HOT FOR A WHILE.

E) When jaws are immersed in blood or other liquid, remove the liquid by suction or with a compress before activating the vessel sealing.

In case of hemorrhage, locate the bleeding source and occlude or clamp the concerned vessel to stop the hemorrhage before starting vessel sealing.

F) The forceps jaws have to be frequently cleaned with a gauze soaked in saline solution, in order to maintain the best contact with the tissues to be sealed.

AUTO-CLEANING: To clean the jaws during intervention, soak them in saline solution and then remove them from the solution without opening the jaws; activate the vessel sealing for a few seconds in order to boil the liquid. Then clean the jaws with a compress.

10. Reprocessing procedure:

Preliminary instructions:

- **Service life:** THERMOCLAMP and THERMOCISION forceps are designed to withstand about one hundred reprocessing cycles. Nevertheless, the effective life is subject to checks before use (see § 4).

- **Safety:** Staff must be protected against contamination by soiled devices. The required protections include at least wearing a blouse, a pair of gloves, a mask and goggles (or a face shield).

- **Tools:** manual cleaning requires a brush (Ref LTA medical: brush Ref 520). Do not use a wire brush or abrasive utensil that may scratch the surface of the device.

- **Use clean utensils:** LTA brushes can be washed in machine at 95 ° C and support steam sterilization up to 135 ° C.

- Washing solutions and rinse water should not be reused. After use, the bins must be cleaned, thoroughly rinsed, drained and dried.

- Time between each reprocessing step: It is important to shorten the time between each step. Too long delays favor the risks of re-contamination and proliferation of microbiological germs. Too long delays for the rinsing step favor chemical attacks on the materials or the deposition of residues.

Preparation at the point of use: - Immediately after use, disconnect the removable cables. - Immerse the devices in enzymatic solution at room temperature (+ 15 ° C to + 25 ° C) for 5 minutes (LAMIDEY NOURY MEDICAL recommends PROLYSTICA 2 X CONCENTRATE 0.4%, from STERIS CORPORATION). Do not leave the immersed devices more than one hour without cleaning. Do not use saline for soaking.

Manual cleaning: in the soaking solution. - Actuate the forceps (close /open) 3 times - Shake the devices for 10 seconds to remove the air bubbles - Carefully brush each device during 2 minutes, insisting on interstices. - Replace the dipping solution with a new freshly prepared solution. Soak for 2 minutes.

Rinse: - Immerse the devices in an osmosis water bath for 20 seconds. Shake the devices to chase the air bubbles. - Renew the rinsing step twice (a total of 3 x 20 seconds of rinsing).

Visual Inspection: - search under illumination and magnification for any trace of dirt or degradation on the devices. Look for traces of corrosion that can be identified by their coloring, crevices, and a rusted surface. Soiled devices will have to go through the cleaning steps. Damaged devices must be removed.

Drying: - Dry the outer surface of the devices with a soft, lint-free cloth.

11. Sterilization - Preparation for bagging: protect all ends of electrodes or sharp forceps with small silicone tubes so as not to puncture the sterilization bag. Do not mount the electrodes on the pencils. Wrap the cables in large loops, without bending.

Steam sterilization: Perform sterilization in accordance with local regulations.

LAMIDEY NOURY MEDICAL recommend a fractionated prevacuum, then steam sterilization at 134 ° C during 13 minutes, drying 20 minutes, cooling 60 minutes. For prion contamination risk, duration of steam sterilization can be increased at 18 minutes.

12. Disposal of accessories -. Disposal of a contaminated device is likely to cause damage to persons and the environment. They must be disposed of according to a specific circuit suitable for contaminated products.

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