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
GE Healthcare

# Revolution™ EVO

Fișa de date a produsului – Rev. 5

**CONFIDENTIAL**  
**SECRET DE AFACERI**

## Cuprins

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GE Healthcare  
9900 Innovation Drive  
Wauwatosa, WI 53226  
SUA

Chalfont St. Giles  
Buckinghamshire  
Marea Britanie

[www.gehealthcare.com](http://www.gehealthcare.com)

imagination at work

SECRET DE AFACERI

CONFIDENTIAL

GE Healthcare

<http://www.gehealthcare.com/>

# EU Revolution EVO EL/EX Gen2

CONFIDENTIAL  
SECRET DE AFACERI

# PRODUCT DETAILS

## I. EU Revolution EVO EL/EX Gen2

| Qty. | Catalog Ref | Description                                                                                                    |
|------|-------------|----------------------------------------------------------------------------------------------------------------|
| 1    | B76242DA    | Revolution EVO system                                                                                          |
| 1    | B77942CB    | 64ch / 128sl Option                                                                                            |
| 1    | B77732CB    | 72KW Option                                                                                                    |
| 1    | B77772CB    | 0.35 Sec Rotation Software Option for Cardiac                                                                  |
| 1    | B78672CB    | ASiR Option                                                                                                    |
| 1    | B76472CB    | VT1700 Patient Table                                                                                           |
| 1    | B76062CB    | Long Cable Collector Kit                                                                                       |
| 1    | B79792CB    | Operator console monitors                                                                                      |
| 1    | B7599ZZ     | Euro misc keyboard kit                                                                                         |
| 1    | B76542CB    | Revolution EVO - European Kit                                                                                  |
| 1    | B76542CB    | Revolution EVO - European Kit                                                                                  |
| 1    | B77922FJ    | Revo EVO Russia Addendum                                                                                       |
| 1    | B75502RE    | Advanced Cardiac Package                                                                                       |
| 1    | B75542RE    | CT Angiography Package                                                                                         |
| 1    | B75682RE    | Head Perfusion Package                                                                                         |
| 1    | B75572RE    | Workflow Package                                                                                               |
| 1    | B78552CA    | Optima Operator Console Desk (Wide Design)                                                                     |
| 1    | B77322CA    | CT chair (no armrest)                                                                                          |
| 1    | B77552DA    | 2 Phase Uninterruptible Power Supply                                                                           |
| 1    | B75352CA    | Table Convenience kit                                                                                          |
| 1    | B78792CB    | UW protocol for Revolution EVO                                                                                 |
| 1    | R23053AC    | Standard service tools L3 - warranty period                                                                    |
| 3    | W90862CT    | Prepaid Tube Watch Option Special POS offering                                                                 |
| 1    | M81521KA    | AW VolumeShare 7 with 32GB of RAM                                                                              |
| 1    | M80281AA    | AW VolumeShare 7 Monitors                                                                                      |
| 1    | M81501PG    | AW HP Workstation Power Cords Kit (EMEA)                                                                       |
| 1    | M80501DV    | Field Engineer Letter                                                                                          |
| 1    | M81521TS    | 3D Suite                                                                                                       |
| 1    | M81521VQ    | Volume Illumination                                                                                            |
| 1    | M81521ED    | Integrated Registration                                                                                        |
| 1    | M81521BP    | OncoQuant                                                                                                      |
| 1    | B79921TA    | CardIQ Function Xpress                                                                                         |
| 1    | B79821RE    | CardIQ Xpress 2.0 Reveal                                                                                       |
| 1    | B79821PE    | CT Comprehensive Neuro Package                                                                                 |
| 1    | B79971JH    | SmartScore 4.0                                                                                                 |
| 1    | B77121BK    | VessellIQ Xpress & AutoBone Xpress                                                                             |
| 1    | M85101BD    | DoseWatch Starter Kit-1 Connection(CT/CVIR)                                                                    |
| 1    | E8007TB     | Cardiac Cable Kit for IPC, Ivy 3150-B, Ivy 7800                                                                |
| 1    | E8007TM     | Ivy Monitor Roll Stand - for 7000 series                                                                       |
| 1    | E80171TD    | Ivy 7800 cardiac monitor                                                                                       |
| 1    | E8004GK     | Fuyo Strap Collector - Six straps, Narrow 6cm, Medium 15cm and Wide 37cm - 54 cm and 106 cm lengths each width |
| 1    | E8004TG     | Coronal headholder                                                                                             |
| 1    | B75342CB    | Axial Headholder, Arm Support                                                                                  |
| 1    | E8016AZ     | Cushion sealed in slicker for CT with GT1700 Table - 1700 Systems (2-pc Set)                                   |
| 1    | E45021BB    | CT MDP CE 125A 400V 50Hz 3 phases                                                                              |
| 1    | E45971NC    | ABB UPS PowerValue 11RT 2kVA, single phase, 7 min runtime                                                      |
| 1    | A82000CT    | 16 credits Welcome pack CT for new customers                                                                   |
| 1    | A11181CT    | Start Me Up 8 days CT                                                                                          |

|   |                    |                                                           |
|---|--------------------|-----------------------------------------------------------|
| 1 | A33331CT           | Make Me an Expert CT                                      |
| 1 | A82016CT           | Additional 16 days of on site application training for CT |
| 1 | W50092CT           | 24 months warranty for Revolution EVO                     |
| 1 | I019875LSI         | Remote Service (Router)                                   |
| 1 | TC_CT_QMS-R_LSI    | Lead Glass                                                |
| 1 | E10261NE-2-RO      | Dual head injector                                        |
| 1 | TC_CT_NONQMS-R_LSI | AW Workstation/Server furniture (Desk and chair)          |

CONFIDENTIAL

SECRET DE AFACERI


## LONG DESCRIPTIONS

All pictures featured are for illustration purposes only, and do not necessarily relate to products or services mentioned in the text.

### B76242DA - Revolution EVO system

Revolution EVO is the next generation Volume CT with Clarity Imaging Chain and advanced iterative reconstruction. Clarity Imaging Chain consists of Clarity Detector, DAS, Performix\*40 Plus X-ray Tube and Iterative reconstruction, and delivers high resolution imaging to meet various customer needs in real clinical situations. Clarity Imaging Chain delivers higher spatial resolution, lower noise, or less-artifact.

#### Clarity Imaging Chain

Revolution EVO's Clarity Imaging Chain consists of Clarity Detector, DAS, Performix\*40 Plus X-ray Tube and ASiR-V reconstruction, to deliver high resolution imaging.

Clarity Imaging Chain provides the following:

- For better performance Volume CT, Clarity Imaging Chain provides enhancement of spatial resolution up to 20% compared with previous GE technology.
- Designed as analog cable free between ASIC and Diode and has a capability to reduce electric noise.
- Designed for less heat generation, up to 90% compared with previous GE technology and all in one DAS / Detector. It has capability to reduce electric noise.
- Designed for less floor-noise, up to 44% compared with previous GE technology and it has capability to reduce electric noise.
- Optimized collimator with ability to reduce scatter noise.
- Performix40\* Plus X-ray tube provides less focus movement.
- A liquid bearing tube that has a capability of less-wear of Tube bearing and is enabled up to 0.35sec rotation speed option with a routine scan. Revolution EVO allows users to utilize helical pitches up to 1.531 and 0.35sec rotation speed option that meets GE's image quality specifications for lower pitch acquisitions. This high pitch and 0.35sec rotation speed enable faster scan times which may allow for shorter breath holds, and may help to avoid sedation, simultaneously (or "as well as") reducing motion artifacts from patient and organ movement. As an example, using this higher pitch, a full-body trauma scan of 1000 mm can be acquired in as little as 6 seconds.

#### Key Features:

- Silent design of Revolution EVO gantry allows significant reduction of audible noise compared with previous GE technology.
- IQ Enhance (IQE) reconstruction reduces helical Artifact Index in thin slice helical scanning. This reduction in artifacts makes it possible to scan at faster helical pitches.
- GE's protocol management is improved with the addition of a workflow improvement feature, which allows easy configuration of back to back Axial or helical scans of the same anatomy at two different X-ray energies (kVps). To further improve registration accuracy, patient immobilization may be utilized. The additionally acquired dual energy data can be post-processed on console or AW workstation using Add/Sub function to gain additional clinical information.
- Adaptive Enhance Level Adjustment (AELA) may improve visual spatial resolution while maintaining pixel noise standard deviation and artifact.

- Organ Dose Modulation provides reduction of radiation dose via X-ray tube current modulation for superficial tissues, such as breasts.
- **AutomA/SmartmA\*** modulates X-ray tube mA to account for specific patient anatomy - based upon data gathered from the scout image. The system predicts the optimal setting for the exam and adjusts mA to these settings.
- Dynamic Z-axis tracking provides automatic and continuous correction of the x-ray beam shape to block unused x-ray at the beginning and end of a helical scan to reduce unnecessary radiation.
- Xstream Display is a multi-purpose LCD display and Xstream Display can show basic patient information on the gantry monitor. The user can confirm patient information in the scan room, enhancing workflow improvement with preset positioning (Default Patient positioning) on gantry display.
- Direct MPR with Auto-Batch feature, affording automatic real-time direct reconstruction and transfer of fully corrected multi-planar images, also allows users to move from routine 2D review to prospective 3D image review of axial, sagittal, coronal, and oblique planes while enabling automated protocol-driven batch reformats to be created and networked to their desired reading location.
- Dose Check provides users with tools to help them manage CT dose in clinical practice and is based on the standard XR-25-2010 published by The Association of Electrical and Medical Imaging Equipment Manufacturers (NEMA).
- Dose Reporting: CTDIvol, DLP, Dose Efficiency displays during scan prescription and provides dose information. The CTDIvol, DLP, and Phantom size used to calculate dose is automatically saved once the user selects End Exam.
- DICOM Structured Dose Report generates a CT Dose Report, which can enable tracking of dose (CTDIvol and DLP) for the patient by the hospital radiation tracking system/RIS/HIS.

### Scan mode: Helical

- Helical Scan Speeds: Full 360° rotational scans: 0.7, 0.8, 0.9, 1.0 second
- Helical Pitch (nominal): 0.516 to 1.531
- Cardiac Pitch: 0.16 to 0.325
- Selectable kV: 80, 100, 120, 140
- Selectable mA: 10 to 560, 5mA increments
- Reconstruction Algorithms: Soft Tissue, Standard, Detail, Chest, Bone, Bone Plus, Lung, Ultra, Edge, Edge Plus

### Scan Mode: Axial & Cine

- Scan Speeds: 0.7, 0.8, 0.9, 1.0, and 2.0 second full scans (360° acquisition).
- Selectable kV: 80, 100, 120, 140
- Selectable mA: 10 to 560, 5mA increments
- Scan Plane Geometry:  $\pm 30^\circ$  gantry tilt,  $0.5^\circ$  increments
- Reconstruction Algorithms: Soft Tissue, Standard, Detail, Chest, Bone, Bone Plus, Lung, Ultra, Edge, Edge Plus

### Image Quality

- 0.28mm high resolution

### System Components:

Gantry Advanced slip ring design continuously rotates the generator, Performix\*40 Plus, Clarity detector and data acquisition system around the patient.

- Aperture: 70 cm
- Maximum SFOV: 50 cm
- Rotational Speeds: 360 degrees in 0.7, 0.8, 0.9, 1.0 seconds

**SECRET DE AFACERI**

- Tilt: +/- 30 degrees, speed 1 degree/sec

Multi-purpose LCD touch screen display with workflow features

Integrated start scan button with countdown timer to indicate when x-ray will turn on.

X-ray Tube: Performix\*40 plus liquid metal bearing tube unit offers an optimized design for exams requiring a number of scans without tube cooling.

- Performix\*40 Plus with 7.0MHU of storage and capability of 72 kw operation provides increased helical performance with greater patient throughput
- Wide range of technique (10 mA to 560 mA, in 5 ma increments) gives technologist and physician flexibility to tailor protocols to specific patient needs for optimizing patient dose.
- Heat storage capacity: 7.0MHU(Performix\*40 Plus)
- Dual Focal Spots:
  - Small Focal Spot: 0.7 (W) x 0.6 (L) Nominal Value; (IEC 60:193)
  - Large Focal Spot: 0.9 (W) x 0.9 (L) Nominal Value; (IEC 60:193)

### High Voltage Generator:

High Frequency on-board generator allows for continuous operation during scan.

72kW based system

- kV: 80, 100, 120, 140
- Max Power: 72kW
- mA: 10 to 600mA, 5mA increments

48kW based system

- kV: 80, 100, 120, 140
- Max Power (Hardware): 72kW
- mA: 10 to 400mA, 5mA increments

### Clarity Detector:

64ch based system

- 54,272 individual elements composed by 64 rows of 0.625mm thickness at isocenter. All data is acquired as thin slice at 0.625mm with the option of thicker slice from image reconstruction or processing. 98% absorption efficiency.

32ch based system

- 54,272 individual elements composed by 64 rows of 0.625mm thickness at isocenter. All data is acquired as thin slice at 1.25mm with the option of thicker slice from image reconstruction or processing. 32x 0.625mm or 32x1.25mm scan mode. 98% absorption efficiency.

### Clarity DAS (Data Acquisition System):

The Clarity DAS dramatically reduces noise and improves image performance.

- 2,460 Hz maximum sample rate.
- 861 - 1968 views per rotation.

### Revolution EVO Operator Console:

- Dual Intel Xeon 4110 CPU, 12 core, max turbo frequency 3.0 GHz.

GE Healthcare

SECRET DE AFACERI

CONFIDENTIAL

- 96 GB RAM.
- 3,000GB Disk (system, image, scan disks) stores up to 460,000 512\*2images and 3520 scan rotations at 64 slice mode or up to 1,500 scan data files, or up to 300 exams.
- Reconstruction speed with Standard reconstruction: Up to 50 frames per second.

Warranty: The published Company warranty in effect on the date of shipment shall apply. The Company reserves the right to make changes.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

Laser alignment devices contained within this product are appropriately labeled according to the requirements of the Center for Devices and Radiological Health.

Asterisk\*: Trademark of General Electric Company

### B77942CB - 64ch / 128sl Option

64ch-128sl option will provide 0.625 mm acquisition for the full 40 mm of detector coverage. Combined with the overlapped reconstruction feature, it enables 128 slices per rotation in axial scanning modes and delivers improved Z-axis visualization performance relative to non-overlapped reconstruction.

### B77732CB - mA Power Option

The mA power option upgrades the maximum allowable mA selection of the on-board high frequency generator by 40% from 400 mA max to 560 mA. More mA can be used to image large patients or at faster rotation times you can maintain the mAs prescribed.

### B77772CB - 0.35 Sec Rotation Software Option

0.35 second gantry rotation - provides the necessary software and hardware to enable a 360-degree rotation in 0.35 seconds to lower temporal resolution while scanning the heart. This software will also enable an additional rotation time of 0.37 sec in a 360-degree rotation. These additional rotation times will enhance the user's ability to acquire cardiac images of the heart with less motion. These faster rotation times are supported by the Performix™ Pro tube/generator combination which will supply the peak power required to maintain IQ with drastically reduced imaging times.

NOTE: This item does not contain any other hardware or software for cardiac acquisition or analysis.

### B78672CB - ASiR

ASiR\* (Adaptive Statistical Iterative Reconstruction) dose reduction technology.

- ASiR reconstruction technology may enable reduction in pixel noise standard deviation (a measurement of image noise). The ASiR reconstruction algorithm may allow for reduced mA in the acquisition of images, thereby reducing the dose required.
- A reconstruction technology that may enable improvement in low contrast detectability.

## B76472CB - VT1700 Patient Table

VT1700 is designed for flexible positioning with 1700mm long scannable range and 500lb (227kg) patient weight capacity.

## B76062CB - Long Cable Collector Kit

Revolution EVO and Optima CT660 Long Cable set

## B79792CB - Operator console monitors

19inch HP color LCD monitor for CT console side use. The collector includes one monitor itself EliteDisplay E190i. The BCAT includes 2pcs collectors.

## B75502RE - Advanced Cardiac Package

Advanced cardiac package is innovative cardiac options package and gives customers access to low dose cardiac acquisition and productive post-processing workflow.

This package is including :

- SmartScore Pro
- ECG trace
- Cardiac enhance filter
- CardIQ SnapShot
- SnapShot Pulse
- ECG Wave on Gantry
- SnapShot Assist
- SnapShot Freeze
- CardIQ Reveal on OC and CardEP on OC.

CardIQ SnapShot is an integrated cardiovascular helical image acquisition option.

The CardIQ SnapShot Software can be used to Acquire ECG Gated CT images of Cardiovascular anatomy with improved temporal resolution to reduce heart motion effects and artifacts. More specifically, with CardIQ SnapShot option, users can acquire cardiac images of patients using the following cardiac imaging techniques: (1) Retrospectively EKG-gated helical scanning method - SnapShot: primarily used for cardiac morphology imaging, with this technique, cardiac images of single or multiple cardiac phases at any given Z-axis location can be acquired and generated. (2) EKG-gated Multi-slice CINE Scan mode: used primarily for coronary artery calcification scoring (CACS) studies or for cardiac morphology Imaging. Once a specific imaging model is selected, helical pitch and/or gantry rotation speed will be automatically selected for optimal scan coverage and image quality.

SnapShot Pulse is an axial cardiac scanning technique that reduces patient significantly compared to a helical acquisition and improves cardiac workflow, without compromising image quality.

SnapShot Pulse uses prospectively triggered axial acquisitions synchronized by the patient heart rate, in which X-rays are turned on only during the required heart phase and turned off completely at all other times. In essence, the technique captures a complete picture of the heart using a series of three to four snap shots taken at precise patient table positions and precisely timed to correspond to a specific phase of the cardiac cycle. SnapShot Pulse helps improve workflow by reducing the size of image set to be reconstructed, reviewed and post processed.

A typical SnapShot Pulse series consists of 280 - 400 images, compared with up to 2,000 to 3,000 images in a typical helical cardiac scan series. Since there is a smaller number of images to reconstruct, SnapShot Pulse takes less time to transfer, process, and review, yet still delivers the same amount of clinical information as a helical cardiac exam.

g Imaging (Helical cardiac acquisition) and 0.35 sec rotation. Not available on GE 32 channel CT systems.

SnapShot Assist offers intelligent motion correction to enhance temporal resolution for coronary artery imaging. Snap Shot Assist provides assistance for the technologist for optimal cardiac protocol selection. Major sub-systems of the scanner have been re-engineered to improve imaging performance and reduce the radiation dose required for diagnostic studies.

**SnapShot Freeze** is an intelligent motion correction algorithm designed to reduce blurring of coronary arteries due to motion artifacts. SnapShot Freeze reduces motion artifacts up to 6X, equivalent to a 0.058s Equivalent Gantry Rotation Speed with Effective Temporal Resolution of 29msec\*. This benefit is delivered by characterizing the vessel motion (path and velocity) to derive the optimal vessel position at the target phase.

SmartScore Pro is CT Operator Console Acquisition Software for Prospective Gating.

The ECG trace provided by the Ivy monitor (option) will be displayed on the CT operator's console with this option. Allowing the user to display the live trace of the patient's heart rate and display the actual location of the window of time when the image is being acquired. It will provide easy access to patient cardiac output status and assist in providing visual feedback for optimum acquisition start.

Cardiac Image Filters provides users the capability to reconstruct filtered images using three steps of noise (pixel noise standard deviation) reduction for helical and axial cardiac imaging, which may allow a reduction of mA while maintaining an acceptable level of image performance.

ECG trace provides users the capability to display the heart rate and ECG waveform based on the data from the ECG equipment on the Xstream Display to review the patient heart rate during cardiac scanning.

CardIQ Xpress Reveal is an integrated post processing image analysis software dedicated for the application of cardiovascular imaging on the Console. The CardIQ Xpress Reveal software option can be used to effectively display, reformat and analyze 2D or 3D cardiac CT images for qualitative or quantitative assessment of heart anatomy and coronary artery vessels from single or multiple cardiac phase image data sets. CardIQ Xpress Reveal is launched via its own link or from Volume Viewer applications. It provides the user with both single and multiple cardiac phase analysis protocols.

The operator has a variety of different 2D, 3D or reformatted protocols to choose from to perform analysis and measurements. They include: display of coronary vessel tree, angiographic view, 2D and 3D rendering of single or multiple coronary artery vessels or grafts, automatic reformation of cross sectional cardiac images into planes along short or long axis of the heart, one-touch cath views for 3D or reformatted images, 3D angiographic view phase registration, plaque density measurements and color mapping of the non-calcified & calcified plaque, IVUS-like views, 3D ejection fraction, 4D aortic and Mitral valve views, relative perfusion, transparency views and beating heart images from single or multiple cardiac phase image data sets.

The CardIQ Xpress Reveal tool can be applied to standard axial or helical CT images. These images can be acquired on GE's multi-slice CT scanners using the cardiac CT SnapShot Pulse, Segment, Burst or Burst Plus imaging acquisition option.

### Clinical Benefits:

Cardiovascular CT imaging using multi-slice CT technology is a new and exciting clinical application which may make significant impacts to cardiovascular disease management as a non-invasive imaging technique. Multi-slice detector CT, which has been quickly adopted by the clinical community, has the advantage of being easy-to-use, reliable and accessible, as compared to other invasive or non-invasive cardiac imaging techniques. One of the critical

components for an effective cardiac CT application is a fully integrated post-processing and analysis tool tailored for cardiac imaging. The CardIQ Xpress Reveal option is designed to provide an easy-to-use and time-effective means for cardiovascular image manipulation. Clinical applications include: imaging of cardiac morphology, coronary artery imaging and assessment of relative perfusion, assessment of plaque, bypass graft patency, post intervention follow-up and functional assessment.

#### CardIQ Xpress Reveal simplifies user workflow by:

- Pre-processing the images & models for quick review of the exam
- Loading images into the auto launch area for real-time review of multiple exams
- Stream-lining protocol selection
- Easy switching from one protocol to the other without exiting the application
- Single click one-touch cath views
- Batch movie output within cardiac reformat
- User defined layouts within vessel analysis for simplified viewing and filming
- Multi-phase load to single phase review

#### The CardIQ Xpress reveal option allows the user to:

- Extract, render and display 2D/3D coronary vascular tree images with automatic vessel tracking & labeling with single click of a protocol. Images can be reviewed in axial, reformat, curved, oblique MPVR, and cross-section views.
- Various measurements of coronary artery vessels to include stenosis, density and length of stenotic areaPlaqID to color code non-calcified and calcified plaque with volume measurements.
- 2D reformat review with predefined views to review all coronary vessels.
- Color enhanced relative perfusion defect pattern recognition for detection of ischemic heart disease with 4 color patterns.
- Automatically render data for streamlined reading to include: 3D rendered heart, angiographic view, tree VR, and ejection fraction
- Reformat standard axial CT images of single or multiple cardiac phases automatically into short, long and two chamber long axis of the heart for easy review.
- Perform functional evaluation of the heart and cine capabilities for multiphase beating heart images with one easy click
- Automatic extraction of the left ventricle with automatic selection of ES and ED for ejection fraction & volume measurements. 4D aortic valve and mitral valve views with one touch Select protocols within the review step area allowing user to select a different protocol without exiting the application
- Pre-defined VR IVUS-like views for virtually determining the different compositions of the plaque
- One touch angiographic view protocol display coronary vessel tree and myocardium with automatic removal of heart chambers for cath comparative view.
- Heart transparency model allowing for full visualization of coronaries in relations to the heart chambers with the ability to fade out the chambers of the heart
- Oblique reformat views in the standard cath angles to provide an easy analysis of the coronary vessels
- Load multi-phase images, review the data and decide which phase or phases will be reviewed for further processing by dropping the non-essential phases
- Phase registration - ability to register images from different cardiac phases into a unique data set. The data set can then be saved as a 3D object and/or used for further analysis.

#### System requirements: VolumeViewer on Console

CardEP for CT Operators Console on the Xstream Platform. CardEP is Integrated Post Processing Image Analysis Software dedicated for Electrophysiology applications on GE's Xstream Operators Console.

The CardEP software option can be used to effectively display, reformat and analyze 2D or 3D Cardiac CT images for qualitative or quantitative assessment of cardiac chambers and veins.

The Operator has a variety of different 2D, 3D or reformatted images to choose from to perform analysis and measurements.

They include: automatic 3D VR atrium views, One touch 3D VR heart, multi-phase image reformation, automatic tracking of either pulmonary veins or coronary sinus with cardiac veins, reformation of cross-sectional pulmonary vein images with the ability to orient images along short or long axis of the heart, one touch EP views, EP navigator views, phase registration protocols and batch movie capabilities.

Along with these capabilities all protocols allow the user to load in multi-phase data for more accurate analysis. One of the critical components for an effective cardiac CT Application is a fully integrated post-processing analysis tool tailored to cardiac images. The CardEP option is designed to provide an easy-to-use and time-effective means for electrophysiologists to improve workflow.

## **B75542RE - CT Angiography Package**

CT Angiography package is innovative post-processing SW package and gives customers productivity improvement with decreased time to first clinically relevant Image and analysis. This package is including AVA Xpress and AutoBone Xpress.

CT AVA is a Highly Automated Software Post-Processing Package for the CT Operator's Console. It is an Additional Tool for the Analysis of 3D Angiography Data Providing a Number of Display, Measurement and Batch Filming/Archive Features to Study User-Selected Vessels Which Include Stenosis Analysis; Pre/Post Stent Planning Procedures and Directional Vessel Tortuosity Visualization.

### **Clinical Benefits**

- **Decreased Operator Dependence:** Currently there is Heavy Operator Dependence to Produce True Vessel Cross Sections and Vessel Profiles. This Software Eliminates the Need for the Operator to Manually Identify the Center of the Vessel.
- **Automated Batch Filming and the Ability to Rotate Around a Vessel,** Reduces the Risk of Overlooking Vascular Structures.
- **Quick AVA - Two click vessel analysis Measurement Tools:** Quantitative Information on User-Selected Vessel Segments, Aids in the Proper Selection of Prosthesis

Distances to Bifurcations or Other Landmarks are Critical for Clinical Decisions Increased Value of Reports: A Single Report Provides a Complete 3D Context; Measurements Cross-References and 3D Views. Consistency in the Format and Style of the Reports Also Help Referring Physicians.

### **Productivity Benefits**

**Decreased Time to First Clinically Relevant Image:** Automatic Centerline Detection Provides a Quick 3D Value Understanding of a Selected Vessel. The Anatomy Becomes Visible Once Two Points Identifying the Section of Interest Have Been Defined.

- **Background Auto-Filming:** Replaces Manual Filming.

AutoBone Express is a Software Package that provides Automatic Segmentation of Bony Structures and Calcified Plaques Optimized for the latest CTA Acquisition Techniques.

### **AutoBone Xpress Clinical Benefits:**

- **Click Segmentation of Bony Structures** to facilitate Vascular Structures Visualization for any Anatomy including Head and Neck CTA

- 1-Click Automatic Segmentation of Calcifications for Abdominal CTA and Run-Off Exams. Side-by-Side display of Vessels in 3D MIP with and without Calcifications provides a Direct Access to Calcified Plaques effect on Vessel Lumen.

Operator Productivity Benefits Include:

- Decreased time to First Clinically Relevant Image Segmenting Automatically Bony Structures and providing a Quick 3D MIP Overview of Vascular Structures.
- Synchronized Viewports enabling Fast confirmation of Results on Reformatted and Native Images.
- AutoSelect Segmentation Tools may be used to Refine Segmentation by Quickly Adding or Removing Structures.
- The resulting Volume Rendered Image can be manipulated to View Vessels Only.
- Transparent Bones can be restored for Landmarks. Calcifications can also be visualized in Transparency to Show Lumen.
- Optimized Layouts for each Anatomy for Fast and Relevant Visualization.

### System Requirements:

VolumeViewer on the Console

## B75682RE - Head Perfusion Package

Head Perfusion package is an innovative options package that allows 80mm wide coverage acquisition and image analysis for head Perfusion study. It is including VolumeShuttle and CT Perfusion 4D Neuro.

VolumeShuttle innovatively provides the 80-mm of coverage necessary for accurate dynamic neuro angiographic and perfusion studies with a single contrast injection. GE's exclusive real-time scan control, system architecture, and fast, smooth table acceleration and deceleration enable the patient to be effortlessly shuttled back and forth between two adjacent axial locations, with minimal inter-scan delay

The GE CT Scanner system uniquely designed to make it all possible - as a result of these key scanner attributes:

- The 40-mm high resolution V-Res detector with micro voxel technology
- Real-time system controls to precisely control table movement and X-ray control

**VolumeShuttle** provides the wider coverage margin needed to allow for patient variability in the Circle of Willis (80mm) and from the basal ganglia to lateral ventricles (>60mm) - all with the existing 40-mm-wide detector and without the multiple contrast injections necessary with today's standard CT systems. VolumeShuttle innovatively provides the 80-mm of Z-coverage necessary for accurate dynamic neuro angiographic and perfusion studies with a single contrast injection. GE's exclusive real-time scan control, system architecture, and fast, smooth table acceleration and deceleration enable the patient to be effortlessly shuttled back and forth between two adjacent axial locations, with minimal inter-scan delay.

The Revolution EVO system uniquely designed to make it all possible - as a result of these key scanner attributes:

- The 40-mm high resolution V-Res detector with micro voxel technology.
- Real-time system controls to precisely control table movement and X-ray control.

This technology works by scanning axially in one location and the moving the patient to an adjacent position in @1 second. Another axial acquisition is performed, followed by a shuttle back to the previous position. This cycle continues for the duration of the exam up to 40 seconds. Each cycle of two acquisitions is approximately 3 seconds. VolumeShuttle provides the wider coverage margin needed to allow for patient variability in the Circle of Willis (80mm) and from the basal ganglia to lateral ventricles (>60mm) - all with the existing 40-mm-wide detector and without the multiple contrast injections necessary with today's standard CT systems.

CT Perfusion 4D Neuro Package is an image analysis software package that allows the evaluation of dynamic CT data following an injection of a compact bolus of contrast material, generating information with regards to changes in image intensity over time. The software provides a quick and reliable assessment of the type and extent of cerebral perfusion disturbances by providing qualitative and quantitative information on various perfusion related parameters, which may be related to acute stroke, brain tumor angiogenesis and treatment thereof. The key perfusion parameters that CT Perfusion, 4D

Neuro Package generates are:

- Regional Blood Volume (BV; ml/100g)
- Regional Blood Flow (BF; ml/min/100g)
- Regional Mean Transit Time (rMTT;s)
- Capillary Permeability Surface Area Product (PS)
- Time of Arrival (IRF T0)
- Transit Time to IRF Peak (Tmax;sec)

The user now has the ability to visualize all the information in true volumetric form. Additional elements of Perfusion 4D include Smart Map, a new algorithm that improves the image quality of the functional maps in the presence of noise. Perfusion 4D also includes a new streamlined workflow for Tissue Classification. Tissue Classification may aid the clinician in determining the status of the tissue based on blood volume and one of blood flow, mean transit time, or Tmax.

Productivity is enhanced through the protocol driven design of the user interface. An example of this is the Brain Stroke Protocol (Automatic) that completes the processing with one touch reducing the time required to process the exam and to enhance repeatability.

System requirements: VolumeViewer on the Console.

## **B75572RE - Workflow Package**

Workflow package is software package for productivity enhancement. This package includes Image check and One-stop scanning mode.

**Image Check** provides 340x340 matrix images for confirming reconstructed image coverage in real time and tracking up to 1800mm length with less than 1 sec delay.

Reconstruction time is up to 55 fps.

Revolution EVO's exceptional one stop scanning mode provides a streamlined workflow on the Xstream Display such as "Patient selection", "Protocol selection" and "Confirm". Pre-scanning can be accomplished in as few as five touches.

**SECRET DE AFACERI**  
**CONFIDENTIAL**

## **B78552CA - Optima Operator Console Desk (Wide Design)**

Optima Desk is the desk designed with ergonomics. This table design enables the efficient use of space while enhancing workflow and technologist comfort. Benefits include:

- Improved ergonomics for technologist
- Wide desktop for working space improvement
- Flexible location of Console hardware

## **B77552DA - 2 Phase Uninterruptible Power Supply**

Vertiv Uninterruptible Power Supply with custom designed cables to interconnect with GE scanners. The UPS Primarily Backs Up the System Computer Functions.

Bridges Short Power Outages and Provides Time for Crossover from Normal Main Power to Emergency Power. Must be Located Within Eight Feet of the PDU.

## **B75352CA - TABLE Convenience kit**

Content: IV pole, IV Tray & Tray Assey - IV Board - IV Strap Wrist

## **B78792CB - UW protocol for Revolution EVO**

Dose-optimized protocols designed, developed and validated by a team of physicians, physicists and technologists at the University of Wisconsin -Madison School of Medicine and Public Health. Tested in a variety of settings including outpatient imaging centers, the American Family Children's Hospital and one of the nation's premier academic medical centers, these protocols are designed to be used "as is", requiring no adjustment by the operator. The protocols cover multiple clinical applications and are designed for patients of varying sizes. Each protocol has been designed to deliver diagnostic quality images while using as little dose as reasonably possible. Included with the protocols is a 300 plus page manual containing information on the proper patient set-up, contrast parameters, scan range instructions, reformatting instructions and technical acquisition parameters for each protocol.

## **R23053AC - Standard service tools L3 - warranty period**

GE Healthcare has reclassified its service tools, diagnostics and documentation into various classes (please refer to the Service Licensing Notification statement at the beginning of this Quotation). The Standard License provides access to service tools used to perform basic level service on the Equipment and is included at no charge for the warranty period.

## **M81521KA - AW VS7 workstation**

AW VolumeShare 7 with Two Flat Panel Monitors and 32GB of RAM.

AW VolumeShare 7 is a multi-modality image review, comparison and post processing workstation built with simplicity and power at its core. Powerful software is optimized to take advantage of state of the art 64 bit technology and multiple cores to ensure leading edge performance.

GE Healthcare

**SECRET DE AFACERI**  
**CONFIDENTIAL**

## AW VolumeShare 7 features include:

### Hardware:

- HP Z820 Workstation
- CPU: 2x Intel Xeon E5-2630 Six Core 2.6 GHz
- CPUs with 15MB Shared L3 Cache each and 1866 MHz Dual Front Side Bus
- RAM: 32GB (8x4GB) DDR3 1866 MHz ECC DIMM
- NVIDIA Quadro NVS 310, 512MB Graphics card
- 1x 300GB SAS 10k rpm Hard Disk for OS and Apps
- 2x 300GB SAS 10k rpm Hard Disks for Data
- 2 x 19" 1280x1024 color monitors

### Software:

- GE Healthcare HELiOS 6 operating system
- Volume Viewer for advanced post-processing
- Demo Exams for training and exploration
- Fast access to information you need through optional RIS integration & priors post-fetch
- Efficient workflow through dynamic load, end review and Key Image Notes features
- Productivity package to pre-process exams and allow up to 8 simultaneous sessions
- Applications usage monitor to track and view usage of your system
- Smart layouts with Volume Viewer General review protocol that optimizes comparison and single exam layouts
- Enhanced multi-modality contouring tool with support for PET SUVs
- Support for external DICOM USB media and preference management tool to exchange preferences across users
- Support for optional, broad suite of multi-modality advanced applications

## M80281AA – AW VolumeShare 7 Monitors

AW VolumeShare 7 Monitors are two high-quality monitors offering bright and high contrast imagery suited to the display of medical images per the AW VolumeShare Indications for Use.

Each provides a 19" 1280x1024 (5:4 aspect ratio) display that complies with international medical and patient safety standards and offers the following specifications:

- Maximum luminance (panel typical) : 330 nit
- DICOM Part 14 calibrated luminance: 215 nit
- Contrast ratio (panel typical) : 900:1
- An ambient light sensor
- Brightness non-uniformity (measured as per DIN6868-157) : +/-25%

## M80501DV - Field Engineer Letter

Field Engineer Letter - Operator Manual on paper for AW

**SECRET DE AFACFRI**

**CONFIDENTIAL**

## M81521ED - Integrated Registration

Full Fusion Package Integrated Registration will be delivered on AW VolumShare 7 or AW Server 3.2

Integrated Registration is designed to provide easy comparison of three dimensional (3D) anatomical images from Computed Tomography (CT), MRI (Magnetic Resonance Imaging), PET (Positron Emission Tomography), Single Photon Emission Computed Tomography (SPECT) and X-Ray Angiography (XA)\*.

It allows registration and fusion between two volumetric acquisitions, which come from either the same or from different acquisition modalities.

Major features and enhancements are:

- Ability to combine any two of the 5 modalities together.
- Automatic propagation of registration across series acquired in the same patient exam (i.e. same frame of reference) and to any series from any loaded exam that have been manually grouped together.
- Full compatibility of the 3 different registration methods: automatic, manual and landmark that can be combined together to provide an optimal result.
- 2D, 3D and hybrid 2D/3D Fusion capabilities.
- Access to Volume Viewer\*\* functionalities including MPR, Slab and oblique reformations, triple oblique easy definition, Volume Rendering, 3D display, distance and ROI measurements (the ROI measurement only work on the rigid registered images, not on the non-rigid registered images), layout management, segmentations, film and save.
- Ability to save registered data as new DICOM series or as Registered DICOM object (except from SPECT saving which is currently a limitation).
- Ability to draw and save contours as RTSS DICOM objects.

Summary of operation:

- User loads DICOM 3 CT, MR, PET, SPECT and/or XA data into a Integrated Registration protocol.
- Registration is performed based on reference and moving series selection.
- User reviews the quality of the registration with visualization tools and validates results.
- Optional: user defines and saves the contours of structures of interest.
- Registration results are saved.

\* For XA modality series, Integrated Registration currently supports only the 3D X Ray Angiography (i.e., 3D X-Ray Angiography images stored as CT Image Storage DICOM objects) images acquired with GE Innova equipment and reconstructed with the Innova3DXR application.

## M81521BS - ONCOQUANT MULTIMODALITY PACKAGE

The OncoQuant Multi-Modality Package contains OncoQuant and Integrated Registration - Full Fusion.

OncoQuant is an oncology workflow enhancement tool that provides multi-modality image and dataset reviews. It provides user-friendly tools to follow lesion size over time, apply study criteria and provide tabulated results to the Oncology team. It is an optional package available on both the GE AW VolumeShare7 and the AW Server 3.2 platforms.

OncoQuant is fully integrated within the standard Volume Viewer 7 protocols and therefore works as a toolset rather than a standalone application. Because of this tight integration, OncoQuant, as a product, benefits from the new Volume Viewer 7 improvements.

#### Major features and enhancements are:

- An integrated kit of Oncology Tools compatible in any standard reading protocol in Volume Viewer 7 to aid routine oncology reads.
- Adaptable Workflow for standard clinical reading to advanced research using tools supporting RECIST 1.0, 1.1 and WHO criteria.
- A Multi-Modality reading platform allowing comparison and correlation of CT, MR, PET/CT, and 3D X-Ray data. Features:
  - Automatic Multi-Modality full coverage and regional registration (CT, MR, PET, 3D X-Ray) if the user has purchased Integrated Registration.
- Dedicated smart PACS-like layout for facilitating oncology review and follow-up studies from the Volume Viewer 5.
- Full access to entire set of 3D visualization tools.
- One consistent contouring tool for all modalities (CT and MR images and PET SUV's).
- Benefit from Lung VCAR algorithms and DCA as a tool inside the routine oncology workflow if the user has purchased Lung VCAR.
- Advanced support for output with new oncology save state.
- New intuitive Summary Table of Findings supporting guided follow-up for standard or more advanced studies like RECIST.
- Export statistical results (for excel) and images to USB and in DICOM to the filmer.

## B79921TA - CARDIQ FUNCTION XPRESS

### CardIQ Function Xpress CardIQ

Function Xpress is an Integrated post processing image analysis software for Cardiovascular CT on GE's Advantage Workstation and AW Server.

CardIQ Function Xpress allows users to Non-Invasively Image the Functional Parameters of the Heart such as ejection fraction (EF) and ventricular volumes. CardIQ Function Xpress uses multi-phase cardiac gated datasets for processing.

The software automatically detects endocardial and epicardial contours for assessment of left ventricular (LV), right ventricular (RV) and left atrial functional parameters.

- Automatically select each chamber of the heart for individual chamber volume analysis. The software automatically selects the LV 97%, LA 87% and RV 96% of the time.
- Automatic end diastolic and end systolic selection for LV, RV and LA ejection fractions >91% of the time.
- Behind the scene processing & loading of function data for real time review of ejection fraction, volume analysis and myocardial analysis.
- One click activation for 4D beating heart
- Automatic selection of epicardium and endocardium for myocardial analysis.
- Bulls eye plots representing wall motion, wall thickness and wall thickening.
- Automatic left atrium volume calculation with exclusion of the pulmonary veins.

**CONFIDENTIAL**  
**SECRET DE AFACERI**

- Single click visual wall motion activation with short axis images in the basal, mid and distal orientation along with a 2-chamber long axis view.
- Flexible reporting tool with graphical representations.
- Display table of key functional parameters for instant visualization.

### System requirements:

AW Workstation with VolumeShare 7 or later or AW Server 3.2 or later.  
Requires AutoLaunch and Pre-Processing.

## B79821RE - CARDIQ XPRESS 2.0 REVEAL

CardIQ Xpress 2.0 Reveal is an integrated post processing image analysis software for Cardiovascular CT on GE's Advantage Workstation.

The optional CardIQ Xpress Reveal software can be used to effectively display, reformat and analyze 2D, 3D, and GSI CT images for qualitative or quantitative assessment of the anatomy of the heart and coronary artery vessels from single or multiple cardiac phase image data sets. When used with CardIQ Function, CardIQ Xpress Reveal can also provide functional assessment including relative perfusion information.

CardIQ Xpress Reveal can be launched directly or from within Volume Viewer applications using gated axial, helical or GSI CT images; including images created using the SnapShot Freeze intelligent motion correction option.

The software includes a variety of different 2D, 3D or reformatted protocols including:

- Display of the coronary vessel tree, angiographic view
- 2D and 3D rendering of single or multiple coronary artery vessels or grafts
- Automatic reformation of cross sectional cardiac images into planes along short or long axis of the heart
- One-touch cath views for 3D or reformatted images
- 3D angiographic view phase registration
- Color mapped plaque density measurements
- IVUS-like views
- 3D ejection fraction
- 4D aortic and Mitral valve views
- Relative perfusion
- Transparency views and beating heart images from single or multiple cardiac phase image data sets.

CardIQ Xpress Reveal combines simplified user workflow with SnapShot Freeze intelligent motion correction imaging.

- Pre-processing the images & models including SnapShot Freeze exams, for faster review
- Loading images into the auto launch area area for real-time review of multiple exams
- Easy switching from one protocol to the other without exiting the application
- Single click one-touch cath views
- Batch movie output within cardiac reformat
- User defined layouts within vessel analysis for simplified viewing and filming
- Multi-phase load to single phase review

**SECRET DE AFACERI**  
**CONFIDENTIAL**

The CardIQ Xpress reveal option allows the user to:

- Rendering and display of 2D/3D coronary vascular tree images with automatic vessel tracking & labeling with single click of a protocol. Images can be reviewed in axial, reformat, curved, oblique MPVR, and cross section views
- Measurements of coronary arteries including: stenosis, density and length of stenotic area
- PlaqID to color code non-calcified and calcified plaque with volume measurements.
- 2D reformat review with predefined views to review all coronary vessels.
- Color enhanced relative perfusion defect pattern recognition for detection of ischemic heart disease with 4 color patterns
- Automatically render data for streamlined reading to include: 3D rendered heart, angiographic view, tree VR, and ejection fraction.
- Reformat standard axial CT images of single or multiple cardiac phases automatically
- into short, long and two chamber long axis of the heart for easy review Perform functional evaluation of the heart and cine capabilities for multiphase beating heart images with one easy click
- Extraction of the left ventricle and automated ejection fraction and volume measurements. Note: CardIQ Function Xpress is needed if myocardial wall motion, mass, wall thickness or chamber volumes for the Right Ventricle, Left Atrium, Right Atrium is needed.
- 4D aortic valve and mitral valve views with one touch
- Ability to select different protocols without exiting the application
- Pre-defined VR IVUS-like views for virtually determining plaque compositions
- One touch angiographic view protocol displays coronary vessel tree and myocardium with automatic removal of heart chambers for cath comparative view
- Heart transparency model allowing for full visualization of coronaries in relations to the heart chambers with the ability to fade out the chambers of the heart
- Oblique reformat views in the standard cath angles for easy analysis of the coronary vessels
- Load multi-phase images, review the data and decide which phase or phases will be reviewed for further processing by dropping the non-essential phases
- Phase registration - ability to register images from different cardiac phases into a unique data set. The data set can then be saved as a 3D object and/or used for further analysis

#### System requirements:

- AW VolumeShare 7 or AW Server 3.2
- Auto Launch and Preprocessing Option

### B79971JH - SmartScore 4.0

SmartScore 4.0 provides a valid baseline test that helps detect calcium plaque burden using minimal radiation dose. It allows to better assess a patient's overall risk for coronary artery disease and communicate results to patients in comprehensive, easy-to-read reports.

Automatic detection reduces or eliminates the need for manual visual detection of lesions, with all scorable lesions above a threshold highlighted in bright green. Lesions are also correctly classified in a 3D space with fewer mouse clicks.

SmartScore 4.0 uses the conventional Agaston-Janowitz calcium scoring technique with a threshold of 130 HU and is adjusted for the appropriate image slice thickness.

The Mass scoring algorithm is quickly being adopted as an objective measuring tool, reported in milligrams, for a more quantifiable calcium calculation. SmartScore 4.0 also computes a Volume score algorithm for calcium quantification in cubic millimeters.

#### Features:

- Fully integrated patient reports
- Personalized reports (from a patient's overall cardiac history to exam-specific images, diagrams, tables and calcium rankings.
- Calcium distribution and population graphs
- PACS function for future retrieval

## B77121BK - VesselIQ XPRESS AND AUTOB

VesselIQ Xpress & AutoBone Xpress VesselIQ Xpress provides an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment plans from a set of CTA images.

There are new features introduced in the VolumeShare 7 release including:

- Auto Abdominal Aorta Vessel tracking which is a completely automated protocol with autobone removal, auto vessel tracking and automatic labeling of the abdominal aorta vasculature.
- Fast Tracking which provides automatic real-time feedback for auto-detected centerlines to speed up vessel tracking.
- New editing tools that allow for flexibility in editing based on the size of the vessel being edited.

This software supports the physician in:

- Assessment of aneurysms with or without thrombus (false lumen) for size and volume measurements with the capability to track the size and volume over time, stenosis analysis, pre/post stent and surgical planning and directional vessel tortuosity visualization.
- Automatic tools for the segmentation of bony structures in the brain and neck and other vascular areas for accurate identification of the vessels, single or double click vessel analysis.
- Sizing the vessel, analyzing calcified and which is a completely automated protocol non-calcified plaque to determine the densities of plaque within a vessel, measure areas of abnormalities within a vessel (like stenosis, plaque, thrombus, dissection or leakage).
- Semi-automated detection and segmentation of thrombus for subsequent measurements within the application.
- Dedicated anatomy based protocols for improved workflow.
- Compare a patient's previous exam to their current exam in order to measure and track any changes over time of their vascular structures.
- After review of the exams, there are multiple ways to film, archive and capture information for future review.

System Requirements: AW VolumeShare 7 or AW Server 3.2

Note: All softwares are Non-Transferable to other hardware and are Non-Returnable

**SECRET DE AFACERI**  
**CONFIDENTIAL**

## E8007TB - Cardiac Cable Kit for IPC, Ivy 3150-B, Ivy 7800

Cardiac Cable Kit for IPC

## E8016AZ - Slicker for 1700 Table

- Increase system uptime by protecting table from spills and particulate contaminants
- Easy to install and comfortable for patients
- Will not interfere with normal operation of CT table
- Clear PVC plastic facilitates faster cleanup of blood and fluids
- Prevents contaminate build up in hard to clean areas
- Thermosealed seams and flaps
- Recommended for trauma centers and sites concerned about exposure to blood and fluid-borne disease

Two pieces/set, CT cushion and cover.

## A82000CT - Welcome pack - 16 Credits

Get the best use of your GE Equipment with the GE CARES program.

The Welcome pack is an option offering you 16 credits to be used at your convenience, after the Start me Up application days.

The Education credits packages are designed to provide flexible training options to be used in promoting learner retention, supporting employee turnover needs and allows for efficient and effective skill development. Credits may be used for Clinical Education on GEHC Diagnostic Imaging products located at Customer's facilities.

Credits may be used for trainings conducted at Customer's facility, via remote training sessions, at GE Healthcare Academies and GE Healthcare partner sites as follows (upon availability):

| Education Service                                                                                       | Number of training credits |
|---------------------------------------------------------------------------------------------------------|----------------------------|
| <i>On-site Education Services</i>                                                                       |                            |
| 1-day onsite training at customer's facility                                                            | 8 credits                  |
| <i>Off-site Education Services</i>                                                                      |                            |
| 1-day classroom for 1 person                                                                            | 4 credits                  |
| 1-day immersion for a radiographer                                                                      | 4 credits                  |
| <i>Online Education Services</i>                                                                        |                            |
| 1-hour remote training                                                                                  | 1 credit                   |
| 1-day remote classroom                                                                                  | 4 credits                  |
| Access to 1 remote only Elevate package for 5 participants                                              | 8 credits                  |
| Access to 1 remote Elevate package for 5 participants with 1-day onsite training at customer's facility | 16 credits                 |

Credits are valid for 12 months from Equipment warranty start, or date of signed agreement if purchased standalone. Unused Credits at the end of this period are forfeited without refund or credit. Additional credits may be available for purchase separately.

## A11181CT - Start Me Up 8 days CT

### Customized on-site Application training


Get the best use of your GE Equipment with the GE CARES program.

The program is divided into 3 steps, it is tailored to adapt to staff's level of expertise to deliver the most adapted training.

We define together in advance your application needs and pre-training for an efficient on-site training experience (8 days).

After your on-site application, we offer you unique access to an expert for live assistance and application training follow-up directly on your equipment.

We support you throughout the life of your imaging equipment by providing privileged access to interactive clinical webinars with experts and on-line courses

**Duration:** 8 days onsite training

**Target attendees:** Radiographers, Technologists, Radiologists, Physicians, Cardiologists

**Course Leader:** GE Clinical Education Specialist

**Certification:** GE Healthcare Education Certificate

**Language:** Local Language when available or English

**Number of Participants:** Maximum 4 people per session

**Expiry date:** Available during the warranty period

## A33331CT - Make Me an Expert CT

### Customized on-site Application training

Get the best use of your GE Equipment with the GE CARES program.

GE Healthcare


CONFIDENTIAL  
SECRET DE AFACERI

The program is divided into 3 steps, it is tailored to adapt to staff's level of expertise to deliver the most adapted training.

We define together in advance your application needs and pre-training for an efficient on-site training experience. After your on-site application, we offer you unique access to an expert for live assistance and application training follow-up directly on your equipment.

We support you throughout the life of your imaging equipment by providing privileged access to interactive clinical webinars with experts and on-line courses

**Target attendees:** Radiographers, Technologists, Radiologists, Physicians, Cardiologists

**Course Leader:** GE Clinical Education Specialist

**Certification:** GE Healthcare Education Certificate

**Language:** Local Language when available or English

**Number of Participants:** Maximum 4 people per session

**Expiry date:** Available during the warranty period

## A82016CT - Additional 16 days of on site application training for CT

At GE Healthcare we believe that continual personal development, learning and improving clinical practice is essential for all Healthcare professionals.

Demands on healthcare professionals' time is continually increasing and the GE Healthcare Clinical Education Department is focused on providing relevant education solutions. Our goal is to provide opportunities to gain the knowledge and skills to optimize equipment usage, clinical practice and patient care.

The Training credits packages are designed to provide flexible training options to be used in promoting learner retention, supporting employee turnover needs and allows for efficient and effective skill development. Credits may be used for Clinical Education on GEHC Diagnostic Imaging products located at Customer's facilities.

Training credits may be used for trainings conducted at Customer's facility (On-site Education Services), via remote training sessions (Online Education Services), at GE Healthcare Academies and GE Healthcare partner sites as follows (upon availability) (Off-site Education Services):

GE Healthcare

SECRET DE AFACERI  
CONFIDENTIAL

23/24

| Education Service                                                                                       | Number of training credits |
|---------------------------------------------------------------------------------------------------------|----------------------------|
| <b>On-site Education Services</b>                                                                       |                            |
| 1-day onsite training at customer's facility                                                            | 8 credits                  |
| <b>Off-site Education Services</b>                                                                      |                            |
| 1-day classroom for 1 person                                                                            | 4 credits                  |
| 1-day immersion for a radiographer                                                                      | 4 credits                  |
| <b>Online Education Services</b>                                                                        |                            |
| 1-hour remote training                                                                                  | 1 credit                   |
| 1-day remote classroom                                                                                  | 4 credits                  |
| Access to 1 remote only Elevate package for 5 participants                                              | 8 credits                  |
| Access to 1 remote Elevate package for 5 participants with 1-day onsite training at customer's facility | 16 credits                 |

This 16-training credits Education Package is valid for 12 months from the product warranty start date when associated with a product purchase; or the effective date of this Agreement for Education Services standalone purchase.

If not spent within this time period, other than because of GEHC's fault, the training credits will expire without refund.

Additional credits may be available for purchase separately.

*GEHC accepts to provide Customer with the Education Services described in the General Terms and Conditions attached hereto based on the training credits package/pool chosen by the Customer in this Quotation. By signing the Quotation, including this long description, the Customer declares having fully read and understood the General Terms and Conditions of Education Services of GEHC as well as this Quotation and fully agrees with and accepts such terms.*

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## Revolution™ EVO

### User Manual

GE Healthcare Japan Corporation does business as GE Healthcare

This manual supports the following configurations:

- Revolution™ EVO EL
- Revolution™ EVO EX
- Revolution™ EVO ES
- Revolution™ EVO
- Revolution™ EVO 48kW

Not all configurations are available in all regions.

This product is certified as a Revolution™ EVO CT Scanner.

Revolution™ EVO  
Operator Manual, English  
5805440-1EN  
Revision: 1  
© 2019 General Electric Company  
All rights reserved.

CONFIDENTIAL  
SECRETE DE AFACERE

## Gantry

### Gantry display

The gantry display provides gantry and table status information.

Figure 7-9: Gantry display with patient information

Table 7-7: Gantry display components

| No. | Component                                                                                                               | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----|-------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.  | <b>Collision indicator</b><br>       | <p>The Collision light illuminates and a picture of CT system with high-lighted collision sensors is posted on the gantry display (Figure 7-10) when there is a possibility that the table, cradle and gantry comes in contact with each other. There are tilt/elevation collision sensors on the front and rear covers at the top of the gantry opening. The collision sensors are active during table elevation and during tilting gantry. The light also illuminates if there is 4.5 kg (10 lb) of resistance during cradle motion in/out or when reaching the travel limits of table, cradle and gantry motion. In the case of interference, you may need to raise or lower the table, check the gantry tilt, or determine if resistance is being caused by any patient restraints, accessories, or by the patient. If you need to tilt the gantry after a collision has been detected, the gantry only tilts in 0.5 degree increments.</p> <p> The picture of CT system with high-lighted collision sensors on the gantry display can be hidden by touching the display.</p> |
| 2.  | <b>Cradle unlocked indicator</b><br> | <p>The Cradle Unlocked Indicator is illuminated, indicating the cradle is unlocked. When the cradle is unlocked, it is "free floating".</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 3.  | <b>Laser light indicator</b>                                                                                            | <p>The Laser Light Indicator is illuminated when the laser alignment</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |


| No. | Component                                                                                                                   | Description                                                                                                                                                                                                                                                                                                                                                        |
|-----|-----------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|     |                                            | light is on. Normally, this light is only on when positioning a patient.<br> <b>CAUTION:</b> For patient safety, it is important to always have patients close their eyes anytime the laser alignment light is on.                                                                |
| 4.  | Lung indicator<br>                         | The lung indicator is not lit on this system.                                                                                                                                                                                                                                                                                                                      |
| 5.  | Heart indicator<br>                        | The heart is lit whenever there is an ECG gating signal seen by the system.                                                                                                                                                                                                                                                                                        |
| 6.  | Vertical height indicator<br> 123.4        | The Vertical Height Indicator displays the vertical height of the table in relation to ISO center.                                                                                                                                                                                                                                                                 |
| 7.  | Horizontal Cradle Position<br> S 12345.6 | The Horizontal Cradle Position displayed is the position of the cradle based on the established anatomic reference of the patient. This reference is established using the internal or external landmark. The number is preceded by an S if the position is superior to the reference point or an I if the position is inferior to the reference point.            |
| 8.  | Gantry tilt indicator<br> I 12.3°        | The Gantry Tilt Indicator shows one of two displays. One display, which is the normal display, indicates the current gantry tilt. If the gantry tilt is superior, the number is preceded with an S. If the gantry tilt is inferior, the number is preceded with an I. This is enlarged during Prescribed Tilt button on the Scan Control Interface Module flashes. |
| 9.  | Patient Name*<br>GE Healthcare                                                                                              | Patient Name is displayed as input on Patient Information screen. <ul style="list-style-type: none"><li>• The system automatically display 2 lines when the name is long.</li><li>• Font size is automatically changed by the number of characters.</li></ul>                                                                                                      |
| 10. | Patient ID*<br>123456                                                                                                       | Patient ID is displayed as input on Patient Information screen. <ul style="list-style-type: none"><li>• Font size is automatically changed by the number of characters.</li></ul>                                                                                                                                                                                  |
| 11. | Age, Birthdate and Gender<br>99M                                                                                            | Age or Birthdate, and Gender are displayed as input on Patient Information screen. <ul style="list-style-type: none"><li>• The patient below one year old, 0 is displayed in Age.</li><li>• Font size is automatically changed by the number of characters.</li></ul>                                                                                              |

SECRET DE AFACERI

| No. | Component                                                                                                                       | Description                                                                                                                                                                             |
|-----|---------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 12. | <b>Exam Description*</b><br><b>Chest-Abdomen</b>                                                                                | Exam Description is displayed as input on Patient Information screen. <ul style="list-style-type: none"> <li>Font size is automatically changed by the number of characters.</li> </ul> |
| 13. | <b>Default Patient Positioning buttons</b><br> | These buttons define the location where the cradle stops with a set pedal. The location can be set by GE Field Engineers.                                                               |
| 14. | <b>Tools button</b><br>                        | There are several functions related to the Display setting.                                                                                                                             |
| 15. | <b>Movie button</b><br>                       | There are some movies that help to explain CT examination procedures to patients.                                                                                                       |

\*Components can Hide/Show during an exam. If you want not to show the Patient Information item during every exam or want to select show components, this can be set by a GE Field Engineer. Contact GE Service to change the selection.

Figure 7-10: A picture of a system with highlighted collision sensor

SECRET DE AFACERI  
CONFIDENTIAL

If the system has Emergency Patient on the Gantry display option, the Gantry display includes 2 buttons.

Figure 7-11: Gantry display with Emergency Patient Option and ECG trace on Gantry display option



Table 7-8: Gantry display component with Emergency Patient Option and ECG trace on Gantry display option

| No. | Component                                                                                                        | Description                             |
|-----|------------------------------------------------------------------------------------------------------------------|-----------------------------------------|
| 1.  | <div>ECG trace</div>          | To display ECG trace on Gantry display. |
| 2.  | <div>Emergency Patient</div>  | To start on emergency exam.             |

Figure 7-12: Enabled ECG trace on Gantry display

CONFIDENTIAL  
SECRET DE AFACERI

### Breath hold lights

Figure 7-19: Breath hold lights


Table 7-13: Breath hold light components

| No. | Description                         |
|-----|-------------------------------------|
| 1.  | Breathe light                       |
| 2.  | Hold breath light                   |
| 3.  | Elapsed breath hold time in seconds |

### Internal components

Table 7-14: Gantry internal components

| Component                              | Description                                                                                                                                                                                                                                 |
|----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Tube and Col-<br>limator               | The anode heat capacity with a Performix™ 40 Plus tube is 7.0 million heat units (MHU) and the cooling rate is 1,070 kilo heat units per minute (KHU/min).<br><br>Performix™ 40 Plus operates at 48kW or 72kW, depending on configuration.  |
| High Fre-<br>quency Gen-<br>erator     | The high frequency generator is composed of an auxiliary box, inverter and high voltage tank. Each tank occupies one cubic foot in the gantry. Together the cathode and the anode provide 72 kilowatts (kW) of the Performix™ 40 Plus tube. |
| Detector/Data<br>Acquisition<br>System | The detector is comprised of a solid scintillator material known as HiLight. HiLight detectors, provides 98 % absorption efficiency. The V-Res™ Detector enables Volume CT with a true 64-channel design with 40mm coverage.                |
| Slip Ring                              | Slip Ring provides data communication path from Detector/DAS to the Digital Interface Processor from there to the system disk.                                                                                                              |

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