

4008 B / 4008 H / 4008 S Hemodialysis Device

Technical Manual

Software version: $\geq 4.5 / \geq 5.3$

Edition: 12A-2016

Part no.: F40003413



**FRESENIUS
MEDICAL CARE**

4008 B / 4008 H /
4008 S

TSC / MA report

FRESENIUS
MEDICAL CAREThese checks must be performed every **24 months** if all of the following requirements have been met:

- Unique assignment of the rotors to the appropriate tube roller pumps
- Software:
4008 B system: from 4.951 to < 5.00 or with 5.201 or higher
4008 H/S systems: from 2.951 to < 3.00 or from 4.311 to < 10.00
- DIP switch PCB LP632 array 2 SW 5 set to OFF.
- DIP switch PCB LP631 array 2 SW 7 set to ON.
- DIP switch PCB LP631 array 2 SW 8 set to ON. (With CDS option only)
- Heater rod made of Titanium

If one of the requirements mentioned above has not been met, the checks have to be performed every **12 months**.

Interval:

- 24 months
- 12 months
(Please check)

Name of technician

Service report no.:

Customer / Customer no.:

Serial number:

Inventory no.:

Operating hours:

- Device type: 4008 B 4008 H 4008 S
- With option: DIASAFE[®] plus OCM ONLINEplus™ SN bibag[®]
- BPM BTM BVM DIASAFE ONLINE HDF 4008 HDF

	No.	Description	Meas. value	✓
	1	Visual inspections		
TSC	1.1	Fuses accessible from the outside comply with the specified values	-	<input type="checkbox"/>
TSC	1.2	Labels and inscriptions are present and legible	-	<input type="checkbox"/>
TSC	1.3	Mechanical condition permits further safe use	-	<input type="checkbox"/>
TSC	1.4	There are no signs of damage or contamination	-	<input type="checkbox"/>
TSC	1.4.1	No signs of damage on the tube roller pump rotors	-	<input type="checkbox"/>
TSC	1.4.2	PCB LP 928 (CAN bus distributor) with covering checked (only 4008 S with E-Code 905 and lower)	-	<input type="checkbox"/>
TSC	1.5	Power cable not damaged	-	<input type="checkbox"/>
MA	1.6	Preventive measures		
MA	1.6.1	Sealing plungers in suction tubes replaced and lubricated; rivets replaced	-	<input type="checkbox"/>
MA	1.6.2	Rubber in the rinse chambers checked for proper function	-	<input type="checkbox"/>
MA	1.6.3	Rinse chamber screws are tight	-	<input type="checkbox"/>
MA	1.6.4	Check valve replaced and, if required, CDS check valves 117/118 and filter 119/120 replaced	-	<input type="checkbox"/>
MA	1.6.5	Pre-UF pump filter, filter downstream of MV 43, filter between rinse chambers and on MV 99, MV 100 replaced	-	<input type="checkbox"/>
MA	1.6.6	Dialysate filter replaced or sieve changed	-	<input type="checkbox"/>
MA	1.6.7	O-rings in dialyzer couplings replaced	-	<input type="checkbox"/>
MA	1.6.8	Sampling valve is functional	-	<input type="checkbox"/>
MA	1.6.9	Fan filter replaced	-	<input type="checkbox"/>
MA	1.6.10	Running band and tube segment of air separation pump replaced	-	<input type="checkbox"/>
MA	1.6.11	MV 84 (yellow label), replaced after 2 years. (Only applicable when using Puristeril 340, Puristeril plus)	-	<input type="checkbox"/>
MA	1.6.12	Connecting piece or equilibration chamber replaced (Only when ONLINEplus™ and DIASAFE [®] plus options are not used)	-	<input type="checkbox"/>
MA	1.6.13	Filter 210 (if present) replaced	-	<input type="checkbox"/>
MA	1.6.14	Disinfectant filter replaced	-	<input type="checkbox"/>
MA	1.6.15	There are no dirty or worn tubes	-	<input type="checkbox"/>
MA	1.6.16	Every 4 years only: Battery replaced (ensure correct polarity!)	-	<input type="checkbox"/>
	2	General checks		
TSC	2.1	Power failure alarm – continuous sound – text displayed: Emergency operation	-	<input type="checkbox"/>

	No.	Description	Meas. value	✓
TSC	2.2	Air separation by air separation pump activated; text displayed if more air must be separated and OD senses blood: Fill program	-	<input type="checkbox"/>
TSC	2.3	Check DIP switches PCB LP 631 (CPU1) DIP switch array 2 SW7 is set to ON. With Central Delivery System: PCB LP 631 (CPU1) DIP switch array 2 SW8 is set to ON.	-	<input type="checkbox"/>
3 Check of the hydraulics				
MA	3.1	Water inlet pressure (reduced) 0.9 bar to 1.4 bar	_____	<input type="checkbox"/>
MA	3.2	Loading pressure 1.25 bar ±0.05 bar From EC495 (4008 B); EC295 (4008 H) and EC275 (4008 S) Loading pressure 1.45 bar ±0.05 bar	_____	<input type="checkbox"/>
MA	3.3	Degassing pressure: -0.81 to -0.85 bar	_____	<input type="checkbox"/>
MA	3.4	Relief pressure of balancing chamber at 800 ml/min Loading pressure: Relief pressure: 1.2 to 1.3 bar 1.9 to 2.1 bar 1.45 ±0.05 bar 2.2 ±0.05 bar	_____	<input type="checkbox"/>
4 Ultrafiltration system and membrane pumps				
TSC	4.1	UF pump, 1 stroke = 1 ml, 60 strokes = 60 ml ±0.5 ml (59.6 g ±0.5 g dialysis water)	_____	<input type="checkbox"/>
MA	4.2	Concentrate pump calibration volume removal / number of strokes	_____	<input type="checkbox"/>
MA	4.3	Bicarbonate pump calibration volume removal / number of strokes	_____	<input type="checkbox"/>
5 Treatment				
MA	5.1	Desired temperature 37 °C ±0.5 °C	_____	<input type="checkbox"/>
MA	5.2	Temperature display 37 °C ±0.5 °C	_____	<input type="checkbox"/>
MA	5.3	Dialysate pressure - Flow off zero point checked - Slope checked	_____ _____	<input type="checkbox"/> <input type="checkbox"/>
TSC	5.4	Conductivity display checked with reference meter If the bibag [®] option is used, connect a bibag [®] ! - CD device - CD ref.	_____ _____	<input type="checkbox"/> <input type="checkbox"/>
6 Extracorporeal components				
MA	6.1	Arterial pressure display checked with reference meter	-	<input type="checkbox"/>
MA	6.2	Venous pressure display checked with reference meter	-	<input type="checkbox"/>
TSC	6.3	Blood pumps: Check of blood pump rate (calibration program: BP rate TEST)	-	<input type="checkbox"/>
TSC	6.4	SN switching pressure checked according to table in the Technical Manual	-	<input type="checkbox"/>
TSC	6.5	Blood pump stop alarm checked	-	<input type="checkbox"/>
TSC	6.6	Venous tube clamp closes after blood alarm	-	<input type="checkbox"/>
TSC	6.7	Pressure of approx. 2 bar in the venous bubble catcher. Pressure must not drop by more than 0.1 bar within 3 minutes	-	<input type="checkbox"/>
7 Options				
7.1 bibag[®]				
MA	7.1.1	bibag [®] connector, O-rings replaced	-	<input type="checkbox"/>
MA	7.1.2	Switching pressure of PSW 134 checked, 130 mbar, + 30 mbar	_____	<input type="checkbox"/>
7.2 DIASAFE				
MA	7.2.1	DIASAFE filter life checked	-	<input type="checkbox"/>
MA	7.2.2	Hydrophobic filter 111 replaced	-	<input type="checkbox"/>
MA	7.2.3	O-rings in the DIASAFE's dialysate couplings checked	-	<input type="checkbox"/>
7.3 DIASAFE^{® plus}				
MA	7.3.1	DIASAFE ^{® plus} filter life checked	-	<input type="checkbox"/>
MA	7.3.2	Hydrophobic filter 111 replaced	-	<input type="checkbox"/>
7.4 4008 HDF				
TSC	7.4.1	2nd UF pump, 1 stroke = 1 ml, 60 strokes = 60 ml ±0.5 ml (59.6 g ±0.5 g dialysis water)	_____	<input type="checkbox"/>
7.5 ONLINE-HDF (and DIASAFE)				
MA	7.5.1	Filter life of DIASAFE and ONLINE filter checked	-	<input type="checkbox"/>
MA	7.5.2	Hydrophobic filter 111 replaced	-	<input type="checkbox"/>
MA	7.5.3	O-rings in the DIASAFE's dialysate couplings checked	-	<input type="checkbox"/>
MA	7.5.4	HDF pump rotor checked (smooth running, wear)	-	<input type="checkbox"/>
MA	7.5.5	Fastening strap for Luer-lock checked	-	<input type="checkbox"/>
TSC	7.5.6	Substitute pump <input type="checkbox"/> with DC motor: Pump adjusted or <input type="checkbox"/> with stepper motor: Pump rate checked (calibration program: HDF-Pump-Rate Check)	Desired / actual: ____ / ____	<input type="checkbox"/>

	No.	Description	Meas. value	✓
TSC	7.5.7	Substitute pump stop – Following blood alarm – After triggering the bypass function – After opening the blood pump door	– – –	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
TSC	7.5.8	Check substitute pump for proper function: – Rinse program, delivery rate: 400 ml/min – Hot rinse program, delivery rate: 150 ml/min – Disinfection program, delivery rate: 400 ml/min	– – –	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	7.6	ONLINEplus™ (and DIASAFE® plus)		
MA	7.6.1	Filter life of DIASAFE® plus and ONLINEplus™ checked	–	<input type="checkbox"/>
MA	7.6.2	Hydrophobic filters 111 and 184 replaced	–	<input type="checkbox"/>
MA	7.6.3	O-rings at substitute port 195 and rinse port 194 replaced. O-rings lubricated with Unisilikon grease before installation	–	<input type="checkbox"/>
TSC	7.6.4	Tube in the tube pinch valve 193 (ONL1) replaced	–	<input type="checkbox"/>
TSC	7.6.5	Check of DIP switch PCB LP632 (CPU2) DIP switch array 2 switch 5 is set to OFF in HPU (hydraulics processing unit)	–	<input type="checkbox"/>
	7.7	OCM		
MA	7.7.1	OCM temperature / conductivity compensation test completed – CD cell 7 – CD cell 110 – OCM PULSE calibration required?	– _____ _____ <input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	7.8	BPM 4008 The Technical Measurement Checks (TMC) additionally stipulated for specific countries form part of the Technical Safety Checks (TSC).		
MA	7.8.1	Labels and inscriptions are present and legible	–	<input type="checkbox"/>
MA	7.8.2	Mechanical condition permits further safe use	–	<input type="checkbox"/>
MA	7.8.2.1	Tube connection correctly connected to device	–	<input type="checkbox"/>
MA	7.8.2.2	Internal blood pressure module, printed circuit boards, cable connections properly attached	–	<input type="checkbox"/>
MA	7.8.2.3	Damaged tubes or cuffs have been replaced	–	<input type="checkbox"/>
MA	7.8.2.4	Indicating elements checked	–	<input type="checkbox"/>
MA	7.8.2.5	Touch panel checked	–	<input type="checkbox"/>
TSC TMC	7.8.3	Leakage test: Pressure leakage rate less than 6 mmHg/min	–	<input type="checkbox"/>
TSC TMC	7.8.4	Calibration: Pressure values Tolerance 250 mmHg ±3 mmHg 200 mmHg ±3 mmHg 150 mmHg ±3 mmHg 100 mmHg ±3 mmHg 050 mmHg ±3 mmHg	_____ _____ _____ _____ _____	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
TSC TMC	7.8.5	Safety valve: Drained at 320 mmHg, ±10 mmHg	–	<input type="checkbox"/>
TSC TMC	7.8.6	Blood pressure measurement performed	–	<input type="checkbox"/>
	7.9	BTM 4008		
		No further Technical Safety Checks and maintenance procedures must be performed		
	7.10	BVM 4008		
		No further Technical Safety Checks and maintenance procedures must be performed		

	No.	Description	Meas. value	✓
	8	Check of the electrical safety According to (DIN) EN 62353:2008, IEC 62353:2007 Measurement points: (see Chapter 3.2.3.1 on page 3-39) For 4008 HDF option, check additional measuring point!		
TSC	8.1	Visual inspections performed according to item 1	-	<input type="checkbox"/>
TSC	8.2	Protective earth resistance, maximum 0.3 Ω (with power cable)	_____ Ω	<input type="checkbox"/>
TSC	8.3	Device leakage current measured <input type="checkbox"/> Device leakage current – differential current measurement according to fig. 5 or <input type="checkbox"/> Device leakage current – direct measurement according to fig. 4 Nominal voltage of power supply: _____ V Device leakage current, mains polarity 1 For line voltage _____ μA _____ V Normalized to nominal voltage (maximum 500 μA, see Additional conditions) _____ μA Device leakage current, mains polarity 2 For line voltage _____ V Normalized to nominal voltage (maximum 500 μA, see Additional conditions) _____ μA		<input type="checkbox"/>
	9	Functional checks		
TSC	9.1	T1 Test completed	-	<input type="checkbox"/>

Test equipment used:

Temperature, conductivity, pressure
(type, serial number): _____

Protective earth resistance, leakage current
(type, serial number): _____

Remarks:

Date:	Signature:	Stamp:
-------	------------	--------

The device is released for its intended use (attach inspection label)	<input type="checkbox"/> Yes <input type="checkbox"/> No	
Date of next check:		
Remarks:		
Date:	Signature:	Stamp:



Warning

A disinfection must be performed on completion of the checks.