INSTRUCTIONS FOR USE

Axium™ Prime Detachable Coil and Instant Detacher (I.D.)

fr	Spirale détachable Axium™ Prime et système de détachement instantané (I.D.)	sk	Oddeliteľná cievka Axium™ Prime a rýchly oddeľovač
de	Ablösbarer Axium™ Prime Coil und Schnellablösesystem (I.D.)	ro	Spirala detașabilă Axium™ Prime și dispozitivul de detașare instantanee (D.I.)
it	Spirale staccabile Axium™ Prime e dispositivo di distacco istantaneo	bg	Отделяща се намотка Axium™ Prime и устройство за незабавно освобождаване (I.D.)
es	Bobina desmontable Axium™ Prime y dispositivo de desprendimiento instantáneo	et	Eraldatav koil Axium™ Prime ja kiireraldaja
sv	Axium™ Prime frånkopplingsbar coil och snabbfrånkopplare (I.D.)	hr	Odvojiva zavojnica i brzi odvajač Axium™ Prime
nl	Loslaatbare Axium™ Prime-coil en loslaatsysteem	lt	„Axium™ Prime“ atkabinamoji spirale ir skubaus atkabinimo priemonė (ID)
pt	Espiral destacável e dispositivo de destacamento instantâneo (DDI) Axium™ Prime	lv	Axium™ Prime atvienojamā spirāle un tūlītējās atvienošanas ierīce
fi	Irrotettava Axium™ Prime -koili ja pikairrotin	sl	Snemljiva tuljava Axium™ Prime in pripomoček za hitro snemanje (I.D.)
da	Axium™ Prime aftagelig coil og lynfrakobler (I.D.)	sr	Axium™ Prime odvojiva spirala i sredstvo za trenutno odvajanje (T.O.)
el	Αποσπώμενο σπείραμα Axium™ Prime και εργαλείο άμεσης απόσπασης (I.D.)	pt-br	Mola destacável Axium™ Prime e separador instantâneo (S.I.)
cs	Odpojiteľná spirála Axium™ Prime a rychloodpojovač (I.D.)	ko	Axium™ Prime 분리형 코일 및 인스턴트 분리기(I.D.)
hu	Axium™ Prime leválasztható spirál és gyorsleválasztó (I.D.)	vi	Cuộn coil nút mạch Axium™ Prime và Instant Detacher (I.D.)
ru	Отсоединяемая спираль Axium™ Prime и система моментального отсоединения I.D. (Instant Detacher)	kk	Axium™ Prime ажырайтын шарғысы және жылдам ажыратқыш (Ж.А.)
pl	Odłączalna spirala Axium™ Prime i przyrząd do szybkiego odłączania Instant Detacher (I.D.)	mk	Одделив намот Axium™ Prime и инстант одделувач (I.D.)
tr	Axium™ Prime Brakılabilir Koil ve ID Anında Koparıcı	id	Koil Lepas Pasang Axium™ Prime dan Pelepas Instan (ID)
no	Axium™ Prime-coil med frigjøringsmekanisme og hurtigfrigjøringsenhet (Instant Detacher, I.D.)	ar	الملف القابل للفصل للفتل Axium™ Prime وجهاز الفصل الفوري (I.D.)

TABLE OF CONTENTS

Axiom™ Prime Detachable Coil and Instant Detacher (I.D.)

English	3
Français	4
Deutsch	6
Italiano	8
Español	10
Svenska	12
Nederlands	14
Português	16
Suomi	18
Dansk	20
Ελληνικά	22
Česky	24
Magyar	26
Русский	28
Polski	30
Türkçe	32
Norsk	34
Slovenčina	36
Română	38
Български	40
Eesti	42
Hrvatski	44
Lietuvių	45
Latviski	47
Slovenščina	49
Srpski	51
Português (Brasil)	53
한국어	55
Tiếng Việt	57
Қазақша	59
Македонски	61
Bahasa Indonesia	63
العربية	67
Symbol Glossary	68

English en

Instructions for Use

Axium™ Prime Detachable Coil and Instant Detacher (I.D.)

CAUTION

- Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.
- This device should be used only by physicians with a thorough understanding of angiography and percutaneous neurointerventional procedures.

DESCRIPTION

The Axium™ Prime Detachable Coil consists of a platinum embolization coil attached to a composite implant delivery pusher with a radiopaque positioning marker and a hand-held Instant Detacher (I.D.) which when activated detaches the coil from the delivery pusher tip. The Instant Detacher (I.D.) is sold separately.

DEVICE COMPATIBILITY

- The following devices are required for use with the Axium™ Prime Detachable Coil:
- Axium™ Prime Detachable Coils should only be delivered through a microcatheter with a minimum inside diameter of 0.0165"–0.017" with two marker bands.
 - Instant Detacher (I.D.).

Other Accessory Products (required to perform a procedure)

- 6–8F Guide Catheter*
- Microcatheter (see above)*
- Guidewires compatible with microcatheter*
- Continuous saline/heparin saline flush set*
- Rotating hemostatic valves (RHV)*
- 3-Way stopcock*
- 1-Way stopcock*
- IV pole*
- Arterial Sheath*

*Not provided as part of the system; chosen based upon physician experience and preference.

INTENDED PURPOSE / INDICATIONS FOR USE

The Axium™ Prime Detachable Coils are intended for the endovascular embolization of intracranial aneurysms. The Axium™ Prime Detachable Coils are also intended for the embolization of other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Access site complications such as vessel occlusion, intimal damage, pseudoaneurysm, arteriovenous fistula, necrosis, pain, infection, inflammation, granuloma, hematoma, nerve damage, or pathological hand cold intolerance
- Delayed neoplasia
- Hemorrhage
- Ischemia
- Neurological deficits including hand dysfunction, stroke, and death
- Thromboembolic episodes
- Complications of radiation exposure such as alopecia, burns ranging in severity from skin reddening to ulcers, cataracts,
- Vascular thrombosis
- Vasospasm
- Vessel perforation

WARNINGS

- The Axium™ Prime Detachable Coil, the dispenser track, and the introducer sheath are supplied in a sterile and non-pyrogenic, unopened and undamaged package. The package should be checked for potential damage. Damaged Axium™ Prime Detachable Coils must not be used, as they may result in patient injury.
- Axium™ Prime Detachable Coils are intended for one use only. The Instant Detacher (I.D.) is supplied sterile and intended for single patient use. After use do not resterilize and/or reuse. Reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn, may result in patient injury, illness or death.
- Do not use if sterile packaging has been compromised or damaged.
- Damaged implant delivery pusher and/or coils may affect coil delivery to, and stability inside, the vessel or aneurysm, possibly resulting in coil migration or stretching.
- Do not rotate the implant delivery pusher during or after delivery of the coil into the aneurysm. Rotating the delivery pusher during or after coil delivery into the aneurysm may result in a stretched coil or premature detachment of the coil from the implant delivery pusher, which could result in coil migration.
- Do not use hemostats or a torquer in an attempt to advance delivery pusher. This may result in a kinked/damaged pusher which may lead to premature detachment.

WARNINGS

- Verify that the distal shaft of the microcatheter is not under stress before the Axium™ Prime Detachable Coil detachment. Axial compression or tensile forces could be stored in the microcatheter causing the tip to move during the Axium™ Prime Detachable Coil delivery. Microcatheter tip movement could cause the aneurysm or vessel to rupture.
- Advancing the delivery pusher beyond the microcatheter tip once the coil has been deployed and detached involves risk of aneurysm or vessel perforation.
- If undesirable movement of the Axium™ Prime Detachable Coil can be seen under fluoroscopy following coil placement and prior to detachment, remove the coil and replace with another more appropriately sized Axium™ Prime Detachable Coil. Movement of the coil may indicate the coil could migrate once it is detached. Angiographic controls should also be performed prior to detachment to ensure that the coil mass is not protruding into the parent vessel.
- High quality, digital subtraction fluoroscopic road mapping is mandatory to achieve safe catheterization of the aneurysm or vessel and correct placement of the first coil. With smaller aneurysms this is a particularly important step.
- If Axium™ Prime Detachable Coil repositioning is necessary, take special care to retract coil under fluoroscopy in a one-to-one motion with the implant pusher. If the coil does not move with a one-to-one motion, or repositioning is difficult, the coil has been stretched and could possibly break. Gently remove and discard both the catheter and coil.
- Due to the delicate nature of the Axium™ Prime Detachable Coil, the tortuous vascular pathways that lead to certain aneurysms and vessels, and the varying morphologies of intracranial aneurysms, a coil may occasionally stretch while being maneuvered. Stretching is a precursor to potential malfunctions such as coil breakage and migration.
- If resistance is encountered while withdrawing an Axium™ Prime Detachable Coil, which is at an acute angle relative to the catheter tip, it is possible to avoid coil stretching or breaking by carefully repositioning the distal tip of the catheter at the ostium of the aneurysm, or just slightly inside the parent artery.
- Take care not to puncture gloves or sterile drape while handling implant delivery pusher.
- Multiple placements of Axium™ Prime Detachable Coils may be required to achieve the desired occlusion of some aneurysms or vessels.
- The long term effect of this product on extravascular tissues has not been established so care should be taken to retain this device in the intravascular space.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.
- The safety and effectiveness of this device for radial neurovasculature access in direct comparison to a transfemoral approach has not been demonstrated. The risks and benefits for radial access against a transfemoral approach should be carefully weighed and considered for each patient.

PRECAUTIONS

- Handle the Axium™ Prime Detachable Coil with care to avoid damage before or during treatment.
- Do not advance the Axium™ Prime Detachable Coil against a noted resistance until the cause of the resistance is cleared by fluoroscopy. This may lead to the destruction of the coil and/or catheter or perforation of the vessel.
- It is essential to confirm catheter compatibility with the Axium™ Prime Detachable Coil. The outer diameter of the Axium™ Prime Detachable Coil should be checked to ensure that the coil will not block the catheter.
- The dispenser track, introducer sheath and the Instant Detacher (I.D.) are not intended to make contact within the patient.
- Do not use the Axium™ Prime Detachable Coil and the Instant Detacher (I.D.) after the expiration date printed on the product label.
- In order to achieve optimal performance of the Axium™ Prime Detachable Coil and to reduce the risk of thromboembolic complications, it is critical that a continuous infusion of appropriate flush solution be maintained.
- Advance and retract Axium™ Prime Detachable Coils slowly and smoothly, especially in tortuous anatomy. Remove the coil if unusual friction or "scratching" is noted. If friction is noted in a second coil, carefully examine both the coil and the catheter for possible damage such as catheter shaft buckling or kinking, or an improperly fused joint.
- If delivery pusher buckling or kinking is noted, grasp the distal most portion of the delivery pusher, distal to the kink, buckle or break and remove from microcatheter.
- Do not advance the coil with force if the coil becomes lodged within or outside the microcatheter. Determine the cause of resistance and remove the system when necessary.
- If resistance is encountered when withdrawing the implant delivery pusher, draw back on the infusion catheter simultaneously until the delivery pusher can be removed without resistance.
- If resistance is noted during coil delivery, remove the system and check for possible damage to the catheter.
- If using radial artery access, perform a screening examination of the radial artery per institutional practices to ensure that radial access is appropriate for the patient.
- Operators should take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.

- Dispose of device in accordance with hospital, administrative, and/or local government policy.

PREPARATIONS FOR USE

- In order to achieve optimal performance of the Axium™ Prime Detachable Coil and to reduce the risk of thromboembolic complication, it is advised that a continuous saline flush be maintained between
 - the arterial sheath and the guiding catheter,
 - the microcatheter and guiding catheter and
 - the microcatheter and the implant delivery pusher and the Axium™ Prime Detachable Coil.
- Place the appropriate guiding catheter following recommended procedures. Connect a rotating haemostatic valve (RHV) to the hub of the guiding catheter. Attach a 3-way stopcock to the side arm of the RHV, then connect a line for the continuous flush.
- Attach a second RHV to the hub of the microcatheter. Attach a 1-way stopcock to the side arm of the RHV, then connect a line for continuous flush.

For Axium™ Prime Detachable Coils:
One drop from the pressure bag every 3–5 seconds is suggested
- Check all fittings so that air is not introduced into guiding catheter or microcatheter during continuous flush.

MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

- MR** Non-clinical testing and electromagnetic/thermal modeling demonstrated that the Axium™ Prime Detachable Coil is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
- Static magnetic field of 1.5-Tesla or 3-Tesla only.
 - Maximum spatial gradient magnetic field of 4000 gauss/cm (40T/m) or less.
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg during operation for the MR system.

Under the scan conditions defined above, the Axium™ Prime Detachable Coil is expected to produce a maximum temperature rise of less than or equal to 1.9°C after 15-minutes of continuous scanning. The Axium™ Prime Detachable Coil can be scanned for 60 minutes of continuous scan time. Allow 30 minutes of cooling before additional scans.

In non-clinical testing, the image artifact caused by the device extends approximately 5 mm from the Axium™ Prime Detachable Coil when imaged with a gradient echo pulse sequence and a 3-Tesla MRI system.

In non-clinical testing, the magnetically induced displacement force and magnetically induced torque were tested and zero displacement and zero torque was measured.

DIRECTIONS FOR USE

- Slowly and simultaneously remove the Axium™ Prime Detachable Coil and introducer sheath from the dispenser track. Inspect proximal implant delivery pusher for irregularities. If irregularities exist, replace with a new Axium™ Prime Detachable Coil.
- Slowly advance the Axium™ Prime Detachable Coil out of the introducer sheath into the palm of your gloved hand and inspect for irregularities of the coil or the detachment zone. Due to potential risks of irregularities, a visual proof should be performed. If irregularities exist, replace with a new Axium™ Prime Detachable Coil (See Figure 1).


Figure 1

- Implant
 - Detachment Zone
 - Coil Alignment Marker
 - Delivery Pusher
 - Hypotube Break Indicator (HBI)
 - Positive Load Indicator
- Gently immerse the Axium™ Prime Detachable Coil and its detachment zone in heparinized saline. Take care not to stretch the coil during this procedure, in order to preserve the coil memory. While still immersed in the heparinized saline, point introducer sheath vertically into saline and gently retract the distal tip of the coil into the introducer sheath.
 - Insert the distal end of the introducer sheath through the rotating hemostatic valve (RHV) and into the hub of the microcatheter until the sheath is firmly seated. Tighten the RHV around the introducer sheath to prevent back flow of blood, but not so tight as to damage the coil during its introduction into the catheter.
 - Transfer the Axium™ Prime Detachable Coil into the microcatheter by advancing the implant pusher in a smooth, continuous manner (1–2 cm strokes). Stop approximately 15 cm from the proximal end of implant pusher. Do not attempt to advance the entire implant pusher into the introducer sheath this may cause kinking in the proximal pusher. Once the flexible portion of the implant pusher has entered the catheter shaft, loosen the RHV and remove the introducer sheath over the implant pusher's proximal end. Once completed, tighten the RHV around the implant pusher. Leaving the introducer sheath in place will interrupt normal infusion of flushing solution and allow back flow of blood into the microcatheter.
 - Visually verify that the flushing solution is infusing normally. Once confirmed, loosen the RHV enough to advance the implant pusher, but not enough to allow back flow of blood into the implant pusher.
 - Advance the Axium™ Prime Detachable Coil under fluoroscopy and position carefully at the desired site. If coil placement is unsatisfactory, slowly withdraw by pulling on the implant pusher, then slowly advance again to reposition the coil. If the coil size is inappropriate, remove and replace with an appropriately sized coil.

1. If delivery pusher buckling or kinking is noted, grasp the distal most portion of the delivery pusher, distal to the kink, buckle or break and remove from microcatheter.

WARNING

Do not use hemostats or a torquer in an attempt to advance delivery pusher. This may result in a kinked/damaged pusher which may lead to premature detachment.

8. Continue to advance the Axiom™ Prime Detachable Coil until the coil alignment marker of the implant pusher is just distal of the proximal marker of the microcatheter (See Figure 2).

1. Advance coil alignment marker just beyond proximal catheter marker band, then retract delivery pusher until coil alignment marker creates a "T" with the proximal catheter marker band under fluoroscopy. This relieves forward stress that may lead to false positive detachment.


Figure 2

1. Coil Alignment Marker
 2. Proximal Catheter Marker Band
 3. Distal Catheter Marker Band
9. Tighten RHV to prevent movement of the implant pusher.
 10. Remove the Instant Detacher (I.D.) from its protective packaging and place it within the sterile field. Instant Detacher (I.D.) is packaged separately as a sterile device for single patient use only.
 11. Confirm again under fluoroscopy that the coil alignment marker of the implant pusher creates a "T" with the proximal marker of the microcatheter.
 12. Verify that the RHV is firmly locked around the implant pusher before attaching the Instant Detacher (I.D.) to ensure that the coil does not move during the connection process. Ensure that the implant pusher is straight between the RHV and the Instant Detacher (I.D.). Straightening this section of the implant pusher optimizes alignment to the Instant Detacher (I.D.).
 13. Hold the proximal end of the implant pusher at the distal end of the load indicator. Advance the Instant Detacher (I.D.) over the proximal end of the implant pusher until the load indicator fully enters the funnel and the pusher is firmly seated in the actuator (See Figure 3).


Figure 3

1. Proximal end of Load Indicator
2. Grip pusher here when advancing Instant Detacher (I.D.) over implant pusher.
3. Advance Instant Detacher (I.D.) over implant pusher.
4. Ready to Detach
5. RHV and Implant pusher are in a straight line and distal end of Load Indicator is fully inserted in the Instant Detacher (I.D.) funnel.

Note: If indicator band is still showing as in Figure 4 below, the Instant Detacher (I.D.) should be advanced further until the implant pusher is fully seated in the funnel as in Figure 3 above.


Figure 4

1. Not Fully Loaded

14. To detach coil, place the Instant Detacher (I.D.) into palm and retract the Thumb-slide back until it stops and clicks, and slowly allow the Thumb-slide to return to its original position. Remove the Instant Detacher (I.D.).

Note: The Instant Detacher (I.D.) can also be removed at the end of stroke if desired. To remove the Instant Detacher (I.D.) at the end of stroke, hold the thumb slide in its most rearward position and remove the Instant Detacher (I.D.) (See Figure 5).


Figure 5

1. Ready to Detach
 2. Distal end of Load indicator fully in Instant Detacher (I.D.) funnel.
15. Successful coil detachment must be verified by fluoroscopic monitoring to ensure that the coil has detached. Slowly pull back the implant pusher while watching the fluoroscopy to make sure that the coil does not move. In the unlikely event the coil moves, repeat steps 12–14. If necessary, advance the implant pusher to re-establish the coil and catheter marker alignment. Verify coil detachment as above.
 16. If you wish to confirm detachment, grasp the positive load indicator between your thumb and forefinger of your left hand and the proximal end of the implant delivery pusher with your thumb and forefinger of your right hand. Gently pull on the proximal end of the implant delivery pusher. If it moves freely from the hypo-tube, the system is detached properly. If it does not, repeat steps 13–15.
Note: If the Coil does not detach after 3 attempts, discard the Instant Detacher (I.D.) and replace with a new Instant Detacher (I.D.).
 17. In the rare event that the coil does not detach and cannot be removed from the implant delivery pusher, use the following steps for detachment.
 1. Grip the hypotube approximately 5 cm distal of the positive load indicator at the hypotube break indicator and bend the implant delivery pusher just distal to the HBI 180 degrees.
 2. Next straighten the pusher back, continue bending and straightening until the pusher tubing opens exposing the release element (Figure 6).


Figure 6

3. Gently separate the proximal and distal ends of the open pusher. Then, under fluoroscopy, pull the proximal portion of the implant delivery pusher approximately 2–3 cm to confirm implant detachment per IFU (Figure 7).


Figure 7

18. Once coil detachment has been detected and fluoroscopically confirmed slowly withdraw the implant pusher from the microcatheter.

WARNING

- If attempt to detach fails, remove coil from treatment area and microcatheter and replace with a new Axiom™ Prime Detachable Coil under fluoroscopy.
- If coil becomes prematurely detached, remove implant pusher and confirm under fluoroscopy:
 - Advance next coil to push remaining tail of prematurely detached coil into treatment area.
 - Remove prematurely detached coil with the appropriate retrieval device.

19. Repeat steps 1–18, if additional coil placements are required.
20. Once the procedure is complete, discard the Instant Detacher (I.D.).

WARNING

- Do not resterilize the Instant Detacher (I.D.). For Single Patient use only.
- The Instant Detacher (I.D.) is intended for a maximum of 25 cycles.

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Micro Therapeutics, Inc.
Medtronic B.V.