INgezim® CSF COMPAC

Prod Ref: 11.CSF.K3

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials with Positive Control	1	1mL	2	1mL
Vials with Negative Control	1	1mL	2	1mL
Vials with Peroxidase Conjugate, ready to use	1	30mL	2	30mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Diluent (DE01-01), ready to use	1	125mL	1	125mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter)

TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (**Blocking Elisa**). We make a brief description of the technique below: The antigen is fixed in a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen absorbed on plate while if the serum sample does not contain specific antibodies, they will not bind the antigen. If we add a specific monoclonal antibody (Mab) against the viral antigen coated to the plate (conjugated with peroxidase), it will compete with the antibodies of the serum. If the serum samples contain specific antibodies, they will not permit the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material from the plate, we can detect the presence or absence of labelled Mab by adding the specific substrate that in presence of the peroxidase will develop a colorimetric reaction. The antigen coated to the plate in our kit consists of a purified recombinant protein from the virus (the E2), which is the major structural protein from the CSFV and the most antigenic one.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

Control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Samples to be used in the assay must be porcine sera. **Do not use highly hemolyzed or contaminated samples. You might obtain false positive results.** To be assayed a 1/2 dilution must be done in the supplied diluent.

This dilution could be done directly on the plate by adding 50 µL of diluent and 50 µL of sera into each well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e., 40 mL of concentrate plus 960 mL of dH2O). When ready this solution remains stable between +2°C and +8°C until the expire date described at the label of the concentrated solution.

Control sera:

Once opened, control sera are stable for one month between +2°C and +8°C. In case that they are not going to be used in this period, we recommend to store them under -20°C.

The control sera need to be diluted 1/2 in diluent prior to be used, in the same way than the sera samples (you can do the dilution directly into the plate by adding 50 µL of diluent plus 50 µL of control sera).

Conjugate: The conjugate is ready to use.

Bring to room temperature and add 100 µL/well.

TEST PROCEDURE

1. All reagents must be allowed to come to room temperature (20°C to 25°C) before use.
2. Add 50 µL of supplied diluent to each well. Add 50 µL of positive control sera to two wells (e.g., A1 and B1), and 50 µL of the negative control sera (e.g., A2 and B2). Add 50 µL of sera samples to test on each remainder wells. We recommend the use of two wells for each control. Seal the plate and **incubate for 90 ± 5 min at 36 ± 1°C**.
3. Wash 5 times as previously described.
4. Add 100 µL of specific conjugate to each well. Seal the plate and **incubate for 15 minutes at 36 ± 1°C**.
5. Wash 5 times as previously described
6. Add 100 µL of substrate to each well. Keep the plate for **10 min at room temperature**.
7. Add 100 µL of stop solution, to each well.
8. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

VALIDATION OF THE TEST:

The test could be considered valid when the OD of the Negative Control (NC) is, at least ,4 times higher than the OD of the Positive Control (PC):

$$\frac{\text{Negative control}}{\text{Positive control}} \geq 4$$

RESULTS INTERPRETATION:

When you are running duplicates, OD of the sample, will be calculated as the arithmetic mean of OD values in both wells.

Calculate the blocking % of each sample as follow:

$$\% = 100 - \left(\frac{\text{OD Sample}}{\text{OD Negativo control}} \right) \times 100$$

- All samples with blocking % higher or equal than 35 must be considered Positive.
- All samples with blocking % lower than 35 must be considered Negative.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

9191.INGE
9175.INGE

Distribuido en
Distributed in

por
by