INgezim® ADV gB

Prod Ref: 11.GB. K3

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KIT COMPOSITION

REAGENT	2 Plates		5 Plates (divisible)		5 Plates (indivisible)		10 Plates	
	Uni.	Vol.	Uni.	Vol.	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5		5		10	
Vials with Positive Control	1	2mL	2	2mL	1	2mL	1	2mL
Vials with Negative Control	1	2mL	2	2mL	1	2mL	1	2mL
Vials with Peroxidase Conjugate, ready to use	1	30mL	2	30mL	2	30mL	1	110mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL	1	125mL	1	250mL
Bottles with Diluent (DE31-01), ready to use	1	15mL	1	15mL	1	15mL	1	30mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL	1	60mL	1	120mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL	1	60mL	1	120mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter)

TECHNICAL BASIS

The gB test uses a monoclonal antibody (Mab) specific for gB glycoprotein of ADV and microtiter plates coated with ADV antigen. After sample addition, if antibodies are absent in the test serum, the gB antigen will remain free and the Mab peroxidase-conjugated specific for gB, will bind. After substrate addition a colour reaction is produced. On the other hand, if gB antibodies are present in the test sample they will bind to the gB antigen of the well and will block the binding of the conjugated Mab. In this case no colour reaction will be produced.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.
Control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Serum samples and controls must be used without previous dilution adding 80 µL/well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

Conjugate:

The conjugate is ready to use.

Controls (+) y (-):

Control sera should be used like samples (80 µL/well).

TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add 20 µL of serum diluent to each well and then, 80 µL of serum (controls and samples) to each well. Gently, mix the plate, seal it and incubate for **1 hour at 37°C or 16-18 hours at 4°C**.
3. Wash 3 times following the described in Information on Washing Steps
4. Add 100 µL of conjugate to each well. Seal the plate and incubate for **30 minutes at 37°C**
5. Wash 5 times following the procedure.
6. Add 100 µL of substrate, to each well. Keep the plate for **10 min at room temperature in a dark place**.
7. Add 100 µL of stop solution to each well.
8. Read the OD at 450 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

TEST VALIDATION:

The test must be considered valid if:

- Positive control OD is Lower than 0.3
- Negative control OD is higher than 0.8

INHIBITION PERCENTAGE (IP) VALUES CALCULATION:

$$IP = 1 - \left(\frac{M}{CN} \right) \times 100$$

Where:

- M = Mean of OD 450 nm of the sample.
- CN = Mean of OD 450 nm of Negative control.

RESULT INTERPRETATION:

- IP higher than 30%: **Positive** samples.
- IP lower than 30%: **Negative** samples.

PRACTICAL EXAMPLE:

- OD 450 CN = 1.360
- OD 450 CP = 0.110
- OD 450 Sample 1 = 0.373
- OD 450 Sample 2 = 1.295

Validation:

- CN = 1.360 (> 0.8)
- CP = 0.110 (< 0.3)
- IP sample 1 = $[1 - (0.373/1.360)] \times 100 = 72.5$
- IP sample 2 = $[1 - (1.295/1.360)] \times 100 = 4.7$

Interpretation:

- Sample 1: Positive to gB antibodies (the animal has been in touch with a vaccinal and/or field ADV strain).
- Sample 2: Negative to gB antibodies (the animal has never been in touch with any of ADV strains. This result will be obtained with non-infected and non vaccinated pigs).

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

Distribuido en
Distributed in

por
by