INSTRUCTIONS FOR USE

Onyx™ Liquid Embolic System

This page is intentionally left blank.

TABLE OF CONTENTS

English.....	5
Français.....	6
Deutsch.....	7
Italiano.....	10
Español.....	12
Svenska.....	14
Nederlands.....	16
Português.....	18
Suomi.....	20
Dansk.....	21
Ελληνικά.....	23
Česky.....	25
Magyar.....	27
Русский.....	29
Polski.....	31
Türkçe.....	33
Norsk.....	34
Slovenčina.....	36
Română.....	38
Български.....	39
한국어.....	41
عربي.....	45
Symbol Glossary.....	46

This page is intentionally left blank.

English

Instructions for Use

EN

Onyx™ Liquid Embolic System

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

Description

Onyx™ LES is a non-adhesive liquid embolic agent comprised of EVOH (ethylene vinyl alcohol) copolymer dissolved in DMSO (dimethyl sulfoxide), and suspended micronized tantalum powder to provide contrast for visualization under fluoroscopy. The Onyx Liquid Embolic System (LES™) consists of a 1.5 ml vial of Onyx™ LES, a 1.5 ml vial of DMSO, two 1 ml Onyx™ LES delivery syringes and one 1 ml DMSO syringe. A DMSO compatible delivery micro catheter that is indicated for use in the neuro vasculature is used to access the embolization site. The Onyx™ LES Syringe-Catheter Interface adapter, which is sold separately, is an accessory to the Onyx™ LES system that provides an interface between the Onyx™ LES delivery syringe and the delivery catheter.

Principle Of Operation

Onyx™ LES is delivered through a micro catheter into the AVM under fluoroscopic control. The DMSO solvent dissipates into the blood and interstitial fluids, causing the EVOH copolymer and suspended tantalum to precipitate in situ into a spongy, coherent embolus. Onyx™ LES immediately forms a skin as the polymeric embolus solidifies from the outside to the inside, while traveling more distally in the vessel. Since Onyx™ LES is non-adhesive, the micro catheter can be left in place while slow, controlled injections are performed. Post embolization angiography can be conducted with the delivery micro catheter in place, enabling the physician to make additional injections through the same micro catheter, if necessary.

Indications For Use

Embolization of lesions in the peripheral and neurovasculature, including arteriovenous malformations and hypervascular tumors.

How Supplied

Onyx™ LES is available in three product formulations, Onyx™-18 LES (6% EVOH), Onyx™-20 LES (6.5% EVOH) and Onyx™-34 LES (8% EVOH):

- Onyx™-18 LES and Onyx™-20 LES: Recommended when feeding pedicle injections will be conducted close to the nidus.
- Onyx™-34 LES: Recommended for embolizing higher flow and larger fistulous components.

Contraindications

Onyx™ LES is not indicated for use with premature infants (<1,500 g) or individuals with significant liver function impairment.

Potential Complications

Potential complications include, but are not limited to:

- Hematoma.
- Arterial thrombosis.
- Ischemic events due to embolic migration, vasospasm, thrombosis.
- Catheter entrapment.
- Catheter rupture.
- Device migration and cast movement.
- Hemorrhagic accidents: vascular rupture – perforation.
- Hemodynamic changes induced by the embolization may result in hemorrhagic complications.
- Hemorrhagic complications related to attempts to remove entrapped catheter.
- These ischemic or hemorrhagic complications may result in various functional neurological deficits and possibly death.

WARNINGS

Due to the possibility of electrical arcing with the tantalum metal in the Onyx™ LES material, use of monopolar electrocautery devices for surgical resection of DAVMs or arteriovenous fistula embolized with Onyx™ LES should be avoided. Bipolar devices should be used with caution.

This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

- Use only the Onyx™ LES syringes to inject DMSO and Onyx™ LES. Other syringes may not be compatible with DMSO.
 - Use only DMSO compatible micro catheters. Other micro catheters may not be compatible with DMSO and their use can result in thromboembolic events due to catheter degradation.
 - Always consider the potential for Onyx™ LES to interact with other embolic agents, e.g., cyanoacrylates, coated coils, particles and/ or embolic spheres.
 - Performing embolization to occlude blood vessels is a high-risk procedure. The procedure should be carried out by a specialist with the appropriate neuro or peripheral interventional training, and a thorough knowledge of the pathology to be treated, angiographic techniques, and super-selective embolization.
 - The patient's anigoanatomy must be amenable to microcatheter tip placement at a location that is distal to arterial branches that potentially may feed a cranial nerve.
 - Adjunctive coil use should be considered if angiography shows that venous drainage of the AVM appears almost simultaneously with arterial opacification.
 - Premature solidification of Onyx™ LES may occur if micro catheter luer contacts saline, blood or contrast of any amount.
 - Failure to continuously mix Onyx™ LES for the required time may result in inadequate suspension of the tantalum, resulting in inadequate fluoroscopic visualization during delivery. Inject Onyx™ LES immediately after mixing. If Onyx™ LES injection is delayed, tantalum settling can occur within the syringe resulting in poor visualization of Onyx™ LES during injection.
 - Using the Onyx™ LES Syringe-Catheter Interface Adapter will reduce labeled dead space of the micro catheter. Failure to make the appropriate allowances may result in unintended embolization.
 - Reference the Onyx Syringe-Catheter Interface Adapter Instructions for Use for a list of compatible micro catheters and dead space information.
- When injecting Onyx™ LES, fluoroscopic visualization should reveal Onyx™ LES traveling through the catheter lumen. It is recommended to obtain fluoroscopic imaging prior to reaching the minimum catheter + adapter deadspace in order to visualize the embolic material before it exits the tip of the catheter.

- Only use thumb pressure to inject Onyx™ LES. Using palm of hand to advance plunger may result in catheter rupture due to overpressurization in the event of catheter occlusion.
- Inject Onyx™ LES and DMSO at a slow, steady rate, not to exceed 0.3 ml/minute. Animal studies have shown that a rapid injection of DMSO into the vasculature may lead to vasospasm and/or angioneurosis.
- Do not allow more than 1 cm of Onyx™ LES to reflux back over catheter tip. Angioarchitecture, vasospasm, excessive Onyx reflux, or prolonged injection time may result in difficult catheter removal and catheter entrapment. Excessive force to remove an entrapped catheter may cause serious intracranial hemorrhage. The long term effects of an entrapped catheter that is left in a patient is unknown, but could potentially include clot formation, infection or catheter migration.
- DO NOT interrupt Onyx™ LES injection for longer than two minutes prior to re-injection. Solidification of Onyx™ LES may occur at the catheter tip resulting in catheter occlusion, and use of excessive pressure to clear the catheter may result in catheter rupture.
- Adequate fluoroscopic visualization must be maintained during Onyx™ LES delivery or non-target vessel embolization may result. If visualization is lost at any time during the embolization procedure, HALT Onyx™ LES delivery until adequate visualization is re-established. Testing has shown that over-pressurization and rupture can occur if 0.05 ml of Onyx™ LES is injected and is not visualized exiting the catheter tip.

- Increased resistance to Onyx™ LES injection – Stop Injection. Do not attempt to clear or overcome resistance by applying increased injection pressure. If this occurs, determine the cause of resistance (for example, Onyx™ LES occlusion in catheter lumen), and replace catheter. Use of excessive pressure may result in catheter rupture and embolization of unintended areas.
- After using a micro catheter with Onyx™ LES, do not attempt to clear the micro catheter or to inject any material through it. Attempts to clear catheter may lead to embolus or embolization of unintended area.
- If Onyx™ LES escapes outside the vascular space, as might occur if the vessel wall is compromised, a subacute inflammatory response to the material may occur.

PRECAUTIONS

Inspect product packaging prior to use. Do not use if sterile barrier is opened or damaged.

Use prior to expiration date.

Review both the catheter and interface adapter Instructions for Use prior to using Onyx™ LES.

Verify that the catheters and accessories used in direct contact with the Onyx™ LES polymer are clean and compatible with the material and do not trigger precipitation or degrade with contact. Refer to the respective Warnings and Directions for Use sections.

Failure to wait a few seconds to retrieve micro catheter after Onyx™ LES injection may result in fragmentation of Onyx™ LES into non-target vessels.

Difficult catheter removal or catheter entrapment may be caused by one or more of the following factors:

- Long catheterization time
- Angio-architecture: very distal arteriovenous malformation fed by afferent, lengthened, small, or tortuous pedicles
- Vasospasm
- Reflux
- Injection time

To reduce the risk of catheter entrapment, carefully select catheter placement and manage reflux to minimize the factors listed above.

Onyx™ LES reflux along the distal tip of the micro catheter: Apart from the risk of ischemic complications due to unintended embolization, significant reflux may result in entrapment of the micro catheter causing difficult removal. The amount of reflux that can be accepted must always be compared to the angio-architecture of the malformation to minimize risk of unintended embolization or difficult catheter removal. In general, minimize the reflux to less than 1 cm along the distal tip of the micro catheter. All other factors may affect this limit.

- Should catheter removal become difficult, the following technique may assist in catheter retrieval:
 - Carefully pull the catheter to assess any resistance to removal.
 - If resistance is felt, remove any "slack" in the catheter.
 - Gently apply traction to the catheter (approximately 3-4 cm of stretch to the catheter).
 - Hold this traction for a few seconds and release. Assess traction on vasculature to minimize risk of hemorrhage.
 - This process can be repeated intermittently until catheter is retrieved.
 - When using micro catheters, do not apply more than 20 cm of traction to catheter, to minimize risk of catheter separation.
- For entrapped catheters:
 - Under some difficult clinical situations, it may be safer to leave a flow-directed catheter in the vascular system, rather than risk rupturing the malformation and, consequently a hemorrhage, by exercising too much traction on an entrapped catheter.
 - This is accomplished by stretching the catheter and cutting the shaft near the entry point of vascular access allowing the catheter to remain in the artery.

- If catheter breaks during removal, distal migration or coiling of the catheter may occur.
- Same day surgical resection should be considered to minimize risk of thrombosis.

Storage

Store Onyx™ LES and DMSO between -20° and 55°C. Prior to use, maintain product temperature between 19° and 24°C. If product freezes due to exposure to colder temperatures, thaw at room temperature before use.

Directions For Use

1. Shake Onyx™ LES at least 20 minutes on an Onyx™ LES mixer¹ at a setting of 8. Continue mixing until ready to inject Onyx™ LES per step 7.
2. Confirm micro catheter placement with injection of contrast agent per institutional procedure.
3. Flush contrast from micro catheter hub with 10 ml of saline. Leave the syringe connected.
4. Filling catheter deadspace: aspirate approximately 0.8 ml of sterile DMSO into an Onyx™ LES 1 ml syringe. Inject the DMSO into the delivery micro catheter in sufficient volume to fill the catheter deadspace. Refer to delivery catheter labeling for dead space volume.
5. When filling the syringe with Onyx™ LES, aspirate approximately 1.1 ml of Onyx™ LES so that the piston head is beyond the 1 ml graduation line on the syringe barrel.

NOTE

This extra volume of Onyx™ LES will be used to fill (prime) the interface adapter.

6. Attach the interface adapter to the syringe.
7. Inject Onyx™ LES through the interface adapter until the distal end of the piston head is even with the 1 ml graduation line on the syringe barrel.
8. Any excess Onyx™ LES on the tip of the interface adapter may be wiped off with a clean, dry particle free cloth.
9. Remove the DMSO syringe from the catheter, overfill and wash the luer hub with the balance of the DMSO.
10. Immediately connect the Onyx™ LES syringe firmly to the micro catheter hub, making sure there is no air in the hub during the connection.
11. Inject Onyx™ LES to displace DMSO. Based on clinical practice, it is recommended that Onyx™ LES be injected at a steady rate of 0.16 ml/minute (0.25 ml/90 second). Do not exceed 0.3 ml/minute.

WARNINGS

Failure to continuously mix Onyx™ LES for the required time (20 minutes) may result in inadequate fluoroscopic visualization during delivery.

Inject Onyx™ LES immediately after mixing. Tantalum settling can occur within the syringe resulting in poor visualization of Onyx™ LES during injection.

Use thumb pressure only to inject Onyx™ LES. Using palm of hand to advance plunger may result in catheter rupture due to over-pressurization in the event of catheter occlusion.

Stop injection if Onyx™ LES is not visualized exiting catheter tip. If the catheter becomes occluded, over-pressurization can occur. During Onyx™ LES injection, continuously verify that Onyx™ LES is exiting the catheter tip. Testing has shown that over-pressurization and rupture can occur if 0.05 ml of Onyx™ LES is injected and is not visualized exiting the catheter tip.

WARNINGS

Stop injection if increased syringe resistance is felt. Increased injection force may be due to catheter occlusion. Continued injection into an occluded catheter will result in catheter rupture.

Do not interrupt Onyx™ LES injection for longer than two minutes prior to re-injection. Solidification of Onyx™ LES may occur at the catheter tip resulting in catheter occlusion, and use of excessive pressure to clear the catheter, may result in catheter rupture.

Stop injection if abnormal dilation of the micro catheter diameter is observed. The increase in micro catheter diameter can signify a blocked catheter. The increase in micro catheter diameter can best be determined by double incidence fluoroscopy at low magnification¹ or with 'road mapping' procedure, examining the largest possible section of the micro catheter.

Adequate fluoroscopic visualization must be maintained during Onyx™ LES delivery or non-target vessel embolization may result. If visualization is lost at any time during the embolization procedure, HALT Onyx™ LES delivery until adequate visualization is re-established.

12. If a second Onyx™ LES syringe is required for an embolization, do not remove the adapter from the catheter. When the new Onyx™ LES syringe is ready, simply remove the empty syringe from the adapter's proximal end and connect the new syringe to the adapter.
13. Upon completion of Onyx™ LES injection, wait a few seconds, slightly aspirate syringe, and then gently pull the catheter to separate the Onyx™ LES cast.

The following is an optional catheter retrieval technique:

- Upon completion of Onyx™ LES injection, wait a few seconds and then slightly aspirate.
- Remove all slack from the catheter by putting a few centimeters of traction on the catheter to create a slight tension in the catheter system.
- Firmly hold the catheter and then pull it using a quick wrist snap motion (from left to right) 10 – 15 cm to remove the catheter from the Onyx™ LES cast.

NOTES

It is not necessary to pull the catheter any more than 20 cm.

This technique has not been qualified for the Rebar catheters.

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Micro Therapeutics, Inc.