When an executed macro entry in the *Log* window is edited in the comment column, the macro color change is reversed in the *Macro* window, and, the color change is visible on all the workstations that are joined to the study.

## 14.3 Macro toolbar

When the edit mode is on, icons to the right of the **Edit Mode** icon become active.

Figure 14-1 Macro toolbar

Item	Name	Description
1	Folder Up	Navigates to the parent folder. This button is disabled at the top-level macro folder.
2	Print	Prints all of the items in a selected folder or all of the items in all selected folders.
3	Edit Mode	Turns edit mode on and off.
4	Insert Folder	Creates a new folder.
5	Insert Item	Inserts a new macro item in a selected folder.
6	Edit Item/Folder	Edits item descriptions or folder descriptions.
7	Delete	Deletes an item in a list or deletes an entire folder and all the contents.
8	Cut	Removes an item from a folder to be placed in another folder or in a different location in the folder.
9	Copy	Copies an item in a folder to be placed in another folder.
10	Paste	Places an item that was cut or copied into another location.
11	Navigation path	Shows the ancestry for the current folder location.

## 14.4 Customize macros

### NOTE

Always build macros in the sequence they will be used. For example, do not build a macro command that says "Zero Transducer" before building the macro command to "Set up Transducer".

**Edit** mode must be on to modify any macro functions. Click the **Edit Mode** icon and the edit mode icons appear.

### 14.4.1 Create macro folders

1. To create a new macro folder, click the **Insert Folder** icon. A text box appears in the **Macro** column.

#### NOTE

To insert the new folder to the **Categories** directory, select the **Categories** folder.

To add the new macro folder in a specific place, highlight the folder below where the new macro folder should appear before clicking the **Insert Folder** icon. For example, to locate the new folder above the pre-case folder, highlight the pre-case folder and click the **New Folder** icon. The new folder will appear above the pre-case folder.

2. Type the name of the new macro folder in the text box. The new macro folder appears in the macro category window.

### 14.4.2 Delete macro folders

1. To delete a macro folder, click **Categories** from the *Macro* window.
2. Select the folder to be deleted and click the **Delete** icon.
3. Confirm the deletion of the folder.

### 14.4.3 Create a new macro

1. Select the name of the desired macro folder in the **Macro** list.
2. Click the **Insert Item** icon in the toolbar. The *Macro Builder* window appears.

Field	Description
Macro name	Enter the macro name or comment or note text. This text box has a 64 character maximum limit.
No Log entry	Check this box to prevent <b>Macro Name</b> text from being entered into the <i>Log</i> window. If this box is checked, the macro name will display in parenthesis, for example, <b>(Patient Arrives)</b> . This macro statement will not be entered into the <i>Log</i> window.
Execute next macro	Check this box to execute the next macro without additional clicks.
Macro commands	Folder list of available commands. Some folders have a plus sign (+) to the left which indicates there are sub-choices.
Command	Actual command actions that can be executed by the macro are displayed in the <i>Command</i> window.
Step	List of the actions performed by the macro.

3. Type the name of the macro in the **Macro name** text box.

#### NOTE

The **Macro name** text box has a 64 character maximum limit. It is possible to type more than 64 characters but only 64 characters will appear in the Macro name when displayed in the *Macro* window.

#### NOTE

When a macro is executed, the color change is visible for all the macros that have the same macro name. It is advised to give your macro a distinct name.

4. To create a simple text macro (for example, a standard comment or note):

CONFIDENTIAL

SECRET DE  
AFACERI

- a. Type the macro name and click **OK**.
  - b. Select or deselect the **No Log entry** and/or **Execute next macro** text boxes as desired.
  - c. Click **OK** to close the *Macro Builder* window. The new macro appears in the macro folder in the *Macro* window.
5. To create a more complex macro:
- a. Choose a parameter from the **Macro commands:** window (on the left). The available macro command displays in the *Command* window (on the right).
  - b. Double-click the macro command in the *Command* window to add it to the macro.

**NOTE**

Some of the parameters have a dialog that displays after the parameter is double-clicked that may require text (for example, **Lead On** prompts for the name of the Lead).

- c. Select or deselect the **No Log entry** and/or **Execute next macro** text boxes as desired (see [the Forms Macro table on page 181](#) for information about these text boxes).
- d. Complete steps a and b to add more macro commands to the macro, if desired.

**NOTE**

Use the **Delete**, **Cut**, **Copy** and **Paste** icons to remove and rearrange the macro steps in the **Step** list as desired.

- e. Click **OK** to close the *Macro Builder* window. The new macro appears in the macro folder in the *Macro* window.

**NOTE**

To add a standard or custom form macro, see [15.3.1 Forms macro on page 196](#).

## 14.5 Edit macros

**NOTE**

Macros should not be edited while studies are in progress anywhere on the Centricity Cardiology INW network. Editing macros while live studies are in progress elsewhere on the network may cause permanent loss of macro folders and/or commands.

1. Click the **Edit Mode** icon in the *Macro* window.
2. Select the macro folder containing the macro to edit.
3. Double-click the macro or click the **Edit Item/Folder** icon. The *Macro Builder* window appears.
4. See [14.4.3 Create a new macro on page 181](#) to make any necessary edits to the macro.
5. Click **OK** to close the *Macro Builder* window.
6. Click **Edit Mode** to exit the edit mode.

## 14.6 Import/export macros

For additional information about importing/exporting custom data, see [5.6 Import and export custom data on page 77](#).

CONFIDENTIAL

SECRET DE  
AFACERI

## 14.7 Macro commands

Channels	ECG	Real-Time		
				<b>Lead Off:</b> Turn off lead I, II, III, AVF, AVL, AVR, or V1 - V6.
				<b>Lead On:</b> Turn on lead I, II, III, AVF, AVL, AVR, or V1 - V6.
		<b>Review</b>		<b>Lead Off:</b> Turn off lead I, II, III, AVF, AVL, AVR, or V1 - V6.
				<b>Lead On:</b> Turn on lead I, II, III, AVF, AVL, AVR, or V1 - V6.
		<p><b>NOTE</b></p> <p>When setting up V leads, labels, such as V1, V2, must be used, not just V.</p>		
	<b>Pressure</b>	<b>Real-Time</b>	P1	<p><b>Change Label:</b> Set or change a specific pressure channel label.</p> <p><b>NOTE</b></p> <p>Labels must be uppercase letters (<b>AO</b>, not <b>ao</b>).</p>
			P1	<b>Mean Off:</b> Turn the mean pressure tracing off.
			P1	<b>Mean On:</b> Turn the mean pressure tracing on.
			P1	<b>Off:</b> Turn the pressure channel off.
			P1	<b>On:</b> Turn the pressure channel on.
			P1	<b>Phasic Off:</b> Turn the phasic tracing off.
			P1	<b>Phasic On:</b> Turn the phasic tracing on.
			P1	<b>Scale:</b> Set or change the pressure scale.
			P1	<b>Set Active:</b> Make the pressure channel the active channel.
			P1	<b>Zero:</b> Zero the individual pressure channel.
		<b>P2 - P4 Real-Time:</b> All the same commands and functions.		
			FFR	<b>Equalize:</b> Equalize P1 and P2 channels.
			FFR	<b>Off:</b> Turn off the FFR channel.
			FFR	<p><b>On:</b> Turn on the FFR channel.</p> <p>If the DFR channel is currently on, this changes to the FFR channel without requiring equalization to be done again.</p>
			DFR	<b>Equalize:</b> Equalize P1 and P2 channels.
			DFR	<b>Off:</b> Turn off the DFR channel.
			DFR	<p><b>On:</b> Turn on the DFR channel.</p> <p>If the FFR channel is currently on, this changes to the DFR channel without requiring equalization to be done again.</p>
			Zero All	<p>Zero all available pressure channels simultaneously.</p> <p><b>NOTE</b></p> <p>There will be a 15 second delay using this command versus selecting zero for each transducer in use.</p>

CONFIDENTIAL

SECRET DE AFACERI

		<b>Review</b>	P1	<b>Change Label:</b> Set or change a specific pressure channel label.
			P1	<b>Mean Off:</b> Turn the mean pressure tracing off.
			P1	<b>Mean On:</b> Turn the mean pressure tracing on.
			P1	<b>Off:</b> Turn the pressure channel off.
			P1	<b>On:</b> Turn the pressure channel on.
			P1	<b>Phasic Off:</b> Turn the phasic tracing off.
			P1	<b>Phasic On:</b> Turn the phasic tracing on.
			P1	<b>Scale:</b> Set or change the pressure scale.
			P1	<b>Set Active:</b> Make the pressure channel the active channel.
		<b>P2 – P4 Review:</b> ALL the same commands and functions.		
		<b>Active</b>		<b>Change Label:</b> Set or change a specific channel label.
				<b>Mean Off:</b> Turn the mean pressure tracing off.
				<b>Mean On:</b> Turn the mean pressure tracing on.
				<b>Phasic Off:</b> Turn the phasic tracing off.
				<b>Phasic On:</b> Turn the phasic tracing on.
				<b>Scale:</b> Set or change the pressure scale.
				<b>Zero:</b> Zero the individual pressure channel.
		<b>NOTE</b> Active means the channel is select in the <i>Real-Time</i> window.		

Measures	Pressure	Real-Time	All Pres-sures	All Pressures: Measure all of the available pressure waveforms.
			DFR	Do DFR measurement.
			FFR	Do FFR measurement.
			P1	<b>P1:</b> Measure the current pressure waveform on this channel.
			P2	<b>P2:</b> Measure the current pressure waveform on this channel.
			P3	<b>P3:</b> Measure the current pressure waveform on this channel.
			P4	<b>P4:</b> Measure the current pressure waveform on this channel.
		<b>Review</b>	All Pres-sures	<b>All Pressures:</b> Measure all of the available pressure waveforms.
			P1	<b>P1:</b> Measure the current pressure waveform on this channel.
			P2	<b>P2:</b> Measure the current pressure waveform on this channel.
			P3	<b>P3:</b> Measure the current pressure waveform on this channel.

CONFIDENTIAL

SECRET DE AFACERI

			P4	<b>P4:</b> Measure the current pressure waveform on this channel.
		<b>Active</b>		Measure the active pressure in the <i>Real-Time</i> window.
		<b>A-wave Analysis On/Off</b>		Prompt user to enter on or off parameter. <b>NOTE</b> Should be used if A-Fib or A-Flutter is present.
	<b>Manual Cardiac Output</b>			Open the <i>Manual Cardiac Output</i> window.
	<b>O2 Saturation</b>			Open the <i>Blood Flow Analysis</i> window and an O <sub>2</sub> sat entry line.
	<b>Peak to Peak</b>			Open the <i>Peak-to-Peak</i> window.
	<b>Pullback</b>			Open the <i>Pullback</i> measurement window.
	<b>Thermodilution Output</b>			Open the <i>Thermodilution Cardiac Output</i> window.
	<b>Valve Analysis</b>			Open the <i>Valve Analysis</i> measurement window.

<b>Results</b>	<b>Report</b>	<b>Generate</b>	Generate specific report formats by name of customization.
		<b>Print</b>	Print a specific report that has been generated.
		<b>Print All</b>	Print all report formats in the <i>Reports</i> window.
	<b>Blood Flow Analysis</b>		Open the <i>Blood Flow Analysis</i> window (Fick, shunt, and O <sub>2</sub> sats).
	<b>Pressure</b>		Open the <i>Pressure</i> results window.
	<b>Resistances</b>		Open the <i>Resistances</i> results window.
	<b>Stroke Work</b>		Open the <i>Stroke Work</i> results window.
	<b>Valve Area</b>		Open the <i>Valve Area</i> results window.

<b>Vitals</b>	<b>Log NBP On/Off (Auto NBP On/Off)</b>	Start the NBP cuff inflation and begin logging vitals per interval set. <b>NOTE</b> Set the <b>Auto NBP</b> interval from the <b>Measurements</b> tab in Set-Up.
	<b>Log Vitals Interval</b>	Prompt user for the number of minutes between vital logging. <b>NOTE</b> Do not use this command if auto-logging NBP.
	<b>Log Vitals On/Off</b>	Turn on and off the auto-logging of vitals.
	<b>Set Vital Pressure</b>	Allow selection of P1, P2, P3, P4, or NBP for blood pressure logging. <b>NOTE</b> For most catheterization labs, this should be set to NBP.
	<b>Start/Stop NBP</b>	Immediately start or stop the NBP cuff. (STAT NBP).

<b>Study</b>	<b>Case Event</b>	Open the <b>Case Event</b> dialog box.
	<b>Complications</b>	Open the <b>Complications</b> master list.

CONFIDENTIAL

SECRET DE AFACERI

	<b>Conscious Sedation</b>	Open the pre/post <i>Conscious Sedation</i> window for documentation.	
	<b>Contrast</b>	Open the <b>Contrast</b> master list.	
	<b>Forms</b>	Open the selected registry form or custom form from the drop-down list.	
	<b>Medications</b>	Open a dialog box with a specific drug selected and allows input of drug amount and time. <b>NOTE</b> The medication must already exist in the medication master list. When setting this macro, enter the medication exactly as it is entered in the master list.	
	<b>Medications List</b>	Open the <b>Medication</b> master list.	
	<b>Patient Info</b>	Open the <i>Patient Information</i> window.	
	<b>Pediatric On/Off</b>	Put the system in pediatric mode for additional pressure site labels.	
	<b>Phase</b>	Open the <b>Phase List</b> for changes during a procedure.	
	<b>Procedures</b>	Open the <b>Procedures</b> master list.	
	<b>Radiology</b>	Open the <b>Radiology</b> master list.	
	<b>Report</b>	Open the <i>Reports</i> window.	
	<b>Save On/Off</b>	Turn save on and off (toggles).	
	<b>ST Segment</b>	Open the <i>ST Segment</i> window.	
	<b>Staff</b>	Open the <b>Staff</b> master list.	
	<b>Supplies</b>	Allow the user to enter a specific supply (for example, 6F Sheath). <b>NOTE</b> The supply must already exist in the supply master list. When setting this macro, enter the supply exactly as it is entered in the master list. When the macro is executed, a dialog box will appear stating <b>Item Found</b> . Click <b>Insert</b> .	
	<b>Supplies List</b>	Open the <b>Supply</b> master list.	
<b>Display</b>	<b>Real-Time</b>	<b>Active Page</b>	Set the active page to 1 or 2.
		<b>Grid Major</b>	Set the invasive blood pressure grid to 3 graticules.
		<b>Grid Minor</b>	Set the invasive blood pressure grid to 11 graticules.
		<b>Grid Subminor</b>	Set the invasive blood pressure grid to 41 graticules.
		<b>Print 12 Lead</b>	Print a real-time 12 lead (7 lead if using five ECG cables). <b>NOTE</b> Active Page: 2 must be set before attempting to acquire a real-time 12 lead.
		<b>Sweep Speed</b>	Allow user to set or change sweep speed in the <i>Real-Time</i> window.
	<b>Clock</b>	<b>Close Stopwatch</b>	Close the <i>Stopwatch</i> window.
		<b>Close Timer</b>	Close the <i>Timer</i> window.
		<b>Reset Stopwatch</b>	Reset the stopwatch (counts forward).

CONFIDENTIAL

SECRET DE  
AFACERI

		<b>Set Stopwatch Title</b>	Allow the user to label the stopwatch.
		<b>Set Timer</b>	Allow the user to set the timer (in minutes).
		<b>Set Timer Title</b>	Allow the user to label the <i>Timer</i> window.
		<b>Start/Stop Stopwatch</b>	Start/stop the stopwatch.
		<b>Start/Stop Timer</b>	Start/stop the timer (counts down).
	<b>Review</b>	<b>Active Page</b>	Allow selection of the active page.
		<b>Grid Major</b>	Set the invasive blood pressure grid to 3 graticules.
		<b>Grid Minor</b>	Set the invasive blood pressure grid to 11 graticules.
		<b>Grid Subminor</b>	Set the invasive blood pressure grid to 41 graticules.
		<b>Print</b>	Perform a quick print of the <i>Review</i> window.
		<b>Print 12 Lead</b>	Print a 12 lead (7 lead if using five ECG cables). <b>NOTE</b> <b>Active Page: 2</b> must be set before attempting to acquire a review 12 lead.
		<b>Sweep Speed</b>	Allow user to set or change the sweep speed in the <i>Review</i> window.
		<b>Update</b>	Save 10 seconds of data and update the <i>Review</i> window.
	<b>Image</b>	<b>Save</b>	Grab an image from the <i>Image 1</i> or <i>Image 2</i> windows (optional feature) and save to the study.
		<b>Save ICE</b>	Grab an image from the <i>Image 3</i> window (optional feature) and save to the study.
	<b>Printer Setup</b>		Opens the <i>Printer Set-Up</i> window to change settings.
<b>Macro Tools</b>	<b>Delay</b>		Allow user to initiate a specific delay in seconds before continuing to run the rest of the macro.
	<b>Note</b>		Open the <i>Note</i> window for a real-time comment to be charted.
	<b>Pause</b>		Create a pause in the macro command. <b>NOTE</b> The user must select the <b>OK</b> box on the <i>Review</i> screen in order to continue the macro.
<b>Intervention</b>	<b>Add Coronary Intervention</b>		Open the <i>Coronary Intervention Event</i> window.
	<b>Add Graft Intervention</b>		Open the <i>Graft Intervention Event</i> window.
	<b>Add Peripheral Intervention</b>		Open the <i>Peripheral Intervention Event</i> window.

CONFIDENTIAL

SECRET DE AFACERI

## 7.3.2 ST Segment and the Log window

An ST Segment Event is placed in the *Log* window with each ST Segment reference obtained.

The ST Segment event contains the time of the event and an automatically generated segment number (example ST Segment: 1) description. Any saved segment may be brought into the *ST Segment* window by selecting the event time. This will also update the signal data in the *Review* window. Double-clicking on the event summary will bring up the *ST Segment Event Description* dialog box, in which the segment description can be changed.

An ST Segment section may be added to any report format. This section will display the time of the ST Segment, as well as the ST Segment values for all surface leads. The report will separate the segment measurements by phase. See [Chapter 16 Reports on page 199](#).

## 7.4 Image Viewer

Use Image Viewer to view and annotate snapshots and images in a study.

During a study, you can annotate a captured image to highlight unique characteristics of the captured data. Annotations can be saved with the image as part of the study. The image and annotations can also be saved to a file or copied to the clipboard.

### 7.4.1 Annotation

#### 7.4.1.1 Add annotation objects

Add lines, rectangles, ellipses, polylines, polygons, pointers, freehand drawings, highlights, text, and notes to the image.

1. Select an annotation type in the annotation toolbar, then left-click and drag.
2. When adding a text or note annotation, press **Ctrl + Enter** to create a new line of text below the current line.
3. To complete adding a text or note annotation, press **Enter** or left-click outside of the object.

#### 7.4.1.2 Select annotation objects

Copy, paste, move, resize, delete, and change properties of annotations.

- To select an individual object, move the cursor over the object until the cursor turns into a circle, then left-click.
- To select multiple objects:
  1. Hold down the shift key.
  2. Move the cursor over an object until the cursor turns into a circle, then left-click.
  3. Repeat the previous step for each object you want to select.
  4. Release the shift key.
  - Alternatively, draw a focus box around the objects by holding down the left-mouse button while dragging.
- To select all objects, right-click in the image and choose **Select All**.

## 6.3.4 FFR and DFR equalization

When using FFR or DFR in a study, the **P1** and **P2** channels must be equalized before you do an FFR or DFR measurement. When the catheter and pressure wire are in place, right-click **FFR/DFR Value Display** to do the equalize function.

- The FFR or DFR display must be turned on. Select either **Real-Time > FFR Display On** or **Real-Time > DFR Display On** to turn on the display or toggle between them..
- To turn on the FFR or DFR display with a macro, see [Chapter 14 Macros on page 179](#) for instructions.

After you select **Equalize**, the **FFR/DFR Value Display** shows that the system is **Equalizing**. This process may fail if the catheter and pressure wire are too far apart. If this happens, the FFR/DFR Value Display again shows **Equalize**. The largest allowable offset between the two channels is 10 mmHg, and the largest allowable time difference is 100 ms. Adjust the position of the catheter and pressure wire and select **Equalize** again.

When equalization is complete, the system calibrates the **P1** and **P2** channels so that mean, amplitude, and time are aligned. The system inserts a note into the study with the mean calibration value, which is called the offset. The system also displays the Pd/Pa ratio in the pressure status area. You can now do FFR or DFR measurements. See [12.3.9 Do an FFR measurement on page 159](#) or [12.3.10 Do a DFR measurement on page 160](#) for more information.

The calibration resets when the FFR or DFR channel is turned off, the study is closed, or if either **P1** or **P2** is zeroed. Calibration is maintained when you toggle between FFR and DFR without turning off the channel. This allows the system to change between FFR and DFR measurements without recalibrating.

- Follow FFR or DFR accessory preparation and checks recommended by the accessory manufacturer. Ensure accessories are not damaged prior to use.
- To equalize FFR or DFR (equalize **P1** and **P2** channels) using a macro instead of right-clicking the **FFR/DFR Value Display**, see [Chapter 14 Macros on page 179](#) for instructions.

## 6.3.5 Shortcut to turn channels on and off

### 6.3.5.1 Turn off a channel

Channels can be turned on or off quickly using the following shortcut:

## 12.3.9 Do an FFR measurement

### CAUTION

#### INACCURATE VALUES

Erroneous FFR or DFR values may occur with irregular heartbeats. Computerized interpretation is only significant when used in conjunction with clinical findings. Ensure values are verified by a qualified healthcare practitioner.

Note the following when doing an FFR measurement:

- Follow FFR accessory preparation and checks recommended by the accessory manufacturer. Ensure accessories are not damaged prior to use.
- Ensure that the beat count is set to the appropriate value.
- Pressure channels must be equalized before you do an FFR measurement. When equalization is completed properly, the FFR value should be close to 1.00.
- Before you do a real-time FFR measurement:
  - Both **P1** and **P2** must be active.
  - The FFR channel must be active.

### NOTE

Use caution when choosing labels. Pa and PA are both listed. When choosing Pa (Pressure artery - for FFR) or PA (Pulmonary Artery), the values will fall into different groups. Be sure to choose the correct abbreviation.

### NOTE

The pressure wire connector must be connected to **P2** for the FFR measurement to be valid.

1. To open the *FFR* window, select **Real-Time > FFR Display On**.
2. Do the FFR equalization as described in [6.3.4 FFR and DFR equalization on page 88](#).
3. To do a real-time DFR measurement, do one of the following:
  - Select **Measurements > FFR/DFR** or **Measurements > FFR**.
  - Execute the macro command **Measures > Pressure > Real-Time > FFR > .**
  - Push **Action + F** on the keyboard.

The *DFR Measurement* dialog opens.

4. Select a **Coronary Artery Segment**, as needed.
5. Select **Start** to initiate the FFR measurement process.


While the FFR measurement is in progress, the *FFR Measurement* dialog remains open and provides the option to complete or cancel the measurement.

If full disclosure saving is not on, saving turns on for the duration of the FFR measurement.

6. Select **Stop** to complete the FFR measurement.

The system inserts an FFR event into the study at the point of the minimum FFR value. This event is visible in the *Log* window and the *Review* window.

Figure 7-1 Holter window

Item	Name	Description
1	<b>Signal</b>	Select the desired signal from the <b>Signal</b> drop-down menu.
2	<b>Speed</b>	You can control how much of the signal is visible from the <b>Speed</b> drop-down menu. A smaller number displays the beats closely together. A larger number displays the beats further apart.
3	<b>Lines</b>	You can select the desired number of lines to display in the <i>Holter</i> window.

## 7.3 ST Segment window

The *ST Segment* window is a non-invasive method of recording changes in the electrical conduction system of the heart throughout an invasive procedure. ST segment information may be saved at any point during a case and compared to previous saved segments. The continuous read-out of the ST segment status allows the user an instant comparison between the saved segments and the patient's current ST status. Mac-Lab has the ability to store ST segment values and print these values in the report.

The *ST Segment* window allows for capturing ECG signals and the associated ST segment measurement from the cath amplifier. This window may be opened side by side with the *Real-Time* window for comparison of reference with real-time ST segments.

To launch the window, place the *Real-Time* screen on page 2. Select **Windows > ST Segment**. Position the cursor to the left of the *Real-Time* window and click to drop the *ST Segment* window in place. The *ST Segment* window appears.

The real-time ST segment value is displayed in the same color as the associated ECG signal.

The reference ST segment is displayed in white.

### 11.1.3 Fick cardiac output

The Fick cardiac output is calculated by entering saturations in the *Blood Flow Analysis* window. The *VO2 Worksheet* (oxygen consumption) window is used to calculate the estimated oxygen consumption value.

1. Press the **Action + J** keys. The *Blood Flow Analysis* window appears.
2. Select the label from the pull-down menu in the **Oximetry Inputs** field or type the abbreviation for the site label in the **Oximetry Input** field.
3. Use the **TAB** key through the saturation line or select the window.
4. Type the **Hb, Sat** and **pO2** values.

**NOTE**

The **pO2** values are shown in units of either mmHg or kPa, depending on the mode previously selected by the user.

5. Click the **Insert Item** icon to add additional blood samples.
6. Finish each saturation line by pressing **Enter**.

**NOTE**

A blood sample from the mixed venous group (MV) and a blood sample from the systemic artery group (SA) or pulmonary artery group (PA) are required for a Fick cardiac output.

#### 11.1.3.1 VO2 worksheet

1. Select **VO2 Worksheet** in the *Blood Flow Analysis* window. The *VO2 Worksheet* window appears.
2. Select **VO2 Formulas**. The *VO2 Formula Selection* window appears.
3. Select the appropriate radio button for the desired method for calculating VO<sub>2</sub>.

- LaFarge: type a heart rate to obtain a VO<sub>2</sub> value. If other formulas are used, heart rate information is not needed.

**NOTE**

K1 is 11.49 for male and 17.04 for female.

- K\*BSA: use the following table

**NOTE**

K \* BSA is generally used for pediatric patients.

**Table 11-1 K\*BSA**

Age	K*BSA
< 3 months	152 * BSA
3 months to 1 year	176 * BSA
1 year to 2 years	175 * BSA
2 years to 4 years	168 * BSA
> 4 year	149 * BSA

4. Select **OK** to return to the **VO2 Worksheet** window.

## A.4 Continuously recorded parameters

The following values are acquired and continuously recorded by the cath amplifier for display on the Mac-Lab system.

Parameter	Units
Continuously recorded parameters	
Heart Rate	BPM
Non-invasive Blood Pressure	mmHg or kPa
Respiration Rate (by thoracic impedance)	breaths per minute
SpO <sub>2</sub>	%
ST segment	60 ms post J
Temperature	degrees Celsius
Continuously recorded parameters from the CO <sub>2</sub> module (optional)	
Inspired CO <sub>2</sub>	mmHg
Expired CO <sub>2</sub>	mmHg
Respiration Rate	breaths per minute

## A.5 Automatically derived formulas

### Computations

The hierarchy in computations requiring cardiac output for valve areas where multiple cardiac outputs have been measured is Manual Entry, Thermodilution and FICK CO. For stroke work and resistance calculations where multiple cardiac outputs have been measured, the hierarchy is Manual Entry, Qp or Qs, Thermal CO and FICK CO.

#### Mitral Diastolic Filling Period: (DFPm)

Formula:  $DFPm = (LVFT * HR) / 1000$

Units: sec/min

Reference: Grossman, 4th ed., p. 153

Left ventricular filling time (LVFT) is in msec

#### Fractional Flow Reserve: (FFR)

Formula:  $Mean Pd / Mean Pa$

Units: None

Reference: Kern, 4th ed., p. 512

### NOTE

The Mac-Lab system always treats pressure channel 2 as Pd and pressure channel 1 as Pa, regardless of current pressure site labels.

#### Tricuspid Diastolic Filling Period: (DFPt)

Formula:  $DFPt = (RVFT * HR) / 1000$

**NOTE**

The current heart rate, SpO<sub>2</sub>, respiration rate, temperature, blood pressure (either invasive or NBP as indicated in the Measurement Configuration (see [6.6.7.2 Measurements on page 114](#)) and optional inspired and/or expired CO<sub>2</sub> are recorded in the *Log* window.

### 6.4.2.19 Auto-log NBP (optional)

To automatically log non-invasive blood pressure measurements with vitals, select **Measurements > Auto-Log NBP**. Set the vitals configuration to **Auto-Log NBP** when using NBP as the vital logging pressure source.

**Auto-Log Vitals** is most useful when you use an invasive blood pressure as the blood pressure channel.

Double-click the vitals event in the *Log* window to open a vitals event window. You can enter a level of consciousness and a level of pain in the vitals event window. Because the auto-log NBP records a set of vitals with each reading, this information may be redundant.

### 6.4.2.20 Vital configuration

Select **Measurements > Vital Configuration** to open the *Vital Configuration* tab in the *System Settings* window. See Vital Configuration.

### 6.4.2.21 Audible indicators

**NOTE**

The following instructions only apply when audible indicators are enabled from the *Vital Configuration* tab (see Vital Configuration).

**WARNING****AUDIBLE INDICATORS MALFUNCTION**

The audible indicators are for reference only. The audible indicators are not designed for use as a patient alarm.

1. Perform one of the following to silence the audible indicator:
  - Click the **Turn On Or Off Audible Indicator** icon in the toolbar.
  - Select **Measurements > Silence Audible Indicator** from the *Real-Time* or *Review* screen.
  - Press **Ctrl + D** on the keyboard.
2. Click **Yes** at the confirmation window.

**NOTE**

A confirmation window is not displayed if the icon is selected from the toolbar.

### 6.4.2.22 Data extraction

All channels of data can be exported from the system for a user-selectable time. The data is extracted in either binary or text.

On large data sets, the data extraction operation may take several minutes to complete. During this time, the application does not respond to other user inputs. If you attempt to interact with either the application or the *Data Extraction* window, the system displays an error message and the data extraction bar stops advancing. If this happens, ignore the message. The operation will complete successfully even if the data extraction bar stops advancing.

# Manual de utilizare Mac-Lab

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SECRET DE  
AFACERI

5222001-1RO  
Revizia 1.0  
Romanian | Limba română  
© 2022 General Electric Company

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