**Classical Swine Fever Virus Antigen Test Kit
(CSFV) / Serum Plus**

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Test With Confidence™

IDEXX

Classical Swine Fever Virus Antigen Test Kit (CSFV) / Serum Plus

For veterinary use only.

Name and Intended Use

IDEXX CSFV Ag Serum Plus is IDEXX's enzyme immunoassay for the detection of CSFV antigens in swine serum, plasma and tissue (preferably tonsils, spleen, mesenteric lymph node or kidneys) samples.

General Information

Classical Swine Fever Virus (CSFV), Bovine Viral Diarrhoea Virus (BVDV), and Border Disease Virus (BDV) are the 3 members of the genus *Pestivirus* within the family *Flaviviridae*. CSFV causes serious losses in the pig industry since it is an obligate pathogen and can cause widespread deaths. Pigs infected with highly virulent CSFV strains may shed a high amount of virus before showing clinical signs of the disease. Animals that survive an acute or sub-acute infection develop antibodies and will no longer spread the virus. Moderately virulent, less pathogenic strains may lead to chronic infection, when pigs excrete the virus continuously or intermittently until death. Fetuses of susceptible pregnant sows, exposed to CSFV, might become infected through the placenta. Congenital infection can result in abortion, fetal mummification, stillborn and/or weak piglets or embryonic malformations. The most frequent outcome of congenital infection with low virulent strains is the birth of persistently infected piglets, spreading the virus without signs of the disease, in the absence of immune response. Pigs may also be infected with BVDV or BDV. Although these infections in pigs are usually mild and self limiting, virus detection in such cases is very rare. Since the IDEXX CSFV Ag Serum Plus test has high sensitivity for all pestiviruses, the detection of BVDV or BDV cannot be excluded. Positive results should therefore be confirmed with a CSFV specific assay.

Descriptions and Principles

A microtitration format has been configured by immobilizing specific monoclonal antibodies on the plates. CSFV Ag of the sample is captured on the plates. After incubation of the test sample in the well, captured CSFV Ag is detected by specific polyclonal antiserum for CSFV (E^{rns}) and a horseradish-peroxidase Conjugate. Next, unbound Conjugate is washed away and a Substrate/chromogen solution is added. In the presence of enzyme, Substrate is converted into a product which reacts with the chromogen to generate a blue color. Upon addition of the Stop Solution, a yellow color is generated. The absorbance at a single wavelength of 450 nm [A(450)] or a dual wavelength of 450 nm and 650 nm [A(450/650)] is measured using a spectrophotometer. The corrected OD value of the samples is calculated by using the absorbance [A(450)] or [A(450/650)] obtained with the test sample and corrected for the absorbance of the Negative Control.

Reagents		Volume
1	Anti-Er ^{ns} mAb Coated Plate	2
2	Positive Control	1 x 1.6 mL
3	Negative Control	1 x 1.6 mL
4	Conjugate	1 x 25 mL
5	Tissue Soaking Buffer Concentrate (2X)	1 x 115 mL
A	TMB Substrate N.12	1 x 20 mL
B	Stop Solution N.3	1 x 20 mL
C	Wash Concentrate (10X)	1 x 125 mL
D	Detection Solution	1 x 15 mL
Other Components: Zip lock bag		1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with single 450 nm filter or dual filters at 450 nm and 650 nm)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent
- Microplate covers (lid, aluminium foil or adhesive) and/or humid chamber
- Centrifuge (capacity 1500 x g)
- Incubator capable of maintaining a temperature of +37°C ($\pm 3^{\circ}\text{C}$)

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Reagents

Wash Solution

The Wash Concentrate (10X) must be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1:10 with distilled/deionized water before use (e.g., 30 mL of Wash Concentrate (10X) plus 270 mL of distilled/deionized water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Tissue soaking buffer

The Tissue Soaking Buffer Concentrate (2X) must be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Tissue Soaking Buffer Concentrate (2X) must be diluted 1:2 with distilled/deionized water before use. The amount of diluted Tissue Soaking Buffer needs to be calculated as required. When prepared under sterile conditions, the diluted Tissue Soaking Buffer can be stored for up to two weeks at 2–8°C.

Preparation of Samples

Serum and Plasma Samples

Fresh or frozen serum and plasma samples can be tested.

Tissues

Use preferably fresh tissues but, if necessary, tissues can be stored frozen. Process one or two tissues from each animal submitted, preferably tonsils, spleen, mesenteric lymph node or kidneys.

- a. With scissors, cut tissue into small pieces (250 mg).
- b. Place the tissue pieces into an appropriate centrifuge tube or an Eppendorf® tube and add 1 mL Tissue Soaking Buffer. Vortex. Leave at 18–26°C for 1–2 hours.
- c. Centrifuge at 1,500 x g for 10 minutes and use 50 µL of the clean supernatant for testing as described in the Test Procedure.

The overnight incubation protocol (12–18 hours at 2–8°C) should be used for tissue samples.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

- 2 Dispense 50 μ L of Detection Solution into each well.

- 3 Dispense 50 μ L of Negative Control (NC) into duplicate wells.

- 4 Dispense 50 μ L of Positive Control (PC) into duplicate wells.

- 5 Dispense 50 μ L of samples into remaining wells.

- 6 Mix the content of the microwells by gently tapping the plate or use a microplate shaker.

- 7 Incubate samples:
 - a. Serum and plasma samples: incubate for 2 hours (\pm 5 min.) at +37°C (\pm 3°C) or 12–18 hours at 2–8°C. With either option, plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.
 - b. Tissue samples: incubate overnight (12–18 hours) at 2–8°C. The plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.

- 8 Remove the solution and wash each well with approximately 300 μ L of Wash Solution 5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 9 Dispense 100 μ L of Conjugate into each well.

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- 13 Incubate 10 minutes (\pm 1 min.) at 18–26°C away from direct light.

14 Dispense 100 μ L of Stop Solution N.3 into each well.

15 Measure and record the absorbance of the samples and controls at 450 nm or using a dual wavelength of 450 nm and 650 nm.

16 Calculations:

Controls

$$NC\bar{x} = \frac{NC1 A(450) + NC2 A(450)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450) + PC2 A(450)}{2}$$

Validity criteria

$$PC\bar{x} - NC\bar{x} \geq 0.150$$

$$NC\bar{x} \leq 0.250$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S-N = \text{Sample A}(450) - NC\bar{x}$$

The presence or absence of CSFV antigen in the sample is determined by corrected OD value (S-N) for each sample.

17 Interpretation:

Negative

Positive

$$S-N \leq 0.300$$

$$S-N > 0.300$$

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: [idexx.com/contact](https://www.idexx.com/contact)

* IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. All other product and company names and logos are trademarks or registered trademarks of their respective holders.

© 2023 IDEXX Laboratories, Inc. All rights reserved.

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DANGER / DANGER / PERIGO / PELIGRO / GEFAHR

H317 / H412 / P261 / P280 / P333+P313

Positive control / Negative control / Conjugate / Tissue Soaking Buffer Concentrate (2X) / Wash concentrate (10X) / Detection Solution – May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapors. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention.

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Stop solution No.3 – Causes severe skin burns and eye damage. May cause an allergic skin reaction. Causes serious eye damage. Do not breath mist/vapors. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. If skin irritation or rash occurs: Get medical advice/attention.

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IDEXX Laboratories, Inc.
One IDEXX Drive
Westbrook, Maine 04092
USA

Manufacturer
IDEXX Switzerland GmbH
Stationsstrasse 12
CH-3097 Liebefeld-Bern
Switzerland

EU-Representative
IDEXX B.V.
P.O. Box 1334
2130 EK Hoofddorp
The Netherlands

idexx.com

IDEXX