

ORIGIN 



Signature
ORTHOPAEDICS

Table of Contents

Introduction	3
Origin Stem Features	4
Surgical Technique	5
<i>Preoperative Planning</i>	5
<i>Femoral Neck Resection</i>	6
<i>Femoral Preparation</i>	7
<i>Compaction Broaching</i>	8
<i>Calcar Reaming</i>	8
<i>Trial Reduction</i>	9
<i>Femoral Component Insertion</i>	10
<i>Femoral Head Impaction</i>	11
Origin Implant Sizing Guide	12
Origin Instrument Tray	13
Origin Instruments	14
Origin Preoperative Templates	17
Origin Implants	18
Appendix: Compaction Broaching Techniques and Tips	20

Introduction

The Signature Orthopaedics Origin Instrument system is an optimised instrument set for implantation of the Signature Orthopaedics cementless hip stem. The Origin instrument set features an adaptable broach handle that quick-connects to several instrument attachments to reduce the overall number of instruments and minimize the size and weight of the instrument tray.

Indications

Signature Orthopaedics' hip replacement range is intended to replace a hip joint where bone stock is sufficient to support the implant. When a surgeon has selected prosthetic replacement as the preferred treatment, the devices are indicated for:

- Non-inflammatory degenerative joint disease including osteoarthritis or avascular necrosis
- Inflammatory joint disease including rheumatoid arthritis (excluding TSI stem)
- Correction of functional deformity including congenital hip dysplasia
- Traumatic injury involving the hip joint including traumatic arthritis or femoral head or neck fracture
- Failed previous hip surgery including internal fixation or joint fusion, reconstruction, hemiarthroplasty, surface replacement, or total replacement.

Contraindications

In general, prosthetic components require adequate bone support for correct fit and function. The use of prosthetic components is therefore contraindicated where any pathological condition may reduce the quantity and or strength of the bone which is supporting the prosthesis. Some contraindications are relative to the extent and severity of conditions and the benefits of prosthetic arthroplasty should be considered based on the patient's overall evaluation and the possibility of alternative treatment. Examples of such conditions include; osteoporosis, osteomalacia, osteogenesis imperfecta, or hypophosphatemia. Other contraindications include:

- Conditions limiting blood supply to the bone or joint.
- Systemic or local infection.
- Previous high dose radiotherapy.
- Psychological or neurological conditions which would restrict the patient's ability or compliance in restricting physical activity.
- Skeletal immaturity
- Conditions or activity which may place excessive load on the components such as; obesity, muscle, tendon & ligament deficiencies, multiple joint disabilities, and Charcot joints.

Origin Stem Features

Cementless Hip Stem

- Proven geometry, material and coating (Ti6Al4V with HA coating).
- Threaded proximal feature aids in positioning and removal.

1. Multiple offsets

2. 12/14 Taper

3. Low-profile lateral shoulder,

Enables easy insertion in reduced insertion techniques, including anterior approach.

4. Collared, non-collared and Coxa Vara Versions

5. Step geometry


Converts hoop stresses to compressive loads.

6. Grooves

Vertical and horizontal grooves provide rotational and axial stability.

7. Hydroxyapatite Plasma Spray Coating


- Thickness 120-190 μ m
- Tensile Strength >15 MPa
- Crystallinity \geq 80%
- HA Content \geq 70%



1

Preoperative Planning

Origin X-Ray templates can be used over anterior/posterior and lateral radiographs to help determine the correct size to restore the patient's anatomy. Templates are 115% magnification.



2


Preoperative Planning

The Logical cup can be used with any surgical approach that the surgeon selects.

- a. Posterior approach
- b. Posterolateral/anterolateral approach
- c. Anterior approach

Note:


Prior to the following steps, complete all steps detailed in 111-12-0003 for the Logical acetabular cup implantation.






Femoral Neck Resection

The osteotomy guide should be used in conjunction with preoperative planning, to determine the level of the femoral neck resection. This can be performed in multiple steps, depending on surgeon preference.



Optional technique:


The femoral head extractor may be used with the T-handle or under power to aid in the removal of the resected head, especially during an anterior approach technique.



4




Femoral Preparation

Enter the femoral canal as laterally as possible with the Box Osteotome to initiate access to the medullary canal. The Canal Reamer may be used as needed to open the natural axis of the femoral canal for broach preparation.

**Optional technique:**

While the Origin system is intended to be a broach-only system, the Origin instrument tray contains instruments for optional use to ensure proper axial alignment along the femoral canal and to induce lateral bias where needed:


- Starter broach induces lateral bias by rasping beneath the greater trochanter
- Canal reamer creates a guide hole for the distal end of the trochanteric reamer
- Trochanteric reamer removes trochanteric bone tissue laterally to ensure neutral alignment with the femoral axis during broaching



5

Compaction Broaching

The broach should run parallel to the posterior cortex following the natural anatomy of the femur. Begin with the smallest broach and increase the size of the broach sequentially until longitudinal and rotational stability is achieved; broaching should then be stopped. Careful preoperative planning is key to help selection of the final broach size. The version will be determined by the natural version of the femur.


**Note:**


Refer to Appendix: **Compaction Broaching Techniques and Tips** for detail on the theory of compaction broaching, why this technique is used for Origin stems, and how to form a bed of compacted bone that will maximize the longevity and stability of Origin stems in many bone types.

6



Calcar Reaming

With the broach in situ, use the Calcar Reamer to achieve a flat resection surface. Slide the reamer over the broach quick connect fitting to maintain the resection angle. Carefully advance the reamer towards the broach face and into the resected edge of the femur until it bottoms out against the broach face.



**Trial Reduction**

With the final broach still in situ, attach the appropriate trial neck and trial head. Reduce the hip and assess what adjustments, if any, are required to provide stability through a full range of motion. Remove the trial head, trial neck and final broach. DO NOT irrigate or dry the femoral canal. This will help to preserve the compacted cancellous bone quality and encourage biological fixation of the stem.



**Instrument identification:**

Trial heads are colour coded based on offset. Refer to Origin Implants Sizing Guide in this surgical technique for more details.




Femoral Component Insertion

When implanting the definitive stem (same size as final broach) in the femoral canal, ensure that it is directed in by hand. This will help avoid changing the version as a precautionary measure. There should be no more than 15-20mm between the resection line and the top of the HA coating on the stem. If the stem does not readily go down this far, the surgeon should broach again. Once the stem is placed, lightly tap the stem impactor to fully seat. DO NOT over-impact as this may lead to splitting of the femur.



Instrument operation:


To connect the Origin stem to the stem positioner, first slide the inner shaft of the positioner through the outer shaft, spinning the strikeplate so the inner shaft threads in and falls through. Turn the threaded tip of the inner shaft into the female threads on the Origin stem until a snug hold is achieved to prevent damage to the threads.



9

Femoral Head Impaction


A final trial reduction is carried out to confirm joint stability and range of motion. Clean and dry the stem taper to remove any particulate matter or debris. Place the femoral head onto the taper and lightly tap it using the head impactor. Ensure that bearing surfaces are clean and finally reduce the hip.




Origin Implant Sizing Guide

Cementless Hip Stems and Femoral Heads


SIZE	OFFSET	LT	L	W	O	α°	F	G	
8	Standard Offset	116	94	38.2	36.3	135	6.3	7.0	
	Standard Offset Collared	115	94	38.2	38.2	135	6.3	7.9	
9	Standard Offset	131	109	40.3	36.6	135	6.5	8.1	
	High Offset	131	109	41.7	45.5	135	6.3	8.1	
	Standard Offset Collared	131	109	38.2	39.5	135	6.3	8.1	
	Cox Vara	121	109	40.3	45.5	125	6.3	8.1	
10	Standard Offset	141	119	40.3	39.5	135	6.6	8.5	
	High Offset	141	119	42.4	46.5	135	6.6	8.5	
	Standard Offset Collared	141	119	38.2	39.5	135	6.6	8.5	
	Cox Vara	141	119	40.3	46.5	125	6.6	8.5	
11	Standard Offset	146	124	40.3	42.0	135	7.5	9.0	
	High Offset	146	124	42.0	47.0	135	7.5	9.0	
	Standard Offset Collared	146	124	38.2	40.0	135	7.5	9.0	
11	Cox Vara	146	124	40.3	47.0	125	7.5	9.0	
	12	Standard Offset	151	129	40.3	41.0	135	8.4	9.5
		High Offset	151	129	42.4	46.0	135	8.4	9.5
Standard Offset Collared		151	129	38.2	41.0	135	8.4	9.5	
12	Cox Vara	151	129	40.3	45.0	125	8.4	9.5	
	13	Standard Offset	156	134	40.3	41.5	135	9.3	9.8
		High Offset	156	134	42.4	48.5	135	9.3	9.8
Standard Offset Collared		156	134	38.2	41.5	135	9.3	9.8	
13	Cox Vara	156	134	40.3	48.5	125	9.3	9.8	
	14	Standard Offset	161	139	40.3	42.0	135	10.0	10.0
		High Offset	161	139	42.4	49.0	135	10.0	10.0
Standard Offset Collared		161	139	38.2	42.0	135	10.0	10.0	
14	Cox Vara	161	139	40.3	49.0	125	10.0	10.0	
	15	Standard Offset	166	144	40.3	43.0	135	10.8	10.0
		High Offset	166	144	42.4	50.0	135	10.8	10.0
Standard Offset Collared		166	144	38.2	43.0	135	10.8	10.0	
15	Cox Vara	166	144	40.3	50.0	125	10.8	10.0	
	16	Standard Offset	171	149	40.3	43.5	135	11.8	10.3
		High Offset	171	149	42.4	50.5	135	11.8	10.3
Standard Offset Collared		171	149	38.2	43.5	135	11.8	10.3	
16	Cox Vara	171	149	40.3	50.5	125	11.8	10.3	
	18	Standard Offset	181	159	40.3	46.5	135	13.4	11.0
		High Offset	181	159	42.4	52.0	135	13.4	11.0
Standard Offset Collared		181	159	38.2	44.5	135	13.4	11.0	
18	Cox Vara	181	159	40.3	52.5	125	13.4	11.0	



		Ø28	Ø32	Ø36	Ø40		
S	CrCo Ceramic	-3.5	+4.0	-4.0	-4.0		
	111-152-011	111-152-611	111-152-021	111-152-621	111-152-031	111-152-631	111-152-041
M	CrCo Ceramic	+0.0	+0.0	+0.0	+0.0		
	111-152-012	111-152-612	111-152-022	111-152-622	111-152-032	111-152-632	111-152-042
L	CrCo Ceramic	+4.0	+4.0	+4.0	+4.0		
	111-152-013	111-152-613	111-152-023	111-152-623	111-152-033	111-152-633	111-152-043
XL	CrCo Ceramic	+8.0	+8.0	+8.0	+8.0		
	111-152-024	111-152-624	111-152-034	111-152-634	111-152-044	111-152-644	




Origin Instrument



Standard Origin Instruments


Straight Broach Handle

112-182-032




Curved Broach Handle

112-182-001




Calcar Reamer

P12-01-0012




Osteotomy Guide

112-212-011




Quick Release Box Chisel

112-212-005




Quick Release Head Impactor

112-212-017



Quick Release Stem Impactor


112-212-006



Standard Origin Instruments


Origin Trial Neck

- 112-182-020 Standard Offset
- 112-182-021 High Offset
- 112-182-022 Coxa Vara




Stem Positioner

- 112-212-015 Assembly
- 112-212-021 Inner



Stepped Entry Reamer

- 112-162-001



IM Drill

- 112-182-087




Tapered Pin Reamer

- 112-182-013



Femoral Head Extractor


- 112-182-117



Origin Instruments


Origin Broaches

111-18-0809	Size 8
111-18-0909	Size 9
111-18-1009	Size 10
111-18-1109	Size 11
111-18-1209	Size 12
111-18-1309	Size 13
111-18-1409	Size 14
111-18-1509	Size 15
111-18-1609	Size 16
111-18-1809	Size 18




Trial Femoral Heads

111-182-040	Ø28mm -3.5mm Green
111-182-041	Ø28mm 0.0mm Yellow
111-182-042	Ø28mm +4.0mm Orange
111-182-017	Ø32mm -3.5mm Green
111-182-018	Ø32mm 0.0mm Yellow
111-182-019	Ø32mm +4.0mm Orange
111-182-020	Ø32mm +8.0mm Red
111-182-021	Ø36mm -3.5mm Green
111-182-022	Ø36mm 0.0mm Yellow
111-182-023	Ø36mm +4.0mm Orange
111-182-024	Ø36mm +8.0mm Red
111-182-043	Ø40mm -3.5mm Green
111-182-044	Ø40mm 0.0mm Yellow
111-182-045	Ø40mm +4.0mm Orange
111-182-046	Ø40mm +8.0mm Red



Origin Preoperative Templates


112-212-108	Size 8 Template
112-212-109	Size 9 Template
112-212-110	Size 10 Template
112-212-111	Size 11 Template
112-212-112	Size 12 Template
112-212-113	Size 13 Template
112-212-114	Size 14 Template
112-212-115	Size 15 Template
112-212-116	Size 16 Template
112-212-118	Size 18 Template



Optional Origin Instruments

Stem Positioner UniBody

112-212-026



Canal Reamer

112-212-023




Trochanteric Canal Reamer

112-212-022




XL Tapered Pin Reamer

112-182-148




Powered Femoral Head Extractor

112-25-0116




Curved Broach Handle

112-25-0044




Starter Broach

112-212-328




Quick Release Stem Inserter

112-212-323




Anteversio Grip

112-182-034



Broach Handle Power Adaptor

112-212-018



Origin Implants

Origin Standard Offset Stem

111-18-0800	Size 8
111-18-0900	Size 9
111-18-1000	Size 10
111-18-1100	Size 11
111-18-1200	Size 12
111-18-1300	Size 13
111-18-1400	Size 14
111-18-1500	Size 15
111-18-1600	Size 16
111-18-1800	Size 18

Origin High Offset Stem



111-18-0901	Size 9
111-18-1001	Size 10
111-18-1101	Size 11
111-18-1201	Size 12
111-18-1301	Size 13
111-18-1401	Size 14
111-18-1501	Size 15
111-18-1601	Size 16
111-18-1801	Size 18

Origin Standard Offset Collared Stem

111-18-0902	Size 9
111-18-1002	Size 10
111-18-1102	Size 11
111-18-1202	Size 12
111-18-1302	Size 13
111-18-1402	Size 14
111-18-1502	Size 15
111-18-1602	Size 16
111-18-1802	Size 18

Origin High Offset Collared


111-18-0903	Size 9
111-18-1003	Size 10
111-18-1103	Size 11
111-18-1203	Size 12
111-18-1303	Size 13
111-18-1403	Size 14
111-18-1503	Size 15
111-18-1603	Size 16
111-18-1803	Size 18



Origin Implants


Origin Coxa Vara Stem

111-18-0904	Size 9
111-18-1004	Size 10
111-18-1104	Size 11
111-18-1204	Size 12
111-18-1304	Size 13
111-18-1404	Size 14
111-18-1504	Size 15
111-18-1604	Size 16
111-18-1804	Size 18



Femoral Heads

111-152-011	Ø28mm -3.5mm CrCo
111-152-012	Ø28mm 0.0mm CrCo
111-152-013	Ø28mm +4.0mm CrCo
111-152-021	Ø32mm -3.5mm CrCo
111-152-022	Ø32mm 0.0mm CrCo
111-152-023	Ø32mm +4.0mm CrCo
111-152-024	Ø32mm +8.0mm CrCo
111-152-031	Ø36mm -3.5mm CrCo
111-152-032	Ø36mm 0.0mm CrCo
111-152-033	Ø36mm +4.0mm CrCo
111-152-034	Ø36mm +8.0mm CrCo
111-152-041	Ø40mm -3.5mm CrCo
111-152-042	Ø40mm 0.0mm CrCo
111-152-043	Ø40mm +4.0mm CrCo
111-152-044	Ø40mm +8.0mm CrCo
111-152-611	Ø28mm -3.5mm Ceramic
111-152-612	Ø28mm 0.0mm Ceramic
111-152-613	Ø28mm +4.0mm Ceramic
111-152-621	Ø32mm -3.5mm Ceramic
111-152-622	Ø32mm 0.0mm Ceramic
111-152-623	Ø32mm +4.0mm Ceramic
111-152-624	Ø32mm +8.0mm Ceramic
111-152-631	Ø36mm -3.5mm Ceramic
111-152-632	Ø36mm 0.0mm Ceramic
111-152-633	Ø36mm +4.0mm Ceramic
111-152-634	Ø36mm +8.0mm Ceramic
111-152-641	Ø40mm -3.5mm Ceramic
111-152-642	Ø40mm 0.0mm Ceramic
111-152-643	Ø40mm +4.0mm Ceramic
111-152-644	Ø40mm +8.0mm Ceramic



Appendix: Compaction Broaching Techniques and Tips

Successful long-term fixation of the Origin cementless hip stem relies on a “sheath” of compacted trabecular bone formed by the compaction broaching system. A successful procedure is equally dependent on the instruments provided in the tray as well as the surgeon’s technique and understanding of the underlying mechanism.

The trabecular sheath is a cavity of dense cancellous bone that is formed line-to-line with the corresponding implant size. By design, the contact stress between compacted trabecular bone and implant is highly controlled down the length of the stem. According to Wolff’s Law, this dense sheath and the trabeculae that support it will strengthen over time in response to applied loads until the stiffness of the sheath closely matches the stiffness of the implant stem. A continuous stiffness gradient at the implant-bone interface minimizes stress shielding and micromotion.

It is important to select an appropriate stem size (and final broach size) to prevent any direct contact between the implant and the true cortical walls. Maintaining axial alignment with the canal from the first broach will make this determination easier. Instruments are included in the Origin case to aid in lateralization where needed. A change in pitch during impaction of the broach is often the first sign of an appropriate amount of compaction. Check for axial stability before checking for rotational stability as this can deform the trabecular sheath.

It is better to stop broaching when stability is achieved with a slightly countersunk broach than to attempt to force an oversized broach into the canal. The HA line on the implant stem will sit below the resection line but the stem will achieve exceptional stability. Consider a longer neck or higher head offset.

Managing Different Femoral Canal Geometries


The Dorr femur type system classifies femurs based on a ratio that relates to the geometry of the femoral canal:

$$R = D_2/D_1$$

see figures at the bottom left


An R-value less than ~0.5 implies a very fluted canal that will more likely bind distally if the canal is not over-reamed. Preoperative templating is especially important for this reason. Refer to the Origin Implant Sizing guide in this technique to help decide what size reamer to use. The three types are as follows.

- Type A $R < 0.5$
- Type B $R = 0.5$ to 0.9
- Type C $R > 0.9$



Type A

Reaming the femoral canal to the distal tip of the definitive stem ensures good compaction of the metaphyseal cancellous bone and prevents binding against the cortical wall, which can lead to stress shielding and an ill-fixed stem.



Types B & C

The Origin broaches preserve the natural anatomy of the femoral canal. Rotational and axial stability and a change in pitch indicate an appropriate amount of compaction broaching.




Signature Orthopaedics
7 Sirius Rd
Lane Cove West, Sydney, 2066
NSW, Australia



Signature
ORTHOPAEDICS



LOGICAL G Ha SERIES CUP



Cup Size (mm)	Head Size (mm)	Optional Head (mm)	Optional Head (mm)
44	28		
46	28		
48	32		
50	32		
52	32	36	
54	32	36	
56	32	36	40
58	32	36	40
60	32	36	40
62	32	36	40
64	32	36	40

Cup Liner

- Clinically proven geometry and material (GUR1020 compression moulded bar stock).
- Available with 10 / 20° hood, or without hood.
- UHMWPE bar stock is cross-linked at 7.5MRads and then remelted. The product is subsequently machined from remelted bar stock and sterilised with ETO.

Cup

- Clinically proven geometry and material (TiAl4V)
- Clinically proven porous coating sprayed with Ha
- Cup available with standard 3-hole and no hole options

HA Coating

- Crystallinity - 68.2%
- HA content - 98.3%
- Ca/P ratio - 1.67
- Density - 99%
- Porosity - <5%
- Pull out strength - 35MPa (on Ti6Al4V)

Porous Coating

- Tensile Strength > 60MPa
- Shear Strength > 30MPa
- Porosity 55-65% and pore size 250-300 microns.



Tender Committee
Orthopaedic Implants Tender
MOLDAVIA

The 3rd of May 2019

Subject : Hip Prosthesis / Range of Femoral Heads / Additional information

Dear Committee members,

In addition to my previous letter dated current April 25th I am pleased to give you hereunder some more information about femoral heads on the market.

Our standard range of CoCr femoral head offers the diameters 28, 32, 36 and 40 with 4 different lengths of neck usually classified as S, M, L and XL. This range is an international standard to cover the requirements of hip arthroplasty and complies with the norms on the market. The entire range is CE marked and also FDA approved.

To allow a wide range of combinations between components of hip prosthesis and then allow a perfect reconstruction of anatomy, the femoral heads could be used with standard stems and lateralized stems.

In addition to this, Signature Orthopaedics propose beyond the standard XLPE liners an entire range of lateralized XLPE liners that could be use in specific cases where surgeon may need this option to restore perfectly the anatomy of patients.

All those products are CE approved and available on demand on a routine basis.

In unlikely case of patient demanding a custom head, like I wrote in my previous letter, Signature Orthopaedics is in position to propose and produce custom head based on surgeon specific. In this case the component could be propose as a custom implant which do no required the CE mark. This is a usual situation for us as we do a lot of custom implants for example for tumor cases; In such situation as Signature Orthopaedics is ISO13485 and complies with the 93/42/EEC annex 2 directive we have agreement to produce custom implants.

With my Best Regards

Olivier PERNOIS

Signature Orthopaedics France

Director

A handwritten signature in black ink, appearing to read 'Olivier Pernois', written over the printed name and title.



Tender Committee
Orthopaedic Implants Tender
MOLDAVIA

The 25th of April 2019

Subject : Hip Prosthesis / Range of Femoral Heads

Dear Committee members,

By the present letter I hereby confirm that in addition to the standard range of metallic femoral heads (made of CoCr), SIGNATURE ORTHOPAEDICS is in position to offer custom heads with other neck length (offsets) to be used with the ARIA stem.

This will be produced in case of special request for specific case received from the user.

The 4 offsets introduced in our leaflet are the main standard used worldwide.

With my Best Regards

Olivier PERNOIS

Signature Orthopaedics France

Director

SIGNATURE ORTHOPAEDICS France SAS

Espace d'Entreprises - L'Arobase

2 Rue Georges Charpak

81100 CASTRES - France

T : +33(0)5 63 73 51 83 - F : +33(0)5 63 73 51 84

RCS CASTRES 815 115 696

N° TVA FR 89 815 115 696



Signature
ORTHOPAEDICS

EC DECLARATION OF CONFORMITY

We

Manufacturer: Signature Orthopaedics Europe Limited
Address: Unit A, IDA Business & Technology Park
Garrycastle, Athlone
N37 DY26 Co. Westmeath
Ireland

hereby declare under our sole responsibility that the following products

Product: Femoral Stems and Heads (refer attached list)
Scope: All production lots

being Class III medical devices according to Commission Directive 2005/50/EC, are placed on the European Market in accordance with the legislation of the United Kingdom transposing the MDD 93/42/EEC as amended by 2007/47/EC, Annex II.

The following harmonised standards have been applied:

ISO 13485:2016 Medical Device – Quality Management System – Requirements for Regulatory Purposes

Notified Body:

TÜV SÜD
Munich, Germany

4 Digit Notified Body Number:

0123

Full Quality Assurance Procedures Certificate: G1 003228 0001 European conformity assessment certificate under Annex II.3 of the Directive 93/42/EEC on Medical Devices

Design Examination Certificate: G7 003228 0007 European conformity assessment certificate under Annex II.4 of the Directive 93/42/EEC on Medical Devices

This Declaration of Conformity is valid until repealed or superseded.

Authorised Signatory:

Signature

Olivier PERNOIS

Name, Position

Attachment: List of Products

10 May 2019

Date



Part Number	Device
111-152-014	Signature -Head 28mm, CoCr size XL
111-152-015	Signature -Head 28mm, CoCr size XXL
111-152-025	Signature -Head 32mm, CoCr size XXL



Product Service

TÜV SÜD Product Service • Octagon House • Concorde Way • Segensworth North • Fareham • PO15 5RL • UK

Choose certainty.
Add value.

To Whom it may concern

Signature Orthopaedics had an ISO 13485:2016 audit booked in November 2018. Due to unforeseen circumstances out of our control, the audit did not go ahead, and the next available audit time has been booked for March 2019.

Two minor NCs were identified in the previous audit in August 2018, and a CAP has been accepted.

There are no concerns regarding the robustness of Signature orthopaedics quality system.

Regards

Paul Jenkins

Operations Manager MHS UK
TUV SUD Product Service
Office: +44 (0)1727 837490
Mobile: +44 (0)7740 235 874

Email address: Paul.Jenkins@tuv-sud.co.uk
www.tuv-sud.co.uk

TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD
 ZERTIFIKAT ♦ CERTIFICATE ♦ 認證證書 ♦ CERTIFICADO ♦ CERTIFICAT




Product Service

Certificate

No. Q1N 002842 0003 Rev. 00

Holder of Certificate: **Signature Orthopaedics Pty Ltd**
 7 Sirius Road
 Lane Cove West NSW 2066
 AUSTRALIA



Scope of Certificate: Design, development and manufacture of cemented and cementless knee and hip systems, bone cement delivery systems, trial implants, spinal implants, osteosynthesis implants, non-resorbable soft tissue fixation implants, radiopaque markers, and related instrumentation including reusable, single use sterile, measuring and those connected to active and non-active medical devices.

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: 75941211

Valid from: 2018-07-20
Valid until: 2019-02-28

Date, 2018-07-20

Stefan Preiß
 Stefan Preiß



TÜV SÜD
 ZERTIFIKAT ♦ CERTIFICATE ♦ 認證證書 ♦ СЕРТИФИКАТ ♦ CERTIFICADO ♦ CERTIFICAT



Certificate

No. Q1N 002842 0003 Rev. 00

Applied Standard(s):

EN ISO 13485:2012 + AC:2012
 Medical devices - Quality management systems
 Requirements for regulatory purposes
 (ISO 13485:2003 + Cor. 1:2009)
 DIN EN ISO 13485:2012
 Upgrade required until 2019-03-31

Facility(ies):

Signature Orthopaedics Pty Ltd
 7 Sirius Road, Lane Cove West NSW 2066, AUSTRALIA

Signature Orthopaedics Europe Limited
 Unit A, IDA Business & Technology Park, Garrycastle, Athlone,
 N37 DY26 Co. Westmeath, IRELAND

-/-

TUV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD
ZERTIFIKAT ◆ CERTIFICATE ◆ 認證書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT



Benannt durch/Designated by
Zentralstelle der Länder
für Gesundheitsschutz
bei Arzneimitteln und
Medizinprodukten
www.zlg.de
ZLG-BS-244.10.08



Product Service

EC Certificate

Full Quality Assurance System
Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)
(Devices in Class IIa, IIb or III)

No. G1 003228 0001 Rev. 00

Manufacturer: Signature Orthopaedics Europe Limited

Unit A
IDA Business & Technology Park
Garrycastle
Athlone
N37 DY26 Co. Westmeath
IRELAND

Facility(ies): Signature Orthopaedics Europe Limited
Unit A, IDA Business & Technology Park, Garrycastle, Athlone,
N37 DY26 Co. Westmeath, IRELAND

Product Category(ies): Cemented and cementless knee
and hip systems, bone cement
delivery system, trials implants,
spinal implants, osteosynthesis implants,
non-resorbable soft tissue fixation implants,
radio opaque markers, and related
instrumentation including reusable,
single use sterile, non active measuring used
during implantation surgery and
those connected to active and
non-active medical devices.

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

Report No.: 75941211_CN_EC

Valid from: 2018-08-06
Valid until: 2020-02-18

Date, 2018-08-06

Stefan Preiß



Wuhu Ruijin Medical Instrument & Devices Co.,Ltd

Add:33rd, East Wanchun Road, Jiujiang District, Wuhu City, Anhui Province, China.

Website:www.whruijin.com

Email:wendy@whruijin.cn

Tel:+86-553-5905318

FAX: +86-553-5223846


Fax:+86-0553-5223846

Medical Electric hip joint drill

Thank you for using our products,we are honored that we can provide good service for your noble work. Before using medical electric drill and saw,please kindly reading the operating manual carefully.

I: Application and Feature:

Medical Electric drill & saw is one of power tools which used in orthopedic operation.Gradually replaced traditional industrial drill and saw,it can be autoclavable up 135°C high temprature and 1.6Mpa high pressure.can be killed some of bacteria and verses which fumigating can not.it also shorten the sterilization, Just under High temperature and high pressure,just 20 minutes,win the valuable time for the patient.the device is optional,and its good looks,small size,light weight,no wires,easy to carry,safe operation.



Picture 1

II: Composition:

1: One Mainframe
2: One charger
3: One battery
4: One drill chuck
5: One sterilizing channel
6: One aluminum alloy box

III: Model Explanation:

Model	Composition
RJ-JZ-I	Handpiece, charger, battery, drill chuck, sterilizing channel, one aluminum alloy box

IV: Main Technical Data:

Hand Piece	Working Voltage	14.4V
	Output Power	$\geq 20W$
	Unit rise in temperature	$\leq 25^{\circ}C$
Charger	Input Voltage	AC 220V, 50Hz
	Output Power	14.4V, 900mA
	Input Power	30VA
Battery	Voltage	14.4V
	Charging time	About 3 hours

V: Technical Data of acetabulum burnishing drill:

Non-load speed	$\geq 300r/min$
Non-load noise	$\leq 63dB (A)$
Diameter run-out	$\leq 0.5mm$



WUHU RUIJIN MEDICAL INSTRUMENT & DEVICE CO., LTD
ADD: 33RD, Wanchun Road, Wuhu Economic District, Wuhu, Anhui Province, China.
PC: 241000 TEL: 86-553-5905318 FAX: 86-553-5223845 WEBSITE: www.whruijin.com/en

Medical Electric Oscillating saw

Application and Feature:

Medical electric swing saw(Picture 1)is used for joint Surgery. Medical Electric Swing saw have been Gradually replaced the Traditional Industrial Drill and saw which only can be fumigating. It can be autoclavable by high temperature and 1.6Mpa High Pressure, can be killed some of bacteria and verses which fumigating can not.it also shorten the sterilization. Just under High temperature and high pressure,just 20 minutes,win the valuable time for the patient.the device is optional,and its good looks, small size, light weight, no wires, easy to carry, safe operation.



Picture 1

Composition:

The surgical power tool is composed of handpiece, charger, battery. The shell is hermetic. The charger is not used for medical instrument.



WUHU RUIJIN MEDICAL INSTRUMENT & DEVICE CO., LTD
 ADD: 33RD, Wanchun Road, Wuhu Economic District, Wuhu, Anhui Province, China.
 PC: 241000 TEL: 86-553-5905318 FAX: 86-553-5223845 WEBSITE: www.whruijin.com/en

Main Technical Feature:


Hand Piece	Working Voltage	14.4V
	Output Power	≥20W
	Unit rise in temperature	≤25°C
Charger	Input Power	AC 220V, 50Hz
	Output Power	14.4V, 700mA
	Input Power	30VA
Battery	Voltage	14.4V
	Charging time	About 3 hours



Technical Data:

Non-load rotating Frequency	≥16000r/min
Non-load noise	≤90dB (A)

Products Classification:

Medical electric saw drill's classification is based on shocking protection, the motor belongs to inner power B type device, charger belongs to II – B type common device; the device can not be used in flammable anaesthesia air & air mix or oxygen mix air.

Marking introduction:  Working condition: Start up from cold status, run 1 minute with load condition, then stop working, after the instrument become cold, then work again.

Marking Introduction: : Class II type instrument;  : B type Instrument

Mark: L: forward R: reversal



Product Service

CERTIFICATE

No. Q2N 16 06 95446 001

Holder of Certificate: Wuhu Ruijin Medical Instrument & Device Co., Ltd.

No. 33, Wanchun Road,
Economic and Technological Development Zone
241000 Wuhu City, Anhui Province
PEOPLE'S REPUBLIC OF CHINA

Facility(ies):

Wuhu Ruijin Medical Instrument & Device Co., Ltd.
No. 33, Wanchun Road, Economic and Technological Development Zone, 241000 Wuhu City, Anhui Province, PEOPLE'S REPUBLIC OF CHINA



Certification Mark:



Scope of Certificate: Production and Distribution of
Medical Electric Saw Drill,
Medical Electric Drill,
Medical Electric Saw,
Saw Blade,
Drill Bit,
Autoclavable Box,
Flexible Reamer

Applied Standard(s):

EN ISO 13485:2012 + AC:2012
Medical devices - Quality management systems -
Requirements for regulatory purposes
(ISO 13485:2003 + Cor. 1:2009)
DIN EN ISO 13485:2012

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system (excluding subclause 7.3), which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: SH16104901

Valid from: 2018-11-08

Valid until: 2019-11-07

Date: 2018-11-08

Sören Frieß
Sören Frieß



Page 1 of 1





Product Service

EC Certificate

Production Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class IIa, IIb or III)

No. G2 16 06 95446 002

Manufacturer: Wuhu Ruijin Medical Instrument & Device Co., Ltd.

No.33, Wanchun Road
Economic and Technological Development Zone
241000 Wuhu City, Anhui Province
PEOPLE'S REPUBLIC OF CHINA



EC-Representative: Wellkang Ltd

Suite B, 29 Harley Street
LONDON
W1G 9QR
UNITED KINGDOM

Product Category(ies): Medical Electric Saw Drill,
Medical Electric Drill,
Medical Electric Saw

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex V. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class IIb and III devices an additional Annex III certificate is mandatory. See also notes overleaf.

Report No.: SH16104901

Valid from: 2016-11-08

Valid until: 2021-11-07



Date: 2016-11-08

Stefan Preis
Stefan Preis

TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

Page 1 of 2