MERCATOR MEDICAL



Medical and protective gloves

Product catalogue



MERCATOR MEDICAL

TABLE OF CONTENTS

Disposable glovesThe meaning of pictograms on the packaging
of diagnostic and protective glovesLatex glovesLatex glovesNitrile glovesSafe Pack systemVinyl glovesHolders – an optimal solutionHousehold and protective glovesHow to put the gloves on?How to take the gloves off?Other products of Mercator MedicalAssortment of medical and protective gloves



DISPOSABLE GLOVES

Single-use gloves are among the most common and popular protective products. Bacteria, mycobacteria, fungi, micro-organisms and viruses can be very easily transferred to your hands from the objects you touch. Observing the basic principles of hand hygiene plus application of gloves selected according to the threat helps to protect the user, the patient and the touched object.

Diagnostic and protective gloves serve both as a medical device and personal protective equipment, and their purpose is to protect our hands, provide the first barrier to bacteria, viruses and other contaminants.

The quality of gloves and their features must comply with very rigorous and detailed specifications indicated in appropriate directives and ordinances end harmonized standards.

Medical gloves – used in the healthcare sector to protect the patient from potential threats.

Conformity assessment under Directive 93/42/EEC and 2007/47/EC on medical devices, as amended.

Class I – Examination gloves. Class I sterile – Sterile examination gloves.

Class IIa – Surgical gloves.

Protective gloves – protect the hands of users against microbiological, mechanical and chemical hazards. Conformity assessment under Ordinance (UE) 2016/425 on Personal Protective Equipment. Category I – for use in low risk environments. Category II – for use in moderate risk environments.

Category III – for use in high risk environments. Type A – the highest protection. Type B – high protection.

Type C – basic protection.



THE MEANING OF PICTOGRAMS ON THE PACKAGING OF DIAGNOSTIC AND PROTECTIVE GLOVES

(2)	single-use product	
NON	non-sterile product	
	the quality of gloves is not guaranteed if the packaging is damaged	
*	protect from sunlight	
Ĵ	protect from moisture	
sc	storage temperatures: 5–35°C	
LATEX	gloves made of natural latex	
LOT	batch code assigned by the manufacturer	
\Box	expiry date	
[m]	date of manufacture	
***	manufacturer, the name and address of the manufacturer, should be entered next to this symbol	,
CE	The mark CONFORMITE EUROPEENNE indicating European conformity. It indicates the product's conformity with all the requirements of the New Approach Directives	
CE 0321	the CE mark with a 4-digit identification number of the notified body indicates that the manufacturing process of the product is subject to independent supervision	
() i	see the instructions for use	
Rï	approved for food contact	
PAP	recyclable – paper	
Ì	the packaging may be treated as municipal waste	
Pictograms	corresponding to the certifications of the EC type exam	inatio
EN 374-1:	003 EN 374-1:2003 EN 374:2003	
?	GLK GLK	
gloves with lin		st

■ THE CCS SYSTEM – COLOUR-CODED SIZE

The CCS system has been designed to help users identify the correct size regardless of the product brand. The coloured bar featured on the unit packaging indicates the specific size of a glove.

Piktogramy	v dodatkowe			
	powdered glc	oves		
	powder-free o	powder-free gloves		
POCIMER	gloves with a	polymer layer in	side	
	gloves with a	skin care layer i	nside	
NITRYL	gloves made	gloves made of nitrile		
VINYL	gloves made	gloves made of polyvinyl chloride		
	additional information inside the packaging			
by weight	indicates the number of items in a packaging by weight			
AQL	AQL (acceptance quality limit) – the glove's level of tightness in accordance with the standards EN 455-2 and EN 374-2			
ASTIM F ANT	protection against viruses and bloodborne pathogens in accordance with the standards ASTM F 1671 and ISO 16604			
ATTA D GOT	protection against cytostatics in accordance with the stan- dard ASTM D 6978			
SHIT FRIENOS	product friendly for the skin			
	electrostatic p EN 1149-1	properties in acc	ordance with th	e standard
Pictograms	correspondir	ng to the certi	fications of th	e EU
ISO 374-1:2016/TYPE A	ISO 374-1-2016/TYPE B	ISO 374-12016/TYPE C	150 374-5:2016	ISO 374-5:2016
TYPE A: examination	TYPE B:	TYPE C:	Protection	Protection



protectio

chemicals

against

bacteria and

funai

protection

against

chemicals

against

chemicals

bacteria, fungi

and viruses



of rubber.





Latex proteins may cause type I allergy, which is why the material is not recommended for people prone to allergic reactions.

ROPERTIES OF THE GLOVES	LATEX
GENERAL - MATERIAL / APPLICATION	••
SKIN-FRIENDLY – CONTAINS NO LATEX PROTEINS CAUSING ALLERGY	•
COMFORTABLE TO USE	•••
STRENGTH/DURABILITY	
TENSILITY/FLEXIBILITY	• • •
PUNCTURE RESISTANCE	• •
TEAR RESISTANCE	• • •
ABRASION RESISTANCE	• •
PROTECTION AGAINST SOLVENTS (ACETONE-BASED)	
PROTECTION AGAINST SOLVENTS (PETROLEUM-BASED)	٠
PROTECTION AGAINST ACIDS AND BASES	• •
PROTECTION AGAINST ALCOHOLS	• •
PROTECTION AGAINST VIRUSES AND BACTERIA (BLOODBORNE PATHOGENS)	•••
BIODEGRADABILITY	•••



Hevea brasiliensis - the tree from which milky latex is extracted and processed in the production

Currently 90% of the world production of natural rubber is located in Asia.

This raw material is widely used in industry, for example for the production of single-use gloves.

Latex gloves provide:

- high flexibility and softness,
- strength and comfortable of use,
- excellent tactile sensitivity, the glove feels like 'second skin'.

Thanks to the natural raw material, latex gloves are biodegradable and environmentally friendly.

● low | ● ● good | ● ● very good

■ EXAMINATION AND PROTECTIVE LATEX GLOVES



try, housework, minor assembly works

comfort® powd		
Glove type	examination, protective, non-sterile	
Colour	cream	
Powdering agent	maize starch (USP)	
External surface	smooth	
AQL	≤ 1.5	
CE classification	 class I – Medical Device category I – Personal Protective Equipment 	
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs	
Distinctive feature	thicker wall – improved barrier properties	
Sizes	XS S M L XL	
Reference number	RD110010 01-05	



Recommended for





Glove type	examination, protective, non-sterile
Colour	blue
Powdering agent	maize starch (USP)
External surface	smooth
AQL	≤ 1.5
CE classification	 class I – Medical Device category I – Personal Protective Equipment
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	colour approved for food contact
Sizes	XS S M L XL
Reference number	RD110040 01-05

LATEX

nursing care, veterinary medicine, beauty and hairdressing salons, food contact, sanitary and cleaning works, furniture indus-

Recommended for

Recommended for



nursing care, veterinary medicine, beauty and hairdressing salons, food contact



Recommended for



intensive care units, outpatient clinics - work with sterile material, dental surgery, catheterization, nursing care, veterinary medicine, beauty salons, food contact

MERCATOR MEDICAL

x ® powde	
Glove type	examination, protective, non-sterile
Colour	cream
Powdering agent	maize starch (USP)
External surface	smooth / textured
AQL	≤ 1.5
CE classification	 class I – Medical Device category I – Personal Protective Equipment
Unit packaging lective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	thinner wall – maximum flexibility
Sizes	XS S M L XL
eference number	RD110100 01-05 / RD110860 01-05

santex[®] sterile



Glove type	examination, protective, sterile
Colour	cream
Powdering agent	maize starch (USP)
External surface	smooth
AQL	≤ 1.5
CE classification	 class I – Medical Device category I – Personal Protective Equipment
Unit packaging lective packaging	1 pair / 50 pairs
Distinctive feature	sterile glove
Sizes	S - XL
Reference number	RC110100 02-05

■ EXAMINATION AND PROTECTIVE LATEX GLOVES



ermagel® coated		
Glove type	examination, protective, non-sterile	
Colour	cream	
Powdering agent	powder free – polymerized	
External surface	textured	
AQL	≤ 1.0	
CE classification	 class I – Medical Device category III – Personal Protective Equipment 	
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs	
Distinctive feature	polymer-coated inner surface – easy donning	
Sizes	XS S M L XL	
Reference number	RD100060 01-05	



Recommended for



Recommended for

 \sim

diagnostics, nursing care, laboratories, contact with selected chemical substances, dentistry, veterinary medicine, beauty and hairdressing salons, food contact, furniture industry, minor assembly works



-	
Glove type	examination, protective, non-sterile
Colour	light yellow
Powdering agent	no
External surface	textured
AQL	≤ 1.5
CE classification	 class I – Medical Device category III – Personal Protective Equipment
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	reduced protein content
Sizes	XS S M L XL
Reference number	RD100050 01-05

LATEX

Recommended for



diagnostics, nursing care, laboratories, contact with selected chemical substances, dentistry, veterinary medicine, beauty and hairdressing salons, food contact, furniture industry, minor assembly works

gh risk HIGH RISK	
examination, protective, non-sterile	
dark blue	
powder free – chlorinated	
textured	
≤ 1.5	
 class I – Medical Device category III – Personal Protective Equipment 	
50 pcs / 10 × 50 pcs	
ambulance [®] high risk: 3 × thicker than the stand- ard glove, long cuff (at least 300 mm) – the highest level of safety	
S M L XL	
RD100110 02-05	



rubber.

Nitrile gloves:

PROPER

SKIN PROTE PROTECT



Nitrile, also known as acrylonitrile butadiene rubber (NBR) is a material composed of organic chemical compounds. The production process of this raw material was modified so that it can compete with natural

Nitrile gloves are an ideal alternative to gloves made of a natural material.

• do not contain latex proteins,

 their content of vulcanization accelerators is limited to a minimum.

Minimized content of chemical accelerators s nificantly reduces the risk of type IV allergy. Minimized content of chemical accelerators sig-

RTIES OF THE GLOVE	NITRYL
GENERAL – MATERIAL / APPLICATION	•••
N-FRIENDLY – CONTAINS NO LATEX PROTEINS CAUSING ALLERGY	•••
COMFORTABLE TO USE	•••
STRENGTH/DURABILITY	• •
TENSILITY/FLEXIBILITY	• •
PUNCTURE RESISTANCE	•••
TEAR RESISTANCE	• •
ABRASION RESISTANCE	•••
ECTION AGAINST SOLVENTS (ACETONE-BASED)	••
TION AGAINST SOLVENTS (PETROLEUM-BASED)	•••
PROTECTION AGAINST ACIDS AND BASES	•••
PROTECTION AGAINST ALCOHOLS	•••
PROTECTION AGAINST VIRUSES AND BACTERIA (BLOODBORNE PATHOGENS)	•••
BIODEGRADABILITY	٠

● low | ● ● good | ● ● very good


■ EXAMINATION AND PROTECTIVE NITRILE GLOVES



itrylex® classic 50		
Glove type	examination, protective, non-sterile	
Colour	white	
Powdering agent	no	
External surface	microtexture with additional texture at finger tips	
AQL	≤ 1,0	
CE classification	 class I – Medical Device category III – Personal Protective Equipment 	
Unit packaging / collective packaging	50 pcs / 20 × 50 pcs	
Distinctive feature	small, light and convenient packaging	
Sizes	XS S M L XL	
Reference number	RD301740 01-05	



Recommended for





nitrylex [®] classic		
Glove type	examination, protective, non-sterile	
Colour	blue / white / violet	
Powdering agent	по	
External surface	fingertip texture	
AQL	≤ 1.0	
CE classification	class I – Medical Device category III – Personal Protective Equipment	
Unit packaging / collective packaging	■ 100 pcs / 10 × 100 pcs ■ 200 pcs / 10 × 200 pcs	
Distinctive feature	a wide range of different colours	
Sizes	XS S M L XL	
Reference number	a'100 RD300190 01-05 RD301690 01-05 RD301430 01-05 a'200 RD300960 01-05 RD301680 01-05 RD301680 01-05 RD300970 01-05	

Recommended for

Recommended for


diagnostics, nursing care, laboratories, cytostatic laboratories, contact with selected chemical substances, dentistry, veterinary medicine, beauty and hairdressing salons, food contact, minor assembly works

isolated patients, transport carts, patient's bed during hospital transport, patient care after chemotherapy, public health nurses

and midwives, motorcycle ER, mobile beauty professionals - hairdressers, beauticians, housework



Recommended for



MERCATOR MEDICAL

nitrylex® grip



Glove type	examination, protective, non-sterile
Colour	blue
Powdering agent	no
External surface	fully textured
AQL	≤ 1.5
CE classification	 class I – Medical Device category III – Personal Protective Equipment
Unit packaging lective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	$2 \times \mbox{thicker}$ than the standard glove, textured on the whole surface
Sizes	XS S M L XL
Reference number	RD300880 01-05

x ® collag	en	PREMIUM	NITRYL
Glove type	examination protection	ve non-sterile	

Glove type	examination, protective, non-sterile
Colour	magenta
Powdering agent	no
External surface	fingertip texture
AQL	≤ 1.0
CE classification	 class I – Medical Device category III – Personal Protective Equipment
Unit packaging lective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	skin care layer inside improves condition of the skin, nice and distinguishing colour
Sizes	XS S M L XL
Reference number	RD300950 01-05

EXAMINATION AND PROTECTIVE NITRILE GLOVES

nursing care, injection, collection and preparation of biological material for analytical research,

cytostatic, contact with selected chemical substances, emergency medical services, minor assembly works

feeding the patient, administering medication, providing wounds, changing the dressing, procedures related to vascular ac-

cess, patient preparation for surgery, dentistry, laboratories, veterinary medicine, beauty salons, contact with food, workshops



rylex [®] complete		
Glove type	examination, protective, non-sterile	
Colour	lavender	
Powdering agent	no	
External surface	microtexture with additional texture at finger tips	
AQL	≤ 1.0	
CE classification	 class I – Medical Device category III – Personal Protective Equipment 	
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs	
Distinctive feature	maximum flexibility and longer cuff (at least 260 mm), well-arranged gloves, two vertical perfo- ration openings thanks to which gloves can be dispensed one by one cuff first	
Sizes	XS S M L XL	
Reference number	RD 301020 01-05	

nitrylex' chemo long MERCATORINEDICAL / col

Recommended for

nitrylex' beFree



nitrylex high risk REACATORIN

nitrylex [®] high ri	isk high risk 🔊 Mitryl
Glove type	examination, protective, non-sterile
Colour	orange
Powdering agent	no
External surface	microtexture with additional texture at finger tips
AQL	≤ 1.5
CE classification	 class I – Medical Device category III – Personal Protective Equipment
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	$2.5 \times$ thicker than the standard glove, longer cuff (at least 270 mm) – the highest degree of safety
Sizes	XS S M L XL
Reference number	RD300010 01-05

Recommended for

Recommended for

 \checkmark



high-risk procedures, central sterile supply departments (C.S.S.D), admission rooms, hospital emergency departments, rectal examination, uniformed services, diagnostics, nursing care, laboratories, cytostatic laboratories, contact with selected chemical substances, medical rescue services, veterinary medicine, food contact, dentistry, car repairs and automotive industry, printing and electronic industry

MERCATORMEDICA / col

Recommended for

	oncology wards, neonatal wards, paediatric wards, geriatric war	
	\sim	care, dental practices, beauty salons, protection of hands of peo

MERCATOR MEDICAL

14

nitrylex[®] chemo long



Glove type	examination, protective, non-sterile
Colour	cornflower blue
Powdering agent	no
External surface	microtexture with additional texture at finger tips
AQL	≤ 1.0
CE classification	 class I – Medical Device
CE Classification	 category III – Personal Protective Equipment
Unit packaging lective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	additional forearm protection, distinctive colour
Sizes	XS S M L XL
Reference number	RD300930 01-05

nitrylex® beFre		
Glove type	examination, protective, non-sterile	
Colour	blueberry	
Powdering agent	no	
External surface	microtexture with additional texture at finger tips	
AQL	≤ 1.0	
CE classification	 class I – Medical Device 	
	 category III – Personal Protective Equipment 	
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs	
	free of rubber accelerators and sulphur (chemical	

Distinctive feature	accelerators) – ideal for people extremely sen- sitive to allergens, product friendly for the skin	
Sizes	XS S M L XL	
Reference number	RD300830 01-05	

ards, cytostatic laboratories – in the non-sterile phase, nursing eople with allergy problems.

■ EXAMINATION AND PROTECTIVE NITRILE GLOVES I



printing and electronic industry, minor assembly works, housework

itrylex® basic	NITRYL
Glove type	examination, protective, non-sterile
Colour	navy blue
Powdering agent	no
External surface	fingertip texture
AQL	≤ 1.5
CE classification	 class I – Medical Device category I – Personal Protective Equipment
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs 200 pcs / 10 × 200 pcs
Distinctive feature	anti-static properties
Sizes	XS S M L XL
Reference number	a'100 RD301050 01-05 a'200 RD300840 01-05



Recommended for





trylex® green		
Glove type	examination, protective, non-sterile	
Colour	mint green	
Powdering agent	no	
External surface	fingertip texture	
AQL	≤ 1.5	
CE classification	 class I – Medical Device category III – Personal Protective Equipment 	
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs	
Distinctive feature	characteristic colour	
Sizes	XS S M L XL	
Reference number	RD301380 01-05	

Recommended for

Recommended for



paediatric wards, beauty and hairdressing salons, dentistry, diagnostics, nursing care, laboratories, contact with selected chemical substances, veterinary medicine, food contact, minor assembly works, housework

nursing care, dentistry, veterinary medicine, beauty and hairdressing salons, food contact, mechanical and automotive industry,



Recommended for

	tattoo artists, beauty parlors, hairdressing salons, mechanical
\checkmark	veterinary medicine, beauty and hairdressing salons, food conta
	housework

16

ex ® pink			
Glove type	examination, protective, non-sterile		
Colour	light pink		
Powdering agent	no		
External surface	fingertip texture		
AQL	≤ 1.5		
CE classification	 class I – Medical Device category III – Personal Protective Equipment 		
Unit packaging lective packaging	100 pcs / 10 × 100 pcs		
Distinctive feature	characteristic colour		
Sizes	XS S M L XL		
Reference number	RD301440 01-05		



nitrylex® black

Glove type examination, protective, non-sterile		
Colour	black	
Powdering agent	no	
External surface fingertip texture		
AQL	QL • <1,5	
CE classification	 class I – Medical Device category I – Personal Protective Equipment 	
Unit packaging lective packaging	100 pcs / 10 × 100 pcs	
Distinctive feature	opaque colour	
Sizes	XS S M L XL	
Reference number	RD301040 01-05	

al workshops, cleaning personnel, nursing care, dentistry, ntact, car repairs and automotive industry, minor assembly works,

SAFE PACK SYSTEM

Innovative Safe Pack dispenser system



Special folding method which allows to dispense a single item. Each time a glove is pulled by the cuff, another glove is prepared, which minimizes the risk of contamination of clean gloves.

Why choose the SAFE PACK system?



Reduces contamination.

This is achieved because each glove is dispensed by cuff from the bottom of the package.

Minimizes the risk of cross-infections.

There is no contact between the user's hand with the working surface of the glove or with the box.

Helps to organize the workplace efficiently. Mounted to the wall, it is always at hand and in one specific location.

Easy disinfection and cleaning of holders. Lightweight design, made of plastic or steel wire covered with white anti-bacterial paint

Focus on hygiene and safety



The SAFE PACK glove dispenser system is recommended by the Polish Association of Epidemiological Nurses Due to the risk of contaminating clean gloves when they are dispensed from the collective packaging, it is recommended that solutions to prevent this should be applied.



Recommended for



diagnostics, nursing care, laboratories, cytostatic laboratories, contact with selected chemical substances, dentistry, veterinary medicine, beauty and hairdressing salons, food contact, minor assembly works

nitrylex




Recommended for



outpatient clinics – work with sterile material, dental surgery, catheterization, medical rescue services, diagnostics, nursing care, especially in ICUs, laboratories, dentistry, veterinary medicine, beauty salons, food contact

nitrylex® one by one



(® sterile STERILE EO		NITRYL	
-----------------------	--	--------	--

Glove type	examination, protective, sterile	
Colour	blue	
Powdering agent no		
External surface	fingertip texture	
AQL	≤ 1.5	
CE classification	 class I – Medical Device category I – Personal Protective Equipment 	
Unit packaging lective packaging	1 pair / 50 pairs	
Distinctive feature	sterile glove	
Sizes	S - XL	
Reference number	RC300150 02-05	



0

and IV allergy.

PROPER

SKIN PROTE PROTECT Р



Vinyl belongs to the group of PVC materials. Thanks to its properties, polyvinyl chloride (PVC), has a broad scope of application in many industries.

Vinyl gloves are skin friendly for users prone to type I

RTIES OF THE GLOVE	VINYL
GENERAL – MATERIAL / APPLICATION	٠
N-FRIENDLY – CONTAINS NO LATEX PROTEINS CAUSING ALLERGY	•••
COMFORTABLE TO USE	٠
STRENGTH/DURABILITY	٠
TENSILITY/FLEXIBILITY	٠
PUNCTURE RESISTANCE	٠
TEAR RESISTANCE	٠
ABRASION RESISTANCE	•
ECTION AGAINST SOLVENTS (ACETONE-BASED)	•
TION AGAINST SOLVENTS (PETROLEUM-BASED)	٠
PROTECTION AGAINST ACIDS AND BASES	٠
PROTECTION AGAINST ALCOHOLS	٠
PROTECTION AGAINST VIRUSES AND BACTERIA (BLOODBORNE PATHOGENS)	•
BIODEGRADABILITY	•

● low | ● ● good | ● ● very good


EXAMINATION AND PROTECTIVE VINYL GLOVES

Turker I	vinylex* powdered	
3		M
e	MERCATOR MEDICAL	ţ=
	MERCENCER CONTRACT	
	-	
	2	No.

vinylex [®] powdered		
Glove type	examination, protective, non-sterile	
Colour	transparent	
Powdering agent	yes	
External surface	smooth	
AQL	≤ 1.5	
CE classification	class I – Medical Device category III – Personal Protective Equipment	
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs	
Distinctive feature	powdered gloves – easier donning	
Sizes	XS S M L XL	
Reference number	RD210150 01-05	



Recommended for





-	
Glove type	examination, protective, non-sterile
Colour	transparent blue
Powdering agent	yes
External surface	smooth
AQL	≤ 1.5
CE classification	 class I – Medical Device category I – Personal Protective Equipment
Unit packaging / collective packaging	100 pcs / 10 × 100 pcs
Distinctive feature	blue colour – approved for contact with food
Sizes	XS S M L XL
Reference number	RD210160 01-05

Recommended for

Recommended for



veterinary medicine, beauty and hairdressing salons, food contact, sanitary and cleaning works

veterinary medicine, beauty and hairdressing salons, food contact, sanitary and cleaning works



vinylex[®] powder-free

Glove type examination, protective, non-sterile	
Colour transparent	
Powdering agent	no
External surface smooth	
AQL	≤ 1.5
CE classification	 class I – Medical Device (Council Directive 93/42/EEC) category III – Personal Protective Equipment (Council Directive 89/686/EEC)
Unit packaging lective packaging	100 pcs / 10 × 100 pcs
istinctive feature no powdering agent – reduced risk of allergic reactions	
Sizes	XS S M L XL
Reference number	RD200180 01-05

HOLDERS – AN OPTIMAL SOLUTION



Designed with the users' needs in mind



Wall holders minimize the risk of contamination of gloves and improve the productivity of work thanks to its better organization.



A wide range of holders types allows the use of **most** examination gloves available on the market according to individual needs.



Metal holders are made of steel wire with a mounting plate and they are covered with white paint with enhanced resistance, which significantly increases their durability.



The optimal design of the wall holders allows easy cleaning and disinfection of their surface.



Wall holders for mounting on flat surfaces with screws or mounting tape without drilling.




Holders suspend the packaging with gloves at the patient's bed (during transport) or at the treatment chair. The system allows access to gloves in the selected place.



Select appropriate wall holders, which allow you to dispense gloves easily and fast during work.

Basic metal brackets

Brackets made of steel wire with a mounting plate; the whole item coated with anti-bacterial paint.



Material	steel	steel
Distinctive feature	 allows to save the working space at the workplace KEUR001 – possibility to attach trans- port hooks. 	 KEUR001 br the bracket to hang on t cart.
Ref. number	KEUR001 / KEUR004	KEUR005

■ SAFE PACK holders

Brackets for gloves dispensed one by one, arranged in the "Safe Pack" bottom dispensing system.



	Safe Pack single / triple holder	Safe Pack single / triple holder	Safe Pack holder for the Modura rail
Material	plastic	steel	steel
Distinctive feature	allows to save the working space,reduces the risk of contamination.	 allows to save the working space, reduces the risk of contamination.	no installation or drilling required,may be hung on medical columns.
Ref. number	KEUR010	KEUR007 / KEUR008	KEUR012

Together, let's ensure a proper place

HOUSEHOLD AND PROTEC-**TIVE GLOVES**

health. Reusable.

weak detergents.

- work.

Gloves coated with cotton lining, which increases the product's durability, mechanical resistance and practical features of the gloves.

Rubber household gloves have an extended cuff -300 mm, which protects not only hands, but forearms as well.





Household and protective gloves protect hands from external factors that are undesirable or harmful to

Protection against harmful effects of cleaners and

• Protection against chafes and minor cuts during

 Increased safety at work in contact with slippery and wet objects.

> Protective gloves are used for household works, but also in light industry or professional services.

HOUSEHOLD AND PROTECTIVE GLOVES





ideall [®] cotton			
Glove type	household, protective, non-sterile		
Colour	yellow		
Internal surface	soft cotton lining		
External surface with projections on the gripping side			
• Category I – Personal Protective Equipmen (Council Directive 89/686/EEC)			
Unit packaging / collective packaging	1 pair / 240 pairs		
Sizes	S M L XL		
Reference number	RG100070 02-05		



1	
w	
28	
M	

ideall® yellow			
Glove type	household, protective, non-sterile		
Colour	yellow		
Internal surface	flock-lined		
External surface	with projections on the gripping side		
CE classification	 Category I – Personal Protective Equipment (Council Directive 89/686/EEC) 		
Unit packaging / collective packaging	1 pair / 240 pairs		
Sizes	S M L XL		
Reference number	RG100040 02-05		



ideall® orange	
Glove type	household, protective, non-sterile
Colour	orange
Internal surface	flock-lined
External surface	with projections on the gripping side
CE classification	 Category I – Personal Protective Equipment (Council Directive 89/686/EEC)
Unit packaging / collective packaging	1 pair / 240 pairs
Sizes	S M L XL
Reference number	RG100470 02-05

Recommended for




sanitary and cleaning works in public places, cleaning works at facilities, utility rooms and at home, professional cleaning services, refurbishment works, gardening, agriculture, food processing and food contact, mechanical works, housework

MERCATOR MEDICAL



REMEMBER! BEFORE TAKING A GLOVE FROM THE BOX, PERFORM APPROPRIATE HAND HYGIENE ACTIVITIES.

■ HOW TO PUT THE GLOVES ON? I








Take another glove from the box with the hand wearing a glove.



MERCATOR MEDICAL









REMEMBER! ONCE YOU HAVE TAKEN THE GLOVES OFF AND DISPOSED THEM AS HEALTHCARE WASTE, REPEAT THE HAND HYGIENE ACTIVITIES.

OTHER PRODUCTS OF MERCATOR MEDICAL



When it comes to protection against infections – trust experts.

Use structured approach to prevent infections by choosing:

- medical gloves: examination and protective, surgical, household gloves,
- caps, masks, aprons and overalls,
- operative sets with drapes and necessary accessories used to secure the operative field during surgical procedures,
- medical dressings.



Nitrile gloves for children

Ask about the product: (54) 416 96 00 / 800 633 422



Mercator Medical S.A.

Company Head Office ul. Heleny Modrzejewskiej 30 31-327 Kraków tel. +48 126 655 400 fax. +48 126 655 415 e-mail: medical@mercatormedical.eu

C

Sales Office and Logistics Center Stary Brześć 25E 87-880 Brześć Kujawski tel/fax (054) 416 96 00

nit

nitrylex¹

e-mail: bok@pl.mercatormedical.eu

