TAVI Analysis

Automated segmentation & pre-procedure planning software for transcatheter aortic valve implantation.

Transcatheter Aortic Valve Implantation (TAVI) or Replacement (TAVR) is a procedure for select patients with severe symptomatic aortic stenosis (narrowing of the aortic valve opening). TAVI/TAVR is performed on a beating heart and a replacement valve is delivered via a catheter.

CT provides important information that is important for successful TAVI/TAVR procedures. CT is used to help determine aortic annulus size, to guide selection of appropriate valve, provide dimensions of the entire aorta, help determine the access path for the catheter and give guidance for C-arm angulation for deployment of the device.

GE's TAVI Analysis software provides a streamlined, guided workflow to enable efficient, consistent work-ups of your TAVI studies with connectivity directly to the interventional suite.

Overview

TAVI Analysis is a post-processing software package for the Advantage Workstation (AW) and AW Server platforms. It is a planning tool used for TAVI/TAVR procedures. It automatically segments the aorta and displays the aortic valve in multiple views for quick and easy measurements of the annulus. TAVI Analysis provides guided workflow and semi-automated tools to help evaluate appropriate access pathways and can communicate directly with the interventional suite.

Highlights

- Automated Aortic segmentation displayed to highlight calcific areas.
- Guided workflow for complete pre-procedure planning measurements.
- Wizard panel for valve measurements used for sizing requirements.
- Perpendicular view demonstrates working angles for valve deployment.
- Auto bone segmentation
- Fast centerline tracking to subclavian and femoral arteries.
- Lumen and curved views showing entire vascular tree.
- Quantification wizard panel.
- Intuitive and easily accessible tools
- Connectivity to the interventional suite

Visit us:

www.gehealthcare.com/aw/applications/

SECRET DE
AFACERI

Features

- One comprehensive workflow for valve measurements and access path measurements.
- Automatic segmentation of transparent volume rendered model of the aorta with calcifications highlighted
- Auto detection of valve plane.
- Intuitive workflow for acquiring the necessary measurements for aortic valve sizing
- One click perpendicular view provides guidance for C-arm angulation for deployment of device.
- Automatic bone segmentation for the chest, abdomen and pelvis.
- Real time vessel tracking provides fast centerline detection based on a single start point.
- Guided vessel tracking tools allow for easy access planning for all transcatheter approaches:
 - Subclavian Artery
 - Femoral Artery
 - Apical Approach
- Peripheral artery measurements automatically sent to summary table to allow easy exporting.
- Save State and send to the Cath lab for sharing access to models and measurements.

Indications for Use

TAVI Analysis for CT

TAVI Analysis is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

TAVI Analysis is a post processing application option for the Advantage Workstation (AW) platform, CT Scanner or PACS stations, which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and directional vessel tortuosity visualization. TAVI Analysis automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel, analyzing calcified and non-calcified plaque to determine the densities of plaque within a coronary artery, measure areas of abnormalities within a vessel.

TAVI Analysis for X-Ray

TAVI Analysis is a software post processing option for the Advantage Workstation (AW) Platform, which can be used in the analysis of the 3D angiography data. It provides a number of display, measurement and batch filming/archive features and will aid Physicians in studying user selected vessels for stenosis analysis, pre & post stent planning, and directional vessel tortuosity visualization.

System Requirements

Minimum platform release:

- AW Workstations VolumeShare 7 and AW Server 3.1 or later.
- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.


SECRET DE
AFACERI

Features Detail

Valve Definition


- The TAVI protocol allows users to select single or multiple cardiac phases.
- Automatic segmentation of aorta and auto detection of valve plane with double oblique tool.
- Ability to define minimum and maximum diameters above the annulus along the length of the aorta.

- Analyze, trace and edit annulus contour. Available annulus measurements include:
 - Minimum diameter
 - Maximum diameter
 - Perimeter
 - Area
 - Equiv. diameter based on perimeter and area

- Aortic measurement tools are also available which include:
 - Diameters above annulus
 - Ostia Heights
 - Sinus diameters
- One click perpendicular view demonstrates working angles for valve deployment in the cath lab.


Aorta Analysis

- Automatic bone segmentation and real time vessel tracking.
- Guided vessel tracking tools allow for easy access planning for all transcatheter approaches:
 - Subclavian Artery
 - Femoral Artery
 - Apical Approach

- Quantification view has dedicated tools capable of doing measurements for stenosis, aneurysms, 3D angles and distance measurements.
- Predefined measurements list provides a template for a standard order and location from which to place all measurements.

- Peripheral artery measurements automatically sent to summary table to allow easy exporting.

Phase	Measurement	Value
CT - phase 495	Min Diameter	22.6 ± 0.5 mm
CT - phase 495	Max Diameter	22.6 ± 0.5 mm
CT - phase 495	Perimeter	86.9 ± 0.7 mm
CT - phase 495	Area	77.1 ± 0.8 mm²
CT - phase 495	Equiv. Diameter	22.6 ± 0.5 mm
CT - phase 495	Min Diameter	22.6 ± 0.5 mm
CT - phase 495	Max Diameter	22.6 ± 0.5 mm
CT - phase 495	Perimeter	86.9 ± 0.7 mm
CT - phase 495	Area	77.1 ± 0.8 mm²
CT - phase 495	Equiv. Diameter	22.6 ± 0.5 mm

- Summary table converts to PDF format for a simple and concise way to share measurements and images with heart team for planning discussions.

Procedure Preparation

- Provides valve measurements, perpendicular view and gantry angulation in the cath lab to prepare for TAVI/TAVR procedure.
- Direct communication with Heart Vision software for easy transition from planning to guidance.
- 3D and calcium overlay VR models to aid in visualization during procedure.

Indications for Use

TAVI Analysis for CT

TAVI Analysis is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

TAVI Analysis is a post processing application option for the Advantage Workstation (AW) platform, CT Scanner or PACS stations, which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and directional vessel tortuosity visualization.

TAVI Analysis automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel, analyzing calcified and non-calcified plaque to determine the densities of plaque within a coronary artery, measure areas of abnormalities within a vessel.

TAVI Analysis for X-Ray

TAVI Analysis is a software post processing option for the Advantage Workstation (AW) Platform, which can be used in the analysis of the 3D angiography data. It provides a number of display, measurement and batch filming/archive features and will aid Physicians in studying user selected vessels for stenosis analysis, pre & post stent planning, and directional vessel tortuosity visualization.

System Requirements

Minimum platform release:

- AW Workstations VolumeShare 7 and AW Server 3.1 or later.
- Color Monitors
- Single or Dual Display Monitors for AW Server in landscape or portrait orientations

Regulatory Compliance

This product complies with the European CE marking regulation following Medical Devices Directive: Directive 93/42/EEC.

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company.
DOC1471986

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HeartVision 2

Structural Heart 3D guidance

Guiding and deploying precisely devices during structural heart procedures may be delicate. Under 2D fluoroscopic imaging, the visualization of structures such as the aortic annulus plane may be challenging. The overlay of the 3D anatomy and enhanced visualization of calcified structures could help you during these complex interventional procedures.

Overview

HeartVision 2 overlays prepared 3D datasets on live fluoroscopy to support guidance and deployment of devices such as the one used during structural heart procedures. The Calcification Visualization Enhancement is a mode that allows the user to enhance moving contrasted structures in the image to improve the visualization.

What's new

- **Intuitive user interface, all at table side.** Control HeartVision 2 at table side through an intuitive user interface. Register your anatomy through a simple guided workflow.
- **Bi-view mode eases the registration.** Bi-view registration allows the user to fuse pre-op CT or MR data with the live fluoroscopic images with an accuracy better than 2.5mm (1).
- **Enhanced calcifications visibility.** By suppressing nonmoving contrasted structures, the Calcification Visualization Enhancement mode allows the user to enhance the Contrast to Noise Ratio of moving contrasted structures by a factor of 2.6 (2).
- **New 3D Contour rendering** Improved intra-aorta visualization on the fused image.

Visit us:

www.gehealthcare.com/heartvision2

SECRET DE
AFACERI

Intuitive user interface

HeartVision 2 can be controlled directly on the AW or at table side through the Central Touchscreen (Fig 1). The user can:

- Launch and close the application
- Use the Bi-View registration mode
- select of the 3D objects to be displayed
- Optimize the display of the 3D Overlay
- Review and record sequences with their 3D overlay

Fig 1

Bi-view registration

- The Bi-View registration mode allows achieve a proper 3D registration between the anatomy seen under fluoroscopy and the 3D model
- It helps the user manage two images acquired under different angulations and then optimize the 3D registration in translation and rotation. An accuracy score (Fig 2) is provided to help increase registration confidence. Correct patient motion at table side or in control room.

Fig 2

Display controls

Optimize overlaid objects display:

- Opacity of the 3D overlay
- 2 rendering modes (Fig 3): Surface rendering or 3D Contours
- Color controls
- Front/ back view
- Displayed Hounsfield Units threshold
- Show/ hide landmarks and planning lines
- Zoom/ roam the displayed images
- Show/ hide ECG display


Fig 3

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Calcification Visualization Enhancement

The Calcification Visualization Enhancement mode can be activated. It enhances moving contrasted structures in the image while suppressing nonmoving contrasted structures (Fig 4).

Without

Fig 4

With

GE imagination at work

HeartVision 2

SECRET DE
AFACER!

Intended Use

HeartVision 2 software application is intended to enable users to load 3D datasets and overlay and register in real time these 3D datasets with radiosopic or radiographic images of the same anatomy in order to support catheter/device guidance during interventional procedures.

Regulatory Compliance

HeartVision 2 complies with the European Medical Devices Directive requirements for CE marking: Directive 93/42/EEC.

In non CE mark countries, HeartVision 2 cannot be placed on the market or put into service until it has obtained all required regulatory authorizations. HeartVision 2 is cleared by the U.S. Food and Drug Administration for commercial availability in the U.S.

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Disclaimer

(1) For systems with Innova-IQ table or Omega V table and with augmented calibration, rigid geometric phantom, 3D volume with voxel size of 0.95mm, frontal plane, no L-arm or table motion after registration, Region of interest of 10cm.


(2) Average value, measured on experimental set-up including a chest phantom and moving aortic calcification phantom, using Calcification Visualization Enhancement feature at level 100%.

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, Innova and Discovery are trademarks of General Electric Company.

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Vision 2

3D live guidance

Guiding and deploying devices through complex anatomies may be delicate. Under 2D fluoroscopic imaging, such procedures can require significant amount of contrast media, prolonged radiation exposure, and lot of time. A 3D view of vessels anatomy could help you perform complex interventional procedures.

Overview

Vision 2 overlays prepared 3D datasets on live fluoroscopy to support localization and guidance of catheters, coils and other devices during interventional procedures.

What's new

- **Intuitive user interface, all at table side.** Control Vision 2 at table side through an intuitive user interface. Register your anatomy through a simple guided workflow
- **Register the 3D anatomy with accuracy.** Vision 2 allows the user to automatically fuse CBCT data over the live fluoroscopic images with an accuracy better than 1.8 mm (1). Vision 2 allows the user to visualize mis-registrations between CBCT data and live fluoroscopic images and correct them at any time in both translation and rotation with an accuracy better than 1 mm (2).
- **Bi-view mode eases the registration.** Bi-view registration allows the user to fuse pre-op CT or MR data with the live fluoroscopic images with an accuracy better than 2.5mm (3).
- **Radiation dose savings during registration.** Vision 2 helps achieve up to 99% radiation dose savings when fusing a pre-operative 3D model over fluoroscopy (4).

SECRET DE
AFACER!

Intuitive user interface

Vision 2 can be controlled directly on the AW or at table side through the Central Touchscreen (Fig 1). The user can:

- Launch and close the application
- Use the Bi-View registration mode
- select of the 3D object to be displayed
- Optimize the display of the 3D Overlay
- Review and record sequences with their 3D overlay

Fig 1

Bi-view registration

- The Bi-View registration mode allows achieve a proper 3D registration between the anatomy seen under fluoroscopy and the 3D model
- It helps the user manage two images (Fig 2) acquired under different angulations and then optimize the 3D registration in translation and rotation. An accuracy score (Fig 3) is provided to help increase registration confidence. Correct patient motion at table side or in control room.


Fig 2


Fig 3

Display controls

Optimize overlaid objects display:

- Opacity of the 3D overlay
- 3 rendering modes (Fig 4): Digitally Reconstructed Radiography, Volume rendering or 2D Contours
- Color controls
- Front/ back view
- Displayed Hounsfield Units threshold
- Show/ hide landmarks and planning lines
- Show/ hide vessel centerlines
- Zoom/ roam the displayed images


Fig 4

Review controls

The review controls allow the user to:

- Review stored Vision sequences
- Review one of the last ten acquired sequences
- Create movies
- Store photos of fused sequences.

Movies and photos are stored on the AW as DICOM objects.

The Auto Store function allows store automatically any new sequence of acquired fused images.

Biplane

On biplane systems, Vision 2 allows the user to automatically switch the 3D model display on either the frontal or lateral views based on pedal press. Frontal or biplane pedal press activates the 3D overlay on the frontal view. Lateral pedal press activates the 3D overlay on lateral view.

GE imagination at work

Vision 2

SECRET DE
AFACENI

Intended Use

Vision 2 software applications is intended to enable users to load 3D datasets and overlay and register in real time these 3D datasets with radioscopy or radiographic images of the same anatomy in order to support catheter/device guidance during interventional procedures.

Regulatory Compliance

Vision 2 complies with the European Medical Devices Directive requirements for CE marking: Directive 93/42/EEC.

In non CE mark countries, Vision 2 cannot be placed on the market or put into service until it has obtained all required regulatory authorizations. Vision 2 is cleared by the U.S. Food and Drug Administration for commercial availability in the U.S.

System Requirements

- Advantage workstation (AW) at the AW 4.7 (VolumeShare 7) level or above and dedicated Volume Viewer products.
- One of the following X-ray system: Innova™ IGS 520, Innova IGS 530, Innova IGS 540, Innova IGS 620, Innova IGS 630, Discovery™ IGS 730 or Discovery IGS 740
- The following options shall be activated as well:
 - The AW video shall be displayed in the exam room either with a

dedicated AW monitor or using the X-ray system's Large Display Monitor.

- A point-to-point Ethernet connection between the X-ray Vascular interventional system and the Advantage Workstation product is required.
- A Central Touchscreen shall be installed with the IGS X-Ray system.

Disclaimer

(1) For systems with Innova-IQ table or Omega V table, rigid geometric phantom, frontal plane, L-arm 0°, table motion inferior to 80mm between 3D acquisition and live fluoroscopy, accuracy obtained for at least 90% of the images acquired on each system.

(2) For systems with Innova-IQ table, Omega V table. Using rigid geometric phantom, CBCT data, frontal plane, L-arm 0°, Region of interest of 10cm.

(3) The accuracy is determined by engineering analysis using rigid phantom data having a 3D voxel size with no dimension greater than 0.95 mm. This accuracy is obtained with 2 fluoroscopic images taken at optimal angulation, with system calibrated with "augmented calibration", without motion of L-arm and table between the 2 sequences and after

validation of the registration.

Accuracy is measured on the frontal plane in Region of interest of 10 cm.

(4) Based on the dose of the procedure step needed using a CBCT acquisition to register pre-operative data vs. a Bi-View registration process. The stated dose reduction does not reflect the entire interventional procedure, rather to a specific step in the procedure.


The dose for the CBCT acquisition is from typical exposure settings (Innova CT 40°/s, 30fps, IQ Standard, Normal, Nominal FOV). The dose from the Bi-View registration process is from two spatially separated, 2-seconds fluoroscopic acquisitions, with typical exposure settings (3.75 fps, IQ Standard, Normal, Max Dose Reduction, Nominal FOV). The dose data for all acquisitions are from the Air Kerma per IEC 60601-2-43 conditions, provided in the interventional X-Ray user manual. In clinical practice, the use of Vision 2 may reduce patient radiation dose depending on the clinical task, patient size, anatomical location and clinical practice.

GE imagination at work

© 2016 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, Innova and Discovery are trademarks of General Electric Company.

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Integrated Registration

Multi-Modality Image Management

In Radiology practice today, multiple complementary imaging modalities are frequently relied on to image the same patient, as each modality provides unique anatomical and physiological information. 3D modalities such as CT and MRI provide information about the relationship of anomalous findings to surrounding structures and tissues. The Integrated Registration application provides you with a streamlined workflow to combine and display images from any two of the five major modalities (CT, MR, PET, SPECT, and X-ray angiography), aimed at improving your ability to provide more clinically relevant diagnostic information, and improving surgical and treatment planning.

Overview

Integrated Registration provides you with the capability to align and fuse two volumetric acquisitions from either the same or different acquisition modalities. With it, you can easily compare 3D anatomical images from CT, MR with PET, SPECT, and X-ray angiography¹ for a comprehensive analysis.

Highlights

- Load and register on the fly with drag and drop.
- Easy registration using automatic, manual and landmark methods.
- Multiple 2D and 3D fusion capabilities.
- Automatic registration propagation across series and one-click access to pre-defined protocols.
- Rigid and deformable registration.
- Drawing and saving of contours as RTSS DICOM objects.
- Access to regional registration on the fly for an easy review of challenging areas.
- Compatible with GE Advantage SIM MD and other standard treatment planning software.

Visit us:

www.gehealthcare.com/aw/applications/integrated-registration/

SECRET DE
AFACERI

Features

- Allows you to load exams and series from CT, MR, PET, SPECT or XA together for registration in the same session.
- Offers you a choice of four automatic registration methods.
- Allows you to customize screen layouts to suit your review needs and provides comparison of a current exam with previous exams to let you evaluate disease progression and treatment efficacy.
- Combines functional and anatomical images from different modalities, providing you with enhanced perspective of the area of interest.
- Allows you to define one or more contours around relevant anatomy and save them as RTSS DICOM objects.
- Provides you the ability to save registered data as new DICOM series or as a registration DICOM object (does not include SPECT).
- Define volume of interest in one model for automatic report to another model for subsequent use in Radiation Therapy (RT) planning.²
- Provide referring physicians with clear, detailed reports.

System Requirements

- Available on AW Server 3.1 and above and recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048).
- Available on AW4.7 and above

Image Requirements

Acquisition position and plane, slices resolution, number, thickness, and spacing may be different in the acquisition of the two exams to be registered.

Each image set should meet basic requirements, however:

- Field of view, matrix size, and display center should be the same for all images inside one series.
- Orientation should be the same for all images in the series.
- Series should include more than one image.
- Tilted acquisitions are not supported by mutual information based automatic algorithms; they can be registered with manual, landmark, or automatic "matching boundaries" registration methods.
- Datasets used as reference and moving should have some matching anatomical location.

Indications for Use

Integrated Registration provides easy means for comparison of three-dimensional (3D) images from Computed Tomography (CT), Magnetic Resonance Imaging (MRI), Emission Tomography (PET or SPECT) and X-Ray Angiography images (XA). To help physicians in diagnostic radiology or therapy planning, Integrated Registration allows 3D registration between volumetric acquisitions that may come from the same acquisition modality or from different acquisition modalities.

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References

¹For XA modality series, Integrated Registration currently supports only 3D X-Ray Angiography images (stored as CT Image Storage DICOM objects) acquired with GE Innova equipment and reconstructed with the Innova3DXR application.

²Requires DICOM RT capability on RT planning system.

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
JB32775XX
GE Medical Systems SCS
283 rue de la Minière
78530 Buc France.

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Vizitați-ne:

www.gehealthcare.com/raul/applications/integrated-registration/

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
SECRET DE
AFACERI

FlightPlan for EVAR

Simplified sizing & planning for EVAR

Planning of Abdominal and Thoracic Endovascular Aneurysm Repair (EVAR) procedures may be delicate specifically when dealing with complex anatomies.


Assessment of patients' eligibility for endovascular treatment of the aneurysm, sizing of the endograft and extraction of the meaningful information to help guide the procedure are essential to the success of the procedure.

Overview

FlightPlan for EVAR is a dedicated application to help plan EVAR procedures. It lets you visualize the vascular anatomy, perform key anatomical measurements, choose the treatment strategy, size the endograft and save key information that can be used during the intervention.

Highlights

- Plan an EVAR in a few minutes on MR or CT Angiography
- 0-click Bone Removal and tracking of the abdominal vessels on CTA
- Perform key aorta sizing measurements in a few clicks
- Analyze the profile of the endograft sealing zones
- Generate and save volumes, C-arm angulations and vessel ostia contours to be used during the procedure

Visit us:

www.gehealthcare.com/ow/applications/oncoquant/

Features summary

Visualization of the anatomy

0-click bone removal, automatic vessel tracking and labeling and centerline extraction of the abdominal aorta and vessels. Manual revision of the centerline can be performed to anticipate potential anatomical deformations happening during the procedure.

Anatomical measurements for endograft sizing

Key anatomical measurements dedicated to EVAR planning can be performed in few clicks and reported in a summary graphic.

To help user to position measurement points:

- The aorta bifurcation and ostium of visceral arteries are automatically detected and marked by arrows in lumen view
- The profile of the vessel diameter is displayed along the lumen view of the vessel

Additional measurements are available for EVAR planning:

- 0-click aneurysm neck analysis
- 1-click assessment of the endograft main body and leg length using rulers automatically positioned
- 1-click assessment for visceral vessels for distance
- 1-click clock position for planning of Fenestrated EVAR

Standard measurements are automatically reported on the graphic together with key anatomical

images.

The report with its key images is then available in pdf or DICOM format.

Preparation for pre-operative 3D image Fusion

Get prepared for pre-operative 3D Image Fusion in few clicks with:

- 0-click creation of the vascular and bone models for EVAR cases
- 0-click detection and saving of key vessel ostia contours
- Saving of the optimal C-arm angulations
- 1-click export of the data to EVAR Vision or saving for later use.

This preparation can be performed following the endograft sizing or as a stand-alone task performed just prior the beginning of the case.

System Requirements

- AW VolumeShare 7 workstation or above
- Volume Viewer
- Volume Viewer Innova
- VesselIQ Xpress
- AutoBone Xpress
- MR VesselIQ Xpress

Indications for Use

FlightPlan for EVAR CT is intended to provide an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment paths from a set of Computed Tomography (CT) Angiographic images.

FlightPlan for EVAR CT is a post processing application option which can be used in the analysis of 2D and 3D CT Angiography images/data derived from DICOM 3.0 compliant CT scans for the purpose of

cardiovascular and vascular disease assessment. This software is designed to support the physician in assessment of stenosis analysis, pre/post stent planning and directional vessel tortuosity visualization.

FlightPlan for EVAR CT automatic visualization tools provide the users with the capabilities to facilitate segmentation of bony structures for accurate identification of the vessels. Once vessels are visualized, tools are available for sizing the vessel, measure areas of abnormalities within a vessel.

FlightPlan for EVAR MR is intended to provide an optimized non-invasive application to facilitate vascular anatomy and pathology analysis from a set of DICOM 3.0 compliant 3D contrast-enhanced Magnetic Resonance Angiographic (MRA) images.

FlightPlan for EVAR MR is a post processing application, which can be used in the analysis of MRA data for the purpose of vascular disease assessment. This software is designed to assist radiologists and other clinicians in the evaluation and assessment of vascular anatomy and disease with the capability to provide a set of tools for visualizing directional vessel tortuosity, for sizing the vessel, measuring areas of abnormalities within a vessel.

Regulatory Compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

GE imagination at work

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, imagination at work, AutoBone, VesselIQ are trademarks of General Electric Company.
All other third party trademarks are the property of their respective owners.
DOC1761261.

SECRET DE
AFACERI

Features details

1. Visualization of the anatomy

0-click bone removal, automatic vessel tracking and labeling and centerline extraction of the abdominal aorta and vessels.


Manual revision of the centerline can be performed to anticipate potential anatomical deformations happening during the procedure after endograft insertion.

3D, lumen and axial views of the aorta and iliac arteries are automatically populated.

2. Anatomical measurements for endograft sizing

Key anatomical measurements dedicated to EVAR planning can be performed in few clicks and reported in a graphic summary. They include diameter, length and angle measurements which can provide additional information for you to determine if patients are eligible for endovascular treatment and to size the endograft.


To help user to position measurement points:


- the aorta bifurcation and ostium of visceral arteries are automatically detected and identified with arrows on the lumen view
- The profile of the vessel diameter is displayed along the lumen view of the vessel

Additional measurements are available for EVAR planning:


- 0-click aneurysm neck analysis with assessment of the neck profile, minimum and maximum diameters as well as percentage of diameter variation along the zone on CTA

- 1-click assessment of the endograft main body and leg length using 2 rulers automatically positioned at the proximal and distal endograft landing zone.

- 1-click assessment of visceral vessels, distance
- 1-click clock position for planning of Fenestrated EVAR


Standard measurements are automatically reported on the graphic together with key anatomical images. Manual measurements can also be saved on a report with associated images. The report with the key images is then available in pdf or DICOM format.

3. Pre-operative preparation of 3D images for Fusion


Get prepared for 3D Image Fusion in few clicks:

- 0-click creation of the vascular and bone models for EVAR cases
- 0-click detection and saving of key vessel ostia contours
- Stored positions of the optimal C-arm angulations
- 1-click export of the data to EVAR Vision or saving for later use.

This preparation can be performed following the endograft sizing or as a stand-alone task performed in the Discovery or Innova control room just prior to the beginning of the case.

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DOC1761261.

SECRET DE AFACERI

VessellQ™ Xpress with AutoBone™ Xpress

Fast and Easy to Use Vessel Analysis with Automated Real-Time Tracking, Auto-Labeling and Simplified User Interface.

Peripheral vascular disease (PVD) is a condition similar to coronary artery disease and carotid artery disease. In PVD, fatty deposits build up in the inner linings of the artery walls. These blockages restrict blood circulation, mainly in arteries leading to the kidneys, stomach, arms, legs and feet.

There are different clinical techniques used to diagnose PVD. Computed Tomography Angiography (CTA) can potentially provide sufficient information for therapeutic decision-making and reduced imaging costs for the management of PVD. One of the advantages of computed tomography (CT) is that it shows both vascular & nonvascular anatomy and also provides 3D information for best view. Also, CT may be more useful for patients with limited or no peripheral access or with stents, who may not be candidates for magnetic resonance angiography (MRA). Complementary to the ability to perform CTA is to perform clinical assessment by a dedicated post-processing tool.

Overview

VessellQ™ Xpress with Autobone Xpress is a post-processing software package for the Advantage Workstation (AW) platform, AW Server, CT scanners and PACS reading stations. It is a tool to analyze CT angiographic data including stenosis analysis, thrombus, pre and post stent planning procedures, and directional vessel tortuosity visualization. Autobone Xpress provides zero-click bone segmentation for head, neck and other anatomy. The automation and ease of use streamlines workflow.

Highlights

- Automated detection of aorta and iliacs with auto-labeling of vessels
- Automated real-time fast tracking for all vasculature
- New contextual wheel menu including easy access to advanced options
- New centerline edition capabilities
- Improved predefined measurements and edit capabilities
- New reporting interface which includes measurements in summary table
- New user interface

Visit us:

<https://www.gehealthcare.com/products/advanced-visualization/all-applications/autobone-vessellq-xpress>

SECRET DE AFACERI

Features

- Anatomy based protocols for guided workflow in the analysis for your peripheral vascular CTA acquisitions.
- Automatic bone removal for efficient review of vasculature.
- Automated abdominal aorta CTA protocol with vessel tracking and vessel labeling.
- **Fast tracking detects entire length of vessel and displays the centerline based on a single start point.**
- Advanced editing tools bringing new flexibility for faster editing based on targeted anatomy.
- **Thrombus detection with measurements and comparison over time.**
- Dynamic AVA provides access to vessel tracking at any time, from any protocol, on any vessel.
- Advanced measurement tools for stenosis detection and color identification of HU value for quantification of different densities.
- Create a queue of cases that are pre-loaded and pre-processed in the background while you read another exam.
- Measurements are organized in the summary table allowing for easy viewing and exporting.
- Fully customizable workflows with adjustable layouts, personalized parameter and settings, custom review steps.

- Save State lets you save and restore the state of the processed images at any stage

System Requirements

- AW Workstation AW 4.7 ext 14 or higher
- AW Server 3.2 Ext. 3.2 or higher

Indications for Use

VessellIQ Xpress is an optional, non-invasive, optimized, post-processing application intended to provide images and tools to analyze vascular anatomy and pathology, aiding physicians in diagnosis and determination of treatment paths, from a set of Computed Tomography (CT) Angiographic images. VessellIQ Xpress is an option for the Advantage Workstation (AW) platform, CT Scanner, and/or PACS, which can be used in the analysis of 2D and 3D CT Angiography images/data for the purpose of cardiovascular and vascular disease assessment. This software-only device is designed to support physician assessment for a wide variety of clinical uses such as stenosis analysis, pre/post stent planning, pre/post valve replacement planning, and directional vessel tortuosity visualization. VessellIQ Xpress' automatic visualization tools provide users the capability to segment bony structures for accurate identification of the vessels. Additional tools enable analysis of the vascular anatomy including the aorta, valves, and branching vessels for: anatomical sizing;

density and volume analysis of segmented vasculature and calcified / non-calcified plaque; and measurements of abnormalities. The TAVI Analysis option for VessellIQ Xpress is a planning tool used for Trans Aortic Valve Implantation (TAVI) procedures. It automatically segments the aorta and displays the aortic valve in multiple views for measurements of anatomic structures commonly needed for aortic annulus replacement planning. TAVI Analysis provides guided workflow and semi-automated tools to aid in evaluation of appropriate access pathways for interventional procedure planning.

AutoBone Xpress option is a software package that is intended to facilitate segmentation of bony structures and calcifications from abdominal and extremity CT Angiography data.

Regulatory Compliance

This product complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices (MDR). This product or its features may not be available in some other countries or regions. Please contact your sales associate.

Rx Only

GE imagination at work

© 2021 General Electric Company.
All rights reserved. Data subject to change.
GE, GE Monogram, Imagination at work, AutoBone Xpress, VessellIQ Xpress are trademarks of General Electric Company.
All other third-party trademarks are the property of their respective owners
Doc1471990

SECRET DE
AFACERI

Features Detail

Anatomy-based protocols


- The program offers you predefined protocols for fast, easy review of all vascular structures:
 - Intra-cerebral
 - Carotids
 - Thoracic Aorta
 - Abdominal Aorta
 - Runoff
 - Spinal
 - Extremity
- Protocols may be customized to meet your requirements and saved for quick access.

AutoBone™ Xpress

- Anatomical driven protocols with automatic bone removal.
- Improve vessel visualization by removing obstructive bony detail


- Automatic calcium removal for vessel lumen visualization

Auto Aorta


- With a single click, the program automatically removes bones and performs centerline tracking of the aorta and iliacs and labels vessels.

- Simply confirm the automated tracking or modify it by adding additional branch points, deleting a single point, or deleting multiple points from any view.


Fast Tracking

- Semi-automatic detection of complete vessel from a single starting point

- Advanced editing with variable interpolation tool for quick centerline edits adjusted for length of vessel.

- Bulls eye adjustment capability on x-section views to perfect tracking of the smallest most tortuous vessels.
- Volume measurements of the thrombus allow you to track change over time. Additional measurements include :
 - Minimum diameter
 - Maximum diameter
 - Mean diameter
 - Cross-sectional area


Vessel Label Data Base

- Each protocol is associated with a database of vessel names
- Customize vessel label database to meet your requirements for every anatomical region.
- A right mouse click provides access to the list of vessel names per anatomy.


Dynamic AVA

- Real-time tracking of the vessel center line and instant visualization of an unfolded view for quick vessel inspection.
- Immediate access to vessel analysis tools.
- Edit contours with a selection of drawing tools with variable interpolation for fast editing.

Advanced Vessel Analysis Tools

- Semi-automated vessel analysis lets you track, extract, visualize, and size vessels.
- Track a vessel and its branches for quick review in curved, oblique longitudinal or cross-sectional views.
- Multi-point deposition capabilities for improved center line tracking of challenging cases.
- When you edit a vessel contour, the program automatically interpolates the edit to neighboring slices.
- Generate and save rotating movies
 - Curved planar reformation and cross-sectional (lumen) views
 - Best L-section reformatted images
- Select vessel segments for quantitative or qualitative analysis.


- Available measurements include:
 - Distance
 - Cross-sectional area

- Single or dual reference point comparison
- Relative percent stenosis / aneurysm
- Color Identification
- 3D Angle
- Thrombus Extraction
- Volume
- Mean, min, max diameter
-
- Move freely between review steps and protocols without leaving the exam.

Stenosis Detection & ColorID

- With the customizable one- or two-point deposit technique, you can automatically quantify degree of stenosis and segment length.
- By assigning customizable color mapping to CT number ranges within an image, you can identify plaque on an axial image, curved reformat, or MPVR image. Up to ten levels of color coding are available. Colors transition smoothly from one density to another.
- Make volume and area measurements to track plaque size over time.

- Custom sizeable 3D ROI lets you apply color over all anatomy.
- Generate multiple ROIs on one vessel or multiple vessels.
- Automatically include each defined plaque, its volume and area into the report.

Quick AVA

- Dynamic AVA, provides access to vessel analysis from any protocol.
- Launch the Quick AVA feature from the main tool bar.
- Perform vessel analysis from any 3D or reformatted image with one- or two-point clicks.
- If you detect a lesion, you can analyze it by dropping a point above and below the section. The vessel analysis tool launches automatically for fast review in curved reformat, cross-sectional, lumen, and MRP view.
- Extend tracking proximally or distally for extended views at any time.


Summary Table & Report

- Measurements are automatically captured in summary table for easy review.

- Simply click on any measurement in the summary table and you will be taken to that measurement location.

- Prepare report with a single click or generate an advanced report with customized options
- Send the report to database, film or PDF.

Save State & Protocol


- Save the current state of tracking and measurements for easy reloading with a single mouse click.
Save state captures:
 - 3D images
 - Vessel tracking
 - Measurements
 - Vessel contours.
- Transfer the Save State to the PACS or other AW workstations for future review.
- Save custom protocols including anatomical labels and measurements.

Summary

VesselIQ Xpress provides user-friendly tools and an intuitive interface to help streamline your vessel analysis workflow and help you generate more comprehensive reports for referring physicians.


Multiple Study Comparison

- With VesselIQ Xpress, compare a previous study with the current study. You also have the flexibility to launch comparison within Quick AVA.

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