



RESPIRATORY
Continuous Positive Airway Pressure


Boussignac CPAP
Just breathe



Value Life

What is **Boussignac** CPAP ?

Boussignac CPAP is a Non Invasive Ventilation (NIV) device generating Continuous Positive Airway Pressure (CPAP).



How does **Boussignac** CPAP work?

- The operating principle of **Boussignac** CPAP is simple as it is based on the creation of a virtual valve.
- The pressure level obtained at this virtual valve depends on the gas flow rate supplied to the device. When you increase the flow rate, you increase the pressure. When you decrease the flow rate, you decrease the pressure.
- **Boussignac** CPAP is an open system and works without a mechanical valve. Thus, it adapts continuously to the medical treatment and to the patient's physiology, with optimum safety concerning the pressures generated.

What are the indications for Boussignac CPAP ?

• Pre-hospital & and Intra-Hospital Emergency Department Treatment

Main indication:

- ▶ Acute Cardiogenic Pulmonary Edema (ACPE): (6) Willem Dieperink et al - BMC Cardiovascular Disorders 2007
 - Application of CPAP: 7.5 to 10 cmH₂O

Other clinical indications:

- ▶ Conscious drowned person: (7) Dottorini M, et al - Chest. 1996
- ▶ Acute severe asthma: Boussignac CPAP can be used in conjunction with nebulization: (9) Laurent Brochard et al - Respiratory Care. Oct. 2011 (10) Service d'Aide Médicale d'Urgence (SAMU) de Lille: Protocole d'utilisation de la CPAP de Boussignac avec nébuliseur. 14 mars 2007

Proven benefits of Boussignac CPAP

- High FiO₂
- Rapid hypoxia correction with an increase of SpO₂
- Rapid correction of ACPE's clinical signs
- Patient's tolerance to an open system

• Intra-hospital treatment

Operating room

- ▶ Bariatric surgery, cardiac surgery, thoraco-abdominal surgery: (12) Wong D., et al - 2011. Can J Anesth
- ▶ Preoperative: pre-oxygenation : (11) Delay J.M., Jaber S. - Presse médicale 2012
- ▶ Extubation with positive pressure
- ▶ Postoperative: immediately after extubation: (13) Neligan P., et al. - 2009. Anesthesiology

Intensive care unit and resuscitation

- ▶ Post-resuscitation: stabilization and weaning after extubation: (14) Dieperink W, et al. - 2008. Respiration
- Fiberoptic bronchoscopy in hypoxemic patients: (16) Jaber S., et al - Am J Resp Crit Care Med Vol.162, 2000


Paediatrics

- ▶ Infant bronchiolitis: (17) Fleming P.F., et al. - 2012. J Paediatr Child Health


Proven benefits of Boussignac CPAP

- Alveolar recruitment
- Prevention of atelectasis
- Improved lung volume
- Reduced risk of reintubation
- Decreased length of in-patient hospital stay
- Proven decreased rate of ventilator associated pneumonia and duration of oxygen requirement without prolonging the hospital stay

What are the features of a good CPAP device?




- **The Work-of-Breathing (WOB) needs to be reduced for the patient: Boussignac CPAP decreases the WOB**
 The **WOB** is the expenditure of energy required to inspire air in the lungs.
 The **Delta P** is the difference between inspiratory pressure and expiratory pressure.



Boussignac CPAP



Competitor A



Competitor B

- Expiratory work - Pressure (Tracheal) [cmH₂O]
- CPAP level - Pressure (Airway) [cmH₂O]
- Inspiratory work - Pressure (Oesophagus) [cmH₂O]

(1) Valero P, Khoury A., et al. Poster: SRLF-2013

- **According to the above schematics, Boussignac CPAP was reported to have the same performances as complex capital equipment respiratory devices and has the advantage of having a stable CPAP level.**

Boussignac CPAP is...

Safe


- **Open System:** if necessary the patient can breathe atmospheric air and thus constantly breathe the required volume of gas
- Spontaneous ventilation is possible even if the gas flow stops
- **No mechanical parts**
- Precise control of pressures generated with specific manometer
- **No risk of barotrauma / volotrauma:** the pressure in the lungs cannot be superior to the one in the CPAP valve thanks to open system
- **No risk of hypoventilation**

Easy to use

- Light
- User-friendly
- **Mucus suctioning is possible without interruption of the treatment**
- Possibility to connect a nebulizer
- Regulation of FiO₂ with an attachable ring

Effective

- **Decreases the WOB**
- Flow rate of gas available for inspiration is **280 to 320 l/min**
- The difference between inspiratory pressure and expiratory pressure (Delta P) is only **1.5 +/- 0.2 cmH₂O (8)**
- Regulation of inspired FiO₂



Boussignac CPAP Features:

- Pressure variation in the airways during respiratory cycle is near 1 cmH₂O, irrespective of the CPAP level access chosen
- Airways remain accessible for any other medical intervention (aspiration, bronchoscopy...etc)
- CPAP level is not altered and medical interventions can be applied without any disruption of the treatment
- Peak flow rate is high, thus reducing the WOB
- In practice: clinical signs of hypoxia are rapidly reduced and O₂ saturation increases rapidly



Boussignac CPAP

What are the features and benefits?



No risk of misconnection



Supplementary port

Allows pressure measurement, EtCO₂ monitoring or additional oxygen or air delivery



FiO₂ ring



Nebulizer

Boussignac CPAP and Nebulization

A specific nebulizer fitted with a separate oxygen extension tube, allows an optimum flow rate (6l/min) for drug particles between 2 to 4 µm ensuring a perfect broncho-alveolar spread.

Integral O₂ tube

With international colour coding for O₂ delivery




Printed tube

Clearly marked tubing for improved safety



Special « star » shaped O₂ tube

Reduces the risk of tube kinking and the associated cessation of gas delivery



Noise reducer

This device reduces the noise by 6 decibels. With the noise reducer Boussignac CPAP now reaches a level of approximately 69 decibels (=to shower running 70 dbls). This accessory is optimum for intra-hospital use.

Product range and ordering information



Boussignac CPAP	
CPAP & manometer connector	code 5570.13

Single use accessories	
Manometer connector	code 5558.053 / 5558.203
E.T. Tube connector	code 555.01
vyconnector (Y connector)	code 884.06
Nebulizer	code 5569.01
FiO ₂ ring (angle)	code 5566.01
FiO ₂ ring (straight)	code 5566.02
Noise Reduction Device	code 5558.91

Mask and Harness	
Mask (size 4) ≈ S	code 5557.45
Mask (size 5) ≈ M	code 5557.55
Mask (size 6) ≈ L	code 5557.65
Silicone fixation harness	code 5559.01
Fabric fixation harness	code 5559.03

Reusable accessories	
O ₂ Flowmeter (Afnor connection)	code 5563.02
O ₂ Flowmeter (BSI connection)	code 5563.41
O ₂ Flowmeter (DIN connection)	code 5563.42
O ₂ Flowmeter (Nordic connection)	code 5563.44
Manometer	code 527.01
Air flowmeter (Afnor connection)	code 5563.01
Air flowmeter (BSI connection)	code 5563.31
Air flowmeter (DIN connection)	code 5563.32
Air flowmeter (Nordic connection)	code 5563.34

Emergency CPAP kit contents

2 x CPAP & 2 x Manometer connector	code 5570.13
1 Mask (size 4) ≈ S	code 5557.45
1 Mask (size 5) ≈ M	code 5557.55
1 Mask (size 6) ≈ L	code 5567.65
1 Harness	code 5559.01
1 O ₂ Flowmeter	
1 Manometer	code 527.01
1 x 20ml syringe	
1 nebulizer	code 5569.01
1 FiO ₂ ring	code 5566.01



Emergency CPAP kit code

code 5562.802 (with Afnor flowmeter)

code 5562.841 (with BSI flowmeter)

code 5562.842 (with DIN flowmeter)

code 5562.700 (without flowmeter)



CPAP set with nebulizer

1 CPAP
1 mask
1 manometer connector
1 nebulizer + T piece

code 5571.303/403/503/603 (scented mask)
code 5571.300/400/500/600 (odourless mask)

1 CPAP
1 mask
1 manometer connector
1 nebulizer + T piece
1 harness

code 5572.303/403/503/603

1 CPAP
1 nebulizer + T piece

code 5577.013

CPAP set without nebulizer

1 CPAP
1 mask
1 manometer connector

code 5561.303/403/503/603

1 CPAP
1 Silicone harness
1 mask
1 manometer connector

code 5562.303/403/503/603

References

- (1) Valero P., Khoury A, et al. Comparison between 3 devices delivering continuous Positive Airway Pressure (CPAP). Poster. SRLF-2013.
- (2) Richard JC, Cordioli RL, Brochard L, et al. Testing the Boussignac CPAP system on an active bench model simulating spontaneous ventilation and comparing its efficacy and resistive properties to other CPAP systems. Laboratory University Hospital, Geneva.
- (3) Templier F, et al. Boussignac continuous positive airway pressure system: practical use in a prehospital medical care unit. 2003. Eur J Emerg Med 10(2): 87-93.
- (4) Mattu A. Lawner B. Management of congestive heart failure. Heart Fail Clin. 2009 : 19-24.
- (5) Leman P, et al. Simple lightweight disposable continuous positive airways pressure mask to effectively treat acute pulmonary edema. Emerg Med Australia 17 (3): 224-230.
- (6) Willem Dieperink, Iwan Van der Horst, et al. Boussignac continuous positive airway pressure for the management of acute cardiogenic pulmonary edema. BMC Cardiovascular Disorders 2007, 7: 40.
- (7) Dottorini M, et al. Nasal-continuous positive airway pressure in the treatment of near-drowning in fresh-water. Chest. 1996 ; 110 : 1122-1124.
- (8) M. Chinellato, A.C. Astolfi, L. Aigle, E. Chinellato. Œdème aigu du poumon d'immersion. Ann. Fr. Med. Urgence (2015) 5 : 187-188.
- (9) Laurent Brochard, Bernard Maitre, et al. Aerosol delivery with the Boussignac CPAP device. Respiratory Care. Oct. 2011.
- (10) Service d'Aide Médicale d'Urgence (SAMU) de Lille: Protocole d'utilisation de la CPAP de Boussignac avec nébuliseur. 14 mars 2007.
- (11) Delay J.M., Jaber S. Respiratory preparation before surgery in patients with chronic respiratory failure. Presse médicale 2012. 41 : 225-233.
- (12) Wong D., et al. A comparison between the Boussignac continuous positive airway pressure mask and the venturi mask in terms of improvement in the PaO₂/FiO₂ ratio in morbidly obese patients undergoing bariatric surgery. 2011. Can J Anesth. 58 (6) : 532-539.
- (13) Neligan P, et al. Continuous positive airway pressure via the Boussignac system immediately after extubation improves lung function in morbidly obese patients with obstructive sleep apnea undergoing laparoscopic bariatric surgery. 2009. Anesthesiology 110 (4) : 878-884.
- (14) Dieperink W, et al. Boussignac continuous positive airway pressure for weaning with tracheostomy tubes. 2008. Respiration 75 (4) : 427-431.
- (15) Belenguer-Muncharaz, A, et al. Non Invasive Ventilation in severe pneumonia due to H1N1 virus. Med Intensiva 2011 35 (8): 470-477.
- (16) Maitre b, jaber s, Maggiore s.m, bergot e, richard j.c, Bakthiari h, housset b, Boussignac g, brochard l, Unité de Pneumologie, Service de Soins Intensifs et Unité INSERM U 492, Hôpital H. Mondor, AP-HP, Université Paris 12 - Créteil (94) France - Continuous Positive Airway Pressure during Fiberoptic Bronchoscopy in Hypoxemic Patients - Am J Resp Crit Care Med Vol.162. pp 1063-1067, 2000.
- (17) Fleming P.F., et al. Use of continuous positive airway pressure during stabilisation and infants with suspected bronchiolitis. 2012. J Paediatr Child Health.

Notes

 ANAESTHESIA EMERGENCY

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
Export customer service: +33 (0)1.39.92.64.15
Fax: +33 (0)1.39.92.64.44 • www.vygon.com


www.vygon.com